



IMMULITE 3gAllergy

Allergen Component Testing

A valuable diagnostic decision-making tool

[siemens-healthineers.com](https://www.siemens-healthineers.com)

Introduction

Component allergen testing provides clinicians with improved specificity for IgE-mediated allergies as compared to allergen extract testing. Following extract testing, component testing allows for differential diagnoses, differentiation between co-sensitization and cross-reactivity, and improvement of inclusion criteria for immunotherapy as well as the follow-up of immunotherapy.¹

Allergen extracts contain both allergenic and non-allergenic elements, offering direction as to whether reflex testing with component allergens should be conducted. Component allergens are purified, isolated allergenic proteins found within the extract. Allergen extracts provide high levels of sensitivity, whereas allergen components offer the advantage of high specificity and help to generate a comprehensive clinical picture. Thus, testing with extracts often precedes component testing, which may be recommended as a follow-up step. Without component testing, there may be not enough specificity in extract testing for a number of allergens to provide for a differential diagnosis.²

Component testing is also advantageous to differentiate between cross-reactivity and co-sensitization—key for both allergen and autoimmune testing. Cross-reactivity occurs when IgE antibodies recognize allergens from one source but react to a similar protein from another source. Co-sensitization involves the presence of IgE toward epitopes that are not shared between allergenic sources.^{1,2}

As with allergen extracts, the quality of a component remains paramount, and the Siemens Healthineers portfolio of 3gAllergy™ specific allergens complies with the acceptable clinical and analytical performance standards. The IMMULITE® 3gAllergy™ component allergen menu, containing some of the most prevalent and reactogenic allergens, provides clinicians with the necessary tools to satisfy comprehensive allergen testing needs. And as an aid in clinical decision making, literature-supported testing algorithms are provided to guide clinicians through the recommended step-by-step testing process.

For this white paper, 11 3gAllergy component allergens were selected with the goal of:

- Briefly explaining component allergens
- Reviewing the testing algorithms based on published literature
- Showing clinical performance

The 11 selected allergens cover a spectrum of component allergens available from Siemens Healthineers:

- Apple Component Allergen, rMal d 1 (*Malus domestica*)
- Apple Component Allergen, rMal d 4 (*Malus domestica*)
- Birch Pollen Component Allergen, nBet v 1 (*Betula verrucosa*)
- Birch Pollen Component Allergen, rBet v 2 (*Betula verrucosa*)
- Cherry Component Allergen, rPru av 1 (*Prunus avium*)
- Cherry Component Allergen, rPru av 3 (*Prunus avium*)
- Cherry Component Allergen, rPru av 4 (*Prunus avium*)
- Dust Mites Component Allergen, nDer p 1 & nDer p 2 (*Dermatophagoides pteronyssinus*)
- Mugwort Pollen Major Allergen, nArt v 1 (*Artemisia vulgaris*)*
- Olive Pollen Component Allergen, nOle e 1 (*Olea europaea*)
- Peach Component Allergen, nPru p 3 (*Prunus persica*)


*Not available for sale in the U.S.

Apple Component Allergen, rMal d 1 (*Malus domestica*)

Mal d 1 is a PR-10 protein associated with oral allergy syndrome (OAS) to apple.³ It is a homologous protein to Bet v 1, which has been identified as a primary pollen sensitizer eliciting specific IgE antibodies.⁴ Although Mal d 1 shares only 57% sequence homology with Bet v 1, 75% of the Mal d 1 tertiary structure binds anti-Bet v 1 antibodies.⁵ Approximately 50–93% of birch pollen-allergic patients develop concomitant OAS reactions to fruits, nuts, and vegetables. Mal d 1-allergic individuals do not experience systemic reactions, as PR-10 proteins are susceptible to degradation by heat and gastric digestion.^{6,7}



Testing algorithm according to published literature^{3,4,6,8}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the rMal d 1-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	50	1	51
Negative	4	116	120
Total	54	117	171

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
92.6% (82.1–97.9%)	99.1% (95.3–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F49 (Reference Method) A464 (Test Method)	145	92% (133/145)	77% (36/47)	99% (97/98)


[†]This data is from Siemens Healthineers verification studies.



Apple Component Allergen, rMal d 4 (*Malus domestica*)

Mal d 4, a member of the profilin family of allergens, is an approximately 14 kD protein associated with oral allergy syndrome (OAS) to apple.⁹ It is a homologous protein to Bet v 2, which has been identified as a minor pollen sensitizer in approximately 10–30% of pollen-allergic individuals.¹⁰ As with other profilins, Mal d 4 is highly thermolabile and rapidly enzyme-degradable, accounting for its inability to elicit systemic reaction and its loss of allergenicity in cooked foods.⁸ Primary sensitization to Mal d 4 develops through pollinosis and cross-reactivity to Bet v 2 and is not presumed to arise directly from apple ingestion without previous sensitization to birch or grass profilin.^{3,8,11}

Testing algorithm according to published information^{3,4,8,10-12}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the rMal d 4-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	15	1	16
Negative	39	116	155
Total	54	117	171

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
27.8% (16.5–41.6%)	99.1% (95.3–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]


	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F49 (Reference Method)	159	79% (126/159)	49% (27/55)	95% (99/104)
A796 (Test Method)				

[†]This data is from Siemens Healthineers verification studies.

Birch Pollen Component Allergen, nBet v 1 (*Betula verrucosa*)

Birch (*Betula verrucosa*) pollen is a major cause of allergy in the northern parts of Europe and America that is attributed to about 22% of allergic individuals suffering from pollinosis.¹³ Six birch pollen allergens have already been identified (Bet v 1, Bet v 2, Bet v 3, Bet v 4, Bet v 6, Bet v 7). Bet v 1, a major allergenic protein of birch pollen, is a 17 kDa protein consisting of several isoallergens¹⁴ and recognized by IgE antibodies from almost all birch pollen-allergic patients.¹⁵ Oral allergy syndrome (OAS) as a result of primary sensitization to Bet v 1 has been previously reported.

Testing algorithm according to published literature^{13,16-18}



Clinical Performance⁶⁰

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the nBet v 1-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	37	5	42
Negative	6	95	101
Total	43	100	143

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
86% (not available)	95% (not available)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F49 (Reference Method) A796 (Test Method)	143	92% (132/143)	79% (42/53)	100% (90/90)


[†]This data is from Siemens Healthineers verification studies.



Birch Pollen Component Allergen, rBet v 2 (*Betula verrucosa*)

Profilin allergens are responsible for multiple pollen and food sensitization with extensive cross-reactivity.^{19,20} The IgE conformational epitopes in profilin are highly conserved, which is key to maintaining its cross-reactive nature in plant allergens. The Bet v 2 profilin from birch pollen may be used to evaluate IgE reactivity in patients with suspected birch allergy or other cross-reactivities between mugwort, grass pollen, celery, carrots, and hazelnut.²¹⁻²³

Testing algorithm according to published literature^{16-18, 24}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the rBet v 2-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥0.10 kU/L)	25	1	26
Negative	45	116	161
Total	70	117	187

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
35.7% (24.6–48.1%)	99.1% (95.3–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
T3 (Reference Method) A89 (Test Method)	156	80% (125/156)	51% (30/59)	98% (95/97)


[†]This data is from Siemens Healthineers verification studies.

Cherry Component Allergen, rPru av 1 (*Prunus avium*)

Pru av 1 is a 17.7 kDa PR-10 protein associated with oral allergy syndrome (OAS) to cherry.^{10,25} It is a homologous protein to Bet v 1, which has been identified as a primary pollen sensitizer eliciting specific IgE antibodies and is considered a major allergen.^{10,25,26} Although Pru av 1 shares only 59–64% sequence homology with Bet v 1, 75% of the tertiary structures of the two proteins are virtually identical, and preincubation of cherry-allergic patient sera with Bet v 1 inhibits binding by Pru av 1.²⁷ Up to 90% of cherry-allergic patients manifest a concomitant allergy to birch pollen, as primary sensitization arises via pollinosis.^{28,29} Pru av 1-allergic individuals typically do not experience systemic reactions, as PR-10 proteins are rapidly degraded by heat and gastric digestion.^{28,29}



Testing algorithm according to published literature^{10,12,25-29}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the rPru av 1-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	33	5	38
Negative	4	112	116
Total	37	117	154

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
89.2% (74.6–97.0%)	95.7% (90.3–98.6%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F242 (Reference Method)	159	89% (142/159)	75% (46/61)	98% (96/98)
A597 (Test Method)				


[†]This data is from Siemens Healthineers verification studies.



Cherry Component Allergen, rPru av 3 (*Prunus avium*)

Plant lipid transfer proteins (LTP) are highly conserved proteins of approximately 10 kDa.³⁰ They are typically associated with more-severe and systemic reactions such as urticaria and anaphylaxis in some populations, but less-severe reactions in the form of oral allergy syndrome (OAS) in others.^{28,29} Because these proteins are so highly conserved, sensitization to LTP from one plant can result in allergic responses to other taxonomically related or unrelated fruits and vegetables.²⁸⁻³⁰ Pru av 3 is an LTP isolated from cherry^{26,31} and may be used to evaluate specific IgE reactivity in patients with suspected cherry allergy.^{26,28-31} Monosensitization to Pru av 3 is rare; studies suggest that peach LTP (Pru p 3) is the likely primary allergic sensitizer, triggering allergic response to other Rosaceae via LTP cross-reactivity.^{28,29}

Testing algorithm according to published literature^{12,26,28-31}




Cherry Component Allergen, rPru av 4 (*Prunus avium*)

Pru av 4, a member of the profilin family of allergens, is an approximately 15 kDa protein associated with oral allergy syndrome (OAS) to cherry.^{28,32} It is a homologous protein to Bet v 2, which has been identified as a minor pollen sensitizer in approximately 10–30% of pollen-allergic individuals.^{10,20,28} Pru av 4 shares over 70% amino acid homology and a very similar tertiary structure with other profilins.²⁰ It is highly thermolabile and rapidly enzyme-degradable, accounting for its inability to elicit systemic reaction and its loss of allergenicity in cooked foods.¹⁰ Primary sensitization to Pru av 4 develops through pollinosis and cross-reactivity to Bet v 2 and is not presumed to arise directly from cherry ingestion without previous sensitization to birch or grass profilin.¹⁰



Testing algorithm according to published literature^{12,20,26,28,29}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the rPru av 4-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	9	1	10
Negative	28	116	144
Total	37	117	154

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
24.3% (11.8–41.2%)	99.1% (95.3–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F242 (Reference Method) A600 (Test Method)	156	80% (125/156)	50% (26/52)	95% (99/104)

[†]This data is from Siemens Healthineers verification studies.



Dust Mites Component Allergen, nDer p 1 and nDer p 2 (*Dermatophagoides pteronyssinus*)

Development of allergic asthma in adults and children is commonly associated with sensitization to house dust mites (HDM).³³ IgE sensitivities are often directed toward major HDM allergens such as Der p 1 (24 kDa) and Der p 2 (15 kDa) from *Dermatophagoides pteronyssinus* (allergen code: D1).^{34,35} Der p 1 and Der p 2 individually are recognized by more than 80% of D1-sensitized patient IgE and have a predictive value for D1 greater than 95%.³⁶ Cross-reactivity has been observed between homologous proteins Der f 1 and Der f 2 from *Dermatophagoides farinae*, to which *D. pteronyssinus* shares 80–90% sequence identity.³⁷

Clinical Performance⁶²

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the nDer p 1- and nDer p 2-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
nDer p 1			
Positive (≥0.10 kU/L)	49	3	52
Negative	7	109	116
Total	56	112	168
nDer p 2			
Positive (≥0.10 kU/L)	51	2	53
Negative	5	110	115
Total	56	112	168

Table 2. Clinical performance: sensitivity and specificity.


nDer p 1		nDer p 2	
Sensitivity (95% CI)	Specificity (95% CI)	Sensitivity (95% CI)	Specificity (95% CI)
88% (79–96%)	97% (94–100%)	91% (84–99%)	98% (96–101%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

nDer p 1 and nDer p 2	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
D1 (Reference Method)	107	99% (106/107)	98% (50/51)	100% (56/56)
A310 and A316 (Test Method)				

Testing algorithm according to published literature^{12,38,39}




[†]This data is from Siemens Healthineers verification studies.

Mugwort Pollen Major Allergen, nArt v 1 (*Artemisia vulgaris*)*

Mugwort (*Artemisia vulgaris*) is one of the main causes of seasonal pollinosis in Europe and is also found throughout the Northern Hemisphere.⁴⁰ Six mugwort pollen allergens have been identified (Art v 1, Art v 2, Art v 3, Art v 4, Art v 5, Art v 6).⁴¹ Art v 1, a glycoprotein, has an approximate molecular weight of 24–28 kDa and appears as a double-band due to its heterogenous glycosylation.⁴⁰ Native Art v 1 is preferred over the recombinant molecule produced by *E. coli* due to only 30–50% recognition by nArt v 1-positive sera.⁴² Art v 1 is a major allergenic protein of mugwort pollen that is recognized by more than 70% of mugwort-sensitized patient IgE.⁴³ A homologous protein in ragweed (Amb a 4) has been identified.⁴⁴



Testing algorithm according to published literature^{12,45}



Clinical Performance⁶⁰

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the nArt v 1-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	20	7	27
Negative	10	93	103
Total	30	100	130

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
67% (50–84%)	93% (88–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
W6 (Reference Method) A753 (Test Method)	142	83% (117/142)	93% (25/27)	80% (92/115)

*Not available for sale in the U.S.


†This data is from Siemens Healthineers verification studies.



Olive Pollen Component Allergen, nOle e 1 (*Olea europaea*)

Olive tree (*Olea europaea*) pollen is a major cause of type 1 seasonal allergy in the Mediterranean and can also be found in areas with a Mediterranean-type climate, including parts of Australia, South Africa, and North America.^{46,47} Ten olive tree pollen allergens have been identified (Ole e 1, Ole e 2, Ole e 3, Ole e 4, Ole e 5, Ole e 6, Ole e 7, Ole e 8, Ole e 9, Ole e 10).⁴⁸ Ole e 1 is a major allergenic protein of olive tree pollen that exists in two main forms: glycosylated (~20 kDa) and non-glycosylated (~18.5 kDa).⁴⁹ Homologous proteins are found in other members of the Oleaceae family, such as ash (Fra e 1), privet (Lig v 1), and lilac (Syr v 1).⁵⁰ More than 80% of patients sensitized to olive pollen have IgE reactivity to the Ole e 1 allergenic molecule, and common symptoms of exposure include asthma, rhinitis, and conjunctivitis.⁵¹⁻⁵³

Testing algorithm according to published literature^{12,53,54}



Clinical Performance⁶⁰

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the nOle e 1-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	37	3	40
Negative	2	97	99
Total	39	100	139

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
95% (not available)	97% (not available)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
T9 (Reference Method) A482 (Test Method)	139	97% (135/139)	98% (39/40)	97% (96/99)


[†]This data is from Siemens Healthineers verification studies.

Peach Component Allergen, nPru p 3 (*Prunus persica*)

Plant lipid transfer proteins (LTP) are highly conserved proteins of approximately 10 kDa.³⁰ They are found in the seeds, stems, flowers, leaves, and pollen of plants and are considered panallergens.³⁰ LTPs are typically associated with more-severe and systemic reactions such as urticaria and anaphylaxis.^{30,55} Because these proteins are so highly conserved, sensitization to LTP from one plant can result in unexpected allergic responses to other taxonomically unrelated fruits or vegetables.^{30,55,56} Pru p 3 is an LTP isolated from peach skin and pulp³⁰ and may be used to evaluate specific IgE reactivity in patients with suspected peach allergy.^{30,55-57} While peach LTP is typically the primary allergic sensitizer, cross-reactivity to LTP of other taxonomically related or unrelated plants may be responsible for peach allergy in some patients.⁵⁷



Testing algorithm according to published literature^{12,55-59}



Clinical Performance⁶¹

Clinical performance was demonstrated by testing serum samples from clinically diagnosed atopic patients and apparently healthy individuals against the nPru p 3-specific allergen. The results were obtained using the IMMULITE 2000/XPi 3gAllergy Specific IgE assay. Overall agreement, sensitivity, and specificity are presented in the tables below.

Table 1. Clinical performance: overall agreement.

	Atopic	Non-atopic	Total
Positive (≥ 0.10 kU/L)	11	0	11
Negative	27	117	144
Total	38	117	155

Table 2. Clinical performance: sensitivity and specificity.

Sensitivity (95% CI)	Specificity (95% CI)
28.9% (15.4–45.9%)	100% (96.9–100%)

Additional clinical performance of the specific allergens was demonstrated in comparison to the whole extract allergen. The results are presented below.

Table 3. Clinical performance of the specific allergens in comparison to the whole extract allergen.[†]

	Number of Samples	Overall Percent Agreement	Positive Percent Agreement	Negative Percent Agreement
F95 (Reference Method) A603 (Test Method)	156	84% (122/146)	49% (23/47)	100% (99/99)

[†]This data is from Siemens Healthineers verification studies.

Conclusion

Allergen component testing plays a critical role, as it allows clinicians to provide a better differential diagnosis than extract testing alone. The era of protein mixture testing has given way to this new frontier of highly specific, single-allergen component testing, proving it to be an aid in the clinical diagnosis of IgE-mediated allergic disorders. Being equipped with an in-house portfolio of rigorously tested and validated component allergens may provide for a more comprehensive and accelerated testing process.

The Siemens Healthineers portfolio of commonly tested component allergens is currently available for testing on the IMMULITE 2000/XPi platform and has undergone rigorous clinical and analytical scrutiny to validate the quality of the tests. Please refer to the IMMULITE 2000/XPi 3gAllergy menu online for the full list of available component allergens.

References:

- Hoffmann K, Hilger C, Santos A, de las Vecillas L, Dramburg S. Molecular Allergy Users Guide 2.0. EAACI 2022. Eds. https://hub.eaaci.org/resources_guidelines/molecular-allergy-users-guide-2-0/
- Ansotegui I, et al. IgE allergy diagnostics and other relevant tests in allergy, a WHO position paper. WAO J. 2020;13:100080.
- Fernández-Rivas M, Bolhaar S, González-Mancebo E, Asero R, van Leeuwen A, Bohle B, et al. Apple allergy across Europe: how allergen sensitization profiles determine the clinical expression of allergies to plant foods. J Allergy Clin Immunol. 2006;118(2):481-8.
- Ebo DG, Bridts CH, Verweij MM, De Knop KJ, Hagendorens MM, De Clerck LS, et al. Sensitization profiles in birch pollen-allergic patients with and without oral allergy syndrome to apple: lessons from multiplexed component-resolved allergy diagnosis. Clin Exp Allergy. 2010 Feb;40(2):339-47.
- Jenkins JA, Griffiths-Jones S, Shewry PR, Breiteneder H, Mills EN. Structural relatedness of plant food allergens with specific reference to cross-reactive allergens: an in silico analysis. J Allergy Clin Immunol. 2005 Jan;115(1):163-70.
- Mauro M, Russello M, Incorvaia C, Gazzola G, Frati F, Moingeon P, et al. Birch-apple syndrome treated with birch pollen immunotherapy. Int Arch Allergy Immunol. 2011;156(4):416-22.
- Wuthrich B, Schindler C, Leuenberger P, Ackermann-Lieblich U. Prevalence of atopy and pollinosis in the adult population of Switzerland (SAPALDIA study). Swiss Study on Air Pollution and Lung Diseases in Adults. Int Arch Allergy Immunol. 1995 Feb;106(2):149-56.
- Andersen MB, Hall S, Dragsted LO. Identification of European allergy patterns to the allergen families PR-10, LTP, and profilin from Rosaceae fruits. Clin Rev Allergy Immunol. 2011;41(1):4-19.
- Oberhuber C, Ma Y, Marsh J, Rigby N, Smole U, Radauer C, et al. Purification and characterisation of relevant natural and recombinant apple allergens. Mol Nutr Food Res. 2008;52(suppl 2):S208-19.
- Breiteneder H, Radauer C. A classification of plant food allergens. J Allergy Clin Immunol. 2004;113(5):821-30.
- Ma Y, Zuidmeer L, Bohle B, Bolhaar ST, Gadermaier G, Gonzalez-Mancebo E, et al. Characterization of recombinant Mal d 4 and its application for component-resolved diagnosis of apple allergy. Clin Exp Allergy. 2006;36(8):1087-96.
- Matricardi PM, et al. EAACI molecular allergology user's guide. Pediatr Allergy Immunol. 2016 May;27 Suppl 23:1-250.
- Horak F, Jager S. Die Erreger des Heufiebers. 19 pp. Verlag Urban & Schwarzenberg Munich-Vienna-Baltimore, 1979.
- Swoboda I, et al. J Biol Chem. 1995;270:2607-13.
- Ipsen H, Lowenstein H. J Allergy Clin Immunol. 1983;72:150-9.
- Menz G, Dolecek C, Schonheit-Kenn U, Ferreira F, Moser M, Schneider T, Suter M, Boltz-Nitulescu G, Ebner C, Kraft D, Valenta R. Serological and skin-test diagnosis of birch pollen allergy with recombinant Bet v 1, the major birch pollen allergen. Clin Exp Allergy. 1996;26(1):50-60.
- Kazemi-Shirazi L, Pauli G, Purohit A, Spitzauer S, Froschl R, Hoffmann-Sommergruber K, Breiteneder H, Scheiner O, Kraft D, Valenta R. Quantitative IgE inhibition experiments with purified recombinant allergens indicate pollen-derived allergens as the sensitizing agents responsible for many forms of plant food allergy. J Allergy Clin Immunol. 2000;105(1 Pt 1):116-25.
- Moverare R, Westritschnig K, Svensson M, Hayek B, Bende M, Pauli G, Sorva R, Haahtela T, Valenta R, Elfman L. Different IgE reactivity profiles in birch pollen-sensitive patients from six European populations revealed by recombinant allergens: an imprint of local sensitization. Int Arch Allergy Immunol. 2002;128(4):325-35.
- Valenta R, Duchene M, Pettenburger K, Sillaber C, Valent P, Bettelheim P, Breitenbach M, Rumpold H, Kraft D, Scheiner O. Identification of profilin as a novel pollen allergen; IgE autoreactivity in sensitized individuals. Science. 1991;253(5019):557-60.
- Scheurer S, Wangorsch A, Nerkamp J, Skov PS, Ballmer-Weber B, Wuthrich B, Hausteiner D, Vieths S. Cross-reactivity within the profilin panallergen family investigated by comparison of recombinant profilins from pear (Pyr c 4), cherry (Pru av 4) and celery (Api g 4) with birch pollen profilin Bet v 2. J Chromatogr B Biomed Sci Appl. 2001;756(1):315-25.
- Valenta R, Sperr WR, Ferreira F, Valent P, Sillaber C, Tejkl M, Duchene M, Ebner C, Lechner K, Kraft D. Induction of specific histamine release from basophil with purified natural and recombinant birch pollen allergens. J Allergy Clin Immunol. 1993;91(1):88-97.
- Van Ree R, Fernandez-Rivas M, Cuevas M, van Wijngaarden M, Aalberse RC. J Allergy Clin Immunol. 1995;95(3):726-34.
- Ferreira F, Hawranek T, Gruber N, Wopfner N, Mari A. Allergy. 2004;59:243-67.
- Matricardi PM, et al. EAACI molecular allergology user's guide. Pediatr Allergy Immunol. 2016;27(suppl 23):1-250.
- Scheurer S, Metzner K, Hausteiner D, Vieths S. Molecular cloning, expression and characterization of Pru a 1, the major cherry allergen. Mol Immunol. 1997;34(8-9):619-29.

26. Scheurer S, Pastorello EA, Wangorsch A, Kästner M, Haustein D, Vieths S. Recombinant allergens Pru av 1 and Pru av 4 and a newly identified lipid transfer protein in the in vitro diagnosis of cherry allergy. *J Allergy Clin Immunol*. 2001;107(4):724-31.
27. Neudecker P, Schweimer K, Nerkamp J, Scheurer S, Vieths S, Sticht H, et al. Allergic cross-reactivity made visible: solution structure of the major cherry allergen Pru av 1. *J Biol Chem*. 2001;276(25):22756-63.
28. Reuter A, Lidholm J, Andersson K, Ostling J, Lundberg M, Scheurer S, et al. A critical assessment of allergen component-based in vitro diagnosis in cherry allergy across Europe. *Clin Exp Allergy*. 2006;36(6):815-23.
29. Ballmer-Weber BK, Scheurer S, Fritsche P, Enrique E, Cistero-Bahima A, Haase T, et al. Component-resolved diagnosis with recombinant allergens in patients with cherry allergy. *J Allergy Clin Immunol*. 2002;110(1):167-73.
30. Asero R, Mistrello G, Roncarolo D, de Vries SC, Gautier MF, Ciurana CL, et al. Lipid transfer protein: a pan-allergen in plant-derived foods that is highly resistant to pepsin digestion. *Int Arch Allergy Immunol*. 2001 Jan-Mar;124(1-3):67-9.
31. Hartz C, Lauer I, del Mar San Miguel Moncin M, Cistero-Bahima A, Foetisch K, Lidholm J, Vieths S, et al. Comparison of IgE-binding capacity, cross-reactivity and biological potency of allergenic non-specific lipid transfer proteins from peach, cherry and hazelnut. *Int Arch Allergy Immunol*. 2010;153(4):335-46.
32. IUIS Allergen Nomenclature Sub-Committee Allergen Nomenclature. Available from: <http://www.allergen.org/viewallergen.php?aid=552>. Accessed June 2012.
33. Vallance G, et al. House dust mite control measures in the treatment of asthma. *Ther Clin Risk Manag*. 2006;2(4):347-54.
34. Wong C, et al. House dust mite Der p 1 elevates the release of inflammatory cytokines and expression of adhesion molecules in co-culture of human eosinophils and bronchial epithelial cells. *Int Immunol*. 2006;8(8):1327-35.
35. Brown A, et al. House dust mite Der p 1 down regulates defenses of the lung by inactivating elastase inhibitors. *Am J Respir Cell Mol Biol*. 2003;29:381-389.
36. Trombone A, et al. Use of a chimeric ELISA to investigate immunoglobulin E antibody responses to Der p 1 and Der p 2 in mite-allergic patients with asthma, wheezing and/or rhinitis. *Clin Exp Allergy*. 2002 Sep;32(9):1323-8.
37. Thomas W, Smith W. Towards defining the full spectrum of important house dust mite allergens. *Clin Exp Allergy*. 1999 Dec;29(12):1583-7.
38. Pajno G, et al. Prevention of new sensitizations in asthmatic children monosensitized to house dust mite by specific immunotherapy. A six year follow-up study. *Clin Exp Allergy*. 2001 Sep;31(9):1392-7.
39. Aki T, et al. Immunochemical characterization of recombinant and native tropomyosins as a new allergen from the house dust mite, *Dermatophagoides farinae*. *J Allergy Clin Immunol*. 1995 Jul;96(1):74-83.
40. Jahn-Schmid B, Kelemen P, Himly M, Bohle B, Fischer G, Ferreira F, et al. The T cell response to Art v 1, the major mugwort pollen allergen, is dominated by one epitope. *J Immunol*. 2002;169:6005-11.
41. Gadermaier G, Wopfner N, Wallner M, Egger M, Didierlaurent A, Regl G, et al. Array-based profiling of ragweed and mugwort pollen allergens. *Allergy*. 2008 Nov;63(11):1543-9.
42. Himly M, Jahn-Schmid B, Dedic A, Kelemen P, Wopfner N, Altmann F, et al. Art v 1, the major allergen of mugwort pollen, is a modular glycoprotein with a defensin-like and hydroxyprolinerich domain. *FASEB J*. 2003 Jan;17(1):106-8.
43. Jimeno L, Duffort O, Serrano C, Barber D, Polo F. Monoclonal antibody-based ELISA to quantify the major allergen of *Artemisia vulgaris* pollen, Art v 1. *Allergy*. 2004 Sep;59(9):995-1001.
44. Leonard R, Wopfner N, Pabst M, Stadlmann J, Petersen B, Duus J, et al. A new allergen from ragweed (*Ambrosia artemisiifolia*) with homology to Art v 1 from Mugwort. *J Biol Chem*. 2010 Aug;285:27192-200.
45. Wopfner N, Gadermaier G, Egger M, Asero R, Ebner C, Jahn-Schmid B, et al. The spectrum of allergens in ragweed and mugwort pollen. *Int Arch Allergy Immunol*. 2005;138:337-46.
46. Liccardi G, D'Amato G. Oleaceae pollinosis: a review. *Int Arch Allergy Immunol*. 1996;111(3):210-7.
47. Wheeler AW. Review. Hypersensitivity to allergens of the pollen from the olive tree (*Olea europaea*). *Clin Exp Allergy*. 1992;22:1052-7.
48. Barral P, et al. A major allergen from pollen defines a novel family of plant proteins and shows intra- and interspecies cross-reactivity. *J Immunol*. 2004 Mar;172(6):3644-5.
49. Batanero E, et al. Glycosylation site for the major allergen from olive tree pollen. Allergenic implications of the carbohydrate moiety. *Mol Immunol*. 1994 Jan;31(1):31-7.
50. Obispo T, et al. The main allergen of *Olea europaea* (Ole e 1) is also present in other species of the Oleaceae family. *Clin Exp Allergy*. 1993;23(4):311-6.
51. van Ree R, et al. A sensitive monoclonal antibody sandwich ELISA for the measurement of the major olive pollen allergen Ole e 1. *Int Arch Allergy Immunol*. 2000 Jul;122(3):224-8.
52. Rodriguez R, et al. Allergenic diversity of the olive pollen. *Allergy*. 2002;57:6-16.
53. Vourdas D, et al. Double-blind, placebo-controlled evaluation of sublingual immunotherapy with standardized olive pollen extract in pediatric patients with allergic rhinoconjunctivitis and mild asthma due to olive pollen sensitization. *Allergy*. 1998;53:662-72.
54. Rodriguez R, et al. The spectrum of olive pollen allergens. *Int Arch Allergy Immunol*. 2001;125:185-95.
55. Metz-Favre C, Pauli G, Bessot JC, De Blay F. Molecular allergology in practice: an unusual case of LTP allergy. *Eur Ann Allergy Clin Immunol*. 2011 Dec;43(6):193-5.
56. Asero R. Lipid transfer protein cross-reactivity assessed in vivo and in vitro in the office: pros and cons. *J Investig Allergol Clin Immunol*. 2011;21(2):129-36.
57. Asero R. Plant food allergies: a suggested approach to allergen-resolved diagnosis in the clinical practice by identifying easily available sensitization markers. *Int Arch Allergy Immunol*. 2005 Sep;138(1):1-11. Epub 2005 Aug 5.
58. Vieira T, Lopes C, Pereira AM, Araújo L, Moreira A, Delgado L. Microarray based IgE detection in poly-sensitized allergic patients with suspected food allergy—an approach in four clinical cases. *Allergol Immunopathol (Madr)*. 2012 May;40(3):172-80.
59. Gamboa PM, Sanz ML, Lombardero M, Barber D, Sánchez-Monje R, Goikoetxea MJ, Antépara I, Ferrer M, Salcedo G. Component-resolved in vitro diagnosis in peach-allergic patients. *J Investig Allergol Clin Immunol*. 2009;19(1):13-20.
60. 510(k) Substantial Equivalence Determination Decision Summary Assay Only Template. [Accessed 3-19-2022]. Available from: https://www.accessdata.fda.gov/cdrh_docs/reviews/K112523.pdf
61. 510(k) Substantial Equivalence Determination Decision Summary Assay Only Template. [3-18-2022]. Available from: https://www.accessdata.fda.gov/cdrh_docs/reviews/K132801.pdf
62. 510(k) Substantial Equivalence Determination Decision Summary Assay Only Template. [Accessed 3-19-2022]. Available from: https://www.accessdata.fda.gov/cdrh_docs/reviews/K093987.pdf

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. By constantly bringing breakthrough innovations to market, we enable healthcare professionals to deliver high-quality care, leading to the best possible outcome for patients.

Our portfolio, spanning from in-vitro and in-vivo diagnostics to image-guided therapy and innovative cancer care, is crucial for clinical decision-making and treatment pathways. With our strengths in patient twinning, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the biggest challenges in healthcare. We will continue to build on these strengths to help fight the world's most threatening diseases, improving the quality of outcomes, and enabling access to care.

We are a team of 66,000 highly dedicated employees across more than 70 countries passionately pushing the boundaries of what's possible in healthcare to help improve people's lives around the world.

3gAllergy, IMMULITE, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Laboratory Diagnostics
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000