

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\20min_protocols\hip\hip\pd_tse_sag

TA: 5:40

PAT: Off

Voxel size: 0.5×0.5×2.5 mm

Rel. SNR: 1.00

SIEMENS: tse

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	T > S40.5 > C0.1
Phase enc. dir.	A >> P
Rotation	17.00 deg
Slice group 3	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	T > S37.5 > C-15.7
Phase enc. dir.	A >> P
Rotation	32.00 deg
Slice group 4	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	T > C-36.9 > S29.4
Phase enc. dir.	H >> F
Rotation	-44.00 deg
Slice group 5	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	C > T-39.4 > S22.0
Phase enc. dir.	A >> P
Rotation	-49.00 deg
Slice group 6	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	C > T-26.2 > S14.7
Phase enc. dir.	H >> F
Rotation	-48.00 deg
Slice group 7	
Slices	1
Dist. factor	50 %
Position	R69.0 P2.2 F0.3
Orientation	C > T-12.0 > S4.9

Phase enc. dir.	R >> L
Rotation	-40.00 deg
Phase oversampling	100 %
FoV read	150 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	2100 ms
TE	22 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	BO1,2

Contrast

MTC	Off
Magn. preparation	None
Flip angle	140 deg
Fat suppr.	None
Water suppr.	None
Restore magn.	On

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

PAT mode	None
Matrix Coil Mode	Auto (CP)

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.	None

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off

System

Body	Off
------	-----

SIEMENS MAGNETOM TrioTim syngo MR B17

BO1	On
BO2	On

Positioning mode	ISO
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R69.0 P2.2 F0.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	212 mm
F >> H	208 mm
A >> P	188 mm

Physio

1st Signal/Mode	None

Dark blood	Off

Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	1
Bandwidth	180 Hz/Px
Flow comp.	No
Allowed delay	30 s
Echo spacing	11.2 ms

Define	Turbo factor
Turbo factor	4
Echo trains per slice	160
RF pulse type	Normal
Gradient mode	Fast