

The Only Fully-Automated Cyclosporine Assays for Enhanced Productivity

The need: To monitor cyclosporine levels in transplant patients

Cyclosporine A (CSA) is one of the primary immunosuppressive agents used with patients who have undergone cardiac, renal, and liver transplants. Toxic side effects, a narrow window of therapeutic efficacy, and the need to maintain consistent drug levels to prevent organ rejection all require frequent CSA monitoring. Although there is no universally established therapeutic range for CSA in whole blood, specific therapeutic ranges have been developed for each organ type and by each transplant institution.

Siemens Healthcare Diagnostics Dimension® CSA and CSAE Assays Meet the Need

The two Dimension cyclosporine assays are the only automated cyclosporine assays on the market, requiring no manual pretreatment steps. This feature will increase your productivity and turn-around time and reduce the potential for errors caused by manual sample splitting.

Add Cyclosporine testing to your routine

With Dimension CSA and CSAE, results are available around the clock, thanks to immediate testing made possible by the direct use of whole blood. Results can be reported rapidly with continuous bar code identification and online result transmission.

Dimension Cyclosporine (CSA) and Cyclosporine Extended Range (CSAE) Assays

Answers for life.



Dimension Cyclosporine (CSA) and Cyclosporine Extended Range (CSAE) Assays

Cyclosporine Flex® Reagent Cartridge

Mnemonic	CSA								
Product Code	DF89								
Packaging	80 Tests/Kit, 20 Tests/Flex								
Calibration	Dimension CSA Calibrator, Cat. No. DC89 Calibrator levels: 0, 80, 180, 330, 500 ng/mL Calibration frequency: 30 days or new reagent lot								
Pretreatment	No manual sample pretreatment – the sample is mixed on board by ultrasound, an aliquot taken, pretreated and analyzed fully automated.								
Quality Control	BioRad or MORE								
Stability	Closed on board (at 2-8° C): 30 days Open well stability: 3 days								
Assay Principle	Magnetic particle immunoassay (ACMIA)								
Type of Measurement	Photometric								
Sample Type	EDTA Whole Blood (200 μL/cup, 5 μL/test)								
Reportable Range	25-500 ng/mL								
Sensitivity Limit	Analytical: 25 ng/mL Functional: Not supplied								
Precision	Control Level	With	in-run	Total Precision					
	(ng/mL)	SD	%CV	SD	%CV				
	108	5.3	4.9	12.2	11.3				
	228	11.1	4.9	15.9	7.0				
	328	17.8	5.4	25.6	7.8				
Method	HPLC	Slope	Int	R	n				
Comparison	All	1.19	-13.7	0.935	661				
Companison	Heart	1.28	-16.0	0.893	111				
	Liver	1.18	-13.3	0.961	201				
	Kidney	1.13	-8.9	0.920	273				
	Abbott TDx*								
	All	0.83	-20.7	0.957	667				
	Heart	0.83	-20.9	0.927	111				
	Liver	0.82	-26.3	0.952	201				
	Kidney	0.83	-18.8	0.966	279				
Specificity	For specific cross-reactivity information,								

Cyclosporine (CSAE) Flex® Reagent Cartridge

•	•		3		9			
Mnemonic	CSAE							
Product Code	DF108							
Packaging	80 Tests/Kit, 20 Tests/Flex							
Calibration	Dimension CSAE Calibrator, Cat. No. DC108A Calibrator levels: 200, 400, 800, 1400, 2000 ng/mL Calibration frequency: 30 days or new reagent lot							
Pretreatment	No manual sample pretreatment – the sample is mixed on board by ultrasound, an aliquot taken, pretreated and analyzed fully automated.							
Quality Control	BioRad or MORE							
Stability	Closed on board (at 2-8° C): 30 days Open well stability: 3 days							
Assay Principle	Magnetic particle immunoassay (ACMIA)							
Type of Measurement	Photometric							
Sample Type	EDTA Whole Blood (200 μL/cup, 3 μL/test							
Reportable Range	350-2000 ng/mL							
Sensitivity Limit	Analytical: Not supplied Functional: <350 ng/mL							
Precision	Control Level							
	(ng/mL)	SD	%CV	SD	%CV			
	488.7	13.99	2.86	28.77	5.89			
	866.1	18.61	2.15	47.71	5.5			
	1301.9	34.41	2.64	68.89	5.29			
Method	HPLC	Slope	Int	R	n			
Comparison	All	1.09	17.1	0.986	140			
	Heart	0.93	93.7	0.983	35			
	Liver	1.00	59.9	0.980	40			
	Kidney	1.10	15.2	0.990	60			
	Abbott TDx®							
	All	1.13	-67.2	0.980	140			
	Heart	0.996	2.52		35			
	Liver	1.11	-61.7	0.970	40			
	Kidney	1.09	-32.3	0.973	60			
Specificity	For specif							

For more information about Cyclosporine Assays, contact your local sales representative or call **1-800-242-3233, option 3**.

see Instructions for Use.

References:

- 1. Hardman, JG, Limbird LE, Malinoff PB, Ruddon RW. <u>The Pharmacological Basis of Therapeutics</u>. 9th ed. New York, NY: McGraw Hill; 1996.
- 2. Kaplan LA, Pesce AJ. <u>Clinical Chemistry-Theory</u>, <u>Analysis</u>, <u>Correlation</u>. 3rd ed. St. Louis, MO: Mosby; 1996.

Siemens Healthcare Diagnostics Inc. 1717 Deerfield Road Deerfield, IL 60015-0778 IISA Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

© 2008 Siemens Healthcare Diagnostics Inc.
Order No. A91DX-0701146-C1-4A00
Printed in USA
Dimension, Flex and all associated marks are
trademarks of Siemens Healthcare Diagnostics
Inc. All other trademarks and brands are

the property of their respective owners.