

SIEMENS

VSIM DICOM Conformance Statement



Revision A • 12/02

58 57 698 OTH 2523 English

Siemens Medical Solutions USA, Inc.
Oncology Care Systems Group
4040 Nelson Avenue
Concord, CA 94520

The software with the VSIM is furnished under a license or nondisclosure agreement. The software may be used or copied only in accordance with the terms of the agreement. It is against the law to copy the software on any medium except as specifically allowed in the license or nondisclosure agreement.

Information in this manual is subject to change without notice. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording, for any purpose without the express written permission of Siemens Medical Solutions USA, Inc.

For reasons of product safety, use only original accessories made by Siemens or supplied by Siemens. The user assumes full responsibility and risk if third-party accessories not supplied by Siemens are used.

VSIM is a trademark of Siemens Medical Solutions USA, Inc.

All other products or services mentioned in this manual are identified by the trademarks or service marks of their respective companies or organizations. Siemens Medical Solutions USA, Inc., disclaims any responsibility for specifying which marks are owned by which companies or organizations.



CAUTION: RESTRICTED USE DEVICE

Federal law restricts use of this device to or on the order of a therapeutic radiologist.

VSIM

DICOM Conformance Statement


© 2002 by Siemens Medical Solutions USA, Inc.

All rights reserved.

58 57 698

Revision A

Software Version 1.0



EEC Declaration of Conformity

Application of Council Directive: Medical Device Directive 93/42/EEC

Standards to which conformity is declared:

EN60601-1, EN60601-1-1, EN60601-1-2, IEC60601-1-4, **IEC60601-2-1**, EN55011, EN61000-4-2, EN61000-4-3, EN61000-4-4, EN61000-4-5, EN61000-4-6

Manufacturer's name: Siemens Medical Solutions USA, Inc., Oncology Care Systems Group

Manufacturer's address: 4040 Nelson Avenue, Concord, CA, USA 94520

Type of equipment: Radiotherapy equipment

Note Any modification of the named device without written authorization by Siemens Medical Solutions USA, Inc., will invalidate this declaration.

Contents

Purpose	1
Structure	1
1 Introduction	5
Purpose	5
Scope	5
Definitions, Abbreviations	5
Definitions	5
Abbreviations	5
References	6
Connectivity and Inter-operability	6
2 Implementation Model Verification	9
Application Data Flow Diagram	9
Functional Definitions of Application Entities	9
Sequencing of real World Activities	9
3 AE Specification Verification	11
Verification AEs Specification	11
Association Establishment Policies	11
Association Initiation Policy	11
Association Acceptance Policy	12
4 Implementation Model Storage	13
Application Data Flow Diagram	13
Functional Definitions of Application Entities	14
Sequencing of real World Activities	14
5 AE Specification Storage	15
Storage AEs Specification	15
.....	18
Association Establishment Policies	18
Association Initiation Policy	18
Association Acceptance Policy	26

6	Implementation Model Storage Commit	37
	Application Data Flow Diagram	37
	Functional Definitions of Application Entities	38
	Sequencing of real World Activities	38
7	AE Specification Storage Commitment	39
	Association Establishment Policies	39
	Association Initiation Policy	40
	Association Acceptance Policy	42
8	Implementation Model Query/Retrieve	45
	Application Data Flow Diagram	45
	Functional Definitions of Application Entities	46
	Sequencing of real World Activities	46
9	AE Specification Query/Retrieve	47
	Query/Retrieve Service AEs Specification	47
	Association Establishment Policies	48
	Association Initiation Policy	49
	Association Acceptance Policy	55
10	Implementation Model Print	71
	Application Data Flow Diagram	71
	Functional Definitions of Application Entities	71
	Sequencing of real World Activities	72
11	AE Specification Print	73
	Print Management AE Specification	73
	Association Establishment Policies	74
	Association Initiation Policy	74
12	Implementation Model Worklist	87
	Application Data Flow Diagram	87
	Functional Definitions of Application Entities	87
	Sequencing of real World Activities	88
13	AE Specification Worklist	89
	Modality Worklist Service AEs Specification	89
	Association Establishment Policies	89
	Association Initiation Policy	90

14	Implementation Model MPPS	99
	Application Data Flow Diagram	99
	Functional Definitions of Application Entities	99
	Sequencing of real World Activities.....	100
15	AE Specification MPPS	101
	Modality Performed Procedure Step AE Specification	101
	Association Establishment Policies	101
	Association Initiation Policy	102
16	Communication Profiles	111
	Supported Communication Stacks	111
	TCP/IP Stack	111
17	Extensions/Specializations/Privatizations	113
	Standard Extended/Specialized/Private SOPs	113
	Private Transfer Syntaxes	113
18	Configuration	115
	AE Title / Presentation Address Mapping	115
	DICOM Verification.....	115
	DICOM Storage AE Title.....	115
	DICOM Query/Retrieve AE Title.....	115
	DICOM Print AE Title.....	115
	Configurable Parameters	116
	Storage, Storage Commitment and Query/Retrieve.....	116
	Print	116
	Modality Worklist.....	117
	Default Parameters.....	117
19	Support of Extended Character Sets.....	119
20	Introduction	123
	Purpose	123
	Scope.....	123
	Definitions, Abbreviations.....	123
	Definitions	123
	Abbreviations.....	123
	References	124
	Connectivity and Inter operability	124

21	Implementation Model	127
	Application Data Flow Diagram.....	127
	Functional definitions of AE's	128
	Sequencing of Real World Activities.....	128
	File Meta Information Options	128
22	AE Specifications.....	129
	DICOM Archive Specification	129
	File Meta Information for the Application Entity.....	130
	Real-World Activities for this Application Entity.....	130
23	Augmented and Private Profiles	135
	Augmented Application Profiles	135
	AUG-GEN-CD.....	135
	AUG-CTMR-xxxxx	135
	AUG-XA1K-CD	135
	Private Application Profiles	135
24	Extensions, Specializations and Privatization of SOP Classes and Transfer Syntaxes	137
	SOP Specific Conformance Statement for Basic Directory	137
	Extension, Specialization for SIEMENS Non-Image Objects	137
25	Configuration	139
	AE Title Mapping	139
	DICOM Media Storage AE Title.....	139
26	Support of Extended Character Sets.....	141
27	Clinical Context.....	145
	Roles and Service Class Options	145
	File Set Creator.....	145
	File Set Reader	146
	File Set Updater.....	146
28	PRI-SYNGO Profiles.....	147
	SOP Classes and Transfer Syntaxes.....	147
	Physical Media and Media Formats.....	156
	Directory Information in DICOMDIR	156
	Basic Directory IOD Specialization.....	157
	Additional Keys	157
	Private Directory Record Keys.....	159

Icon Images	160
Other Parameters	160
Multiframe JPEG Format.....	160
Attribute Value Precedence.....	160
A Appendix	163
SIEMENS Private Non-Image IOD	163
Siemens Non-Image IOD - E-R Model	163
Siemens Non-Image IOD - Module Table	165
Siemens Non-Image IOD - Modules.....	165
Siemens Standard Extended Modules	169
CSA Image Header Module	169
.....	170
CSA Series Header Module.....	170
.....	170
MEDCOM Header Module	170
.....	173
MEDCOM OOG Module	173
syngo Report Data.....	174
.....	175
syngo Report Info.....	175
Registry of DICOM Data Elements.....	176
Standard Extensions of all SOP Classes.....	177
Image Type	179
Patient Position.....	182
DICOM Print SCU - detailed status displays.....	182
Common Status Information	183
Additional Status Information - AGFA Printers	192
.....	192
Additional Status Information - Kodak PACS Link (formerly Imation)	192
Additional Status Information - Kodak 190I.....	193
Additional Status Information - Kodak 2180/1120	194
Additional Status Information - Codonics	194
Additional DICOM Execution Status Information	197
RT IOD specific implementation details	197



List of Figures

Figure 1. Application Data Flow Diagram - Verification SCU 9

Figure 2. Application Data Flow Diagram - Storage SCU 13

Figure 3. Application Data Flow Diagram - Storage SCP 13

Figure 4. Application Data Flow Diagram Storage Commitment SCU 37

Figure 5. Application Data Flow Diagram Storage Commitment SCP 38

Figure 6. Application Data Flow Diagram - Query/Retrieve SCU 45

Figure 7. Application Data Flow Diagram - Query/Retrieve SCP 46

Figure 8. Application Data Flow Diagram - Print SCU 71

Figure 9. Application data flow diagram - Basic Worklist SCU 87

Figure 10. Application data flow diagram - MPPS SCU 99

Figure 11. DICOM Information Model 164

List of Tables

Table 1. Initiation presentation context - Verification	12
Table 2. SOP Classes as Storage SCU	15
Table 3. Private SOP Classes as Storage SCU	16
Table 4. SOP Classes as Storage SCP	16
Table 5. Private SOP Classes as Storage SCP	18
Table 6. Initiation presentation context - Storage	19
Table 7. Acceptable presentation contexts - Storage	27
Table 8. Storage Commitment SOP Class as an SCU and SCP	39
Table 9. Proposed presentation contexts - request Storage commitment	41
Table 10. Proposed presentation contexts - send Storage commitment response	41
Table 11. Accepted presentation contexts - Receive Storage commitment request	42
Table 12. Accepted presentation contexts - Update Flags (receive commit response)	43
Table 13. Query/Retrieve SOP Classes as an SCU	47
Table 14. Query/Retrieve SOP Classes as an SCP	47
Table 15. Patient Root and Study Root query attributes	50
Table 16. C-FIND response status	52
Table 17. C-MOVE response status	54
Table 18. Patient level attributes, Patient Root Information Model	57
Table 19. Study level attributes, Patient Root Information Model	58
Table 20. Series level attributes, Patient Root Information Model	59
Table 21. Study level attributes, Study Root Information Model	60
Table 22. Image level attributes, Patient Root Information Model	60
Table 23. Series level attributes, Study Root Information Model	62
Table 24. Image level attributes, Study Root Information Model	63
Table 25. Patient instance level, Patient Study Only Information Model	63
Table 26. Study level attributes, Patient Study Only Information Model	64
Table 27. C-FIND SCP return status	65
Table 28. Presentation Contexts - Get SCP	66
Table 29. C-GET SCP return status	67
Table 30. Presentation Contexts - Move SCP	68
Table 31. C-MOVE return status	69
Table 32. Basic Grayscale Print Management Meta SOP Classes	73
Table 33. Basic Color Print Management Meta SOP Classes	73

Table 34. Presentation Contexts.....	75
Table 35. Used Basic Film Session N-CREATE-RQ attributes	77
Table 36. Attributes of N-DELETE-RQ on Basic Film Session SOP Class	77
Table 37. Basic Film Session SOP status	78
Table 38. Used Film Box N-CREATE-RQ attributes	79
Table 39. Attributes of N-DELETE-RQ on Basic Film Box SOP Class.....	80
Table 40. Basic Film Box SOP status.....	80
Table 41. Mandatory Basic Grayscale Image Box N-SET attributes.....	81
Table 42. Basic Grayscale Image Box SOP status	81
Table 43. Mandatory Basic Color Image Box N-SET attributes.....	82
Table 44. Basic Color Image Box SOP status	82
Table 45. Used Presentation LUT N-CREATE-RQ attributes	83
Table 46. Attributes of N-CREATE-RSP on Presentation LUT SOP Class.....	83
Table 47. Presentation LUT SOP status	84
Table 48. Used Printer N-EVENT report attributes.....	84
Table 49. Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes	84
Table 50. Used Print Job N-EVENT report attributes	85
Table 51. SOP Classes as an SCU	89
Table 52. Proposed presentation contexts	90
Table 53. Supported Broad Worklist Query Search Key Attributes	91
Table 54. Basic Worklist C-FIND-RSP Return Key Attributes	92
Table 55. Patient based “narrow query” Search Key Attributes	96
Table 56. C-FIND Response Status	97
Table 57. SOP Classes as an SCU	101
Table 58. Proposed presentation contexts	102
Table 59. Performed Procedure Step N-CREATE Attributes	103
Table 60. MPPS N-CREATE Response Status.....	106
Table 61. Performed Procedure Step N-SET Attributes.....	107
Table 62. MPPS N-CREATE Response Status.....	109
Table 63. Application profiles, Activities, and Roles for DICOM Archive	129
Table 64. STD-GEN-CD Supported SOP Classes for FSR role.....	131
Table 65. PRI-SYNGO Universal VSim Class Profiles.....	136
Table 66. Basic Directory Extension for Non-Image Objects	137
Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes	147
Table 68. PRI-SYNGO Additional DICOMDIR Keys.....	157
Table 69. PRI-SYNGO Additional PRIVATE Record Keys.....	159
Table 70. CSA Non-Image IOD Modules	164
Table 71. CSA Non-Image Module	165
Table 72. CSA Image IOD Modules.....	168
Table 73. CSA Image Header Module	169
Table 74. CSA Series Header Module	169

Table 75. MEDCOM Header Module	170
Table 76. MEDCOM History Information.....	172
Table 77. MEDCOM OOG Module.....	173
Table 78. syngo Report Data Module.....	174
Table 79. Standard Extensions of all SOP Classes	178
Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class.....	183
Table 81. Printer Status Infos: Additional Agfa printer status infos	192
Table 82. Printer Status Infos: Additional Kodak infos for Pacs Link (formerly Imation cameras).....	192
Table 83. Printer Status Infos: Additional Kodak infos for Kodak 190.....	193
Table 84. Printer Status Infos: Additional Kodak infos for 2180/1120.....	194
Table 85. Printer Status Infos: Additional Codonics infos	194
Table 86. Additional DICOM Execution Status Infos	197
Table 87. RT Series Table	199
Table 88. RT Plan IOD module Table	199
Table 89. RT Structureset IOD Module Table.....	200
Table 90. RT Image IOD Module Table.....	201
Table 91. RT Structure Set Module	202
Table 92. ROI Contour Module	207
Table 93. ROI Observations Module	208
Table 94. RT General Plan Module	209
Table 95. RT Patient Setup Module	210
Table 96. RT Beam Module	211
Table 97. RT Image Module	215

Introduction

Purpose

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

The DICOM Conformance Statement (DCS) is supplied here as a complete document for all Services including Offline Media Storage Services. The VSim DICOM network implementation acts as SCU and SCP for the DICOM Storage and Storage Commitment, as SCU and SCP for the Query/Retrieve service and as SCU for the DICOM Basic Print, Modality Worklist and Modality Performed Procedure Step service.

Structure

The VSim Conformance Statement is subdivided into three Parts which relate to individual documents needed to declare Conformance according to the requirements of "Part 2 - Conformance" of the DICOM Standard.

Those parts are:

- **DICOM Conformance Statement** for Network related Services
- A privately defined *syngo* **Application Profile**
- **Offline Media Conformance Statement** to support local archive media.

A general appendix follows these sections.

Network Conformance Statement

This part contains the Conformance declaration for the DICOM Network Services:

- Storage - User/Provider (includes Verification - User/Provider)
- Storage Commitment - User/Provider
- Query/Retrieve - User/Provider
- Basic Grayscale/Color Print - User
- Basic Worklist - User
- Modality Performed Procedure Step - User

Introduction

Purpose

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

The application described in this conformance statement are part of the SIEMENS VSim based on *syngo*® software¹. The VSim DICOM network implementation acts as SCU and SCP for the DICOM Storage, Storage Commitment and Query/Retrieve services and as SCU for the DICOM Print, DICOM Basic Worklist and Modality Performed Procedure Step Services. Verification is supported in SCU (only via Service environment) and SCP role.

Scope

This DICOM Conformance Statement refers to Siemens OCS VSim product using *syngo*® software applications.

The VSim target Operating System is Microsoft Windows 2000.

Definitions, Abbreviations

Definitions

DICOM	Digital Imaging and Communications in Medicine
DIMSE	DICOM Message Service Element
DIMSE-C	DICOM Message Service Element with Composite information objects

Abbreviations

ACR	American College of Radiology
AE	DICOM Application Entity
ASCII	American Standard Code for Information Interchange
Med-WS	Multimodality-Workstation
CSE	Customer Service Engineer
DB	Database

1. *syngo* is a registered trademark of Siemens AG

DCS	DICOM Conformance Statement
DSA	Digital Subtraction Angiography
IIDC	Image-intensifier Distortion Correction
IOD	DICOM Information Object Definition
ISO	International Standard Organization
NEMA	National Electrical Manufacturers Association
O	Optional Key Attribute
PDU	DICOM Protocol Data Unit
R	Required Key Attribute
RIS	Radiology Information System
RWA	Real-World Activity
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM server)
SOP	DICOM Service-Object Pair
U	Unique Key Attribute

References

[1] Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-15, 2000

Connectivity and Inter-operability

This Conformance Statement by itself does not guarantee successful inter-operability of Siemens equipment with non-Siemens equipment. The user (user’s agent) should be aware of the following issues:

- Inter-operability
Inter-operability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee inter-operability of Siemens equipment with non-Siemens equipment. It is the user’s responsibility to analyze thoroughly the application requirements and to specify a solution that integrates Siemens equipment with non-Siemens equipment.
- Validation
Siemens equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.

Where Siemens equipment is linked to non-Siemens equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation test will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.

- New versions of the DICOM Standard

The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. Siemens is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, Siemens reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-Siemens provider linking to Siemens equipment, also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into Siemens equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

Implementation Model Verification

The Siemens VSim DICOM Service Tool application requests Verification to verify the ability of a foreign DICOM application on a remote node to respond to DICOM messages.

Responding to Verification requests from remote nodes is handled by the Storage SCP application.

Application Data Flow Diagram

The VSim DICOM network implementation acts as SCU for the C-ECHO DICOM network service. The VSim target Operating System is Microsoft Windows 2000.




Figure 1. Application Data Flow Diagram - Verification SCU

Functional Definitions of Application Entities

The Siemens VSim DICOM Service Tool application opens an association when a "Test" of a remote application is requested during a configuration session. This can be done when entering new data for remote application configuration or to verify existing configuration data.

Sequencing of real World Activities

Newly entered data have to be saved first, before a "test" of these data is possible.

AE Specification Verification

Verification AEs Specification

Association Establishment Policies

General

The Siemens VSim DICOM Service Tool application attempts to open an association for verification request whenever the “Test” function is activated during network configuration of a remote DICOM application.

Number of Associations

The Siemens VSim DICOM Service Tool application initiates one association at a time to request verification.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens VSim DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.5.9.20000101

and an Implementation Version Name of

- “SIEMENS_SWFVB10A”

Association Initiation Policy

The Siemens VSim DICOM Service Tool application attempts to initiate a new association for

- DIMSE C-ECHO

service operations.

Associated Real-World Activity - Verification

Associated Real-World Activity -Request Verification “Test”

The associated Real-World activity is a C-ECHO request initiated by Service and Configuration SW environment whenever a “Test” is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than “Success” this will be indicated in the service environment and the association is closed.

Proposed Presentation Contexts

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Table 1. Initiation presentation context - Verification

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotia- tion
Name	UID	Name List	UID List		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1	SCU	None

SOP Specific Conformance - Verification SCU

The Application conforms to the definitions of the Verification SCU in accordance to the DICOM Standard.

Association Acceptance Policy

The Verification SCP is part of the Storage SCP - see section on page 26

Implementation Model Storage

The Siemens VSim DICOM Application Entity both originates associations for Storage of DICOM Composite Information Objects in Remote Application Entities and receives association requests for Storage from Remote Application Entities.

Application Data Flow Diagram

The VSim DICOM network implementation acts as SCU and SCP for the C-STORE DICOM network service and as SCP for the C-ECHO DICOM network service. The VSim target Operating System is Microsoft Windows 2000.




Figure 2. Application Data Flow Diagram - Storage SCU




Figure 3. Application Data Flow Diagram - Storage SCP

Functional Definitions of Application Entities

The Storage SCU is invoked by the job control interface that is responsible for processing network archival tasks. The job consist of data describing the composite image objects selected for storage and the destination. An association is negotiated with the destination application entity and the image data is transferred using the C-STORE DIMSE-Service. Status of the transfer is reported to the job control interface.

The Storage SCP component of the Siemens VSim DICOM application is operating as background server process. It is existing when the machine is powered on and waits for Storage association requests. Upon accepting an association with a negotiated Presentation Context it starts to receive the Composite Image Objects and imports them to local database. Verification requests will be processed and responded by Storage SCP component too.

Sequencing of real World Activities

not applicable.

AE Specification Storage

Storage AEs Specification

The VSim Storage service class user/service class provider applications use one AE when initiating/receiving associations to/from remote DICOM nodes.

SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as an SCU.

Table 2. SOP Classes as Storage SCU

SOP Class Name	SOP Class UID
CR Computed Radiography Object Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Information Object Storage	1.2.840.10008.5.1.4.1.1.2
DX Digital X-Ray Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1
DX Digital X-Ray Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1
MG Digital Mammography Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
MG Digital Mammography Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2
MR Image Information Object Storage	1.2.840.10008.5.1.4.1.1.4
NM (Nuclear Medicine) Image Information Object Storage	1.2.840.10008.5.1.4.1.1.20
PET Image Storage	1.2.840.10008.5.1.4.1.1.128
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
SC (Secondary Capture) Image Information Object Storage	1.2.840.10008.5.1.4.1.1.7
US Multi-frame Image Information Object Storage	1.2.840.10008.5.1.4.1.1.3.1

Table 2. SOP Classes as Storage SCU

SOP Class Name	SOP Class UID
US Image Information Object Storage	1.2.840.10008.5.1.4.1.1.6.1
X-Ray Angiographic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.1
X-Ray Radiofluoroscopic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.2
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1

SIEMENS VSim DICOM products provide Private Conformance to the following DICOM V3.0 conform private SOP Classes as an SCU:

Table 3. Private SOP Classes as Storage SCU

SOP Class Name	SOP Class UID
CSA Non-Image	1.3.12.2.1107.5.9.1

SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as SCP

:

Table 4. SOP Classes as Storage SCP

SOP Class Name	SOP Class UID
CR Computed Radiography Object Storage	1.2.840.10008.5.1.4.1.1.1
CT Image Information Object Storage	1.2.840.10008.5.1.4.1.1.2
DX Digital X-Ray Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1
DX Digital X-Ray Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1

Table 4. SOP Classes as Storage SCP

SOP Class Name	SOP Class UID
MG Digital Mammography Storage - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
MG Digital Mammography Storage - For Presentation	1.2.840.10008.5.1.4.1.1.1.2
MR Image Information Object Storage	1.2.840.10008.5.1.4.1.1.4
NM (Nuclear Medicine) Image Information Object Storage	1.2.840.10008.5.1.4.1.1.20
PET Image Storage	1.2.840.10008.5.1.4.1.1.128
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5
SC (Secondary Capture) Image Information Object Storage	1.2.840.10008.5.1.4.1.1.7
US Multi-frame Image Information Object Storage	1.2.840.10008.5.1.4.1.1.3.1
US Image Information Object Storage	1.2.840.10008.5.1.4.1.1.6.1
X-Ray Angiographic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.1
X-Ray Radiofluoroscopic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.2
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1
Verification	1.2.840.10008.1.1
Ultrasound Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6

Table 4. SOP Classes as Storage SCP

SOP Class Name	SOP Class UID
Ultrasound Multi-frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3

SIEMENS VSim DICOM products provide Private Conformance to the following DICOM V3.0 conform private SOP Classes as an SCP:

Table 5. Private SOP Classes as Storage SCP

SOP Class Name	SOP Class UID
CSA Non-Image	1.3.12.2.1107.5.9.1

Association Establishment Policies

General

The existence of a job queue entry with network destination or an internal trigger from processing a retrieve request will activate the DICOM Storage Application. An association request is sent to the destination AE and upon successful negotiation of a Presentation Context the transfer is started.

The default PDU size used is 28 KB.

Number of Associations

The Siemens VSim DICOM application initiates several associations at a time, one for each destination to which a transfer request is being processed in the active job queue list. There is no limit for the number of initiated associations.

The Siemens VSim DICOM application is able to accept multiple associations at a time. It can handle in parallel 10 associations, and a configurable number of association requests can be queued until one of the 10 associations is released.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Association Initiation Policy

If a job with network destination gets active in the job list or a retrieve sub-operation is processed, the Siemens VSim DICOM application attempts to initiate a new association for

- DIMSE C-STORE

service operations.

Associated Real-World Activity

Associated Real-World Activity -Send DICOM Objects to a Network destination

The associated Real-World activity is a C-STORE request initiated by an internal daemon process triggered by a job with network destination or the processing of an external C-MOVE retrieve request. If the process successfully establishes an association to a remote Application Entity, it will transfer each image one after another via the open association. If the C-STORE Response from the remote Application contains a status other than “Success” or “Warning” the association is aborted.

Proposed Presentation Contexts - Send Images

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table

:

Table 6. Initiation presentation context - Storage

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
CR Image Storage	1.2.840.10008.5.1.4.1.1.1	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None

Table 6. Initiation presentation context - Storage

DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1.1 1.2.840.10008.5.1.4.1.1.1.1	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2.1 1.2.840.10008.5.1.4.1.1.1.2	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
NM Image Storage	1.2.840.10008.5.1.4.1.1.20	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
PET Image Storage	1.2.840.10008.5.1.4.1.1.128	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None

Table 6. Initiation presentation context - Storage

RT Structure Set Storage	1.2.840.10008.5.1.4 .1.1.481.3	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
RT Dose Storage	1.2.840.10008.5.1.4 .1.1.481.2	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
RT Image Storage	1.2.840.10008.5.1.4 .1.1.481.1	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
RT Plan Storage	1.2.840.10008.5.1.4 .1.1.481.5	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
SC Image Storage	1.2.840.10008.5.1.4 .1.1.7	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
US Multi-frame Image Storage	1.2.840.10008.5.1.4 .1.1.3.1	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
US Image Storage	1.2.840.10008.5.1.4 .1.1.6.1	JPEG Lossy Extended *1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline *1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None

Table 6. Initiation presentation context - Storage

X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4 .1.1.12.1	JPEG Lossy Extended * 1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline * 1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
X-Ray Radiofluorosc opic Image Storage	1.2.840.10008.5.1.4 .1.1.12.2	JPEG Lossy Extended * 1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline * 1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
Waveform Storage	1.2.840.10008.5.1.4 .1.1.9.1.1 1.2.840.10008.5.1.4 .1.1.9.1.2 1.2.840.10008.5.1.4 .1.1.9.1.3 1.2.840.10008.5.1.4 .1.1.9.2.1 1.2.840.10008.5.1.4 .1.1.9.3.1 1.2.840.10008.5.1.4 .1.1.9.4.1	JPEG Lossy Extended * 1 (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline * 1 (Process 1) compression Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCU	None
CSA Non- Image Storage	1.3.12.2.1107.5.9.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1	SCU	None

***1** The Transfer Syntax used is strongly influenced by the fact of “how was the accepted Transfer Syntax at the time when the Instance was received”. e.g. the Instances received with JPEG Lossy Transfer Syntaxes will not be converted and can only be sent out with the same Transfer Syntax.

Note Uncompressed transfer syntaxes are proposed together in a single presentation context for each abstract syntax.

The “MOVE destinations” must be configured as Storage destinations. This would include the configuration of Transfer Syntax capabilities.

Not all the listed transfer syntaxes will be proposed all the time. For some abstract syntaxes only a list of uncompressed (UC) transfer syntaxes (one or more) will be proposed, for other abstract syntaxes also JPEG Lossless (LL) syntax will be proposed and/or a list of JPEG Lossy (LY) transfer syntaxes. The contents of this lists is configurable, e.g. UC could be configured to contain only Implicit Little Endian for instance.

Depending on the real world activity initiating the C-STORE, we have the following behaviors:

- if the C-STORE is initiated by a user, a configuration parameter called QualityFactor(Q) will be used to decide which transfer syntax lists will be proposed. Q can take values between 0 and 100. If Q=0, only UC will be proposed. If Q = 100, UC and LL will be proposed. Else UC and LY will be proposed.
- if the C-STORE is initiated by the C-MOVE SCP, there is another configuration parameter called Compression Types Supported (CTS) which will be used to decide what transfer syntaxes are proposed. CTS can take integer values. If CTS=0 or CTS > 3, UC will be proposed. If CTS=1, UC and LY will be proposed. If CTS = 2, UC and LL will be proposed. If CTS = 3, UC, LL and LY will be proposed.

SOP Specific Conformance to Storage SOP classes

The VSim will create SC IOD type images when performing special function that create Derived Images. The SC IOD will be a Standard Extended SC Storage SOP Class. The 3D viewing application is able to create further 2D (SC) images to document results from 3D processing.

The VSim (DICOM) application will not change private attributes as long as no modification is done. During a “Save Image” operation all private attributes not defined within the VSim DICOM application will be removed when the new object instance is created.

For association and DIMSE level time-outs, please refer to Configuration section of this document.

- Optional Attributes

Data Dictionary of DICOM Type 2 and 3 IOD Attributes

Please see the related Image Object definition tables in the Annex for a list of all DICOM IOD attributes of type 2 and 3 which are encoded in the VSim applications.

- Specialized Information Object Definitions

The DICOM images created by Siemens VSim DICOM application conform to the DICOM IOD definitions (Standard extended IODs). But they will contain additional private elements which have to be discarded by a DICOM system when modifying the image.

The DICOM nodes are responsible for data consistency when modifying images. All unknown private attributes have to be removed upon modification!

Data Dictionary of applied private IOD Attributes

Please see “Siemens Standard Extended Modules” in the Appendix for a list of possible private IOD attributes.

- Image Pixel Attribute Description for Grayscale Images

The Siemens VSim DICOM application supports the Monochrome2 Photometric Interpretation with the unsigned integer 16 bit gray scale pixel and graphic overlay format. The lower 12 bits are used for pixel and the higher 4 bits are used for the graphic overlay:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = “MONOCHROME1”
- photometric interpretation (attribute 0028,0004) = “MONOCHROME2”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8, 16
- bits stored (attribute 0028,0101) = 8, 10, 12
- high bit (attribute 0028,0102) = 7, 9, 11

Overlay plane

- overlay type (attribute 60xx, 0040) = “G”
- overlay bits allocated (attribute 60xx, 0100) = 16
- overlay bit position (attribute 60xx, 0102) = 12, 13, 14, 15

Overlay plane

- overlay type (attribute 60xx, 0040) = “G”
- bits allocated (attribute 60xx, 0100) = 1
- bit position (attribute 60xx, 0102) = 0
- overlay data (attribute 60xx, 3000) = supported.

The Siemens VSim DICOM application sends also the MONOCHROME1 and MONOCHROME2 photometric interpretation pixel format with binary 2’s complement integer and 16 bits allocated. Possible values:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = “MONOCHROME1”
- photometric interpretation (attribute 0028,0004) = “MONOCHROME2”
- pixel representation (attribute 0028, 0103) = 1
- bits allocated (attribute 0028, 0100) = 16
- bits stored (attribute 0028,0101) = 16
- high bit (attribute 0028,0102) = 15

Overlay plane

- overlay type (attribute 60xx, 0040) = “G”
- overlay bits allocated (attribute 60xx, 0100) = 1
- overlay bit position (attribute 60xx, 0102) = 0
- overlay data (attribute 60xx, 3000) = supported.
- Image Pixel Attribute Description for Color Images

The Siemens VSim DICOM application supports the RGB color image description with the unsigned integer 24 bit color image plane pixel format:

- samples per pixel (attribute 0028, 0002) = 3
- photometric interpretation (attribute 0028,0004) = “RGB”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8
- bits stored (attribute 0028,0101) = 8
- high bit (attribute 0028,0102) = 7
- planar configuration (attribute 0028,0006) = 0.

The Siemens VSim DICOM application supports the “Palette Color” color image description with unsigned integer and 2’s complement pixel format:

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = “PALETTE COLOR”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8, 16
- bits stored (attribute 0028,0101) = 8, 16
- high bit (attribute 0028,0102) = 7, 15

Both 8bit and 16bit palettes are supported - but no Segmented Palette Color LUTs.

The Siemens VSim DICOM application supports the YBR_FULL color image description with the unsigned integer 24 bit color image plane pixel format:

- samples per pixel (attribute 0028, 0002) = 3
- photometric interpretation (attribute 0028,0004) = “YBR_FULL”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8
- bits stored (attribute 0028,0101) = 8
- high bit (attribute 0028,0102) = 7
- planar configuration (attribute 0028,0006) = 0.

Association Acceptance Policy

The Siemens VSim DICOM application attempts to accept a new association for

- DIMSE C-ECHO
- DIMSE C-STORE

service operations. Any Information Objects transmitted on that association will be checked on conformance and stored in database if check was successful.

Associated Real-World Activity - Receive

- Associated Real-World Activity - Receiving Image Objects from a Remote Node

The daemon receiving process will accept an association and will receive any images transmitted on that association and will store the images on disk in the own data base if the conformance check is performed successfully.

- Presentation Context Table

The Siemens VSim DICOM application will accept Presentation Contexts as shown in the following table:

Table 7. Acceptable presentation contexts - Storage

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
CR Image Storage	1.2.840.10008.5.1.4.1.1.1	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
DX Image Storage	1.2.840.10008.5.1.4.1.1.1.1.1 1.2.840.10008.5.1.4.1.1.1.1	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None

Table 7. Acceptable presentation contexts - Storage

MG Image Storage	1.2.840.10008.5.1.4.1.1.1.2.1 1.2.840.10008.5.1.4.1.1.1.2	JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	SCP	None
		JPEG Lossless Process 14 (selection value 1)	1.2.840.10008.1.2.4.70		
		JPEG Lossy Baseline (Process 1) compression	1.2.840.10008.1.2.4.50		
		RLE Lossless	1.2.840.10008.1.2.5		
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Explicit VR Big Endian	1.2.840.10008.1.2.2	SCP	None
		Implicit VR Little Endian	1.2.840.10008.1.2		
		JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51		
		JPEG Lossless Process 14 (selection value 1)	1.2.840.10008.1.2.4.70		
		JPEG Lossy Baseline (Process 1) compression	1.2.840.10008.1.2.4.50		
NM Image Storage	1.2.840.10008.5.1.4.1.1.20	RLE Lossless	1.2.840.10008.1.2.5	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		Implicit VR Little Endian	1.2.840.10008.1.2		
		JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51		
PET Image Storage	1.2.840.10008.5.1.4.1.1.128	JPEG Lossless Process 14 (selection value 1)	1.2.840.10008.1.2.4.70	SCP	None
		JPEG Lossy Baseline (Process 1) compression	1.2.840.10008.1.2.4.50		
		RLE Lossless	1.2.840.10008.1.2.5		
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		Implicit VR Little Endian	1.2.840.10008.1.2		
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

Table 7. Acceptable presentation contexts - Storage

RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
SC Image Storage	1.2.840.10008.5.1.4.1.1.7	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
US Multi-frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
US Image Storage	1.2.840.10008.5.1.4.1.1.6.1	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None

Table 7. Acceptable presentation contexts - Storage

Ultrasound Image Storage (Retired) ^a	1.2.840.10008.5.1.4.1.1.6	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
Ultrasound Multi-frame Image Storage (retired) ^a	1.2.840.10008.5.1.4.1.1.3	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	JPEG Lossy Extended (Process 2 & 4) JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression RLE Lossless Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian	1.2.840.10008.1.2.4.51 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.5 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2	SCP	None

Table 7. Acceptable presentation contexts - Storage

Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2.4.51	SCP	None
	1.2.840.10008.5.1.4.1.1.9.1.2	JPEG Lossless Process 14 (selection value 1)	1.2.840.10008.1.2.4.70		
	1.2.840.10008.5.1.4.1.1.9.1.3	JPEG Lossy Baseline (Process 1) compression	1.2.840.10008.1.2.4.50		
	1.2.840.10008.5.1.4.1.1.9.2.1	RLE Lossless	1.2.840.10008.1.2.5		
	1.2.840.10008.5.1.4.1.1.9.3.1	Explicit VR Little Endian	1.2.840.10008.1.2.1		
	1.2.840.10008.5.1.4.1.1.9.3.1	Explicit VR Big Endian	1.2.840.10008.1.2.2		
	1.2.840.10008.5.1.4.1.1.9.4.1	Implicit VR Little Endian	1.2.840.10008.1.2		
CSA Non-Image Storage	1.3.12.2.1107.5.9.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

a. US Retired and US Multiframe Retired images are converted to US Images/US Multiframe images before storing them into the local database. The conversion creates new images, which implies new UIDs.

- **SOP Specific Conformance Statement - Receiving Images**

The Siemens VSim DICOM application conforms to the Full Storage Class at Level 2.

Upon successful receiving a C-STORE-RQ, the Siemens VSim DICOM receiver performs a quick plausibility test on the received image and of the own resources. If this test succeeds, it returns the status SUCCESS, otherwise one of the following status codes is returned:

- **Refused (A700):**
This error status indicates a lack of Resources (e.g. not enough disk space) on the VSim modality.
- **Error (A900 or C000):**
An error occurred while processing the image which makes it impossible to proceed. The image will not be stored and the association aborted.

ATTENTION!

Only after sending the response, the image will be saved into the database. If during this operation an error occurs, the association will be aborted. This implies that a C-STORE-RSP with status SUCCESS does not mean that the image was successfully stored into the database.

In order to be sure that the sent images were successfully stored in the database, the sending application should use Storage Commitment.

If an image instance is received that is identified by an SOP Instance UID that is already used by an Instance stored in database then the actual received image will be discarded. The existing Instance is not superseded.

The following sections will differentiate the attribute contents required for Image Viewing. The Siemens VSim DICOM application supports more formats for Storage of images.

Image Pixel Attribute Acceptance Criterion
for Grayscale Images - Viewing

The VSim accepts the MONOCHROME1 and MONOCHROME2 photometric interpretation pixel format and graphic overlay with unsigned integer and 8 or 16 bits allocated. Accepted values:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = “MONOCHROME1”
- photometric interpretation (attribute 0028,0004) = “MONOCHROME2”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8, 16
- bits stored (attribute 0028,0101) = 8, 10, 12
- high bit (attribute 0028,0102) = 7, 9, 11
- only aspect ratio 1:1 is allowed

Overlay plane

- overlay type (attribute 60xx, 0040) = “G”
- bits allocated (attribute 60xx, 0100) = 16
- bit position (attribute 60xx, 0102) = 12, 13, 14, 15 (only bits above high bit permitted)

Graphic Overlay will be shifted to fill Overlay Planes from Bit 12 and consecutive.

Overlay plane

- overlay type (attribute 60xx, 0040) = “G”
- bits allocated (attribute 60xx, 0100) = 1
- bit position (attribute 60xx, 0102) = 0
- overlay data (attribute 60xx, 3000) = supported

The Siemens VSIM Image Viewing application accepts also the MONOCHROME1 and MONOCHROME2 photometric interpretation pixel format with binary 2's complement integer and 16 bits allocated. Accepted values:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = "MONOCHROME1"
- photometric interpretation (attribute 0028,0004) = "MONOCHROME2"
- pixel representation (attribute 0028, 0103) = 1 (signed)
- bits allocated (attribute 0028, 0100) = 16
- bits stored (attribute 0028,0101) = 16
- high bit (attribute 0028,0102) = 15
- only aspect ratio 1:1 is allowed

Overlay plane

- overlay type (attribute 60xx, 0040) = "G"
- bits allocated (attribute 60xx, 0100) = 1
- bit position (attribute 60xx, 0102) = 0
- overlay data (attribute 60xx, 3000) = supported

For MOD LUT, both the linear LUT (Rescale Slope/Intercept) and the MOD LUT SQ are supported and considered when pixel data is displayed. However there are two limitations. The MOD LUT SQ will be ignored in the following cases:

- 8-Bit signed pixels
- the pixel format is changed by the MOD LUT (e.g. 8bit -> 16bit)

If the MOD LUT SQ contains multiple LUTs, then only the first one is used.

For VOI LUT, both the linear LUT (Window Center/Width) and the VOI LUT SQ are supported (VOI LUT SQ with 8 or 16 bit LUT data)

But if both, a VOI LUT SQ and a linear MOD LUT, are specified within one image, then the value for Rescale Slope is restricted to 1.

If the VOI LUT SQ contains multiple LUTs, then only the first one is used.

Only Rectangular and Circular Shutter Shape is supported in this version. Image with other Shutter Shapes will be displayed w/o shutter.

Image Pixel Attribute Acceptance Criterion for Color Images - Viewing

The Siemens VSim Multi-Modality Viewing application supports the RGB color image description with the unsigned integer 24 bit color image plane pixel format. Accepted values:

- samples per pixel (attribute 0028, 0002) = 3
- photometric interpretation (attribute 0028,0004) = “RGB”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8
- bits stored (attribute 0028,0101) = 8
- high bit (attribute 0028,0102) = 7
- planar configuration (attribute 0028,0006) = 0 (pixel interleave).

The Siemens VSim Multi-Modality Viewing application supports the “Palette Color” image description with the unsigned integer and 2’s complement pixel format.

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = “PALETTE COLOR”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8 and bits stored (attribute 0028,0101) = 8
- bits allocated (attribute 0028, 0100) = 16 and bits stored (attribute 0028,0101) = 16
- high bit (attribute 0028,0102) = 7,15

Both 8-Bit and 16-Bit palettes are supported, but NO Segmented Palette Color LUTs.

The Siemens VSim Multi-Modality Viewing application supports the YBR_FULL image description with the unsigned integer and 2’s complement pixel format.

- samples per pixel (attribute 0028, 0002) = 3
- photometric interpretation (attribute 0028,0004) = “YBR_FULL”
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8 and bits stored (attribute 0028,0101) = 8
- high bit (attribute 0028,0102) = 7

If Siemens VSim software is making any persistent changes on a YBR image, the resulting new image will be saved with Photometric Interpretation = “RGB”.

- Presentation Context Acceptance Criterion

The Siemens VSim DICOM application will accept any number of verification or storage SOP classes that are listed above. There is no limit on the number of presentation contexts accepted, except for the DICOM limit. In the event that the Siemens VSim DICOM application runs out of resources, it will reject the association request.

- Transfer Syntax Selection Policies

The Siemens VSim DICOM application currently supports

- the Implicit VR Little Endian, the Explicit VR Little Endian and Explicit VR Big Endian Transfer Syntaxes
- the JPEG Lossless Non-hierarchical Transfer Syntax
- the JPEG Baseline and JPEG Extended Transfer Syntaxes (JPEG Lossy).
- the RLE Lossless Transfer Syntax

Any proposed presentation context which includes one of these transfer syntaxes will be accepted. Any proposed presentation context that does not include one of these transfer syntaxes will be rejected.

The order of preference in accepting Transfer Syntaxes within Presentation Contexts or Presentation Contexts with single Transfer Syntaxes is:

- 1 JPEG Lossy Extended
- 2 JPEG Lossless non-hierarchical
- 3 JPEG Lossy Baseline
- 4 RLE Lossless
- 5 Explicit VR Little Endian
- 6 Explicit VR Big Endian
- 7 Implicit VR Little Endian

With Implicit VR Little Endian Transfer Syntax the VSim DICOM application will remove any Private Attributes not known to the application. Decision on removal of a Private Element is done if there is NO entry in the dictionary of attributes of the VSim DICOM application.

With RLE Lossless Transfer Syntax the VSim DICOM application will decompress the image before storing it into the database.

Therefore any Explicit VR Transfer Syntax shall preferably be used by the Storage SCU's when sending Composite Image Instances to the VSim DICOM application.

Implementation Model Storage Commit

The Storage Commitment service class defines an application-level class of service which facilitates the commitment to storage. It performs an additional task of commitment of composite objects apart from the network based storage of images as defined by the Storage Service class. The VSim DICOM implementation supports the Storage Commitment Push Model as SCU and SCP.

Application Data Flow Diagram

The VSim DICOM network implementation acts as SCU/SCP for the Storage Commitment Push Model Service using the Storage Commitment Service Class. The VSim target Operating System is Microsoft Windows 2000.




Figure 4. Application Data Flow Diagram Storage Commitment SCU




Figure 5. Application Data Flow Diagram Storage Commitment SCP

Functional Definitions of Application Entities

With each successfully completed send job, the VSim DICOM Application will create a Storage Commitment Push Model Identifier from the SOP Instances sent. Then an a Storage Commit Request is triggered. Depending on configuration, the VSim DICOM application will keep the association open for responses with a configurable time-out, or closes the association and expects responses on a different association that has to be established by the remote Storage Commitment SCP.

The commitment status derived from the related trigger response will be indicated in the related Status Flags of the related entity. It is possible to create triggers (“auto rules”) from this event.

The Transaction UUIDs of the pending commitment request are kept “open” for a configurable amount in time (default: 1h). If the “open time” for a pending commitment request has elapsed w/o a related response from the provider, the Transaction UUID is removed and the related entities are indicated as “commit failed”.

In any case, commitment will only be requested for previously and successfully sent images.

The Storage Commitment SCP is running in background and is ready to receive request when the system is started.

Sequencing of real World Activities

The Storage Commitment trigger is automatically derived from the successful completion of a Send Job.

AE Specification Storage Commitment

SIEMENS VSim DICOM application provides Standard Conformance to the following DICOM V3.0 SOP Class as an SCU and SCP:

Table 8. Storage Commitment SOP Class as an SCU and SCP

SOP Class Name	SOP Class UID
Storage Commitment Push Model	1.2.840.10008.1.20.1

Association Establishment Policies

General

With a Send Job successfully completed, the DICOM application will generate an Storage Commitment Identifier which references to all Instances of the processed job. The Commit Request is then sent over a single opened association. The VSim will wait for Status responses of the Storage Commitment Request. If the Provider accepts the Storage Commitment with Success Status, the generated Transaction UID, together with study identification data and a time-stamp, is kept. Depending on configuration, the association is closed when the configured time-out has elapsed or a response was received before. If the association is closed before a response was received, the response is then expected on a different association. Multiple Storage Commitment Requests can be pending.

The default PDU size used will be 28 KB.

Number of Associations

The Siemens VSim DICOM application initiates several associations at a time, one for each storage commitment request being processed. There is no limit for the number of initiated associations.

The Siemens VSim DICOM application is able to accept more associations at a time. It can handle in parallel 10 associations, and a configurable number of association requests can be queued until one of the 10 associations is released.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens VSim DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- “SIEMENS_SWFVB10A”

Association Initiation Policy

The VSim DICOM Application Entity acts as a Service Class User (SCU) for the

- Storage Commitment Push Model Service Class (to request commitment for storage of instances previously sent).

To do so, the VSim will issue a

- N-ACTION DIMSE to request commitment or a
- N-EVENT-REPORT DIMSE to respond to a received storage commitment request and the association was closed by the remote system prior to response.

Real World Activity

- Associated Real-World Activity - Job Completed

The VSim Storage Commitment application sends the commit request (N-ACTION-RQ) message and waits for acceptance of this request (N-ACTION-RSP). After receiving this, the transaction is marked as “waiting”.

Depending on a configuration value, the association will then be closed or kept open. In the first case, there is another configurable value giving the number of minutes (by default 60) to wait for the corresponding commit response (N_EVENT_REPORT). In the second case, this time is the (also configurable) time-out for the association. For both cases, if the commit response (N_EVENT_REPORT) does not arrive during the configured time, the transaction will be marked as failed.

If the commit response (N_EVENT_REPORT) received has the status of “complete - failure exists”, the transaction is marked as failed, else the transaction is marked as “completed”; In both cases, a message is shown to the user.

- Proposed Presentation Contexts - Job Completed

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Table 9. Proposed presentation contexts - request Storage commitment

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP Specific Conformance Statement- Job Completed

Storage Commitment is supported for all the SOP class UIDs as mentioned in 'Acceptable presentation contexts - Storage' in the Storage SCP section of this document.

The Referenced Study Component Sequence is not supported.

Storage Media File-Set ID and UID Attributes will not be supported in the commitment request (N-ACTION primitive) invoked by the Storage Commitment SCU.

- Associated Real-World Activity - Send Commit Response

Acting as an Storage Commitment Provider, the VSim Storage Commitment AE received an Storage Commitment request, carried out the request, and is ready to send back the response, but the association is not open anymore. In this case it will by itself initiate an association to send the storage commitment response (N_EVENT_REPORT) to the SCU.

- Proposed Presentation Contexts - Send Commitment Response

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Table 10. Proposed presentation contexts - send Storage commitment response

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP Specific Conformance Statement - Send Commitment Response

Storage Media File-Set ID and UID Attributes will not be supported in the N-EVENT-REPORT primitive invoked by the Storage Commitment SCP.

Association Acceptance Policy

The VSim DICOM Application Entity acts as a Service Class Provider (SCP) for the

- Storage Commitment Push Model Service Class (to proof commitment for storage of instances previously received).

To do so, the VSim attempts to accept a

- N-ACTION DIMSE to receive a commitment request for the instance included or a
- N-EVENT-REPORT DIMSE to receive a storage commitment response from a previous request and the SCP behavior requires a different association than the commit request.

Associated Real-World Activity - Commit SCP

- Associated Real-World Activity - Receive Commit Request

When receiving an Storage Commitment request the VSim DICOM application will perform the necessary steps to check the received list Instances against the local database or, if configured, check the Instances with the attached archive system.

- Accepted Presentation Contexts - Receive Commit Request

The Siemens VSim DICOM application will accept Presentation Contexts as shown in the following table:

Table 11. Accepted presentation contexts - Receive Storage commitment request

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotia- tion
Name	UID	Name List	UID List		
Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP-specific Conformance Statement - Receive Commit Request

The VSim Storage Commitment DICOM Application can be configured to run on an archive system.

If the Storage Commitment Application is running on an archive system, it will interact with this archive system in order to commit the storage of images and will send back to the SCU the result of the operation.

If not running on an archive node, the VSim Storage Commitment AE will return success for images that are stored in the local database and failure for images that are not. However, the committed images can later be deleted by the user on the SCP side!

Associated Real-World Activity - Commit SCU

- Associated Real-World Activity - Update Flags
The VSim Storage Commitment DICOM Application has sent a Storage Commitment Request and, being configured to receive response on a separate association, has closed the association, and now it gets an association request from the Storage Commitment SCP that want to send the results. Consecutive to start-up, the VSim DICOM application will await Storage commitment Notification triggers. Any incoming Notification will be checked for validity, that is, if the related Transaction UID is still part of the Pending Request Queue.

If the Notification is valid, the Notification Identifier is evaluated and the related Instances marked with the related status. The over-all Commit Status of the higher Information Entities is derived from propagation of the States of all Image entities included in a study.

The Status Flags directly affected by Storage Commitment results and indicated in the different entities of the Patient Browser list can be one of

- “AC” or “SC” - Successful Commitment, A means archived to configured Archive destination, whereas S means sent to any other destination
- “A” or “Sf” - Commitment failed.
- “A?” or “S?” - Commitment request is sent, response is pending.

In case of failure the user has to repeat the transfer of images to the Archive destination. Another Storage Commitment will be performed after sending is completed successfully.

- Accepted Presentation Contexts - Update Flags

The Siemens VSim DICOM application accepts Presentation Contexts shown in the table:

Table 12. Accepted presentation contexts - Update Flags (receive commit response)

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotia- tion
Name	UID	Name List	UID List		
Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP-specific Conformance Statement - Update Flags
If the Commitment response (N_EVENT_REPORT) received has the status of “complete - failure exists”, the transaction is marked as failed, else the transaction is marked as “completed”; In both cases, a message is shown to the user.

The related status flags are set for the committed images in the local database.

The VSim DICOM application will NOT support the Storage Media File Set ID attributes.

Implementation Model Query/Retrieve

The query/retrieve service class defines an application-level class of services which facilitates the management of images and patient data against the well defined information model of DICOM and allows a DICOM AE to retrieve images from a remote DICOM node or to request a remote DICOM AE to initiate a transfer of images to another DICOM AE. The VSim DICOM query/retrieve application supports the query/retrieve services to act as SCU and SCP.

Application Data Flow Diagram

The VSim DICOM network implementation acts as SCU and SCP for the query/retrieve network service. The VSim target Operating System is Microsoft Windows 2000.




Figure 6. Application Data Flow Diagram - Query/Retrieve SCU




Figure 7. Application Data Flow Diagram - Query/Retrieve SCP

Functional Definitions of Application Entities

The VSim DICOM query/retrieve SCU requests the remote query/retrieve SCP to perform a search and match to the keys specified in the request in order to display the results in the VSim user interface. Depending on user action (Import) the VSim DICOM SCU sends a C-MOVE DIMSE service to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query/Retrieve SCP) to the VSim Storage SCP.

The VSim DICOM query/retrieve SCP responds to C-FIND DIMSE services from remote SCU and depending on further remote request a C-GET or a C-MOVE involves the Siemens VSim DICOM query/retrieve SCP application to initiate a C-STORE sub-operation (by triggering and parametrizing the own Storage SCU) to send image objects to a remote Storage SCP.

All components of the Siemens DICOM query/retrieve SCP application are operating as background server processes. They are existing when the machine is powered on and then respond to queries based on the records stored in its database.

Sequencing of real World Activities

Retrieve of images is only possible if results from a previous “Search...” operation exist and those entities can be selected for “Import”.

AE Specification Query/Retrieve

Query/Retrieve Service AEs Specification

The Query/Retrieve SCU requests that the remote SCP performs a match of all keys specified in the request, against the information in its database and the identified images will be moved over a different (C-MOVE) storage association.

The Query/Retrieve SCP responds to queries based on the records in its database and images will be sent to the requesting SCU or to a different storage destination.

SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as SCU:

Table 13. Query/Retrieve SOP Classes as an SCU

SOP Class Name	SOP Class UID
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2

SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as an SCP:

Table 14. Query/Retrieve SOP Classes as an SCP

SOP Class Name	SOP Class UID
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2
Patient Root Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.1.3
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2
Study Root Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.2.3

Table 14. Query/Retrieve SOP Classes as an SCP

SOP Class Name	SOP Class UID
Patient/Study Only Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.3.1
Patient/Study Only Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.3.2
Patient/Study Only Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.3.3

Note See also the Storage DICOM Conformance Statement section of the Siemens VSim DICOM application to compare for conformance of the C-STORE sub-operation generated by the C-GET or C-MOVE DIMSE services. Furthermore compare the supported Storage Service SOP classes described in the Storage DICOM Conformance Statement of the Modality to which the images shall be transferred to.

Association Establishment Policies

General

With the “Search...” function the query data are input and the DICOM query/retrieve application is started. An query request will be sent out to one remote node that can be selected from a list of configured Query Providers and the response data will be displayed for the user. Upon request (Import), the retrieval of selected items is initiated.

The default PDU size used will be 28 KB.

Number of Associations

The Siemens VSim DICOM application initiates one association for each query request being processed to a remote node. There is no limit for the number of initiated associations.

The Siemens VSim DICOM application is able to accept more associations at a time. It can handle in parallel 10 associations, and a configurable number of association requests can be queued until one of the 10 associations is released.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens VSim DICOM software provides a single Implementation Class UID of

- “1.3.12.2.1107.5.9.20000101”

and an Implementation Version Name of

- modify “SIEMENS_SWFVB10A”

Association Initiation Policy

The query user interface will request the query-data from user and triggers one C-FIND request to the selected remote node. The response data will be displayed in the query UI for further data navigation.

When requesting Import of related items the browser requests the retrieve application to send a C-MOVE request to the related remote node. Images will then be received by the Storage SCP as described in the related section.

Real World Activity - Find SCU

- Associated Real-World Activity - “Search”

The associated Real-World activity is to fill out a query form with search data and pass it as query to the network application which issues a C-FIND over a previously build association. The remote SCP will respond with related data-entries that will be passed to a browser application. When data transfer is finished the association is closed.

- Proposed Presentation Contexts - “Search”

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotia- tion
Name	UID	Name List	UID List		
Patient Root Query/ Retrieve Find	1.2.840.10008.5.1.4 .1.2.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/ Retrieve Find	1.2.840.10008.5.1.4 .1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

It is configurable which of the two query models (or both) are to be used by the VSim DICOM Query SCU application.

- SOP Specific Conformance Statement - Find SCU
The Siemens VSim DICOM Query/Retrieve SCU supports hierarchical queries with all mandatory search keys. The interactive querying of attributes on IMAGE level is not supported by the Query SCU, hence retrieval of individual Objects is possible. The following table describes the search keys for the different query models that the SCU supports. Matching is either wildcard, which means that the user can supply a string containing wildcards, or universal, which means that the attribute is requested as return value.

Table 15. Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Patient Level^a					
Patient Name	(0010,0010)	R	Wildcard ^b	enter value	yes
Patient ID	(0010,0020)	U	Wildcard ^b	enter value	yes
Patient's Birth date	(0010,0030)	O	universal (Null)	-	yes
Patient's Sex	(0010,0040)	O	universal (Null)	-	yes
Number of Patient related Studies	(0020,1200)	O	universal (Null)	-	yes
Number of Patient related Series	(0020,1202)	O	universal (Null)	-	yes
Number of Patient related Instances	(0020,1204))	universal (Null)	-	yes
Study Level					
Patient Name ^c	(0010,0010)	R	Wildcard ^b	enter value	yes
Patient ID	(0010,0020)	R	Wildcard ^b	enter value	yes
Study Instance UID	(0020,000D)	U	single value	select from list	yes
Study ID	(0020,0010)	R	universal (Null)	-	yes
Study Date	(0008,0020)	R	universal (Null)	-	yes

Table 15. Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Study Time	(0008,0030)	R	universal (Null)	-	yes
Accession Number	(0008,0050)	R	universal (Null)	-	yes
Study Description	(0008,1030)	O	universal (Null)	-	yes
Referring Physician's Name	(0008,0090)	O	universal (Null)	-	yes
Name of Physician Reading Study	(0008,1060)	O	universal (Null)	-	yes
Modalities in Study	(0008,0061)	O	universal (Null)	-	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	yes
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	yes
Number of Study related Series	(0020,1206)	O	universal (Null)	-	yes
Number of Study related Instances	(0020,1208)	O	universal (Null)	-	yes
Series Level					
Series Instance UID	(0020,000E)	U	single value	select from list	yes
Series Number	(0020,0011)	R	universal (Null)	-	yes
Modality	(0008,0060)	R	universal (Null)	-	yes
Series Date	(0008,0021)	O	universal (Null)	-	yes
Series Time	(0008,0031)	O	universal (Null)	-	yes
Series Description	(0008,103E)	O	universal (Null)	-	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	yes
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	yes
Protocol Name	(0018,1030)	O	universal (Null)	-	yes
Perf. Procedure Step Start Date	(0040,0244)	O	universal (Null)	-	yes
Perf. Procedure Step Start Time	(0040,0245)	O	universal (Null)	-	yes
Number of Series related Instances	(0020,1209)	O	universal (Null)	-	yes

Table 15. Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Image Level					
SOP Instance UID	(0008,0018)	U	single value	-	yes
Image Number	(0020,0013)	R	universal (Null)	-	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	yes
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	yes
Instance Date	(0008,0023)	O	universal (Null)	-	yes
Instance Time	(0008,0033)	O	universal (Null)	-	yes
Number of Frames	(0028,0008)	O	universal (Null)	-	yes
Image Comments	(0020,4000)	O	universal (Null)	-	yes

- a.Patient Root Information Model only
- b.Always a '*' is appended to the user-supplied string
- c.Study Root Information Model only

The syngo Query / Retrieve application is using only a subset of the search attributes supported by study transfer. StudyTransfer API supports all DICOM query models: Patient_Root, Study_Root, Patient/Study_Only.

The StudyTransfer Find SCU is supporting the same attributes as the StudyTransfer Find SCP-

The Find SCU interprets following status codes:

Table 16. C-FIND response status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources	A700	(0000,0902)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	Cxxx	(0000,0901) (0000,0902)

Table 16. C-FIND response status

Service Status	Meaning	Protocol Codes	Related Fields
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

Real World Activity - Move SCU

- Associated Real-World Activity - Move SCU “Import”

When selecting a data entry in the Query UI and activating the “Import” function, a retrieval request is passed to the archival application which issues a C-MOVE service according to the Patient Root or Study Root query model. (The Storage Service Class Conformance Statement of the SCP will describe the C-STORE service which is generated by the C-MOVE service.)

The transferred image data are processed as described in the storage class SCP descriptions.

The possibility to request the remote C-MOVE provider (remote application that responded to the C-FIND) to move data to an application entity other than the C-MOVE SCU (the VSim DICOM application) is NOT USED.

C-MOVE operation on Patient Level is not supported by the Query UI.

- Proposed Presentation Contexts - Move SCU
The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.2.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Note C-MOVE extended negotiation will not be supported by the SCU.

- SOP Specific Conformance Statement - Move SCU

At association establishment time the C-MOVE presentation context shall be negotiated. The C-STORE sub-operations must be done on a different association to transfer images to the own Storage Service Class SCP.

The Move SCU interprets following status codes:

Table 17. C-MOVE response status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
Failed	Unable to process	Cxxx	(0000,0901) (0000,0902)

Table 17. C-MOVE response status

Service Status	Meaning	Protocol Codes	Related Fields
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures or Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

Association Acceptance Policy

The Siemens VSim DICOM application will accept associations for the following DIMSE-C operations as SCP:

- C-FIND
- C-GET
- C-MOVE
- C-FIND-CANCEL
- C-GET-CANCEL
- C-MOVE-CANCEL

Extended negotiation - which is relational retrieve - is NOT supported for the above listed services. The VSim DICOM application does support multiple C-FIND requests over the same association, while multiple C-MOVE or C-GET operations are NOT supported over the same association.

Real World Activity - Find SCP

- Associated Real-World Activity - Find request to SCP

The associated Real-World activity is to respond query requests to an SCU with the query model Patient Root, Study Root and Patient/Study Only. Relational retrieve operation is

NOT supported. With a C-FIND-CANCEL request the running query can be canceled at any time.

Multiple C-FIND requests over the same association are supported.

- Accepted Presentation Contexts - Find SCP

The Siemens VSim DICOM application will accept Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Patient/Study Only Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.3.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Note C-FIND extended negotiation will not be supported by the SCP.
The order of preference for accepting Transfer Syntaxes is: 1. Explicit VR Little Endian, 2. Explicit VR Big Endian, 3. Implicit VR Little Endian.

- SOP Specific Conformance Statement - Find SCP

The Siemens VSim DICOM Query/Retrieve SCP supports hierarchical queries with all mandatory and optional search keys.

The query attribute contents will be treated case-sensitive.

With wildcard queries the symbol “?” is treated as “*” by the C-FIND SCP application. As a consequence the query string of “?abc*” will be processed as “*abc*”.

If the value for the patient-level unique key “Patient ID” is not known, it may be returned with zero length. The attribute “Image Comments” will not be include in the C-FIND-RSP, if it is not set in the DB, even if it was requested as return key in the related C-FIND-RQ.

Usage of Storage Media File-Set ID, Retrieve AE Title with C-FIND-RSP message:

- The Storage Media File-Set ID - if existent - can be returned at Study/Series/Image Level. Only on Image Level, the values of ONLINE, NEARLINE or OFFLINE are returned to indicate the Storage Location of the related Instance.
- The Retrieve AE Title - if existent - can only be returned at Image Level (for Patient Root and Study Root models) or Study Level (for Patient/Study Only model).

Relational Queries are not supported.

A remote DICOM AE can cancel the running query by sending a C-FIND-CANCEL. Matches are possibly continuing (more C-FIND responses with status PENDING) until the cancel operation has completed.

The supported attributes on the various query levels of the three supported information models are list in the tables of the following sections.

Patient Root Information Model

Table 18. Patient level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal
Patient ID	(0010,0020)	U	single value, wildcard, universal
Patient's Birth Date	(0010,0030)	O	single value, range, universal
Patient's Birth Time	(0010,0032)	O	single value, range, universal
Patient's Sex	(0010,0040)	O	single value, wildcard, universal
Ethnic Group	(0010,2160)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal
Number of Patient related Instances	(0020,1204)	O	universal

Table 19. Study level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Study Instance UID	(0020,000D)	U	single value, list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,1060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal
Number of Study related Instances	(0020,1208)	O	universal

Table 20. Series level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Series Instance UID	(0020,000E)	U	single value, list of UID
Series Number	(0020,0011)	R	single value, universal
Modality	(0008,0060)	R	single value, wildcard, universal
Laterality	(0020,0060)	O	single value, wildcard, universal
Body Part Examined	(0018,0015)	O	single value, wildcard, universal
Patient Position	(0018,5100)	O	single value, wildcard, universal
Smallest Pixel Value in Series	(0028,0108)	O	single value, universal
Largest Pixel Value in Series	(0028,0109)	O	single value, universal
Protocol Name	(0018,1030)	O	single value, wildcard, universal
Series Date	(0008,0021)	O	single value, range, universal
Series Time	(0008,0031)	O	single value, range, universal
Series Description	(0008,103E)	O	single value, wildcard, universal
Operators Name	(0008,1070)	O	single value, wildcard, universal
Performing Physician's Name	(0008,1050)	O	single value, wildcard, universal
Perf. Procedure Step Start Date	(0040,0244)	O	universal
Perf. Procedure Step Start Time	(0040,0245)	O	universal
Number of Series related Instances	(0020,1209)	O	universal

Table 21. Image level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
SOP Instance UID	(0008,0018)	U	single value, list of UID
Image Number	(0020,0013)	R	single value, universal
Image Date	(0008,0023)	O	single value, range, universal
Image Time	(0008,0033)	O	single value, range, universal
Modality	(0008,0060)	O	single value, wildcard, universal
Image Comments	(0020,4000)	O	universal

Study Root Information Model

Table 22. Study level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal
Patient ID	(0010,0020)	R	single value, wildcard, universal
Patient’s Birth Date	(0010,0030)	O	single value, range, universal
Patient’s Birth Time	(0010,0032)	O	single value, range, universal
Patient’s Sex	(0010,0040)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal

Table 22. Study level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Number of Patient related Instances	(0020,1204)	O	universal
Study Instance UID	(0020,000D)	U	single value, list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,8060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal

Table 22. Study level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Number of Study related Instances	(0020,1208)	O	universal

Table 23. Series level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Series Instance UID	(0020,000E)	U	single value, list of UID
Series Number	(0020,0011)	R	single value, universal
Modality	(0008,0060)	R	single value, wildcard, universal
Laterality	(0020,0060)	O	single value, wildcard, universal
Body Part Examined	(0018,0015)	O	single value, wildcard, universal
Patient Position	(0018,5100)	O	single value, wildcard, universal
Smallest Pixel Value in Series	(0028,0108)	O	single value, universal
Largest Pixel Value in Series	(0028,0109)	O	single value, universal
Protocol Name	(0018,1030)	O	single value, wildcard, universal
Series Date	(0008,0021)	O	single value, range, universal
Series Time	(0008,0031)	O	single value, range, universal
Series Description	(0008,103E)	O	single value, wildcard, universal
Operators Name	(0008,1070)	O	single value, wildcard, universal

Table 23. Series level attributes, Study Root Information Model

Performing Physician's Name	(0008,1050)	O	single value, wildcard, universal
Performed Procedure Step Start Date	(0040,0244)	O	universal
Performed Procedure Step Start Time	(0040,0245)	O	universal
Number of Series related Instances	(0020,1209)	O	universal

Table 24. Image level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
SOP Instance UID	(0008,0018)	U	single value, list of UID
Image Number	(0020,0013)	R	single value, universal
Image Date	(0008,0023)	O	single value, range, universal
Image Time	(0008,0033)	O	single value, range, universal
Modality	(0008,0060)	O	single value, wildcard, universal
Image Comments	(0020,4000)	O	universal

Patient Study Only Information models

Table 25. Patient instance level, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal

Table 25. Patient instance level, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Patient ID	(0010,0020)	U	single value, wildcard, universal
Patient’s Birth Date	(0010,0030)	O	single value, range, universal
Patient’s Birth Time	(0010,0032)	O	single value, range, universal
Patient’s Sex	(0010,0040)	O	single value, wildcard, universal
Ethnic Group	(0010,2160)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal
Number of Patient related Instances	(0020,1204)	O	universal

Table 26. Study level attributes, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Study Instance UID	(0020,000D)	U	single value, list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal

Table 26. Study level attributes, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,8060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal
Number of Study related Instances	(0020,1208)	O	universal

The Find SCP returns following status codes:

Table 27. C-FIND SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources	A700	(0000,0902)

Table 27. C-FIND SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	C001	(0000,0901) (0000,0902)
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

Real World Activity - Get SCP

- Associated Real-World Activity - Get request to SCP
The associated Real-World activity is to respond to retrieve requests initiated from a foreign SCU. The SCP supports the query model Patient Root, Study Root and Patient/Study Only. The Storage Service Class Conformance Statement describes the C-STORE service which is generated by the C-GET service. Relational retrieve operation is NOT supported.

Multiple C-GET requests over the same association are NOT supported.

- Accepted Presentation Contexts - Get SCP
The Siemens VSim DICOM application will accept Presentation Contexts as shown in the following table:

Table 28. Presentation Contexts - Get SCP

Presentation Context Table				
Abstract Syntax		Transfer Syntax		Role
Name	UID	Name List	UID List	
				Extended Negotiation

Table 28. Presentation Contexts - Get SCP

Patient Root Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.1.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.2.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Patient/Study Only Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.3.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Note C-GET extended negotiation will not be supported by the SCP.

The order of preference for accepting Transfer Syntaxes is: 1. Explicit VR Little Endian, 2. Explicit VR Big Endian, 3. Implicit VR Little Endian.

- SOP Specific Conformance Statement - Get SCP

At association establishment time the C-GET presentation context must be negotiated along with the C-STORE sub-operations which must be accomplished on the same association as the C-GET operation. Relational retrieve operation is NOT supported.

All unique keys have to be supplied according to the selected Query/Retrieve Level. The related tables in the C-FIND SCP section will give information about “U” marked key attributes. The Get SCP returns following status codes:

Table 29. C-GET SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	C001	(0000,0901) (0000,0902)

Table 29. C-GET SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures or Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

- Real World Activity - Move SCP

Associated Real-World Activity - Move request to SCP

The associated Real-World activity is to respond to retrieve requests to an SCU. The SCP supports the query model Patient Root, Study Root and Patient/Study Only. The Storage Service Class Conformance Statement describes the C-STORE service which is generated by the C-MOVE service. Relational retrieve operation is NOT supported.

Multiple C-MOVE requests over the same association are NOT supported.

- Accepted Presentation Contexts - Move SCP

The Siemens VSim DICOM application will accept Presentation Contexts as shown in the following table:

Table 30. Presentation Contexts - Move SCP

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		

Table 30. Presentation Contexts - Move SCP

Patient Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.2.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Patient/Study Only Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.3.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Note C-MOVE extended negotiation will not be supported by the SCP.

The order of preference for accepting Transfer Syntaxes is: 1. Explicit VR Little Endian, 2. Explicit VR Big Endian, 3. Implicit VR Little Endian.

- SOP Specific Conformance Statement - Move SCP

At association establishment time the C-MOVE presentation context shall be negotiated. The C-STORE sub-operations is done on a different association, specified in the C-MOVE request, to transfer images to a remote SCP of the Storage Service Class. Relational retrieve operation is NOT supported.

All unique keys have to be supplied according to the selected Query/Retrieve Level. The related tables in the C-FIND SCP section will give information about “U” marked key attributes. Multiple C-MOVE requests over the same association are NOT supported.

The Move SCP returns following status codes:

Table 31. C-MOVE return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

Table 31. C-MOVE return status

Service Status	Meaning	Protocol Codes	Related Fields
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	C001	(0000,0901) (0000,0902)
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures or Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

Implementation Model Print

The Print Management Service Classes define an application-level class of services which facilitate the printing of images on a hardcopy medium. The print management SCU and print management SCP are peer DICOM print management application entities. The *syngo* DICOM print application supports the print management DIMSE services to act as SCU.

Application Data Flow Diagram

The *syngo* DICOM network implementation acts as SCU for the print management network service. The VSim target Operating System is Microsoft Windows 2000.




Figure 8. Application Data Flow Diagram - Print SCU

Functional Definitions of Application Entities

The Print SCU is invoked by the user interface to setup film-sheet layout and whenever an image is ready to be printed on film. The Print SCU will hold and maintain all data needed to compile a complete film-sheet from the data (images, layout, configuration) received. Whenever a film-sheet is ready to print the related data is used to supply the Information to the SOP Classes of the Print Management Service Class. A queue is maintained in order to intermediately store several film-sheets in case of resource problems on printer. The SCU will only supply and require the mandatory SOP Classes of the Print Management Service Class.

Sequencing of real World Activities

Not applicable.

AE Specification Print

Print Management AE Specification

The *syngo* print management SCU (HCS) invokes print management DIMSE services to transfer images from the local AE to the remote SCP AE to print images with defined layout on a selected network-based DICOM hardcopy printer. This is done in a “full-page” print mode.

SIEMENS *syngo* DICOM products provide Standard Conformance to the following DICOM V3.0 Print Management Meta SOP Classes as an SCU:

Table 32. Basic Grayscale Print Management Meta SOP Classes

SOP Class Name	SOP Class UID
Basic Grayscale Print Management Meta SOP Class	1.2.840.10008.5.1.1.9
- Basic Film Session SOP Class	1.2.840.10008.5.1.1.1
- Basic Film Box SOP Class	1.2.840.10008.5.1.1.2
- Basic Grayscale Image Box SOP Class	1.2.840.10008.5.1.1.4
- Printer SOP Class	1.2.840.10008.5.1.1.16
Print Job SOP Class	1.2.840.10008.5.1.1.14
Presentation LUT SOP Class	1.2.840.10008.5.1.1.23

Table 33. Basic Color Print Management Meta SOP Classes

SOP Class Name	SOP Class UID
Basic Color Print Management Meta SOP Class	1.2.840.10008.5.1.1.18
- Basic Film Session SOP Class	1.2.840.10008.5.1.1.1
- Basic Film Box SOP Class	1.2.840.10008.5.1.1.2
- Basic Color Image Box SOP Class	1.2.840.10008.5.1.1.4.1

Table 33. Basic Color Print Management Meta SOP Classes

SOP Class Name	SOP Class UID
- Printer SOP Class	1.2.840.10008.5.1.1.16
Print Job SOP Class	1.2.840.10008.5.1.1.14

Association Establishment Policies

General

Whenever a film is completely set up and printed by command or automatism, the job is prepared for processing. As soon as the queue is ready to process, the job is activated and worked according to the processing data. The related Print application will initiate an association to the print destination and process the printing of the related information.

The default PDU size used will be 28 KB.

Number of Associations

The Siemens *syngo* DICOM application initiates one association at a time for each different print device configured.

Asynchronous Nature

The Siemens *syngo* DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens *syngo* DICOM software provides a single Implementation Class UID of

- 1.2.12.2.1107.5.9.20010101
- and an Implementation Version Name of
- SIEMENS_SWFVC20A

Association Initiation Policy

Triggered by the Print job queue the Print Management SCU establishes an association by using the DICOM association services. With the help of the N-GET request for the Printer SOP Class the Status is determined before printing.

With no problem encountered with the N-CREATE/N-SET Services for the related Basic Print SOP Classes the film sheet is set up for printing and the image(s) is (are) transferred to the printer device.

After the last film is printed from queue, the Print application will leave open the association for another 60 seconds. If a new film job is ready for printing within this time-limit, the job will be immediately processed over the still open association. If there is no new job, the association is closed if the time-out elapsed. This is done to optimize automated printing.

During the “idle-time” (no open association to printer) the Print application will issue a cyclic camera status request (using N-GET of Printer SOP Class) every 5 minutes.

Real-World Activity - Print

- Associated Real-World Activity - Printing a Printer Job Queue Entry

Whenever a film sheet is prepared by the user it is forwarded to the Printer Job queue. As soon as the associated Printer device is available the job is activated and association is set up.

The film sheet is internally processed, converted to a Standard/1-1 page and then the page image is sent. Status is controlled by awaiting any N-EVENT message all through the transfer until the last image or film-sheet is sent.

If the response from the remote application contains a status other than Success or Warning the association is aborted.

- Proposed Presentation Contexts

The Siemens *syngo* DICOM application will propose Presentation Contexts as shown in the following table:

Table 34. Presentation Contexts

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Basic Grayscale Print Management Meta SOP class	1.2.840.10008.5.1.1.9	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic Color Print Management Meta SOP class	1.2.840.10008.5.1.1.18	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Table 34. Presentation Contexts

Basic film session SOP class	1.2.840.10008.5.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic film box SOP class	1.2.840.10008.5.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic grayscale image box SOP class	1.2.840.10008.5.1.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic color image box SOP class	1.2.840.10008.5.1.1.4.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Printer SOP class	1.2.840.10008.5.1.1.16	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Print Job SOP class	1.2.840.10008.5.1.1.14	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Presentation LUT SOP class	1.2.840.10008.5.1.1.23	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP Specific Conformance - Meta SOP Classes

The Siemens *syngo* DICOM print management SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class and the Basic Color Print Management Meta SOP Class.

The application uses a setting platform to define the properties of the connected DICOM SCP, e.g.:

- maximum number of print jobs in the queue

- maximum number of print copies
- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP
- lookup table definition.

The printing is only suspended in the case of a failure return status of the SCP.

Basic film session SOP class

The Basic Film Session information object definition describes all the user defined parameters which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes and that are printed on one hardcopy printer.

The *syngo* DICOM print management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the following attributes:

Table 35. Used Basic Film Session N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Number of Copies	(2000,0010)	U	1
Medium Type	(2000,0030)	U	BLUE FILM CLEAR FILM PAPER
Film Destination	(2000,0040)	U	MAGAZINE PROCESSOR

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session - see below:

Table 36. Attributes of N-DELETE-RQ on Basic Film Session SOP Class

Attribute Name	Tag	Source of Information
Requested SOP Instance UID	(0000,1000) -> (0000,1001)	Affected SOP Instance UID of N-CREATE-RSP on Basic Film Session

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

The Basic Film Session SOP class interprets the following status codes (from N_CREATE_RSP, N_DELETE_RSP messages):

Table 37. Basic Film Session SOP status

Service Status	Meaning	Protocol Codes
Failure	Film session SOP instances hierarchy does not contain film box SOP instances	C600
	Unable to create print job, print queue is full	C601
	Image size is larger than images box size	C603
Warning	Memory allocation not supported	B600
	Film session printing is not supported	B601
	Film box does not contain image box (empty page)	B602
Success	Film belonging to the film session are accepted for printing	0000

Basic Film Box SOP class

The Basic Film Box information object definition describes all the user defined parameter of one film of the film session. The Basic Film Box information description defines the presentation parameters which are common for all images on a given sheet of film.

The Basic Film Box refers to one or more Image Boxes.

Supported Service Elements as SCU are:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP class N-CREATE-RQ message uses the following attributes (the actual values for each attribute depend on DICOM printer configuration within the VSim DICOM print management SCU)

Table 38. Used Film Box N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Image Display Format	(2010,0010)	M	STANDARD\1,1
Referenced Film Session Sequence	(2010,0500)	M	
>Referenced SOP Class UID	(0008,1150)	M	1.2.840.10008.5.1.1.1
>Referenced SOP Instance UID	(0008,1155)	M	
Film Orientation	(2010,0040)	M	PORTRAIT
Film Size ID	(2010,0050)	M	8INx10IN, 10INx12IN, 10INx14IN, 11INx14IN, 14INx14IN, 14INx17IN, 24CMx24CM, 24CMx30CM
Magnification Type	(2010,0060)	M	BILINEAR, CUBIC, NONE, REPLICATE
Border Density	(2010, 0100)	U	BLACK, WHITE
MaxDensity	(2010,0130)	U	>0
MinDensity	(2010,0120)	U	>0 but <50
Illumination	(2010,015E)	U	>0 Required if Presentation LUT is present
Reflective Ambient Light	(2010,0160)	U	>0 Required if Presentation LUT is present
Referenced Presentation LUT Sequence	(2050,0500)	U	

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally to be further used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the *syngo* DICOM print manager will issue a N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Box - see below:

Table 39. Attributes of N-DELETE-RQ on Basic Film Box SOP Class

Attribute Name	Tag	Source of Information
Requested SOP Instance UID	(0000,1000) -> (0000,1001)	Affected SOP Instance UID of N-CREATE-RSP on Basic Film Box

The Basic Film Box SOP class interprets the following status codes:

Table 40. Basic Film Box SOP status

Service Status	Meaning	Protocol Codes
Failure	Unable to create print job; print queue is full	C602
	Image size is larger than image box size	C603
Warning	Film box does not contain image box (empty page)	B603
	Requested MinDensity or MaxDensity outside of Printer’s operating range	B605
Success	Film accepted for printing	0000

Basic Grayscale Image Box SOP Class

The Basic Grayscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the following attributes:

Table 41. Mandatory Basic Grayscale Image Box N-SET attributes

Attribute name	Tag	Usage SCU	Supported Values
Image Position	(2020,0010)	M	1
Basic Grayscale Image Sequence	(2020,0110)	M	
>Samples Per Pixel	(0028,0002)	M	1
>Photometric Interpretation	(0028,0004)	M	MONOCHROM E2
>Rows	(0028,0010)	M	
>Columns	(0028,0011)	M	
>Pixel Aspect Ratio	(0028,0034)	M	
>Bits Allocated	(0028,0100)	M	8
>Bits Stored	(0028,0101)	M	8
>High Bit	(0028,0102)	M	7
>Pixel Representation	(0028,0103)	M	0
>Pixel Data	(7FE0,0010)	M	

The Grayscale Image Box SOP class interpret the following status codes:

Table 42. Basic Grayscale Image Box SOP status

Service Status	Meaning	Protocol Codes
Warning	Requested MinDensity or MaxDensity outside of Printer's operating range	B605
Failure	Image contains more pixel than printer can print in Image Box	C603
	Insufficient memory in printer to store the image	C605
Success		0000

Basic Color Image Box SOP Class

The Basic Color Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Color Image Box SOP Class uses only the N-SET-RQ with the following attributes:

Table 43. Mandatory Basic Color Image Box N-SET attributes

Attribute name	Tag	Usage SCU	Supported Values
Image Position	(2020,0010)	M	1
Basic Color Image Sequence	(2020,0111)	M	
>Samples Per Pixel	(0028,0002)	M	3
>Photometric Interpretation	(0028,0004)	M	RGB
>Planar Configuration	(0028,0006)	M	0
>Rows	(0028,0010)	M	
>Columns	(0028,0011)	M	
>Pixel Aspect Ratio	(0028,0034)	M	
>Bits Allocated	(0028,0100)	M	8
>Bits Stored	(0028,0101)	M	8
>High Bit	(0028,0102)	M	7
>Pixel Representation	(0028,0103)	M	0
>Pixel Data	(7FE0,0010)	M	

The Color Image Box SOP class interpret the following status codes:

Table 44. Basic Color Image Box SOP status

Service Status	Meaning	Protocol Codes
Warning	Image size larger than image box size	B604

Table 44. Basic Color Image Box SOP status

Service Status	Meaning	Protocol Codes
Failure	Image contains more pixel than printer can print in Image Box	C603
	Insufficient memory in printer to store the image	C605
Success		0000

Presentation LUT SOP Class

The objective of the Presentation LUT is to realize image hardcopy printing tailored for specific modalities, applications and user preferences.

The output of the Presentation LUT is Presentation Values (P-Values). P-Values are approximately related to human perceptual response. They are intended to facilitate common input for hardcopy. P-Values are intended to be independent of the specific class or characteristics of the hardcopy device.

The Presentation LUT SOP Class uses only the N-CREATE-RQ with the following attributes:

Table 45. Used Presentation LUT N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Presentation LUT Shape	(2050,0020)	U	IDENTITY

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests on the Basic Film Box (N-CREATE-RQ) and on the Presentation LUT (N-DELETE-RQ) - see below:

Table 46. Attributes of N-CREATE-RSP on Presentation LUT SOP Class

Attribute Name	Tag	Source of Information
Requested SOP Instance UID	(0000,1000) -> (0000,1001)	Affected SOP Instance UID of N-CREATE-RSP on Presentation LUT

The Presentation LUT SOP class interprets the following status codes:

Table 47. Presentation LUT SOP status

Service Status	Meaning	Protocol Codes
Warning	Requested MinDensity or MaxDensity outside of HCD’s operating range. HCD will use its respective minimum or maximum density value instead.	B605
Success	Presentation LUT successfully created	0000

Printer SOP Class

The Printer SOP Class is the possibility to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The SCU uses the mandatory N-EVENT Report DIMSE service to monitor the changes of the printer status in an asynchronous way.

It can directly ask the Printer (SCP) for its status or receive Events from the Printer asynchronously:

- N-GET as SCU
- N-EVENT-REPORT as SCU

In both cases the following information is supported:

Table 48. Used Printer N-EVENT report attributes

Event type name	Event	Attributes	Tag	Usage SCU
Normal	1			
Warning	2	Printer Status Info	(2110,0020)	U
Failure	3	Printer Status Info	(2110,0020)	U

Table 49. Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes

Attribute name	Tag	Usage SCP	supported values
Printer Status	(2110,0010)	M	NORMAL, FAILURE, WARNING
Printer Status Info	(2110,0020)	M	see tables in the Annex for details

Note For a detailed description on how *syngo* reacts on different printer status messages, please refer to the Annex section “DICOM Print SCU – detailed status displays”.

Print Job SOP Class

The Print Job SOP Class is the possibility to monitor the execution of the print process.

The *syngo* DICOM Print Management application supports the optional N-EVENT-REPORT DIMSE Service to receive the changes of the Print Job Status in an asynchronous way.

It can receive Events from the Print SCP asynchronously:

- N-EVENT-REPORT

The following information is supported:

Table 50. Used Print Job N-EVENT report attributes

Event type name	Event	Attributes	Tag	Usage SCU
Normal	1	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U
Printing	2	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U

Table 50. Used Print Job N-EVENT report attributes

Event type name	Event	Attributes	Tag	Usage SCU
Done	3	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U
Failure	4	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U

Note For a detailed description on how *syngo* reacts on different printer status messages, please refer to the Annex section “DICOM Print SCU – detailed status displays”.

Implementation Model Worklist

The Basic Worklist Management Service class defines an application-level class of service which facilitates the transfer of worklists from the information system to the imaging modality. The worklist is queried by the AE and supplies the SCU with the scheduled tasks which have to be performed on the modality. The VSim DICOM worklist application supports the worklist service as SCU.

Application Data Flow Diagram

The VSim DICOM network implementation acts as SCU for the Basic Worklist Service using the Modality Worklist SOP Class. The VSim target Operating System is Microsoft Windows 2000.

:




Figure 9. Application data flow diagram - Basic Worklist SCU

Functional Definitions of Application Entities

The worklist SCU (“broad query”) is invoked from the patient browser user interface or by timer to request the worklist from a remote Information System (Modality Worklist Class SCP). This is done to perform a match to the internal worklist query keys specified in the C-Find DIMSE service issued for the Modality Worklist Model.

The worklist SCP responses to the C-FIND query and scheduled imaging service requests (scheduled procedure steps) and patient demographic information will be downloaded from the information system to the VSim. All information retrieved will be hold in the scheduling database for usage during Patient registration procedure.

Furthermore, in Patient Registration dialog, it is possible to update/complete the entered data by a “Get Worklist” function. Some of the entered data will then be used as matching criteria (“narrow query”) for the issue worklist query. With the response data then the Patient Registration dialog fields are populated according availability within the worklist response identifier.

Sequencing of real World Activities

The “narrow” (interactive) Worklist Query requires that sufficient matching keys or a unique matching key are/is entered before the query is issued. Only then a single response can be expected to complete the registration dialog.

AE Specification Worklist

Modality Worklist Service AEs Specification

The Modality worklist SCU (patient registration in conjunction with the network application) requests that the remote SCP performs a match of all keys specified in the query against the information in its worklist database.

SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Class as an SCU:

Table 51. SOP Classes as an SCU

SOP Class Name	SOP Class UID
Modality Worklist Information Model - FIND	1.2.840.10008.5.1.4.31

Association Establishment Policies

General

It is possible to configure a cyclic update of the modality scheduler database through a background worklist request with date/time and modality information.

In addition the user can request worklist update with “Update Worklist”. No duplicate entries will be added in the Scheduler DB. Entries are uniquely identified by the Study Instance UID (0020,000D) for the Requested Procedure and the SPS ID (0040,009) in the SPS Sequence (0040,0100).

An interactive patient based worklist query can be issued with search criteria entered during patient registration.

The default PDU size used will be 28 KB.

Number of Associations

The Siemens VSim DICOM application initiates one association at a time to query worklist entry data.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens VSim DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- “SIEMENS_SWFVC20A”

Association Initiation Policy

The network application will cyclically query the worklist and by request of patient registration interface. Ever then it establishes an association by using the DICOM association services. During association establishment the negotiation of SOP classes to exchange the capabilities of the SCU and the SCP is not supported.

The following DIMSE-C operation is supported as SCU:

- C-FIND

Real World Activity

- Associated Real-World Activity - Update (Query) Worklist

A network application will perform worklist queries with the C-FIND request at regular intervals. In addition it can be triggered by immediate request. The received worklist items will be compared with the contents of the local scheduler database. New items will be inserted into scheduler database.

- Proposed Presentation Contexts - Update Worklist

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Table 52. Proposed presentation contexts

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Modality Worklist Information Model- FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP Specific Conformance Statement- Update Worklist

Search Key Attributes of the Worklist C-FIND

The Siemens VSim DICOM worklist SCU supports “broad worklist queries” with all required search keys. The following tables describe the “broad query” search keys that the SCU supports.

Table 53. Supported Broad Worklist Query Search Key Attributes

Attribute name	Tag	Matching Key Type	Query Value
Scheduled Procedure Step			
Scheduled Procedure Step Sequence	(0040,0100)	R	
>Scheduled Station AE Title	(0040,0001)	R	<own AET> or “*” ^a
>Scheduled Procedure Step Start Date	(0040,0002)	R	<act. Date>-<act. Date> or range from UI ^b
>Scheduled Procedure Step Start Time	(0040,0003)	R	<zero length> or range from UI ^b
>Modality	(0008,0060)	R	“*” or <own Modality> ^a

a.This depends on user configuration (Options->Configuration->Patient Registration) if the “own AET” or “Modality” is provided. Use the “HIS/RIS” tabcard for confirmation.
b.It depends on user configuration (Options->Configuration->Patient Registration) if the actual Date with a full time range or an interactive input dialog for date/time specification is used.

Return Key Attributes of the Worklist C-FIND-RSP

The Siemens VSim DICOM worklist SCU supports worklist queries with return key attributes of all types. The following tables describe the return keys that the SCU supports.

An “x” in the **UI** column will indicate that the attribute is visualized when browsing the Worklist results with Patient Browser and/or during Patient Registration. The Patient Browser display is additionally influenced by the related Browser configuration.

An “x” in the **IOD** column will indicate that the related attribute is included into the SOP Instances of the IOD’s created during processing of this worklist request.

Table 54. Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
SOP Common				
Specific Character Set	(0008,0005)	1C	-	x
Scheduled Procedure Step				
Scheduled Procedure Step Sequence	(0040,0100)	1		
>Scheduled Station AE Title	(0040,0001)	1	x	
>Scheduled Procedure Step Start Date	(0040,0002)	1	x	
>Scheduled Procedure Step Start Time	(0040,0003)	1	x	
>Scheduled Procedure Step End Date	(0040,0004)	3	-	
>Scheduled Procedure Step End Time	(0040,0005)	3	-	
>Modality	(0008,0060)	1	x	
>Scheduled Performing Physician's Name	(0040,0006)	1	x	x ^a
>Scheduled Procedure Step Description	(0040,0005)	1C	x	x
>Scheduled Station Name	(0040,0010)	2	x	
>Scheduled Procedure Step Location	(0040,0011)	2	x	
>Scheduled Action Item Code Sequence	(0040,0008)	1C	-	
>>Code Value	(0008,0100)	1C	x	
>>Coding Scheme Designator	(0008,0102)	1C	x	
>> Coding Scheme Version	(0008,0103)	3	x	
>>Code Meaning	(0008,0104)	3	x	
>Pre-Medication	(0040,0012)	2C	x	
>Scheduled Procedure Step ID	(0040,0009)	1	x	x
>Requested Contrast Agent	(0032,1070)	2C	x	x
> Scheduled Procedure Step Status	(0040,0020)	3	x	

Table 54. Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
>Comments on the Scheduled Procedure Step	(0040,0400)	3	-	
Requested Procedure				
Requested Procedure ID	(0040,1001)	1	x	x
Requested Procedure Description	(0032,1060)	1C	x	x
Requested Procedure Code Sequence	(0032,1064)	1C	-	
>Code Value	(0008,0100)	1C	x	
>Code Scheme Designator	(0008,0102)	1C	x	
> Coding Scheme Version	(0008,0103)	3	x	
>Code Meaning	(0008,0104)	3	x	
Study Instance UID	(0020,000D)	1	-	x
Referenced Study Sequence	(0008,1110)	2	-	
>Referenced SOP Class UID	(0008,1150)	1C	-	
>Referenced SOP Instance UID	(0008,1155)	1C	-	
Requested Procedure Priority	(0040,1003)	2	x	
Patient Transport Arrangements	(0040,1004)	2	-	
Reason for the Requested Procedure	(0040,1002)	3	-	
Confidentiality Code	(0040,1008)	3	-	
Reporting Priority	(0040,1009)	3	-	
Names of intended Recipients of Results	(0040,1010)	3	-	
Requested Procedure Comments	(0040,1400)	3	x	
Requested Procedure Location	(0040,1005)	3	-	
Imaging Service Request				
Accession Number	(0008,0050)	2	x	x
Requesting Physician	(0032,1032)	2	x	x

Table 54. Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
Referring Physician’s Name	(0008,0090)	2	x	x
Reason for Imaging Service Request	(0040,2001)	3	-	
Imaging Service Request Comments	(0040,2400)	3	x	
Requesting Service	(0032,1033)	3	x	
Issuing Date of Imaging Service Request	(0040,2004)	3	-	
Issuing Time of Imaging Service Request	(0040,2005)	3	-	
Placer Order Number / Imaging Service Request	(0040,2006)	3	-	
Filler Order Number / Imaging Service Request	(0040,2007)	3	-	
Order entered by ...	(0040,2008)	3	-	
Order Enterer’s location	(0040,2009)	3	-	
Order Callback Number	(0040,2010)	3	-	
Visit Identification				
Admission ID	(0038,0010)	2	x	
Issuer of Admission ID	(0038,0011)	3	-	
Visit Status				
Current Patient Location	(0038,0300)	2	x	
Visit Relationship				
Referenced Patient Sequence	(0008,1120)	2	-	
>Referenced SOP Class UID	(0008,1150)	2	-	
>Referenced SOP Instance UID	(0008,1155)	2	-	
Visit Admission				
Institution Name	(0008,0080)	3	x	
Admitting Diagnosis Description	(0008,1080)	3	x	
Patient Identification				

Table 54. Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
Patient's Name	(0010,0010)	1	x	x
Patient ID	(0010,0020)	1	x	x
Patient Demographic				
Patients Birth Date	(0010,0030)	2	x	x
Patient's Sex	(0010,0040)	2	x	x
Patient's Weight	(0010,1030)	2	x	x
Confidential constraint on patient data	(0040,3001)	2	x	x
Patient's Address	(0010,1040)	3	x	
Military Rank	(0010,1080)	3	x	
Ethnic Group	(0010,2160)	3	x	
Patient Comment	(0010,4000)	3	x	
Patient Medical				
Patient State	(0038,0500)	2	x	x
Pregnancy Status	(0010,21C0)	2	x	
Medical Alerts	(0010,2000)	2	x	x
Contrast Allergies	(0010,2110)	2	x	x
Special Needs	(0038,0050)	2	x	x
Smoking Status	(0010,21A0)	3	x	
Last Menstrual Date	(0010,21D0)	3	x	
Additional Patient History	(0010,21B0)	3	x	

a. "Scheduled Performing Physician's Name" is taken as default for "Performing Physicians's Name"

- Associated Real-World Activity - Get Worklist

With “Get Worklist” the contents of certain input fields of the Patient Registration UI are used to form a worklist request identifier. With the response data the Patient Registration dialog input is completed. The response data are additionally placed in the scheduler database.

- Proposed Presentation Contexts - Get Worklist

This RWE will propose the same Presentation Contexts as with “Update Worklist”. Please see table in section •.

- SOP Specific Conformance Statement - Get Worklist

Search Key Attributes of the Worklist C-FIND

The Siemens VSim DICOM worklist SCU supports “narrow worklist queries” with all required search keys. The following tables describe the “narrow query” search keys that the SCU supports.

Table 55. Patient based “narrow query” Search Key Attributes

Attribute name	Tags	Matching Key Type	Query Value
Scheduled Procedure Step			
Scheduled Procedure Step Sequence	[0040,0100]	R	
>Scheduled Performing Physician’s Name	[0040,0006]	R	input from UI or <zero length>
Requested Procedure			
Requested Procedure ID	[0040,1001]	O	input from UI or <zero length>
Imaging Service Request			
Accession Number	[0008,0050]	O	input from UI or <zero length>
Referring Physician’s Name	[0008,0090]	O	input from UI or <zero length>
Visit Status			
Current Patient Location	[0038,0300]	O	input from UI or <zero length>

Table 55. Patient based “narrow query” Search Key Attributes

Attribute name	Tags	Matching Key Type	Query Value
Patient Identification			
Patient’s Name	(0010,0010)	R	input from UI or <zero length>
Patient ID	(0010,0020)	R	input from UI or <zero length>

Return Key Attributes of the Worklist C-FIND

Please see list for “Update Worklist” RWE.

Status Codes of the Worklist C-FIND

The worklist SCU interprets following status codes:

Table 56. C-FIND Response Status

Service Status	Meaning	Status Codes (0000,0900)	Related Fields
Refused	Out of Resources	A700	(0000,0902)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	Cxxx	(0000,0901) (0000,0902)
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None

Table 56. C-FIND Response Status

Service Status	Meaning	Status Codes (0000,0900)	Related Fields
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

Implementation Model MPPS

The Modality Performed Procedure Step Service class defines an application-level class of service which facilitates the transfer of procedure, billing and radiation dose information from the imaging modality to the information system. The Performed Procedure Step is created and set by the AE and supplies the SCP with the information about a real-world procedure which is performed on the modality. The VSim DICOM Modality Performed Procedure Step application supports the MPPS service as SCU.

Application Data Flow Diagram

The VSim DICOM network implementation is a Windows 2000 application and acts as SCU for the Basic Worklist Service using the Modality Performed Procedure Step SOP Class.

:




Figure 10. Application data flow diagram - MPPS SCU

Functional Definitions of Application Entities

With registering a Patient (i.e. a Scheduled Procedure Step from Worklist), the VSim DICOM application will create a MPPS Instance and communicate it to the MPPS SCP.

Furthermore a manual update can be performed with the VSim MPPS user interface. Only there it is possible to set the state of the MPPS to “Completed” or “Discontinued”. If done so, the DICOM application will no longer allow updates on the related MPPS Instance.

The VSim will not only allow a “1:1 -relationship” of Scheduled Procedure Steps and Performed Procedure Steps, but also supports the “simple group-case”, “complex group-case” and “append case” from the respective IHE-scenarios.

The VSim will support creation of “unscheduled cases” by allowing MPPS Instances to be communicated for locally registered Patients.

Sequencing of real World Activities

An MPPS Update is only possible, if the MPPS Instance was created with Patient Registration before.

AE Specification MPPS

Modality Performed Procedure Step AE Specification

The Modality Performed Procedure Step SCU (Patient Registration and MPPS UI) provide information about a performed real-world Procedure to a remote SCP (Information System).
SIEMENS VSim DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Class as an SCU:

Table 57. SOP Classes as an SCU

SOP Class Name	SOP Class UID
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3

Association Establishment Policies

General

The creation of MPPS Instance is done automatically by VSim whenever a patient is registered for image acquisition through the Patient Registration dialog.
Further updates on the MPPS data can be done interactively from the related MPPS user interface. The MPPS “Complete” or “Discontinued” states can only be set from user interface.
The default PDU size used will be 28 KB.

Number of Associations

The Siemens VSim DICOM application initiates one association at a time to create or set MPPS instance.

Asynchronous Nature

The Siemens VSim DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

Implementation Identifying Information

The Siemens VSim DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- “SIEMENS_SWFVC20A”

Association Initiation Policy

The VSim DICOM Application Entity acts as a Service Class User (SCU) for the

- Modality Performed Procedure Step Service Class (to notify a RIS about status of a procedure while it is performed).

To do so, the VSim will issue a

- N-CREATE DIMSE according to the CREATE Modality Performed Procedure Step SOP Instance operation or a
- N-SET DIMSE to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.

Real World Activity

- Associated Real-World Activity - Patient registered

A patient is registered by the Patient Registration “Exam” action. From this event the trigger to create a MPPS Instance is derived. The related Instance is then immediately communicated to the configured (R)IS system. An association is established and the MPPS Instance is sent.

- Proposed Presentation Contexts - Patient registered

The Siemens VSim DICOM application will propose Presentation Contexts as shown in the following table:

Table 58. Proposed presentation contexts

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

- SOP Specific Conformance Statement- Patient registered

Attributes used for the Performed Procedure Step N-CREATE

The Siemens VSim DICOM Modality Performed Procedure Step SCU informs the remote SCP when the examination of a scheduled procedure step will be performed (i.e. the patient is registered). The N-CREATE message is sent when the examination is started with successful registration of the patient data. The following table describes the supported attributes of a N-CREATE message.

Table 59. Performed Procedure Step N-CREATE Attributes

Attribute name	Tag	Required Type	Value
SOP Common			
Specific Character Set	(0008,0005)	1C	from MWL or created
Performed Procedure Step Relationship			
Scheduled Step Attribute Sequence	(0040,0270)	1	
>Study Instance UID	(0020,000D)	1	from MWL or created
>Referenced Study Sequence	(0008,1110)	2	from MWL or <zero length>
>>Referenced SOP Class UID	(0008,1150)	1C	
>>Referenced SOP Instance UID	(0008,1155)	1C	
>Accession Number	(0008,0050)	2	from MWL or user input
>Placer Order Number/Imaging Service Request	(0040,2016)	3	from MWL or <zero length>
>Filler Order Number/Imaging Service Request	(0040,2017)	3	from MWL or <zero length>
>Requested Procedure ID	(0040,0001)	2	from MWL or user input
>Requested Procedure Description	(0032,1060)	2	from MWL or <zero length>
>Scheduled Procedure Step ID	(0040,0009)	2	from MWL or <zero length>
>Scheduled Procedure Step Description	(0040,0007)	2	from MWL or <zero length>
>Scheduled Action Item Sequence	(0040,0008)	2	from MWL or <zero length>
>>Code Value	(0008,0100)	1C	
>>Coding Scheme Designator	(0008,0102)	1C	

Table 59. Performed Procedure Step N-CREATE Attributes

Attribute name	Tag	Required Type	Value
>>Coding Scheme Version	(0008,0103)	3	
>>Code Meaning	(0008,0104)	3	
Patient's Name	(0010,0010)	2	from MWL or user input
Patient ID	(0010,0020)	2	from MWL or user input or created
Patient's Birth Date	(0010,0030)	2	from MWL or user input
Patient's Sex	(0010,0040)	2	from MWL or user input
Referenced Patient Sequence	(0008,1120)	2	from MWL or <zero length>
>Referenced SOP Class UID	(0008,1150)	1C	
>Referenced SOP Instance UID	(0008,1155)	1C	
Performed Procedure Step Informations			
Performed Procedure Step ID	(0040,0253)	1	from SPS ID or created
Performed Station AE Title	(0040,02410	1	own AE Title
Performed Station Name	(0040,0242)	2	own hostname
Performed Location	(0040,0243)	2	from SPS Location or <zero length>
Performed Procedure Step Start Date	(0040,0244)	1	created
Performed Procedure Step Start Time	(0040,0245)	1	created
Performed Procedure Step Status	(0040,0252)	1	“IN PROGRESS”
Performed Procedure Step Description	(0040,0254)	2	from SPS Description or <zero length>
Performed Procedure Type Description	(0040,0255)	2	<zero length>

Table 59. Performed Procedure Step N-CREATE Attributes

Attribute name	Tag	Required Type	Value
Procedure Code Sequence	(0008,1032)	2	from Requested Procedure Code or <zero length>
>Code Value	(0008,0100)	1C	
>Coding Scheme Designator	(0008,0100)	1C	
>Coding Scheme Version	(0008,0100)	3	
>Code Meaning	(0008,0100)	3	
Performed Procedure Step End Date	(0040,0250)	2	<zero length>
Performed Procedure Step End Time	(0040,0251)	2	<zero length>
Image Acquisition Results			
Modality	(0008,0060)	1	<modify/insert value>
Study ID	(0020,0010)	2	from Requested Procedure ID or created
Performed Action Item Code Sequence	(0040,0260)	2	from Scheduled Action ITeM Sequence or <zero length>
>Code Value	(0008,0100)	1C	
>Coding Scheme Designator	(0008,0100)	1C	
>Coding Scheme Version	(0008,0100)	3	
>Code Meaning	(0008,0100)	3	
Performed Series Sequence	(0040,0340)	2	
>Performing Physicians's Name	(0008,1050)	2C	from MWL or user input
>Operator's Name	(0008,1070)	2C	user input
>Series Instance UID	(0020,000E)	1C	created
>Series Description	(0008,103E)	2C	<zero length>
>Retrieve AE Title	(0008,0054)	2C	<zero length>

Table 59. Performed Procedure Step N-CREATE Attributes

Attribute name	Tag	Required Type	Value
>Referenced Image Sequence	(0008,1140)	2C	<zero length>
>Referenced Standalone SOP Instance Sequence	(0040,0220)	2C	<zero length>

Status Codes of the Performed PProcedure Step N-CREATE

The Performed Procedure Step SCU interprets following status codes:

Table 60. MPPS N-CREATE Response Status

Service Status	Meaning	Status Codes (0000,0900)
Failed	Processing Failure	0110
	No such attribute	0105
	Invalid attribute value	0106
	Duplicate SOP Instance	0111
	No such SOP Instance	0112
	No such SOP Class	0118
	Class Instance conflict	0119
	Missing attribute	0120
	Missing attribute value	0121
	Resource limitation	0213
Success	Matching is complete - No final Identifier is supplied	0000

- Associated Real-World Activity - MPPS UI-Update
- With the MPPS UI the status of the MPPS Instance can be set to “COMPLETED” or “DISCONTINUED”. There is no cyclic update during performance of the procedure.

- Proposed Presentation Contexts - MPPS UI-Update

This RWE will propose the same Presentation Contexts as with “Patient registered”. Please see table in section •.

- SOP Specific Conformance Statement - MPPS UI-Update

Attributes used for the Performed Procedure Step N-SET

The Siemens VSim DICOM Modality Performed Procedure Step SCU informs the remote SCP about the performed examination and its status. The N-SET message is sent only per ended examination (finished status “COMPLETED” or incomplete status “DISCONTINUED”). The following table describes the supported attributes of a N-SET message.

Table 61. Performed Procedure Step N-SET Attributes

Attribute name	Tag	Required Type	Value
Performed Procedure Step Informations			
Performed Procedure Step Status	(0040,0252)	3	“COMPLETED” or “DISCONTINUED”
Performed Procedure Step Description	(0040,0254)	3	from SPS Description or user input
Performed Procedure Type Description	(0040,0255)	3	user input
Procedure Code Sequence	(0008,1032)	3	from Requested Procedure Code or <modify and add specific data>
>Code Value	(0008,0100)	1C	
>Coding Scheme Designator	(0008,0100)	1C	
>Coding Scheme Version	(0008,0100)	3	
>Code Meaning	(0008,0100)	3	
Performed Procedure Step End Date	(0040,0250)	1	created
Performed Procedure Step End Time	(0040,0251)	1	created
Image Acquisition Results			

Table 61. Performed Procedure Step N-SET Attributes

Attribute name	Tag	Required Type	Value
Performed Action Item Code Sequence	(0040,0260)	3	from Scheduled Action Item Sequence or <modify and add specific data>
>Code Value	(0008,0100)	1C	
>Coding Scheme Designator	(0008,0100)	1C	
>Coding Scheme Version	(0008,0100)	3	
>Code Meaning	(0008,0100)	3	
Performed Series Sequence	(0040,0340)	2	
>Performing Physicians’s Name	(0008,1050)	2C	from MWL or user input
>Protocol Name	(0018,1030)	1C	from related SOP Instance
>Operator’s Name	(0008,1070)	2C	user input
>Series Instance UID	(0020,000E)	1C	from related SOP Instance
>Series Description	(0008,103E)	2C	from related SOP Instance
>Retrieve AE Title	(0008,0054)	2C	from Storage Commitment response or <zero length>
>Referenced Image Sequence	(0008,1140)	2C	Series related SOP Instances as items
>>Referenced SOP Class UID	(0008,1150)	1C	
>>Referenced SOP Instance UID	(0008,1155)	1C	
>Referenced Standalone SOP Instance Sequence	(0040,0220)	2C	<zero length>
All other attributes from Radiation Dose Module		3	
All other attributes from Billing and Material Code Module		3	

Status Codes of the Performed PROcedure Step N-SET

The Performed Procedure Step SCU interprets following status codes:

Table 62. MPPS N-CREATE Response Status

Service Status	Meaning	Status Codes (0000,0900)
Failed	Processing Failure: Performed procedure Step Object may no longer be updated	0110
	No such attribute	0105
	Invalid attribute value	0106
	No such SOP Instance	0112
	Invalid object instance	0117
	No such SOP class	0118
	Class Instance conflict	0119
	Missing attribute value	0121
	Resource limitation	0213
Success	Matching is complete - No final Identifier is supplied	0000

Communication Profiles

Supported Communication Stacks

The Siemens VSim DICOM application provide DICOM V3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard. The VSim target Operating System is Microsoft Windows 2000.

TCP/IP Stack

The Siemens VSim DICOM application uses the TCP/IP stack from the target operating system upon which it executes. It uses the MergeCOM-3 subroutine library from Merge Technologies Inc. that is based on a Berkeley socket interface.

API

The Siemens VSim DICOM application uses the MergeCOM library that is based on a TCP/IP socket interface.

Physical Media Support

The Siemens VSim DICOM application is indifferent to the physical medium over which TCP/IP executes; it inherits this from target operating system upon which it executes.

Extensions/Specializations/ Privatizations

Standard Extended/Specialized/Private SOPs

Please refer to the Appendix for further information on these topics. A detailed overview is given there.

Private Transfer Syntaxes

Not applicable.

Configuration

AE Title / Presentation Address Mapping

To ensure unique identification within the network the hostname should be used as part of the AE Titles (see examples below, hostname = name1). The string can be up to 16 characters long and must not contain any extended characters, only 7-bit ASCII characters (excluding Control Characters) are allowed according to DICOM Standard.

DICOM Verification

The Verification Service uses the AE configuration of the DICOM Service that is checked with the C-ECHO message. e.g. Verification will use the Storage AE, if initiated to check the configuration of a remote DICOM node.

DICOM Storage AE Title

The DICOM Storage application provides the application entity title which can be configured via Service UI:

e.g. STU_NAME1

The port number is set to the fixed value of
104

DICOM Query/Retrieve AE Title

The DICOM Query/Retrieve application uses the same application entity title as the DICOM Storage AE.

DICOM Print AE Title

The DICOM Print application provides the application entity title:

e.g. PRI_NAME1 (*No input of AETs starting with a numeric char is possible*)

The port number is set to the fixed value of
108

Configurable Parameters

The Application Entity Titles, host names and port numbers for the remote AE are configured using the Siemens VSim Service/Installation Tool. For each AET the list of services supported can be configured.

Storage, Storage Commitment and Query/Retrieve

The Siemens VSim Service/Installation Tool can be used to set the AET's, port-numbers, host-names, IP-addresses and capabilities for the remote nodes (SCP's). The service user can select transfer syntaxes, compression modes and query models for each SCP separately.

- a quality factor which determines the proposed transfer syntax in case that an user has initiated the C-STORE. By convention, 0 means: Only Uncompressed Transfer Syntax(es) are proposed, 100 means: Lossless Transfer Syntax is proposed, and any other value between 1 and 99 means that an JPEG Lossy Transfer Syntax is proposed. One Uncompressed Transfer Syntax will be proposed in any case. This parameter is general for all destination nodes.
- a “compression type supported” which determines the proposed transfer syntax in case that the C-STORE was initiated as a suboperation of an incoming C-MOVE-RQ. By convention, 0 means: Only Uncompressed Transfer Syntax(es) are proposed, 1 means: Lossless Transfer Syntax is proposed, and 2 means that an JPEG Lossy Transfer Syntax is proposed. One uncompressed transfer syntax will be proposed in any case. This parameter can be set for each configured destination node.

Additional configurable parameters for Storage Commitment are:

When acting as SCU:

- flag to indicate whether the association will be kept open to receive the response or to close the association and be prepared to receive the response on another association.
- time-out which defines how long the association of N-ACTION is kept to receive a N-EVENT-REPORT on the same association.
- time-out for N-EVENT-REPORT (applicability of transaction UID)
(default 1 h)

When acting as SCP:

- flag to indicate if an archive system is installed

Print

The Siemens VSim Service/Installation Tool can be used to configure the SCP (DICOM printer). The AET, host-name, IP-address and port-number can be set.

Modality Worklist

The Service application can be used to set the AETs, port numbers, host names, IP addresses, capabilities and time-outs for the remote nodes (SCPs)

Additional configurable parameters for Modality Worklist Query are:

- Query Waiting time - the time to wait for the C-FIND-RSP after sending the C-FIND-RQ (default 20 sec.)
- Max Query Match Number - the maximum number of entries accepted in one worklist (default is 200)
- Query Interval: the time between two C-FIND-RQ to the Hospital Information system (default is 60 min.)
- Broad Worklist Query behavior: two values are defined:
- Set the AE Title search attribute to the own AE Title, and the Modality search attribute to “*”.
- Set the Modality search attribute to the own modality and the AE Title search attribute to “*”.

Default Parameters

This installation tool also uses some default parameters:

- max PDU size set to 28672 Bytes (28 kB)
- time-out for accepting/rejecting an association request: 60 s
- time-out for responding to an association open/close request: 60 s
- time-out for accepting a message over network: 60 s
- time-out for waiting for data between TCP/IP-packets: 60 s

The time-outs for waiting for a Service Request/Response message from the remote node are as follows:

- for Storage SCP/SCU: 600 s
- for Storage Commitment SCU:
time-out for Response to N-ACTION: 600 s
- for Query/Retrieve SCP/SCU: 600 s
- for Print Management SCU:
- time-out for Response to N-SET-RQ: 240 s
- time-out for Response to other Requests: 60 s

Support of Extended Character Sets

The Siemens VSim DICOM application supports the ISO 8859 Latin 1 (ISO-IR 100) character set.

When there is a mismatch between the SCS tags (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM con-form:

Try to import with ISO_IR 100. If ISO_IR 100 fails, convert each illegal character to a '?'.



Application Profile Conformance Statement

This chapter will contain the Conformance Statement to all "Offline Media Application Profiles (inc. private extensions)" supported by the VSim archive options.

Those application profiles supported shall be:

- Standard Application Profiles
- Augmented Application Profiles
- *syngo* private Application Profile Class

Purpose

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

The applications described in this conformance statement is the SIEMENS OCS VSim, based on *syngo*® software. The VSim DICOM offline media storage service implementation acts as FSC,FSU and/or FSR for the specified application profiles and the related SOP Class instances.

Scope

This DICOM Conformance Statement refers to SIEMENS VSim software.

Definitions, Abbreviations

Definitions

DICOM	Digital Imaging and Communications in Medicine
DIMSE	DICOM Message Service Element
DIMSE-C	DICOM Message Service Element with Composite information objects

Abbreviations

ACR	American College of Radiology
AE	DICOM Application Entity
ASCII	American Standard Code for Information Interchange
DB	Database
DCS	DICOM Conformance Statement
DVD	Digital Versatile Disk
FSC	File Set Creator
FSR	File Set Reader
FSU	File Set Updater
IOD	DICOM Information Object Definition
ISO	International Standard Organization

Med-WS	Multimodality Workstation
MOD	Magneto-Optical Disk
NEMA	National Electrical Manufacturers Association
O	Optional Key Attribute
PDU	DICOM Protocol Data Unit
R	Required Key Attribute
RAM	Random Access Memory
RWA	Real-World Activity
U	Unique Key Attribute

References

[2] Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-15, 2000

Connectivity and Inter operability


This Conformance Statement by itself does not guarantee successful inter operability of SIEMENS equipment with non-SIEMENS equipment. The user (user’s agent) should be aware of the following issues:

- **Inter operability**
Inter operability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networks environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee inter operability of SIEMENS equipment with non-SIEMENS equipment. It is the user’s responsibility to analyze thoroughly the application requirements and to specify a solution that integrates SIEMENS equipment with non-SIEMENS equipment.
- **Validation**
SIEMENS equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.
Where SIEMENS equipment is linked to non-SIEMENS equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation test will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user’s agent) to specify the appropriate test suite and to carry out the additional validation tests.

- **New versions of the DICOM Standard**

The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. SIEMENS is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, SIEMENS reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-SIEMENS provider linking to SIEMENS equipment, also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into SIEMENS equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

Application Data Flow Diagram



The DICOM archive application will serve as an interface to the CD-R and (optional) MOD offline medium device. It serves interfaces to include the offline media directory into the browser and to copy SOP instances to a medium or retrieve SOP Instances from medium into local storage.

The DICOM Archive application will support the 120mm CD-R medium, the 130mm 2.3GB R/W MOD, the 130mm 4.1 GB R/W MOD, the 120 mm 4.7 GB DVD-RAM and 1.44 MB diskette media.

The FSU role will update new SOP Instances only to media with pre-existing File-sets conforming to the Application Profiles supported.

The contents of the DICOMDIR will be temporarily stored in Archive-Database.

Functional definitions of AE's

The VSim DICOM offline media storage application consist of the DICOM Archive application entity serving all interfaces to access offline media. The DICOM Archive application is capable of

- 1 creating a new File-set onto an unwritten medium (Export to...).
- 2 updating an existing File-set by writing new SOP Instances onto the medium (Export to...).
- 3 importing SOP Instances from the medium onto local storage (Import)
- 4 reading the File-sets DICOMDIR information into temporary database and pass it to display applications.

Sequencing of Real World Activities

The DICOM Archive application will not perform updates before the Directory information of the DICOMDIR is completely read.

When performing updates, the SOP instances are checked for existence before updating. Duplicate instances will be avoided.

File Meta Information Options

The Implementation Class UID is:

- 1.3.12.2.1107.5.9.200010101

and an Implementation Version Name of

- "SIEMENS_SWFVB10A"

DICOM Archive Specification

The DICOM Archive provides Standard conformance to Media Storage Service Class (Interchange Option). In addition Augmented conformance is provided to store extra data attributes important for the full feature support of the VSim product SW. Details are listed in Table 63..

Table 63. Application profiles, Activities, and Roles for DICOM Archive

Application Profiles Supported	Real World Activity	Role	SC Option
PRI-SYNGO-CD	Browse Directory Information	FSR	Interchange
PRI-SYNGO-MOD23 (optional)	Import into local Storage	FSR	Interchange
PRI-SYNGO-MOD41(optional)			
PRI-SYNGO-DVD-RAM47	Export to local archive media	FSC, FSU	Interchange
AUG-GEN-CD			
AUG-CTMR-MOD650 *1	Browse Directory Information	FSR	Interchange
AUG-CTMR-MOD12 *1			
AUG-CTMR-MOD23 *1			
AUG-CTMR-CD *1			
AUG-XA1K-CD *1			
STD-GEN-CD			
STD-CTMR-MOD650	Import into local Storage	FSR	Interchange
STD-CTMR-MOD12			
STD-CTMR-MOD23			
STD-CTMR-CD			
<modify STD-XABC-CD>			
<modify STD-XA1K-CD>			
TD-US-zz-yF-xxxxxx *2	Export to local archive media	FSC, FSU	Interchange
STD-WVFM-GEN-FD			

*1 - With no Private SOP Class used, the PRI-SYNGO-CD Profile definitions are appropriate to describe the Augmentation of the related STD Profiles.

*2 - All combinations of the following values for xx, yF and xxxxxx are supported: yF={SF|MF}, xx={ID|SC|CC}, xxxxxx={FLOP|MOD128|MOD230|MOD540|MOD650|MOD12|MOD23|DVD-RAM47|CDR}

On VSim based Workstation Products the Private Extended *syngo* Profile (PRI-SYNGO-CD or PRI-SYNGO-DVD-RAM47 or optional PRI-SYNGO-MOD23/PRI-SYNGO-MOD41) will be preferably used by the system. The General Purpose Interchange Profile (STD-GEN-CD), Basic Cardiac Profile (STD-XABC-CD) and will be supported with read capability of the related media.

File Meta Information for the Application Entity

The Source Application Entity Title is set by configuration. See Chapter "Configuration" for details.

Real-World Activities for this Application Entity

Real-World Activity: Browse Directory Information

The DICOM Archive application acts as FSR using the interchange option when requested to read the media directory.

The DICOM archive application will read the DICOMDIR and insert those directory entries, that are valid for the application profiles supported, into a local database. The database can then be used for browsing media contents.

Note IconImageSQ is also supported in DICOMDIR. But only those IconImages with BitsAllocated (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the Browser.

- Application Profiles for the RWA: Browse Directory Information

See Table 63. for the Application Profiles listed that invoke this Application Entity for the Browse Directory Information RWA.

Real-World Activity: Import into local Storage

The DICOM Archive application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the local storage.

The SOP Instance selected from the media directory will be copied into the local storage. Only SOP Instances, that are valid for the application profile supported and are listed as supported by the Storage SCP Conformance section (Network DCS,), can be retrieved from media storage. This is due to the fact that the Browse Directory Information will filter all SOP Instances not matching the Application profiles supported.

During operation no "Attribute Value Precedence" is applied to the SOP Instances. Detached Patient Management is not supported.

For media conforming to the STD-GEN-CD Profile the following SOP classes will be supported as an FSR:

Table 64. STD-GEN-CD Supported SOP Classes for FSR role

Information Object Definition	Service Object Pair Class UID	Transfer Syntax and UID
CR Image	1.2.840.10008.5.1.4.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
DX Image- For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
DX Image- For Presentation	1.2.840.10008.5.1.4.1.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
MG Image- For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
MG Image- For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
NM Image	1.2.840.10008.5.1.4.1.1.20	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
PET Image	1.2.840.10008.5.1.4.1.1.128	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
RT Dose	1.2.840.10008.5.1.4.1.1.481.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
RT Image	1.2.840.10008.5.1.4.1.1.481.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1

Table 64. STD-GEN-CD Supported SOP Classes for FSR role

Information Object Definition	Service Object Pair Class UID	Transfer Syntax and UID
RT Plan	1.2.840.10008.5.1.4.1.1.481.5	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
RT Structure Set	1.2.840.10008.5.1.4.1.1.481.3	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Ultrasound Multi-Frame Image (Retired)	1.2.840.10008.5.1.4.1.1.3	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Ultrasound Multi-Frame Image	1.2.840.10008.5.1.4.1.1.3.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
12-lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1

Table 64. STD-GEN-CD Supported SOP Classes for FSR role

Information Object Definition	Service Object Pair Class UID	Transfer Syntax and UID
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
CSA Non-Image	1.3.12.2.1107.5.9.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1

- Application Profiles for the RWA: Import into local Storage

See Table 63. for the Application Profiles listed that invoke this Application Entity for the Import into Local Storage RWA.

Real-World Activity: Export to local Archive Media

The DICOM Archive application acts as FSU (for media with existing DICOM file-set) or FSC (media not initialized) using the interchange option when requested to copy SOP Instances from the local storage to local Archive medium.

The DICOM Archive application will receive a list of SOP Instances to be copied to the local archive medium. According to the state of the medium inserted (new medium, Medium with DICOM file-set) the validity of the SOP Instances according to the applicable profile is checked. Only valid SOP Instances are accepted.

When the DICOM archive application is requested to copy SOP Instances, the preferred application profile according configuration (AUG-XA1K-CD or PRI-SYNGO-xxx) will be used to validate and copy the referred SOP Instances. When creating a new file-set no Descriptor File will be allocated and the related ID is not used.

The DICOM archive application will not close the medium.

- Application Profiles for the RWA: Export to local Archive Media

See Table 63. for the Application Profiles listed that invoke this Application Entity for the Export to local Archive Media RWA.

Augmented Application Profiles

AUG-GEN-CD

With no private Siemens Non-Images stored onto Medium, the definitions of the PRI-SYNGO-CD Profile are applicable to denote the augmentations for the STD-GEN-CD Standard Profile.

Storage of Private Information Objects will only be supported with reference to a Private Application Profile (see next section).

AUG-CTMR-xxxxx

With no private Siemens Non-Images stored onto Medium, the definitions of the PRI-SYNGO-CD Profile are applicable to denote the augmentations for the STD-CTMR-MOD650, STD-CTMR-MOD12, STD-CTMR-MOD23 and STD-CTMR-CDR Standard Profiles.

Storage of Private Information Objects will only be supported with reference to a Private Application Profile (see next section).

AUG-XA1K-CD

With no private Siemens Non-Images stored onto Medium, the definitions of the PRI-SYNGO-CD Profile are applicable to denote the augmentations for the STD-XA1K-CD Standard Profile.

Storage of Private Information Objects will only be supported with reference to a Private Application Profile (see next section).

Private Application Profiles

This document defines an Application Profile Class for "*syngo*® speaking¹" modalities or applications.

The identifier for this class shall be PRI-SYNGO. This class is intended to be used for interchange of extended and private Information Objects via CD-R or re-writable magneto-optical disk (MOD) offline media between dedicated acquisition or workstation modalities build from a common *syngo* architecture.

1.

syngo is registered trademark of Siemens AG.

The specific application profiles in this class are shown in Table 65.

Table 65. PRI-SYNGO Universal VSim Class Profiles

Application Profile	Identifier	Description
"syngo speaking" System on CD-R	PRI-SYNGO-CD	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Siemens Non-Image IOD).
"syngo speaking" System on 2.3 GB MOD	PRI-SYNGO-MOD23	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Siemens Non-Image IOD).
"syngo speaking" System on 4.1 GB MOD	PRI-SYNGO-MOD41	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Siemens Non-Image IOD).
"syngo speaking" System on Floppy Disk	PRI-SYNGO-FD	Handles interchange of Waveform SOP isntances and privately defined SOP Instances (Siemens Non-Image IOD).
"syngo speaking" System on 4.7 GB DVD-RAM	PRI-SYNGO-DVD-RAM47	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Siemens Non-Image IOD).

Equipment claiming conformance for this *syngo* Application Profile Class shall make a clear statement on handling of the private defined SOP Instances.

Extensions, Specializations and Privatization of SOP Classes and Transfer Syntaxes

The SOP Classes listed refer in majority to those created by the equipment to which this conformance Statement is related to. For SOP classes not listed in this section, please refer to the Storage section of the DICOM Conformance Statement of the product. This will include all SOP Instances that can be received and displayed and therefor will be included into offline media storage even though these SOP Instances are not created by the equipment serving the Media Storage Service.

SOP Specific Conformance Statement for Basic Directory

Extension, Specialization for SIEMENS Non-Image Objects

According to the PRI-SYNGO Application Profile Class the usage of the Private Creator UIDs and further optional keys for the Directory Records referring to SIEMENS Non-Image Objects are listed in the following tables.

Table 66. Basic Directory Extension for Non-Image Objects

Attribute	Tag	Value used
Private Record UID	[0004,1432]	1.3.12.2.1107.5.9.1
SOP Class UID	[0008,0016]	1.3.12.2.1107.5.9.1

For those "Non-Images" no Icon Image Sequence will be generated.

AE Title Mapping

DICOM Media Storage AE Title

The DICOM Storage application provides the application entity title:

CsaImageManager

Support of Extended Character Sets

The Siemens VSim DICOM application supports the ISO 8859 Latin 1 (ISO-IR 100) character set.

When there is a mismatch between the SCS tags (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM con-form:

Try to import with ISO_IR 100. If ISO_IR 100 fails, convert each illegal character to a '?'.



VSim Offline Media Application Profile

This contains a VSim specific Application Profile Class.

Structure of this Application Profile is defined in Part 11 of the 2000 DICOM Standard.

It is needed to describe the requirements for Offline Media Storage of the private IOD (Non-Image IOD).

This application profile facilitates the interchange of original acquired and derived images and private data related to them. Typical media interchange would be from in-lab acquisition equipment to dedicated workstations and archive systems with specific extensions to handle the private data objects (in both directions).

Additionally, images (from MR, CT, US) used to prepare procedures, multi-modality images (e.g. integrated US) and images derived from primary diagnostic images, such as annotations, quantitative analysis images, reference images, screen capture images may be interchanged via this profile.

Roles and Service Class Options

This Application Profile uses the Media Storage Service Class defined in PS 3.4 with the Interchange Option.

The Application Entity shall support one or more of the roles of File Set Creator (FSC), File Set Reader (FSR), and File Set Updater (FSU), defined in PS 3.10.

File Set Creator

The Application Entity acting as a File-Set Creator generates a File Set under the PRI-SYNGO Application Profiles.

File Set Creators shall be able to generate the Basic Directory SOP Class in the DICOMDIR file with all the subsidiary Directory Records related to the Image SOP Classes and Private SOP Classes stored in the File Set.

In case of PRI-SYNGO-CD Profile, the FSC shall offer the ability to either finalize the disc at the completion of the most recent write session (no additional information can be subsequently added to the disc) or to allow multi-session (additional information may be subsequently added to the disc).

Note A multiple volume (a logical volume that can cross multiple physical media) is not supported by this Application Profile Class. If a set of Files, e.g., a Study, cannot be written entirely on one CD-R, the FSC will create multiple independent DICOM File-Set such that each File-Set can reside on a single CD-R medium controlled by its individual DICOMDIR file. The user of the FSC can opt to use written labels on the discs to reflect that there is more than one disc for this set of files (e.g., a Study).

File Set Reader

The role of the File Set Reader shall be used by Application Entities which receive the transferred File Set.

File Set Readers shall be able to read all the defined SOP Instances files defined for the specific Application Profiles to which a conformance claim is made, using all the defined Transfer Syntaxes.

File Set Updater

The role of the File Set Updater shall be used by Application Entities which receive a transferred File Set and update it by the addition of processed information.

File Set Updaters shall be able to read and update the DICOMDIR file. File-Set Updaters do not have to read the image/private information objects. File-Set Updaters shall be able to generate any of the SOP Instances files defined for the specific Application Profiles to which a conformance claim is made, and to read and update the DICOMDIR file.

In case of PRI-SYNGO-CD Profile, the FSU shall offer the ability to either finalize a disc at the completion of the most recent write session (no additional information can be subsequently added to the disc) or to allow multi-session (additional information may be subsequently added to the disc).

Note (for CD-R) If the disc has not been finalized, the File-Set Updater will be able to update information assuming there is enough space on the disc to write a new DICOMDIR file, the information, and the fundamental CD-R control structures. CD-R control structures are the structures that inherent to the CD-R standards; see PS 3.12

SOP Classes and Transfer Syntaxes

This Application Profiles are based on the Media Storage Service Class with the Interchange Option.

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
Basic Directory	1.2.840.10008.1.3.10	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
CR Image	1.2.840.10008.5.1.4.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
CR Image	1.2.840.10008.5.1.4.1.1.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
CR Image	1.2.840.10008.5.1.4.1.1.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
CR Image	1.2.840.10008.5.1.4.1.1.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
CR Image	1.2.840.10008.5.1.4.1.1.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
CT Image	1.2.840.10008.5.1.4.1.1.2	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
CT Image	1.2.840.10008.5.1.4.1.1.2	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
CT Image	1.2.840.10008.5.1.4.1.1.2	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
DX Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
DX Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
DX Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
DX Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
DX Image - For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
DX Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
DX Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
DX Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
DX Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
DX Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
MG Image - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
MG Image - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
MG Image - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
MG Image - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
MG Image - For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
MG Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
MG Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
MG Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
MG Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
MG Image - For Presentation	1.2.840.10008.5.1.4.1.1.1.2	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
MR Image	1.2.840.10008.5.1.4.1.1.4	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
MR Image	1.2.840.10008.5.1.4.1.1.4	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
MR Image	1.2.840.10008.5.1.4.1.1.4	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
NM Image	1.2.840.10008.5.1.4.1.1.20	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
NM Image	1.2.840.10008.5.1.4.1.1.20	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
NM Image	1.2.840.10008.5.1.4.1.1.20	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
NM Image	1.2.840.10008.5.1.4.1.1.20	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
NM Image	1.2.840.10008.5.1.4.1.1.20	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
PET Image	1.2.840.10008.5.1.4.1.1.128	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
PET Image	1.2.840.10008.5.1.4.1.1.128	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
PET Image	1.2.840.10008.5.1.4.1.1.128	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
PET Image	1.2.840.10008.5.1.4.1.1.128	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
PET Image	1.2.840.10008.5.1.4.1.1.128	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
RT Dose	1.2.840.10008.5.1.4.1.1.481.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
RT Dose	1.2.840.10008.5.1.4.1.1.481.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
RT Image	1.2.840.10008.5.1.4.1.1.481.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
RT Image	1.2.840.10008.5.1.4.1.1.481.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
RT Image	1.2.840.10008.5.1.4.1.1.481.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
RT Image	1.2.840.10008.5.1.4.1.1.481.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
RT Image	1.2.840.10008.5.1.4.1.1.481.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
RT Plan	1.2.840.10008.5.1.4.1.1.481.5	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
RT Plan	1.2.840.10008.5.1.4.1.1.481.5	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
RT Structure Set	1.2.840.10008.5.1.4.1.1.481.3	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
RT Structure Set	1.2.840.10008.5.1.4.1.1.481.3	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	No	Yes	No
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	No	Yes	No
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	No	Yes	No
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	No	Yes	No
Ultrasound Image (Retired)	1.2.840.10008.5.1.4.1.1.6	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
Ultrasound Multiframe Image (Retired)	1.2.840.10008.5.1.4.1.1.3	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	No	Yes	No
Ultrasound Multiframe Image (Retired)	1.2.840.10008.5.1.4.1.1.3	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	No	Yes	No
Ultrasound Multiframe Image (Retired)	1.2.840.10008.5.1.4.1.1.3	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	No	Yes	No
Ultrasound Multiframe Image (Retired)	1.2.840.10008.5.1.4.1.1.3	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	No	Yes	No
Ultrasound Multiframe Image (Retired)	1.2.840.10008.5.1.4.1.1.3	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	Yes	Yes	Yes
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	Yes	Yes	Yes
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	RLE Lossless (decompression) 1.2.840.10008.1.2.5	No	Yes	No
Waveform SOP Classes	1.2.840.10008.5.1.4.1.1.9.1.1 1.2.840.10008.5.1.4.1.1.9.1.2 1.2.840.10008.5.1.4.1.1.9.1.3 1.2.840.10008.5.1.4.1.1.9.2.1 1.2.840.10008.5.1.4.1.1.9.3.1 1.2.840.10008.5.1.4.1.1.9.4.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes

Table 67. PRI-SYNGO SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FS R	FS U
Waveform SOP Classes	1.2.840.10008.5.1.4.1.1.9.1.1 1.2.840.10008.5.1.4.1.1.9.1.2 1.2.840.10008.5.1.4.1.1.9.1.3 1.2.840.10008.5.1.4.1.1.9.2.1 1.2.840.10008.5.1.4.1.1.9.3.1 1.2.840.10008.5.1.4.1.1.9.4.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	Yes	Yes	Yes
CSA Non-Image	1.3.12.2.1107.5.9.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	Yes	Yes	Yes
Detached Patient Management	1.2.840.10008.3.1.2.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	No	No	No

Physical Media and Media Formats

The PRI-SYNGO-CD Profile requires the 120mm CD-R physical media with the ISO/IEC 9660 Media Format, as defined in PS3.12.

The PRI-SYNGO-MOD23 Profile requires the 130mm 2.3 GB R/W MOD physical medium with the PCDOS Media Format, as defined in PS3.12.

The PRI-SYNGO-MOD41 Profile requires the 130mm 4.1 GB R/W MOD physical medium with the PCDOS Media Format, as approved by DICOM Suppl. 62.

The PRI-SYNGO-FD Profile requires the 1.44 MB diskette physical medium with the PCDOS Media Format, as defined in PS3.12.

The PRI-SYNGO-DVD-RAM47 Profile requires the 120 mm 4.7 GB DVD-RAM physical medium with the UDF - Universal Disk Format version 1.5, as defined in PS3.12.

Directory Information in DICOMDIR

Conforming Application Entities shall include in the DICOMDIR File the Basic Directory IOD containing Directory Records at the Patient and subsidiary levels appropriate to the SOP Classes in the File-set. All DICOM files in the File-set incorporating SOP Instances defined for the specific Application profile shall be referenced by Directory Records.

Note DICOMDIRs with no directory information are not allowed by this Application Profile. Privately defined IODs will be referenced by "PRIVATE" Directory Records.

Basic Directory IOD Specialization

This Application Profile makes use of optional attributes of the Basic Directory IOD to support recognition of Patient’s Storage Service request results in spanning multiple volumes (file sets). Therefore the File Set Descriptor File can be used and is then referenced by optional Basic Directory IOD attributes. If existent, the specified Descriptor File may be used by FSR applications. Any FSU, FSC shall make a clear Statement if the Descriptor File mechanism is used according to the specialization defined in this Application Profile.

The Descriptor Files shall have the following contents:

One single Line without any control-characters and according to the Basic Character-Set having the following defined text:

"MULTIVOLUME: xx of yy"

xx, yy are replaced by the actual Number of the volume (xx) and the Total Number of Volumes in the set (yy).

If used, the Descriptor File shall have the File ID "README" and reside in same directory level as the DICOMDIR. It is referenced by the attribute [0004,1141] File-set Descriptor File ID having the defined content of "README".

Additional Keys

File-set Creators and Updaters are required to generate the mandatory elements specified in PS 3.3, Annex F of the DICOM Standard. Table 68.<Italic>PRI-SYNGO Additional DICOMDIR Keys specifies the additional associated keys. At each directory record level other additional data elements can be added, but it is not required that File Set Readers be able to use them as keys. Refer to the Basic Directory IOD in PS 3.3.

Table 68. PRI-SYNGO Additional DICOMDIR Keys

Key Attribute	Tag	Directory Record Level	Type	Notes
Date of Birth	[0010,0030]	PATIENT	2C	required, if present in SOP Instance
Patient’s Sex	[0010,0040]	PATIENT	2C	required, if present in SOP Instance
Series Date	[0008,0021]	SERIES	3	
Series Time	[0008,0031]	SERIES	3	
Institute Name	[0008,0080]	SERIES	2C	required, if present in SOP Instance

Table 68. PRI-SYNGO Additional DICOMDIR Keys

Key Attribute	Tag	Directory Record Level	Type	Notes
Institution Address	[0008,0081]	SERIES	2C	required, if present in SOP Instance
Series Description	[0008,103E]	SERIES	3	
Performing Physician’s Name	[0008,1050]	SERIES	2C	required, if present in SOP Instance
Image Type	[0008,0008]	IMAGE	1C	required, if present in SOP Instance
SOP Class UID	[0008,0016]	IMAGE	3	
SOP Instance UID	[0008,0018]	IMAGE	3	
Image Date	[0008,0023]	IMAGE	3	
Image Time	[0008,0033]	IMAGE	3	
Referenced Image Sequence	[0008,1140]	IMAGE	1C	required, if present in SOP Instance
>Referenced SOP Class UID	[0008,1150]	IMAGE		
>Referenced SOP Instance UID	[0008,1155]	IMAGE		
Image Position (Patient)	[0020,0032]	IMAGE	2C	required, if present in SOP Instance
Image Orientation (Patient)	[0020,0037]	IMAGE	2C	required, if present in SOP Instance
Frame of Reference UID	[0020,0052]	IMAGE	2C	required, if present in SOP Instance
Rows	[0028,0010]	IMAGE	3	
Columns	[0028,0011]	IMAGE	3	
Pixel Spacing	[0028,0030]	IMAGE	1C	required, if present in SOP Instance
Calibration Image	[0050,0004]	IMAGE	2C	required, if present in SOP Instance

Table 68. PRI-SYNGO Additional DICOMDIR Keys

Key Attribute	Tag	Directory Record Level	Type	Notes
Icon Image Sequence	[0088,0200]	IMAGE	3	required for Image SOP Classes
>Samples per Pixel	[0028,0002]			1
>Photometric Interpretation	[0028,0004]			MONOCHROME2
>Rows	[0028,0010]			128 for XA, 64 for others
>Columns	[0028,0011]			128 for XA, 64 for others
>Bits Allocated	[0028,0100]			8
>Bits Stored	[0028,0101]			8
>High Bit	[0028,0102]			7
>Pixel Representation	[0028,0103]			0
>Pixel Data	[7FE0,0010]			Icon Image
Curve Number	[0020,0024]	CURVE	1C	required, if present in SOP Instance

Private Directory Record Keys

Private Directory Records are supported by this Application Profile Class at the following Level - IMAGE. The PRIVATE Directory Records will have required elements in addition to the mandatory elements specified in PS 3.3.

Table 69. will list the additional required keys for PRIVATE Directory Records.

Table 69. PRI-SYNGO Additional PRIVATE Record Keys

Key Attribute	Tag	Type	Notes
Private Record UID	(0004,1432)	1	see Conformance Statement
SOP Class UID	[0008,0016]	1C	required, if present in SOP Instance
SOP Instance UID	[0008,0018]	1C	required, if present in SOP Instance
Image Type	[0008,0008]	3	

Table 69. PRI-SYNGO Additional PRIVATE Record Keys

Key Attribute	Tag	Type	Notes
Acquisition Date	[0008,0022]	3	
Acquisition Time	[0008,0032]	3	
Acquisition Number	[0020,0012]	3	
CSA Data Type	[0029,xx08]	1	private owner code = SIEMENS CSA NON-IMAGE
CSA Data Version	[0029,xx09]	3	private owner code = SIEMENS CSA NON-IMAGE

Icon Images

Directory Records of type SERIES or IMAGE may include Icon Images. The Icon Image pixel data shall be as specified in PS 3.3 "Icon Image Key Definition", and restricted such that Bits Allocated (0028,0100) and Bits Stored (0028,0101) shall be equal 8, and Rows (0028,0010) and Columns (0028,0011) shall be equal to 128 for XA Images and 64 for all other Images. The Photometric Interpretation (0028,0004) shall always be restricted to "MONOCHROME2".

PRIVATE Directory Records will not contain Icon Image information.

Other Parameters

This section defines other parameters common to all specific Application Profiles in the PRI-SYNGO class which need to be specified in order to ensure interoperable media interchange.

Multiframe JPEG Format

The JPEG encoding of pixel data shall use Interchange Format (with table specification) for all frames.

Attribute Value Precedence

The values of attributes contained in a Detached Patient Management SOP Instance referenced by a DICOMDIR PATIENT Directory Record shall take precedence over the values of those attributes contained in the SOP Instance referenced by a subsidiary Directory Record. The DICOMDIR Directory Records shall have key attribute values in accordance with this precedence.

Note This allows patient identification and demographic information to be updated without changing the composite Image IOD files. The DICOMDIR file thus is critical in establishing the link between the updated information and the image. As an example, at the time an

Image file was written, the patient's name therein was incorrect, or inconsistent with the Hospital Information System records. Subsequently, a Detached Patient management file with the corrected name is added to the file-set. The FSR should use the name from the Patientfile rather than in the image File.

A

Appendix

SIEMENS Private Non-Image IOD

For encoding binary data-streams not representing image data, Siemens has created a private "Non-Image IOD" according to the rules governed by the DICOM Standard. The following section will roll-out the definition of this Private IOD. It can be communicated with Network Storage Service and Offline Media Storage Services.


The Siemens "Non-Image IOD" is identified by a private Non-Image Storage SOP Class UID of

"1.3.12.2.1107.5.9.1"

Siemens Non-Image IOD - E-R Model

The E-R model in A.1.2 depicts those components of the DICOM Information Model which directly refer to the Siemens Non-Image IOD. The Frame of Reference IE, Overlay IE, Modality Lookup-Table IE, VOI Lookup-Table IE and Curve IE are not components of the Siemens Non-Image IOD.

Figure 11. DICOM Information Model



Siemens Non-Image IOD - Module Table

Table 70. CSA Non-Image IOD Modules

IE	Module	Reference	Usage
Patient	Patient	[1] PS3.3 C.7.1.1	M

Table 70. CSA Non-Image IOD Modules

Study	General Study	[1] PS3.3 C.7.2.1	M
	Patient Study	[1] PS3.3 C.7.2.2	U
Series	General Series	[1] PS3.3 C.7.3.1	M
Equipment	General Equipment	[1] PS3.3 C.7.5.1	U
CSA	CSA Image Header		U
	CSA Series Header		U
	MEDCOM Header		U
	CSA Non-Image		M
	SOP Common	[1] PS3.3 C.12.1	M

Siemens Non-Image IOD - Modules

CSA Non-Image Module

The table in this section contains private IOD Attributes that describe CSA Non-Images..

Table 71. CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
Image Type	(0008,0008)	-	3	Image identification characteristics.
Acquisition Date	(0008,0022)	-	3	The date the acquisition of data that resulted in this data set started.
Acquisition Time	(0008,0032)	-	3	
Conversion Type	(0008,0064)	-	3	Describes the kind of image conversion. Defined Terms: DV = Digitized Video, DI = Digital Interface, DF = Digitized Film, WSD = Workstation.

Table 71. CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
Referenced Image Sequence	(0008,1140)	-	3	A sequence which provides reference to a set of Image SOP Class/Instance identifying other images significantly related to this data set. Encoded as sequence of items: (0008,1150) and (0008,1155).
Derivation Description	(0008,2111)	-	3	A text description of how this data set was derived.
Source Image Sequence	(0008,2112)	-	3	A Sequence which identifies the set of Image SOP Class/Instance pairs of the Images which were used to derive this data set. Zero or more Items may be included in this Sequence. Encoded as sequence of items: (0008,1150) and (0008,1155).
Patient Position	(0018,5100)	-	3	Patient position descriptor relative to the equipment.
Acquisition Number	(0020,0012)	-	3	A number identifying the single continuous gathering of data over a period of time which resulted in this data set.
Image Number	(0020,0013)	-	3	A number that identifies this data set.
Frame of Reference UID	(0020,0052)	-	3	Uniquely identifies the frame of reference for a Series.
Image Comments	(0020,4000)	-	3	User-defined comments about the image.
Quality Control Image	(0028,0300)	-	3	Indicates whether or not this image is a quality control or phantom image. If this Attribute is absent, then the image may or may not be a quality control or phantom image. Enumerated Values: YES, NO.

Table 71. CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
Burned In Annotation	(0028,0301)	-	3	Indicates whether or not image contains sufficient burned in annotation to identify the patient and date the image was acquired. If this Attribute is absent, then the image may or may not contain burned in annotation. Enumerated Values: YES, NO.
Lossy Image Compression	(0028,2110)	-	3	Specifies whether an Image has undergone lossy compression. Enumerated Values: 00 = Image has NOT been subjected to lossy compression, 01 = Image has been subjected to lossy compression.
Lossy Image Compression Ratio	(0028,2112)	-	3	Describes the approximate lossy compression ratio(s) that have been applied to this image. May be multi valued if successive lossy compression steps have been applied.
CSA Data Type	(0029,xx08)	SIEMENS CSA NON-IMAGE	1	CSA Data identification characteristics. Defined Terms: BSR REPORT = BSR Study Report Data 3D EDITOR 3D FLY PATH = Fly Through Data 3D FLY VRT = Fly Through Data 3D FUSION MATRIX = Fusion Data RAW DATA NUM 4 = NUMARIS/4 Raw Data RAW DATA SOM 5 = SOMARIS/5 Raw Data RT3D CONFIG = InSpaceIS Data SPEC NUM 4 = NUMARIS/4 Spectroscopy
CSA Data Version	(0029,xx09)	SIEMENS CSA NON-IMAGE	3	Version of CSA Data Info (0029,xx10) format and CSA Non-Image Data (7FE1,xx10) format.

Table 71. CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
CSA Data Info	(0029,xx10)	SIEMENS CSA NON-IMAGE	3	Information to describe the CSA Data (7FE1,xx10).
CSA Data	(7FE1,xx10)	SIEMENS CSA NON-IMAGE	2	Binary data as byte stream.

Siemens Standard Extended Modules

Table 72. CSA Image IOD Modules

IE	Module	Reference	Usage	Note
Image or Document	CSA Image Header		U	private GG information
	CSA Series Header		U	
	MEDCOM Header		U	private <i>syngo</i> information
	MEDCOM OOG		U	if object graphics is attached to image

CSA Image Header Module

The table in this section contains private IOD Attributes that describe the CSA Image Header..

Table 73. CSA Image Header Module

Attribute Name	Tag	Owner	Type	Notes
CSA Image Header Type	(0029,xx08)	SIEMENS CSA HEADER	1	CSA Image Header identification characteristics. Defined Terms: NUM 4 = NUMARIS/4 SOM 5 = SOMARIS/5
CSA Image Header Version	(0029,xx09)	SIEMENS CSA HEADER	3	Version of CSA Image Header Info (0029,xx10) format.
CSA Image Header Info	(0029,xx10)	SIEMENS CSA HEADER	3	Manufacturer model dependent information.

CSA Series Header Module

The table in this section contains private IOD Attributes that describe the CSA Series Header..

Table 74. CSA Series Header Module

Attribute Name	Tag	Owner	Type	Notes
CSA Series Header Type	(0029,xx18)	SIEMENS CSA HEADER	1	CSA Series Header identification characteristics. Defined Terms: NUM 4 = NUMARIS/4
CSA Series Header Version	(0029,xx19)	SIEMENS CSA HEADER	3	Version of CSA Series Header Info (0029,xx20) format.
CSA Series Header Info	(0029,xx20)	SIEMENS CSA HEADER	3	Manufacturer model dependent information.

MEDCOM Header Module

The table in this section contains private IOD Attributes that describe MEDCOM Header.

Table 75. MEDCOM Header Module

Attribute Name	Tag	Owner	Type	Notes
MedCom Header Type	(0029,xx08)	SIEMENS MEDCOM HEADER	1C	MedCom Header identification characteristics. Defined Terms: MEDCOM 1 Required if MedCom Header Info (0029,xx10) present.
MedCom Header Version	(0029,xx09)	SIEMENS MEDCOM HEADER	2C	Version of MedCom Header Info (0029,xx10) format. Required if MEDCOM Header Info (0029,xx10) present.
MedCom Header Info	(0029,xx10)	SIEMENS MEDCOM HEADER	3	Manufacturer model dependent information. The value of the attribute MedCom Header Info (0029,xx10) can be build up in each user defined format.
MedCom History Information	(0029,xx20)	SIEMENS MEDCOM HEADER	3	MedCom defined Patient Registration history information. See .
PMTF Information 1	(0029,xx31)	SIEMENS MEDCOM HEADER	3	Transformation Information
PMTF Information 2	(0029,xx32)	SIEMENS MEDCOM HEADER	3	Transformation Information
PMTF Information 3	(0029,xx33)	SIEMENS MEDCOM HEADER	3	Transformation Information
PMTF Information 4	(0029,xx34)	SIEMENS MEDCOM HEADER	3	Transformation Information

Table 75. MEDCOM Header Module

Attribute Name	Tag	Owner	Type	Notes
PMTF Information 5	(0029,xx35)	SIEMENS MEDCOM HEADER	3	Transformation Information
Application Header Sequence	(0029,xx40)	SIEMENS MEDCOM HEADER	3	Sequence of Application Header items. Zero or more items are possible.
>Application Header Type	(0029,xx41)	SIEMENS MEDCOM HEADER	1C	Application Header identification characteristics. Required, if Sequence is sent.
>Application Header ID	(0029,xx42)	SIEMENS MEDCOM HEADER	3	Identification of an application header
>Application Header Version	(0029,xx43)	SIEMENS MEDCOM HEADER	3	Version of CSA Series Header Info (0029,xx44) format.
>Application Header Info	(0029,xx44)	SIEMENS MEDCOM HEADER	3	Application dependent information.
Workflow Control Flags	(0029,xx50)	SIEMENS MEDCOM HEADER	3	Eight free definable flags.
Archive Management Flag Keep Online	(0029,xx51)	SIEMENS MEDCOM HEADER	3	Flag to control remote archive management system to keep the image always online (also when already archived). Enumerated Values: 00 = remote control not required 01 = keep image online
Archive Management Flag Do Not Archive	(0029,xx52)	SIEMENS MEDCOM HEADER	3	Flag to control remote archive management system not to archive the related image. Enumerated Values: 00 = remote control not required 01 = don't archive image

Table 75. MEDCOM Header Module

Attribute Name	Tag	Owner	Type	Notes
Image Location Status	(0029,xx53)	SIEMENS MEDCOM HEADER	3	Image location status to control retrieving. Defined Terms: ONLINE = retrieving has to be done as usual, NEARLINE = move request to SCP and delay according to value of Estimated Retrieve Time (0029,xx54), OFFLINE = invoking a retrieve operation initiates an operator request, INVALID = invoking a retrieve operation would always result in an error.
Estimated Retrieve Time	(0029,xx54)		3	Estimated retrieve time in seconds. A value less then zero (< 0) indicates location is OFFLINE or INVALID.
Data Size of Retrieved Images	(0029,xx55)		3	Data size of images in MByte.

MEDCOM History Information

The value of the attribute MEDCOM History Information (0029,xx20) is defined in the following way:

Table 76. MEDCOM History Information

Part	Name	Type	Bytes	Notes
header	Identifier	string	32	always "CSA HISTORY" e.g. "V1.10"
	Version	string	32	
n items	Class Name	string	64	
	Modification String	string	1024	

MEDCOM OOG Module

The table in this section contains private IOD Attributes that describe MEDCOM Object Oriented Graphics (OOG). This module is used whenever object graphics is drawn on the image and need to be stored as graphic object properties. Given the condition that the module contents was not removed by other modalities, the graphic objects remain re-animatable if such an image was transferred and is then retrieved back.

Table 77. MEDCOM OOG Module

Attribute Name	Tag	Owner	Type	Notes
MedCom OOG Type	(0029,xx08))	SIEMENS MEDCOM OOG	1	MEDCOM Object Oriented Graphics (OOG) identification characteristics. Defined Terms: MEDCOM OOG 1 MEDCOM OOG 2
MedCom OOG Version	(0029,xx09))	SIEMENS MEDCOM OOG	3	Version of MEDCOM OOG Info (0029,xx10) format.
MedCom OOG Info	(0029,xx10))	SIEMENS MEDCOM OOG	3	MEDCOM Object Oriented Graphics (OOG) data.

The graphics objects are also fully drawn in the Image Overlay Plane for compatibility with other products which do not support the MedCom OOG module. Any system not supporting the MedCom OOG module shall remove the OOG module and it's contents when modifying the image overlay plane content.

syngo **Report Data**

The module contains private IOD Attributes that describe *syngo* reports. This module is used when *syngo* report data are added to DICOM SR and DICOM SC objects.

Table 78. *syngo* Report Data Module

Attribute Name	Tag	Private Creator	Type	Notes
<i>syngo</i> Report Type	(0029,xx08))	SIEMENS CSA REPORT	1	<i>syngo</i> report characteristics, e.g. report creating application. Defined Terms: CT_LUNGCARE MR_ARGUS This attribute value will be used to identify the corresponding application during generic extension dll management. A restricted character set is used: only A-Z and underscore are supported.
<i>syngo</i> Report Version	(0029,xx09))	SIEMENS CSA REPORT	3	Version of <i>syngo</i> Report Data (0029,xx10) format.
<i>syngo</i> Report Data	(0029,xx10))	SIEMENS CSA ENVELOPE	3	A representation of DICOM SR attributeContent Sequence (0040,A730). This includes the document relationship and document content. This data will typically be represented using an XML encoding according to a Siemens private scheme.
<i>syngo</i> Report Presentation	(0029,xx11))	SIEMENS CSA ENVELOPE	3	A representation of the recommended presentation for the <i>syngo</i> Report Data (0029,xx10). This presentation will typically be encoded in XSLT.

Table 78. syngo Report Data Module

Attribute Name	Tag	Private Creator	Type	Notes
SR Variant	(0029,xx15)	SIEMENS CSA REPORT	3	DICOM SR variant. Enumerated Values: 0 = Basic Text SR (1.2.840.10008.5.1.4.1.1.88.11 1 = Enahanced SR (1.2.840.10008.5.1.4.1.1.88.22 2 = Comprehensive SR (1.2.840.10008.5.1.4.1.1.88.33 3 = Mammography Cad SR (1.2.840.10008.5.1.4.1.1.88.50 4 = Key Object Selection Document (1.2.840.10008.5.1.4.1.1.88.59
SC SOP Instance UID	(0029,xx17)	SIEMENS CSA REPORT	3	DICOM SOP Instance UID of <i>syngo</i> based SC Image representing the <i>syngo</i> report object. This UID will be used to identify the resulting SC object after SR to SC conversion.

syngo **Report Info**

The module *syngo* Report Info contains all DICOM SR attributes expect the Contents Sequence (0040,A730). This module is only used during SR to SC conversion.

Registry of DICOM Data Elements

<modify this table and add your own attribute, if needed>

Tag	Private Owner Code	Name	VR	VM
(0029,xx08)	SIEMENS CSA NON-IMAGE	CSA Data Type	CS	1
(0029,xx09)	SIEMENS CSA NON-IMAGE	CSA Data Version	LO	1
(0029,xx10)	SIEMENS CSA NON-IMAGE	CSA Data Info	OB	1
(0029,xx08)	SIEMENS CSA HEADER	CSA Image Header Type	CS	1

(0029,xx09)	SIEMENS CSA HEADER	CSA Image Header Version	LO	1
(0029,xx10)	SIEMENS CSA HEADER	CSA Image Header Info	OB	1
(0029,xx18)	SIEMENS CSA HEADER	CSA Series Header Type	CS	1
(0029,xx19)	SIEMENS CSA HEADER	CSA Series Header Version	LO	1
(0029,xx20)	SIEMENS CSA HEADER	CSA Series Header Info	OB	1
(0029,xx08)	SIEMENS CSA REPORT	<i>syngo</i> Report Type	CS	1
(0029,xx09)	SIEMENS CSA REPORT	<i>syngo</i> Report Version	LO	1
(0029,xx15)	SIEMENS CSA REPORT	SR Variant	US	1
(0029,xx17)	SIEMENS CSA REPORT	SC SOP Instance UID	UI	1
(0029,xx10)	SIEMENS CSA ENVELOPE	<i>syngo</i> Report Data	OB	1
(0029,xx11)	SIEMENS CSA ENVELOPE	<i>syngo</i> Report Presentation	OB	1
(0029,xx08)	SIEMENS MEDCOM HEADER	MedCom Header Type	CS	1
(0029,xx09)	SIEMENS MEDCOM HEADER	MedCom Header Version	LO	1
(0029,xx10)	SIEMENS MEDCOM HEADER	MedCom Header Info	OB	1
(0029,xx20)	SIEMENS MEDCOM HEADER	MedCom History Information	OB	1
(0029,xx31)	SIEMENS MEDCOM HEADER	PMTF Information 1	LO	1
(0029,xx32)	SIEMENS MEDCOM HEADER	PMTF Information 2	UL	1
(0029,xx33)	SIEMENS MEDCOM HEADER	PMTF Information 3	UL	1
(0029,xx34)	SIEMENS MEDCOM HEADER	PMTF Information 4	CS	1
(0029,xx35)	SIEMENS MEDCOM HEADER	PMTF Information 5	UL	1
(0029,xx40)	SIEMENS MEDCOM HEADER	Application Header Sequence	SQ	1
(0029,xx41)	SIEMENS MEDCOM HEADER	Application Header Type	CS	1
(0029,xx42)	SIEMENS MEDCOM HEADER	Application Header ID	LO	1
(0029,xx43)	SIEMENS MEDCOM HEADER	Application Header Version	LO	1
(0029,xx44)	SIEMENS MEDCOM HEADER	Application Header Info	OB	1

(0029,xx50)	SIEMENS MEDCOM HEADER	Workflow Control Flags	LO	8
(0029,xx51)	SIEMENS MEDCOM HEADER	Arch. Mgmnt Flag Keep Online	CS	1
(0029,xx52)	SIEMENS MEDCOM HEADER	Arch. Mgmnt Flag Do not Archive	CS	1
(0029,xx53)	SIEMENS MEDCOM HEADER	Image Location Status	CS	1
(0029,xx54)	SIEMENS MEDCOM HEADER	Estimated Retrieve Time	DS	1
(0029,xx55)	SIEMENS MEDCOM HEADER	Data Size of Retrieved Images	DS	1
(0029,xx08)	SIEMENS MEDCOM OOG	MEDCOM OOG Type	CS	1
(0029,xx09)	SIEMENS MEDCOM OOG	MEDCOM OOG Version	LO	1
(0029,xx10)	SIEMENS MEDCOM OOG	MEDCOM OOG Info	OB	1
(7FE1,xx10)	SIEMENS CSA NON-IMAGE	CSA Data	OB	1

Standard Extensions of all SOP Classes

The following tables list the data dictionary of all DICOM IOD attributes where the DICOM standard definitions are extended:

Note: Additional standard and private tags in SQ which are not defined by the IOD module are not supported. These tags are ignored.

Table 79. Standard Extensions of all SOP Classes

Attribute Name	Tag	Private Creator	Type	Notes
Image Type	(0008,0008 -)		1	see additional Defined Terms: Defined Terms for value 3: OTHER Defined Terms for value 4: CSA 3D EDITOR CSA 3D FLY PATH CSA 3D FLY VRT CSA 3D FUSION CSA AVERAGE CSA BLACK IMAGE CSA RESAMPLED CSA MIP CSA MPR CSA MPR CURVED CSA MPR THICK CSA SSD CSA SUBTRACT CT_SOM4 * ECAT ACF ECAT NORMAL ECAT 3D SINO ECAT 3D SINO FLT SHS *
Patient Position	(0018,5100 -)		2C	see Some of the additional terms are: HLS HLP FLS FLP HLDL HLDR FLDL FLDR

All SOP classes may contain additional type 3 attributes which DICOM standard defines in a different DICOM IOD or DICOM SOP class (attributes from Normalized SOP classes).

This is the case for example for

- Rescale Slope (0028,1053)
- Rescale Intercept (0028,1052)

which are also used in the MR IOD.

Image Type

The Image Type (0008,0008) attribute identifies important image identification characteristics. These characteristics are:

- 1 Pixel Data Characteristics:
 - is the image an ORIGINAL Image; an image whose pixel values are based on original or source data, or
 - is the image a DERIVED Image; an image whose pixel values have been derived in some manner from the pixel value of one or more other images.
- 2 Patient Examination Characteristics:
 - is the image a PRIMARY Image; an image created as a direct result of the Patient examination, or
 - is the image a SECONDARY Image; an image created after the initial Patient examination.
- 3 Modality Specific Characteristics (SOP Specific Characteristics).
- 4 Implementation specific identifiers; other implementation specific identifiers shall be documented in an implementation's conformance claim.

The Image Type attribute is multi-valued and shall be provided in the following manner:

- Value 1 shall identify the Pixel Data Characteristics; Enumerated Values for the Pixel Data Characteristics are:
 - ORIGINAL = identifies an Original Image
 - DERIVED = identifies a Derived Image
- Value 2 shall identify the Patient Examination Characteristics; Enumerated Values for the Patient Examination Characteristics are:
 - PRIMARY = identifies a Primary Image
 - SECONDARY = identifies a Secondary Image
- Value 3 shall identify any Image IOD specific specialization, the following terms are defined in addition to the DICOM standard definitions:
 - OTHER = is also used for converted non-Axial and non-Localizer CT images
 - MPR = for 3D MPR images
 - PROJECTION IMAGE = for 3D MIP and SSD images
- Value 4 which are implementation specific, the following terms are defined in addition to the DICOM standard definitions:

- original *syngo* generated data set types:

CSA 3D EDITOR = object created by 3D Editor

CSA 3D FLY PATH = object created by Fly Through Path

CSA 3D FLY VRT = object created by Fly Through Volume Rendering Technique

CSA 3D FUSION = object created by Fusion

CSA AVERAGE = image was created by Average

CSA BLACK IMAGE = SC Image with black pixels, only graphics information is of interest

CSA MIP = image created by Maximum Intensity Projection

CSA MIP THIN = image created by Maximum Intensity Projection

CSA MPR = image created by Multi Planar Reconstruction

CSA MPR CURVED = image created by Multi Planar Reconstruction

CSA MPR THICK = image created by Multi Planar Reconstruction

CSA MPR THIN = image created by Multi Planar Reconstruction

CSA RESAMPLED = derived image created by zooming or panning original image

CSA REPORT = *syngo* reporting (documentation of diagnosis)

CSA RESULT = *syngo* reporting (postprocessing results)

CSA SSD = SC Image as Shaded Surface Display

CSA SUBTRACT = image was created by Subtraction

ECAT ACF = CTI PET Attenuation Correction

ECAT NORMAL = CTI PET Normalization

ECAT 3D SINO = CTI PET 3D Sinogram Short

ECAT 3D SINO FLT = CTI PET 3D Sinogram Float

- converted images

CT_SOM4 NONE = converted SOMARIS image

CT_SOM4 CONV = converted SOMARIS Convolution Kernel file

CT_SOM4 DART = converted SOMARIS Dental Artificial image

CT_SOM4 DEVA = converted SOMARIS Dental Evaluation image

CT_SOM4 DGRA = converted SOMARIS Dental Graphics image

CT_SOM4 DMEA = converted SOMARIS Dynamic Measurement image

CT_SOM4 DPAN = converted SOMARIS Dental Panorama image

CT_SOM4 DPAR = converted SOMARIS Dental Paraxial image

CT_SOM4 EBT = converted SOMARIS Evolution image
CT_SOM4 HIS = converted SOMARIS Histogram Graphics image
CT_SOM4 HISC = converted SOMARIS Histogram Graphics image
CT_SOM4 MUL = converted SOMARIS Multiscan image
CT_SOM4 OEVA = converted SOMARIS Osteo Evaluation image
CT_SOM4 OTOM = converted SOMARIS Osteo Tomogram image
CT_SOM4 OTOP = converted SOMARIS Osteo Topogram image
CT_SOM4 PLOT = converted SOMARIS Plot image
CT_SOM4 QUAL = converted SOMARIS Quality image
CT_SOM4 R2D = converted SOMARIS 2D Rebuild image
CT_SOM4 R3D = converted SOMARIS 3D Rebuild image
CT_SOM4 R3DE = converted SOMARIS 3D Rebuild image
CT_SOM4 RMAX = converted SOMARIS Maximum Intensity Projection image
CT_SOM4 RMIN = converted SOMARIS Minimum Intensity Projection image
CT_SOM4 ROT = converted SOMARIS Rotation Mode image
CT_SOM4 RRAD = converted SOMARIS Radiographic Projection image
CT_SOM4 RVIT = converted SOMARIS Vessel Image Tool image
CT_SOM4 RVRT = converted SOMARIS Volumetric Rendering image
CT_SOM4 SAVE = converted SOMARIS Evolution Screen Save image
CT_SOM4 SCAN = converted SOMARIS Standard Mode image
CT_SOM4 SEQ = converted SOMARIS Sequence Mode image
CT_SOM4 SER = converted SOMARIS Serial Mode image
CT_SOM4 SIN = converted SOMARIS Sinogram image
CT_SOM4 SINC = converted SOMARIS Sinogram image
CT_SOM4 SPI = converted SOMARIS Spiral Mode image
CT_SOM4 STA = converted SOMARIS Static Mode image
CT_SOM4 TAB = converted SOMARIS Correction Table image
CT_SOM4 TOP = converted SOMARIS Topogram image
CT_SOM4 GTOP = converted SOMARIS Topo Graphics image
CT_SOM4 PEVG = converted SOMARIS Pulmo Evaluation image
CT_SOM4 PEVI = converted SOMARIS Pulmo Evaluation image

CT_SOM4 PUL = converted SOMARIS Pulmo Respiration image
CT_SOM4 PROT = converted SOMARIS Protocol image
CT_SOM4 TEXT = converted SOMARIS Text image
CT_SOM4 ICD = converted SOMARIS Interventional Cine image
SHS DENT = converted MagicView Dental Tomogram image
SHS DPAN = converted MagicView Dental Panorama image
SHS DPAR = converted MagicView Dental Paraxial image
SHS 3D_CURVED = converted MagicView image
SHS 3D_MIP = converted MagicView Maximum Intensity Projection image
SHS 3D_MPR = converted MagicView Multi Planar Reconstruction image
SHS 3D_SSD = converted MagicView Shaded Surface Display image
SHS 3D_VRT = converted MagicView Volumetric Rendering image

Patient Position

The Patient Position attribute (0018,5100) defines the patient position relative to the equipment.

The Defined Terms for this value were extended for the Siemens products. Here the patient is not positioned HeadFirst/FeetFirst when facing the front of the imaging equipment but HeadLeft or FeetLeft.

the new values are:

- HLS (Head left - Supine)
- HLP (Head left - Prone)
- FLS (Feet left - Supine)
- FLP (Feet left - Prone)
- HLDL (Head left - Decubitus left)
- HLDR (Head left - Decubitus right)
- FLDL (Feet left - Decubitus left)
- FLDR (Feet left - Decubitus right)

DICOM Print SCU - detailed status displays

The following tables document the behavior of the VSim DICOM Print AE in response to messages received for the printer SOP class and the print job SOP class.

Common Status Information

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
NORMAL	Camera is ready.	Camera is ready.	<None>/idle
BAD RECEIVE MGZ	There is a problem with the film receive magazine. Films from the printer cannot be transported into the magazine.	Problem with receive magazine.	<None>/interact
BAD SUPPLY MGZ	There is a problem with a film supply magazine. Films from this magazine cannot be transported into the printer.	Problem with supply magazine.	<None>/interact
CALIBRATING	Printer is performing self calibration, it is expected to be available for normal operation shortly.	Self calibration. Please wait.	<None>/idle
CALIBRATION ERROR	An error in the printer calibration has been detected, quality of processed films may not be optimal.	Problem in calibration. Film quality may not be optimal.	<None>/interact
CHECK CHEMISTRY	A problem with the processor chemicals has been detected, quality of processed films may not be optimal.	Problem with chemistry. Film quality may not be optimal.	<None>/interact
CHECK SORTER	There is an error in the film sorter	Error in film sorter.	<None>/interact
CHEMICALS EMPTY	There are no processing chemicals in the processor, films will not be printed and processed until the processor is back to normal.	Camera chemistry empty. Please check.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngo</i>/camera symbol
CHEMICALS LOW	The chemical level in the processor is low, if not corrected, it will probably shut down soon.	Camera chemistry low. Please check.	<None>/interact
COVER OPEN	One or more printer or processor covers, drawers, doors are open.	Camera cover, drawer or door open.	<None>/interact
ELEC CONFIG ERR	Printer configured improperly for this job.	Camera configured improperly for this job. Queue stopped.	Queue for this camera will be STOPPED/Queue stopped
ELEC DOWN	Printer is not operating due to some unspecified electrical hardware problem.	Camera electrical hardware problem.	<None>/interact
ELEC SW ERROR	Printer not operating for some unspecified software error.	Camera software problem. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
EMPTY 8x10	The 8x10 inch film supply magazine is empty.	8x10 film supply empty.	<None>/interact
EMPTY 8x10 BLUE	The 8x10 inch blue film supply magazine is empty.	8x10 blue film supply empty.	<None>/interact
EMPTY 8x10 CLR	The 8x10 inch clear film supply magazine is empty.	8x10 clear film supply empty.	<None>/interact
EMPTY 8x10 PAPER	The 8x10 inch paper supply magazine is empty.	8x10 paper supply empty.	<None>/interact
EMPTY 10x12	The 10x12 inch film supply magazine is empty.	10x12 film supply empty.	<None>/interact
EMPTY 10x12 BLUE	The 10x12 inch blue film supply magazine is empty.	10x12 blue film supply empty.	<None>/interact
EMPTY 10x12 CLR	The 10x12 inch clear film supply magazine is empty.	10x12 clear film supply empty.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
EMPTY 10x12 PAPR	The 10x12 inch paper supply magazine is empty.	10x12 paper supply empty.	<None>/interact
EMPTY 10x14	The 10x14 inch film supply magazine is empty.	10x14 film supply empty.	<None>/interact
EMPTY 10x14 BLUE	The 10x14 inch blue film supply magazine is empty.	10x14 blue film supply empty.	<None>/interact
EMPTY 10x14 CLR	The 10x14 inch clear film supply magazine is empty.	10x14 clear film supply empty.	<None>/interact
EMPTY 10x14 PAPR	The 10x14 inch paper supply magazine is empty.	10x14 paper supply empty.	<None>/interact
EMPTY 11x14	The 11x14 inch film supply magazine is empty.	11x14 film supply empty.	<None>/interact
EMPTY 11x14 BLUE	The 11x14 inch blue film supply magazine is empty.	11x14 blue film supply empty.	<None>/interact
EMPTY 11x14 CLR	The 11x14 inch clear film supply magazine is empty.	11x14 clear film supply empty.	<None>/interact
EMPTY 11x14 PAPR	The 11x14 inch paper supply magazine is empty.	11x14 paper supply empty.	<None>/interact
EMPTY 14x14	The 14x14 inch film supply magazine is empty.	14x14 film supply empty.	<None>/interact
EMPTY 14x14 BLUE	The 14x14 inch blue film supply magazine is empty.	14x14 blue film supply empty.	<None>/interact
EMPTY 14x14 CLR	The 14x14 inch clear film supply magazine is empty.	14x14 clear film supply empty.	<None>/interact
EMPTY 14x14 PAPR	The 14x14 inch paper supply magazine is empty.	14x14 paper supply empty.	<None>/interact
EMPTY 14x17	The 14x17 inch film supply magazine is empty.	14x17 film supply empty.	<None>/interact
EMPTY 14x17 BLUE	The 14x17 inch blue film supply magazine is empty.	14x17 blue film supply empty.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
EMPTY 14x17 CLR	The 14x17 inch clear film supply magazine is empty.	14x17 clear film supply empty.	<None>/interact
EMPTY 14x17 PAPR	The 14x17 inch paper supply magazine is empty.	14x17 paper supply empty.	<None>/interact
EMPTY 24x24	The 24x24 inch film supply magazine is empty.	24x24 film supply empty.	<None>/interact
EMPTY 24x24 BLUE	The 24x24 inch blue film supply magazine is empty.	24x24 blue film supply empty.	<None>/interact
EMPTY 24x24 CLR	The 24x24 inch clear film supply magazine is empty.	24x24 clear film supply empty.	<None>/interact
EMPTY 24x24 PAPR	The 24x24 inch paper supply magazine is empty.	24x24 paper supply empty.	<None>/interact
EMPTY 24x30	The 24x30 inch film supply magazine is empty.	24x30 film supply empty.	<None>/interact
EMPTY 24x30 BLUE	The 24x30 inch blue film supply magazine is empty.	24x30 blue film supply empty.	<None>/interact
EMPTY 24x30 CLR	The 24x30 inch clear film supply magazine is empty.	24x30 clear film supply empty.	<None>/interact
EMPTY 24x30 PAPR	The 24x30 inch paper supply magazine is empty.	24x30 paper supply empty.	<None>/interact
EMPTY A4 PAPR	The A4 paper supply magazine is empty.	A4 paper supply empty.	<None>/interact
EMPTY A4 TRANS	The A4 transparency supply magazine is empty.	A4 transparency supply empty.	<None>/interact
EXPOSURE FAILURE	The exposure device has failed due to some unspecified reason.	Exposure device has failed.	<None>/interact
FILM JAM	A film transport error has occurred and a film is jammed in the printer or processor.	Film jam.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
FILM TRANSP ERR	There is a malfunction with the film transport, there may or may not be a film jam.	Film transport problem.	<None>/interact
FINISHER EMPTY	The finisher is empty.	Finisher is empty.	<None>/interact
FINISHER ERROR	The finisher is not operating due to some unspecified reason.	Finisher problem.	<None>/interact
FINISHER LOW	The finisher is low on supplies	Finisher low.	<None>/interact
LOW 8x10	The 8x10 inch film supply magazine is low.	8x10 film supply low.	<None>/interact
LOW 8x10 BLUE	The 8x10 inch blue film supply magazine is low.	8x10 blue film supply low.	<None>/interact
LOW 8x10 CLR	The 8x10 inch clear film supply magazine is low.	8x10 clear film supply low.	<None>/interact
LOW 8x10 PAPR	The 8x10 inch paper supply magazine is low.	8x10 paper supply low.	<None>/interact
LOW 10x12	The 10x12 inch film supply magazine is low.	10x12 film supply low.	<None>/interact
LOW 10x12 BLUE	The 10x12 inch blue film supply magazine is low.	10x12 blue film supply low.	<None>/interact
LOW 10x12 CLR	The 10x12 inch clear film supply magazine is low.	10x12 clear film supply low.	<None>/interact
LOW 10x12 PAPR	The 10x12 inch paper supply magazine is low.	10x12 paper supply low.	<None>/interact
LOW 10x14	The 10x14 inch film supply magazine is low.	10x14 film supply low.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
LOW 10x14 BLUE	The 10x14 inch blue film supply magazine is low.	10x14 blue film supply low.	<None>/interact
LOW 10x14 CLR	The 10x14 inch clear film supply magazine is low.	10x14 clear film supply low.	<None>/interact
LOW 10x14 PAPR	The 10x14 inch paper supply magazine is low.	10x14 paper supply low.	<None>/interact
LOW 11x14	The 11x14 inch film supply magazine is low.	11x14 film supply low.	<None>/interact
LOW 11x14 BLUE	The 11x14 inch blue film supply magazine is low.	11x14 blue film supply low.	<None>/interact
LOW 11x14 CLR	The 11x14 inch clear film supply magazine is low.	11x14 clear film supply low.	<None>/interact
LOW 11x14 PAPR	The 11x14 inch paper supply magazine is low.	11x14 paper supply low.	<None>/interact
LOW 14x14	The 14x14 inch film supply magazine is low.	14x14 film supply low.	<None>/interact
LOW 14x14 BLUE	The 14x14 inch blue film supply magazine is low.	14x14 blue film supply low.	<None>/interact
LOW 14x14 CLR	The 14x14 inch clear film supply magazine is low.	14x14 clear film supply low.	<None>/interact
LOW 14x14 PAPR	The 14x14 inch paper supply magazine is low.	14x14 paper supply low.	<None>/interact
LOW 14x17	The 14x17 inch film supply magazine is low.	14x17 film supply low.	<None>/interact
LOW 14x17 BLUE	The 14x17 inch blue film supply magazine is low.	14x17 blue film supply low.	<None>/interact
LOW 14x17 CLR	The 14x17 inch clear film supply magazine is low.	14x17 clear film supply low.	<None>/interact
LOW 14x17 PAPR	The 14x17 inch paper supply magazine is low.	14x17 paper supply low.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
LOW 24x24	The 24x24 inch film supply magazine is low.	24x24 film supply low.	<None>/interact
LOW 24x24 BLUE	The 24x24 inch blue film supply magazine is low.	24x24 blue film supply low.	<None>/interact
LOW 24x24 CLR	The 24x24 inch clear film supply magazine is low.	24x24 clear film supply low.	<None>/interact
LOW 24x24 PAPR	The 24x24 inch paper supply magazine is low.	24x24 paper supply low.	<None>/interact
LOW 24x30	The 24x30 inch film supply magazine is low.	24x30 film supply low.	<None>/interact
LOW 24x30 BLUE	The 24x30 inch blue film supply magazine is low.	24x30 blue film supply low.	<None>/interact
LOW 24x30 CLR	The 24x30 inch clear film supply magazine is low.	24x30 clear film supply low.	<None>/interact
LOW 24x30 PAPR	The 24x30 inch paper supply magazine is low.	24x30 paper supply low.	<None>/interact
LOW A4 PAPR	The A4 paper supply magazine is low.	A4 paper supply low.	<None>/interact
LOW A4 TRANS	The A4 transparency supply magazine is low.	A4 transparency supply low.	<None>/interact
NO RECEIVE MGZ	The film receive magazine no available.	Film receiver not available.	<None>/interact
NO RIBBON	The ribbon cartridge needs to be replaced.	Replace ribbon cartridge.	<None>/interact
NO SUPPLY MGZ	The film supply magazine specified for this job is not available.	Film supply not available.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
CHECK PRINTER	The printer is not ready at this time, operator intervention is required to make the printer available.	Check camera.	<None>/interact
CHECK PROC	The processor is not ready at this time, operator intervention is required to make the printer available.	Check processor.	<None>/interact
PRINTER DOWN	The printer is not operating due to some unspecified reason.	Camera down.	<None>/interact
PRINTER INIT	The printer is not ready at this time, it is expected to become available without intervention, For example, it may be in a normal warm-up state.	Camera initializing.	<None>/idle
PRINTER OFFLINE	The printer has been disabled by an operator or service person.	Camera off-line.	<None>/interact
PROC DOWN	The processor is not operating due to some unspecified reason.	Processor down.	<None>/interact
PROC INIT	The processor is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state.	Processor initializing.	<None>/idle
PROC OVERFLOW FL	Processor chemicals are approaching the overflow full mark.	Processor chemicals overflow.	<None>/interact
PROC OVERFLOW HI	Processor chemicals have reached the overflow full mark.	Processor chemicals near overflow.	<None>/interact

Table 80. Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for <i>syngo</i> /camera symbol
QUEUED	Print job in Queue	-	<None>/idle
RECEIVER FULL	The Film receive magazine is full.	Receiver full.	<None>/interact
REQ MED NOT INST	The requested film, paper, or other media supply magazine is installed in the printer, but may be available with operator intervention.	Install media supply.	<None>/interact
REQ MED NOT AVAI	The requested film, paper, or other media requested is not available on this printer.	Media supply not available on this camera. Queue stopped. Change camera.	Queue for this camera will be STOPPED/queue stopped
RIBBON ERROR	There is an unspecified problem with the print ribbon.	Error with print ribbon.	<None>/interact
SUPPLY EMPTY	The printer is out of film.	Camera out of film.	<None>/interact
SUPPLY LOW	The film supply is low.	Film supply low.	<None>/interact
UNKNOWN	There is an unspecified problem.	Unspecified problem with camera.	<None>/interact

Additional Status Information - AGFA Printers

Table 81. Printer Status Infos: Additional Agfa printer status infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
WARMING UP	Printer is in the warm-up stage. Spooling of print jobs to disk is still possible.	Camera is warming up.	<None>idle
OFFLINE	Printer is switched off-line. Spooling of print jobs to disk is still possible.	Camera is switched off-line.	<None>/interact
NONE	General printer warning, no specific information is available. Spooling of print jobs to disk is still possible.	-	<None>/idle
-	-	-	-

Additional Status Information - Kodak PACS Link (formerly Imation)

Table 82. Printer Status Infos: Additional Kodak infos for Pacs Link (formerly Imation cameras)

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
SUPPLY MGZ ERR	The supply magazine has an error.	Film supply has an error.	<None>/interact
-	-	-	-

Additional Status Information - Kodak 190I

Table 83. Printer Status Infos: Additional Kodak infos for Kodak 190

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngol</i> camera symbol
PRINTER STOPPED	The printer has stopped	Camera has stopped.	<None>/interact
FATAL ERROR	Fatal error.	Fatal error. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
-	-	-	-

Additional Status Information - Kodak 2180/1120

Table 84. Printer Status Infos: Additional Kodak infos for 2180/1120

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngo</i> / camera symbol
PRINTER NOT RDY	Printer not ready.	Camera not ready.	<None>/interact
CHECK PROCESSOR	Check processor.	Check processor.	<None>/interact
NO TONER	No toner.	No toner.	<None>/interact
FATAL	Fatal error.	Fatal error. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
-	-	-	-

Additional Status Information - Codonics

Table 85. Printer Status Infos: Additional Codonics infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngo</i> / camera symbol
STANDARD	Printer is ready.	Camera is ready.	<None>/Normal
LOAD A-SIZE	Load A-size media.	Load A-size media.	<None>/interact
LOAD A-DVPAPER	Load A-size black and white paper.	Load A-size black and white paper.	<None>/interact
LOAD A-CVPAPER	Load A-size color paper.	Load A-size color paper.	<None>/interact
LOAD A-CVTRANS	Load A-size transparencies.	Load A-size transparencies.	<None>/interact
LOAD A4-SIZE	Load A4-size media.	Load A4-size media.	<None>/interact
LOAD A4-DVPAPER	Load A4-size black and white paper.	Load A4-size black and white paper.	<None>/interact
LOAD A4-CVPAPER	Load A4-size color paper.	Load A4-size color paper.	<None>/interact

Table 85. Printer Status Infos: Additional Codonics infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
LOAD A4-CVTRANS	Load A4-size transparencies.	Load A4-size transparencies.	<None>/interact
LOAD LA-SIZE	Load LA-size media.	Load LA-size media.	<None>/interact
LOAD LA-DVPAPER	Load LA-size black and white paper.	Load LA-size black and white paper.	<None>/interact
LOAD LA-CVPAPER	Load LA-size color paper.	Load LA-size color paper.	<None>/interact
LOAD LA-CVTRANS	Load LA-size transparencies.	Load LA-size transparencies.	<None>/interact
LOAD LA4-SIZE	Load LA4-size media.	Load LA4-size media.	<None>/interact
LOAD LA4-DVPAPER	Load LA4-size black and white paper.	Load LA4-size black and white paper.	<None>/interact
LOAD LA4-CVPAPER	Load LA4-size color paper.	Load LA4-size color paper.	<None>/interact
LOAD LA4-CVTRANS	Load LA4-size transparencies.	Load LA4-size transparencies.	<None>/interact
LOAD XLA-SIZE	Load XLA-size media.	Load XLA-size media.	<None>/interact
LOAD XLA-DVPAPER	Load XLA-size black and white paper.	Load XLA-size black and white paper.	<None>/interact
LOAD XLA-CVPAPER	Load XLA-size color paper.	Load XLA-size color paper.	<None>/interact
LOAD XLA-CVTRANS	Load XLA-size transparencies.	Load XLA-size transparencies.	<None>/interact
LOAD XLA4-SIZE	Load XLA4-size media.	Load XLA4-size media.	<None>/interact
LOAD XLA4-DVPAPE	Load XLA4-size black and white paper.	Load XLA4-size black and white paper.	<None>/interact

Table 85. Printer Status Infos: Additional Codonics infos


Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
LOAD XLA4-CVPAPE	Load XLA4-size color paper.	Load XLA4-size color paper.	<None>/interact
LOAD XLA4-CVTRAN	Load XLA4-size transparencies.	Load XLA4-size transparencies.	<None>/interact
LOAD XLW-SIZE	Load XLW-size media.	Load XLW-size media.	<None>/interact
LOAD XLW-DVPAPER	Load XLW-size black and white paper.	Load XLW-size black and white paper.	<None>/interact
LOAD XLW-CVPAPER	Load XLW-size color paper.	Load XLW-size color paper.	<None>/interact
LOAD 8X10-SIZE	Load 8x10 media.	Load 8x10 media.	<None>/interact
LOAD 8X10-DVFILM	Load 8x10 black and white film.	Load 8x10 black and white film.	<None>/interact
SUPPLY MISSING	The film supply magazine specified for this job is not available.	Film supply not available.	<None>/interact
RIBBON MISSING	Ribbon is missing.	Ribbon is missing.	<None>/interact
RIBBON EMPTY	Ribbon is empty.	Ribbon is empty.	<None>/interact
TOP COVER OPEN	Top cover of printer is open.	Top cover of camera is open.	<None>/interact
-	-	-	-

Additional DICOM Execution Status Information

Table 86. Additional DICOM Execution Status Infos

Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
INVALID PAGE DES	The specified page layout cannot be printed or other page description errors have been detected.	Film Job cannot be printed on this camera. Queue stopped. Please redirect film job.	Queue for this camera will be STOPPED/queue stopped
INSUFFIC MEMORY	There is not enough memory available to complete this job.	Not enough memory available in camera. Queue stopped. Please continue queue or change camera.	Queue for this camera will be STOPPED/queue stopped
NONE	General printer warning, no specific information is available. Spooling of print jobs to disk is still possible.	-	<None>/idle

RT IOD specific implementation details



DICOM INFORMATION MODEL - RADIO THERAPY

Table 87. RT Series Table

Attribute Name	Tag	Type	Notes
Modality	(0008,0060)	1	RTIMAGE if RT Image IOD, RTSTRUCT if RT Structure Set IOD, RTPLAN if RT Plan IOD. RTDOSE & RTRECORD are not supported.
Series Instance UID	(0020,000E)	1	
Series Number	(0020,0011)	2	
Series Description	(0008,103E)	3	

Table 88. RT Plan IOD module Table

IE	Module	Reference	DICOM Usage	Notes
Patient	Patient	C.7.1.1	M	
Study	General Study	C.7.2.1	M	
	Patient Study	C.7.2.2	U	Not supported
Series	RT Series	C.8.8.1	M	
Equipment	General Equipment	C.7.5.1	M	
Plan	RT General Plan	C.8.8.9	M	
	RT Prescription	C.8.8.10	U	Not supported
	RT Tolerance Tables	C.8.8.11	U	Not supported
	RT Patient Setup	C.8.8.12	U	
	RT Fraction Scheme	C.8.8.13	U	Not supported
	RT Beams	C.8.8.14	C	
	RT Brachy Application Setups	C.8.8.15	C	Not supported
	Approval	C.8.8.16	U	Not supported
	Audio	C.10.3	U	Not supported
	SOP Common	C.12.1	M	

Table 89. RT Structureset IOD Module Table

IE	Module	Reference	DICOM Usage	Notes
Patient	Patient	C.7.1.1	M	
Study	General Study	C.7.2.1	M	
	Patient Study	C.7.2.2	U	Not supported
Series	RT Series	C.8.8.1	M	
Equipment	General Equipment	C.7.5.1	M	
Structure Set	Structure Set	C.8.8.5	M	
	ROI Contour	C.8.8.6	M	
	RT ROI Observations	C.8.8.8	M	
	Approval	C.8.8.16	U	Not supported
	SOP Common	C.12.1	M	
	Audio	C.10.3	U	Not supported

Table 90. RT Image IOD Module Table

IE	Module	Reference	DICOM Usage	Notes
Patient	Patient	C.7.1.1	M	
Study	General	C.7.2.1	M	
	Patient Study	C.7.2.2	U	Not supported
Series	RT Series	C.8.8.1	M	
Frame Of Reference	Frame Of Reference	C.7.4.1	U	
Equipment	General Equipment	C.7.5.1	M	
Image	General Image	C.7.6.1	M	
	Image Pixel	C.7.6.3	M	
	Contrast/Bolus	C.7.6.4	C	Not supported
	Cine	C.7.6.5	C	Not supported
	Multi Frame	C.7.6.6	C	Not supported
	RT Image	C.8.8.2	M	
	Modality LUT	C.11.1	U	Not supported
	VOI LUT	C.11.2	U	Not supported
	Approval	C.8.8.16	U	Not supported
	Curve	C.10.2	U	Not supported
	Audio	C.10.3	U	Not supported
	SOP Common	C.12.1	M	

Table 91. RT Structure Set Module

>>Referenced Instance UID	SOP	(0008,1155)	1C	
>>RT Referenced Series Sequence		(3006,0014)	1C	
>>>Series Instance UID		(0020,000E)	1C	
>>>Contour Image Sequence		(3006,0016)	1C	
>>>>Referenced Class UID	SOP	(0008,1150)	1C	
>>>>Referenced Instance UID	SOP	(0008,1155)	1C	
Structure Set ROI Sequence		(3006,0020)	3	
>ROI Number		(3006,0022)	1C	
>Referenced Frame of Reference UID		(3006,0024)	1C	
>ROI Name		(3006,0026)	2C	VSIM Structure Name for the ROI
>ROI Description		(3006,0028)	3	Not supported
>ROI Volume		(3006,002C)	3	Not supported
>ROI Algorithm	Generation	(3006,0036)	2C	
>ROI Description	Generation	(3006,0038)	3	Not supported

Table 92. ROI Contour Module

Attribute Name	Tag	Type	Note
ROI Contour Sequence	(3006,0039)	1	
>Referenced ROI Number	(3006,0084)	1	
>ROI Display Color	(3006,002A)	3	Display Color as three (RGB) values
>Contour Sequence	(3006,0040)	3	
>>Contour Number	(3006,0048)	3	
>>Attached Contours	(3006,0049)	3	Not supported
>>Contour Image Sequence	(3006,0016)	3	
>>>Referenced SOP Class UID	(0008,1150)	1C	
>>>Referenced SOP Instance UID	(0008,1155)	1C	
>>>Referenced Frame Number	(0008,1160)	1C	
>>Contour Geometric Type	(3006,0042)	1C	
>>Contour Slab Thickness	(3006,0044)	3	Not supported
>>Contour Offset Vector	(3006,0045)	3	Not supported
>>Number of Contour Points	(3006,0046)	1C	
>>Contour Data	(3006,0050)	1C	

Table 93. ROI Observations Module

Attribute Name	Tag	Type	Note
RT ROI Observations Sequence	(3006,0080)	1	
>Observation Number	(3006,0082)	1	
>Referenced ROI Number	(3006,0084)	1	
>ROI Observation Label	(3006,0085)	3	VSim Structure Name for the associated ROI
>ROI Observation Description	(3006,0088)	3	Not supported
>RT Related ROI Sequence	(3006,0030)	3	Not supported
>>Referenced ROI Number	(3006,0084)	1C	Not supported
>>RT ROI Relationship	(3006,0033)	3	Not supported
>RT ROI Identification Code Sequence	(3006,0086)	3	Not supported
>Related RT ROI Observations Sequence	(3006,00A0)	3	Not supported
>>Observation Number	(3006,0082)	1C	Not supported
>RT ROI Interpreted Type	(3006,00A4)	2	
>ROI Interpreter	(3006,00A6)	2	
>Material ID	(300A,00E1)	3	Not supported
>ROI Physical Properties Sequence	(3006,00B0)	3	
>>ROI Physical Property	(3006,00B2)	1C	
>>ROI Physical Property Value	(3006,00B4)	1C	

Table 94. RT General Plan Module

Attribute Name	Tag	Type	Note
RT Plan Label	(300A,0002)	1	VSim Plan ID
RT Plan Name	(300A,0003)	3	Not supported
RT Plan Description	(300A,0004)	3	Not supported
Instance Number	(0020,0013)	3	
Operators' Name	(0008,1070)	2	
RT Plan Date	(300A,0006)	2	
RT Plan Time	(300A,0007)	2	
Treatment Protocols	(300A,0009)	3	Not supported
Treatment Intent	(300A,000A)	3	Not supported
Treatment Sites	(300A,000B)	3	Not supported
RT Plan Geometry	(300A,000C)	1	"PATIENT"
Referenced Structure Set Sequence	(300C,0060)	1C	
>Referenced SOP Class UID	(0008,1150)	1C	
>Referenced SOP Instance UID	(0008,1155)	1C	
Referenced Dose Sequence	(300C,0080)	3	Not supported
>Referenced SOP Class UID	(0008,1150)	1C	Not supported
>Referenced SOP Instance UID	(0008,1155)	1C	Not supported
Referenced RT Plan Sequence	(300C,0002)	3	Not supported
>Referenced SOP Class UID	(0008,1150)	1C	Not supported
>Referenced SOP Instance UID	(0008,1155)	1C	Not supported
>RT Plan Relationship	(300A,0055)	1C	Not supported

Table 95. RT Patient Setup Module

Attribute Name	Tag	Type	Note
Patient Setup Sequence	(300A,0180)	1	
>Patient Setup Number	(300A,0182)	1	
>Patient Position	(0018,5100)	1C	Same patient position as specified in the associated Image Series
>Patient Additional Position	(300A,0184)	1C	
>Fixation Device Sequence	(300A,0190)	3	Not supported
Fixation Device Type	(300A,0192)	1C	Not supported
>>Fixation Device Label	(300A,0194)	2C	Not supported
>>Fixation Device Description	(300A,0196)	3	Not supported
>>Fixation Device Position	(300A,0198)	3	Not supported
>Shielding Device Sequence	(300A,01A0)	3	Not supported
>>Shielding Device Type	(300A,01A2)	1C	Not supported
>>Shielding Device Label	(300A,01A4)	2C	Not supported
>>Shielding Device Description	(300A,01A6)	3	Not supported
>>Shielding Device Position	(300A,01A8)	3	Not supported
>Setup Technique	(300A,01B0)	3	ISOCENTRIC or FIXED SSD
>Setup Technique Description	(300A,01B2)	3	
>Setup Device Sequence	(300A,01B4)	3	Present if "Relative" beam is present.
>>Setup Device Type	(300A,01B6)	1C	LASER_POINTER
>>Setup Device Label	(300A,01B8)	2C	Empty
>>Setup Device Description	(300A,01BA)	3	Not supported.
>>Setup Device Parameter	(300A,01BC)	2C	
>>Setup Reference Description	(300A,01D0)	3	Isocenter name with which relative beam is designed.
>Table Top Vertical Setup displacement	(300A,01D2)	3	0
>Table Top Longitudinal Setup displacement	(300A,01D4)	3	0
>Table Top Lateral Setup Displacement	(300A,01D6)	3	0

Table 96. RT Beam Module

>Referenced Dose Sequence	(300C,0080)	3	Not supported
>>Referenced SOP Class UID	(0008,1150)	1C	Not supported
>>Referenced SOP Instance UID	(0008,1155)	1C	Not supported
>Number of Wedges	(300A,00D0)	1	
>Wedge Sequence	(300A,00D1)	1C	Not supported
>>Wedge Number	(300A,00D2)	1C	Not supported
>>Wedge Type	(300A,00D3)	2C	Not supported
>>Wedge ID	(300A,00D4)	3	Not supported
>>Wedge Angle	(300A,00D5)	2C	Not supported
>>Wedge Factor	(300A,00D6)	2C	Not supported
>>Wedge Orientation	(300A,00D8)	2C	Not supported
>>Source to Wedge Tray Distance	(300A,00DA)	3	Not supported
>Number of Compensators	(300A,00E0)	1	
>Total Compensator Tray Factor	(300A,00E2)	3	Not supported
>Compensator Sequence	(300A,00E3)	1C	Not supported
>>Compensator Number	(300A,00E4)	1C	Not supported
>>Compensator Type	(300A,00EE)	3	Not supported
>>Material ID	(300A,00E1)	2C	Not supported
>>Compensator ID	(300A,00E5)	3	Not supported
>>Source to Compensator Tray Distance	(300A,00E6)	2C	Not supported
>>Compensator Rows	(300A,00E7)	1C	Not supported
>>Compensator Columns	(300A,00E8)	1C	Not supported
>>Compensator Pixel Spacing	(300A,00E9)	1C	Not supported
>>Compensator Position	(300A,00EA)	1C	Not supported
>>Compensator Transmission Data	(300A,00EB)	1C	Not supported
>>Compensator Thickness Data	(300A,00EC)	1C	Not supported
>Number of Boli	(300A,00ED)	1	
>Referenced Bolus Sequence	(300C,00B0)	1C	Not supported
>>Referenced ROI Number	(3006,0084)	1C	Not supported
>Number of Blocks	(300A,00F0)	1	
>Total Block Tray Factor	(300A,00F2)	3	Not supported
>Block Sequence	(300A,00F4)	1C	
>>Block Tray ID	(300A,00F5)	3	Not supported
>>Source to Block Tray Distance	(300A,00F6)	2C	
>>Block Type	(300A,00F8)	1C	
>>Block Divergence	(300A,00FA)	2C	
>>Block Number	(300A,00FC)	1C	
>>Block Name	(300A,00FE)	3	Not supported
>>Material ID	(300A,00E1)	2C	
>>Block Thickness	(300A,0100)	2C	
>>Block Transmission	(300A,0102)	2C	
>>Block Number of Points	(300A,0104)	2C	
>>Block Data	(300A,0106)	2C	
>Applicator Sequence	(300A,0107)	3	Not supported
>>Applicator ID	(300A,0108)	1C	Not supported
>>Applicator Type	(300A,0109)	1C	Not supported
>>Applicator Description	(300A,010A)	3	Not supported
>Final Cumulative Meterset Weight	(300A,010E)	1C	
>Number of Control Points	(300A,0110)	1	

>Control Point Sequence	(300A,0111)	1	
>>Control Point Index	(300A,0112)	1C	
>>Cumulative Meterset Weight	(300A,0134)	2C	
>>Referenced Dose Reference Sequence	(300C,0050)	3	Not supported
>>>Referenced Dose Reference Number	(300C,0051)	1C	Not supported
>>>Cumulative Dose Reference Coefficient	(300A,010C)	2C	Not supported
>>Nominal Beam Energy	(300A,0114)	3	Not supported
>>Dose Rate Set	(300A,0115)	3	Not supported
>>Wedge Position Sequence	(300A,0116)	3	Not supported
>>>Referenced Wedge Number	(300C,00C0)	1C	Not supported
>>>Wedge Position	(300A,0118)	1C	Not supported
>>Beam Limiting Device Position Sequence	(300A,011A)	1C	
>>>RT Beam Limiting Device Type	(300A,00B8)	1C	
>>>Leaf/Jaw Positions	(300A,011C)	1C	
>>Gantry Angle	(300A,011E)	1C	
>>Gantry Rotation Direction	(300A,011F)	1C	
>>Beam Limiting Device Angle	(300A,0120)	1C	
>>Beam Limiting Device Rotation Direction	(300A,0121)	1C	
>>Patient Support Angle	(300A,0122)	1C	
>>Patient Support Rotation Direction	(300A,0123)	1C	
>>Table Top Eccentric Axis Distance	(300A,0124)	3	
>>Table Top Eccentric Angle	(300A,0125)	1C	
>>Table Top Eccentric Rotation Direction	(300A,0126)	1C	
>>Table Top Vertical Position	(300A,0128)	2C	For “Relative” beams Control Point(CP) 1 will have table offset values w.r.t. ISOCENTER as defined in Patient Setup Module and the table values are undefined in CP0. For “Absolute” beams, CP 0 holds table values. Pls. Refer to C.8.8.14.5 & C.8.8.14.6 of Part 3.
>>Table Top Longitudinal Position	(300A,0129)	2C	
>>Table Top Lateral Position	(300A,012A)	2C	
>>Isocenter Position	(300A,012C)	2C	
>>Surface Entry Point	(300A,012E)	3	
>>Source to Surface Distance	(300A,0130)	3	

Table 97. RT Image Module

>KVP	(0018,0060)	2C	Not supported
>X-Ray Tube Current	(0018,1151)	2C	Not supported
>Exposure Time	(0018,1150)	2C	Not supported
>Meterset Exposure	(3002,0032)	2C	Not supported
>Beam Limiting Device Sequence	(300A,00B6)	3	Not supported
>>RT Beam Limiting Device Type	(300A,00B8)	1C	Not supported
>>Source to Beam Limiting Device Distance	(300A,00BA)	3	Not supported
>>Number of Leaf/Jaw Pairs	(300A,00BC)	1C	Not supported
>>Leaf Position Boundaries	(300A,00BE)	2C	Not supported
>>Leaf/Jaw Positions	(300A,011C)	1C	Not supported
>Applicator Sequence	(300A,0107)	3	Not supported
>>Applicator ID	(300A,0108)	1C	Not supported
>>Applicator Type	(300A,0109)	1C	Not supported
>>Applicator Description	(300A,010A)	3	Not supported
>Number of Blocks	(300A,00F0)	1C	Not supported
>Block Sequence	(300A,00F4)	2C	Not supported
>>Block Tray ID	(300A,00F5)	3	Not supported
>>Source to Block Tray Distance	(300A,00F6)	2C	Not supported
>>Block Type	(300A,00F8)	1C	Not supported
>>Block Divergence	(300A,00FA)	2C	Not supported
>>Block Number	(300A,00FC)	1C	Not supported
>>Block Name	(300A,00FE)	3	Not supported
>>Material ID	(300A,00E1)	2C	Not supported
>>Block Thickness	(300A,0100)	3	Not supported
>>Block Number of Points	(300A,0104)	2C	Not supported
>>Block Data	(300A,0106)	2C	Not supported
Gantry Angle	(300A,011E)	3	
Beam Limiting Device Angle	(300A,0120)	3	
Patient Support Angle	(300A,0122)	3	
Table Top Eccentric Axis Distance	(300A,0124)	3	
Table Top Eccentric Angle	(300A,0125)	3	
Table Top Vertical Position	(300A,0128)	3	
Table Top Longitudinal Position	(300A,0129)	3	
Table Top Lateral Position	(300A,012A)	3	

The tag or sequence which is supported by VSim is left blank. If the sequence is not supported, all other tags under the sequence are not supported by VSim.