

AXIOM Artis VB11, VB20, VB21

AX

DICOM Conformance Statement

Revision: 01

Status: Released

Last Change: 28-Jul-2004

Copyright by SIEMENS AG Medical Solutions

- - - -

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

Table of Contents

Table of Contents	3
List of Figures	8
List of Tables	9
Network Conformance Statement	11
1 Introduction	12
1.1 Purpose	12
1.2 Scope	12
1.3 Definitions, Abbreviations	12
1.3.1 Definitions	12
1.3.2 Abbreviations	12
1.4 References	13
1.5 Connectivity and Interoperability	13
2 Implementation Model Verification	15
2.1 Application Data Flow Diagram	15
2.2 Functional Definitions of Application Entities	15
2.3 Sequencing of real World Activities	15
3 Application Entity Specification Verification	16
3.1 Verification AEs Specification	16
3.1.1 Association Establishment Policies	16
3.1.2 Association Initiation Policy	16
3.1.3 Association Acceptance Policy	17
4 Implementation Model Storage	18
4.1 Application Data Flow Diagram	18
4.2 Functional Definitions of Application Entities	18
4.3 Sequencing of real World Activities	19
5 Application Entity Specification Storage	20
5.1 Storage AEs Specification	20
5.1.1 Association Establishment Policies	21
5.1.2 Association Initiation Policy	21
5.1.3 Association Acceptance Policy	26
6 Implementation Model Storage Commit	33
6.1 Application Data Flow Diagram	33
6.2 Functional Definitions of Application Entities	33
6.3 Sequencing of real World Activities	34

7	AE Specification Storage Commitment	35
7.1	Storage Commitment AE Specification	35
7.1.1	Association Establishment Policies	35
7.1.2	Association Initiation Policy	36
8	Implementation Model Query/Retrieve	39
8.1	Application Data Flow Diagram	39
8.2	Functional Definitions of Application Entities	40
8.3	Sequencing of real World Activities	40
9	AE Specification Query/Retrieve	41
9.1	Query/Retrieve Service AEs Specification	41
9.1.1	Association Establishment Policies	42
9.1.2	Association Initiation Policy	42
9.1.3	Association Acceptance Policy	48
10	Implementation Model Print	61
10.1	Application Data Flow Diagram	61
10.2	Functional Definitions of Application Entities	61
10.3	Sequencing of real World Activities	61
11	Application Entity Specification Print	62
11.1	Print Management AE Specification	62
11.1.1	Association Establishment Policies	62
11.1.2	Association Initiation Policy	63
12	Implementation Model Worklist	71
12.1	Application Data Flow Diagram	71
12.2	Functional Definitions of Application Entities	71
12.3	Sequencing of real World Activities	72
13	AE Specification Worklist	73
13.1	Modality Worklist Service AEs Specification	73
13.1.1	Association Establishment Policies	73
13.1.2	Association Initiation Policy	74
14	Communication Profiles	81
14.1	Supported Communication Stacks	81
14.1.1	TCP/IP Stack	81
15	Extensions/Specializations/Private Extensions	82
15.1	Standard Extended/Specialized/Private SOPs	82
15.2	Private Transfer Syntaxes	82

16 Configuration	83
16.1 AE Title / Presentation Address Mapping	83
16.1.1 DICOM Verification	83
16.1.2 DICOM Storage AE Title	83
16.1.3 DICOM Query/Retrieve AE Title	83
16.1.4 DICOM Print AE Title	83
16.1.5 DICOM Modality Worklist AE Title	84
16.2 Configurable Parameters	84
16.2.1 Storage and Query/Retrieve	84
16.2.2 Print	84
16.2.3 Worklist	85
16.3 Default Parameters	85
17 Support of Extended Character Sets	86
M.A.S. Offline Media Application Profile	87
1 Class and Profile Identification	88
2 Clinical Context	88
2.1 Roles and Service Class Options	89
2.1.1 File Set Creator	89
2.1.2 File Set Reader	89
2.1.3 File Set Updater	89
3 PRI-XAMAS Profiles	91
3.1 SOP Classes and Transfer Syntaxes	91
3.2 Physical Media and Media Formats	94
3.3 Directory Information in DICOMDIR	94
3.3.1 Basic Directory IOD Specialization	94
3.3.2 Additional Keys	95
3.3.3 Private Directory Record Keys	96
3.3.4 Icon Images	97
3.4 Other Parameters	97
3.4.1 Multiframe JPEG Format	97
3.4.2 Attribute Value Precedence	97
Application Profile Conformance Statement	99
1 Introduction	100
1.1 Purpose	100
1.2 Scope	100
1.3 Definitions, Abbreviations	100
1.3.1 Definitions	100
1.3.2 Abbreviations	100

1.4	References	101
1.5	Connectivity and Interoperability	101
2	Implementation Model	103
2.1	Application Data Flow Diagram	103
2.2	Functional definitions of AE's	103
2.3	Sequencing of Real World Activities	104
2.4	File Meta Information Options.	104
3	AE Specifications	105
3.1	DICOM Archive Specification.	105
3.1.1	File Meta Information for the Application Entity	105
3.1.2	Real-World Activities for this Application Entity	105
4	Augmented and Private Profiles	108
4.1	Augmented Application Profiles.	108
4.2	Private Application Profiles	108
4.2.1	PRI-XAMAS-CD	108
5	Extensions, Specializations and Privatizations of SOP Classes and Transfer Syntaxes	109
5.1	SOP Specific Conformance Statement for Basic Directory	109
5.1.1	Extension, Specialization for SIEMENS Non-Image Objects	109
6	Configuration	110
6.1	AE Title Mapping	110
6.1.1	DICOM Media Storage AE Title	110
7	Support of Extended Character Sets	110
A	Appendix	111
A.1	SIEMENS Private Non-Image IOD	111
A.1.1	Siemens Non-Image IOD - E-R Model.	111
A.1.2	Siemens Non-Image IOD - Module Table	112
A.1.3	Siemens Non-Image IOD - Modules	113
A.2	Siemens Standard Extended Modules	115
A.2.1	CSA Image Header Module	115
A.2.2	CSA Series Header Module.	116
A.2.3	MEDCOM Header Module.	116
A.2.4	MEDCOM OOG Module.	118
A.2.5	Edge Enhancement Module.	119
A.2.6	Angio Viewing Module.	119
A.2.7	Angio Acquisition Data Module	121
A.2.8	Angio Quantitation Module.	123
A.2.9	Original Image Info Module	124

A.2.10 SOP Common Module - Image Type Extensions	125
A.3 Registry of DICOM Data Elements	126
A.4 Extended XA IOD	130
A.5 Private Non-Image IOD	136
A.6 Quant Report Images - derived XA IOD	138
A.7 CT/MR image converted to XA IOD (VA22x and later)	141
A.8 Sorting Order AXIOM Artis (VA10C/VA20F/VA22x) with syngo VA47Z/VA49B 148	
A.8.1 Identification of Images	149
A.8.2 Handling of Performed Procedure Step ID.	151
A.8.3 The sorting algorithms of Artis Viewer	151
A.8.4 The sorting algorithms of Artis Browser.	154
A.9 DICOM Print SCU - detailed status displays	155
A.9.1 Common Status Information	155
A.9.2 Additional Status Information - AGFA Printers	161
A.9.3 Additional Status Information - Kodak PACS Link (formerly Imation). . .	161
A.9.4 Additional Status Information - Kodak 1901.	162
A.9.5 Additional Status Information - Kodak 2180/1120	163
A.9.6 Additional Status Information - Codonics	163
A.9.7 Additional DICOM Execution Status Information.	165

List of Figures

Figure 1:	Application Data Flow Diagram - Verification SCU	15
Figure 2:	Application Data Flow Diagram - Storage SCU.....	18
Figure 3:	Application Data Flow Diagram - Storage SCP	18
Figure 4:	Application Data Flow Diagram Storage Commitment SCU	33
Figure 5:	Application Data Flow Diagram - Query/Retrieve SCU	39
Figure 6:	Application Data Flow Diagram - Query/Retrieve SCP	40
Figure 7:	Application Data Flow Diagram - Print SCU	61
Figure 8:	Application data flow diagram - Basic Worklist SCU:	71
Figure 1:	DICOM Information Model.....	112

List of Tables

Table 1:	Siemens AXIOM Artis Products	12
Table 2:	Initiation presentation context - Verification	17
Table 3:	SOP Classes as an Storage SCU	20
Table 4:	Private SOP Classes as an Storage SCU	20
Table 5:	SOP Classes as Storage SCP	20
Table 6:	PrivateSOP Classes as an Storage SCP	20
Table 7:	Initiation presentation context - Storage	22
Table 8:	Acceptable presentation contexts - Storage	27
Table 9:	Applied Defaults Dictionary of DICOM type 2 and 3 IOD Attributes	28
Table 10:	Storage Commitment SOP Class as an SCU	35
Table 11:	Proposed presentation contexts - request Storage commitment	36
Table 12:	Query/Retrieve SOP Classes as an SCU	41
Table 13:	Query/Retrieve SOP Classes as an SCP	41
Table 14:	Patient Root and Study Root query attributes	43
Table 15:	C-FIND response status	46
Table 16:	C-MOVE response status	47
Table 17:	Patient level attributes, Patient Root Information Model	50
Table 18:	Study level attributes, Patient Root Information Model	50
Table 19:	Series level attributes, Patient Root Information Model	51
Table 21:	Study level attributes, Study Root Information Model	52
Table 20:	Image level attributes, Patient Root Information Model	52
Table 22:	Series level attributes, Study Root Information Model	53
Table 23:	Image level attributes, Study Root Information Model	54
Table 24:	Patient instance level, Patient Study Only Information Model	54
Table 25:	Study level attributes, Patient Study Only Information Model	55
Table 26:	C-FIND SCP return status	56
Table 27:	C-GET SCP return status	57
Table 28:	C-MOVE return status	59
Table 29:	Basic Grayscale Print Management Meta SOP Classes	62
Table 30:	Used Basic Film Session N-CREATE-RQ attributes	65
Table 31:	Attributes of N-DELETE-RQ on Basic Film Session SOP Class	65
Table 32:	Basic Film Session SOP status	66
Table 33:	Used Film Box N-CREATE-RQ attributes	66
Table 34:	Attributes of N-DELETE-RQ on Basic Film Box SOP Class	67
Tabelle 35:	Basic Film Box SOP status	67
Table 36:	Mandatory Basic Grayscale Image Box N-SET attributes	68
Table 37:	Basic Grayscale Image Box SOP status	68
Table 38:	Used Printer N-EVENT report attributes	69
Table 39:	Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes	69
Table 40:	Used Print Job N-EVENT report attributes	70
Table 41:	SOP Classes as an SCU	73
Table 42:	Proposed presentation contexts	74
Table 43:	Supported Broad Worklist Query Search Key Attributes	75
Table 44:	Basic Worklsit C-FIND-RSP Return Key Attributes	75

Table 45:	Patient based "narrow query" Search Key Attributes	79
Table 46:	C-FIND Response Status	80
Table 1:	PRI-XAMAS Universal Angio/Radio Class Profiles.....	88
Table 2:	PRI-XAMAS SOP Classes and Transfer Syntaxes.....	91
Table 3:	PRI-XAMAS Additional DICOMDIR Keys	95
Table 4:	PRI-XAMAS Additional PRIVATE Record Keys	96
Table 1:	SIEMENS AXIOM Artis based products	100
Table 2:	100
Table 3:	Application profiles, Activities, and Roles for DICOM Archive	105
Table 4:	STD-GEN-CD Supported SOP Classes for FSR role	106
Table 5:	Basic Directory Extension for Non-Image Objects.....	109
Table 1:	CSA Non-Image IOD Modules.....	112
Table 2:	CSA Non-Image Module	113
Table 3:	CSA Image IOD Modules	115
Table 4:	CSA Image Header Module	115
Table 5:	CSA Series Header Module	116
Table 6:	MEDCOM Header Module.....	116
Table 7:	MEDCOM History Information.....	118
Table 8:	MEDCOM OOG Module.....	118
Table 9:	Edge Enhancement Module	119
Table 10:	Angio Viewing Module.....	119
Table 11:	Angio Viewing Module.....	121
Table 12:	Angio Viewing Module.....	123
Table 13:	Original Image Info Module	124
Table 14:	Image Type Extentions.....	125
Table 15:	overview of supplied attributes - X-Ray Angiographic image (Original, Derived). 130	
Table 16:	Data Dictionary of private IOD Definition - Study-Report	136
Table 17:	overview of supplied attributes - X-Ray Angiographic image (Quant Report Image) 138	
Table 18:	overview of supplied attributes - XA image (derived) from CT image	142
Table 19:	overview of supplied attributes - XA image (derived) from MR image	145
Table 20:	Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class	155
Table 21:	Printer Status Infos: Additional Agfa printer status infos	161
Table 22:	Printer Status Infos: Additional Kodak infos for Pacs Link (formerly Imation cam- eras)	161
Table 23:	Printer Status Infos: Additional Kodak infos for Kodak 190.....	162
Table 24:	Printer Status Infos: Additional Kodak infos for 2180/1120	163
Table 25:	Printer Status Infos: Additional Codonics infos.....	163
Table 26:	Additional DICOM Execution Status Infos	165

Network Conformance Statement

This part contains the Conformance declaration for the DICOM Network Services:

- Storage - User/Provider
- Storage - Commitment (Push Model) User
- Query/Retrieve - User/Provider
- Basic Greyscale Print - User
- Modality Worklist - User

1 Introduction

1.1 Purpose

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

The applications described in this conformance statement are part of the SIEMENS AXIOM Artis based on syngo® software. The AXIOM Artis DICOM network implementation acts as SCU and SCP for the DICOM Storage and Query/Retrieve services and as SCU for the DICOM Print, DICOM BasicWorklist and DICOM Storage Commitment Services.

Verification is supported in SCU (only via Service environment) and SCP role.

1.2 Scope

This DICOM Conformance Statement refers to SIEMENS AXIOM Artis based products using AXIOM Artis software. The following table relates software names to SIEMENS AXIOM Artis products.

Table 1: Siemens AXIOM Artis Products

Software Name	SIEMENS AXIOM Product
VB11x	AXIOM Artis FC,FA,BC,BA,MP,TA

1.3 Definitions, Abbreviations

1.3.1 Definitions

DICOM	Digital Imaging and Communications in Medicine
DIMSE	DICOM Message Service Element
DIMSE-C	DICOM Message Service Element with Composite information objects

1.3.2 Abbreviations

ACR	American College of Radiology
AE	DICOM Application Entity
ASCII	American Standard Code for Information Interchange
DCS	DICOM Conformance Statement
IOD	DICOM Information Object Definition
ISO	International Standard Organization
MAS or M.A.S.	Modular Angio System
NEMA	National Electrical Manufacturers Association
O	Optional Key Attribute

PDU	DICOM Protocol Data Unit
R	Required Key Attribute
RIS	Radiology Information System
RWA	Real-World Activity
SCP	DICOM Service Class Provider (DICOM server)
SCU	DICOM Service Class User (DICOM client)
SOP	DICOM Service-Object Pair
U	Unique Key Attribute

1.4 References

- [1] Digital Imaging and Communications in Medicine (DICOM), NEMA PS 3.1-16

1.5 Connectivity and Interoperability

This Conformance Statement by itself does not guarantee successful interoperability of SIEMENS equipment with non-SIEMENS equipment. The user (user's agent) should be aware of the following issues:

- Interoperability

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of SIEMENS equipment with non-SIEMENS equipment. It is the user's responsibility to analyse thoroughly the application requirements and to specify a solution that integrates SIEMENS equipment with non-SIEMENS equipment.

- Validation

SIEMENS equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.

Where SIEMENS equipment is linked to non-SIEMENS equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation test will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.

- New versions of the DICOM Standard

the DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. SIEMENS is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, SIEMENS reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-SIEMENS provider linking to SIEMENS equipment,

also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into SIEMENS equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

2 Implementation Model Verification


The Siemens AXIOM Artis DICOM Service Tool application requests Verification to verify the ability of an foreign DICOM application on a remote node to respond to DICOM messages.

Responding to Verification requests from remote nodes is handled by the Storage SCP application.

2.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows XP application and acts as SCU for the C-ECHO DICOM network service.

Figure 1: Application Data Flow Diagram - Verification SCU



2.2 Functional Definitions of Application Entities

The Siemens AXIOM Artis DICOM Service Tool application opens an association when a "Test" of a remote application is requested during a configuration session. This can be done when entering new data for remote application configuration or to verify existing configuration data.

2.3 Sequencing of real World Activities

Newly entered data have to be saved first, before a "test" of these data is possible.

3 Application Entity Specification Verification

3.1 Verification AEs Specification

3.1.1 Association Establishment Policies

3.1.1.1 General

The Siemens AXIOM Artis DICOM Service Tool application attempts to open an association for verification request whenever the "Test" function is activated during network configuration of a remote DICOM application.

3.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM Service Tool application initiates one association at a time to request verification.

3.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

3.1.1.4 Implementation Identifying Information

The Siemens <product> DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20010101

and an Implementation Version Name of

- "SIEMENS_SWFVC20B"

3.1.2 Association Initiation Policy

The Siemens AXIOM Artis DICOM Service Tool application attempts to initiate a new association for

- DIMSE C-ECHO

service operations.

3.1.2.1 Associated Real-World Activity

3.1.2.1.1 Associated Real-World Activity -Request Verification

The associated Real-World activity is a C-ECHO request initiated by Service and Configuration SW environment whenever a "test" is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than "Success" this will be indicated in the service environment and the association is closed.

3.1.2.1.2 Proposed Presentation Contexts

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Table 2: Initiation presentation context - Verification

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotia- tion
Name	UID	Name List	UID List		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None

3.1.2.1.3 SOP Specific Conformance - Verification SCU

None.

3.1.3 Association Acceptance Policy

The Verification SCP is part of the Storage SCP - see section 5.1.3 on page 26

4 Implementation Model Storage

The Siemens AXIOM Artis DICOM Application Entity both originates associations for Storage of DICOM Composite Information Objects in Remote Application Entities and receives association requests for Storage from Remote Application Entities.

4.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows XP application and acts as SCU and SCP for the C-STORE DICOM network service and as SCP for the C-ECHO DICOM network service.

Figure 2: Application Data Flow Diagram - Storage SCU





Figure 3: Application Data Flow Diagram - Storage SCP



4.2 Functional Definitions of Application Entities

The Storage SCU is invoked by the job control interface that is responsible for processing network archival tasks. The job consists of data describing the composite image objects selected for storage and the destination. An association is negotiated with the destination application entity and the image data is transferred using the C-STORE DIMSE-Service. Status of the transfer

is reported to the job control interface. Sending can also be activated by the "Send 512 to..." interface of the Viewer. It can be configured that sending in 512 format is the default mode of operation for the viewer's send functionality.

The Storage SCP component of the Siemens AXIOM Artis DICOM application is operating as background server process. It is existing when the machine is powered on and waits for Storage association requests. Upon accepting an association with a negotiated Presentation Context it starts to receive the Composite Image Objects and imports them to local database. Verification requests will be processed and responded by Storage SCP component too.

4.3 Sequencing of real World Activities

not applicable.

5 Application Entity Specification Storage

5.1 Storage AEs Specification

The AXIOM Artis Storage service class user/service class provider applications use one AE when initiating/receiving associations to/from remote DICOM nodes.

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as an SCU:

Table 3: SOP Classes as an Storage SCU

SOP Class Name	SOP Class UID
X-Ray Angiographic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.1

SIEMENS AXIOM Artis DICOM products provide Private Conformance to the following DICOM V3.0 conform private SOP Classes as an SCU:

Table 4: Private SOP Classes as an Storage SCU

SOP Class Name	SOP Class UID
CSA Non-Image	1.3.12.2.1107.5.9.1

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as SCP:

Table 5: SOP Classes as Storage SCP

SOP Class Name	SOP Class UID
MR Image Information Object Storage	1.2.840.10008.5.1.4.1.1.4
CT Image Information Object Storage	1.2.840.10008.5.1.4.1.1.2
X-Ray Angiographic Image Information Object Storage	1.2.840.10008.5.1.4.1.1.12.1
Verification	1.2.840.10008.1.1

SIEMENS AXIOM Artis DICOM products provide Private Conformance to the following DICOM V3.0 conform private SOP Classes as an SCP:

Table 6: PrivateSOP Classes as an Storage SCP

SOP Class Name	SOP Class UID
CSA Non-Image	1.3.12.2.1107.5.9.1

5.1.1 Association Establishment Policies

5.1.1.1 General

The existence of a job queue entry with network destination will activate the DICOM Storage Application. An association request is sent to the destination AE and upon successful negotiation of a Presentation Context the transfer is started.

The default PDU size used will be 28 KB.

5.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM application initiates several associations at a time, one for each destination to which a transfer request being processed in the active job queue list.

The Siemens AXIOM Artis DICOM application is able to accept multiple associations at a time. It can handle up to ten associations in parallel.

The number of simultaneous DICOM associations can be configured via the Service UI. The dialog can be found in "Configuration/DICOM/General".

5.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

5.1.1.4 Implementation Identifying Information

The Siemens AXIOM Artis DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20010101

and an Implementation Version Name of

- "SIEMENS_SWFVC20B"

5.1.2 Association Initiation Policy

If a job with network destination gets active in the job list the Siemens AXIOM Artis DICOM application attempts to initiate a new association for

- DIMSE C-STORE

service operations.

5.1.2.1 Associated Real-World Activity

5.1.2.1.1 Associated Real-World Activity -Send Image Objects to a Network destination

The associated Real-World activity is a C-STORE request initiated by an internal daemon process triggered by a job with network destination. If the process successfully establishes an association to a remote Application Entity, it will transfer each image one after another via the open association. If the C-STORE Response from the remote Application contains a status other than "Success" or "Warning" the association is aborted.

5.1.2.1.2 Proposed Presentation Contexts

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Table 7: Initiation presentation context - Storage

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCU	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCU	None
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCU	None

Table 7: Initiation presentation context - Storage

CSA Non-Image Storage	1.3.12.2.1107.5.9.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1	SCU	None
-----------------------	---------------------	--	---	-----	------

Note

JPEG compression Transfer Syntaxes are supported only for monochrome images (Photometric Interpretation = "MONOCHROME1" or "MONOCHROME2")

5.1.2.1.3 SOP Specific Conformance to Storage SOP classes

The AXIOM Artis System itself will create XA IOD type images from image acquisition and postprocessing applications. Furthermore a private "Non-Image IOD" will be used to store numerical data and textual report results (e.g. Study Report). The XA IOD will be a "Standard Extended XA Storage" SOP Class. For certain imported/received IOD Types a conversion to a new XA image will be performed in order to enable viewing of those Instances on the AXIOM Artis based System.

A AXIOM Artis System together with its DICOM application will not change private attributes as long as no modification is done. During a "Store Monitor" or "Store Reference Image" operation all private attributes not defined within the AXIOM Artis DICOM application will be removed when the new object instance is created.

For association and DIMSE level time-outs, please refer to the Configuration section of this document.

5.1.2.1.3.1 Optional Attributes

Data Dictionary of DICOM Type 2 and 3 IOD Attributes

Please see "A.4 Table 15:overview of supplied attributes - X-Ray Angiographic image (Original, Derived)" in the Appendix for a list of all DICOM IOD attributes of type 2 and 3 which are encoded in the AXIOM Artis application.

5.1.2.1.3.2 Specialized Information Object Definitions

The DICOM images created by Siemens AXIOM Artis DICOM application conform to the DICOM IOD definitions (Standard extended IODs). But they will contain additional private elements which have to be discarded by a DICOM system when modifying the image.

The DICOM nodes are responsible for data consistency when modifying images. All unknown private attributes have to be removed upon modification!

Data Dictionary of applied private IOD Attributes

Please see "A.2 Siemens Standard Extended Modules" in the Appendix for definition of the private attributes and the IOD tables in the Appendix for usage in the related IOD encoding.

5.1.2.1.3.3 Image Pixel Attribute Description for Grayscale Images

The Siemens AXIOM Artis DICOM application supports the Monochrome2 Photometric Interpretation with the unsigned integer 16 bit grayscale pixel and graphic overlay format. The lower 12 bits are used for pixel and the higher 4 bits are used for the graphic overlay:

Pixel plane

- + samples per pixel (attribute 0028, 0002) = 1
- + photometric interpretation (attribute 0028,0004) = "MONOCHROME2"
- + pixel representation (attribute 0028, 0103) = 0
- + bits allocated (attribute 0028, 0100) = 8, 16
- + bits stored (attribute 0028,0101) = 8, 10, 12
- + high bit (attribute 0028,0102) = 7, 9, 11

Overlay plane

- + overlay type (attribute 60xx, 0040) = "G"
- + bits allocated (attribute 60xx, 0100) = 1
- + bit position (attribute 60xx, 0102) = 0
- + overlay data (attribute 60xx, 3000) = supported (with "bits allocated = 8")

5.1.2.1.3.4 Private Information Object Definitions

To fulfill all application requirements the AXIOM Artis DICOM implementation will use private IOD's to store Data currently not defined in the DICOM Standard according to the DICOM information model. The privately defined IOD will contain all references to identify the Patient/Study/Series/Images to which the related information is belonging to.

This private Information will include

- CARE-Graph
- Complete Study Report

All IOD used to store this private Information will be based on various Instances of the SIEMENS Non-Image IOD.

Please see "A.1:SIEMENS Private Non-Image IOD" and "A.5:Private Non-Image IOD" in the Appendix for a detailed overview of the private IOD definition and the IOD tables for value encoding.

5.1.2.1.4 Associated Real-World Activity -Send Image Objects 512 to a Network destination

From the Viewing application it is possible to issue special send requests for the job queue. The images in this case are "pre-processed" to be compliant to the Basic Cardiac Image definitions (512² Matrix with 8-Bits allocated). Consecutive to the pre-processing stage the transfer is processed as with the normal "Send to..." RWA.

5.1.2.1.5 Proposed Presentation Contexts - "Send 512 to ..."

The proposed Presentation Contexts for this RWA are identical to those described in section 5.1.2.1.2.

5.1.2.1.6 SOP Specific Conformance - "Send 512 to..."

The AXIOM Artis System itself will create XA IOD type images from image acquisition and postprocessing applications. Furthermore a private "Non-Image IOD" will be used to store numerical data and textual report results (e.g. Study Report). The XA IOD will be a "Standard Extended XA Storage" SOP Class. Please see the Appendix for detailed descriptions.

The "Send 512 to..." will allow conversion of any image that is displayed in the AXIOM Artis viewing application to an image that conforms to the Basic Cardiac format. Conversion is limited to reduction of Matrix Size and Pixel Depth and is restricted to be applied to "image-type" instances only. The converted image is a temporary copy and therefore a new Instance UID is generated.

Quant Report images will not be downsized and are sent "as is".

The SOP Instance UID of the converted images will be created in a reproducible way by appending a ".512" to the original UID.

The conversion will generate a Pixel Matrix output that is described by the following Attribute values:

- [0029,0010] Rows = **512**
- [0028,0011] Columns = **512**
- [0028,0100] Bits allocated = **8**
- [0028,0101] Bits Stored = **8**
- [0028,0102] High Bit = **7**

In order to preserve integrity of the data with respect to the original Matrix Size and Pixel Depth, the following attributes will be adapted accordingly for the resulting image:

- [0018,1602] Shutter Left Vertical Edge
- [0018,1604] Shutter Right Vertical Edge
- [0018,1606] Shutter Upper Horizontal Edge
- [0018,1608] Shutter Lower Horizontal Edge
- [0018,1610] Center of Circular Shutter
- [0018,1612] Radius of Circular Shutter
- [0018,1702] Collimator Left Vertical Edge
- [0018,1704] Collimator Right Vertical Edge
- [0018,1706] Collimator Upper Horizontal Edge
- [0018,1708] Collimator Lower Horizontal Edge
- [0018,1710] Center of Circular Collimator
- [0018,1712] Radius of Circular Collimator

- [0018,1720] Vertices of the Polygonal Collimator
- [0018,1164] Imager Pixel Spacing
- [0028,1050] Window Center
- [0028,1051] Window Width
- [0028,3002] MOD LUT Descriptor
- [0028,3006] MOD LUT Data
- [6000,0010] Overlay Rows
- [6000,0011] Overlay Columns
- [6000,3000] Overlay Data
- private data attributes

5.1.3 Association Acceptance Policy

The Siemens AXIOM Artis DICOM application attempts to accept a new association for

- DIMSE C-ECHO
- DIMSE C-STORE

service operations. Any Information Objects transmitted on that association will be checked on conformance and stored in database if check was successful.

5.1.3.1 Associated Real-World Activity

5.1.3.1.1 Associated Real-World Activity - Receiving Image Objects from a remote Node

The daemon receiving process will accept an association and will receive any images transmitted on that association and will store the images on disk in the own data base if the conformance check is performed successfully.

5.1.3.1.2 Presentation Context Table

The Siemens AXIOM Artis DICOM application will accept Presentation Contexts as shown in the following table:

The AXIOM Artis DICOM storage provider application will accept CT and MR IOD.

The AXIOM Artis DICOM storage provider will only accept MONOCHROME_{EX} encoded images. Color encoded images (RGB or Palette Color) are not supported.

Table 8: Acceptable presentation contexts - Storage

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCP	None
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCP	None
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian JPEG Lossless Process 14 (selection value 1) JPEG Lossy Baseline (Process 1) compression JPEG Lossy Extended (Process 2 & 4)	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.4.70 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	SCP	None
CSA Non-Image Storage	1.3.12.2.1107.5.9.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1	SCP	None
Verification	1.2.840.10008.1.1	Implicit VR Little Endian Explicit VR Big Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2 1.2.840.10008.1.2.1	SCP	None

Note

JPEG compression Transfer Syntaxes are supported only for monochrome images (Photometric Interpretation = "MONOCHROME1" or "MONOCHROME2")

5.1.3.1.3 SOP Specific Conformance Statement

The Siemens AXIOM Artis DICOM application conforms to the Full Storage Service Class at Level 2. In the event of a successful C-STORE operation, the image has successfully been written on disk in the Siemens AXIOM Artis image format. For Private Attributes of VR=SQ only a nesting level of one is supported. This means that Private Sequences containing another Sequence will be removed from the image header during Storage.

The AXIOM Artis DICOM receiver returns the status Success upon successful operation otherwise one of the following status codes is returned and the association is aborted:

- Refused (0xA700):
This error status indicates a lack of Resources (e.g. not enough disk space) on the AXIOM Artis modality.
- Error (0xA900 or 0xC000):
An error occurred while processing the image which makes it impossible to proceed. The image will not be stored and the association is aborted.
- Processing Error (0x0110)
An error occurred while processing the image, which makes it impossible to proceed.

Attention! Only after sending the response, the image will be saved into the database. If during this operation an error occurs, the association will be aborted. This implies that a C-STORE-RSP with status SUCCESS does not necessarily mean that the image was successfully stored into the database.

If an image instance is received that is identified by an SOP Instance UID that is already used by an Instance stored in database then the actual received image will be discarded. The existing Instance is not superseded.

5.1.3.1.3.1 Optional Attributes acceptance criterion - applied defaults

Upon receiving IOD other than created by the AXIOM Artis DICOM application a check for conformity is done and with successful completion the IOD is accepted for storage. Hereby defaults will be set when optional attributes are encountered with "no value" coding. Please see the following tables for details.

Attributes for minimal acceptance criterion will be listed. Any defaults applicable are denoted.

Table 9: Applied Defaults Dictionary of DICOM type 2 and 3 IOD Attributes

Attribute Name	Tag	Default Value
Study Date	[0008,0020]	<unknown, zero length>
Image Date	[0008,0023]	<unknown, zero length>
Study Time	[0008,0030]	<unknown, zero length>
Image Time	[0008,0033]	<unknown, zero length>

Table 9: *Applied Defaults Dictionary of DICOM type 2 and 3 IOD Attributes*

Attribute Name	Tag	Default Value
Accession Number	[0008,0050]	<unknown, zero length>
Manufacturer	[0008,0070]	<unknown, zero length>
Referring Physician's Name	[0008,0090]	<unknown, zero length>
Study Description	[0008,1030]	<unknown, zero length>
Patient's Name	[0010,0010]	<unknown, zero length>
Patient ID	[0010,0020]	<unknown, zero length>
Patient's Birth Date	[0010,0030]	<unknown, zero length>
Patient's Sex	[0010,0040]	<unknown, zero length>
Patient Comments	[0010,4000]	<unknown, zero length>
Contrast/Bolus Agent	[0018,0010]	<unknown, zero length>
KVP	[0018,0060]	<unknown, zero length>
Exposure Time	[0018,1150]	<unknown, zero length>
X-Ray Tube Current	[0018,1151]	<unknown, zero length>
Exposure	[0018,1152]	<unknown, zero length>
Study ID	[0020,0010]	"_"
Series Number	[0020,0011]	<unknown, zero length>
Image Number	[0020,0013]	<unknown, zero length>
Patient Orientation	[0020,0020]	<unknown, zero length>
Window Center	[0028,1050]	Open window, e.g. bits stored=12, window center = 2048
Window Width	[0028,1051]	Open window, e.g. bits stored=12, window width = 4095
Recommended Viewing Mode	[0028,1090]	<unknown, zero length>
Representative Frame Number	[0028,6010]	<unknown, zero length>
Calibration Image	[0050,0004]	<unknown, zero length>

5.1.3.1.3.2 Image Pixel Attribute Acceptance Criterion for Grayscale Images

The Siemens AXIOM Artis Image Viewing application accepts the monochrome 1 and monochrome 2 photometric interpretation pixel format and graphic overlay with unsigned integer and 8 or 16 bits allocated. Accepted values:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = "MONOCHROME1", "MONOCHROME2"
- pixel representation (attribute 0028, 0103) = 0
- bits allocated (attribute 0028, 0100) = 8, 16
- bits stored (attribute 0028,0101) = 8, 10, 12
- high bit (attribute 0028,0102) = 7, 9, 11
- Only aspect ratio 1:1 is supported

Overlay plane

- + overlay type (attribute 60xx, 0040) = "G"
- + bits allocated (attribute 60xx, 0100) = 16
- + bit position (attribute 60xx, 0102) = 12, 13, 14, 15

Graphic Overlay will be shifted to fill Overlay Planes from Bit 12 and consecutive.

Overlay plane

- + overlay type (attribute 60xx, 0040) = "G"
- + bits allocated (attribute 60xx, 0100) = 1
- + bit position (attribute 60xx, 0102) = 0
- + overlay data (attribute 60xx, 3000) = supported

The Siemens AXIOM Artis Image Viewing application accepts also the MONOCHROME1 and MONOCHROME2 photometric interpretation pixel format with binary 2's complement integer and 16 bits allocated. Accepted values:

Pixel plane

- samples per pixel (attribute 0028, 0002) = 1
- photometric interpretation (attribute 0028,0004) = "MONOCHROME1", "MONOCHROME2"
- pixel representation (attribute 0028, 0103) = 1
- bits allocated (attribute 0028, 0100) = 16
- bits stored (attribute 0028,0101) = 16
- high bit (attribute 0028,0102) = 15

Overlay plane

- + overlay type (attribute 60xx, 0040) = "G"
- + bits allocated (attribute 60xx, 0100) = 1
- + bit position (attribute 60xx, 0102) = 0

- + overlay data (attribute 60xx, 3000) = supported.

If the MOD LUT SQ contains multiple LUTs, then only the first one is used.

For VOI LUT, only the linear LUT (Window Center/Width) is supported. The contents VOI LUT SQ are ignored.

With multi-frame images the Overlay Graphics is only supported as replication of the 1st frame's overlay data.

5.1.3.1.4 Presentation Context Acceptance Criterion

The Siemens AXIOM Artis DICOM application will accept any number of verification or storage SOP classes that are listed above. There is no limit on the number of presentation contexts accepted, except for the DICOM limit. In the event that the Siemens AXIOM Artis DICOM application runs out of resources, it will reject the association request.

5.1.3.1.5 Transfer Syntax Selection Policies

The Siemens AXIOM Artis DICOM application currently supports

- the Implicit VR Little Endian, the Explicit VR Little Endian and Explicit VR Big Endian Transfer Syntaxes
- the JPEG Lossless Non-hierarchical Transfer Syntax
- the JPEG Baseline and JPEG Extended Transfer Syntaxes (JPEG Lossy).

Any proposed presentation context which includes one of these transfer syntaxes will be accepted. Any proposed presentation context that does not include one of these transfer syntaxes will be rejected.

With Implicit VR Little Endian Transfer Syntax the AXIOM Artis DICOM application will remove any Private Attributes not known to the application. Decision on removal of a Private Element is done if there is NO entry in the dictionary of attributes of the AXIOM Artis DICOM application.

Therefore any Explicit VR Transfer Syntax shall preferably be used by the Storage SCU's when sending Composite Image Instances to the AXIOM Artis DICOM application.

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.


6 Implementation Model Storage Commit

The Storage Commitment service class defines an application-level class of service which facilitates the commitment to storage. It performs an additional task of commitment of composite objects apart from the network based storage of images as defined by the Storage Service class. The AXIOM Artis DICOM implementation supports the Storage Commitment Push Model as SCU.

6.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows XP application and acts as SCU for the Storage Commitment Push Model Service using the Storage Commitment Service Class.

Figure 4: Application Data Flow Diagram Storage Commitment SCU



6.2 Functional Definitions of Application Entities

With each successfully completed send job, the AXIOM Artis DICOM Application will create a Storage Commitment Push Model Identifier from the SOP Instances sent. Then an a Storage Commit Request is triggered. Depending on configuration, the AXIOM Artis DICOM application will keep the association open for responses with a configurable time-out, or closes the association and expects responses on a different association that has to be established by the remote Storage Commitment SCP.

The commitment status derived from the related trigger response will be indicated in the related Status Flags of the related entity. It is possible to create triggers ("auto rules") from this event.

The Transaction UUIDs of the pending commitment request are kept "open" for a configurable amount in time (default: 1h). If the "open time" for a pending commitment request has elapsed w/o a related response from the provider, the Transaction UUID is removed and the related entities are indicated as "commit failed".

In any case, commitment will only be requested for previously and successfully sent images.

6.3 Sequencing of real World Activities

The Storage Commitment trigger is automatically derived from the successful completion of a Send Job.

7 AE Specification Storage Commitment

7.1 Storage CommitmentAE Specification

SIEMENS AXIOM Artis DICOM application provides Standard Conformance to the following DICOM V3.0 SOP Class as an SCU:

Table 10: Storage Commitment SOP Class as an SCU

SOP Class Name	SOP Class UID
Storage Commitment Push Model	1.2.840.10008.1.20.1

7.1.1 Association Establishment Policies

7.1.1.1 General

With a Send Job successfully completed, the DICOM application will generate a Storage Commitment Identifier which references to all Instances of the processed job. The Commit Request is then sent over a single opened association. The AXIOM Artis will wait for Status responses of the Storage Commitment Request. If the Provider accepts the Storage Commitment with Success Status, the generated Transaction UID, together with study identification data and a timestamp, is kept. Depending on configuration, the association is closed when the configured timeout has elapsed or a response was received before. If the association is closed before a response was received, the response is then expected on a different association. Multiple Storage Commitment Requests can be pending.

The default PDU size used will be 28 KB.

7.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM application initiates several associations at a time, one for each storage commitment request being processed.

7.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

7.1.1.4 Implementation Identifying Information

The Siemens AXIOM Artis DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- "SIEMENS_SWFVX47A"

7.1.2 Association Initiation Policy

The AXIOM Artis DICOM Application Entity acts as a Service Class User (SCU) for the

- Storage Commitment Push Model Service Class (to request commitment for storage of instances previously sent).

To do so, the AXIOM Artis will issue a

- N-ACTION DIMSE to request commitment or a
- N-EVENT-REPORT DIMSE to respond to a received storage commitment request and the association was closed by the remote system prior to response.

7.1.2.1 Real World Activity

7.1.2.1.1 Associated Real-World Activity - Job Completed

The AXIOM Artis Storage Commitment application sends the commit request (N-ACTION-RQ) message and waits for acceptance of this request (N-ACTION-RSP). After receiving this, the transaction is marked as "waiting".

Depending on a configuration value, the association will then be closed or kept open. In the first case, there is another configurable value giving the number of minutes (by default 60) to wait for the corresponding commit response (N_EVENT_REPORT). In the second case, this time is the (also configurable) time-out for the association. For both cases, if the commit response (N_EVENT_REPORT) does not arrive during the configured time, the transaction will be marked as failed. The AXIOM Artis does not resend objects from a failed Storage Commitment result in any case.

If the commit response (N_EVENT_REPORT) received has the status of "complete - failure exists", the transaction is marked as failed, else the transaction is marked as "completed"; In both cases, a message is shown to the user.

7.1.2.1.2 Proposed Presentation Contexts - Job Completed

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Table 11: Proposed presentation contexts - request Storage commitment

Presentation Context Table				
Abstract Syntax		Transfer Syntax		Role
Name	UID	Name List	UID List	
				Extended Negotiation

Storage Commitment Push Model	1.2.840.10008.1.20.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

7.1.2.1.3 SOP Specific Conformance Statement- Job Completed

Storage Commitment is supported for all the SOP class UIDs as mentioned in 'Acceptable presentation contexts - Storage' in the Storage SCP section of this document.

The Referenced Study Component Sequence is not supported.

Storage Media File-Set ID and UID Attributes will not be supported in the commitmant request (N-ACTION primitive) invoked by the Storage Commitment SCU.

7.1.2.1.4 Associated Real-World Activity - Update Flags

The AXIOM Artis Storage Commitment DICOM Application has sent a Storage Commitment Request and, being configured to receive response on a separate association, has closed the association, and now it gets an association request from the Storage Commitment SCP that want to send the results. Consecutive to start-up, the AXIOM Artis DICOM application will await Storage commitment Notification triggers. Any incoming Notification will be checked for validity, that is, if the related Transaction UID is still part of the Pending Request Queue.

If the Notification is valid, the Notification Identifier is evaluated and the related Instances marked with the related status. The over-all Commit Status of the higher Information Entities is derived from propagation of the States of all Image entities included in a study.

The Status Flags directly affected by Storage Commitment results and indicated in the different entities of the Patient Browser list can be one of

- "AC" or "SC" - Successful Commitment, A means archived to configured Archive destination, whereas S means sent to any other destination
- "Af" or "Sf" - Commitment failed.
- "A?" or "S?" - Commitment request is sent, response is pending.

In case of failure the user has to repeat the transfer of images to the Archive destination. Another Storage Commitment will be performed after sending is completed successfully.

7.1.2.1.5 SOP-specific Conformance Statement - Update Flags

If the Commitment response (N_EVENT_REPORT) received has the status of "complete - failure exists", the transaction is marked as failed, else the transaction is marked as "completed"; In both cases, a message is shown to the user.

The related status flags are set for the committed images in the local database.

For images sent in downsized format, the flag will be set for the original images still available on the system.

The AXIOM Artis DICOM application will NOT support the Storage Media File Set ID attributes.

8 Implementation Model Query/Retrieve

The query/retrieve service class defines an application-level class of services which facilitates the management of images and patient data against the well defined information model of DICOM and allows a DICOM AE to retrieve images from a remote DICOM node or to request a remote DICOM AE to initiate a transfer of images to another DICOM AE. The AXIOM Artis DICOM query/retrieve application supports the query/retrieve services to act as SCU and SCP.

8.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows XP application and acts as SCU and SCP for the query/retrieve network service.

Figure 5: Application Data Flow Diagram - Query/Retrieve SCU





Figure 6: Application Data Flow Diagram - Query/Retrieve SCP

8.2 Functional Definitions of Application Entities

The AXIOM Artis DICOM query/retrieve SCU requests the remote query/retrieve SCP to perform a search in its database (storage of SOP instances) and match to the keys specified in the request in order to display the results in the AXIOM Artis user interface. Depending on user action (Import) the AXIOM Artis DICOM SCU sends a C-MOVE DIMSE service to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query/Retrieve SCP) to the AXIOM Artis Storage SCP.

The AXIOM Artis DICOM query/retrieve SCP responds to C-FIND DIMSE services from remote SCU and depending on further remote request a C-GET or a C-MOVE involves the Siemens AXIOM Artis DICOM query/retrieve SCP application to initiate a C-STORE sub operation (by triggering and parametrizing the own Storage SCU) to send image objects to a remote Storage SCP.

All components of the Siemens SCP DICOM query/retrieve SCP application are operating as background server processes. They are existing when the machine is powered on and then respond to queries based on the records stored in its database.

8.3 Sequencing of real World Activities

Retrieve of images is only possible if results from a previous "Search..." operation exist and those entities can be selected for "Import".

9 AE Specification Query/Retrieve

9.1 Query/Retrieve Service AEs Specification

The Query/Retrieve SCU request that the remote SCP performs a match of all keys specified in the request, against the information in its database and the identified images will be moved over a different (C-MOVE) storage association.

The Query/Retrieve SCP responds to queries based on the records based on its database and images will be send to the requesting SCU or to a different storage destination.

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as SCU and SCP::

Table 12: Query/Retrieve SOP Classes as an SCU

SOP Class Name	SOP Class UID
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Classes as an SCP:

Table 13: Query/Retrieve SOP Classes as an SCP

SOP Class Name	SOP Class UID
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1
Patient Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.1.2
Patient Root Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.1.3
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1
Study Root Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.2.2
Study Root Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.2.3
Patient/Study Only Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.3.1
Patient/Study Only Query/Retrieve Information Model - MOVE	1.2.840.10008.5.1.4.1.2.3.2
Patient/Study Only Query/Retrieve Information Model - GET	1.2.840.10008.5.1.4.1.2.3.3

Note

See also the Storage DICOM Conformance Statement of the Siemens AXIOM Artis DICOM application to compare for conformance of the C-STORE sub-operation generated by the C-GET DIMSE service with the Storage Service SOP classes described in the Storage DICOM

Conformance Statement of the Modality to which the images shall be transferred to.

9.1.1 Association Establishment Policies

9.1.1.1 General

With the "Search..." function the query data are input and the DICOM query/retrieve application is started. An query request will be sent out to one or more remote nodes and the response data will be displayed for the user. Upon request (Import), the retrieval of selected items is initiated.

The default PDU size used will be 28 KB.

9.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM application initiates one association for each query request being processed to a remote node. The max. number of active associations is configured.

The Siemens AXIOM Artis DICOM application is able to accept multiple associations at a time. It can handle up to ten associations in parallel.

9.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

9.1.1.4 Implementation Identifying Information

The Siemens AXIOM Artis DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.5.9.20000101

and an Implementation Version Name of

- "SIEMENS_SWFVX47A"

9.1.2 Association Initiation Policy

The query user interface will request the query-data from the user and triggers one C-FIND request to the selected remote node. The response data will be displayed in the query UI for further data navigation.

When requesting Import of related items the browser requests the retrieve application to send a C-MOVE request to the related remote node. Images will then be received by the Storage SCP as described in the related section.

The query/retrieve operation can be canceled at any time by a respective C-FIND-CANCEL or C-MOVE-CANCEL.

9.1.2.1 Real World Activity - Find SCU

9.1.2.1.1 Associated Real-World Activity - Find SCU (Search...)

The associated Real-World activity is to fill out a query form with search data and start the search query. The DICOM application creates a query identifier and issues a C-FIND over a previously build association. The remote SCP will respond with related data-entries that will be passed to a browser application. When data transfer is finished the association is closed.

9.1.2.1.2 Proposed Presentation Contexts - Find SCU

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

It is configurable which of the two query models (or both) are to be used by the AXIOM Artis DICOM Query SCU application. If both Abstract Syntaxes are configured, The C-FIND SCU will use the Patient Root Model only for C-FIND requests on PATIENT level. For all other levels it will use the STUDY root model.

9.1.2.1.3 SOP Specific Conformance Statement - Find SCU

The Siemens AXIOM Artis DICOM Query/Retrieve SCU supports hierarchical queries with all mandatory search keys. The following table describes the search keys for the different query models that the SCU supports

Table 14: Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Patient Level^a					

Table 14: Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Patient Name	(0010,0010)	R	Wildcard ^b	enter value	yes
Patient ID	(0010,0020)	U / R	Wildcard ^b	enter value	yes
Patient's Birth date	(0010,0030)	O	universal (Null)	enter value	yes
Patient's Sex	(0010,0040)	O	universal (Null)	enter value	yes
Number of Patient related Studies	(0020,1200)	O	universal (Null)	-	yes ^c
Number of Patient related Series	(0020,1202)	O	universal (Null)	-	no
Number of Patient related Instances	(0020,1204)	O	universal (Null)	-	no
Study Level					
Patient Name	(0010,0010)	R	Wildcard ^b	enter value	yes
Patient ID	(0010,0020)	U / R	Wildcard ^b	enter value	yes
Patient's Birth date ^d	(0010,0030)	O	universal (Null)	enter value	yes
Patient's Sex ^d	(0010,0040)	O	universal (Null)	enter value	yes
Study Instance UID	(0020,000D)	U	single value	-	yes
Study ID	(0020,0010)	R	universal (Null)	enter value	yes
Study Date	(0008,0020)	R	universal (Null)	enter-value ^e	yes
Study Time	(0008,0030)	R	universal (Null)	-	yes
Accession Number	(0008,0050)	R	universal (Null)	-	yes
Study Description	(0008,1030)	O	universal (Null)	-	yes
Referring Physician's Name	(0008,0090)	O	universal (Null)	-	yes
Name of Physician Reading Study	(0008,1060)	O	universal (Null)	-	yes
Modalities in Study	(0008,0061)	O	universal (Null)	-	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	no
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	no
Number of Study related Series	(0020,1206)	O	universal (Null)	-	yes ^f
Number of Study related Instances	(0020,1208)	O	universal (Null)	-	no
Series Level					

Table 14: Patient Root and Study Root query attributes

Attribute name	Tag	Type	Matching	User input	return value display
Series Instance UID	(0020,000E)	U	single value	-	yes
Series Number	(0020,0011)	R	universal (Null)	-	yes
Modality	(0008,0060)	R	universal (Null)	enter value	yes
Series Date	(0008,0021)	O	universal (Null)	-	yes
Series Time	(0008,0031)	O	universal (Null)	-	yes
Series Description	(0008,103E)	O	universal (Null)	enter value	yes
Body Part Examined	(0018,0015)	O	universal (null)	enter value	yes
Performing Physician's Name	(0008,1050)	O	universal (Null)	enter value	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	yes
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	yes
Protocol Name	(0018,1030)	O	universal (Null)	-	yes
Perf. Procedure Step Start Date	(0040,0244)	O	universal (Null)	-	yes
Perf. Procedure Step Start Time	(0040,0245)	O	universal (Null)	-	yes
Request Attributes Sequence	(0040,0275)	O	universal (Null)	-	yes
> Requested Procedure ID	(0040,1001)	O	universal (Null)	-	yes
> Scheduled Procedure Step ID	(0040,0009)	O	universal (Null)	-	yes
Number of Series related Instances	(0020,1209)	O	universal (Null)	-	yes
Image Level					
SOP Instance UID	(0008,0018)	U	single value	-	no
Image Number	(0020,0013)	R	universal (Null)	-	yes
Storage Media File-Set ID	(0008,0130)	O	universal (Null)	-	no
Retrieve AE Title	(0008,0054)	O	universal (Null)	-	no
Instance Date	(0008,0023)	O	universal (Null)	-	no
Instance Time	(0008,0033)	O	universal (Null)	-	no
Number of Frames	(0028,0008)	O	universal (Null)	-	yes
Image Comments	(0020,4000)	O	universal (Null)	-	no

- a. Patient Root Information Model only
- b. Always a '*' is appended to the user-supplied string
- c. Implicitly visualized in the UI if no study and series search attributes have been entered
- d. Study Root Information Model only
- e. Date range also possible

f. Implicitely if no series search attributes have been entered

The Find SCU interprets following status codes:

Table 15: C-FIND response status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources	A700	(0000,0902)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	Cxxx	(0000,0901) (0000,0902)
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

9.1.2.2 Real World Activity - Move SCU

9.1.2.2.1 Associated Real-World Activity - Move SCU (Import...)

When selecting a data entry in the Query UI and activating the "Import" function, a retrieval request is passed to the DICOM application which issues a C-MOVE service according to the Patient Root or Study Root query model. (The Storage Service Class Conformance Statement describes the C-STORE service which is generated by the C-MOVE service.)

The transferred image data are processed as described in the storage class SCP descriptions.

The possibility to request the remote C-MOVE provider (remote application that responded to the C-FIND) to move data to an application entity other than the C-MOVE SCU (the AXIOM Artis DICOM application) is NOT USED.

9.1.2.2.2 Proposed Presentation Contexts - Move SCU

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table			
Abstract Syntax	Transfer Syntax	Role	Extended

Name	UID	Name List	UID List		Negotia- tion
Patient Root Query/Re- trieve Move	1.2.840.10008.5.1. 4.1.2.1.2	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None
Study Root Query/Re- trieve Move	1.2.840.10008.5.1. 4.1.2.2.2	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None

9.1.2.2.3 SOP Specific Conformance Statement - Move SCU

At association establishment time the C-MOVE presentation context shall be negotiated. The C-STORE sub-operations must be done on a different association to transfer images to the own Storage Service Class SCP.

The Move SCU interprets following status codes:

Table 16: C-MOVE response status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
Failed	Unable to process	Cxxx	(0000,0901) (0000,0902)
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures or Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

9.1.3 Association Acceptance Policy

The Siemens AXIOM Artis DICOM application will accept associations for the following DIMSE-C operations as SCP:

- C-FIND
- C-GET
- C-MOVE
- C-FIND-CANCEL
- C-GET-CANCEL
- C-MOVE-CANCEL

Extended negotiation - which is relational retrieve - is NOT supported for the above listed services. The AXIOM Artis DICOM application does support multiple C-FIND requests over the same association, while multiple C-MOVE or C-GET operations are not supported over the same association.

9.1.3.1 Real World Activity - Find SCP

9.1.3.1.1 Associated Real-World Activity - Find request to SCP

The associated Real-World activity is to respond query requests to an SCU with the query model Patient Root, Study Root and Patient/Study Only. Relational retrieve operation is NOT supported. With a C-FIND-CANCEL request the running query can be canceled at any time.

Multiple C-FIND requests over the same association are supported.

9.1.3.1.2 Proposed Presentation Contexts - Find SCP

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.1.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCP	None
Study Root Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.2.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCP	None
Patient/Study Only Query/Retrieve Find	1.2.840.10008.5.1.4.1.2.3.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCP	None

Note

C-FIND extended negotiation will not be supported by the SCP.

The order of preference for accepting Transfer Syntaxes is: 1. Explicit VR Little Endian, 2. Explicit VR Big Endian, 3. Implicit VR Little Endian.

9.1.3.1.3 SOP Specific Conformance Statement - Find SCP

The Siemens AXIOM Artis DICOM Query/Retrieve SCP supports hierarchical queries with all mandatory and optional search keys.

The query attribute contents will be treated case-sensitive.

With wildcard queries the symbol "?" is treated as "*" by the C-FIND SCP application. As a consequence the query string of "?abc*" will be processed as "*abc*".

If the value for the patient-level unique key "Patient ID" is not known, it may be returned with zero length. The attribute "Image Comments" will not be included in the C-FIND-RSP, if it is not set in the DB, even if it was requested as return key in the related C-FIND-RQ.

Usage of Storage Media File-Set ID, Retrieve AE Title with C-FIND-RSP message:

- The Storage Media File-Set ID - if existent - can be returned at Study/Series/Image Level. Only on Image Level, the values of ONLINE, NEARLINE or OFFLINE are returned to indicate the Storage Location of the related Instance.
- The Retrieve AE Title - if existent - can only be returned at Image Level (for Patient Root and Study Root models) or Study Level (for Patient/Study Only model).

Relational Queries are not supported.

A remote DICOM AE can cancel the running query by sending a C-FIND-CANCEL. Matches are possibly continuing (more C-FIND responses with status PENDING) until the cancel operation has completed.

The supported attributes on the various query levels of the three supported information models are listed in the tables of the following sections.

9.1.3.1.3.1 Patient Root Information Model**Table 17:** Patient level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal
Patient ID	(0010,0020)	U	single value, wildcard, universal
Patient's Birth Date	(0010,0030)	O	single value, range, universal
Patient's Birth Time	(0010,0032)	O	single value, range, universal
Patient's Sex	(0010,0040)	O	single value, wildcard, universal
Ethnic Group	(0010,2160)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal
Number of Patient related Instances	(0020,1204)	O	universal

Table 18: Study level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Study Instance UID	(0020,000D)	U	single value, list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal

Table 18: Study level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,1060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal
Number of Study related Instances	(0020,1208)	O	universal

Table 19: Series level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Series Instance UID	(0020,000E)	U	single value, list of UID
Series Number	(0020,0011)	R	single value, universal
Modality	(0008,0060)	R	single value, wildcard, universal
Laterality	(0020,0060)	O	single value, wildcard, universal
Body Part Examined	(0018,0015)	O	single value, wildcard, universal
Patient Position	(0018,5100)	O	single value, wildcard, universal
Smallest Pixel Value in Series	(0028,0108)	O	single value, universal
Largest Pixel Value in Series	(0028,0109)	O	single value, universal
Protocol Name	(0018,1030)	O	single value, wildcard, universal
Series Date	(0008,0021)	O	single value, range, universal
Series Time	(0008,0031)	O	single value, range, universal
Series Description	(0008,103E)	O	single value, wildcard, universal
Operators Name	(0008,1070)	O	single value, wildcard, universal
Performing Physician's Name	(0008,1050)	O	single value, wildcard, universal
Perf. Procedure Step Start Date	(0040,0244)	O	universal
Perf. Procedure Step Start Time	(0040,0245)	O	universal

Table 19: Series level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
Number of Series related Instances	(0020,1209)	O	universal

Table 20: Image level attributes, Patient Root Information Model

Attribute name	Tag	Type	Matching
SOP Instance UID	(0008,0018)	U	single value, list of UID
Image Number	(0020,0013)	R	single value, universal
Image Date	(0008,0023)	O	single value, range, universal
Image Time	(0008,0033)	O	single value, range, universal
Modality	(0008,0060)	O	single value, wildcard, universal
Image Comments	(0020,4000)	O	universal

9.1.3.1.3.2 Study Root Information Model

Table 21: Study level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal
Patient ID	(0010,0020)	R	single value, wildcard, universal
Patient's Birth Date	(0010,0030)	O	single value, range, universal
Patient's Birth Time	(0010,0032)	O	single value, range, universal
Patient's Sex	(0010,0040)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal
Number of Patient related Instances	(0020,1204)	O	universal
Study Instance UID	(0020,000D)	U	single value, list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal

Attribute name	Tag	Type	Matching
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,8060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal
Number of Study related Instances	(0020,1208)	O	universal

Table 22: Series level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
Series Instance UID	(0020,000E)	U	single value, list of UID
Series Number	(0020,0011)	R	single value, universal
Modality	(0008,0060)	R	single value, wildcard, universal
Laterality	(0020,0060)	O	single value, wildcard, universal
Body Part Examined	(0018,0015)	O	single value, wildcard, universal
Patient Position	(0018,5100)	O	single value, wildcard, universal
Smallest Pixel Value in Series	(0028,0108)	O	single value, universal
Largest Pixel Value in Series	(0028,0109)	O	single value, universal

Protocol Name	(0018,1030)	O	single value, wildcard, universal
Series Date	(0008,0021)	O	single value, range, universal
Series Time	(0008,0031)	O	single value, range, universal
Series Description	(0008,103E)	O	single value, wildcard, universal
Operators Name	(0008,1070)	O	single value, wildcard, universal
Performing Physician's Name	(0008,1050)	O	single value, wildcard, universal
Performed Procedure Step Start Date	(0040,0244)	O	universal
Performed Procedure Step Start Time	(0040,0245)	O	universal
Number of Series related Instances	(0020,1209)	O	universal

Table 23: Image level attributes, Study Root Information Model

Attribute name	Tag	Type	Matching
SOP Instance UID	(0008,0018)	U	single value,list of UID
Image Number	(0020,0013)	R	single value, universal
Image Date	(0008,0023)	O	single value, range, universal
Image Time	(0008,0033)	O	single value, range, universal
Modality	(0008,0060)	O	single value, wildcard, universal
Image Comments	(0020,4000)	O	universal

9.1.3.1.3.3 Patient Study Only Information models

Table 24: Patient instance level, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Patient Name	(0010,0010)	R	single value, wildcard, universal
Patient ID	(0010,0020)	U	single value, wildcard, universal
Patient's Birth Date	(0010,0030)	O	single value, range, universal
Patient's Birth Time	(0010,0032)	O	single value, range, universal

Attribute name	Tag	Type	Matching
Patient's Sex	(0010,0040)	O	single value, wildcard, universal
Ethnic Group	(0010,2160)	O	single value, wildcard, universal
Patient Comments	(0010,4000)	O	wildcard, universal
Number of Patient related Studies	(0020,1200)	O	universal
Number of Patient related Series	(0020,1202)	O	universal
Number of Patient related Instances	(0020,1204)	O	universal

Table 25: Study level attributes, Patient Study Only Information Model

Attribute name	Tag	Type	Matching
Study Instance UID	(0020,000D)	U	single value,list of UIDs
Study ID	(0020,0010)	R	single value, wildcard, universal
Study Date	(0008,0020)	R	single value, range, universal
Study Time	(0008,0030)	R	single value, range, universal
Accession Number	(0008,0050)	R	single value, wildcard, universal
Referring Physician's Name	(0008,0090)	O	single value, wildcard, universal
Study Description	(0008,1030)	O	single value, wildcard, universal
Admitting Diagnoses Description	(0008,1080)	O	single value, wildcard, universal
Patient's Age	(0010,1010)	O	single value, wildcard, universal
Patient's Size	(0010,1020)	O	single value, universal
Patient's Weight	(0010,1030)	O	single value, universal
Occupation	(0010,2180)	O	single value, wildcard, universal
Additional Patient History	(0010,21B0)	O	wildcard, universal
Name of Physician reading Study	(0008,8060)	O	single value, wildcard, universal
Modalities in Study	(0008,0061)	O	multiple values, universal
Number of Study related Series	(0020,1206)	O	universal
Number of Study related Instances	(0020,1208)	O	universal

The Find SCP returns following status codes:

Table 26: C-FIND SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources	A700	(0000,0902)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	C001	(0000,0901) (0000,0902)
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

9.1.3.2 Real World Activity - Get SCP

9.1.3.2.1 Associated Real-World Activity - Get request to SCP

The associated Real-World activity is to respond to retrieve requests initiated from a foreign SCU. The SCP supports the query model Patient Root, Study Root and Patient/Study Only. The Storage Service Class Conformance Statement describes the C-STORE service which is generated by the C-GET service. Relational retrieve operation is NOT supported.

Multiple C-GET requests over the same association are NOT supported.

9.1.3.2.2 Accepted Presentation Contexts - Get SCP

The Siemens AXIOM Artis DICOM application will accept Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		

Patient Root Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.1.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.2.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Patient/Study Only Query/Retrieve Get	1.2.840.10008.5.1.4.1.2.3.3	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

9.1.3.2.3 SOP Specific Conformance Statement - Get SCP

At association establishment time the C-GET presentation context must be negotiated along with the C-STORE sub-operations which must be accomplished on the same association as the C-GET operation. Relational retrieve operation is NOT supported.

The Get SCP returns following status codes:

Table 27: C-GET SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
Failed	Unable to process	C001	(0000,0901) (0000,0902)
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures or Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

Table 27: C-GET SCP return status

Service Status	Meaning	Protocol Codes	Related Fields
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)

9.1.3.3 Real World Activity - Move SCP

9.1.3.3.1 Associated Real-World Activity - Move request to SCP

The associated Real-World activity is to respond to retrieve requests to an SCU. The SCP supports the query model Patient Root, Study Root and Patient/Study Only. The Storage Service Class Conformance Statement describes the C-STORE service which is generated by the C-MOVE service. Relational retrieve operation is NOT supported.

Multiple C-MOVE requests over the same association are NOT supported.

9.1.3.3.2 Accepted Presentation Contexts - Move SCP

The Siemens AXIOM Artis DICOM application will accept Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Study Root Query/Retrieve Move	1.2.840.10008.5.1.4.1.2.2.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Patient/Study Only Query/ Retrieve Move	1.2.840.10008.5.1. 4.1.2.3.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

9.1.3.3.3 SOP Specific Conformance Statement - Move SCP

At association establishment time the C-MOVE presentation context shall be negotiated. The C-STORE sub-operations is done on a different association, specified in the C-MOVE request, to transfer images to a remote SCP of the Storage Service Class. Relational retrieve operation is NOT supported.

The Move SCP returns following status codes:

Table 28: C-MOVE return status

Service Status	Meaning	Protocol Codes	Related Fields
Refused	Out of Resources - Unable to calculate number of matches	A701	(0000,0902)
	Out of Resources - Unable to perform sub operations	A702	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
Failed	Unable to process	C001	(0000,0901) (0000,0902)
Cancel	Sub-operations terminated due to Cancel Indication	FE00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Warning	Sub-operations Complete - One or more Failures of Warnings	B000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Pending	Sub-operations are continuing	FF00	(0000,1020) (0000,1021) (0000,1022) (0000,1023)
Success	Sub-operations Complete - No Failures or Warning	0000	(0000,1020) (0000,1021) (0000,1022) (0000,1023)


10 Implementation Model Print

The Print Management Service Classes define an application-level class of services which facilitate the printing of images on a hardcopy medium. The print management SCU and print management SCP are peer DICOM print management application entities. The AXIOM Artis DICOM print application supports the print management DIMSE services to act as SCU.

10.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows XP application and acts as SCU for the print management network service.

Figure 7: Application Data Flow Diagram - Print SCU



Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

10.2 Functional Definitions of Application Entities

The Print SCU is invoked by the user interface to setup film-sheet layout and whenever an image is ready to be printed on film. The Print SCU will hold and maintain all data needed to compile a complete film-sheet from the data (images, layout, configuration) received. Whenever a film-sheet is ready to print, the related data is used to supply the Information to the SOP Classes of the Print Management Service Class. A queue is maintained in order to intermediately store several film-sheets in case of resource problems on printer. The SCU will only supply and require the mandatory SOP Classes of the Print Management Service Class.

10.3 Sequencing of real World Activities

Not applicable.

11 Application Entity Specification Print

11.1 Print Management AE Specification

The AXIOM Artis print management SCU (HCS) invokes print management DIMSE services to transfer images from the local AE to the remote SCP AE to print images with defined layout on a selected network-based DICOM hardcopy printer. This is done in a "full-page" print mode.

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 Print Management Meta SOP Classes as an SCU:

Table 29: Basic Grayscale Print Management Meta SOP Classes

SOP Class Name	SOP Class UID
Basic Greyscale Print Management MetaSOP Class	1.2.840.10008.5.1.1.9
- Basic Film Session SOP Class	1.2.840.10008.5.1.1.1
- Basic Film Box SOP Class	1.2.840.10008.5.1.1.2
- Basic Grayscale Image Box SOP Class	1.2.840.10008.5.1.1.4
- Printer SOP Class	1.2.840.10008.5.1.1.16
Print Job SOP Class	1.2.840.10008.5.1.1.14

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

11.1.1 Association Establishment Policies

11.1.1.1 General

Whenever a film is completely set up and printed by command or automatism, the job is prepared for processing. As soon as the queue is ready to process, the job is activated and processed according to the processing data. The related Print application will initiate an association to the print destination and process the printing of the related information. The association will be closed, if there is no further film to print within a fixed time-frame (1 min.). With no current association open, the Print manager will perform a cyclic check of the camera status (5 min. interval).

The default PDU size used will be 28 KB.

11.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM application initiates one association at a time for each different print device configured.

11.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

11.1.1.4 Implementation Identifying Information

The Siemens AXIOM Artis DICOM software provides a single Implementation Class UID of

- 1.2.12.2.1107.5.5.9.20000101

and an Implementation Version Name of

- "SIEMENS_SWFVX47A"

11.1.2 Association Initiation Policy

Triggered by the Print Job queue the Print Management SCU establishes an association by using the DICOM association services. With the help of the N-GET request for the Printer SOP Class the Status is determined before printing.

If no problem is encountered, then with the N-CREATE/N-SET Services for the related Basic Print SOP Classes the film sheet is set up for printing and the image(s) is(are) transferred to the printer device.

After the last film is printed from queue, the Print application will leave open the association for another 60 seconds. If a new film job is ready for printing within this time-limit, the job will be immediately processed over the still open association. If there is no new job, the association is closed if the time-out elapsed. This done to optimize automated printing.

During the "idle-time" (no open association to printer) the Print application will issue a cyclic camera status request (using N-GET of Printer SOP Class) every 5 minutes.

11.1.2.1 Real-World Activity - Print

11.1.2.1.1 Associated Real-World Activity - Printing a Printer Job Queue Entry

Whenever a film sheet is prepared by the user, it is forward to the Printer Job queue. As soon as the associated Printer device is available, the job is activated and association is set up.

The film sheet is internally processed, converted to a Standard/1-1 page and then the page image is sent. Status is controlled by awaiting any N-EVENT message all through the transfer until the last image or film sheet is sent.

If the response from the remote application contains a status other than Success or Warning, the association is aborted.

11.1.2.1.2 Proposed Presentation Contexts

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Basic Gray-scale Print Management Meta SOP class	1.2.840.10008.5.1.1.9	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic film session SOP class	1.2.840.10008.5.1.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic film box SOP class	1.2.840.10008.5.1.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic gray-scale image box SOP class	1.2.840.10008.5.1.1.4	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Printer SOP class	1.2.840.10008.5.1.1.16	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Print Job SOP class	1.2.840.10008.5.1.1.14	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

11.1.2.1.3 SOP Specific Conformance - Meta SOP Classes

The Siemens AXIOM Artis DICOM print management SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class.

The application uses a setting platform to define the properties of the connected DICOM SCP, e.g.:

- maximum number of print jobs in the queue
- maximum number of print copies

- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP
- lookup table definition.

The printing is only suspended in the case of a failure return status of the SCP.

11.1.2.1.3.1 Basic film session SOP class

The Basic Film Session information object definition describes all the user defined parameters which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes and that are printed on one hardcopy printer.

The AXIOM Artis DICOM print management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the following attributes:

Table 30: Used Basic Film Session N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Number of Copies	(2000,0010)	U	
Medium Type	(2000,0030)	U	BLUE FILM CLEAR FILM PAPER
Film Destination	(2000,0040)	U	MAGAZINE PROCESSOR

The number of copies sent to the DICOM printer is always 1, the job is sent n times for n copies.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session - see below:

Table 31: Attributes of N-DELETE-RQ on Basic Film Session SOP Class

Attribute Name	Tag	Source of Information
Requested SOP Instance UID	(000,1000) -> (0000,1001)	Affected SOP Instance UID of N-CREATE-RSP on Basic Film Session

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

The Basic Film Session SOP class interprets the following status codes (from N-CREATE-RSP, N-DELETE-RSP messages):

Table 32: Basic Film Session SOP status

Service Status	Meaning	Protocol Codes
Failure	Film session SOP instances hierarchy does not contain film box SOP instances	C600
	Unable to create print job, print queue is full	C601
	Image size is larger than images box size	C603
Warning	Memory allocation not supported	B600
	Film session printing is not supported	B601
	Film box does not contain image box (empty page)	B602
Success	Film belonging to the film session are accepted for printing	0000

11.1.2.1.3.2 Basic Film Box SOP class

The Basic Film Box information object definition describes all the user defined parameter of one film of the film session. The Basic Film Box information description defines the presentation parameters which are common for all images on a given sheet of film.

The Basic Film Box refers to one or more Image Boxes.

Supported Service Elements as SCU are:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP class N-CREATE-RQ message uses the following attributes (the actual values for each attribute depend on DICOM printer configuration within the AXIOM Artis DICOM print management SCU)

Table 33: Used Film Box N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Image Display Format	(2010,0010)	M	STANDARD\1,1
Referenced Film Session Sequence	(2010,0500)	M	
>Referenced SOP Class UID	(0008,1150)	M	1.2.840.10008.5.1.1.1
>Referenced SOP Instance UID	(0008,1155)	M	
Film Orientation	(2010,0040)	M	PORTRAIT

Table 33: Used Film Box N-CREATE-RQ attributes

Attribute name	Tag	Usage SCU	Supported Values
Film Size ID	(2010,0050)	M	8INx10IN, 10INx12IN, 10INx14IN, 11INx14IN, 14INx14IN, 14INx17IN, 24CMx24CM, 24CMx30CM
Magnification Type	(2010,0060)	M	BILINEAR, CUBIC, NONE, REPLICATE
Border Density	(2010,0100)	U	BLACK, WHITE
Max Density	(2010,0130)	U	>0
Min Density	(2010,0120)	U	50>value>0
Illumination	(2010,015E)	U	>0 Required if Presentation LUT is present
Reflective Ambient Light	(2010,0160)	U	>0 Required if Presentation LUT is present
Referenced Presentation LUT Sequence	(2050,0500)	U	

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally to be further used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the AXIOM Artis DICOM print manager will issue a N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Box - see below:

Table 34: Attributes of N-DELETE-RQ on Basic Film Box SOP Class

Attribute Name	Tag	Source of Information
Requested SOP Instance UID	(0000,1000) -> (0000,1001)	Affected SOP Instance UID of N-CREATE-RSP on Basic Film Box

The Basic Film Box SOP class interprets the following status codes:

Tabelle 35: Basic Film Box SOP status

Service Status	Meaning	Protocol Codes
Failure	Unable to create print job; print queue is full	C602
	Image size is larger than image box size	C603

Table 35: Basic Film Box SOP status

Service Status	Meaning	Protocol Codes
Warning	Film box does not contain image box (empty page)	B603
	Requested MinDensity or MaxDensity outside of Printer's operating range	B605
Success	Film accepted for printing	0000

11.1.2.1.3.3 Basic Greyscale Image Box SOP Class

The Basic Greyscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the following attributes:

Table 36: Mandatory Basic Grayscale Image Box N-SET attributes

Attribute name	Tag	Usage SCU	Supported Values
Image Position	(2020,0010)	M	1
Basic Grayscale Image Sequence	(2020,0110)	M	
>Samples Per Pixel	(0028,0002)	M	1
>Photometric Interpretation	(0028,0004)	M	MONOCHROME2
>Rows	(0028,0010)	M	
>Columns	(0028,0011)	M	
>Pixel Aspect Ratio	(0028,0034)	M	
>Bits Allocated	(0028,0100)	M	8
>Bits Stored	(0028,0101)	M	8
>High Bit	(0028,0102)	M	7
>Pixel Representation	(0028,0103)	M	0
>Pixel Data	(7FE0,0010)	M	

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

The Grayscale Image Box SOP class interpret the following status codes:

Table 37: Basic Grayscale Image Box SOP status

Service Status	Meaning	Protocol Codes
Warning	Requested MinDensity or MaxDensity outside of Printer's operating range	B605
Failure	Image contains more pixel than printer can print in Image Box	C603
	Insufficient memory in printer to store the image	C605

Table 37: Basic Grayscale Image Box SOP status

Service Status	Meaning	Protocol Codes
Success		0000

11.1.2.1.3.4 Printer SOP Class

The Printer SOP Class is the possibility to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The SCU uses the mandatory N-EVENT Report DIMSE service to monitor the changes of the printer status in an asynchronous way.

It can directly ask the Printer (SCP) for its status or receive Events from the Printer asynchronously:

- N-GET as SCU
- N-EVENT-REPORT as SCU

In both cases the following information is supported:

Table 38: Used Printer N-EVENT report attributes

Event type name	Event	Attributes	Tag	Usage SCU
Normal	1			
Warning	2	Printer Status Info	(2110,0020)	U
Failure	3	Printer Status Info	(2110,0020)	U

Table 39: Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes

Attribute name	Tag	Usage SCP	supported values
Printer Status	(2110,0010)	M	NORMAL, FAILURE, WARNING
Printer Status Info	(2110,0020)	M	see tables in the Annex for details

Note

For a detailed description on how AXIOM Artis reacts on different printer status messages, please refer to the annex section "DICOM Print SCU - detailed status displays".

11.1.2.1.3.5 Print Job SOP Class

The optional Print Job SOP Class is the possibility to monitor the execution of the print process.

The AXIOM Artis DICOM Print Management application supports the mandatory N-EVENT-REPORT DIMSE Service of this class to receive the changes of the Print Job Status in an asynchronous way.

It can receive Events from the Print SCP asynchronously:

- N-EVENT-REPORT

Note

The AXIOM Artis does not support receiving N-EVENT-REPORT from a camera during a print session. To support this behavior, in many cases this can be configured in the camera.

The following information is supported:

Table 40: Used Print Job N-EVENT report attributes

Event type name	Event	Attributes	Tag	Usage SCU
Normal	1	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U
Printing	2	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U
Done	3	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U
Failure	4	Execution Status Info	(2100,0030)	U
		Print Job ID	(2100,0010)	- (Print Queue Management SOP Class not supported)
		Film Session Label	(2000,0050)	U
		Printer Name	(2110,0030)	U


12 Implementation Model Worklist

The Basic Worklist Management Service class defines an application-level class of service which facilitates the transfer of worklists from the information system to the imaging modality. The worklist is queried by the AE and supplies the SCU with the scheduled tasks which have to be performed on the modality. The AXIOM Artis DICOM worklist application supports the worklist service as SCU.

12.1 Application Data Flow Diagram

The AXIOM Artis DICOM network implementation is a Windows NT application and acts as SCU for the Basic Worklist Service using the Modality Worklist SOP Class.

Figure 8: Application data flow diagram - Basic Worklist SCU:



Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

12.2 Functional Definitions of Application Entities

The worklist SCU ("broad query") is invoked from the patient browser user interface or by timer to request the worklist from an remote Information System (Modality Worklist Class SCP) to perform a match to the internal worklist query keys specified in the C-Find DIMSE service issued for the Modality Worklist Model.

The worklist SCP responds to the C-FIND query and scheduled imaging service requests (scheduled procedure steps) and patient demographic information will be downloaded from the information system to the AXIOM Artis modality. All information retrieved will be hold in the scheduling database for usage during Patient registration procedure.

Furthermore, during Patient Registration dialog, it is possible to update/complete the entered data by a "Get Worklist" function. Some of the entered data will then be used as matching criteria ("narrow query") for the issue worklist query. With the response data then the Patient Registration dialog fields are populated according availability within the worklist response identifier.

12.3 Sequencing of real World Activities

The "narrow" (interactive) Worklist Query requires that sufficient matching keys or a unique matching key are/is entered before the query is issued. Only then a single response can be expected to complete the registration dialog.

13 AE Specification Worklist

13.1 Modality Worklist Service AEs Specification

The Modality worklist SCU (patient registration in conjunction with the network application) requests that the remote SCP performs a match of all keys specified in the query against the information in its worklist database.

SIEMENS AXIOM Artis DICOM products provide Standard Conformance to the following DICOM V3.0 SOP Class as an SCU:

Table 41: SOP Classes as an SCU

SOP Class Name	SOP Class UID
Modality Worklist Information Model - FIND	1.2.840.10008.5.1.4.31

13.1.1 Association Establishment Policies

13.1.1.1 General

There will be a cyclic update of the modality scheduler database by background worklist request with date/time and modality information.

In addition the user can request worklist update with "Update Worklist". The modality will not insert duplicate entries identified by unique identifiers.

An interactive worklist query can be issued with search criteria entered during patient registration.

The default PDU size used will be 28 KB.

13.1.1.2 Number of Associations

The Siemens AXIOM Artis DICOM application initiates one association at a time to query worklist entry data.

13.1.1.3 Asynchronous Nature

The Siemens AXIOM Artis DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

13.1.1.4 Implementation Identifying Information

The Siemens AXIOM Artis DICOM software provides a single Implementation Class UID of

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- "SIEMENS_SWFVX47A"

13.1.2 Association Initiation Policy

The network application will cyclically query the worklist and by request of patient registration interface. Ever then it establishes an association by using the DICOM association services. During association establishment the negotiation of SOP classes to exchange the capabilities of the SCU and the SCP is not supported.

The following DIMSE-C operation is supported as SCU:

- C-FIND

13.1.2.1 Real World Activity

13.1.2.1.1 Associated Real-World Activity - Update (Query) Worklist

A network application will perform worklist queries with the C-FIND request at regular intervals. In addition it can be triggered by immediate request. The received worklist items will be compared with the contents of the local scheduler database. New items will be inserted into scheduler database.

13.1.2.1.2 Proposed Presentation Contexts - Update Worklist

The Siemens AXIOM Artis DICOM application will propose Presentation Contexts as shown in the following table:

Table 42: *Proposed presentation contexts*

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Extended Negotiation
Name	UID	Name List	UID List		
Modality Worklist Information Model- FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

13.1.2.1.3 SOP Specific Conformance Statement- Update Worklist

Search Key Attributes of the Worklist C-FIND

The Siemens AXIOM Artis DICOM worklist SCU supports "broad worklist queries" with all required search keys. The following tables describe the "broad query" search keys that the SCU supports.

Table 43: Supported Broad Worklist Query Search Key Attributes

Attribute name	Tag	Matching Key Type	Query Value
Scheduled Procedure Step			
Scheduled Procedure Step Sequence	(0040,0100)	R	
>Scheduled Station AE Title	(0040,0001)	R	<own AET> or "*" ^a
>Scheduled Procedure Step Start Date	(0040,0002)	R	<act. Date>-<act. Date> or range from UI ^b
>Scheduled Procedure Step Start Time	(0040,0003)	R	<zero length> or range from UI ^b
>Modality	(0008,0060)	R	"*" or <own Modality> ^a
>Scheduled Performing Physician's Name	(0040,0006)	R	NULL

a. This depends on user configuration (Options->Configuration->Patient Registration) if the "own AET" or "Modality" is provided. Use the "Search" tabcard and (un-)check the "Worklist for Local/Modality" item.

b. It depends on user configuration (Options->Configuration->Patient Registration) if the actual Date with a full time range or an interactive input dialog for date/time specification is used.

Return Key Attributes of the Worklist C-FIND-RSP

The Siemens AXIOM Artis DICOM worklist SCU supports worklist queries with return key attributes of all types. The following tables describe the return keys that the SCU supports.

An "x" in the **UI** column will indicate that the attribute is visualized when browsing the Worklist results with Patient Browser and/or during Patient Registration. The Patient Browser display is additionally influenced by the related Browser configuration.

An "x" in the **IOD** column will indicate that the related attribute is included into the SOP Instances of the IOD's created during processing of this worklist request.

Table 44: Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
SOP Common				
Specific Character Set	(0008,0005)	1C	-	x
Scheduled Procedure Step				
Scheduled Procedure Step Sequence	(0040,0100)	1		
>Scheduled Station AE Title	(0040,0001)	1	x	
>Scheduled Procedure Step Start Date	(0040,0002)	1	x	
>Scheduled Procedure Step Start Time	(0040,0003)	1	x	

Table 44: Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
>Scheduled Procedure Step End Date	(0040,0004)	3	-	
>Scheduled Procedure Step End Time	(0040,0005)	3	-	
>Modality	(0008,0060)	1	x	
>Scheduled Performing Physician's Name	(0040,0006)	1	x	x ^a
>Scheduled Procedure Step Description	(0040,0007)	1C	x	x
>Scheduled Station Name	(0040,0010)	2	x	
>Scheduled Procedure Step Location	(0040,0011)	2	x	
>Scheduled Protocol Code Sequence	(0040,0008)	1C	-	
>>Code Value	(0008,0100)	1C	x	
>>Coding Scheme Designator	(0008,0102)	1C	x	
>> Coding Scheme Version	(0008,0103)	3	x	
>>Code Meaning	(0008,0104)	3	x	
>Pre-Medication	(0040,0012)	2C	x	
>Scheduled Procedure Step ID	(0040,0009)	1	x	x
>Requested Contrast Agent	(0032,1070)	2C	x	x
> Scheduled Procedure Step Status	(0040,0020)	3	x	
>Comments on the Scheduled Procedure Step	(0040,0400)	3	-	
Requested Procedure				
Requested Procedure ID	(0040,1001)	1	x	x
Requested Procedure Description	(0032,1060)	1C	x	x
Requested Procedure Code Sequence	(0032,1064)	1C	-	
>Code Value	(0008,0100)	1C	x	
>Code Scheme Designator	(0008,0102)	1C	x	
> Coding Scheme Version	(0008,0103)	3	x	
>Code Meaning	(0008,0104)	3	x	
Study Instance UID	(0020,000D)	1	-	x
Referenced Study Sequence	(0008,1110)	2	-	
>Referenced SOP Class UID	(0008,1150)	1C	-	
>Referenced SOP Instance UID	(0008,1155)	1C	-	
Requested Procedure Priority	(0040,1003)	2	x	
Patient Transport Arrangements	(0040,1004)	2	-	
Reason for the Requested Procedure	(0040,1002)	3	-	

Table 44: Basic Worksit C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
Confidentiality Code	(0040,1008)	3	-	
Reporting Priority	(0040,1009)	3	-	
Names of intended Recipients of Results	(0040,1010)	3	-	
Requested Procedure Comments	(0040,1400)	3	x	
Requested Procedure Location	(0040,1005)	3	-	
Imaging Service Request				
Accession Number	(0008,0050)	2	x	x
Requesting Physician	(0032,1032)	2	x	x
Referring Physician's Name	(0008,0090)	2	x	x
Reason for Imaging Service Request	(0040,2001)	3	-	
Imaging Service Request Comments	(0040,2400)	3	x	
Requesting Service	(0032,1033)	3	x	
Issuing Date of Imaging Service Request	(0040,2004)	3	-	
Issuing Time of Imaging Service Request	(0040,2005)	3	-	
Placer Order Number / Imaging Service Request	(0040,2016)	3	-	
Filler Order Number / Imaging Service Request	(0040,2017)	3	-	
Order entered by ...	(0040,2008)	3	-	
Order Enterer's location	(0040,2009)	3	-	
Order Callback Number	(0040,2010)	3	-	
Visit Identification				
Admission ID	(0038,0010)	2	x	
Issuer of Admission ID	(0038,0011)	3	-	
Institution Name	(0008,0080)	3	x	x
Institution Address	(0008,0081)	3		x
Visit Status				
Current Patient Location	(0038,0300)	2	x	
Visit Relationship				
Referenced Patient Sequence	(0008,1120)	2	-	
>Referenced SOP Class UID	(0008,1150)	2	-	
>Referenced SOP Instance UID	(0008,1155)	2	-	
Visit Admission				
Admitting Diagnoses Description	(0008,1080)	3	x	

Table 44: Basic Worklist C-FIND-RSP Return Key Attributes

Attribute name	Tag	Return Key Type	UI	IOD
Patient Identification				
Patient's Name	(0010,0010)	1	x	x
Patient ID	(0010,0020)	1	x	x
Other Patient IDs	(0010,1000)	3	x	x
Patient Patient Names	(0010,1001)	3	x	x
Patient Demographic				
Patients Birth Date	(0010,0030)	2	x	x
Patient's Sex	(0010,0040)	2	x	x
Patient's Size	(0010,1020)	2	x	x
Patient's Weight	(0010,1030)	2	x	x
Confidential constraint on patient data	(0040,3001)	2	x	x
Patient's Address	(0010,1040)	3	x	x
Military Rank	(0010,1080)	3	x	x
Ethnic Group	(0010,2160)	3	x	x
Patient Comment	(0010,4000)	3	x	x
Patient Medical				
Patient State	(0038,0500)	2	x	x
Pregnancy Status	(0010,21C0)	2	x	x
Medical Alerts	(0010,2000)	2	x	x
Contrast Allergies	(0010,2110)	2	x	x
Special Needs	(0038,0050)	2	x	x
Smoking Status	(0010,21A0)	3	x	x
Last Menstrual Date	(0010,21D0)	3	x	x
Additional Patient History	(0010,21B0)	3	x	x

a. "Scheduled Performing Physician's Name" is taken as default for "Performing Physicians's Name"

13.1.2.1.4 Associated Real-World Activity - Get Worklist

With "Get Worklist" the contents of certain input fields of the Patient Registration UI are used to form a worklist request identifier. With the response data the Patient Registration dialog input is completed. The response data are additionally placed in the scheduler database.

13.1.2.1.5 Proposed Presentation Contexts - Get Worklist

This RWE will propose the same Presentation Contexts as with "Update Worklist". Please see table in section 13.1.2.1.2.

13.1.2.1.6 SOP Specific Conformance Statement - Get Worklist

Search Key Attributes of the Worklist C-FIND

The Siemens AXIOM Artis DICOM worklist SCU supports "narrow worklist queries" with all required search keys. The following tables describe the "narrow query" search keys that the SCU supports.

Table 45: Patient based "narrow query" Search Key Attributes

Attribute name	Tags	Matching Key Type	Query Value
Scheduled Procedure Step			
Scheduled Procedure Step Sequence	[0040,0100]	R	
>Scheduled Performing Physician's Name	[0040,0006]	R	input from UI or <zero length>
Requested Procedure			
Requested Procedure ID	[0040,1001]	O	input from UI or <zero length>
Imaging Service Request			
Accession Number	[0008,0050]	O	input from UI or <zero length>
Referring Physician's Name	[0008,0090]	O	input from UI or <zero length>
Visit Status			
Current Patient Location	[0038,0300]	O	input from UI or <zero length>
Patient Identification			
Patient's Name	(0010,0010)	R	input from UI or <zero length>
Patient ID	(0010,0020)	R	input from UI or <zero length>

Return Key Attributes of the Worklist C-FIND

Please see list for "Update Worklist" RWE.

Status Codes of the Worklist C-FIND

The worklist SCU interprets following status codes:

Table 46: C-FIND Response Status

Service Status	Meaning	Status Codes (0000,0900)	Related Fields
Refused	Out of Resources	A700	(0000,0902)
Failed	Identifier does not match SOP Class	A900	(0000,0901) (0000,0902)
	Unable to process	Cxxx	(0000,0901) (0000,0902)
Cancel	Matching terminated due to Cancel request	FE00	None
Success	Matching is complete - No final Identifier is supplied	0000	None
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this identifier	FF01	Identifier

14 Communication Profiles

14.1 Supported Communication Stacks

The Siemens AXIOM Artis DICOM application provide DICOM TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard.

14.1.1 TCP/IP Stack

The Siemens AXIOM Artis DICOM application uses the TCP/IP stack from the Windows NT system upon which it executes. It uses the MergeCOM-3 subroutine library from Merge Technologies Inc. that is based on a Berkeley socket interface.

14.1.1.1 API

The Siemens AXIOM Artis DICOM application uses the MergeCOM library that is based on a TCP/IP socket interface.

14.1.1.2 Physical Media Support

The Siemens AXIOM Artis DICOM application is indifferent to the physical medium over which TCP/IP executes; it inherits this from the Windows XP system upon which it executes.

15 Extensions/Specializations/Privatizations

15.1 Standard Extended/Specialized/Private SOPs

Please refer to Appendix A for all informations on these topics. A detailed overview is given there.

15.2 Private Transfer Syntaxes

Not applicable.

16 Configuration

16.1 AE Title / Presentation Address Mapping

To ensure unique identification within the network the hostname should be used as part of the AE Titles (see examples below, hostname = masca1). The string can be up to 16 characters long and must not contain any extended characters, only 7-bit ASCII characters (excluding Control Characters) are allowed according to DICOM Standard.

Note

The current implementation of AXIOM Artis does not support "Spaces" and special characters (e.g. & < > ") in the AE title string.

16.1.1 DICOM Verification

The Verification Service uses the AE configuration of the DICOM Service that is checked with the C-ECHO message. e.g. Verification will use the Storage AE, if initiated to check the configuration of a remote DICOM node.

16.1.2 DICOM Storage AE Title

The DICOM Storage application provides the application entity title which can be configured via Service UI in "Configuration/DICOM/General":

e.g. STU_MASCA1

The port number is set to the fixed value of

104

16.1.3 DICOM Query/Retrieve AE Title

The DICOM Query/Retrieve application uses the same application entity title as the DICOM Storage AE.

16.1.4 DICOM Print AE Title

The DICOM Print application provides the application entity title which can be configured via Service UI in "Configuration/DICOM/General":

e.g. PRI_MASCA1

The port number is set to the fixed value of

108

16.1.5 DICOM Modality Worklist AE Title

The DICOM Modality Worklist application provides the application entity title which can be configured via Service UI in "Configuration/DICOM/General" :

e.g. HRI_MASCA1

16.2 Configurable Parameters

The Application Entity Titles, host names and port numbers are configured using the Siemens AXIOM Artis Service/Installation Tool. For each AET the list of services supported can be configured.

16.2.1 Storage and Query/Retrieve

The Siemens AXIOM Artis Service/Installation Tool can be used to set the AET's, port-numbers, host-names, IP-addresses and capabilities for the remote nodes (SCP's). The service user can select transfer syntaxes, compression modes and query models for each SCP separately.

The sending of the private SOP Class (Exam Report) can be disabled.

For sending, it can be configured for each destination if downsizing (512, 8 bit) shall be used. Downsizing can be configured depending on image type.

Additional configurable parameters for Storage Commitment are:

When acting as SCU:

- flag to indicate whether the association will be kept open to receive the response or to close the association and be prepared to receive the response on another association.
- time-out which defines how long the association of N-ACTION is kept to receive a N-EVENT-REPORT on the same association. The same value is used to wait for a N-EVENT-REPORT on another association. (default 1h)

When acting as SCP:

- flag to indicate if an archive system is installed

16.2.2 Print

The Siemens AXIOM Artis Service/Installation Tool can be used to configure the SCP (DICOM printer). The AET, host-name, IP-address and port-number can be set.

These parameters have defaults as per configuration and can be modified:

Default camera (yes/no), pixel size, additional or changed film-sheet formats (e.g. 14x14inch, 14x17inch, ...), list with mapping pixel size to each film-sheet format, minimal density, stored printed film jobs, media type, film destination.

16.2.3 Worklist

The Siemens AXIOM Artis Service/Installation Tool can be used to configure the SCP (Worklist Provider). The AET, host-name, IP-address, port-number and time-outs can be set.

Additional configurable parameters for Basic Worklist Service are:

- *Query Waiting Time* - the time to wait for the C-FIND-RSP after sending the C-FIND-RQ (default 20s)
- *Max Query Match Number* - the maximum number of entries accepted as worklist response (default is 100)
- *Query Interval* - the time between two cyclic C-FIND-RQ to the HIS (default is 60 min)
- Broad Worklist Query behaviour: two values are defined:
 - Set the AE Title search attribute to the own AE Title, and the Modality search attribute to "*".
 - Set the Modality search attribute to the own modality and the AE Title search attribute to "*".

16.3 Default Parameters

The installation tool also uses some default parameters:

- max PDU size set to 28672 Bytes (28 kB)
- time-out for accepting/rejecting an association request: 240 s
- time-out for responding to an association open/close request: 240 s
- time-out for accepting a message over network: 240 s
- time-out for waiting for data between TCP/IP-packets: 240 s

The time-outs for waiting for a Service Request/Response message from the remote node are as follows:

- for Storage SCP/SCU: 600 s
- for Storage Commitment SCU:
time-out for response to N-ACTION: 1h
- for Query/Retrieve SCP/SCU: 600 s
- for Basic Worklist SCU; configurable, see section 16.2 on page 84
- for Print Management SCU:
 - time-out for Response to N-SET-RQ: 60 s
 - time-out for Response to other Requests: 60 s

17 Support of Extended Character Sets

The Siemens AXIOM Artis DICOM application supports the ISO 8859 Latin 1 (ISO-IR 100) character set.

Also the Japanese language character sets JIS X 0201 (ISO-IR 13 Japanese katakana and ISO-IR 14 Japanese romaji), JIS X 0208 (ISO-IR 87 Japanese kanji) and JIS X 0212 (IOS-IR 159 Supplementary Japanese kanji) are supported.

When there is a mismatch between the SCS tags (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM con-form:

Try to import with ISO_IR 100. If ISO_IR 100 fails, convert each illegal character to a '?'.

M.A.S. Offline Media Application Profile

Contains a AXIOM Artis specific Application Profile since DICOM does not provide a "General Angio/X-Ray Application Profile" for Storage on "Offline Media".

Structure of this Application Profile is defined in Part 11 of the DICOM Standard.

It is needed to describe the requirements for Offline Media Storage of the private IOD (MAS Non-Image IOD).

The profile is shared with LEONARDO Workstation and therefore support of MOD Media is included, even if not used with AXIOM Artis based systems.

1 Class and Profile Identification

This document defines an Application Profile Class for Universal X-Ray Modular Angiographic / Radiographic applications.

The identifier for this class shall be PRI-XAMAS. This class is intended to be used for interchange of extended and private Information Objects via CD-R or re-writable magneto-optical disk (MOD) offline media between dedicated Angiography, Neuro-Radiology, Radiography applications build from a common architecture.

The specific application profiles in this class are shown in Table 1:

Table 1: PRI-XAMAS Universal Angio/Radio Class Profiles

Application Profile	Identifier	Description
Universal X-Ray Modular Angio/ Radio System on CD-R	PRI-XAMAS-CD	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Dose Reports, Object Graphics, Quantification Results and more).
Universal X-Ray Modular Angio/ Radio System on 2.3 GB MOD	PRI-XAMAS-MOD23	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Dose Reports, Object Graphics, Quantification Results and more).
Universal X-Ray Modular Angio/ Radio System on 4.1 GB MOD ^a	PRI-XAMAS-MOD41	Handles interchange of Composite SOP Instances and privately defined SOP Instances (Dose Reports, Object Graphics, Quantification Results and more).

a. Definition of this profile is based on DICOM Part12 changes introduced by accepted Supplement 44.

Equipment claiming conformance for this Universal X-Ray Modular Angio/Radio System shall make a clear statement on handling of the private defined SOP Instances.

2 Clinical Context

This application profile facilitates the interchange of X-Ray acquisition and derived images and private data related to them. Typical media interchange would be from in-lab acquisition equipment to dedicated angio workstations and archive systems with specific extensions to handle the private data objects (in both directions).

Additionally, images (from MR, CT, US) used to prepare interventional procedures, multi-modality images (e.g. integrated US) and images derived from primary X-Ray images, such as annotations, quantitative analysis images, reference images, screen capture images may be interchanged via this profile.

2.1 Roles and Service Class Options

This Application Profile uses the Media Storage Service Class defined in PS 3.4 with the Interchange Option.

The Application Entity shall support one or more of the roles of File Set Creator (FSC), File Set Reader (FSR), and File Set Updater (FSU), defined in PS 3.10.

2.1.1 File Set Creator

The Application Entity acting as a File-Set Creator generates a File Set under the PRI-XAMAS Application Profiles.

File Set Creators shall be able to generate the Basic Directory SOP Class in the DICOMDIR file with all the subsidiary Directory Records related to the Image SOP Classes and Private SOP Classes stored in the File Set.

In case of PRI-XAMAS-CD Profile, the FSC shall offer the ability to either finalize the disc at the completion of the most recent write session (no additional information can be subsequently added to the disc) or to allow multi-session (additional information may be subsequently added to the disc).

Note

A multiple volume (a logical volume that can cross multiple physical media) is not supported by this Application Profile Class. If a set of Files, e.g., a Study, cannot be written entirely on one CD-R, the FSC will create multiple independent DICOM File-Set such that each File-Set can reside on a single CD-R medium controlled by its individual DICOMDIR file. The user of the FSC can opt to use written labels on the discs to reflect that there is more than one disc for this set of files (e.g., a Study).

2.1.2 File Set Reader

The role of the File Set Reader shall be used by Application Entities which receive the transferred File Set.

File Set Readers shall be able to read all the defined SOP Instances files defined for the specific Application Profiles to which a conformance claim is made, using all the defined Transfer Syntaxes.

2.1.3 File Set Updater

The role of the File Set Updater shall be used by Application Entities which receive a transferred File Set and update it by the addition of processed information.

File Set Updaters shall be able to read and update the DICOMDIR file. File-Set Updaters do not have to read the image/private information objects. File-Set Updaters shall be able to generate any of the SOP Instances files defined for the specific Application Profiles to which a conformance claim is made, and to read and update the DICOMDIR file.

In case of PRI-XAMAS-CD Profile, the FSU shall offer the ability to either finalize a disc at the completion of the most recent write session (no additional information can be subsequently added to the disc) or to allow multi-session (additional information may be subsequently added to the disc).

Note (for CD-R)

If the disc has not been finalized, the File-Set Updater will be able to update information assuming there is enough space on the disc to write a new DICOMDIR file, the information, and the fundamental CD-R control structures. CD-R control structures are the structures that inherent to the CD-R standards; see PS 3.12

3 PRI-XAMAS Profiles

3.1 SOP Classes and Transfer Syntaxes

This Application Profiles are based on the Media Storage Service Class with the Interchange Option.

Table 2: PRI-XAMAS SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FSR	FSU
Basic Directory	1.2.840.10008.1.3.10	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	M
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	M	M	O
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	O	M	O
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	M	M	O
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	O	M	O
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
X-Ray Radiofluoroscopic Image	1.2.840.10008.5.1.4.1.1.12.2	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O

Table 2: PRI-XAMAS SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FSR	FSU
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
CR Image	1.2.840.10008.5.1.4.1.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
CR Image	1.2.840.10008.5.1.4.1.1.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
CR Image	1.2.840.10008.5.1.4.1.1.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
CR Image	1.2.840.10008.5.1.4.1.1.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
CT Image	1.2.840.10008.5.1.4.1.1.2	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
CT Image	1.2.840.10008.5.1.4.1.1.2	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
MR Image	1.2.840.10008.5.1.4.1.1.4	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
MR Image	1.2.840.10008.5.1.4.1.1.4	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

Table 2: PRI-XAMAS SOP Classes and Transfer Syntaxes

Information Object Definition	Service Object Pair Class UID	Transfer Syntax UID	FSC	FSR	FSU
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
Ultrasound Image	1.2.840.10008.5.1.4.1.1.6.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
Ultrasound Multiframe Image	1.2.840.10008.5.1.4.1.1.3.1	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
NM Image	1.2.840.10008.5.1.4.1.1.20	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	M	M	O
NM Image	1.2.840.10008.5.1.4.1.1.20	JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70	O	M	O
NM Image	1.2.840.10008.5.1.4.1.1.20	Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2	O	M	O
NM Image	1.2.840.10008.5.1.4.1.1.20	JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51	O	O	O
CSA Non-Image	1.3.12.2.1107.5.9.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	O	M	O
Detached Patient Management	1.2.840.10008.3.1.2.1.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1	O	see DCS	O

FSC, FSR, FSU - denote the Requirements for those roles.

O - Optional

M - Mandatory

3.2 Physical Media and Media Formats

The PRI-XAMAS-CD Profile requires the 120mm CD-R physical media with the ISO/IEC 9660 Media Format, as defined in PS3.12.

The PRI-XAMAS-MOD23 Profile requires the 130mm 2.3 GB R/W MOD physical medium with the PC DOS Media Format, as defined in PS3.12.

The PRI-XAMAS-MOD41 Profile requires the 130mm 4.1 GB R/W MOD physical medium with the PC DOS Media Format, as defined in PS3.12.

3.3 Directory Information in DICOMDIR

Conforming Application Entities shall include in the DICOMDIR File the Basic Directory IOD containing Directory Records at the Patient and subsidiary levels appropriate to the SOP Classes in the File-set. All DICOM files in the File-set incorporating SOP Instances defined for the specific Application profile shall be referenced by Directory Records.

Note

DICOMDIRs with no directory information are not allowed by this Application Profile
Privately defined IODs will be referenced by "PRIVATE" Directory Records.

3.3.1 Basic Directory IOD Specialization

This Application Profile makes use of optional attributes of the Basic Directory IOD to support recognition of Patient's Storage Service request results in spanning multiple volumes (file sets). Therefore the File Set Descriptor File can be used and is then referenced by optional Basic Directory IOD attributes. If existent, the specified Descriptor File may be used by FSR applications. Any FSU, FSC shall make a clear Statement if the Descriptor File mechanism is used according to the specialization defined in this Application Profile.

The Descriptor Files shall have the following contents:

One single Line without any control-characters and according to the Basic Character-Set having the following defined text:

"MULTIVOLUME: xx of yy"

xx, yy are replaced by the actual Number of the volume (xx) and the Total Number of Volumes in the set (yy).

If used, the Descriptor File shall have the File ID "README" and reside in same directory level as the DicomDIR. It is referenced by the attribute [0004,1141] File-set Descriptor File ID having the defined content of "README".

3.3.2 Additional Keys

File-set Creators and Updaters are required to generate the mandatory elements specified in PS 3.3, Annex F of the DICOM Standard. Table 3: *PRI-XAMAS Additional DicomDIR Keys* specifies the additional associated keys. At each directory record level other additional data elements can be added, but it is not required that File Set Readers be able to use them as keys. Refer to the Basic Directory IOD in PS 3.3.

Table 3: PRI-XAMAS Additional DicomDIR Keys

Key Attribute	Tag	Directory Record Level	Type	Notes
Date of Birth	[0010,0030]	PATIENT	2C	required, if present in SOP Instance
Patient's Sex	[0010,0040]	PATIENT	2C	required, if present in SOP Instance
Series Date	[0008,0021]	SERIES	3	
Series Time	[0008,0031]	SERIES	3	
Institute Name	[0008,0080]	SERIES	2C	required, if present in SOP Instance
Institution Address	[0008,0081]	SERIES	2C	required, if present in SOP Instance
Series Description	[0008,103E]	SERIES	3	
Performing Physician's Name	[0008,1050]	SERIES	2C	required, if present in SOP Instance
Image Type	[0008,0008]	IMAGE	1C	required, if present in SOP Instance
SOP Class UID	[0008,0016]	IMAGE	3	
SOP Instance UID	[0008,0018]	IMAGE	3	
Image Date	[0008,0023]	IMAGE	3	
Image Time	[0008,0033]	IMAGE	3	
Referenced Image Sequence	[0008,1140]	IMAGE	1C	required, if present in SOP Instance
>Referenced SOP Class UID	[0008,1150]	IMAGE		
>Referenced SOP Instance UID	[0008,1155]	IMAGE		
Image Position (Patient)	[0020,0032]	IMAGE	2C	required, if present in SOP Instance
Image Orientation (Patient)	[0020,0037]	IMAGE	2C	required, if present in SOP Instance
Frame of Reference UID	[0020,0052]	IMAGE	2C	required, if present in SOP Instance
Rows	[0028,0010]	IMAGE	3	

Table 3: PRI-XAMAS Additional DICOMDIR Keys

Key Attribute	Tag	Directory Record Level	Type	Notes
Columns	[0028,0011]	IMAGE	3	
Pixel Spacing	[0028,0030]	IMAGE	1C	required, if present in SOP Instance
Calibration Image	[0050,0004]	IMAGE	2C	required, if present in SOP Instance
Icon Image Sequence	[0088,0200]	IMAGE	3	required for XA SOP Class
>Samples per Pixel	[0028,0002]			1
>Photometric Interpretation	[0028,0004]			MONOCHROME2
>Rows	[0028,0010]			128 for XA, 64 for others
>Columns	[0028,0011]			128 for XA, 64 for others
>Bits Allocated	[0028,0100]			8
>Bits Stored	[0028,0101]			8
>High Bit	[0028,0102]			7
>Pixel Representation	[0028,0103]			0
>Pixel Data	[7FE0,0010]			Icon Image

3.3.3 Private Directory Record Keys

Private Directory Records are supported by this Application Profile Class at the following Level - IMAGE. The PRIVATE Directory Records will have required elements in addition to the mandatory elements specified in PS 3.3.

Table 4: will list the additional required keys for PRIVATE Directory Records.

Table 4: PRI-XAMAS Additional PRIVATE Record Keys

Key Attribute	Tag	Type	Notes
Private Record UID	(0004,1432)	1	see Conformance Statement
SOP Class UID	[0008,0016]	1C	required, if present in SOP Instance
SOP Instance UID	[0008,0018]	1C	required, if present in SOP Instance
Image Type	[0008,0008]	3	
Acquisition Date	[0008,0022]	3	
Acquisition Time	[0008,0032]	3	
Acquisition Number	[0020,0012]	3	
CSA Data Type	[0029,xx08]	1	private owner code = SIEMENS CSA NON-IMAGE
CSA Data Version	[0029,xx09]	3	private owner code = SIEMENS CSA NON-IMAGE

3.3.4 Icon Images

Directory Records of type SERIES or IMAGE may include Icon Images. The Icon Image pixel data shall be as specified in PS 3.3 "Icon Image Key Definition", and restricted such that Bits Allocated (0028,0100) and Bits Stored (0028,0101) shall be equal 8, and Rows (0028,0010) and Columns (0028,0011) shall be equal to 128 for XA Images and 64 for all other Images. The Photometric Interpretation (0028,0004) shall always be restricted to "MONOCHROME2".

PRIVATE Directory Records will not contain Icon Image information.

3.4 Other Parameters

This section defines other parameters common to all specific Application Profiles in the PRI-XAMAS class which need to be specified in order to ensure interoperable media interchange.

3.4.1 Multiframe JPEG Format

The JPEG encoding of pixel data shall use Interchange Format (with table specification) for all frames.

3.4.2 Attribute Value Precedence

The values of attributes contained in a Detached Patient Management SOP Instance referenced by a DICOMDIR PATIENT Directory Record shall take precedence over the values of those attributes contained in the SOP Instance referenced by a subsidiary Directory Record. The DICOMDIR Directory Records shall have key attribute values in accordance with this precedence.

Note

This allows patient identification and demographic information to be updated without changing the composite Image IOD files. The DICOMDIR file thus is critical in establishing the link between the updated information and the image. As an example, at the time an Image file was written, the patient's name therein was incorrect, or inconsistent with the Hospital Information System records. Subsequently, a Detached Patient management file with the corrected name is added to the file-set. The FSR should use the name from the Patient file rather than in the image File.

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

Application Profile Conformance Statement

This chapter will contain the Conformance Statement to all "Offline Media Application Profiles (incl. private extensions)" supported by the M.A.S. archive options.

Those application profiles supported shall be:

- 1024 Extended Cardiac
- General Purpose CDR
- M.A.S. private Application Profile

1 Introduction

This Conformance Statement is a "Development Level Version".

That means there are open points and unspecified terms (e.g. xxxxx) mentioned that will be solved during the development phases. The release for product "Product Level Version" will no longer contain any of these above mentioned open points or unspecific terms.

1.1 Purpose

This DICOM Conformance Statement is written according to part PS 3.2 of [1].

The applications described in this conformance statement are the SIEMENS AXIOM Artis based products using the AXIOM Artis software. The AXIOM Artis DICOM offline media storage service implementation acts as FSC, FSU and/or FSR for the specified application profiles and the related SOP Class instances.

1.2 Scope

This DICOM Conformance Statement refers to SIEMENS AXIOM Artis based products using the AXIOM Artis software. The following table relates AXIOM Artis software to SIEMENS products.

Table 1: SIEMENS AXIOM Artis based products

Table 2:

Software Name	SIEMENS AXIOM Artis based products
VB11	AXIOM Artis

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

1.3 Definitions, Abbreviations

1.3.1 Definitions

DICOM	Digital Imaging and Communications in Medicine
DIMSE	DICOM Message Service Element
DIMSE-C	DICOM Message Service Element with Composite information objects

1.3.2 Abbreviations

ACR	American College of Radiology
AE	DICOM Application Entity
ASCII	American Standard Code for Information Interchange

CDR, CD-R	Compact Disk Recordable
DB	Database
DCS	DICOM Conformance Statement
FSC	File Set Creator
FSR	File Set Reader
FSU	File Set Updater
IOD	DICOM Information Object Definition
ISO	International Standard Organization
MAS or M.A.S.	Modular Angio System
NEMA	National Electrical Manufacturers Association
O	Optional Key Attribute
PDU	DICOM Protocol Data Unit
R	Required Key Attribute
RWA	Real-World Activity
U	Unique Key Attribute

1.4 References

- [1] Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-16

1.5 Connectivity and Interoperability

This Conformance Statement by itself does not guarantee successful interoperability of SIEMENS equipment with non-SIEMENS equipment. The user (user's agent) should be aware of the following issues:

- **Interoperability**

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of SIEMENS equipment with non-SIEMENS equipment. It is the user's responsibility to analyze thoroughly the application requirements and to specify a solution that integrates SIEMENS equipment with non-SIEMENS equipment.

- **Validation**

SIEMENS equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.


Where SIEMENS equipment is linked to non-SIEMENS equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation test will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.

- **New versions of the DICOM Standard**

the DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. SIEMENS is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, SIEMENS reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-SIEMENS provider linking to SIEMENS equipment, also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into SIEMENS equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

2 Implementation Model

2.1 Application Data Flow Diagram



Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

The DICOM archive application will serve as an interface to the CD-R offline medium device. It serves interfaces to include the offline media directory into the browser and to copy SOP instances to a medium or retrieve SOP Instances from medium into local storage.

The DICOM Archive application will support the 120 mm CD-R medium, the 130mm 2.3GB R/W MOD and the 130mm 4.1GB R/W MOD media.

The FSU role will update new SOP Instances only to media with pre-existing File-sets conforming to the Application Profiles supported.

The contents of the DICOMDIR will be temporarily stored in Archiv-Database.

2.2 Functional definitions of AE's

The AXIOM Artis DICOM offline media storage application consist of the DICOM Archive application entity serving all interfaces to access offline media. The DICOM Archive application is capable of

1. creating a new File-set onto an unwritten medium (Export to...).
2. updating an existing File-set by writing new SOP Instances onto the medium (Export to...).
3. importing SOP Instances from the medium onto local storage

-
4. reading the File-sets DICOMDIR information into temporary database and pass it to display applications.

2.3 Sequencing of Real World Activities

The DICOM Archive application will not perform updates before the Directory information of the DICOMDIR is completely read.

When performing updates, the SOP instances are checked for existence before updating. Duplicate instances will be avoided.

2.4 File Meta Information Options

The Implementation Class UID is:

- 1.3.12.2.1107.5.9.20000101

and an Implementation Version Name of

- "SIEMENS_SWFVX47A"

3 AE Specifications

3.1 DICOM Archive Specification

The DICOM Archive provides Standard conformance to Media Storage Service Class (Interchange Option). In addition Augmented conformance is provided to store extra data attributes important for the full feature support of the AXIOM Artis product SW. Details are listed in Table 3:.

Table 3: Application profiles, Activities, and Roles for DICOM Archive

Application Profiles Supported	Real World Activity	Role	SC Option
PRI-XAMAS-CD AUG-XA1K-CD *1	Browse Directory Information	FSR	Interchange
	Import into local Storage	FSR	Interchange
	Export to local archive media	FSC,FSU	Interchange
STD-GEN-CD STD-XABC-CD STD-XA1K-CD	Browse Directory Information	FSR	Interchange
	Import into local Storage	FSR	Interchange

*1 With no Private SOP Class used, the PRI-XAMAS-CD profile definitions are appropriate to describe the Augmentation of the STD-XA1K-CD profile.

On AXIOM Artis based Cardiac Systems the Private Extended MAS Profile (PRI-XAMAS-CD) will be preferably used by the system. The General Purpose Interchange Profile (STD-GEN-CD), Basic Cardiac Profile (STD-XABC-CD) and 1024 X-Ray Angiographic Profile (STD-XA1K-CD) will be supported with read capability of the related media.

3.1.1 File Meta Information for the Application Entity

The Source Application Entity Title is set by configuration. See Chapter "Configuration" for details.

3.1.2 Real-World Activities for this Application Entity

3.1.2.1 Real-World Activity: Browse Directory Information

The DICOM Archive application acts as FSR using the interchange option when requested to read the media directory.

The DICOM archive application will read the DICOMDIR and insert those directory entries, that are valid for the application profiles supported, into a local database. The database can then be used for browsing media contents.

Note

IconImageSQ is also supported in DICOMDIR. But only those IconImages with BitsAllocated (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the browser.

3.1.2.1.1 Application Profiles for the RWA: Browse Directory Information

See Table 3: for the Application Profiles listed that invoke this Application Entity for the Browse Directory Information RWA.

3.1.2.2 Real-World Activity: Import into local Storage

The DICOM Archive application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the local storage.

The SOP Instance selected from the media directory will be copied into the local storage. Only SOP Instances, that are valid for the application profile supported and are listed as supported by the Storage SCP Conformance section (Network DCS, 5.1.3), can be retrieved from media storage. This is due to the fact that the Browse Directory Information will filter all SOP Instances not matching the Application Profiles supported.

During operation no "Attribute Value Precedence" is applied to the SOP Instances. Detached Patient Management is not supported.

For media conforming to the STD-GEN-CD Profile the following SOP classes will be supported as an FSR:

Table 4: STD-GEN-CD Supported SOP Classes for FSR role

Information Object Definition	Service Object Pair Class UID	Transfer Syntax and UID
X-Ray Angiographic Image	1.2.840.10008.5.1.4.1.1.12.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
CT Image	1.2.840.10008.5.1.4.1.1.2	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
MR Image	1.2.840.10008.5.1.4.1.1.4	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1
CSA Non-Image	1.3.12.2.1107.5.9.1	Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1

3.1.2.2.1 Application Profiles for the RWA: Import into local Storage

See Table 3: for the Application Profiles listed that invoke this Application Entity for the Import into Local Storage RWA.

3.1.2.3 Real-World Activity: Export to local Archive Media

The DICOM Archive application acts as FSU (for media with existing DICOM file-set) or FSC (media not initialized) using the interchange option when requested to copy SOP Instances from the local storage to local Archive medium.

The DICOM Archive application will receive a list of SOP Instances to be copied to the local archive medium. According to the state of the medium inserted (new medium, Medium with DICOM file-set) the validity of the SOP Instances according to the applicable profile is checked. Only valid SOP Instances are accepted.

When the DICOM archive application is requested to copy SOP Instances the preferred application profile according configuration (AUG-XA1K-CD or PRI-XAMAS-CD) will be used to validate and copy the referred SOP Instances. When creating a new file-set, no Descriptor File will be allocated and the related ID is not used.

The DICOM archive application will not close the medium.

With the "Export 512 to..." function of the Viewer or if configured as default transfer mode for the viewer, a copy of images in Cardiac Format (512x512, 8Bit) can be written onto medium. Please refer to the Storage section "Send 512 to..." RWA description to learn more about the type of conversion that is performed on the Instances.

3.1.2.3.1 Application Profiles for the RWA: Export to local Archive

See Table 3: for the Application Profiles listed that invoke this Application Entity for the Export to local Archive RWA.

4 Augmented and Private Profiles

4.1 Augmented Application Profiles

With no private Siemens Non-Images stored onto Medium, the definitions of the PRI-XAMAS-CD Profile are applicable to denote the augmentations for the STD-XA1K-CD Standard Profile.

Storage of Private Information Objects will only be supported with reference to a Private Application Profile (see next section).

4.2 Private Application Profiles

4.2.1 PRI-XAMAS-CD

The DICOM Archive application is conforming to the PRI-XAMAS-CD application profile provided in a separate document. Please refer to next section for detailed definition of SOP classes supported in accordance to this Application Profile.

5 Extensions, Specializations and Privatizations of SOP Classes and Transfer Syntaxes

The SOP Classes listed refer in majority to those created by the equipment to which this conformance Statement is related to. For SOP classes not listed in this section, please refer to the Storage section of the DICOM Conformance Statement of the product. This will include all SOP Instances that can be received and displayed and therefor will be included into offline media storage even though these SOP Instances are not created by the equipment serving the Media Storage Service.

5.1 SOP Specific Conformance Statement for Basic Directory

The Basic Directory Header will not contain Information to indicate the current medium is part of a sequence of multiple media holding the Objects of one patient/study in total.

5.1.1 Extension, Specialization for SIEMENS Non-Image Objects

According to the PRI-XAMAS Application Profile Class the usage of the Private Creator UIDs and further optional keys for the Directory Records referring to SIEMENS Non-Image Objects are listed in the following tables.

Table 5: Basic Directory Extension for Non-Image Objects

Attribute	Tag	Value used
Private Record UID	[0004,1432]	1.3.12.2.1107.5.9.1
SOP Class UID	[0008,0016]	1.3.12.2.1107.5.9.1

For those "Non-Images" no Icon Image Sequence will be generated.

6 Configuration

6.1 AE Title Mapping

6.1.1 DICOM Media Storage AE Title

The DICOM Storage application provides the application entity title:

CsaImageManager

7 Support of Extended Character Sets

The Siemens AXIOM Artis DICOM archive application supports the ISO 8859 Latin 1 (ISO-IR 100) character set.

Also the Japanese language character sets JIS X 0201 (ISO-IR 13 Japanese katakana and ISO-IR 14 Japanese romaji), JIS X 0208 (ISO-IR 87 Japanese kanji) and JIS X 0212 (ISO-IR 159 Supplementary Japanese kanji) are supported.

In case of a mismatch between the SCS tag (0008,0005) and the characters in an IOD received by the system, the following measures are taken to make the characters DICOM conform:

Try to import using ISO_IR 100. If this fails, convert each nonmatching character to a "?".

A Appendix

A.1 SIEMENS Private Non-Image IOD

For encoding binary data-streams not representing image data, Siemens has created a private "Non-Image IOD" according to the rules governed by the DICOM Standard. The following section will roll-out the definition of this Private IOD. It can be communicated with Network Storage Service and Offline Media Storage Services.


The Siemens "Non-Image IOD" is identified by a private Non-Image Storage SOP Class UID of

"1.3.12.2.1107.5.9.1"

A.1.1 Siemens Non-Image IOD - E-R Model

The E-R model in A.1.2 depicts those components of the DICOM Information Model which directly refer to the Siemens Non-Image IOD. The Frame of Reference IE, Overlay IE, Modality Lookup-Table IE, VOI Lookup-Table IE and Curve IE are not components of the Siemens Non-Image IOD.

Figure 1: DICOM Information Model

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.

A.1.2 Siemens Non-Image IOD - Module Table

Table 1: CSA Non-Image IOD Modules

IE	Module	Reference	Usage
Patient	Patient	[1] PS3.3 C.7.1.1	M
Study	General Study	[1] PS3.3 C.7.2.1	M
	Patient Study	[1] PS3.3 C.7.2.2	U
Series	General Series	[1] PS3.3 C.7.3.1	M
Equipment	General Equipment	[1] PS3.3 C.7.5.1	U
CSA	CSA Image Header	A.2.1	U
	CSA Series Header	A.2.2	U
	MEDCOM Header	A.2.3	U
	CSA Non-Image	A.1.3.1	M
	SOP Common	[1] PS3.3 C.12.1	M

A.1.3 Siemens Non-Image IOD - Modules

A.1.3.1 CSA Non-Image Module

The table in this section contains private IOD Attributes that describe CSA Non-Images..

Table 2: CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
Image Type	(0008,0008)	-	3	Image identification characteristics.
Acquisition Date	(0008,0022)	-	3	The date the acquisition of data that resulted in this data set started.
Acquisition Time	(0008,0032)	-	3	The time the acquisition of data that resulted in this data set started.
Conversion Type	(0008,0064)	-	3	Describes the kind of image conversion. Defined Terms: DV = Digitized Video, DI = Digital Interface, DF = Digitized Film, WSD = Workstation.
Referenced Image Sequence	(0008,1140)	-	3	A sequence which provides reference to a set of Image SOP Class/Instance identifying other images significantly related to this data set. Encoded as sequence of items: (0008,1150) and (0008,1155).
Derivation Description	(0008,2111)	-	3	A text description of how this data set was derived.
Source Image Sequence	(0008,2112)	-	3	A Sequence which identifies the set of Image SOP Class/Instance pairs of the Images which were used to derive this data set. Zero or more Items may be included in this Sequence. Encoded as sequence of items: (0008,1150) and (0008,1155).
Patient Position	(0018,5100)	-	3	Patient position descriptor relative to the equipment.
Acquisition Number	(0020,0012)	-	3	A number identifying the single continuous gathering of data over a period of time which resulted in this data set.
Image Number	(0020,0013)	-	3	A number that identifies this data set.
Frame of Reference UID	(0020,0052)	-	3	Uniquely identifies the frame of reference for a Series.
Image Comments	(0020,4000)	-	3	User-defined comments about the image.

Table 2: CSA Non-Image Module

Attribute Name	Tag	Owner	Type	Notes
Quality Control Image	(0028,0300)	-	3	Indicates whether or not this image is a quality control or phantom image. If this Attribute is absent, then the image may or may not be a quality control or phantom image. Enumerated Values: YES, NO.
Burned In Annotation	(0028,0301)	-	3	Indicates whether or not image contains sufficient burned in annotation to identify the patient and date the image was acquired. If this Attribute is absent, then the image may or may not contain burned in annotation. Enumerated Values: YES, NO.
Lossy Image Compression	(0028,2110)	-	3	Specifies whether an Image has undergone lossy compression. Enumerated Values: 00 = Image has NOT been subjected to lossy compression, 01 = Image has been subjected to lossy compression.
Lossy Image Compression Ratio	(0028,2112)	-	3	Describes the approximate lossy compression ratio(s) that have been applied to this image. May be multi valued if successive lossy compression steps have been applied.
CSA Data Type	(0029,xx08)	SIEMENS CSA NON-IMAGE	1	CSA Data identification characteristics. Defined Terms: BSR REPORT = Study Report Data RT3D CONFIG = InSpaceIS Data
CSA Data Version	(0029,xx09)	SIEMENS CSA NON-IMAGE	3	Version of CSA Data Info (0029,xx10) format and CSA Non-Image Data (7FE1,xx10) format.
CSA Data Info	(0029,xx10)	SIEMENS CSA NON-IMAGE	3	Information to describe the CSA Data (7FE1,xx10).
CSA Data	(7FE1,xx10)	SIEMENS CSA NON-IMAGE	2	Binary data as byte stream.

A.2 Siemens Standard Extended Modules

Table 3: CSA Image IOD Modules

IE	Module	Reference	Usage	Note
Image	CSA Image Header	A.2.1	U	private GG information
	CSA Series Header	A.2.2	U	
	MEDCOM Header	A.2.3	U	private <i>syngo</i> information
	MEDCOM OOG	A.2.4	U	if object graphics is attached to image
	Edge Enhancement	A.2.5	U	private Filter Information
	Angio Viewing	A.2.6	U	private Viewing information
	Angio Acquisition	A.2.7	U	additional private information about image acquisition
	Angio Quantitation	A.2.8	U	if image is a Quant Report
	Original Image Info	A.2.9	U	if derived image

A.2.1 CSA Image Header Module

The table in this section contains private IOD Attributes that describe the CSA Image Header.

Table 4: CSA Image Header Module

Attribute Name	Tag	Owner	Type	Notes
CSA Image Header Type	(0029,xx08)	SIEMENS CSA HEADER	1	CSA Image Header identification characteristics. Defined Terms: "NUM 4" = NUMARIS/4 "SOM 5" = SOMARIS/5
CSA Image Header Version	(0029,xx09)	SIEMENS CSA HEADER	3	Version of CSA Image Header Info (0029,xx10) format.
CSA Image Header Info	(0029,xx10)	SIEMENS CSA HEADER	3	Manufacturer model dependent information.

A.2.2 CSA Series Header Module

The table in this section contains private IOD Attributes that describe the CSA Series Header.

Table 5: CSA Series Header Module

Attribute Name	Tag	Owner	Type	Notes
CSA Series Header Type	(0029,xx18)	SIEMENS CSA HEADER	1	CSA Series Header identification characteristics. Defined Terms: "NUM 4" = NUMARIS/4
CSA Series Header Version	(0029,xx19)	SIEMENS CSA HEADER	3	Version of CSA Series Header Info (0029,xx20) format.
CSA Series Header Info	(0029,xx20)	SIEMENS CSA HEADER	3	Manufacturer model dependent information.

A.2.3 MEDCOM Header Module

The table in this section contains private IOD Attributes that describe MEDCOM Header.

Table 6: MEDCOM Header Module

Attribute Name	Tag	Owner	Type	Notes
MedCom Header Type	(0029,xx08)	SIEMENS MED-COM HEADER	1C	MedCom Header identification characteristics. Defined Terms: MEDCOM 1 Required if MedCom Header Info (0029,xx10) present.
MedCom Header Version	(0029,xx09)	SIEMENS MED-COM HEADER	2C	Version of MedCom Header Info (0029,xx10) format. Required if MEDCOM Header Info (0029,xx10) present.
MedCom Header Info	(0029,xx10)	SIEMENS MED-COM HEADER	3	Manufacturer model dependent information. The value of the attribute MedCom Header Info (0029,xx10) can be build up in any user defined format.
MedCom History Information	(0029,xx20)	SIEMENS MED-COM HEADER	3	MedCom defined Patient Registration history information. See A.2.3.1.
PMTF Information 1	(0029,xx31)	SIEMENS MED-COM HEADER	3	Transformation Information
PMTF Information 2	(0029,xx32)	SIEMENS MED-COM HEADER	3	Transformation Information
PMTF Information 3	(0029,xx33)	SIEMENS MED-COM HEADER	3	Transformation Information
PMTF Information 4	(0029,xx34)	SIEMENS MED-COM HEADER	3	Transformation Information

Table 6: MEDCOM Header Module

Attribute Name	Tag	Owner	Type	Notes
PMTF Information 5	(0029,xx35)	SIEMENS MED-COM HEADER	3	Transformation Information
Application Header Sequence	(0029,xx40)	SIEMENS MED-COM HEADER	3	Sequence of Application Header items. Zero or more items are possible.
>Application Header Type	(0029,xx41)	SIEMENS MED-COM HEADER	1C	Application Header identification characteristics. Required, if Sequence is sent.
>Application Header ID	(0029,xx42)	SIEMENS MED-COM HEADER	3	Identification of an application header
>Application Header Version	(0029,xx43)	SIEMENS MED-COM HEADER	3	Version of CSA Series Header Info (0029,xx44) format.
>Application Header Info	(0029,xx44)	SIEMENS MED-COM HEADER	3	Application dependent information.
Workflow Control Flags	(0029,xx50)	SIEMENS MED-COM HEADER	3	Eight free definable flags.
Archive Management Flag Keep Online	(0029,xx51)	SIEMENS MED-COM HEADER	3	Flag to control remote archive management system to keep the image always online (also when already archived). Enumerated Values: 00 = remote control not required 01 = keep image online
Archive Management Flag Do Not Archive	(0029,xx52)	SIEMENS MED-COM HEADER	3	Flag to control remote archive management system not to archive the related image. Enumerated Values: 00 = remote control not required 01 = don't archive image
Image Location Status	(0029,xx53)	SIEMENS MED-COM HEADER	3	Image location status to control retrieving. Defined Terms: ONLINE = retrieving has to be done as usual, NEARLINE = move request to SCP and delay according to value of Estimated Retrieve Time (0029,xx54), OFFLINE = invoking a retrieve operation initiates an operator request, INVALID = invoking a retrieve operation would always result in an error.
Estimated Retrieve Time	(0029,xx54)	SIEMENS MED-COM HEADER	3	Estimated retrieve time in seconds. A value less then zero (< 0) indicates location is OFFLINE or INVALID.
Data Size of Retrieved Images	(0029,xx55)	SIEMENS MED-COM HEADER	3	Data size of images in MByte.

A.2.3.1 MEDCOM History Information

The value of the attribute MEDCOM History Information (0029,xx20) is defined in the following way:

Table 7: MEDCOM History Information

Part	Name	Type	Bytes	Notes
header	Identifier	string	32	always "CSA HISTORY"
	Version	string	32	e.g. "V1.10"
>n items	Class Name	string	64	
	Modification String	string	1024	

A.2.4 MEDCOM OOG Module

The table in this section contains private IOD Attributes that describe MEDCOM Object Oriented Graphics (OOG). This module is used whenever object graphics is drawn on the image and need to be stored as graphic object properties. Given the condition that the module contents was not removed by other modalities, the graphic objects remain re-animatable if such an image was transferred and is then retrieved back.

Table 8: MEDCOM OOG Module

Attribute Name	Tag	Owner	Type	Notes
MedCom OOG Type	(0029,xx08)	SIEMENS MED-COM OOG	1	MEDCOM Object Oriented Graphics (OOG) identification characteristics. Defined Terms: MEDCOM OOG 1 MEDCOM OOG 2
MedCom OOG Version	(0029,xx09)	SIEMENS MED-COM OOG	3	Version of MEDCOM OOG Info (0029,xx10) format.
MedCom OOG Info	(0029,xx10)	SIEMENS MED-COM OOG	3	MEDCOM Object Oriented Graphics (OOG) data.

The graphics objects are also fully drawn in the Image Overlay Plane for compatibility with other products which do not support the MedCom OOG module. Any system not supporting the MedCom OOG module shall remove the OOG module and it's contents when modifying the image overlay plane content.

A.2.5 Edge Enhancement Module

The table in this section contains private IOD Attributes that describe Edge Enhancement extensions due to the Dynamic Cardio Review definition.

Table 9: Edge Enhancement Module

Attribute Name	Tag	Owner	Type	Notes
Standard Edge Enhancement Sequence	[0029,xx00]	CARDIO-D.R. 1.0	3	Standard formula according Dynaview Extensions
>Convolution Kernel Size	[0029,xx01]	CARDIO-D.R. 1.0	1C	x-/y-size value pair. Each value shall be greater or equal to 3. Required if Sequence is present
>Convolution Kernel Coefficients	[0029,xx02]	CARDIO-D.R. 1.0	1C	Row-by-row list of the kernel Coefficients. Required if Sequence is present
>Edge Enhancement Gain	[0029,xx03]	CARDIO-D.R. 1.0	1C	Applied Filter gain Factor. Range is from 0 to 100 Percent. Required if Sequence is present

A.2.6 Angio Viewing Module

The table in this section contains private IOD Attributes that describe additional Attributes for advanced Angio Viewing features.

Table 10: Angio Viewing Module

Attribute Name	Tag	Owner	Type	Notes
Review Mode	[0019,xx00]	SIEMENS SMS-AX VIEW 1.0	3	Special Modes for Angio Review. Defined Terms are 1 = REV_MAXFILL, 2 = REV_LOOP, 3 = REV_SCROLL, 4 = REV_STEREO_LOOP
Anatomical Background Percent	[0019,xx01]	SIEMENS SMS-AX VIEW 1.0	3	Percentage of Mix between Subtracted Image Result and Native Mask. Range is from 0 to 100.
Number of Phases	[0019,xx02]	SIEMENS SMS-AX VIEW 1.0	3	1-4 (1 or # of "Variable Frame Rate" acq phases)
Apply Anatomical Background	[0019,xx03]	SIEMENS SMS-AX VIEW 1.0	3	boolean
Pixel Shift Array	[0019,xx04]	SIEMENS SMS-AX VIEW 1.0	3	4 * Number of Frames (0028,0008)
Brightness	[0019,xx05]	SIEMENS SMS-AX VIEW 1.0	3	SUB windowing

Table 10: Angio Viewing Module

Attribute Name	Tag	Owner	Type	Notes
Contrast	[0019,xx06]	SIEMENS SMS-AX VIEW 1.0	3	SUB windowing
Enabled Shutters	[0019,xx07]	SIEMENS SMS-AX VIEW 1.0	3	visualize Shutter
Native Edge Enhancement Percent Gain	[0019,xx08]	SIEMENS SMS-AX VIEW 1.0	3	Percent gain for native display of images.
Native Edge Enhancement LUT Index	[0019,xx09]	SIEMENS SMS-AX VIEW 1.0	3	
Native Edge Enhancement Kernel Size	[0019,xx0A]	SIEMENS SMS-AX VIEW 1.0	3	
Subtracted Edge Enhancement Percent Gain	[0019,xx0B]	SIEMENS SMS-AX VIEW 1.0	3	Percent gain for subtracted display of images.
Subtracted Edge Enhancement LUT Index	[0019,xx0C]	SIEMENS SMS-AX VIEW 1.0	3	
Subtracted Edge Enhancement Kernel Size	[0019,xx0D]	SIEMENS SMS-AX VIEW 1.0	3	
Fade Percent	[0019,xx0E]	SIEMENS SMS-AX VIEW 1.0	3	
Flipped before Lateral-ity Applied	[0019,xx0F]	SIEMENS SMS-AX VIEW 1.0	3	
Apply Fade	[0019,xx10]	SIEMENS SMS-AX VIEW 1.0	3	
Zoom	[0019,xx12]	SIEMENS SMS-AX VIEW 1.0	3	
Pan X	[0019,xx13]	SIEMENS SMS-AX VIEW 1.0	3	
Pan Y	[0019,xx14]	SIEMENS SMS-AX VIEW 1.0	3	
Native Edge Enhancement Adverse Percent Gain	[0019,xx15]	SIEMENS SMS-AX VIEW 1.0	3	
Subtracted Edge Enhancement Adverse Percent Gain	[0019,xx16]	SIEMENS SMS-AX VIEW 1.0	3	
Invert Flag	[0019,xx17]	SIEMENS SMS-AX VIEW 1.0	3	
Quant 1K Overlay	[0019,xx1A]	SIEMENS SMS-AX VIEW 1.0	3	Only Store Monitor and Store Reference Images in QUANT (128 kByte)

A.2.7 Angio Acquisition Data Module

The table in this section contains private IOD Attributes that describe additional Attributes for saving specific Angio Acquisition data.

Table 11: Angio Viewing Module

Attribute Name	Tag	Owner	Type	Notes
Acquisition Type	[0021,xx00]	SIEMENS SMS-AX ACQ 1.0	3	Technical Type of Acquisition performed to get image result. Defined Terms are 1 = fixed frame rate, 2 = variable frame rate (manually trigrd), 3 = variable frame rate (time triggered), 4 = peri manual to head, 5 = peri slow speed to head, 6 = peri fast speed to head, 7 = peri manual to feet, 8 = peri slow speed to feet, 9 = peri fast speed to feet, 10 = dynavision manual inject., 11 = dynavision automatic, 12 = direct technique, 13 = continuous fluoro, 14 = pulsed fluoro, 15 = ECG triggered fluoro
Acquisition Mode	[0021,xx01]	SIEMENS SMS-AX ACQ 1.0	3	Technical Mode of Acquisition performed to get image result. Defined Terms are: 0 = no mode specified, 1 = Digital Radiographie, 2 = DSA, 3 = Peri-DSA, 4 = DR stepping, 5 = DR Dynavision (nat), 6 = Dynavision (sub), 7 = Card, 8 = Direct technique, 20 = Roadmap Fluoro, 21 = Normal Fluoro, 22 = Card Fluoro
Foot Switch Index	[0021,xx02]	SIEMENS SMS-AX ACQ 1.0	3	
Acquisition Room	[0021,xx03]	SIEMENS SMS-AX ACQ 1.0	3	
Current Time Product	[0021,xx04]	SIEMENS SMS-AX ACQ 1.0	3	(uAs)
Dose	[0021,xx05]	SIEMENS SMS-AX ACQ 1.0	3	(uGy/pulse)
Skin Dose Percent	[0021,xx06]	SIEMENS SMS-AX ACQ 1.0	3	
Skin Dose Accumulation	[0021,xx07]	SIEMENS SMS-AX ACQ 1.0	3	(mGy)

Table 11: Angio Viewing Module

Attribute Name	Tag	Owner	Type	Notes
Skin Dose Rate	[0021,xx08]	SIEMENS SMS-AX ACQ 1.0	3	(mGy/min)
Copper Filter	[0021,xx0A]	SIEMENS SMS-AX ACQ 1.0	3	(in 0.1 mm)
Measuring Field	[0021,xx0B]	SIEMENS SMS-AX ACQ 1.0	3	Decimal values are defined as follows: 1 = left 2 = middle (rectangular) 3 = left and middle 4 = right 5 = left and right 6 = right and middle 7 = left and right and middle 8 = central (approx. round) 9 = left and central 12 = right and central 13 = left and right and central
Post Blanking Circle	[0021,xx0C]	SIEMENS SMS-AX ACQ 1.0	3	left edge, upper edge, radius
Dyna Angles	[0021,xx0D]	SIEMENS SMS-AX ACQ 1.0	3	May be used in previous versions. Data is now encoded in Standard attributes.
Total Steps	[0021,xx0E]	SIEMENS SMS-AX ACQ 1.0	3	Number of Images which make up a Peri or DR-Step acquisition
Dyna X-Ray Info	[0021,xx11]	SIEMENS SMS-AX ACQ 1.0	3	Multiplicity of (3 * Number of Frames). (kV value, uA value, ms value) Only with Dyna Images
Modality LUT Input Gamma	[0021,xx10]	SIEMENS SMS-AX ACQ 1.0	3	
Modality LUT Output Gamma	[0021,xx11]	SIEMENS SMS-AX ACQ 1.0	3	
SH_STPAR	[0021,xx12]	SIEMENS SMS-AX ACQ 1.0	3	1st Byte = version. Technical data-structure with Gantry information.
Acquisition Zoom	[0021,xx13]	SIEMENS SMS-AX ACQ 1.0	3	0 = Off, 1 = On. Not to be mixed up with digital zoom from Viewing module.
Dyna Angulation Step Width	[0021,xx14]	SIEMENS SMS-AX ACQ 1.0	3	Only with Dyna Images
Harmonization	[0021,xx15]	SIEMENS SMS-AX ACQ 1.0	3	
DR Single Flag	[0021,xx16]	SIEMENS SMS-AX ACQ 1.0	3	1 = DR SINGLE, else single Frame DR acquisition
Source to Isocenter	[0021,xx17]	SIEMENS SMS-AX ACQ 1.0	3	(mm)
Pressure Data	[0021,xx18]	SIEMENS SMS-AX ACQ 1.0	3	

Table 11: *Angio Viewing Module*

Attribute Name	Tag	Owner	Type	Notes
ECG Index Array	[0021,xx19]	SIEMENS SMS-AX ACQ 1.0	3	
FD Flag	[0021,xx1A]	SIEMENS SMS-AX ACQ 1.0	3	
SH_ZOOM	[0021,xx1B]	SIEMENS SMS-AX ACQ 1.0	3	1st Byte = version. Technical data-structure with Gantry information.
SH_COLPAR	[0021,xx1C]	SIEMENS SMS-AX ACQ 1.0	3	1st Byte = version. Technical data-structure with Gantry information.
K-Factor	[0021,xx1D]	SIEMENS SMS-AX ACQ 1.0	3	
EVE	[0021,xx1E]	SIEMENS SMS-AX ACQ 1.0	3	
Total Scene Tim	[0021,xx1F]	SIEMENS SMS-AX ACQ 1.0	3	
Restore Flag	[0021,xx20]	SIEMENS SMS-AX ACQ 1.0	3	
Stand Movement Flag	[0021,xx21]	SIEMENS SMS-AX ACQ 1.0	3	
FD Rows	[0021,xx22]	SIEMENS SMS-AX ACQ 1.0	3	
FD Columns	[0021,xx23]	SIEMENS SMS-AX ACQ 1.0	3	
Table Movement Flag	[0021,xx24]	SIEMENS SMS-AX ACQ 1.0	3	

Copyright © Siemens AG 2004. All rights reserved.
 Alle Rechte vorbehalten.

A.2.8 Angio Quantitation Module

The table in this section contains private IOD Attributes that describe additional Attributes for advanced Angio Quantitation and Calibration Results features.

Table 12: *Angio Viewing Module*

Attribute Name	Tag	Owner	Type	Notes
Horizontal Calibration Pixel Size	[0023,xx00]	SIEMENS SMS-AX QUANT 1.0	3	(in mm)
Vertical Calibration Pixel Size	[0023,xx01]	SIEMENS SMS-AX QUANT 1.0	3	(in mm)
Calibration Object	[0023,xx02]	SIEMENS SMS-AX QUANT 1.0	3	
Calibration Object Size	[0023,xx03]	SIEMENS SMS-AX QUANT 1.0	3	

Table 12: Angio Viewing Module

Attribute Name	Tag	Owner	Type	Notes
Calibration Method	[0023,xx04]	SIEMENS SMS-AX QUANT 1.0	3	
Filename	[0023,xx05]	SIEMENS SMS-AX QUANT 1.0	3	
Frame Number	[0023,xx06]	SIEMENS SMS-AX QUANT 1.0	3	
Calibration Factor Multiplicity	[0023,xx07]	SIEMENS SMS-AX QUANT 1.0	3	Multiplicity Horizontal followed by Multiplicity for Vertical

A.2.9 Original Image Info Module

The table in this section contains private IOD attributes that describe additional original image data for derived images (e.g. Store Monitor Image).

Table 13: Original Image Info Module

Attribute Name	Tag	Owner	Type	Notes
View Native	[0025,xx00]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Series Number	[0025,xx01]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Image Number	[0025,xx02]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Win Center	[0025,xx03]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Win Width	[0025,xx04]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Win Brightness	[0025,xx05]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Win Contrast	[0025,xx06]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Frame Number	[0025,xx07]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Mask Frame Number	[0025,xx08]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Opac	[0025,xx09]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Number of Frames	[0025,xx0A]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Scene Duration	[0025,xx0B]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	

Table 13: Original Image Info Module

Attribute Name	Tag	Owner	Type	Notes
Identifier LOID	[0025,xx0C]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original Scene VFR Info	[0025,xx0D]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	Number of phases, then followed by n pairs (LastFrameNumber, then Frame Rate)
Original Frame ECG Position	[0025,xx0E]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Original ECG 1st Frame Offset	[0025,xx0F]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	<retired, no longer used>
Zoom Flag	[0025,xx10]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
<i>reserved</i>	[0025,xx11]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	-	
Number of Mask Frames	[0025,xx12]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	
Number of Fill Frames	[0025,xx13]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	3	

A.2.10 SOP Common Module - Image Type Extentions

Additional values for the Image type attribute are used to designate the purpose and contents of the SOP Instance created by the AXIOM Artis system. Please the following table for details:

Table 14: Image Type Extentions

Type of Scene/Image	Image Type
Single Plane Scene	ORIGINAL\PRIMARY\SINGLE PLANE\SINGLE A (or B)
Biplane Scene	ORIGINAL\PRIMARY\BIPLANE A (or B)
Reference Image	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A(B)\REFIMAGE
Reference Image from CT	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A\REFIMAGE\CT REF
Reference Image from MR	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A\REFIMAGE\MR REF
Dynamic Reference Image (Dynamap)	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A(B)\DYNAMAP
Store Monitor	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A(B)\STORE MONITOR
Secondary Capture	DERIVED\SECONDARY\SINGLE PLANE\PLANE A(B)\SC
Perivision - Mask Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\PERI\MASK\ ⁿ
Perivision - Fill Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\PERI\FILL\ ⁿ
Dynavision - Mask Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\MASK
Dynavision - Fill Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\FILL

Table 14: Image Type Extentions

Type of Scene/Image	Image Type
Dynavision - Injection Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\INJECTION
Dynavision - Washout Phase	ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\WASHOUT
Store Fluoro Loop	DERIVED\PRIMARY\SINGLE PLANE\SINGLE A(or B)\STORE FLUORO
Store Biplan Fluoro Loop	DERIVED\PRIMARY\BIPLANE A(or B)\STORE FLUORO
Quant Result Image	DERIVED\SECONDARY\<3rd to nth value from Original>\QUANT
Study Report (Non-Image Object)	ORIGINAL\PRIMARY\SINGLE PLANE\STUDY REPORT

¹ the term "n" is replaced by the related Step Number where the Image was acquired. Step Numbers start with "1".

A.3 Registry of DICOM Data Elements

Tag	Private Owner Code	Name	VR	VM
(0019,xx00)	SIEMENS SMS-AX VIEW 1.0	Review Mode	US	1
(0019,xx01)	SIEMENS SMS-AX VIEW 1.0	Anatomical Background Percent	US	1
(0019,xx02)	SIEMENS SMS-AX VIEW 1.0	Number of Phases	US	1
(0019,xx03)	SIEMENS SMS-AX VIEW 1.0	Apply Anatomical Background	US	1
(0019,xx04)	SIEMENS SMS-AX VIEW 1.0	Pixelshift Array	SS	4-4n
(0019,xx05)	SIEMENS SMS-AX VIEW 1.0	Brightness	US	1
(0019,xx06)	SIEMENS SMS-AX VIEW 1.0	Contrast	US	1
(0019,xx07)	SIEMENS SMS-AX VIEW 1.0	Enabled Shutters	US	1
(0019,xx08)	SIEMENS SMS-AX VIEW 1.0	Native Edge Enh. Percent Gain	US	1
(0019,xx09)	SIEMENS SMS-AX VIEW 1.0	Native Edge Enh. LUT Index	SS	1
(0019,xx0A)	SIEMENS SMS-AX VIEW 1.0	Native Edge Enh. Kernel Size	SS	1
(0019,xx0B)	SIEMENS SMS-AX VIEW 1.0	Subtr. Edge Enh. Percent Gain	US	1
(0019,xx0C)	SIEMENS SMS-AX VIEW 1.0	Subtr. Edge Enh. LUT Index	SS	1
(0019,xx0D)	SIEMENS SMS-AX VIEW 1.0	Subtr. Edge Enh. Kernel Size	SS	1
(0019,xx0E)	SIEMENS SMS-AX VIEW 1.0	Fade Percent	US	1
(0019,xx0F)	SIEMENS SMS-AX VIEW 1.0	Flipped Befor Laterality Applied	US	1
(0019,xx10)	SIEMENS SMS-AX VIEW 1.0	Apply Fade	DS	1
(0019,xx12)	SIEMENS SMS-AX VIEW 1.0	Zoom	US	1
(0019,xx13)	SIEMENS SMS-AX VIEW 1.0	Pan X	SS	1
(0019,xx14)	SIEMENS SMS-AX VIEW 1.0	Pan Y	SS	1
(0019,xx15)	SIEMENS SMS-AX VIEW 1.0	Native Edge Enh. Adv Perc. Gain	SS	1
(0019,xx16)	SIEMENS SMS-AX VIEW 1.0	Subtr. Edge Enh. Adv Perc. Gain	SS	1
(0019,xx17)	SIEMENS SMS-AX VIEW 1.0	Invert Flag	US	1
(0019,xx1A)	SIEMENS SMS-AX VIEW 1.0	Quant 1K Overlay (for downsize)	OB	1

Tag	Private Owner Code	Name	VR	VM
(0021,xx00)	SIEMENS SMS-AX ACQ 1.0	Acquisition Type	US	1
(0021,xx01)	SIEMENS SMS-AX ACQ 1.0	Acquisition Mode	US	1
(0021,xx02)	SIEMENS SMS-AX ACQ 1.0	Footswitch Index	US	1
(0021,xx03)	SIEMENS SMS-AX ACQ 1.0	Acquisition Room	US	1
(0021,xx04)	SIEMENS SMS-AX ACQ 1.0	Current Time Product	SL	1
(0021,xx05)	SIEMENS SMS-AX ACQ 1.0	Dose	SL	1
(0021,xx06)	SIEMENS SMS-AX ACQ 1.0	Skin Dose Percent	SL	1
(0021,xx07)	SIEMENS SMS-AX ACQ 1.0	Skin Dose Accumulation	SL	1
(0021,xx08)	SIEMENS SMS-AX ACQ 1.0	Skin Dose Rate	SL	1
(0021,xx0A)	SIEMENS SMS-AX ACQ 1.0	Copper Filter	UL	1
(0021,xx0B)	SIEMENS SMS-AX ACQ 1.0	Measuring Field	US	1
(0021,xx0C)	SIEMENS SMS-AX ACQ 1.0	Post Blanking Circle	SS	3
(0021,xx0D)	SIEMENS SMS-AX ACQ 1.0	Dyna Angles	SS	2-2n
(0021,xx0E)	SIEMENS SMS-AX ACQ 1.0	Total Steps	SS	1
(0021,xx0F)	SIEMENS SMS-AX ACQ 1.0	Dyna X-Ray Info	SL	3-3n
(0021,xx10)	SIEMENS SMS-AX ACQ 1.0	Modality LUT Input Gamma	US	1
(0021,xx11)	SIEMENS SMS-AX ACQ 1.0	Modality LUT Output Gamma	US	1
(0021,xx12)	SIEMENS SMS-AX ACQ 1.0	SH_STPAR	OB	1-n
(0021,xx13)	SIEMENS SMS-AX ACQ 1.0	Acquisition Zoom	US	1
(0021,xx14)	SIEMENS SMS-AX ACQ 1.0	Dyna Angulation Step Width	SS	1
(0021,xx15)	SIEMENS SMS-AX ACQ 1.0	Harmonization	US	1
(0021,xx16)	SIEMENS SMS-AX ACQ 1.0	DR Single Flag	US	1
(0021,xx17)	SIEMENS SMS-AX ACQ 1.0	Source to Isocenter	SL	1
(0021,xx18)	SIEMENS SMS-AX ACQ 1.0	Pressure Data	US	1-n
(0021,xx19)	SIEMENS SMS-AX ACQ 1.0	ECG Index Array	SL	1-n
(0021,xx1A)	SIEMENS SMS-AX ACQ 1.0	FD Flag	US	1
(0021,xx1B)	SIEMENS SMS-AX ACQ 1.0	SH_ZOOM	OB	1-n
(0021,xx1C)	SIEMENS SMS-AX ACQ 1.0	SH_COLPAR	OB	1-n
(0021,xx1D)	SIEMENS SMS-AX ACQ 1.0	K-Factor	US	1
(0021,xx1E)	SIEMENS SMS-AX ACQ 1.0	EVE	US	1
(0021,xx1F)	SIEMENS SMS-AX ACQ 1.0	Total Scene Time	SL	1
(0021,xx20)	SIEMENS SMS-AX ACQ 1.0	Restore Flag	US	1
(0021,xx21)	SIEMENS SMS-AX ACQ 1.0	Stand Movement Flag	US	1
(0021,xx22)	SIEMENS SMS-AX ACQ 1.0	FD Rows	US	1
(0021,xx23)	SIEMENS SMS-AX ACQ 1.0	FD Columns	US	1
(0021,xx24)	SIEMENS SMS-AX ACQ 1.0	Table Movement Flag	US	1
Tag	Private Owner Code	Name	VR	VM
(0023,xx00)	SIEMENS SMS-AX QUANT 1.0	Horizontal Calibration Pixel Size	DS	2

(0023,xx01)	SIEMENS SMS-AX QUANT 1.0	Vertical Calibration Pixel Size	DS	2
(0023,xx02)	SIEMENS SMS-AX QUANT 1.0	Calibration Object	LO	1
(0023,xx03)	SIEMENS SMS-AX QUANT 1.0	Calibration Object Size	DS	1
(0023,xx04)	SIEMENS SMS-AX QUANT 1.0	Calibration Method	LO	1
(0023,xx05)	SIEMENS SMS-AX QUANT 1.0	Filename	ST	1
(0023,xx06)	SIEMENS SMS-AX QUANT 1.0	Frame Number	IS	1
(0023,xx07)	SIEMENS SMS-AX QUANT 1.0	Calibration Factor Multiplicity	IS	2

Tag	Private Owner Code	Name	VR	VM
(0025,xx00)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	View Native	US	1
(0025,xx01)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Series Number	US	1
(0025,xx02)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Image Number	US	1
(0025,xx03)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Win Center	US	1
(0025,xx04)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Win Width	US	1
(0025,xx05)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Win Brightness	US	1
(0025,xx06)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Win Contrast	US	1
(0025,xx07)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Frame Number	US	1
(0025,xx08)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Mask Frame Number	US	1
(0025,xx09)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Opac	US	1
(0025,xx0A)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Number of Frames	US	1
(0025,xx0B)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Scene Duration	DS	1
(0025,xx0C)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Identifier LOID	LO	1
(0025,xx0D)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Scene VFR Info	SS	1-n
(0025,xx0E)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original Frame ECG Position	SS	1
(0025,xx0F)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Original ECG 1st Frame Offset (retired)	SS	1
(0025,xx10)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Zoom Flag	SS	1
(0025,xx11)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	<i>reserved</i>	-	-
(0025,xx12)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Number of Mask Frames	US	1

(0025,xx13)	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0	Number of Fill Frames	US	1
-------------	--	-----------------------	----	---

Tag	Private Owner Code	Name	VR	VM
(0029,xx08)	SIEMENS CSA HEADER	CSA Image Header Type	CS	1
(0029,xx09)	SIEMENS CSA HEADER	CSA Image Header Version	LO	1
(0029,xx10)	SIEMENS CSA HEADER	CSA Image Header Info	OB	1
(0029,xx18)	SIEMENS CSA HEADER	CSA Series Header Type	CS	1
(0029,xx19)	SIEMENS CSA HEADER	CSA Series Header Version	LO	1
(0029,xx20)	SIEMENS CSA HEADER	CSA Series Header Info	OB	1
(0029,xx08)	SIEMENS CSA NON-IMAGE	CSA Data Type	CS	1
(0029,xx09)	SIEMENS CSA NON-IMAGE	CSA Data Version	LO	1
(0029,xx10)	SIEMENS CSA NON-IMAGE	CSA Data Info	OB	1
(0029,xx08)	SIEMENS MEDCOM HEADER	MedCom Header Type	CS	1
(0029,xx09)	SIEMENS MEDCOM HEADER	MedCom Header Version	LO	1
(0029,xx10)	SIEMENS MEDCOM HEADER	MedCom Header Info	OB	1
(0029,xx20)	SIEMENS MEDCOM HEADER	MedCom History Information	OB	1
(0029,xx31)	SIEMENS MEDCOM HEADER	PMTF Information 1	LO	1
(0029,xx32)	SIEMENS MEDCOM HEADER	PMTF Information 2	UL	1
(0029,xx33)	SIEMENS MEDCOM HEADER	PMTF Information 3	UL	1
(0029,xx34)	SIEMENS MEDCOM HEADER	PMTF Information 4	CS	1
(0029,xx08)	SIEMENS MEDCOM OOG	MEDCOM OOG Type	CS	1
(0029,xx09)	SIEMENS MEDCOM OOG	MEDCOM OOG Version	LO	1
(0029,xx10)	SIEMENS MEDCOM OOG	MEDCOM OOG Info	OB	1
(0029,xx00)	CARDIO-D.R. 1.0	Standard Edge Enhancement Sequence	SQ	1
(0029,xx01)	CARDIO-D.R. 1.0	Convolution Kernel Size	US	2
(0029,xx02)	CARDIO-D.R. 1.0	Convolution Kernel Coefficients	US	1-n
(0029,xx03)	CARDIO-D.R. 1.0	Edge Enhancement Gain	FL	1

Tag	Private Owner Code	Name	VR	VM
(7FE1,xx10)	SIEMENS CSA NON-IMAGE	CSA Data	OB	1

Note

Please be informed that some of the Private Owner Codes contain double-spaces in the name definitions. The following terms (only double-spaces marked) are defined:

SIEMENS SMS-AX<spc><spc>VIEW 1.0, SIEMENS SMS-AX<spc><spc>ACQ 1.0, SIEMENS SMS-AX<spc><spc>QUANT 1.0, SIEMENS SMS-AX<spc><spc>ORIGINAL IMAGE INFO 1.0

All spaces not specially marked, are single spaces.

A.4 Extended XA IOD

The AXIOM Artis system will create images during acquisition and with postprocessing applications. Those will be encoded as XA Standard Extended SOP Class. Images created during postprocessing will be marked as derived. Please see the following table for a complete overview of supplied Type 1/2/3 Standard and additional Private Attributes:

The Standard DICOM Modules will be used to encode ECG data.

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Specific Character Set	[0008,0005]	from Configuration
Image Type	[0008,0008]	see "A.2.10 SOP Common Module - Image Type Extensions"
SOP Class UID	[0008,0016]	1.2.840.10008.5.1.4.1.1.12.1
SOP Instance UID	[0008,0018]	
Study Date	[0008,0020]	<yyyymmdd>
Series Date	[0008,0021]	<yyyymmdd>
Acquisition Date	[0008,0022]	Original: same as Content Date, Derived: Date of Original Acquisition
Content Date	[0008,0023]	<yyyymmdd> (Date of Creation)
Study Time	[0008,0030]	<hhmmss>
Series Time	[0008,0031]	<hhmmss>
Acquisition Time	[0008,0032]	Original: Same as Content Time, Derived: Time or Original Acquisition
Content Time	[0008,0033]	<hhmmss> (Time of Creation)
Accession Number	[0008,0050]	RIS or "Accession No." input
Modality	[0008,0060]	XA
Manufacturer	[0008,0070]	Siemens
Institution Name	[0008,0080]	RIS or "Institution Name" input
Institution Address	[0008,0081]	RIS
Referring Physician's Name	[0008,0090]	RIS or input
Station Name	[0008,1010]	from Configuration (VB22x and later)
Study Description	[0008,1030]	"Study" input
Series Description	[0008,103E]	"Organ Program" input or "Organ Program" + (REF)/ (SM)/ (FLM) or one of : STUDY REPORT, CARE GRAPH, or QUANT, DYNAMAP
Physicians of Record	[0008,1048]	from RIS [0040,1010] Names of Intended Recipients of Results

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Performing Physicians' Name	[0008,1050]	"Performing Physician 1" \ "Performing Physician 2" input
Operator's Name	[0008,1070]	"Operator 1" \ "Operator 2" input
Admitting Diagnoses Description	[0008,1080]	"Admitting Diagnoses" input
Manufacturer's Model Name	[0008,1090]	AXIOM-Artis
Referenced Study Sequence	[0008,1110]	from RIS
>Referenced SOP Class UID	[0008,1150]	
>Referenced SOP Instance UID	[0008,1155]	
Referenced Patient Sequence	[0008,1120]	from RIS
>Referenced SOP Class UID	[0008,1150]	
>Referenced SOP Instance UID	[0008,1155]	
Referenced Image Sequence	[0008,1140]	Reference to related Plane with BIPLANE A/B Instances
>Referenced SOP Class UID	[0008,1150]	
>Referenced SOP Instance UID	[0008,1155]	
Derivation Description	[0008,2111]	Notes about transformation steps
Start Trim	[0008,2142]	<1st frame to display>
Stop Trim	[0008,2143]	<last frame to display>
Recomm. Display Frame Rate	[0008,2144]	(in f/s)
Patient's Name	[0010,0010]	RIS or "Patient Name" input
Patient ID	[0010,0020]	RIS or "Patient ID" input
Patient's Birth Date	[0010,0030]	RIS or checked "DoB" input
Patient's Sex	[0010,0040]	RIS or input (M or F or O/unknown)
Other Patient IDs	[0010,1000]	from RIS
Other Patient Names	[0010,1001]	from RIS
Patient's Age	[0010,1010]	calculated from "DoB" input
Patient's Size	[0010,1020]	(in meters)
Patient's Weight	[0010,1030]	(in kilograms)
Patient's Address	[0010,1040]	"Address" input
Military Rank	[0010,1080]	from RIS
Medical Alerts	[0010,2000]	from RIS
Contrast Allergies	[0010,2110]	from RIS
Ethnic Group	[0010,2160]	from RIS

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Smoking Status	[0010,21A0]	from RIS
Additional Patient History	[0010,21B0]	from RIS
Pregnancy Status	[0010,21C0]	from RIS
Last Menstrual Date	[0010,21D0]	from RIS
Patient Comments	[0010,4000]	"Additional Info" input
Contrast/Bolus Agent	[0018,0010]	
Cine Rate	[0018,0040]	<acquired frame rate> (in f/s)
KVP	[0018,0060]	<peak KV used> (KV)
Device Serial Number	[0018,1000]	<modality serial number>
Software Versions	[0018,1020]	
Protocol Name	[0018,1030]	"Organ Program" input
Contrast/Bolus Ingredient	[0018,1048]	
Frame Time	[0018,1063]	(msec/frame)
Frame Time Vector	[0018,1065]	<number of Frame values>(msec)
Distance Source to Detector	[0018,1110]	(mm)
Distance Source to Patient	[0018,1111]	(mm)
Estimated Radiographic Magnification Factor	[0018,1114]	<Ratio of SID/SOD>
Exposure Time	[0018,1150]	<duration of x-Ray exposure> (msec)
X-Ray Tube Current	[0018,1151]	(mA)
Average Pulse Width	[0018,1154]	(msec)
Radiation Setting	[0018,1155]	SC GR
Radiation Mode	[0018,115A]	CONTINUOUS PULSED
Image Area Dose Product	[0018,115E]	(dGy*cm*cm)
Intensifier Size	[0018,1162]	<(zoomed) diameter> (mm)
Imager Pixel Spacing	[0018,1164]	<row space, col space>(mm)
Focal Spot	[0018,1190]	
Positioner Motion	[0018,1500]	DYNAMIC STATIC
Positioner Primary Angle	[0018,1510]	(degrees); 0 with DYNAMIC
Positioner Secondary Angle	[0018,1511]	(degrees); 0 with DYNAMIC
Positioner Primary Angle Increment	[0018,1520]	Only with DYNAMIC (absolute values)
Positioner Secondary Angle Increment	[0018,1521]	Only with DYNAMIC (absolute values)
Shutter Shape	[0018,1600]	RECTANGULAR and/or CIRCULAR and/or POLY- GONAL

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Shutter Left Vertical Edge	[0018,1602]	<column number left edge>
Shutter Right Vertical Edge	[0018,1604]	<column number right edge>
Shutter Upper Horizontal Edge	[0018,1606]	<row number upper edge>
Shutter Lower Horizontal Edge	[0018,1608]	<row number lower edge>
Center of Circular Shutter	[0018,1610]	<row number>,<column number>
Radius of Circular Shutter	[0018,1612]	<number of pixels in row direction>
Vertices of the Polygonal Shutter	[0018,1620]	<list of row/column coordinate pairs>
Collimator Shape	[0018,1700]	RECTANGULAR and/or CIRCULAR and/or POLY- GONAL
Collimator Left Vertical Edge	[0018,1702]	<column number left edge>
Collimator Right Vertical Edge	[0018,1704]	<column number right edge>
Collimator Upper Horizontal Edge	[0018,1706]	<row number upper edge>
Collimator Lower Horizontal Edge	[0018,1708]	<row number lower edge>
Center of Circular Collimator	[0018,1710]	<row number>,<column number>
Radius of Circular Collimator	[0018,1712]	<number of pixels in row direction>
Vertices of the Polygonal Collimator	[0018,1720]	<list of row/column coordinate pairs>
Patient Position	[0018,5100]	mandatory "Patient Position" input
Private Creator	[0019,00xx]	SIEMENS SMS-AX VIEW 1.0
<i>Attributes according to "A.2.6 Angio Viewing Module"</i>		
Study Instance UID	[0020,000D]	RIS or system generated
Series Instance UID	[0020,000E]	
Study ID	[0020,0010]	
Series Number	[0020,0011]	
Acquisition Number	[0020,0012]	
Instance Number	[0020,0013]	
Patient Orientation	[0020,0020]	calculated from "Patient Position" input and from Gan- try coordinate data.
Laterality	[0020,0060]	input via L or R marker annotation, else absent
Images in Acquisition	[0020,1002]	
Image Comments	[0020,4000]	not encoded or "SM" or "REF"
Private Creator	[0021,00xx]	SIEMENS SMS-AX ACQ 1.0
<i>Attributes according to "A.2.7 Angio Acquisition Data Module"</i>		
Private Creator	[0023,00xx]	SIEMENS SMS-AX QUANT 1.0

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
<i>Attributes according to "A.2.8 Angio Quantitation Module" if image was calibrated</i>		
Private Creator	[0025,00xx]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0
<i>Attributes according to "A.2.9 Original Image Info Module"</i>		
Samples per Pixel	[0028,0002]	1
Photometric Interpretation	[0028,0004]	MONOCHROME2
Number of Frames	[0028,0008]	
Frame Increment Pointer	[0028,0009]	(0018,1063) or (0018,1065) for VFR
Rows	[0028,0010]	up to 1024
Columns	[0028,0011]	up to 1024
Bits Allocated	[0028,0100]	8 16
Bits Stored	[0028,0101]	8 12
High Bit	[0028,0102]	7 11
Pixel Representation	[0028,0103]	0
Pixel Intensity Relationship	[0028,1040]	LIN LOG DISP
Window Center	[0028,1050]	<NAT value>
Window Width	[0028,1051]	<NAT value>
Recommended Viewing Mode	[0028,1090]	SUB NAT
Lossy Image Compression	[0028,2110]	
Modality LUT Sequence	[0028,3000]	(if [0028,1040] = LOG)
>LUT Descriptor	[0028,3002]	<num of LUT entries>, <first pixel val mapped>, <Entry bits alloc>
>Modality LUT Type	[0028,3004]	US
>LUT Data	[0028,3006]	<array of data, accord. descriptor>
Representative Frame Number	[0028,6010]	<e.g. frame number of Icon Image>
Mask Subtraction Sequence	[0028,6100]	
>Mask Operation	[0028,6101]	AVG_SUB
>Mask Frame Numbers	[0028,6110]	<numbers> (only for AVG_SUB)
Private Creator	[0029,00xx]	CARDIO-D.R. 1.0
Standard Edge Enhancement Sequence	[0029,xx00]	(if applied) contains two items
>Private Creator	[0029,00xx]	CARDIO-D.R. 1.0
>Convolution Kernel Size	[0029,xx01]	<value Pair x-,y-size>
>Convolution kernel Coefficients	[0029,xx02]	<row by row, as specified by size>
>Edge Enhancement Gain	[0029,xx03]	(in %)

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Private Creator	[0029,00xx]	SIEMENS MEDCOM HEADER
<i>Attributes according to "A.2.3 MEDCOM Header Module"</i>		
Requesting Physician	[0032,1032]	from RIS
Requesting Service	[0032,1033]	from RIS
Requested Procedure Description	[0032,1060]	from RIS
Requested Procedure Code Sequence	[0032,1064]	from RIS
>Code Value	[0008,0100]	from RIS
>Coding Scheme Designator	[0008,0102]	from RIS
>Coding Scheme Version	[0008,0103]	from RIS
>Code Meaning	[0008,0104]	from RIS
Requested Contrast Agent	[0032,1070]	from RIS
Study Comments	[0032,4000]	"Exam Comment" input
Special Needs	[0038,0050]	from RIS
Patient State	[0038,0500]	from RIS
Performed Procedure Step Start Date	[0040,0244]	supplied, even if MPPS SOP class is not supported
Performed Procedure Step Start Time	[0040,0245]	supplied, even if MPPS SOP class is not supported
Performed Procedure Step ID	[0040,0253]	supplied, even if MPPS SOP class is not supported, "XAyyyymmddhhmmss" is set with 1st Image acquired
Request Attributes Sequence	[0040,0275]	only if SOP is done under Scheduled Procedure Step, else Sequence is not encoded
>Scheduled Procedure Step Description	[0040,0007]	from RIS
>Scheduled Procedure Step ID	[0040,0009]	from RIS
>Requested Procedure ID	[0040,1001]	from RIS or "Request ID"input
Confidentiality Constraint on Patient Data Description	[0040,3001]	from RIS
Curve Dimensions	[5000,0005]	2
Number of Points	[5000,0010]	<number of data points>
Type of Data	[5000,0020]	ECG
Axis Units	[5000,0030]	DPPS\NONE (DPPS = data points per second)
Data Value Presentation	[5000,0103]	0000h
Curve Data Descriptor	[5000,0110]	0\1
Coordinate Start Value	[5000,0112]	0
Coordinate Step Value	[5000,0114]	<sampling rate>

Table 15: overview of supplied attributes - X-Ray Angiographic image (Original, Derived)

Attribute Name	Tag	Value
Curve Data	[5000,3000]	
Overlay Rows	[60xx,0010]	
Overlay Columns	[60xx,0011]	
Number of Frames in Overlay	[60xx,0015]	<number>
Overlay Description	[60xx,0022]	
Overlay Type	[60xx,0040]	G
Origin	[60xx,0050]	1\1
Image Frame Origin	[60xx,0051]	<1-based number>
Overlay Bits Allocated	[60xx,0100]	1 or 16 (if uncompressed)
Bit Position	[60xx,0102]	0 or 12 (if uncompressed)
Overlay Data	[6000,3000]	<contains Overlay> (not with uncompressed)
Pixel Data	[7FE0,0010]	

A.5 Private Non-Image IOD

The AXIOM Artis system will create numerical data that cannot be correlated to an individual image instance and therefore need to be stored in separate instance(s). This is necessary to correlate the information in the right level of the DICOM data model hierarchy. Since there is no fitting DICOM SOP Class definition, SIEMENS has created a private "Non-Image IOD" to contain numerical data heaps to be managed within a DICOM structure. Please see previous chapters of the Appendix for IOD definition and the following tables for detailed encoding of the different "Non-Image SOP Class Instances".

Table 16: Data Dictionary of private IOD Definition - Study-Report

Attribute Name	Tag	Value
Specific Character Set	[0008,0005]	from Configuration
Image Type	[0008,0008]	ORIGINAL\PRIMARY\SINGLE PLANE\\STUDY REPORT
SOP Class UID	[0008,0016]	1.3.12.2.1107.5.9.1
SOP Instance UID	[0008,0018]	
Study Date	[0008,0020]	<yyyymmdd>
Series Date	[0008,0021]	<yyyymmdd>
Study Time	[0008,0030]	<hhmmss>

Table 16: Data Dictionary of private IOD Definition - Study-Report

Attribute Name	Tag	Value
Series Time	[0008,0031]	<hhmmss>
Accession Number	[0008,0050]	RIS or "Accession Number" input
Modality	[0008,0060]	XA
Manufacturer	[0008,0070]	Siemens
Institution Name	[0008,0080]	
Referring Physician's Name	[0008,0090]	RIS or input
Study Description	[0008,1030]	"Study" input
Series Description	[0008,103E]	STUDY REPORT
Performing Physicians' Name	[0008,1050]	"Performing Physician 1" \ "Performing Physician 2" input
Operator's Name	[0008,1070]	"Operator 1" \ "Operator 2" input
Admitting Diagnoses Description	[0008,1080]	RIS or "Admitting Diagnoses" input
Manufacturer's Model Name	[0008,1090]	<Product Name>
Patient's Name	[0010,0010]	RIS or "Patient Name" input
Patient ID	[0010,0020]	RIS or "Patient ID" input
Patient's Birth Date	[0010,0030]	RIS or checked input
Patient's Sex	[0010,0040]	RIS or input (M or F or O/unknown)
Patient's Age	[0010,1010]	
Patient's Size	[0010,1020]	(in meters)
Patient's Weight	[0010,1030]	(in kilograms)
Patient Comments	[0010,4000]	"Additional Info" input
Protocol Name	[0018,1030]	STUDY REPORT
Study Instance UID	[0020,000D]	from RIS or system generated
Series Instance UID	[0020,000E]	
Study ID	[0020,0010]	
Series Number	[0020,0011]	
Private Creator	[0029,00xx]	SIEMENS CSA NON-IMAGE
Data Type	[0029,xx08]	BSR REPORT
Data Version	[0029,xx09]	1.0
Requested Procedure Description	[0032,1060]	from RIS
Study Comments	[0032,4000]	"Exam Comment" input
Private Creator	[7FE1,00xx]	SIEMENS CSA NON-IMAGE

Table 16: Data Dictionary of private IOD Definition - Study-Report

Attribute Name	Tag	Value
Data	[7FE1,xx10]	<Study Report data> (encoded as 8-Bit ASCII data, with <LF> as line separator)

A.6 Quant Report Images - derived XA IOD

The AXIOM Artis System will create result images form performing Qunatitative Analysis Functions. To ensure image interchange, the resulting reports are converted to XA multi-frames with dedicated header data.

Table 17: overview of supplied attributes - X-Ray Angiographic image (Quant Report Image)

Attribute Name	Tag	Value
Specific Character Set	[0008,0005]	from Configuration
Image Type	[0008,0008]	DERIVED\SECONDARY\<3rd to nth value from Original>\QUANT
SOP Class UID	[0008,0016]	1.2.840.10008.5.1.4.1.1.12.1
SOP Instance UID	[0008,0018]	
Study Date	[0008,0020]	<yyyymmdd>
Series Date	[0008,0021]	<yyyymmdd> (Quant Series)
Acquisition Date	[0008,0022]	Date of Original Acquisition
Content Date	[0008,0023]	<yyyymmdd> (Date of Creation)
Study Time	[0008,0030]	<hhmmss>
Series Time	[0008,0031]	<hhmmss> (Quant Series)
Acquisition Time	[0008,0032]	Time of Original Acquisition
Content Time	[0008,0033]	<hhmmss> (Time of Creation)
Accession Number	[0008,0050]	from Original
Modality	[0008,0060]	XA
Manufacturer	[0008,0070]	Siemens
Institution Name	[0008,0080]	from Original
Referring Physician's Name	[0008,0090]	from Original
Station Name	[0008,1010]	from Configuration (VB22x and later)
Study Description	[0008,1030]	from Original
Series Description	[0008,103E]	input from Analysis

Table 17: overview of supplied attributes - X-Ray Angiographic image (Quant Report Image)

Attribute Name	Tag	Value
Performing Physicians' Name	[0008,1050]	from Original
Operator's Name	[0008,1070]	from Original
Admitting Diagnoses Description	[0008,1080]	from Original
Manufacturer's Model Name	[0008,1090]	<Product Name>
Referenced Patient Sequence	[0008,1120]	from Original
>Referenced SOP Class UID	[0008,1150]	
>Referenced SOP Instance UID	[0008,1155]	
Start Trim	[0008,2142]	<1st frame to display>
Stop Trim	[0008,2143]	<last frame to display>
Recomm. Display Frame Rate	[0008,2144]	1
Patient's Name	[0010,0010]	from Original
Patient ID	[0010,0020]	from Original
Patient's Birth Date	[0010,0030]	from Original
Patient's Sex	[0010,0040]	from Original
Patient's Age	[0010,1010]	from Original
Patient's Size	[0010,1020]	(in meters)
Patient's Weight	[0010,1030]	(in kilograms)
Medical Alerts	[0010,2000]	from Original
Contrast Allergies	[0010,2110]	from Original
Pregnancy Status	[0010,21C0]	from Original
Patient Comments	[0010,4000]	from Original
Contrast/Bolus Agent	[0018,0010]	from Original
Cine Rate	[0018,0040]	1
KVP	[0018,0060]	from Original
Device Serial Number	[0018,1000]	<modality serial number>
Software Versions	[0018,1020]	
Protocol Name	[0018,1030]	from Original
Contrast/Bolus Ingredient	[0018,1048]	from Original
Frame Time	[0018,1063]	1000
Distance Source to Detector	[0018,1110]	from Original
Distance Source to Patient	[0018,1111]	from Original
Exposure Time	[0018,1150]	<duration of x-Ray exposure> (msec)
X-Ray Tube Current	[0018,1151]	from Original

Table 17: overview of supplied attributes - X-Ray Angiographic image (Quant Report Image)

Attribute Name	Tag	Value
Average Pulse Width	[0018,1154]	from Original
Radiation Setting	[0018,1155]	from Original
Radiation Mode	[0018,115A]	from Original
Image Area Dose Product	[0018,115E]	from Original
Intensifier Size	[0018,1162]	from Original
Focal Spot	[0018,1190]	from Original
Positioner Motion	[0018,1500]	from Original
Positioner Primary Angle	[0018,1510]	from Original
Positioner Secondary Angle	[0018,1511]	from Original
Positioner Primary Angle Increment	[0018,1520]	from Original (Only with DYNAMIC)
Positioner Secondary Angle Increment	[0018,1521]	from Original (Only with DYNAMIC)
Patient Position	[0018,5100]	from Original
Private Creator	[0019,00xx]	SIEMENS SMS-AX VIEW 1.0
<i>Attributes according to "A.2.6 Angio Viewing Module"</i>		
Study Instance UID	[0020,000D]	from Original
Series Instance UID	[0020,000E]	
Study ID	[0020,0010]	from Original
Series Number	[0020,0011]	
Acquisition Number	[0020,0012]	from Original
Instance Number	[0020,0013]	
Patient Orientation	[0020,0020]	from Original
Images in Acquisition	[0020,1002]	
Private Creator	[0021,00xx]	SIEMENS SMS-AX ACQ 1.0
<i>Attributes according to "A.2.7 Angio Acquisition Data Module"</i>		
Private Creator	[0023,00xx]	SIEMENS SMS-AX QUANT 1.0
<i>Attributes according to "A.2.8 Angio Quantitation Module" if image was calibrated</i>		
Private Creator	[0025,00xx]	SIEMENS SMS-AX ORIGINAL IMAGE INFO 1.0
<i>Attributes according to "A.2.9 Original Image Info Module"</i>		
Samples per Pixel	[0028,0002]	1
Photometric Interpretation	[0028,0004]	MONOCHROME2
Number of Frames	[0028,0008]	1 to 4
Frame Increment Pointer	[0028,0009]	(0018,1063)
Rows	[0028,0010]	1024

Table 17: overview of supplied attributes - X-Ray Angiographic image (Quant Report Image)

Attribute Name	Tag	Value
Columns	[0028,0011]	1024
Bits Allocated	[0028,0100]	16
Bits Stored	[0028,0101]	12
High Bit	[0028,0102]	11
Pixel Representation	[0028,0103]	0
Pixel Intensity Relationship	[0028,1040]	LIN
Window Center	[0028,1050]	2047
Window Width	[0028,1051]	4095
Lossy Image Compression	[0028,2110]	from Original
Representative Frame Number	[0028,6010]	1
Study Status ID	[0032,000A]	from Original
Scheduled Study Start Date	[0032,1000]	from Original
Scheduled Study Start Time	[0032,1001]	from Original
Scheduled Study Location	[0032,1020]	from Original
Scheduled Study Location AE Title(s)	[0032,1021]	from Original
Requesting Physician	[0032,1032]	from Original
Requested Procedure Description	[0032,1060]	from Original
Requested Contrast Agent	[0032,1070]	from Original
Study Comments	[0032,4000]	from Original
Performed Procedure Step Start Date	[0040,0244]	from Original
Performed Procedure Step Start Time	[0040,0245]	from Original
Performed Procedure Step ID	[0040,0253]	from Original
Request Attributes Sequence	[0040,0275]	from Original
>Scheduled Procedure Step Description	[0040,0007]	from Original
>Scheduled Procedure Step ID	[0040,0009]	from Original
>Requested Procedure ID	[0040,1001]	from Original
Confidentiality Constraint on Patient Data Description	[0040,3001]	from Original
Pixel Data	[7FE0,0010]	

A.7 CT/MR image converted to XA IOD (VA22x and

later)

The Artis system will create reference images from CT or MR images. Those will be encoded as derived XA Standard Extended SOP Class. Please see the following tables for a complete overview of supplied Type 1/2/3 Standard and additional Private Attributes:

Table 18: overview of supplied attributes - XA image (derived) from CT image

Attribute Name	Tag	Value
Specific Character Set	[0008,0005]	from Configuration
Image Type	[0008,0008]	see "A.2.10 SOP Common Module - Image Type Extensions" (Ref Image)
SOP Class UID	[0008,0016]	1.2.840.10008.5.1.4.1.1.12.1
SOP Instance UID	[0008,0018]	<new UID>
Study Date	[0008,0020]	from Original or changed by CorRea (Drag&Drop)
Series Date	[0008,0021]	<yyyymmdd>
Acquisition Date	[0008,0022]	Date of Original Acquisition, if provided by CT
Content Date	[0008,0023]	<yyyymmdd> (Date of Creation)
Study Time	[0008,0030]	from Original or changed by CorRea (Drag&Drop)
Series Time	[0008,0031]	<hhmmss>
Acquisition Time	[0008,0032]	Time of Original Acquisition, if provided by CT
Content Time	[0008,0033]	<hhmmss> (Time of Creation)
Accession Number	[0008,0050]	from Original or changed by CorRea (Drag&Drop)
Modality	[0008,0060]	XA
Manufacturer	[0008,0070]	Siemens
Institution Name	[0008,0080]	from Original (equipment)
Referring Physician's Name	[0008,0090]	from Original or changed by CorRea (Drag&Drop)
Station Name	[0008,1010]	Artis Station Name (VB22x and later)
Study Description	[0008,1030]	from Original or changed by CorRea (Drag&Drop)
Series Description	[0008,103E]	from Original
Institutional Department Name	[0008,1040]	from Original (equipment)
Performing Physicians' Name	[0008,1050]	from Original (series)
Operator's Name	[0008,1070]	from Original (series)
Admitting Diagnoses Description	[0008,1080]	from Original or changed by CorRea (Drag&Drop)
Manufacturer's Model Name	[0008,1090]	AXIOM-Artis
Source Image Sequence	[0008,2112]	Reference to original CT image
>Referenced SOP Class UID	[0008,1150]	1.2.840.10008.5.1.4.1.1.2 (CT image)

Table 18: overview of supplied attributes - XA image (derived) from CT image

Attribute Name	Tag	Value
>Referenced SOP Instance UID	[0008,1155]	from Original [0008,0018]
Patient's Name	[0010,0010]	from Original or changed by CorRea (Drag&Drop)
Patient ID	[0010,0020]	from Original or changed by CorRea (Drag&Drop)
Patient's Birth Date	[0010,0030]	from Original or changed by CorRea (Drag&Drop)
Patient's Sex	[0010,0040]	from Original or changed by CorRea (Drag&Drop)
Patient's Age	[0010,1010]	from Original or changed by CorRea (Drag&Drop)
Patient's Size	[0010,1020]	from Original or changed by CorRea (Drag&Drop)
Patient's Weight	[0010,1030]	from Original or changed by CorRea (Drag&Drop)
Medical Alerts	[0010,2000]	from Original
Contrast Allergies	[0010,2110]	from Original
Pregnancy Status	[0010,21C0]	from Original
Patient Comments	[0010,4000]	from Original or changed by CorRea (Drag&Drop)
Contrast/Bolus Agent	[0018,0010]	from Original (Contrast/bolus)
KVP	[0018,0060]	from Original (image)
Device Serial Number	[0018,1000]	<Artis Serial Number>
Software Versions	[0018,1020]	<Artis Software Version>
Protocol Name	[0018,1030]	from Original "Series Description"
Contrast/Bolus Ingredient	[0018,1048]	from Original (Contrast/bolus)
Estimated Radiographic Magnification Factor	[0018,1114]	1.0
Exposure Time	[0018,1150]	from Original (CT image)
X-Ray Tube Current	[0018,1151]	from Original (CT image)
Radiation Setting	[0018,1155]	GR
Image Area Dose Product	[0018,115E]	<zero length>
Imager Pixel Spacing	[0018,1164]	values taken from original "[0028,0030] Pixel Spacing" (Image plane)
Positioner Primary Angle	[0018,1510]	<zero length>
Positioner Secondary Angle	[0018,1511]	<zero length>
Private Creator	[0019,00xx]	SIEMENS SMS-AX VIEW 1.0
<i>Attributes according to "A.2.6 Angio Viewing Module" (created, if required)</i>		
Study Instance UID	[0020,000D]	from Original or changed by CorRea (Drag&Drop)
Series Instance UID	[0020,000E]	<new UID>
Study ID	[0020,0010]	from Original or changed by CorRea (Drag&Drop)
Series Number	[0020,0011]	"1000-based" counter

Table 18: overview of supplied attributes - XA image (derived) from CT image

Attribute Name	Tag	Value
Acquisition Number	[0020,0012]	from original (image)
Instance Number	[0020,0013]	"2000-based" counter
Patient Orientation	[0020,0020]	from Original (series)
Images in Acquisition	[0020,1002]	
Image Comments	[0020,4000]	from Original (Image)
Private Creator	[0021,00xx]	SIEMENS SMS-AX ACQ 1.0
<i>Attributes according to "A.2.7 Angio Acquisition Data Module" (created, if required)</i>		
Samples per Pixel	[0028,0002]	1
Photometric Interpretation	[0028,0004]	MONOCHROME2
Rows	[0028,0010]	from Original (Image) <i>only 512 supported</i>
Columns	[0028,0011]	from Original (Image) <i>only 512 supported</i>
Bits Allocated	[0028,0100]	16
Bits Stored	[0028,0101]	12 (converted)
High Bit	[0028,0102]	11 (converted)
Pixel Representation	[0028,0103]	0 (converted to unsigned)
Pixel Intensity Relationship	[0028,1040]	LIN or DISP
Window Center	[0028,1050]	calculated from Original Window Center and Rescale Intercept (VOI LUT, MOD LUT)
Window Width	[0028,1051]	from Original (VOI LUT)
Lossy Image Compression	[0028,2110]	from Original (Image)
Study Status ID	[0032,000A]	from Original
Scheduled Study Start Date	[0032,1000]	from Original
Scheduled Study Start Time	[0032,1001]	from Original
Scheduled Study Location	[0032,1020]	from Original
Scheduled Study Location AE Title(s)	[0032,1021]	from Original
Requesting Physician	[0032,1032]	from Original
Requesting Service	[0032,1033]	from Original
Requested Procedure Description	[0032,1060]	from Original
Requested Contrast Agent	[0032,1070]	from Original
Study Comments	[0032,4000]	from Original or changed by CorRea (Drag&Drop)
Performed Procedure Step Start Date	[0040,0244]	from Original (Series), if existent
Performed Procedure Step Start Time	[0040,0245]	from Original (Series), if existent
Performed Procedure Step ID	[0040,0253]	from Original (Series), if existent

Table 18: overview of supplied attributes - XA image (derived) from CT image

Attribute Name	Tag	Value
Request Attributes Sequence	[0040,0275]	from Original (Series), if existent
>item contents as provided		as set from Original (Series)
Overlay Rows	[60xx,0010]	same as [0028,0010]
Overlay Columns	[60xx,0011]	same as [0028,0011]
Overlay Description	[60xx,0022]	from Original (Overlay)
Overlay Type	[60xx,0040]	from Original (Overlay)
Origin	[60xx,0050]	1\from Original (Overlay)
Overlay Bits Allocated	[60xx,0100]	from Original (Overlay)
Bit Position	[60xx,0102]	from Original (Overlay)
Overlay Data	[6000,3000]	from Original (Overlay)
Pixel Data	[7FE0,0010]	possibly converted from original pixel data

Table 19: overview of supplied attributes - XA image (derived) from MR image

Attribute Name	Tag	Value
Specific Character Set	[0008,0005]	from Configuration
Image Type	[0008,0008]	see "A.2.10 SOP Common Module - Image Type Extensions" (Ref Image)
SOP Class UID	[0008,0016]	1.2.840.10008.5.1.4.1.1.12.1
SOP Instance UID	[0008,0018]	<new UID>
Study Date	[0008,0020]	from Original or changed by CorRea (Drag&Drop)
Series Date	[0008,0021]	<yyyymmdd>
Acquisition Date	[0008,0022]	Date of Original Acquisition, if provided by MR
Content Date	[0008,0023]	<yyyymmdd> (Date of Creation)
Study Time	[0008,0030]	from Original or changed by CorRea (Drag&Drop)
Series Time	[0008,0031]	<hhmmss>
Acquisition Time	[0008,0032]	Time of Original Acquisition, if provided by CT
Content Time	[0008,0033]	<hhmmss> (Time of Creation)
Accession Number	[0008,0050]	from Original or changed by CorRea (Drag&Drop)
Modality	[0008,0060]	XA
Manufacturer	[0008,0070]	Siemens
Institution Name	[0008,0080]	from Original (equipment)

Table 19: overview of supplied attributes - XA image (derived) from MR image

Attribute Name	Tag	Value
Referring Physician's Name	[0008,0090]	from Original or changed by CorRea (Drag&Drop)
Station Name	[0008,1010]	Artis Station Name (VB22x and later)
Study Description	[0008,1030]	from Original or changed by CorRea (Drag&Drop)
Series Description	[0008,103E]	from Original
Institutional Department Name	[0008,1040]	from Original (equipment)
Performing Physicians' Name	[0008,1050]	from Original (series)
Operator's Name	[0008,1070]	from Original (series)
Admitting Diagnoses Description	[0008,1080]	from Original or changed by CorRea (Drag&Drop)
Manufacturer's Model Name	[0008,1090]	AXIOM-Artis
Source Image Sequence	[0008,2112]	Reference to original MR image
>Referenced SOP Class UID	[0008,1150]	1.2.840.10008.5.1.4.1.1.4 (MR image)
>Referenced SOP Instance UID	[0008,1155]	from Original [0008,0018]
Patient's Name	[0010,0010]	from Original or changed by CorRea (Drag&Drop)
Patient ID	[0010,0020]	from Original or changed by CorRea (Drag&Drop)
Patient's Birth Date	[0010,0030]	from Original or changed by CorRea (Drag&Drop)
Patient's Sex	[0010,0040]	from Original or changed by CorRea (Drag&Drop)
Patient's Age	[0010,1010]	from Original or changed by CorRea (Drag&Drop)
Patient's Size	[0010,1020]	from Original or changed by CorRea (Drag&Drop)
Patient's Weight	[0010,1030]	from Original or changed by CorRea (Drag&Drop)
Medical Alerts	[0010,2000]	from Original
Contrast Allergies	[0010,2110]	from Original
Pregnancy Status	[0010,21C0]	from Original
Patient Comments	[0010,4000]	from Original or changed by CorRea (Drag&Drop)
Contrast/Bolus Agent	[0018,0010]	from Original (Contrast/bolus)
KVP	[0018,0060]	<zero length>
Device Serial Number	[0018,1000]	<Artis Serial Number>
Software Versions	[0018,1020]	<Artis Software Version>
Protocol Name	[0018,1030]	from Original "Series Description"
Contrast/Bolus Ingredient	[0018,1048]	from Original (Contrast/bolus)
Estimated Radiographic Magnification Factor	[0018,1114]	1.0
Exposure Time	[0018,1150]	from Original (CT image)
X-Ray Tube Current	[0018,1151]	from Original (CT image)
Radiation Setting	[0018,1155]	SC (even if created w/o radiation!)

Table 19: overview of supplied attributes - XA image (derived) from **MR** image

Attribute Name	Tag	Value
Image Area Dose Product	[0018,115E]	<zero length>
Imager Pixel Spacing	[0018,1164]	values taken from original "[0028,0030] Pixel Spacing" (Image plane)
Positioner Primary Angle	[0018,1510]	<zero length>
Positioner Secondary Angle	[0018,1511]	<zero length>
Private Creator	[0019,00xx]	SIEMENS SMS-AX VIEW 1.0
<i>Attributes according to "A.2.6 Angio Viewing Module" (created, if required)</i>		
Study Instance UID	[0020,000D]	from Original or changed by CorRea (Drag&Drop)
Series Instance UID	[0020,000E]	<new UID>
Study ID	[0020,0010]	from Original or changed by CorRea (Drag&Drop)
Series Number	[0020,0011]	"1000-based" counter
Acquisition Number	[0020,0012]	from original (image)
Instance Number	[0020,0013]	"2000-based" counter
Patient Orientation	[0020,0020]	from Original (series)
Images in Acquisition	[0020,1002]	
Image Comments	[0020,4000]	from Original (Image)
Private Creator	[0021,00xx]	SIEMENS SMS-AX ACQ 1.0
<i>Attributes according to "A.2.7 Angio Acquisition Data Module" (created, if required)</i>		
Samples per Pixel	[0028,0002]	1
Photometric Interpretation	[0028,0004]	MONOCHROME2
Rows	[0028,0010]	from Original (Image) <i>only 512 supported</i>
Columns	[0028,0011]	from Original (Image) <i>only 512 supported</i>
Bits Allocated	[0028,0100]	16
Bits Stored	[0028,0101]	12 (converted)
High Bit	[0028,0102]	11 (converted)
Pixel Representation	[0028,0103]	0 (converted to unsigned)
Pixel Intensity Relationship	[0028,1040]	LIN or DISP
Window Center	[0028,1050]	calculated from Original Window Center and Rescale Intercept (VOI LUT, MOD LUT)
Window Width	[0028,1051]	from Original (VOI LUT)
Lossy Image Compression	[0028,2110]	from Original (Image)
Study Status ID	[0032,000A]	from Original
Scheduled Study Start Date	[0032,1000]	from Original

Table 19: overview of supplied attributes - XA image (derived) from MR image

Attribute Name	Tag	Value
Scheduled Study Start Time	[0032,1001]	from Original
Scheduled Study Location	[0032,1020]	from Original
Scheduled Study Location AE Title(s)	[0032,1021]	from Original
Requesting Physician	[0032,1032]	from Original
Requested Procedure Description	[0032,1060]	from Original
Requested Contrast Agent	[0032,1070]	from Original
Study Comments	[0032,4000]	from Original or changed by CorRea (Drag&Drop)
Performed Procedure Step Start Date	[0040,0244]	from Original (Series), if existent
Performed Procedure Step Start Time	[0040,0245]	from Original (Series), if existent
Performed Procedure Step ID	[0040,0253]	from Original (Series), if existent
Request Attributes Sequence	[0040,0275]	from Original (Series), if existent
>item contents as provided		as set from Original (Series)
Overlay Rows	[60xx,0010]	same as [0028,0010]
Overlay Columns	[60xx,0011]	same as [0028,0011]
Overlay Description	[60xx,0022]	from Original (Overlay)
Overlay Type	[60xx,0040]	from Original (Overlay)
Origin	[60xx,0050]	1\from Original (Overlay)
Overlay Bits Allocated	[60xx,0100]	from Original (Overlay)
Bit Position	[60xx,0102]	from Original (Overlay)
Overlay Data	[6000,3000]	from Original (Overlay)
Pixel Data	[7FE0,0010]	possibly converted from original pixel data

 Copyright © Siemens AG 2004. All rights reserved.
 Alle Rechte vorbehalten.

A.8 Sorting Order AXIOM Artis (VA10C/VA20F/VA22x) with syngo VA47Z/VA49B

This section discusses the sorting algorithms for Artis within the patient browser and directory overview.

A.8.1 Identification of Images

Original acquired images are identified by the following attribute contents:

Attribute Name	Tag	Value
Image Type	[0008,0008]	ORIGINAL\PRIMARY\SINGLE PLANE\SINGLE A ORIGINAL\PRIMARY\SINGLE PLANE\SINGLE B ORIGINAL\PRIMARY\BIPLANE A ORIGINAL\PRIMARY\BIPLANE B ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\PERI\MASK/n ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\PERI\FILL/n ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\MASK ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\INJECTION ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\FILL ORIGINAL\PRIMARY\SINGLE PLANE\PLANE A\DYNA\WASHOUT DERIVED\PRIMARY\SINGLE PLANE\SINGLE A\STORE FLUORO DERIVED\PRIMARY\SINGLE PLANE\SINGLE B\STORE FLUORO DERIVED\PRIMARY\BIPLANE A\STORE FLUORO DERIVED\PRIMARY\BIPLANE B\STORE FLUORO
Series Time	[0008,0031]	Same as Acquisition Time.
Acquisition Time	[0008,0032]	Time the acquisition of this image was performed
Content Time	[0008,0033]	Same as Acquisition Time (rule for acquired images)
Series Number	[0020,0011]	Same as Acquisition Number
Acquisition Number	[0020,0012]	Counts number of acquisitions. All images created during performance of this acquisition have the same number.
Instance Number	[0020,0013]	"1-based" counter for images. Indicates the creation order within this acquisition.

Store Monitor images, derived from original images or Fluoro LIH, are identified by the following attribute contents:

Attribute Name	Tag	Value
Image Type	[0008,0008]	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A\STORE MONITOR DERIVED\SECONDARY\SINGLE PLANE\SINGLE B\STORE MONITOR
Series Time	[0008,0031]	Same as Acquisition Time. Store Monitor is in same series.
Acquisition Time	[0008,0032]	Time the acquisition of the original image was performed
Content Time	[0008,0033]	Time this image was created (Store Monitor was performed)
Series Number	[0020,0011]	Same Number as Acquisition Number. (Reference to original Image)
Acquisition Number	[0020,0012]	Counts number of acquisitions. Refers to the acquisition number of the original image where this image was derived from.
Instance Number	[0020,0013]	"300-based" counter for images. Indicates the creation order of Store Monitor image within this series.

Store Reference images, derived from original images or Fluoro LIH, are identified by the following attribute contents:

Attribute Name	Tag	Value
Image Type	[0008,0008]	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A\REFIMAGE DERIVED\SECONDARY\SINGLE PLANE\SINGLE B\REFIMAGE
Series Time	[0008,0031]	Same as Acquisition Time. Store Reference image is in same series.
Acquisition Time	[0008,0032]	Time the acquisition of the original image was performed
Content Time	[0008,0033]	Time this image was created (Store Reference image was performed)
Series Number	[0020,0011]	Same Number as Acquisition Number. (Reference to original Image)
Acquisition Number	[0020,0012]	Counts number of acquisitions. Refers to the acquisition number of the original image where this image was derived from.
Instance Number	[0020,0013]	"100-based" counter for images. Indicates the creation order of Store References image within this series.

Quant Report images, derived from original images or Fluoro LIH, are identified by the following attribute contents:

Attribute Name	Tag	Value
Image Type	[0008,0008]	DERIVED\SECONDARY\SINGLE PLANE\SINGLE A\QUANT DERIVED\SECONDARY\SINGLE PLANE\SINGLE B\QUANT
Series Time	[0008,0031]	Time this Quant image was created (Quant was performed)
Acquisition Time	[0008,0032]	Time the acquisition of the original image was performed
Content Time	[0008,0033]	Time this image was created (Quant was performed)
Series Number	[0020,0011]	Same Number as Acquisition Number. (Reference to original Image)
Acquisition Number	[0020,0012]	Counts number of acquisitions. Refers to the acquisition number of the original image where this image was derived from.
Instance Number	[0020,0013]	"1-based" counter for images. Quant Reports are stored in new series.

A Study Report non-image, created once per study, is identified by the following attribute contents:

Attribute Name	Tag	Value
Image Type	[0008,0008]	ORIGINAL\PRIMARY\SINGLE PLANE\STUDY REPORT
SOP Class UID	[0008,0016]	1.3.12.2.1107.5.9.1
Series Time	[0008,0031]	Study Report is allocated with Patient Registration.
Acquisition Time	[0008,0032]	Not set
Content Time	[0008,0033]	Not set
Series Number	[0020,0011]	Fixed value of "999"
Acquisition Number	[0020,0012]	Not set

Attribute Name	Tag	Value
Instance Number	[0020,0013]	Not set

A.8.2 Handling of Performed Procedure Step ID

Each image contains the attribute "Performed Procedure Step ID" (0040,0235), based on a "YYMMDDHHMMSS" encoding. This date and time is based on the time when the first image is acquired. The "Performed Procedure Step ID" stays the same for all acquired or derived images as long as the patient is re-registered. A re-registered patient with a new study or new series within the existing study will get a newly assigned "Performed Procedure Step ID".

A.8.3 The sorting algorithms of Artis Viewer

AXIOM Artis has implemented sorting algorithms to display DICOM objects in the directory overview in a specific order. Sorting/Grouping in the Viewer Directory Overview occurs according the following rules:

1. Sort by Series Number
2. If Series Number matches: Group on Series Number, Sort by Content Number
3. If Series Number and Content Number each match: Group on Series Number, Sort by Content Time

In the case, that the user has selected "Scenes + Ref Images" or "Scenes + Store Monitor", the reference images or the store monitor images follow directly the image they had been derived from.

Example 1a:

Patient "Test 1" with one Study, three Series with one acquired image, three Store Monitor Images and one Reference Image in each Series. "Scenes" is selected under the "Dir Overview" button

Patient "Test 1"

Study 1 Series 1 Image 1	Study 1 Series 2 Image 1	Study 1 Series 3 Image 1	

	Acquired Image
	Store Monitor Image
	Store Reference Image

Example 1b:

Patient "Test 1" with one Study, three Series with one acquired image, three Store Monitor Images and one Reference Image in each Series. "Scenes + Ref Images" is selected under the "Dir Overview" button

Patient "Test 1"

Study 1 Series 1 Image 1	Study 1 Series 1 Image 100	Study 1 Series 2 Image 1	Study 1 Series 2 Image 100
Study 1 Series 3 Image 1	Study 1 Series 3 Image 100		

	Acquired Image
	Store Monitor Image
	Store Reference Image

Example 1c:

Patient "Test 1" with one Study, three Series with one acquired image, three Store Monitor Images and one Reference Image in each Series. "Scenes + Store Monitor" is selected under the "Dir Overview" button

Patient "Test 1"

Study 1 Series 1 Image 1	Study 1 Series 1 Image 300	Study 1 Series 1 Image 301	Study 1 Series 1 Image 302
Study 1 Series 2 Image 1	Study 1 Series 2 Image 300	Study 1 Series 2 Image 301	Study 1 Series 2 Image 302
Study 1 Series 3 Image 1	Study 1 Series 3 Image 300	Study 1 Series 3 Image 301	Study 1 Series 3 Image 302

	Acquired Image
	Store Monitor Image
	Store Reference Image

Example 1d:

For Scene #2 of Patient "Test 1" now a Quant Report is performed. According to Example 1c, the Overview now sorts as follows:

Patient "Test 1"

Study 1 Series 1 Image 1	Study 1 Series 1 Image 300	Study 1 Series 1 Image 301	Study 1 Series 1 Image 302
--------------------------------	----------------------------------	----------------------------------	----------------------------------

Study 1 Series 2 Image 1	Study 1 Series 2(Quant) Image 1(Quant) Content Time >	Study 1 Series 2 Image 300	Study 1 Series 2 Image 301
Study 1 Series 2 Image 302	Study 1 Series 3 Image 1	Study 1 Series 3 Image 300	Study 1 Series 3 Image 301
Study 1 Series 3 Image 302			

	Acquired Image
	Store Monitor Image
	Store Reference Image
	Quant Report Image

A.8.4 The sorting algorithms of Artis Browser

Besides the default sorting criteria, Artis supports the following sorting algorithms for the patient browser.

The default syngo sorting settings are as follows

- Artis displays images - within a Series - with the order they had been acquired
- Artis displays Series - within a Study - with the order they had been acquired
- Artis displays Studies - within a Patient - with the order they had been acquired

Therefore the Quant Report Series will not necessarily be displayed adjacent to the original scene, if created later at the end of all scenes.

A.9 DICOM Print SCU - detailed status displays

The following tables document the behavior of the LEONARDO DICOM Print AE in response to messages received for the printer SOP class and the print job SOP class.

A.9.1 Common Status Information

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
NORMAL	Camera is ready.	Camera is ready.	<None>/idle
BAD RECEIVE MGZ	There is a problem with the film receive magazine. Films from the printer cannot be transported into the magazine.	Problem with receive magazine.	<None>/interact
BAD SUPPLY MGZ	There is a problem with a film supply magazine. Films from this magazine cannot be transported into the printer.	Problem with supply magazine.	<None>/interact
CALIBRATING	Printer is performing self calibration, it is expected to be available for normal operation shortly.	Self calibration. Please wait.	<None>/idle
CALIBRATION ERR	An error in the printer calibration has been detected, quality of processed films may not be optimal.	Problem in calibration. Film quality may not be optimal.	<None>/interact
CHECK CHEMISTRY	A problem with the processor chemicals has been detected, quality of processed films may not be optimal.	Problem with chemistry. Film quality may not be optimal.	<None>/interact
CHECK SORTER	There is an error in the film sorter	Error in film sorter.	<None>/interact
CHEMICALS EMPTY	There are no processing chemicals in the processor, films will not be printed and processed until the processor is back to normal.	Camera chemistry empty. Please check.	<None>/interact
CHEMICALS LOW	The chemical level in the processor is low, if not corrected, it will probably shut down soon.	Camera chemistry low. Please check.	<None>/interact
COVER OPEN	One or more printer or processor covers, drawers, doors are open.	Camera cover, drawer or door open.	<None>/interact
ELEC CONFIG ERR	Printer configured improperly for this job.	Camera configured improperly for this job. Queue stopped.	Queue for this camera will be STOPPED/Queue stopped

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
ELEC DOWN	Printer is not operating due to some unspecified electrical hardware problem.	Camera electrical hardware problem.	<None>/interact
ELEC SW ERROR	Printer not operating for some unspecified software error.	Camera software problem. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
EMPTY 8x10	The 8x10 inch film supply magazine is empty.	8x10 film supply empty.	<None>/interact
EMPTY 8x10 BLUE	The 8x10 inch blue film supply magazine is empty.	8x10 blue film supply empty.	<None>/interact
EMPTY 8x10 CLR	The 8x10 inch clear film supply magazine is empty.	8x10 clear film supply empty.	<None>/interact
EMPTY 8x10 PAPR	The 8x10 inch paper supply magazine is empty.	8x10 paper supply empty.	<None>/interact
EMPTY 10x12	The 10x12 inch film supply magazine is empty.	10x12 film supply empty.	<None>/interact
EMPTY 10x12 BLUE	The 10x12 inch blue film supply magazine is empty.	10x12 blue film supply empty.	<None>/interact
EMPTY 10x12 CLR	The 10x12 inch clear film supply magazine is empty.	10x12 clear film supply empty.	<None>/interact
EMPTY 10x12 PAPR	The 10x12 inch paper supply magazine is empty.	10x12 paper supply empty.	<None>/interact
EMPTY 10x14	The 10x14 inch film supply magazine is empty.	10x14 film supply empty.	<None>/interact
EMPTY 10x14 BLUE	The 10x14 inch blue film supply magazine is empty.	10x14 blue film supply empty.	<None>/interact
EMPTY 10x14 CLR	The 10x14 inch clear film supply magazine is empty.	10x14 clear film supply empty.	<None>/interact
EMPTY 10x14 PAPR	The 10x14 inch paper supply magazine is empty.	10x14 paper supply empty.	<None>/interact
EMPTY 11x14	The 11x14 inch film supply magazine is empty.	11x14 film supply empty.	<None>/interact
EMPTY 11x14 BLUE	The 11x14 inch blue film supply magazine is empty.	11x14 blue film supply empty.	<None>/interact
EMPTY 11x14 CLR	The 11x14 inch clear film supply magazine is empty.	11x14 clear film supply empty.	<None>/interact
EMPTY 11x14 PAPR	The 11x14 inch paper supply magazine is empty.	11x14 paper supply empty.	<None>/interact
EMPTY 14x14	The 14x14 inch film supply magazine is empty.	14x14 film supply empty.	<None>/interact

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
EMPTY 14x14 BLUE	The 14x14 inch blue film supply magazine is empty.	14x14 blue film supply empty.	<None>/interact
EMPTY 14x14 CLR	The 14x14 inch clear film supply magazine is empty.	14x14 clear film supply empty.	<None>/interact
EMPTY 14x14 PAPER	The 14x14 inch paper supply magazine is empty.	14x14 paper supply empty.	<None>/interact
EMPTY 14x17	The 14x17 inch film supply magazine is empty.	14x17 film supply empty.	<None>/interact
EMPTY 14x17 BLUE	The 14x17 inch blue film supply magazine is empty.	14x17 blue film supply empty.	<None>/interact
EMPTY 14x17 CLR	The 14x17 inch clear film supply magazine is empty.	14x17 clear film supply empty.	<None>/interact
EMPTY 14x17 PAPER	The 14x17 inch paper supply magazine is empty.	14x17 paper supply empty.	<None>/interact
EMPTY 24x24	The 24x24 inch film supply magazine is empty.	24x24 film supply empty.	<None>/interact
EMPTY 24x24 BLUE	The 24x24 inch blue film supply magazine is empty.	24x24 blue film supply empty.	<None>/interact
EMPTY 24x24 CLR	The 24x24 inch clear film supply magazine is empty.	24x24 clear film supply empty.	<None>/interact
EMPTY 24x24 PAPER	The 24x24 inch paper supply magazine is empty.	24x24 paper supply empty.	<None>/interact
EMPTY 24x30	The 24x30 inch film supply magazine is empty.	24x30 film supply empty.	<None>/interact
EMPTY 24x30 BLUE	The 24x30 inch blue film supply magazine is empty.	24x30 blue film supply empty.	<None>/interact
EMPTY 24x30 CLR	The 24x30 inch clear film supply magazine is empty.	24x30 clear film supply empty.	<None>/interact
EMPTY 24x30 PAPER	The 24x30 inch paper supply magazine is empty.	24x30 paper supply empty.	<None>/interact
EMPTY A4 PAPER	The A4 paper supply magazine is empty.	A4 paper supply empty.	<None>/interact
EMPTY A4 TRANS	The A4 transparency supply magazine is empty.	A4 transparency supply empty.	<None>/interact
EXPOSURE FAILURE	The exposure device has failed due to some unspecified reason.	Exposure device has failed.	<None>/interact
FILM JAM	A film transport error has occurred and a film is jammed in the printer or processor.	Film jam.	<None>/interact

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
FILM TRANSP ERR	There is a malfunction with the film transport, there may or may not be a film jam.	Film transport problem.	<None>/interact
FINISHER EMPTY	The finisher is empty.	Finisher is empty.	<None>/interact
FINISHER ERROR	The finisher is not operating due to some unspecified reason.	Finisher problem.	<None>/interact
FINISHER LOW	The finisher is low on supplies	Finisher low.	<None>/interact
LOW 8x10	The 8x10 inch film supply magazine is low.	8x10 film supply low.	<None>/interact
LOW 8x10 BLUE	The 8x10 inch blue film supply magazine is low.	8x10 blue film supply low.	<None>/interact
LOW 8x10 CLR	The 8x10 inch clear film supply magazine is low.	8x10 clear film supply low.	<None>/interact
LOW 8x10 PAPR	The 8x10 inch paper supply magazine is low.	8x10 paper supply low.	<None>/interact
LOW 10x12	The 10x12 inch film supply magazine is low.	10x12 film supply low.	<None>/interact
LOW 10x12 BLUE	The 10x12 inch blue film supply magazine is low.	10x12 blue film supply low.	<None>/interact
LOW 10x12 CLR	The 10x12 inch clear film supply magazine is low.	10x12 clear film supply low.	<None>/interact
LOW 10x12 PAPR	The 10x12 inch paper supply magazine is low.	10x12 paper supply low.	<None>/interact
LOW 10x14	The 10x14 inch film supply magazine is low.	10x14 film supply low.	<None>/interact
LOW 10x14 BLUE	The 10x14 inch blue film supply magazine is low.	10x14 blue film supply low.	<None>/interact
LOW 10x14 CLR	The 10x14 inch clear film supply magazine is low.	10x14 clear film supply low.	<None>/interact
LOW 10x14 PAPR	The 10x14 inch paper supply magazine is low.	10x14 paper supply low.	<None>/interact
LOW 11x14	The 11x14 inch film supply magazine is low.	11x14 film supply low.	<None>/interact
LOW 11x14 BLUE	The 11x14 inch blue film supply magazine is low.	11x14 blue film supply low.	<None>/interact
LOW 11x14 CLR	The 11x14 inch clear film supply magazine is low.	11x14 clear film supply low.	<None>/interact

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/ camera symbol
LOW 11x14 PAPR	The 11x14 inch paper supply magazine is low.	11x14 paper supply low.	<None>/interact
LOW 14x14	The 14x14 inch film supply magazine is low.	14x14 film supply low.	<None>/interact
LOW 14x14 BLUE	The 14x14 inch blue film supply magazine is low.	14x14 blue film supply low.	<None>/interact
LOW 14x14 CLR	The 14x14 inch clear film supply magazine is low.	14x14 clear film supply low.	<None>/interact
LOW 14x14 PAPR	The 14x14 inch paper supply magazine is low.	14x14 paper supply low.	<None>/interact
LOW 14x17	The 14x17 inch film supply magazine is low.	14x17 film supply low.	<None>/interact
LOW 14x17 BLUE	The 14x17 inch blue film supply magazine is low.	14x17 blue film supply low.	<None>/interact
LOW 14x17 CLR	The 14x17 inch clear film supply magazine is low.	14x17 clear film supply low.	<None>/interact
LOW 14x17 PAPR	The 14x17 inch paper supply magazine is low.	14x17 paper supply low.	<None>/interact
LOW 24x24	The 24x24 inch film supply magazine is low.	24x24 film supply low.	<None>/interact
LOW 24x24 BLUE	The 24x24 inch blue film supply magazine is low.	24x24 blue film supply low.	<None>/interact
LOW 24x24 CLR	The 24x24 inch clear film supply magazine is low.	24x24 clear film supply low.	<None>/interact
LOW 24x24 PAPR	The 24x24 inch paper supply magazine is low.	24x24 paper supply low.	<None>/interact
LOW 24x30	The 24x30 inch film supply magazine is low.	24x30 film supply low.	<None>/interact
LOW 24x30 BLUE	The 24x30 inch blue film supply magazine is low.	24x30 blue film supply low.	<None>/interact
LOW 24x30 CLR	The 24x30 inch clear film supply magazine is low.	24x30 clear film supply low.	<None>/interact
LOW 24x30 PAPR	The 24x30 inch paper supply magazine is low.	24x30 paper supply low.	<None>/interact
LOW A4 PAPR	The A4 paper supply magazine is low.	A4 paper supply low.	<None>/interact
LOW A4 TRANS	The A4 transparency supply magazine is low.	A4 transparency supply low.	<None>/interact
NO RECEIVE MGZ	The film receive magazine no available.	Film receiver not available.	<None>/interact

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
NO RIBBON	The ribbon cartridge needs to be replaced.	Replace ribbon cartridge.	<None>/interact
NO SUPPLY MGZ	The film supply magazine specified for this job is not available.	Film supply not available.	<None>/interact
CHECK PRINTER	The printer is not ready at this time, operator intervention is required to make the printer available.	Check camera.	<None>/interact
CHECK PROC	The processor is not ready at this time, operator intervention is required to make the printer available.	Check processor.	<None>/interact
PRINTER DOWN	The printer is not operating due to some unspecified reason.	Camera down.	<None>/interact
PRINTER INIT	The printer is not ready at this time, it is expected to become available without intervention, For example, it may be in a normal warm-up state.	Camera initializing.	<None>/idle
PRINTER OFFLINE	The printer has been disabled by an operator or service person.	Camera off-line.	<None>/interact
PROC DOWN	The processor is not operating due to some unspecified reason.	Processor down.	<None>/interact
PROC INIT	The processor is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state.	Processor initializing.	<None>/idle
PROC OVER-FLOW FL	Processor chemicals are approaching the overflow full mark.	Processor chemicals overflow.	<None>/interact
PROC OVER-FLOW HI	Processor chemicals have reached the overflow full mark.	Processor chemicals near overflow.	<None>/interact
QUEUED	Print job in Queue	-	<None>/idle
RECEIVER FULL	The Film receive magazine is full.	Receiver full.	<None>/interact
REQ MED NOT INST	The requested film, paper, or other media supply magazine is installed in the printer, but may be available with operator intervention.	Install media supply.	<None>/interact
REQ MED NOT AVAI	The requested film, paper, or other media requested is not available on this printer.	Media supply not available on this camera. Queue stopped. Change camera.	Queue for this camera will be STOPPED/queue stopped

Table 20: Printer Status Infos within Printer SOP Class/Execution Status Infos within Print Job SOP Class

Printer Status info/Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
RIBBON ERROR	There is an unspecified problem with the print ribbon.	Error with print ribbon.	<None>/interact
SUPPLY EMPTY	The printer is out of film.	Camera out of film.	<None>/interact
SUPPLY LOW	The film supply is low.	Film supply low.	<None>/interact
UNKNOWN	There is an unspecified problem.	Unspecified problem with camera.	<None>/interact

A.9.2 Additional Status Information - AGFA Printers

Table 21: Printer Status Infos: Additional Agfa printer status infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
WARMING UP	Printer is in the warm-up stage. Spooling of print jobs to disk is still possible.	Camera is warming up.	<None>idle
OFFLINE	Printer is switched off-line. Spooling of print jobs to disk is still possible.	Camera is switched off-line.	<None>/interact
NONE	General printer warning, no specific information is available. Spooling of print jobs to disk is still possible.	-	<None>/idle
-	-	-	-

A.9.3 Additional Status Information - Kodak PACS Link (formerly Imation)

Table 22: Printer Status Infos: Additional Kodak infos for Pacs Link (formerly Imation cameras)

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
SUPPLY MGZ ERR	The supply magazine has an error.	Film supply has an error.	<None>/interact
-	-	-	-

A.9.4 Additional Status Information - Kodak 190I

Table 23: Printer Status Infos: Additional Kodak infos for Kodak 190

Printer Status info	Description	Message string visible in the Status Bar	Other action for <i>syngo</i> / camera symbol
PRINTER STOPPED	The printer has stopped	Camera has stopped.	<None>/interact
FATAL ERROR	Fatal error.	Fatal error. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
-	-	-	-

A.9.5 Additional Status Information - Kodak 2180/1120

Table 24: Printer Status Infos: Additional Kodak infos for 2180/1120

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
PRINTER NOT RDY	Printer not ready.	Camera not ready.	<None>/interact
CHECK PROCESSOR	Check processor.	Check processor.	<None>/interact
NO TONER	No toner.	No toner.	<None>/interact
FATAL	Fatal error.	Fatal error. Queue stopped.	Queue for this camera will be STOPPED/queue stopped
-	-	-	-

A.9.6 Additional Status Information - Codonics

Table 25: Printer Status Infos: Additional Codonics infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
STANDARD	Printer is ready.	Camera is ready.	<None>/Normal
LOAD A-SIZE	Load A-size media.	Load A-size media.	<None>/interact
LOAD A-DVPA-PER	Load A-size black and white paper.	Load A-size black and white paper.	<None>/interact
LOAD A-CVPA-PER	Load A-size color paper.	Load A-size color paper.	<None>/interact
LOAD A-CVTRANS	Load A-size transparencies.	Load A-size transparencies.	<None>/interact
LOAD A4-SIZE	Load A4-size media.	Load A4-size media.	<None>/interact
LOAD A4-DVPA-PER	Load A4-size black and white paper.	Load A4-size black and white paper.	<None>/interact
LOAD A4-CVPA-PER	Load A4-size color paper.	Load A4-size color paper.	<None>/interact
LOAD A4-CVTRANS	Load A4-size transparencies.	Load A4-size transparencies.	<None>/interact
LOAD LA-SIZE	Load LA-size media.	Load LA-size media.	<None>/interact
LOAD LA-DVPA-PER	Load LA-size black and white paper.	Load LA-size black and white paper.	<None>/interact
LOAD LA-CVPA-PER	Load LA-size color paper.	Load LA-size color paper.	<None>/interact
LOAD LA-CVTRANS	Load LA-size transparencies.	Load LA-size transparencies.	<None>/interact

Table 25: Printer Status Infos: Additional Codonics infos

Printer Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
LOAD LA4-SIZE	Load LA4-size media.	Load LA4-size media.	<None>/interact
LOAD LA4-DVPAPER	Load LA4-size black and white paper.	Load LA4-size black and white paper.	<None>/interact
LOAD LA4-CVPAPER	Load LA4-size color paper.	Load LA4-size color paper.	<None>/interact
LOAD LA4-CVTRANS	Load LA4-size transparencies.	Load LA4-size transparencies.	<None>/interact
LOAD XLA-SIZE	Load XLA-size media.	Load XLA-size media.	<None>/interact
LOAD XLA-DVPAPER	Load XLA-size black and white paper.	Load XLA-size black and white paper.	<None>/interact
LOAD XLA-CVPAPER	Load XLA-size color paper.	Load XLA-size color paper.	<None>/interact
LOAD XLA-CVTRANS	Load XLA-size transparencies.	Load XLA-size transparencies.	<None>/interact
LOAD XLA4-SIZE	Load XLA4-size media.	Load XLA4-size media.	<None>/interact
LOAD XLA4-DVPAPE	Load XLA4-size black and white paper.	Load XLA4-size black and white paper.	<None>/interact
LOAD XLA4-CVPAPE	Load XLA4-size color paper.	Load XLA4-size color paper.	<None>/interact
LOAD XLA4-CVTRAN	Load XLA4-size transparencies.	Load XLA4-size transparencies.	<None>/interact
LOAD XLW-SIZE	Load XLW-size media.	Load XLW-size media.	<None>/interact
LOAD XLW-DVPAPER	Load XLW-size black and white paper.	Load XLW-size black and white paper.	<None>/interact
LOAD XLW-CVPAPER	Load XLW-size color paper.	Load XLW-size color paper.	<None>/interact
LOAD 8X10-SIZE	Load 8x10 media.	Load 8x10 media.	<None>/interact
LOAD 8X10-DVFILM	Load 8x10 black and white film.	Load 8x10 black and white film.	<None>/interact
SUPPLY MISSING	The film supply magazine specified for this job is not available.	Film supply not available.	<None>/interact
RIBBON MISSING	Ribbon is missing.	Ribbon is missing.	<None>/interact
RIBBON EMPTY	Ribbon is empty.	Ribbon is empty.	<None>/interact
TOP COVER OPEN	Top cover of printer is open.	Top cover of camera is open.	<None>/interact
-	-	-	-

A.9.7 Additional DICOM Execution Status Information

Table 26: Additional DICOM Execution Status Infos

Execution Status info	Description	Message string visible in the Status Bar	Other action for syngo/camera symbol
INVALID PAGE DES	The specified page layout cannot be printed or other page description errors have been detected.	Film Job cannot be printed on this camera. Queue stopped. Please redirect film job.	Queue for this camera will be STOPPED/queue stopped
INSUFFIC MEMORY	There is not enough memory available to complete this job.	Not enough memory available in camera. Queue stopped. Please continue queue or change camera.	Queue for this camera will be STOPPED/queue stopped
NONE	General printer warning, no specific information is available. Spooling of print jobs to disk is still possible.	-	<None>/idle

Copyright © Siemens AG 2004. All rights reserved.

Alle Rechte vorbehalten.

Copyright © Siemens AG 2004. All rights reserved.
Alle Rechte vorbehalten.