



**SIEMENS**



Conventional



TrueForm

[www.siemens.com/mri](http://www.siemens.com/mri)

# TrueForm Magnet and Gradient Design

Imaging the true form of the anatomy.

Answers for life.

# TrueForm Magnet and Gradient design

TrueForm® design, pioneered and introduced by Siemens in 2007, is the result of numerous technological innovations and addresses the need for field homogeneity. One important aspect of TrueForm design is TrueForm Magnet and Gradient design, which enables imaging of large volumes with Open Bore and short magnet design – without compromising spatial resolution or speed.

## TrueForm Key Benefits

### Increased homogeneity

- Large volume coverage
- Cylindrically optimized homogeneity volume corresponding better to the true form of the human body

### Higher image quality

- Reduction of the unusable edges on images
- Better fat saturation coverage

### Optimized Field of View


- Ability to use a FoV of up to 50 cm x 50 cm x 45 cm
- Supporting a full range of applications

Image without distortions along the edges due to TrueForm Magnet and Gradient design (MAGNETOM Aera, T2 TSE, 45 cm x 45 cm FoV).



# Excellent clinical performance

TrueForm Magnet and Gradient design has established itself in the clinical routine with over 1,000 installations worldwide as of September 2011. The images below show the elimination of commonly known distortions along the edges.



Images showing the effect of TrueForm Magnet and Gradient design without conventional distortions along the edges even in more challenging measurements using fat saturation techniques.

A: MAGNETOM Skyra,  
T1 3D VIBE FatSat,  
40 cm x 40 cm FoV  
*University Medical Center  
Mannheim, Germany*

B: MAGNETOM Skyra,  
T1 3D VIBE FatSat MIP,  
30 cm x 38 cm FoV  
*University Medical Center  
Mannheim, Germany*




What our customers say:

TrueForm Design leads to increased homogeneity and image quality for all organs.”

*Markus Lentschig,  
MR and PET/CT Imaging Center  
Bremen Mitte, Germany*

# Imaging the true form of the anatomy



## Conventional

Comparison of conventional (left) and TrueForm (right) Magnet design showing the increased possible imaging volume due to the optimized homogeneity.



## TrueForm

TrueForm Magnet and Gradient design is optimized for a cylindrically shaped volume rather than a spherical one, as shown above.

TrueForm Magnet design is an innovation that produces a cylindrically optimized homogeneity volume instead of the conventional spherical or elliptical volume. A cylinder corresponds better to the true form of the human body. TrueForm Gradient design also creates a cylindrical shape for the gradient linearity volume.

The two combined result in better image quality by reducing the unusable edges in the images as well as better fat saturation covered in a scan.

TrueForm design reduces the overlap needed between steps for extended Field of View (FoV) exams, reducing the number of steps needed for a given scanning range.

TrueForm Magnet and Gradient design is optimized for a cylindrically shaped volume rather than a spherical one. An ideal cylinder has 1.5 times more volume than a sphere with the same axis lengths, which means a greater homogeneity volume compared to conventional magnets with identical "nominal" specifications.


For MAGNETOM Aera, MAGNETOM Skyra and MAGNETOM Verio, TrueForm Magnet and Gradient design is standard providing the ability to use a FoV up to 50 cm x 50 cm x 45 cm depending on the application.

---

*TrueForm Magnet and Gradient design is standard on MAGNETOM Skyra 3T ...*



*... MAGNETOM Verio 3T ...*



*... and MAGNETOM Aera 1.5T systems.*



On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States. All devices listed herein may not be licensed according to Canadian Medical Devices Regulations. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced

Please find fitting accessories:  
[www.siemens.com/medical-accessories](http://www.siemens.com/medical-accessories)

#### **Global Business Unit**

Siemens AG  
Healthcare Sector  
Magnetic Resonance  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0

#### **Local Contact Information**

##### **In the USA**

Siemens Medical Solutions USA, Inc.  
51 Valley Stream Parkway  
Malvern, PA 19355  
Phone: +1 888-826-9702  
Phone: +1 610-448-4500  
Fax: +1 610-448-2254

##### **In China**

Siemens Medical Park, Shanghai  
278, Zhouzhu Road  
SIMZ, Nanhui District  
Shanghai, 201318, P.R. China  
Phone: +86-21-38895000  
Fax: +86-10-28895001

##### **In Japan**

Siemens-Asahi  
Medical Technologies Ltd.  
Takanawa Park Tower 14F  
20-14, Higashi-Gotanda 3-chome  
Shinagawa-ku, Tokyo 141-8644  
Phone: +81 3 5423 8411

##### **In Asia**

Siemens Pte Ltd  
Healthcare Sector  
Regional Headquarters  
The Siemens Center  
60 MacPherson Road, Singapore  
348615  
Phone: +65 6490-6000  
Fax: +65 6490-6001

#### **Global Siemens Headquarters**

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

#### **Global Siemens Healthcare Headquarters**

Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

#### **Legal Manufacturer**

Siemens AG  
Wittelsbacherplatz 2  
DE-80333 Muenchen  
Germany