# Case Report: Cortical Dysplasia

Michael Kean: Micheal Ditchfield, M.D.

Children's MRI Centre, Royal Children's Hospital, Parkville, Victoria, Australia

### **Patient history**

12-year-old girl presented after inconclusive scan at an outside institution. The patient was scheduled for a standard Epilepsy scan on our 1.5T MAGNETOM Avanto using a standard Head Matrix coil. The patient was then referred to the 3T MAGNETOM Trio and examined using the 32-channel Head Matrix coil.

### Sequence details

The images obtained are from a 1.5T MAGNETOM Avanto with software version syngo MR 13 (MPRAGE 0.7 mm isotropic voxels) and from a 3T MAGNETOM Trio with software version syngo MR B15 (MPRAGE 0.7 mm isotropic, 0.9 mm isotropic voxels using Water Excitation (WE), DarkFluid 3D SPACE and lipid suppressed 135 ms TE single voxel MR spectroscopy – MRS)

## Image findings


The patient presented with migraines and 1.5T MR demonstrated a left frontal lesion possibly a cortical dysplasia or developmental tumor (DNET). 3T MR demonstrates a focally thickened cortical ribbon with blurring of the grey/ white matter interface. The MR appearances are that of a cortical malformation.

#### Discussion

The images obtained using the 1.5T MAGNETOM Avanto clearly show the lesion but the increased signal-to-noise ratio (SNR) of the 32-channel Head Matrix coil and 3T (MPRAGE and 3D SPACE) provide the radiologist with a level of anatomical detail that enables a diagnosis with a higher degree of sensitivity and specificity.



1 MPRAGE image with 0.7 mm isotropic voxel, acquired using a 1.5T MAGNETOM Avanto system, demonstrating the left frontal lesion.



2 Corresponding MPRAGE image acquired using a 3T MAGNETOM Trio and the 32-channel Head Matrix coil. The focally thickened cortical ribbon and the blurring of the grey/white matter interface are clearly visible.

