

The purpose of this document is to provide an introductory walkthrough of one aspect of the PETLINK™ Guideline. Here we present examples of various and representative bit fields as they might be used – fields which are currently active in support of the Siemens mCT 32-bit bin-address packet format. The goal is to help the first-time user gain a better, general understanding of the packet format in use.

A companion example 32-bit list-mode file is provided: **ea_32ba_walkthrough.I32**

Two companion example C-code files are provided:

ml_ea_32ba_walkthrough_1.c

[Generates the example *.I32 file.]

lmsw32ba_verbose_1.c

[Reports contents of *.I32 file verbosely.]

See also the PETLINK Guideline file.

Overview of list-mode file content:

The list-mode file generated shows examples of various packets with fields being incremented or decremented. Three types of tag packets are used – i.e. elapsed time, horizontal bed position, and lost event tally. [The elapsed time and lost event tally packets are incrementing except for the last maximum value packet examples. The horiz. bed position packets are decrementing except for the last, most-negative value packet.] Two types of event packet fields are shown. The 30-bit BA field is incremented up from zero for a few packets and then incremented up to the maximum value for a few packets. The single-bit Prompt field is shown first with Prompt = 0 (Delayed) examples followed by Prompt = 1 (Prompt) examples as the BA field increments.

Appendix 1: Here is the content of the list-mode generating C-code file,

ml_ea_32ba_walkthrough_1.c:

```
//file: ml_ea_32ba_walkthrough_1.c
// make list-mode file - software list mode generator
// Eagle 32-bit bin-address packet field walkthrough
// 9-Oct-2012 wfj

#include <stdio.h>
#include <ctype.h>
#include <math.h>
#include <stdlib.h>
#include <iostream.h>

FILE *streamt;

main (int argc, char **argv)
{

    static int i,j;
    static int e1, e2;
    static int b1, b2;
    static int let1, let2;
    static unsigned __int64 il;

    static char *out_file;

    static int qb;
    static int prompt; // 0 to 1
    static unsigned int ewl;
    static unsigned int ms = 0; // milliseconds
    static int hbp = 0; // horizontal bed position
    static unsigned int lost = 0; // lost event tally quantity
    static unsigned int ba;
```

```

// get a filename from the command line
if (argc < 2)
{
    printf ("usage: %s out_file_name\n",argv[0]);
    exit(1);
}
else
    out_file = argv[1];

    streamt = fopen (out_file, "wb");
    if ( streamt == NULL) {
        printf ("No file opened %s\n",out_file);
        exit (1);
    }

// Load Output File

// Incrementing Elapsed Time Tag Packet
    el = 0x80000000 | (ms++ & 0xffffffff);
    j = fwrite (&el, sizeof(qb), 1, streamt);

// Decrementing Horizontal Bed Position Tag Packet
    bl = 0xc4000000 | (hbp-- & 0xffff);
    j = fwrite (&bl, sizeof(qb), 1, streamt);

// Incrementing Lost Event Tally Tag Packet - Type 7 (GIM) & NonFunctional as a Loss
Tally
    letl = 0xbc000000 | (lost++ & 0xffff);
    j = fwrite (&letl, sizeof(qb), 1, streamt);

// Packets to Increment Up from Zero a Portion of the ba Field with Prompt bit cleared
(Delayed)
    ba = 0; prompt = 0;
    for (ba=0; ba <= 15; ba++) {
        ewl = 0x00000000 | ((prompt & 1) << 30) | (ba & 0x3fffffff);
        j = fwrite (&ewl, sizeof(qb), 1, streamt);
        il = il + 1;
    }

// Incrementing Elapsed Time Tag Packet
    el = 0x80000000 | (ms++ & 0xffffffff);
    j = fwrite (&el, sizeof(qb), 1, streamt);

// Decrementing Horizontal Bed Position Tag Packet
    bl = 0xc4000000 | (hbp-- & 0xffff);
    j = fwrite (&bl, sizeof(qb), 1, streamt);

// Incrementing Lost Event Tally Tag Packet - Type 7 (GIM) & NonFunctional as a Loss
Tally
    letl = 0xbc000000 | (lost++ & 0xffff);
    j = fwrite (&letl, sizeof(qb), 1, streamt);

// Packets to Increment Up to Maximum a Portion of the ba Field with Prompt bit set
    ba = 0; prompt = 1;
    for (ba=0x3fffffff; ba <= 0x3fffffff; ba++) {
        ewl = 0x00000000 | ((prompt & 1) << 30) | (ba & 0x3fffffff);
        j = fwrite (&ewl, sizeof(qb), 1, streamt);
        il = il + 1;
    }

// Maximum Value for Elapsed Time Tag Packet
    ms = 0xffffffff; // all bits in 29-bit field set to one
    el = 0x80000000 | (ms & 0xffffffff);
    j = fwrite (&el, sizeof(qb), 1, streamt);

// Minimum Value (most negative) Horizontal Bed Position Tag Packet
    hbp = 0x80000;
    bl = 0xc4000000 | (hbp & 0xffff);
    j = fwrite (&bl, sizeof(qb), 1, streamt);

```

```

// Maximum Value Lost Event Tally Tag Packet - Type 7 (GIM) & NonFunctional as a Loss
Tally
    lost = 0xffffffff; // all bits in 20-bit field set to one
    let1 = 0xbc000000 | (lost & 0xffffffff);
    j = fwrite (&let1, sizeof(qb), 1, streamt);

    printf (" number of 64-bit packets output: %I64d file size: %I64x\n",i1,8*i1);
    fclose (streamt);
    exit(0);
}

```

Appendix 2: Here is the content of the list-mode verbose reporting C-code file, lmsw32ba_verbose_1.c:

```
// file: lmsw32ba_verbose_1.c
// Verbose Listing for 32-Bit Bin-Address Packets

// 8-Oct-2012 wfj

#include <stdio.h>
#include <ctype.h>
#include <math.h>
#include <stdlib.h>
#include <iostream.h>

FILE *streami;

main (int argc, char **argv)
{
    static __int64 i,j,k,i1,j1,j2;
    static __int64 i64_i1;

    static char *in_file;

    static int qb;
    static unsigned int ba;
    static int prompt;
    static unsigned int ewl;
    static int pkt;
    static int tag;

    static int ms;           // milliseconds
    static int hbp;          // horizontal bed position value
    static int hbp_se;       // sign-extended horiz. bed position value
    static int lost;         // lost event tally value
    static int tof_se;       // sign-extended tof value

    // get a filename from the command line

    if (argc < 1) {
        printf ("usage: %s in_file_name\n",argv[0]);
        exit(1);
    }
    else { in_file = argv[1];}

    streami = fopen (in_file, "rb");
    if ( streami == NULL) {
        printf ("No file opened %s\n",in_file);
        exit (1);
    }

    k = 0; i64_i1 = 0;
```

```

while ((i = fread (&ewl, sizeof(qb), 1, streami) ) != 0) {

    pkt = pkt + 1;

    il = il + 1;

    tag = ((ewl>>31)&1)!=0;

    prompt = 0;

    if(!tag) {
        prompt = (ewl>>30) & 1;
        ba = ewl & 0x3fffffff;
        printf(" EVENT: pkt_cnt: %I64d ewl(h): %8x prompt: %1d ba(h): %8x
ba(d): %10d\n",
                                                    il, ewl,
prompt, ba, ba);
    }

    if ( tag){
        if ( ((ewl >> 28) & 0xe) == 8 ) { // Elapsed Time Tag Packet
            ms = ewl & 0x1fffffff; // Extract 29-bit millisecond field
            printf(" TAG32_ElapsedTime:    pkt_cnt: %I64d ewl(h): %8x
ms(h): %7x ms(d): %9d\n",
                                                    il,
ewl, ms, ms);
        }
        if ( ((ewl >> 24) & 0xff) == 0xc4 ) { // Horizontal Bed Position
Tag Packet
            hbp = ewl & 0xffff; // Extract 20-bit bed position field
            hbp_se = hbp; // Assume hbp is Zero
            if (((hbp >> 19) & 1) == 1) hbp_se = hbp | 0xffff0000; //
Need to Sign Extend
            printf(" TAG32_HorizBedPos:    pkt_cnt: %I64d ewl(h): %8x
hbp(h): %6x hbp_se(d): %8d\n",
                                                    il,
ewl, hbp, hbp_se);
        }
        if ( ((ewl >> 24) & 0xfc) == 0xbc ) { // Lost Event Tally Tag
Packet - Type 7 (GIM)
            lost = ewl & 0xffff; // Extract 20-bit lost tally field
            printf(" TAG32_LostEventTally: pkt_cnt: %I64d ewl(h): %8x
lost(d): %7d\n",
                                                    il,
ewl, lost);
        }
    }

    fclose (streami);

    exit(0);
}

```

Appendix 3: Here is a text report generated by that “verbose” C-code for the example *.l32 file:

```

TAG32_ElapsedTime:      pkt_cnt: 1 ewl(h): 80000000 ms(h):          0 ms(d):          0
TAG32_HorizBedPos:      pkt_cnt: 2 ewl(h): c4000000 hbp(h):          0 hbp_se(d):          0
TAG32_LostEventTally:   pkt_cnt: 3 ewl(h): bc000000 lost(d):          0
EVENT: pkt_cnt: 4 ewl(h):          0 prompt: 0 ba(h):          0 ba(d):          0
EVENT: pkt_cnt: 5 ewl(h):          1 prompt: 0 ba(h):          1 ba(d):          1
EVENT: pkt_cnt: 6 ewl(h):          2 prompt: 0 ba(h):          2 ba(d):          2
EVENT: pkt_cnt: 7 ewl(h):          3 prompt: 0 ba(h):          3 ba(d):          3
EVENT: pkt_cnt: 8 ewl(h):          4 prompt: 0 ba(h):          4 ba(d):          4
EVENT: pkt_cnt: 9 ewl(h):          5 prompt: 0 ba(h):          5 ba(d):          5
EVENT: pkt_cnt: 10 ewl(h):         6 prompt: 0 ba(h):          6 ba(d):          6
EVENT: pkt_cnt: 11 ewl(h):         7 prompt: 0 ba(h):          7 ba(d):          7
EVENT: pkt_cnt: 12 ewl(h):         8 prompt: 0 ba(h):          8 ba(d):          8
EVENT: pkt_cnt: 13 ewl(h):         9 prompt: 0 ba(h):          9 ba(d):          9
EVENT: pkt_cnt: 14 ewl(h):        a prompt: 0 ba(h):         a ba(d):         10
EVENT: pkt_cnt: 15 ewl(h):        b prompt: 0 ba(h):         b ba(d):         11
EVENT: pkt_cnt: 16 ewl(h):        c prompt: 0 ba(h):         c ba(d):         12
EVENT: pkt_cnt: 17 ewl(h):        d prompt: 0 ba(h):         d ba(d):         13
EVENT: pkt_cnt: 18 ewl(h):        e prompt: 0 ba(h):         e ba(d):         14
EVENT: pkt_cnt: 19 ewl(h):        f prompt: 0 ba(h):         f ba(d):         15
TAG32_ElapsedTime:      pkt_cnt: 20 ewl(h): 80000001 ms(h):          1 ms(d):          1
TAG32_HorizBedPos:      pkt_cnt: 21 ewl(h): c40fffff hbp(h):        ffffff hbp_se(d):        -1
TAG32_LostEventTally:   pkt_cnt: 22 ewl(h): bc000001 lost(d):          1
EVENT: pkt_cnt: 23 ewl(h): 7fffffff0 prompt: 1 ba(h): 3fffffff0 ba(d): 1073741808
EVENT: pkt_cnt: 24 ewl(h): 7fffffff1 prompt: 1 ba(h): 3fffffff1 ba(d): 1073741809
EVENT: pkt_cnt: 25 ewl(h): 7fffffff2 prompt: 1 ba(h): 3fffffff2 ba(d): 1073741810
EVENT: pkt_cnt: 26 ewl(h): 7fffffff3 prompt: 1 ba(h): 3fffffff3 ba(d): 1073741811
EVENT: pkt_cnt: 27 ewl(h): 7fffffff4 prompt: 1 ba(h): 3fffffff4 ba(d): 1073741812
EVENT: pkt_cnt: 28 ewl(h): 7fffffff5 prompt: 1 ba(h): 3fffffff5 ba(d): 1073741813
EVENT: pkt_cnt: 29 ewl(h): 7fffffff6 prompt: 1 ba(h): 3fffffff6 ba(d): 1073741814
EVENT: pkt_cnt: 30 ewl(h): 7fffffff7 prompt: 1 ba(h): 3fffffff7 ba(d): 1073741815
EVENT: pkt_cnt: 31 ewl(h): 7fffffff8 prompt: 1 ba(h): 3fffffff8 ba(d): 1073741816
EVENT: pkt_cnt: 32 ewl(h): 7fffffff9 prompt: 1 ba(h): 3fffffff9 ba(d): 1073741817
EVENT: pkt_cnt: 33 ewl(h): 7fffffffA prompt: 1 ba(h): 3fffffffA ba(d): 1073741818
EVENT: pkt_cnt: 34 ewl(h): 7fffffffB prompt: 1 ba(h): 3fffffffB ba(d): 1073741819
EVENT: pkt_cnt: 35 ewl(h): 7fffffffC prompt: 1 ba(h): 3fffffffC ba(d): 1073741820
EVENT: pkt_cnt: 36 ewl(h): 7fffffffD prompt: 1 ba(h): 3fffffffD ba(d): 1073741821
EVENT: pkt_cnt: 37 ewl(h): 7fffffffE prompt: 1 ba(h): 3fffffffE ba(d): 1073741822
EVENT: pkt_cnt: 38 ewl(h): 7fffffffF prompt: 1 ba(h): 3fffffffF ba(d): 1073741823
TAG32_ElapsedTime:      pkt_cnt: 39 ewl(h): 9fffffff ms(h): 1fffffff ms(d): 536870911
TAG32_HorizBedPos:      pkt_cnt: 40 ewl(h): c4080000 hbp(h):      80000 hbp_se(d): -524288
TAG32_LostEventTally:   pkt_cnt: 41 ewl(h): bc0fffff lost(d): 1048575
```

Appendix 4: Here is an example of what the V file viewer utility can show of the *.l32 file:

Getting a copy of Fileviewer.exe: <http://www.fileviewer.com/>

Setting up Fileviewer (Version 12) for 32-bit packet displayed per line:

View>>Hex Mode (selected)
View>>Hex Formats>> Dword (selected)
View>>Hex Formats>> Flip Ends (selected)
View>>Hex Formats>> Set Hex Line Length (4)

What the V utility shows for this *.l32 file:

	03020100
00000000	80000000
00000004	C4000000
00000008	BC000000
0000000C	00000000
00000010	00000001
00000014	00000002
00000018	00000003
0000001C	00000004
00000020	00000005
00000024	00000006
00000028	00000007
0000002C	00000008
00000030	00000009
00000034	0000000A
00000038	0000000B
0000003C	0000000C
00000040	0000000D
00000044	0000000E
00000048	0000000F
0000004C	80000001
00000050	C40FFFFFFF
00000054	BC000001
00000058	7FFFFFF0
0000005C	7FFFFFF1
00000060	7FFFFFF2
00000064	7FFFFFF3
00000068	7FFFFFF4
0000006C	7FFFFFF5
00000070	7FFFFFF6
00000074	7FFFFFF7
00000078	7FFFFFF8
0000007C	7FFFFFF9
00000080	7FFFFFFA
00000084	7FFFFFFB
00000088	7FFFFFFC
0000008C	7FFFFFFD
00000090	7FFFFFFE
00000094	7FFFFFFF
00000098	9FFFFFFF
0000009C	C4080000
000000A0	BC0FFFFFFF