

SIEMENS

Eckert & Ziegler ^{11}C Modular-Lab Chemistry Solution

www.siemens.com/mi

Answers for life.

Introduction

The Eckert and Ziegler (E&Z) ^{11}C Modular-Lab Chemistry Solution comprises Methyl Iodide/Methyl Triflate (MeI/MeOTf) system, Methylation system and HPLC system for the synthesis and purification of ^{11}C -labeled radiopharmaceuticals. Consumable products include PharmTracer cassettes for use with the Methylation system.

The table below shows a list of ^{11}C -labeled intermediates and compounds that can be synthesized using the Modular-Lab Chemistry Solution. The systems that are required for the synthesis of each compound are indicated.


The entire E&Z ^{11}C Modular-Lab Chemistry Solution can either be installed in a single hot cell or in two mini-cells with the Modular-Lab MeI/MeOTf in one mini-cell and the Modular-Lab Methylation system and High Performance Liquid Chromatography (HPLC) system in a second mini-cell.

Biomarker*	MeI/MeOTf system	Methylation cassette system	HPLC system
^{11}C -Methyl Iodide	✓		
^{11}C -Methyl Triflate	✓		
^{11}C -Methionine	✓	✓	
^{11}C -Choline	✓	✓	
^{11}C -Acetate	✓	✓	
^{11}C -Raclopride	✓	✓	✓

* These agents are not approved by the US Food and Drug Administration for routine clinical use. Check with other regulatory authorities for marketing status outside the US.


Background

The PET radioisotope ^{11}C is characterized by its relatively short half life (~20 minutes), which provides the advantage of being used for the synthesis of many short-lived imaging biomarkers, but also has a long enough half-life to support the synthesis of radiopharmaceuticals that require multiple steps. The process flow for the synthesis of various ^{11}C -labeled intermediates and compounds are depicted graphically in the schematic. The first step in the process involves the production of $^{11}\text{C}\text{-CO}_2$ in a PET cyclotron using $^{14}\text{N}\text{-N}_2$ gas as the starting material and the $^{14}\text{N}(\text{p}, \alpha)^{11}\text{C}$ reaction.



In the next step, $^{11}\text{C}\text{-Methyl Iodide}$ or $^{11}\text{C}\text{-Methyl Triflate}$ is synthesized from $^{11}\text{C}\text{-CO}_2$ produced by the cyclotron using the E&Z Modular-Lab Mel/MeOTf system. $^{11}\text{C}\text{-Methyl Iodide}$ and $^{11}\text{C}\text{-Methyl Triflate}$ are the most important and frequently used secondary ^{11}C -labeling precursor used for the synthesis of a variety of compounds. Their versatility allows the synthesis of a broad variety of biologically relevant PET radiopharmaceuticals including $^{11}\text{C}\text{-Methionine}$, $^{11}\text{C}\text{-Choline}$ and $^{11}\text{C}\text{-Raclopride}$ using the E&Z Modular-Lab Methylation system in the third step of the process. The crude $^{11}\text{C}\text{-Raclopride}$ produced is purified further in the final step of the process using the E&Z Modular-Lab HPLC system.

The E&Z Modular-Lab Methylation system can also be used for the synthesis of $^{11}\text{C}\text{-Acetate}$ directly from $^{11}\text{C}\text{-CO}_2$ ($^{11}\text{C}\text{-CO}_2$ is trapped by Mel/MeOTf system) produced by a PET cyclotron.



E&Z Modular-Lab MeI/MeOTf System

Product Description

^{11}C -Methyl Iodide and ^{11}C -Methyl Triflate can be easily produced with the fully automated Modular-Lab synthesis device. The synthesis occurs via a gas-phase* approach that results in the production of labeled compounds with high radiochemical yield and specific activity.

The ^{11}C - CO_2 produced by the cyclotron is reversibly adsorbed to molecular sieves. After controlled release via heating, the adsorbed ^{11}C - CO_2 is converted into ^{11}C -Methane by hydrogenation of ^{11}C - CO_2 in a palladium oven. The resulting ^{11}C -Methane is then adsorbed and processed to ^{11}C -Methyl Bromide, which can be converted to ^{11}C -Methyl Iodide using a sodium iodide oven or ^{11}C -Methyl Triflate using a silver triflate oven. The resulting products, ^{11}C -Methyl Iodide or ^{11}C -Methyl Triflate, are then available for tracer synthesis.

Key features

- Fully automated synthesis process with no user intervention necessary
- Fully automated cleaning routine after each process to help ensure minimum chemical contamination of the system
- Traceability of the complete process, including documentation of all process parameters and functions
- Short synthesis and pre-preparation time
- Turnaround time between successive batches of ^{11}C -Methyl Iodide: 20 to 30 minutes
- Computer for system control
- Modular-Lab Software Professional – GMP, GAMP5 and 21 CFR part 11 compliant
- Trapping of ^{11}C - CO_2 and controlled release for ^{11}C -Acetate

Benefits

- Higher yields and higher specific activity compared to the liquid-phase method
- Multiple batches (>6 runs per day) without need for apparatus replacement or setup
- No liquid nitrogen necessary for the process
- No air and moisture sensitive chemicals necessary

Onsite installation and training are provided by Eckert & Ziegler.

*The gas-phase ^{11}C -Methyl Iodide setup was developed and validated in close cooperation with Bruce H. Mock Ph.D., Associate Professor of Radiology at Indiana University (USA).



E&Z Modular-Lab Methylation System

Product Description

The Modular-Lab Methylation system that is based upon the Modular-Lab PharmTracer synthesis technology allows for the easy and fully automated synthesis of ^{11}C -Methionine, ^{11}C -Choline, ^{11}C -Acetate and ^{11}C -Raclopride, as well as other radiopharmaceuticals, using disposable single-use cassettes. These cassettes guarantee a sterile fluid path and eliminate cross-contamination.

Key features

Modular-Lab Methylation system features

- Modular-Lab PharmTracer system – 4 fold cassette
- Heater module (Max 220°C + with air cooling)
- Syringe module with holder for 20 ml BD syringe
- Activity detector plug (activity max 5-Ci)
- Operates off the same computer as Mel/MeOTf system


Onsite installation and training are provided by Eckert & Ziegler.

Modular-Lab PharmTracer cassettes

- Cassettes are assembled under GMP-compliant clean room conditions, sterilized with gamma-radiation and double vacuum-packed.
- All consumables used are chemical resistant and have been tested for their suitability with the specific syntheses. All materials are free of animal-derived ingredients.
- A shelf life of 12 months can be guaranteed. Due to the cassettes one-time use, no reconditioning and cleaning or sanitation routines are necessary.
- Sterile cassettes are available for ^{11}C -Choline, ^{11}C -Methionine, ^{11}C -Raclopride and ^{11}C -Acetate

A vacuum-packed cassette unit consists of three components:

- Bag A with sterile cassette
- Bag B with sterile vials
- Bag C with sterile accessories (syringes, needles, etc.)



Bag A



Bag B



Bag C



E&Z Modular-Lab HPLC System

Product Description

The Modular-Lab HPLC system is designed to purify the reaction product of a radiopharmaceutical synthesis by semi-preparative HPLC (needed for ^{11}C -Raclopride). The system comprises integrated parts that include injection valve, a sample loop, a fluid detector and an activity detector.

Key features

- Knauer HPLC pump (10 ml pump head)
- Knauer HPLC UV-Detector (fixed wavelength)
- Analytical Flow Cell
- Operates off the same computer as Mel/MeOTf system

Onsite installation and training are provided by Eckert & Ziegler.

E&Z ¹¹C Modular-Lab Chemistry Solution

System Specifications

Compound	Synthesis Time (minutes)	Decay Corrected Yield	Uncorrected Yield	Radiochemical Purity
¹¹ C-Methyl Iodide*	15	>67%	>40%	>95%
¹¹ C-Methyl Triflate*	15	>60%	>30%	>95%
¹¹ C-L-Methionine†	25	>35%	>21%	>95%
¹¹ C-Choline†	25	>61%	>26%	>98%
¹¹ C-Raclopride	35 + 20 min HPLC purification	>10%	>5%	>85%
¹¹ C-Acetate*	17	>57%	>27%	>95%

Modular-Lab MeI/MeOTf System

Dimensions

Product Dimensions 19.7 in (50 cm) wide × 20.1 in (51 cm) high × 19.7 in (50 cm) deep

External support equipment

Electrical Cabinet 11.8 in (30 cm) wide × 15.7 in (40 cm) high × 7.9 in (20 cm) deep;
11.8 in (30 cm) clearance at front, 5.9 in (15 cm) clearance above

Heat Exchanger 10.2 in (26 cm) wide × 14.6 in (37 cm) high × 14.2 in (36 cm) deep;
5.9 in (15 cm) clearance at back; should not be installed more than
1 meter height difference from system

Power Supply 115 or 230 VAC, 1000 W

Modular-Lab Methylation System and Modular-Lab HPLC System

Dimensions

Product Dimensions 15.7 in (40 cm) wide × 18.5 in (47 cm) high × 14.6 in (37 cm) deep

External support equipment (only with HPLC System)

HPLC Pump 4.7 in (12 cm) wide × 3.9 in (10 cm) high × 8.7 in (22 cm) deep

UV Detector 8.9 in (22.6 cm) wide × 16.1 in (41 cm) high × 5.3 in (13.5 cm) deep, 6 kg

Power Supply 115 or 230 VAC, 1000 W

All modules are supplied with a one year manufacturer's warranty from Eckert & Ziegler.

Service and support are provided by Eckert & Ziegler.

Extended service contracts and additional consumables such as the synthesis cassettes can be purchased from Eckert & Ziegler.

Trademarks and service marks used in this material are property of Siemens Medical Solutions USA, Inc. or Siemens AG. All other company, brand, product and service names may be trademarks or registered trademarks of their respective holders.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

All rights reserved
© 2012 Siemens Medical Solutions USA, Inc
MI-405.TM.JV PDF only 05.2012

Global Business Unit

Siemens Medical Solutions USA, Inc.
Molecular Imaging
2501 N. Barrington Road
Hoffman Estates, IL 60192-2061
USA
Telephone: +1847 304 7700
www.siemens.com/mi

Legal manufacturer

Eckert & Ziegler Eurotope GmbH
Robert-Rössle-Strasse 10
13125 Berlin
Germany
Telephone: +49-30-94 10 84 0
Fax: +49-30-94 10 84 112

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Munich
Germany

Global Siemens Headquarters Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Telephone: +49 9131 84-0
www.siemens.com/healthcare