

CRT Implantation Without Lead Apron

Given the duration of less than six minutes, the average fluoroscopy time for a complete MediGuide Technology-controlled CRT implantation is far lower than conventional exposure.

By Hildegard Kaulen
Photos: Andréas Lang

Hindricks likes to chart new territories. He is considered one of the pioneers of cardiac catheter ablation and manages one of the worldwide leading centers for rhythmology. Each year, his department performs 1,200 atrial fibrillation catheter ablations and around 500 CRT implantations. Recently, Hindricks started to break new ground. About nine months ago, he began implanting biventricular pace-makers using MediGuide™ Technology integrated into his Artis® zee system – and, thus far, the Heart Center in Leipzig, Germany, is the only clinic in the world doing it. Hindricks is aiming to significantly reduce radiation exposure for patients and interventionalists, and seeks to further improve the quality of interventions. If the electrodes of the biventricular pace-maker are optimally positioned, cardiac status can improve and morbidity and

mortality decrease almost without side effects – in fact for many years. To date, Hindricks and his team have completed fifty MediGuide Technology-supported CRT implantations. They would have completed more but the supply of necessary equipment like sheaths and guidewires for non-fluoroscopic catheter navigation wasn't in full swing until recently. "Cardiac resynchronization therapy is the most important innovation in heart failure treatment in the past 10 years," Hindricks explains. "All clinical studies proved that patients live longer and better with it, since the technique leads to an improvement in cardiac function. And, there is a clear medical need for more CRT implantations." In Germany, 202 treatments per million inhabitants occur each year, which ranks as number two in the European care statistics. The majority of eligible patients, however, do

Cardiac Resynchronization Therapy (CRT)

CRT is indicated and effective for:

- Advanced drug-refractory heart failure
- Severe LV systolic dysfunction
- Ventricular dyssynchrony

CRT improves:

- Symptoms of heart failure
- Exercise capacity
- Ventricular function and structure
- Heart failure-related hospitalizations
- Mortality risk



not necessarily receive this treatment. According to Hindricks, this is due to a lack of awareness among referring physicians. In many countries, the technique is not covered by insurers, and electrophysiological (EP) centers, as well as appropriately trained physicians, are not available. "We must overcome existing roadblocks and develop better screening programs," he comments. "The left ventricular electrode must be placed in a safe and suitable position. If complications arise, they normally are related to the intervention and include dislocations or infections," he adds.

How Does the MediGuide Technology Work?

MediGuide Technology is not a new mapping system but a cardiovascular treatment platform for fluoroscopy-free navigation. It can be used whenever devices are navigated under fluoroscopic control inside the body. The position of the device is determined by magnetic tracking, and projected in real-time to a pre-recorded cine loop. "The characteristics of MediGuide Technology are its built-in compensation algorithms," Hindricks

explains. "The devices are depicted independent from primary and secondary organ movements. This provides high precision over time, resulting in 4D catheter tracking. Moreover, the fluoroscopy system and the electromagnetic sensor field do not work independently, but are fully registered." Compensation of respiratory and cardiac movement is enabled by a reference sensor in the patient's sternum. This sensor acts as a reference point for the calculation of the catheter position. The electromagnetic field transmitter is mounted inside the fluoroscopy detector. Conventional CRT implantations are associated with considerable radiation exposure. This is caused by the anatomically and technically challenging implantation of the left ventricular electrode, a process that usually requires continuous fluoroscopy. During this procedure, the interventionist stands near the left shoulder of the patient and probes the coronary sinus by means of left anterior oblique (LAO) projections, thus being exposed to considerable scattered radiation. So, they would also benefit from a reduction in the duration of fluoroscopy procedures. Moreover, the C-arm would not interfere

with the interventionist's work, as the patient is no longer exposed to fluoroscopy during the intervention. Consequently, the implanter is able to navigate in an LAO projection without a C-arm occupying this position. "The catheter can be visualized on pre-recorded cine loops even during prolonged periods of examination. We can, therefore, have fluoroscopy performed only at the beginning and at the end of a CRT implantation," Hindricks says. "I dream of an electrophysiological laboratory that doesn't use lead aprons. I think that we will be able to implement this vision in the very near future."

Study Proves Decrease of Radiation Exposure

There is substantive data to back-up Hindricks' vision. During the 50 MediGuide Technology-supported CRT implantations completed, radiation exposure could be reduced significantly compared to the conventional procedure. The average fluoroscopy period for LV electrode positioning was less than three minutes, and for the entire implantation less than six minutes, without any increase in the rate of complications. "I am absolutely convinced that we can reduce the fluoroscopy duration even further," the electrophysiologist says. "We still are in an early stage and use fluoroscopy during the intervention for safety-related reasons. In case of an uncomplicated anatomy, already today a catheter can be navigated safe and reliably without fluoroscopy. In case of stenoses and anomalies of the coronary sinus, we still need it as a safety precaution. As soon as improved algorithms of navigation are available, this will no longer be necessary." So which MediGuide Technology-supported devices are already available for Hindricks? During CRT implantation, he uses a CS sheath and a subselector, both equipped with a MediGuide Technology sensor. The electrophysiologist reaches the target coronary vein by means of a guidewire of 0.14 mm diameter with the MediGuide Technology sensor, which is used to place the electrode in a suitable tributary of the main cardiac vein. The electrode itself is not equipped with a sensor. "This will not be changed in the near future," he says.





“I dream of an electrophysiological laboratory without lead apron. I think that we will be able to implement this vision in the very near future.”

Professor Gerhard Hindricks, MD
Heart Center of Leipzig University, Leipzig, Germany

“The electrode shall retain its proven design.” Additionally, sensor-equipped diagnostic EP catheters are available as needed.

While the radiation exposure of MediGuide Technology-supported CRT implantations is considerably lower than in conventional procedures, the time expenditure is actually higher. Hindricks and his team need about 20 additional minutes (140 minutes compared to the average CRT implantation time of 120 minutes). “The increased time still reflects our study situation,” he comments. “We enter new territories and are the first and only clinic in the world utilizing this procedure in CRT. With this in mind, the additional time expenditure is remarkably low. Nevertheless, intervention times are crucial for the hospitals. They should at least be comparable, if not even better than for the conventional procedure. I know that we have to improve, and I am convinced we will succeed.”

Short Training Curve

Will other interventionalists adapt to MediGuide Technology? “The system offers one of the shortest learning curves

you can imagine,” Hindricks says. “We had colleagues who were familiar with it already after 10 minutes time, which was not only due to the fact that they were excellent electrophysiologists. Due to the projection of the catheter onto a conventional fluoroscopic image, we stay within an environment that is well-known to all users. They obtain additional imaging quality, but no basically new views. This helps them to quickly become accustomed to the MediGuide Technology system.”

While CRT facilitates a longer and better life for most patients, one quarter of the patients equipped with such a device are non-responders. Hence, Hindricks also uses the MediGuide Technology to understand the reasons for treatment failure. “Non-responders presumably require an individualized therapy,” the electrophysiologist says, “not standard therapy.” Therefore, he measures the phase shift between the electrically and mechanically activated myocardial areas. This shift reduces the heart’s pumping efficiency. Since MediGuide Technology depicts this time delay at high temporal resolution,

Hindricks and his team are able to measure the late points and use them to determine the position of the left ventricular electrode in an appropriate coronary vein, which is optimal for the patient. For this approach, they are currently collecting data. First results are promising.

Author: Hildegard Kaulen, PhD, is a molecular biologist. Following positions at Rockefeller University in New York and Harvard Medical School in Boston, she has worked as a freelance scientific journalist for notable daily newspapers and scientific magazines since the mid-nineties.

The outcomes achieved by the Siemens customers described herein were achieved in the customer’s unique setting. Since there is no “typical” hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that others will achieve the same results.

Contact
bernhard.fahn@siemens.com

Heartbeat – IMPRINT
© 2013 by Siemens AG,
Berlin and Munich
All Rights Reserved

Publisher:
Siemens AG
Healthcare Sector
Customer Solutions Division
Henkestrasse 127, 91052 Erlangen, Germany

Chief Editor: Okan Ekinçi (responsible for content)

Siemens AG, Healthcare Sector, Henkestrasse 127,
91052 Erlangen, Germany
email: editor.heartbeat.healthcare@siemens.com

Project Management: Ingrid Schuhmann, Siemens Healthcare

Contributors to this issue: Janina Beilner, Tanja Berbalk,
Philip Buckle, Julie Clark, Anne-Eloise Cournut,
Marco de Andrade, Monika Demuth, Annette Deubel,
Lars Drueppel, Catherine Eby, Bernhard Fahn,
Laurence Friedland, Friedrich Fuchs, Nancy Gunther,
Hanno Herrmann, Yuan Jiang, Carsten Lauer,
Susan Mortifoglio, Andrea McMurray, Alessandro Ortisi,
Michael Pickart, Doris Pischitz, Stephanie Scharff,
Cornelia Schaub, Gudrun Schubmann, Dietrich Till
(all Siemens Healthcare).

PrePress: Clemens Ulrich, Siemens Healthcare
Production: Michael Brumme, Siemens Healthcare

Includes text and photographic contributions from Primafila
AG, CH-Zurich / independent Medien-Design, D-Munich

The entire editorial staff here at Siemens Healthcare extends their
appreciation to all the experts, radiologists, scholars, physicians and
technicians, who donated their time and energy – without payment –
in order to share their expertise with the readers of Heartbeat.

Design and Editorial Consulting: Carmen Weith,
Cavok Consulting s.a.r.l, 67470 Wintzenbach, France

Printer: G. Peschke Druckerei GmbH, Schatzbogen 35,
81829 Munich, Germany

Note in accordance with § 33 Para.1 of the German Federal
Data Protection Law: Dispatch is made using an address
file which is maintained with the aid of an automated data
processing system.

We remind our readers that when printed, X-ray films never
disclose all the information content of the original. Artifacts in
CT, MR, ultrasound, and DSA images are recognizable by their
typical features and are generally distinguishable from existing
pathology. As referenced below, healthcare practitioners are
expected to utilize their own learning, training, and expertise
in evaluating images.

Partial reproduction in printed form of individual contributions
is permitted, provided the customary bibliographical data, such
as author's name and title of the contribution as well as date
and pages of Heartbeat, are named.

The editors request that two copies be sent to their attention.
The consent of the authors and editors is required for the
complete reprint of an article. Manuscripts submitted without
prior agreement as well as suggestions, proposals, and infor-
mation are always welcome; they will be carefully assessed
and submitted to the editorial conference for attention.

Find more about sustainable cardiovascular care under:
www.siemens.com/cardiology



DISCLAIMERS: Practice of Medicine: "The information presented in this magazine is for illustration only and is not intended to be relied upon by the reader for instruction as to the practice of medicine. Healthcare practitioners reading this information are reminded that they must use their own learning, training, and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens Healthcare to be used for any purpose in that regard." Contrast Agents: "The drugs and doses mentioned herein are consistent with the approved labeling for uses and/or indications of the drug. The treating physician bears the sole responsibility for the diagnosis and treatment of patients, including drugs and doses prescribed in connection with such use. The Operating Instructions must always be strictly followed when operating your Siemens system. The source for the technical data are the corresponding data sheets." Trademarks: "All trademarks mentioned in this document are property of their respective owners." Results: "The outcomes achieved by the Siemens customers described herein were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that others will achieve the same results."

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
www.siemens.com/healthcare

www.siemens.com/healthcare-magazine

A91CC-0003-11M1-7600 | Printed in Germany
CC CRM 1550 081310. | © 08.2013, Siemens AG

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Local Contact Information

Asia/Pacific:

Siemens Medical Solutions
Asia Pacific Headquarters
The Siemens Center
60 MacPherson Road
Singapore 348615
Telephone: +65 9622-2026

Canada:

Siemens Canada Limited
Healthcare Sector
1550 Appleby Lane
Burlington, ON L7L 6X7
Canada
Telephone + 1 905 315-6868

Europe/Africa/Middle East:

Siemens AG, Healthcare Sector
Henkestr. 127,
91052 Erlangen
Germany
Telephone: +49 9131 84-0

Latin America:

Siemens S.A., Medical Solutions
Avenida de Pte. Julio A. Roca No
516, Piso 7
C1067ABN Buenos Aires
Argentina
Telephone: +54 11 4340-8400

USA:

Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355-1406
USA
Telephone: +1 888 826-9702