

Supported by*

- syngo Dyna3E
- syngo DynaCT
- syngo DynaCT 360
- syngo InSpace 3D/3D Fusion
- syngo iPilot enhanced
- syngo iGuide Toolbox
- syngo iGuide
- syngo DynaPBV Neuro
- syngo Dyna PBV Body
- syngo Embolization Guidance
- syngo iFlow
- syngo Advanced Roadmap
- syngo iDentify

Courtesy of

PD Dr. B. Meyer, Institute for Diagnostic and Interventional Radiology, Medical School Hannover, Germany

System & Software

Artis Q ceiling VD10 syngo X Workplace VC10

*This list of applications is not complete. Not all applications available for all software versions.

Case Description

Patient History

Patient with history of recurrent pulmonary embolism and chronic thromboembolic pulmonary hypertension (CTEPH). Pulmonary angiogram and *syngo* DynaCT were acquired for diagnostic work-up.

Diagnosis

CTEPH with perfusion defects and web stenoses mainly in the right lung.

Treatment

Diagnostic work-up.

General Comments

Pulmonary angiograms are still the established gold standard for diagnostic work-up in patients with CTEPH. *syngo* DynaCT can provide additional information about chronic embolism by depicting web stenoses and occlusions.

Tips and Tricks

Selective *syngo* DynaCT imaging can be performed in case of insufficient contrast using a central catheter position.

syngo DynaCT can provide additional information about chronic embolism by depicting web stenoses and occlusions.

syngo DynaCT of the pulmonary arteries

Acquisition Protocol	6sDCT Body
Injection Protocol	
Contrast Media (CM):	300 mg lodine/ml
Dilution (CM/Saline):	70%/30%
Injection Volume:	60 ml
Power Injector Used:	Yes
Injection Rate:	8 ml/s
Duration of Injection:	7.5 s
X-ray Delay:	1.5 s
Catheter Position:	Central Catheter Position


Reconstruction Protocol	DynaCT Body Nat Fill HU
VOI Size:	Full
Slice Matrix:	512x512
Kernel Type:	HU
Image Characteristics:	Normal
Reconstruction Mode:	Nat Fill
Viewing Preset:	DynaCT Body
Secondary Recon	No



Thick MIP 23 mm Web stenosis in lower right lobe (arrow)



Thick MIP 23 mm Small Web stenosis in left lobe (arrow)



Thin MPR

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.

All rights reserved.

Siemens Healthcare Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany

Phone: +49 9131 84-0 siemens.com/healthcare

PDF only. | © Siemens Healthcare GmbH, 2016