

# Cardiac Diffusion Tensor MRI Using Simultaneous Multi-Slice Acquisition with a Blipped-CAIPIRINHA Readout

Choukri Mekkaoui<sup>1</sup>; Timothy G. Reese<sup>1</sup>; Marcel P. Jackowski<sup>2</sup>; Himanshu Bhat<sup>3</sup>; David E. Sosnovik<sup>1,4</sup>

<sup>1</sup> Athinoula A. Martinos Center for Biomedical Imaging, Department of Radiology, Massachusetts General Hospital, Harvard Medical School, Boston, MA, USA

<sup>2</sup> Department of Computer Science, Institute of Mathematics and Statistics, University of São Paulo, São Paulo, Brazil

<sup>3</sup> Siemens Healthcare, Charlestown, MA, USA

<sup>4</sup> Cardiovascular Research Center, Cardiology Division, Massachusetts General Hospital, Harvard Medical School, Boston, MA, USA

## Background


Heart muscle is highly anisotropic with an intricate microstructure, well suited to characterization with diffusion tensor imaging (DTI) [1]. The most widely used measure of fiber organization in the myocardium is the helix angle (HA), simply defined as the inclination of the myofiber out of the local short-axis plane. Myofibers in the subendocardium have a positive HA, while those in the subepicardium have a negative HA [1, 2]. These myofibers are further arranged into laminar sheets, which slide against each other allowing the myocardium to thicken during systole [3, 4]. Alterations of this microstructure due to heart disease affect its mechanical efficiency and also could contribute to arrhythmias [5, 6]. These microstructural changes can precede symptoms and thus a non-invasive evaluation could be of significant clinical value.


The motion of the heart is five orders of magnitude greater than the self-diffusion of water. Therefore, approaches that are sensitive to the microscopic diffusion of water, but not to cardiac motion and strain, are needed for successful *in vivo* imaging [7-9]. One approach to enable *in vivo* DTI uses a diffusion-encoded stimulated echo (STE) sequence (Fig. 1), which can be implemented on most clinical scanners [10, 11]. The diffusion-encoded STE sequence is played out over two successive heartbeats. The first and second 90° excitation pulses are applied in the first heart-

beat, and the third excitation pulse in the second heartbeat. The diffusion-encoding gradients are monopolar and placed immediately after the first and third excitation pulses. The main appeal of the STE approach is its conceptual immunity to cardiac motion. Ideally, each monopolar diffusion gradient occurs at exactly the same time in two sequential R-R intervals. Hence, not only is the phase due to the diffusion-encoding gradient unwound, but the influence of cardiac motion on the phase of the magnetization also is unwound.

The use of a STE sequence, however, introduces a high degree of inefficiency into the acquisition due to the dual-gated acquisition. Additionally, the STE is half the amplitude of a spin-echo. Consequently most inves-

tigators have used ~8 averages per slice in order to achieve sufficient signal-to-noise (SNR) during diffusion-encoded STE acquisitions, taking 5-7 minutes per slice [12]. The inefficiency of the STE approach frequently requires the anatomical coverage of the acquisition to be compromised. For instance, only 3 short-axis slices can be imaged in ~15 minutes, which covers only 25% of the myocardium [12]. New approaches to improve anatomical coverage and reduce scan time are thus sorely needed. The development of simultaneous multi-slice (SMS) acquisition using a blipped Controlled Aliasing in Parallel Imaging (blipped-CAIPIRINHA) readout holds great promise [13, 14], and could play a key role in facilitating the more widespread use of cardiac DTI.





**2** (2A) Simultaneous acquisition of 3 slices with rate 3 SMS. The gap between the slices is 500% of slice thickness. The tensor field is represented by supertoroids color-coded by HA. Supertoroid fields resulting from no SMS (2B), rate 2 (2C), and rate 3 (2D) SMS acquisitions of the same mid-ventricular slice are consistent with the transmurality change in HA from positive in the subendocardium to negative in the subepicardium.

## Implementation

The technical details of SMS excitation using blipped-CAIPIRINHA have been described in detail elsewhere in this volume. The technique has been used extensively in the brain [13, 14], and preliminary experience with it in the heart appears promising [15]. In the current article, we describe our experience with this technique for cardiac DTI in healthy volunteers.


Breath-hold DTI was performed on a clinical 3T scanner (MAGNETOM Skyra, Siemens Healthcare, Erlangen, Germany) with a 34-element receive coil<sup>1</sup> (18 anterior and 16 posterior elements). Images were acquired with a diffusion-encoded STE sequence, which was volume-selected in the phase-encode axis using a slab selective radiofrequency (RF) pulse. Acquisition parameters included: FOV 360 x 180 mm, resolution 2.5 x 2.5 x 8 mm<sup>3</sup>, in-plane GRAPPA rate 2, TE 34 ms, b-value 500 s/mm<sup>2</sup>, 10 diffusion-encoding directions, and 8 averages. Twelve short-axis slices were acquired in the systolic sweet spot of the cardiac cycle to mitigate strain effects [11, 16]. Imaging was performed with no SMS, and rates 2 and 3 SMS. HA was derived from the diffusion tensor, which was estimated from the diffusion-weighted images. Fiber tracts were constructed by integrating the primary eigenvector field into streamlines using an adaptive 5<sup>th</sup> order Runge-Kutta approach [5].

## Results and impact

With no SMS, 96 breath-holds were required to cover the entire LV. Using rate 2 SMS, this was reduced to 48 breath-holds, and with rate 3 SMS to 32. With rate 3 SMS, the acquisition time was approximately 20 minutes for whole-heart coverage. Image quality was well preserved using both rates 2 and 3 SMS. This is demonstrated in Figure 2, where the diffusion tensor in each voxel is represented by the supertoroidal model [17]. The glyphs are parameterized by the magnitude and orientation derived from the diffusion tensor and color-coded by HA. The transmural evolution in HA from positive in the subendocardium to negative in the subepicardium is well resolved in all 3 slices with rate 3 SMS. Glyph fields using rates 2 and 3 SMS of a mid-ventricular slice compare favorably with those acquired with no SMS, and are consistent with expected transmural evolution in HA.

Tractography of the heart has previously been performed over a small anatomical range (3-5 slices) or with very large slice gaps. Meaningful tractography requires the entire heart to be imaged without any slice gaps.

<sup>1</sup> The product is still under development and not commercially available yet. Its future availability cannot be ensured.



**3** Tractography of the entire LV, color-coded by HA, of the same subject imaged with no SMS, and rates 2 and 3 SMS. Tracts obtained using rates 2 and 3 SMS compare favorably and are in agreement with those obtained with no SMS.

With no SMS, this takes over 60 minutes to acquire. However, as shown in Figure 3, fiber tracts of the entire LV were successfully obtained with rates 2 and 3 SMS, and are qualitatively comparable with those obtained with no SMS.

## Discussion

DTI of the heart has the potential to improve the understanding, diagnosis and management of a range of cardiovascular diseases. However, the main limitation is long scan times. With current techniques, the acquisition of 3 short-axis slices takes ~20 minutes. We demonstrated that using SMS, scan time was reduced by 3-fold. The simultaneous acquisition of 3 slices (basal, medial, and apical), as shown in Figure 2, takes ~5 minutes. While imaging only 3 slices yields limited coverage of the LV, the utility of this approach has been demonstrated in first-pass perfusion studies of the heart [18]. The additional acquisition of DTI images in the same 3 short-axis slices would add little time to a clinical study and could be of substantial value.

DTI of the entire heart, while more demanding, could provide a unique way to evaluate myocardial microstructure. SMS combined with further technical advances may facilitate the clinical translation of whole-heart DTI, enabling the reliable characterization of myocardial structure in a wide range of patients with cardiac diseases.

## References

- 1 Streeter DD, Jr., Spotnitz HM, Patel DP, Ross J, Jr., Sonnenblick EH. Fiber orientation in the canine left ventricle during diastole and systole. *Circ Res*. 1969;24(3):339-47.
- 2 Scollan DF, Holmes A, Winslow R, Forder J. Histological validation of myocardial microstructure obtained from diffusion tensor magnetic resonance imaging. *Am J Physiol*. 1998;275(6 Pt 2):H2308-18.
- 3 LeGrice IJ, Smaill BH, Chai LZ, Edgar SG, Gavin JB, Hunter PJ. Laminar structure of the heart: ventricular myocyte arrangement and connective tissue architecture in the dog. *Am J Physiol*. 1995;269(2 Pt 2):H571-82.
- 4 Dou J, Tseng WY, Reese TG, Wedeen VJ. Combined diffusion and strain MRI reveals structure and function of human myocardial laminar sheets in vivo. *Magn Reson Med*. 2003;50(1):107-13.
- 5 Mekkaoui C, Huang S, Chen HH, Dai G, Reese TG, Kostis WJ et al. Fiber architecture in remodeled myocardium revealed with a quantitative diffusion CMR tractography framework and histological validation. *J Cardiovasc Magn Reson*. 2012;14:70.
- 6 Trayanova NA. Whole-heart modeling: applications to cardiac electrophysiology and electromechanics. *Circ Res*. 2011;108(1):113-28.
- 7 Sosnovik DE, Wang R, Dai G, Reese TG, Wedeen VJ. Diffusion MR tractography of the heart. *J Cardiovasc Magn Reson*. 2009;11:47.
- 8 Nguyen C, Fan Z, Sharif B, He Y, Dharmakumar R, Berman DS et al. In vivo three-dimensional high resolution cardiac diffusion-weighted MRI: a motion compensated diffusion-prepared balanced steady-state free precession approach. *Magn Reson Med*. 2014;72(5):1257-67. doi:10.1002/mrm.25038.
- 9 Gamper U, Boesiger P, Kozierke S. Diffusion imaging of the in vivo heart using spin echoes--considerations on bulk motion sensitivity. *Magn Reson Med*. 2007;57(2):331-7. doi:10.1002/mrm.21127.
- 10 Reese TG, Weisskoff RM, Smith RN, Rosen BR, Dinsmore RE, Wedeen VJ. Imaging myocardial fiber architecture in vivo with magnetic resonance. *Magn Reson Med*. 1995;34(6):786-91.
- 11 Tseng WY, Reese TG, Weisskoff RM, Wedeen VJ. Cardiac diffusion tensor MRI in vivo without strain correction. *Magn Reson Med*. 1999;42(2):393-403.
- 12 Nielles-Vallespin S, Mekkaoui C, Gatehouse P, Reese TG, Keegan J, Ferreira PF et al. In vivo diffusion tensor MRI of the human heart: reproducibility of breath-hold and navigator-based approaches. *Magn Reson Med*. 2013;70(2):454-65.
- 13 Setsompop K, Cohen-Adad J, Gagoski BA, Raji T, Yendiki A, Keil B et al. Improving diffusion MRI using simultaneous multi-slice echo planar imaging. *Neuroimage*. 2012;63(1):569-80. doi:10.1016/j.neuroimage.2012.06.033.
- 14 Setsompop K, Gagoski BA, Polimeni JR, Witzel T, Wedeen VJ, Wald LL. Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty. *Magn Reson Med*. 2012;67(5):1210-24. doi:10.1002/mrm.23097.
- 15 Lau AZ, Tunncliffe EM, Frost R, Koopmans PJ, Tyler DJ, Robson MD. Accelerated human cardiac diffusion tensor imaging using simultaneous multislice imaging. *Magn Reson Med*. 2015;73(3):995-1004. doi:10.1002/mrm.25200.
- 16 Stoeck CT, Kalinowska A, von Deuster C, Harmer J, Chan RW, Niemann M et al. Dual-phase cardiac diffusion tensor imaging with strain correction. *PLoS One*. 2014;9(9):e107159. doi:10.1371/journal.pone.0107159.
- 17 Mekkaoui C, Chen IY, Chen HH, Kostis WJ, Pereira F, Jackowski MP et al. Differential response of the left and right ventricles to pressure overload revealed with diffusion tensor MRI tractography of the heart in vivo. *Journal of Cardiovascular Magnetic Resonance*. 2015;17(Suppl 1):O3-O. doi:10.1186/1532-429X-17-S1-O3.
- 18 Stab D, Wech T, Breuer FA, Weng AM, Ritter CO, Hahn D et al. High resolution myocardial first-pass perfusion imaging with extended anatomic coverage. *J Magn Reson Imaging*. 2014;39(6):1575-87. doi:10.1002/jmri.24303.



Timothy G.  
Reese



Choukri  
Mekkaoui

## Contact

Choukri Mekkaoui  
Athinoula A. Martinos Center for  
Biomedical Imaging  
149 13<sup>th</sup> Street  
Charlestown, MA 02129, USA  
Phone: +1 617 724-3407  
Fax: 617 726-7422  
mekkaoui@nmr.mgh.harvard.edu