

# Rapid High Spatial Resolution Diffusion MRI at 7 Tesla Using Simultaneous Multi-Slice Acquisition

Cornelius Eichner<sup>1,2</sup>; Robin M. Heidemann<sup>3</sup>

<sup>1</sup> Max Planck Institute for Human Cognitive and Brain Sciences, Leipzig, Germany

<sup>2</sup> A. A. Martinos Center for Biomedical Imaging, MGH, Harvard Medical School, Boston, MA, USA

<sup>3</sup> Siemens Healthcare, Erlangen, Germany

## Abstract

High quality diffusion-weighted MRI (dMRI) data can be obtained with very high isotropic spatial resolution at ultra-high magnetic field strength such as 7 Tesla (T). Due to the high resolution it is necessary to acquire a large number of imaging slices for whole brain coverage, which results in long acquisition times (TA) of more than an hour. Obviously, Simultaneous multi-slice (SMS) acquisition technology is a prerequisite to significantly reduce the extensive acquisition times of these studies. However, the large energy deposition into the subject caused by the employed Multiband (MB) RF-pulses limit the efficiency of SMS methods. This can be addressed by replacing the

MB pulses by Power Independent of Number of Slices (PINS) pulses<sup>1</sup>. However, this comes at the expense of a reduced bandwidth and increased off-resonance dependency, which degrades the image quality. With a new RF pulse design, given the name MultiPINS<sup>1</sup>, the RF energy is further reduced and/or the pulse length is shortened. This is achieved by combining MB with PINS RF pulses.

*In vivo* dMRI results were recorded with the MultiPINS approach at high spatial resolutions at 7T showing a 3-fold scan time reduction.


## Introduction

Diffusion MRI (dMRI) is an essential tool in neuroscience to study the structural connectivity of the human brain *in vivo*. Due to the complex structure of the brain, it is necessary to acquire data with high isotropic spatial resolution. Sub-millimeter isotropic resolution dMRI of the human brain *in vivo* is feasible at 7T [1]. However, due to the high isotropic resolution it is necessary to acquire a large number of imaging slices, which results in long acquisition times (TA) of about an hour and more – a major limitation of this technique.


<sup>1</sup> The product is still under development and not commercially available yet. Its future availability cannot be ensured.

1


A) Multiband (MB) Pulse



B) PINS Pulse



C) MultiPINS Pulse



1 Example SMS RF pulses with resulting slice profiles; energy reductions that resulted from pulse settings used in this work.

Acquiring multiple slices simultaneously, and unfolding them using information from multi-channel coil arrays [2], can address this speed problem and shorten TA. The CAIPIRINHA approach [3] was developed to reduce the g-factor noise for SMS imaging and has been recently adapted to EPI acquisitions [4] (blipped-CAIPIRINHA). However, the high-energy deposition of Multiband (MB) pulses (Fig. 1A) that are typically used for SMS imaging limit the acquisition speed of SMS methods at 7T, due to SAR/power constraints – especially if high flip angles are employed [5]. For moderate slice acceleration factors (e.g. MB = 2), SMS dMRI data can be acquired also at 7T [6]. In the case of higher slice acceleration factors, the RF power will limit the acquisition speed at 7T.

Recently Norris et al. showed that a periodic slice excitation pattern, suitable for SMS acquisition, can be created without significant increase in power deposition by multiplying a single-slice RF pulse with a Dirac comb function, to end up in a Power Independent of Number of Slices (PINS) pulse [7] (Fig. 1B). This approach has been successfully applied to 7T for structural and functional spin-echo experiments [8] as well as RF power consuming sequences such as Turbo Spin-Echo [9]. In this study, we combine ZOOPPA<sup>1</sup>, an outer volume suppression (OVS) technique [1], for diffusion MRI at 7T with blipped-CAIPIRINHA [4] and PINS pulses [7] to obtain high spatial and angular resolution dMRI with significantly reduced TA. Furthermore, by combining MB and PINS RF excitation, we created a new 'Multi-PINS' RF pulse type<sup>1</sup> to further reduce power deposition for SMS excitation (Fig. 1C). For this new MultiPINS pulse, a MB RF pulse is first reshaped to follow the excitation *k*-space traversal of

the blipped gradient waveform of the PINS pulse. The reshaped MB pulse can then be mixed directly with the PINS pulse to create a MultiPINS pulse with suitable excitation characteristic for SMS imaging (Fig. 1C). To minimize SAR, an optimal mixing ratio between MB and PINS (0% being pure PINS and 100% being pure MB pulse along PINS gradient trajectory) can be easily determined empirically prior to acquisition.


## Methods

The high-resolution diffusion MRI data were acquired on a 7T whole-body MR scanner<sup>2</sup> (Siemens Healthcare, Erlangen, Germany) equipped with a 32-element head coil and a gradient system achieving a maximum amplitude of 70 mT/m with a slew-rate of 200 T/m/s. A Stejskal-Tanner diffusion-weighted EPI sequence [10] was modified to employ ZOOPPA OVS and SMS blipped-CAIPIRINHA. SMS is used to accelerate the acquisition by 3 folds, while ZOOPPA is used to reduce the imaging volume-of-interest and associated image distortion artifact. To reduce SAR/power deposition of the SMS method, a PINS/MultiPINS pulse was utilized for RF refocusing. Energy calculations were performed for the MultiPINS SMS pulses, to find out optimal energy settings with short

pulse durations. *In vivo* diffusion-weighted images with 99 slices of 1 mm and 75 slices of 1.4 mm isotropic resolutions were acquired at 60 diffusion directions with a b-value of 1000 s/mm<sup>2</sup> and 7 interspersed b0 non-diffusion-weighted images (for motion correction). To increase the signal-to-noise ratio of the 1 mm isotropic resolution dataset, 4 averages were recorded. In the case of the 1 mm isotropic dataset, multiple fiber orientations were modeled with constrained spherical deconvolution followed by streamline fiber tracking with MRtrix (<http://www.brain.org.au/software/mrtrix/>). The 1.4 mm isotropic data set was acquired as a single average, resulting in a total acquisition time of about 3 minutes. Data were corrected for motion and eddy currents distortion artifacts with FSL and registered to a structural scan using Freesurfer. A color-coded FA map was calculated using the diffusion toolkit (<http://www.trackvis.org/>).

## Results and discussion

We have shown that dMRI data at 7T can be acquired in a significantly shortened acquisition time. We recorded a 1 mm isotropic resolution with 60 diffusion directions in just 6 minutes, by applying blipped-CAIPIRINHA and ZOOPPA OVS



2 Streamline fiber tracking of 100 000 fibers (5 mm slab) of coronal and axial brain slices. Four times averaged 7T dMRI data with 1 mm isotropic resolution.

<sup>1</sup> The product is still under development and not commercially available yet. Its future availability cannot be ensured.

<sup>2</sup> MAGNETOM 7T is ongoing research. All data shown are acquired using a non-commercial system under institutional review board permission. MAGNETOM 7T is still under development and not commercially available yet. Its future availability cannot be ensured.

(see Table 1) (Fig. 2). Furthermore, a whole brain dataset with 1.4 mm isotropic resolution and 60 diffusion directions was recorded in only 3:30 minutes (Fig. 3). We used PINS/MultiPINS pulses for refocusing to reduce SAR, and thus to gain the full benefit of SMS imaging at 7T.

## Conclusion

For ultra-high-resolution dMRI (1 mm isotropic or better) at 7T, the acquisition time becomes the major hindrance for a broad use of this application. We employed the blipped-CAIPIRINHA SMS technique in conjunction with PINS refocusing pulses as well as with a newly developed hybrid MB/PINS RF pulse to record dMRI data at 7T. The application of low power RF pulses in this SMS sequence enables the acquisition of high-resolution dMRI data at 7T in a significantly reduced scan time of 6:15 min for a 1 mm resolution dMRI scan and 3:30 min for a 1.4 mm isotropic resolution dMRI scan. The slice acceleration enables high-resolution dMRI acquisition at 7T within a time-frame short enough for clinical use, as well as combined *in vivo* anatomical, functional and diffusion studies at the same ultra-high-resolution level within a single scan session.


## References

- 1 R.M. Heidemann, A. Anwender, T. Feiweier, T. R. Knösche, R. Turner, "K-space and q-space: Combining ultra-high spatial and angular resolution in diffusion imaging using ZOOPPA at 7T," *NeuroImage*, 60(2), 967–978, 2012.
- 2 D.J. Larkman, J.V. Hajnal, A.H. Herlihy, G.A. Coutts, I.R. Young, G. Ehnholm, "Use of multicoil arrays for separation of signal from multiple slices simultaneously excited," *JMRI*, 13(2), 313–7, 2001.
- 3 F.A. Breuer, M. Blaimer, R.M. Heidemann, M.F. Müller, M.A. Griswold, P.M. Jakob, "Controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) for multi-slice imaging," *Magn Reson Med*, 53(3), 684–91, 2005.
- 4 K. Setsompop, B.A. Gagoski, J.R. Polimeni, T. Witzel, V.J. Wedeen, L.L. Wald, "Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty," *Magn Reson Med*, 67(5), 1210–24, 2012.
- 5 C. Eichner, "Combining ZOOPPA and blipped CAIPIRINHA for diffusion weighted imaging at 7T" In Proceedings of the 29th Annual Meeting ESMRMB, Lisbon, Vol. 4, p. 45111, 2012.
- 6 A.T. Vu, "High resolution whole brain diffusion imaging at 7 T for the Human Connectome Project", *NeuroImage* (2015), <http://dx.doi.org/10.1016/j.neuroimage.2015.08.004>.
- 7 D.G. Norris, P.J. Koopmans, R. Boyacioglu, M. Barth, "Power independent of number of slices (PINS) radiofrequency pulses for low-power simultaneous multislice excitation" *Magn Reson Med*, 66(5), 1234–40, 2011.
- 8 P.J. Koopmans, R. Boyacioglu, M. Barth, D.G. Norris, "Whole brain, high resolution spin-echo resting state fMRI using PINS multiplexing at 7 T," *NeuroImage*, 62(3), 1939–46, 2012.
- 9 D.G. Norris, R. Boyacioglu, J. Schulz, M. Barth, P.J. Koopmans, "Application of PINS radiofrequency pulses to reduce power deposition in RARE/turbo spin echo imaging of the human head," *Magn Reson Med*, 00(July), 2013.
- 10 E.O. Stejskal, J.E. Tanner, "Spin diffusion measurements: spin echoes in the presence of a time-dependent field gradient," *The Journal of Chemical Physics*, 42, 288, 1965.

Table 1

Resolution	1 mm w/o SMS	1 mm MB=3	1.4 mm MB=3
TR (ms)	13600	5000	3000
TE (ms)	68	64	69
Slices	99	99	75
FOV (mm <sup>2</sup> )	180 x 125	180 x 120	180 x 180
TA (min)	4 x 15:52	4 x 6:15	1x 3:30

Comparable diffusion-weighted protocols for slice accelerated and non-accelerated with SMS blipped-CAIPIRINHA. 60 diffusion directions,  $b = 1000 \text{ s/mm}^2$  with different slice acceleration factors (SMS).



3 Axial (3A) and coronal (3B) views of color coded FA diffusion data recorded at 1.4 mm isotropic resolution in 3:30 min.



## Contact

Cornelius Eichner  
Max Planck Institute for Human Cognitive and Brain Sciences  
Stephanstraße 1A  
04103 Leipzig  
Germany  
Phone: +49 341 9940-2428  
[ceichner@cbs.mpg.de](mailto:ceichner@cbs.mpg.de)