



SIEMENS



# Automatic QC System

One QC Batch for the Lifetime of Your Instruments

# AUTOMATIC QC (AQC) OVERVIEW

## Fully automatic quality control

- No operator intervention required
- Three levels of QC in each AQC cartridge

## Schedule is user-programmable

- Up to three QC levels per shift

## No value assignment sheet

- Value assignments programmed into cartridge
- Eliminates QC batch changeover issues

AQC material follows the same fluidic path as patient samples.



## TRADITIONAL QUALITY CONTROL


- Quality-control material is manufactured in batches.
- Each batch is tested and assigned values for all the analytes.
- The assigned values will show batch-to-batch differences.



## TRADITIONAL QUALITY CONTROL

Managing an instrument population requires knowledge of the assigned values for **every** QC batch at **each** level and **when** the QC batch changeovers occurred.


# TRADITIONAL QUALITY CONTROL



## TRADITIONAL QUALITY CONTROL

If you want to get an accurate picture of long-term performance, you must **normalize** the QC data.

# TRADITIONAL QUALITY CONTROL



# SIEMENS AUTOMATIC QUALITY CONTROL

**Unique** approach eliminates matching QC batch changes to QC data, **reduces waste**, and makes any QC program **easier to administer**.



## Siemens AQC Process

# 1

We mix each batch of reagent to **precise** analyte concentrations using **accurate, state-of-the-art process** equipment.

This produces AQC reagent with analyte values very **close** to the target value.



**1. Mix**

2. CALIBRATE

3. VALUE ASSIGN

4. CALCULATE



## Siemens AQC Process

# 2

Our blood-gas instrument pool is calibrated with reagents that are traceable to national standards.



1. MIX

**2. CALIBRATE**

3. VALUE ASSIGN

4. CALCULATE



## Siemens AQC Process

**3** We value assign **each** batch by **comparing** it to a master batch using our blood-gas instrument pool for a total of **144** measurements providing **over 2000** data points.



1. MIX

2. CALIBRATE

**3. VALUE ASSIGN**

4. CALCULATE

## Siemens AQC Process

**4** We calculate the offset of the value assignment to the batch target and program that offset **into** the AQC cartridge memory chip.



1. MIX

2. CALIBRATE

3. VALUE ASSIGN

**4. CALCULATE**

## SIEMENS AQC IN USE



1

When a customer installs an AQC cartridge, the offset values for all the analytes are read by the RAPIDPoint® 500 Blood Gas System.

2

When an AQC is analyzed, the measured values are compared to the assigned value, the offset is subtracted, and the AQC result reported.

Measured Value

Assigned Value



AQC Result

3

The outcome is the same as the normalizing exercise we just saw, but it is performed for **every** result.

## SIEMENS' UNIQUE APPROACH TO AQC

- Eliminates worry about matching QC batch changes to QC data.
- Reduces waste for institutions that throw away the remnants of a batch to ensure all their analyzers are using the same batch.
- Makes any QC program easier to administer.

It's like having only **one batch of QC** for the **lifetime** of your instruments.

