



The outcomes achieved by the Siemens customers described herein were achieved in the customers' unique settings. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that others will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features, and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications, and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

**Local Contact Information**  
Siemens Medical Solutions USA, Inc.  
51 Valley Stream Parkway  
Malvern, PA 19355-1406  
USA  
Phone: +1-888-826-9702  
www.usa.siemens.com/healthcare

**Global Business Unit**  
Siemens AG  
Medical Solutions  
Computed Tomography &  
Radiation Oncology  
Siemensstr. 1  
DE-91301 Forchheim  
Germany  
Phone: +49 9191 18-0  
www.siemens.com/healthcare

**Global Siemens Headquarters**  
Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

**Global Siemens Healthcare Headquarters**  
Siemens AG  
Healthcare  
Henkestrasse 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
www.siemens.com/healthcare

**Legal Manufacturer**  
Siemens AG  
Wittelsbacherplatz 2  
DE-80333 Muenchen  
Germany

Order No. A9111M-CDV-15665-P1-4A00 | Printed in USA 04-2015 | © 2015 Siemens Medical Solutions USA, Inc. | All rights reserved

[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

SIEMENS



Stellar Detector  
Benefits for the  
SOMATOM®  
Definition Flash  
and Definition AS+

[www.usa.siemens.com/CT](http://www.usa.siemens.com/CT)

# Stellar Detector

Excellent image quality at lower dose for most complex CT exams

## Stellar Detector— The First Fully Integrated Detector

The Stellar Detector in the SOMATOM® Definition Flash or Definition AS+ generates ultra-thin slices with very high spatial resolution for an extremely high level of visual detail. It provides ultra-low-signal imaging, which is invaluable, especially for obese and pediatric patients or Dual Energy applications.

Answers for life.

**Key benefits include:**

- High diagnostic image quality at lower dose for all patients and for most complex CT exams.
- Up to 60% dose reduction with SAFIRE\* technology.
- Generates ultra-thin slices to routinely deliver high spatial resolution images that visualize even the finest image details for procedures such as enhanced calcified plaque and stent analysis.
- Increased signal-to-noise ratio on obese patients.
- 0.3-mm spatial resolution.


\*In clinical practice, the use of SAFIRE may reduce CT patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task. The following test method was used to determine a 54-60% dose reduction when using the SAFIRE reconstruction software. Noise, CT numbers, homogeneity, low-contrast resolution, and high-contrast resolution were assessed in a Gammex 438 phantom. Low-dose data reconstructed with SAFIRE showed the same image quality compared to full-dose data based on this test. Data on file.

## Increased Sharpness and Reduced Blooming




## Less Noise for Improved Image Quality

## Improved Grey-White Matter Differentiation

## High-Resolution Imaging and Low kV CT



Courtesy of German Heart Center/Munich, Germany



Detector A: Excellent 80 kV series with HiDynamics

Detector B: 140 kV images with Selective Photon Shield

Courtesy of German Heart Center/Munich, Germany

The Stellar Detector provides high-quality imaging at a temporal resolution of 75 ms for sound calcified lesion evaluation and exclusion of in-stent restenosis in stents. The full electronic integration of the Stellar Detector elements minimizes cross-talk between neighboring detector rows. This significantly reduces slice blurring, resulting in a more precise slice profile, and generates 0.5-mm slices at a spatial resolution of 0.30 mm without increasing dose.

Adaptive Cardio Sequence  
temp resolution: 75 ms  
collimation: 128 x 0.6 mm  
spatial resolution: 0.30 mm  
scan length: 138 mm  
rotation time: 0.28 s  
120 kV, 296 mAs/rotation  
heart rate: 69 bpm  
CTDIvol: 11.0  
DLP: 152 mGycm  
eff. dose: 2.1 mSv

The full electronic integration of the Stellar Detector elements means that no analog electronic components or cables and connectors are mounted on the elements. This enables the Stellar Detector to minimize electronic noise, which significantly increases the signal-to-noise ratio for improved utilization of low signals at the detector, and supports improved image quality in all clinical examinations.

Flash Spiral\* Scanning  
collimation: 64 x 0.6 mm  
spatial resolution: 0.30 mm  
scan time: 2 s  
scan length: 500 mm  
rotation time: 0.28 s  
120 kV, 72 eff. mAs  
CTDIvol: 6 mGy  
DLP: 300 mGycm  
eff. dose: 4.4 mSv

\*Flash Spiral is only available on the SOMATOM Definition Flash.

The Stellar Detector improves visualization of bleeding, brain injury, and skull fractures in patients with head injuries with excellent grey-white matter differentiation for increased diagnostic certainty—while X-CARE protects the patient's eyes. It supports decision-making with unprecedented image sharpness without increasing image noise through Edge Technology. Additionally, TrueSignal Technology increases signal-to-noise ratio because of its more efficient X-ray quanta processing. Grey-white matter differentiation may also be used to visualize the early signs of a stroke.

X-CARE Acquisition  
collimation: 64 x 0.6 mm  
scan time: 8 s  
scan length: 90 mm  
rotation time: 1.0 s  
100 kV, 580 eff. mAs  
CTDIvol: 54 mGy  
DLP: 608 mGycm  
eff. dose: 1.3 mSv

The Stellar Detector with TrueSignal Technology makes HiDynamics imaging possible. Its full electronic integration significantly extends the dynamic bandwidth, which increases the image detail level especially for low-contrast objects in close proximity to high-contrast objects. This is important for low-kV scans, such as in the 80 kV data-set of Dual Energy scans or in 70 kV pediatric examinations.

Dual Energy Acquisition  
collimation: 64 x 0.6 mm  
spatial resolution: 0.30 mm  
scan time: 24 s  
scan length: 975 mm  
rotation time: 0.33 s  
80/140 kV, 121/53 eff. mAs  
CTDIvol: 4.3 mGy  
DLP: 419 mGycm  
eff. dose: 1.6 mSv