# Free-breathing DCE-MRI of the Kidney using GRASP

Philipp Riffel<sup>1</sup>, Kai Tobias Block<sup>2</sup>

<sup>1</sup> Institute of Clinical Radiology and Nuclear Medicine, University Medical Center Mannheim, University of Heidelberg, Germany

#### Challenges in abdominal DCE-MRI

Dynamic Contrast-Enhanced Magnetic Resonance Imaging (DCE-MRI) is a relevant part of standard renal MRI examination and indispensable for the evaluation of renal pathologies, such as renal masses [1-4]. It also enables the measurement of renal perfusion as well as glomerular filtration rate (GFR), since gadolinium-based contrast agents are excreted almost only by the kidneys [5-11]. However DCE-MRI exams remain challenging and the failure rates are undesirably high. Especially in sick patients unable to hold their breath adequately, respiratory-motion artifacts compromise achievable image quality. Another challenge is the tradeoff between spatial and temporal resolution of DCE-MRI data for measuring perfusion parameters [12, 13]. On the one hand, a high spatial resolution without motion artifacts is beneficial for a diagnostic image

quality of morphological images. On the other hand, a high temporal resolution (less than 4 s per volume) is desirable for an accurate calculation of perfusion parameters [13].

## Solution approaches


For acceleration of image acquisition, several parallel imaging techniques, such as generalized autocalibrating partially parallel acquisition (GRAPPA) or controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) are available. However, the inverse correlation of the acceleration factor and signal-to-noise ratio (SNR) leads to a loss in SNR. Additionally, several techniques have been introduced to reduce respiratory-motion artifacts and to manage the tradeoff between spatial and temporal resolution such as keyhole [14], k-t broad-use linear acquisition speed-up technique (k-t BLAST) [15] and time-resolved angiography with stochastic trajectories (TWIST) [16].

#### The GRASP approach

The recent introduction of a new approach, golden-angle radial sparse parallel MRI (GRASP)1 [17, 18], could potentially address all of the above mentioned issues. By using non-Cartesian acquisition techniques, specifically radial k-space sampling, motion-induced ghosting artifacts can be reduced. The radial sampling is performed using the golden-angle radial sampling scheme for a flexible spatiotemporal resolution [18].

To improve the temporal resolution and to suppress undersampling artifacts, iterative reconstruction with through-time total variation regularization is applied. A 2016 study showed that multiphase MRI of the kidney is a suitable candidate for the application of GRASP [17], probably because of the high amount of spatiotemporal data correlations, which is a premise for the application of that technique.

<sup>1</sup> WIP, GRASP is currently under development and is not for sale in the US and in other countries. Its future availability cannot be ensured.



<sup>&</sup>lt;sup>2</sup> Center for Advanced Imaging Innovation and Research (CAl<sup>2</sup>R), Department of Radiology, New York University School of Medicine, New York, NY, USA

## Improved workflow and increased efficiency

As data are acquired continuously with a native phase prior to contrast medium (CM) application (Fig. 1), timing and synchronization errors are no longer possible. Acquiring all information using a single continuous scan instead of several individual scans simplifies workflow even further.

# Full flexibility by retrospective choice of temporal resolution

With GRASP, data are acquired continuously during free-breathing. The sampling scheme uses the golden angle as radial increment, which allows that basically any number of radial spokes can be combined into an image providing the possibility to retrospectively decide which temporal resolution should be reconstructed. For morphological images, a poorer temporal resolution (e.g. 55 spokes per image for a temporal resolution of 6.2 seconds) can be reconstructed. For perfusion assessment, images can be reconstructed from the same data set with 13 spokes per image, resulting in a temporal resolution of 1.5 seconds (Fig. 2).

# Role model for the combination of morphological and functional imaging


One major strength of MRI is the possibility of acquiring both morphologic and functional information. Perfusion imaging thus represents a cornerstone of functional imaging and is critical for achieving greater diagnostic specificity in terms of lesion characterization and assessment of organ function. The GRASP technique is an excellent example of the combination of morphological and functional imaging and shows considerable potential for overcoming some of the inherent


limitations of conventional DCE-MRI. In particular, the technique delivers an overall good image quality and high diagnostic confidence of morphological images (Fig. 3) [17]. Due to use of radial k-space sampling, motion-induced ghosting artifacts are minimal although data is acquired during free-breathing.

Additionally, the motion insensitive technique has been shown to be a


valuable tool to achieve valid and reproducible perfusion parameter maps (Fig. 4) [17].


In conclusion, this technique has the potential to become an integral part of renal MRI. The combination of morphological imaging and functional imaging helps to further entrench MRI as a qualitative and quantitative imaging modality.







64-year-old female patient who underwent an examination of the kidneys after contrast administration. Arterial (3A) and renal (3B) phase of the GRASP reconstruction showing good image quality without respiratory motion artifacts.









48-year-old male patient with histologically proven clear cell RCC. T2w HASTE sequence (4A), arterial (4B) and renal (4C) phase of the GRASP reconstruction, estimated renal plasma flow (RPF) parameter map (4D) showing hypervascularity of the lesion.

#### References

- Scialpi M, Brunese L, Piscioli I, et al. Dynamic contrast-enhanced MR imaging for differentiation of renal cell carcinoma subtypes: myth or reality? Radiology. 2009;252:929; author reply 30-1. doi:10.1148/radiol.2523090695.
- 2. Sun MR, Ngo L, Genega EM, et al. Renal cell carcinoma: dynamic contrastenhanced MR imaging for differentiation of tumor subtypes--correlation with pathologic findings. Radiology. 2009;250:793-802. doi:10.1148/ radiol.2503080995.
- 3. Ho VB, Allen SF, Hood MN, et al. Renal masses: quantitative assessment of enhancement with dynamic MR imaging. Radiology. 2002;224:695-700. doi:10.1148/radiol.2243011048.
- 4. Notohamiprodjo M, Sourbron S, Staehler M, et al. Measuring perfusion and permeability in renal cell carcinoma with dynamic contrast-enhanced MRI:

- a pilot study. J Magn Reson Imaging. 2010;31:490-501. doi:10.1002/ jmri.22028.
- 5. Kang SK, Huang WC, Wong S, et al. Dynamic contrast-enhanced magnetic resonance imaging measurement of renal function in patients undergoing partial nephrectomy: preliminary experience. Invest Radiol. 2013;48:687-92. doi:10.1097/ RLI.0b013e3182909e7b.
- 6. Wentland AL, Sadowski EA, Djamali A, et al. Quantitative MR measures of intrarenal perfusion in the assessment of transplanted kidneys: initial experience. Acad Radiol. 2009;16:1077-85. doi:10.1016/j.acra.2009.03.020.
- 7. Ludemann L, Nafz B, Elsner F, et al. Absolute quantification of regional renal blood flow in swine by dynamic contrast-enhanced magnetic resonance imaging using a blood pool contrast agent. Invest Radiol. 2009;44:125-34. doi:10.1097/RLI.0b013e318193598c.

- 8. Michaely HJ, Herrmann KA, Nael K, et al. Functional renal imaging: nonvascular renal disease. Abdom Imaging. 2007;32:1-16. doi:10.1007/ s00261-005-8004-0.
- 9. Grenier N, Hauger O, Cimpean A, et al. Update of renal imaging. Semin Nucl Med. 2006;36:3-15. doi:10.1053/j. semnuclmed.2005.08.001.
- 10. Vallee JP. Lazevras F. Khan HG. et al. Absolute renal blood flow quantification by dynamic MRI and Gd-DTPA. Eur Radiol. 2000;10:1245-52. doi:10.1007/ s003300000434.
- 11. Sourbron SP, Michaely HJ, Reiser MF, et al. MRI-measurement of perfusion and glomerular filtration in the human kidney with a separable compartment model. Invest Radiol. 2008;43:40-8. doi:10.1097/RLI.0b013e31815597c5.
- 12 Kuhl CK Schild HH Morakkahati N Dynamic bilateral contrast-enhanced MR imaging of the breast: trade-off between spatial and temporal resolution. Radiology. 2005;236:789-800. doi:10.1148/radiol.2363040811.
- 13. Michaely HJ, Sourbron SP, Buettner C, et al. Temporal constraints in renal perfusion imaging with a 2-compartment model. Invest Radiol. 2008;43:120-8. doi:10.1097/RLI.0b013e3181583b0c.
- 14. van Vaals JJ, Brummer ME, Dixon WT, et al. "Keyhole" method for accelerating imaging of contrast agent uptake. J Magn Reson Imaging. 1993;3:671-5.
- 15. Tsao J, Boesiger P, Pruessmann KP. k-t BLAST and k-t SENSE: dynamic MRI with high frame rate exploiting spatiotemporal correlations. Magn Reson Med. 2003;50:1031-42. doi:10.1002/ mrm.10611.
- 16. Herrmann KH, Baltzer PA, Dietzel M, et al. Resolving arterial phase and temporal enhancement characteristics in DCE MRM at high spatial resolution with TWIST acquisition. J Magn Reson Imaging. 2011;34:973-82. doi:10.1002/jmri.22689.
- 17. Riffel P, Zoellner FG, Budjan J, et al. "One-Stop Shop": Free-Breathing Dynamic Contrast-Enhanced Magnetic Resonance Imaging of the Kidney Using Iterative Reconstruction and Continuous Golden-Angle Radial Sampling. Invest Radiol. 2016. doi:10.1097/ RLI.0000000000000299.
- 18. Chandarana H, Feng L, Block TK, et al. Free-breathing contrast-enhanced multiphase MRI of the liver using a combination of compressed sensing, parallel imaging, and golden-angle radial sampling. Invest Radiol. 2013;48:10-6. doi:10.1097/ RLI.0b013e318271869c.

#### Contact

Philipp Riffel, M.D. Institute of Clinical Radiology and Nuclear Medicine University Medical Center Mannheim Medical Faculty Mannheim - University Heidelberg Theodor-Kutzer-Ufer 1-3 68167 Mannheim Germany

Phone: +49 621 383 2067 Fax: +49 621 383 3817 Philipp.Riffel@umm.de

