MAGNETOM Flash


The Magazine of MRI 翻訳版

Issue Number 3/2015 | Simultaneous Multi-Slice Supplement

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
SIEMENS

The MAGNETOM World

Your portal to talks, articles and case studies

siemens.com/magnetom-world

Siemens' global MRI community offers peer-to-peer support and information. Radiologists, physicists, cardiologists, and technologists, have all contributed with publications, presentations, training documents, case studies, and more – all freely available to you via this unique network.

Don't miss the >300 lectures and presentations by international and renowned experts on all aspects of MRI that will allow you to be exposed to new ideas and alternative approaches.

The centerpiece of the MAGNETOM World Internet platform consists of MAGNETOM users' results. Here you will find case reports, articles and application tips allowing you to optimize your daily work.

Put the advantages of the MAGNETOM World to work for you!

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18 膝関節MRIにおけるSimultaneous Multi-Slice併用TSEの評価

Xiaona Li; Zhigang Peng; Panli Zuo; Dingxin Wang; Jianling Cui

Cover image courtesy of Choukri Mekkaoui

(Athinoula A. Martinos Center for Biomedical Imaging 149 13th Street Charlzestown, MA 02129, USA)

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¹ Institute of Radiology, University of Würzburg, Würzburg, Germany
² Comprehensive Heart Failure Center (CHFC), Würzburg, Germany
³ Research Center Magnetic Resonance Bavaria (MRB), Würzburg, Germany
⁴ Siemens Healthcare, Erlangen, Germany

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References

- Atkinson DJ, Burnstein D, Edelman RR. First-pass cardiac perfusion: evaluation with ultrafast MR imaging. *Radiology* 1990; 174:757–762.
- Wilke N, Jerosch-Herold M, Wang Y, Yimei H, Christensen BV, Stillman E, Ugurbil K, McDonald K, Wilson RF. Myocardial Perfusion Reserve: Assessment with Multisection, Quantitative, First-Pass MR Imaging. *Radiology* 1997; 204:373–384.
- Rieber J, Huber A, Erhard I, Mueller S, Schweyer M, Koenig A, Schiele TM, Theisen K, Siebert U, Schoenberg SO, Reiser M, Klauss V. Cardiac magnetic resonance perfusion imaging for the functional assessment of coronary artery disease: a comparison with coronary angiography and fractional flow reserve. *Eur Heart J* 2006; 27:1465–1471.
- Schwitzer J, Nanz D, Kneifel S, Bertschinger K, Büchi M, Knüsel PR, Marincek B, Lüscher TF, Schulthess GK. Assessment of Myocardial Perfusion in Coronary Artery Disease by Magnetic Resonance. *Circulation* 2001; 103:2230–2235.
- Ritter CO, del Savio K, Brackertz A, Beer M, Hahn D, Köstler H. High-resolution MRI for the quantitative evaluation of subendocardial and subepicardial perfusion under pharmacological stress and at rest. *RoFo* 2007; 179:945–952.
- Strach K, Meyer C, Thomas D, Naehle CP, Schmitz C, Litt H, Bernstein A, Cheng B, Schild H, Sommer T. High-resolution myocardial perfusion imaging at 3 T: comparison to 1.5 T in healthy volunteers. *Eur Radiol* 2007; 17:1829–1835.
- Jung B, Honal M, Hennig J, Markl M. k-t-space accelerated myocardial perfusion. *J Magn Reson Imag* 2008; 28:1080–1085.
- Köstler H, Sandstede JJW, Lipke C, Landschütz W, Beer M, Hahn D. Auto-SENSE perfusion imaging of the whole human heart. *J Magn Reson Imag* 2003; 18:702–708.
- Pruessmann KP, Weiger M, Scheidegger MB, Boesiger P. SENSE: sensitivity encoding for fast MRI. *Magn Reson Med* 1999; 42:952–962.
- Kellman P, Derbyshire JA, Agyeman KO, McVeigh ER, Arai AE. Extended coverage first-pass perfusion imaging using slice-interleaved TSENSE. *Magn Reson Med* 2004; 51:200–204.
- Stäb D, Ritter CO, Breuer FA, Weng AM, Hahn D, Köstler H. CAIPIRINHA accelerated SSFP imaging. *Magn Reson Med* 2011; 65:157–164.
- Breuer FA, Blaimer M, Heidemann RM, Mueller MF, Griswold MA, Jakob PM. Controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) for multi-slice imaging. *Magn Reson Med* 2005; 53:684–691.
- Breuer F, Blaimer M, Griswold M, Jakob P. Controlled Aliasing in Parallel Imaging Results in Higher Acceleration (CAIPIRINHA). *Magnetom Flash* 2012; 49:135–142.
- Griswold MA, Jakob PM, Heidemann RM, Nittka M, Jellus V, Wang J, Kiefer B, Haase A. Generalized autocalibrating partially parallel acquisitions (GRAPPA). *Magn Reson Med* 2002; 47:1202–1210.
- Breuer FA, Kannengiesser SAR, Blaimer M, Seiberlich N, Jakob PM, Griswold MA. General formulation for quantitative G-factor calculation in GRAPPA reconstructions. *Magn Reson Med* 2009; 62:739–746.
- Adluru G, Awate SP, Tasdizen T, Whitaker RT, Dibella EVR. Temporally constrained reconstruction of dynamic cardiac perfusion MRI. *Magn Reson Med* 2007; 57: 1027–1036.
- Otazo R, Kim D, Axel L, Sodickson DK. Combination of compressed sensing and parallel imaging for highly accelerated first-pass cardiac perfusion MRI. *Magn Reson Med* 2010; 64:767–776.
- Ge L, Kino A, Griswold M, Mistretta C, Carr JC, Li D. Myocardial perfusion MRI with sliding-window conjugate-gradient HYPR. *Magn Reson Med* 2009; 62:835–839.
- Plein S, Kozerke S, Suerder D, Luescher TF, Greenwood JP, Boesiger P, Schwitzer J. High spatial resolution myocardial perfusion cardiac magnetic resonance for the detection of coronary artery disease. *Eur Heart J* 2008; 29:2148–2155.

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¹ Athinoula A. Martinos Center for Biomedical Imaging, Department of Radiology, Massachusetts General Hospital, Harvard Medical School, Boston, MA, USA

² Department of Computer Science, Institute of Mathematics and Statistics, University of São Paulo, São Paulo, Brazil

³ Siemens Healthcare, Charlestown, MA, USA

⁴ Cardiovascular Research Center, Cardiology Division, Massachusetts General Hospital, Harvard Medical School, Boston, MA, USA

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References

- Hennig, J., Nauerth, A. & Friedburg, H., 1986. RARE imaging: a fast imaging method for clinical MR. *Magnetic Resonance in Medicine* : 3(6), pp.823-833.
- Prenger, E.C. et al., 1994. Comparison of T2 weighted spin echo and fast spin echo techniques in the evaluation of myelination. : *JMRI*, 4(2), pp.179-184.
- Patola, W.B. et al., 2001. A comparison of conventional spin-echo and fast spin-echo in the detection of multiple sclerosis. : *JMRI*, 13(5), pp.657-667.
- Pruessmann, K.P. et al., 1999. SENSE: sensitivity encoding for fast MRI. *Magnetic Resonance in Medicine* : 42(5), pp.952-962.
- Pruessmann, K.P. et al., 2001. Advances in sensitivity encoding with arbitrary k-space trajectories. *Magn Reson in Med* : 46(4), pp.638-651.
- Griswold MA, Jakob PM, Heidemann RM, et al. Generalized autocalibrating partially parallel acquisitions (GRAPPA). *Magn Reson Med*. 2002;47:1202-1210
- Breuer, F.A. et al., 2005. Controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) for multi-slice imaging. *Magnetic Resonance in Medicine* : 53(3), pp.684-691.
- Eichner, C. et al., 2013. Slice accelerated gradient-echo spin-echo dynamic susceptibility contrast imaging with blipped CAIPI for increased slice coverage. *Magnetic Resonance in Medicine* 72(3), pp.770-778.
- Gagoski, B.A. et al., 2015. RARE/turbo spin echo imaging with simultaneous multislice Wave-CAIPI. *Magnetic resonance in medicine* : 73(3), pp.929-938.
- Bilgic, B. et al., 2014. Wave-CAIPI for highly accelerated 3D imaging. *Magnetic Resonance in Medicine* : 73(6), pp.2152-2162.
- Wang, D 2014 Multiband Slice Accelerated TSE: Application in Brain Imaging. *Proc Int Soc Mag Reson Med* 22 , 4317

Contact

Michael Kean, FSMRT
Chief MRI Technologist
The Royal Children's Hospital
Flemington Road
Parkville Victoria 3052
Australia
michael.kean@rch.org.au

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Xiaona Li¹; Zhigang Peng¹; Panli Zuo²; Dingxin Wang³; Jianling Cui¹

¹ Department of Radiology, Hebei 3rd Hospital, Province Biomechanical Key Laboratory of Orthopedics, Shijiazhuang, China

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References

- 1 Xu J, Moeller S, Auerbach EJ et al (2013) Evaluation of slice accelerations using multiband echo planar imaging at 3 T. *Neuroimage* 83:991-1001.
- 2 Moeller S, Yacoub E, Olman CA et al (2010) Multiband multislice GE-EPI at 7 tesla, with 16-fold acceleration using partial parallel imaging with application to high spatial and temporal whole-brain fMRI. *Magn Reson Med* 63:1144-1153.
- 3 Setsompop K, Gagoski BA, Polimeni JR, Witzel T, Wedeen VJ, Wald LL (2012) Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty. *Magn Reson Med* 67:1210-1224.
- 4 Sotiropoulos SN, Moeller S, Jbabdi S et al (2013) Effects of image reconstruction on fiber orientation mapping from multi-channel diffusion MRI: reducing the noise floor using SENSE. *Magn Reson Med* 70:1682-1689.
- 5 Feinberg DA, Moeller S, Smith SM et al (2010) Multiplexed echo planar imaging for sub-second whole brain FMRI and fast diffusion imaging. *PLoS One* 5:e15710.
- 6 Smith SM, Miller KL, Moeller S et al (2012) Temporally-independent functional modes of spontaneous brain activity. *Proc Natl Acad Sci U S A* 109:3131-3136.
- 7 Breuer, FA, Blaimer M, Heidemann RM et al., (2005) Controlled aliasing in parallel imaging results in higher acceleration (CAIPIRINHA) for multi-slice imaging. *Magnetic Resonance in Medicine* 53(3): 684–691.
- 8 Norris DG, Boyacioglu R, Schulz J, Barth M, Koopmans PJ (2014) Application of PINS radiofrequency pulses to reduce power deposition in RARE/turbo spin echo imaging of the human head. *Magn Reson Med* 71:44-49.
- 9 Wang D, Kollasch P, Li X et al., "Multiband Slice Accelerated TSE: Clinical Applications in Brain imaging", *Proceedings of the International Society for Magnetic Resonance in Medicine* (22) 2014, 4317.
- 10 Wang D, Padua A, Ellermann J et al., "Multiband Slice Accelerated TSE for High Resolution Knee Imaging", *Proceedings of the International Society for Magnetic Resonance in Medicine* (22) 2014, 1216.

Contact

Jianling Cui
 Department of Radiology
 The Third Hospital of Hebei Medical University
 Hebei Province Biomechanical Key Laboratory
 of Orthopedics
 Ziqiang Road 139
 Shijiazhuang, Hebei 050051
 China