



SIEMENS

DICOM Conformance Statement

syngo® MR XA10A

siemens.com/dicom

1 Conformance Statement Overview

Syngo® MR conforms to the DICOM Standard and supports the network services as described in Table 1: Network Services and the media services as described in Table 2 - Media Services.

Table 1: Network Services

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|--|-------------------------------|-----------------------|------|---------------------------|---------|
| Verification | | | | | |
| Verification | 1.2.840.10008.1.1 | Yes | | Yes | |
| SOP Classes managed by syngo® MR | | | | | |
| | | Create | Send | Store | Display |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | No | Yes | Yes | Yes |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | No | Yes | Yes | Yes |
| Digital X-Ray Image Storage – For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | No | Yes | Yes | Yes |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2 | No | Yes | Yes | Yes |
| Digital Mammography X-Ray Image Storage – For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | No | Yes | Yes | Yes |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | No | Yes | Yes | Yes |
| Enhanced CT Image Storage | 1.2.840.10008.5.1.4.1.1.2.1 | No | Yes | Yes | Yes |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | No | Yes | Yes | Yes |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | Yes | Yes | Yes | Yes |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | Yes | Yes | Yes | Yes |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | Yes | Yes | Yes | Yes |
| Enhanced MR Color Image Storage | 1.2.840.10008.5.1.4.1.1.4.3 | Yes | Yes | Yes | Yes |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | No | Yes | Yes | Yes |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Yes | Yes | Yes | Yes |
| Multi-frame Single Bit Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.1 | No | Yes | Yes | Yes |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2 | No | Yes | Yes | Yes |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3 | Yes | Yes | Yes | Yes |
| Multi-frame True Color Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | Yes | Yes | Yes | Yes |
| 12-lead ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.1 | No | Yes | Yes | No |
| General ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.2 | No | Yes | Yes | No |

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|--|--------------------------------|-----------------------|-----|---------------------------|-----|
| Ambulatory ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.3 | No | Yes | Yes | No |
| Hemodynamic Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.2.1 | No | Yes | Yes | No |
| Cardiac Electrophysiology Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.3.1 | No | Yes | Yes | No |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 | Yes | Yes | Yes | Yes |
| Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.2 | No | Yes | Yes | No |
| Pseudo-Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.3 | No | Yes | Yes | No |
| Blending Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.4 | No | Yes | Yes | No |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | No | Yes | Yes | Yes |
| Enhanced XA Image Storage | 1.2.840.10008.5.1.4.1.1.12.1.1 | No | Yes | Yes | Yes |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | No | Yes | Yes | Yes |
| Enhanced XRF Image Storage | 1.2.840.10008.5.1.4.1.1.12.2.1 | No | Yes | Yes | Yes |
| X-Ray 3D Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.1 | No | Yes | Yes | Yes |
| Breast Tomosynthesis Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.3 | No | Yes | Yes | Yes |
| Raw Data Storage | 1.2.840.10008.5.1.4.1.1.66 | Yes | Yes | Yes | No |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | No | Yes | Yes | No |
| Spatial Fiducials Storage | 1.2.840.10008.5.1.4.1.1.66.2 | No | Yes | Yes | No |
| Deformable Spatial Registration SOP Class | 1.2.840.10008.5.1.4.1.1.66.3 | No | Yes | Yes | No |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | No | Yes | Yes | Yes |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | No | Yes | Yes | No |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | No | Yes | Yes | No |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | No | Yes | Yes | No |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 | Yes | Yes | Yes | No |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | Yes | Yes | Yes | No |
| Procedure Log Storage Storage | 1.2.840.10008.5.1.4.1.1.88.40 | No | Yes | Yes | No |
| Mammography CAD SR Storage | 1.2.840.10008.5.1.4.1.1.88.50 | No | Yes | Yes | No |
| Key Object Selection Document Storage | 1.2.840.10008.5.1.4.1.1.88.59 | No | Yes | Yes | No |
| X-Ray Radiation Dose SR Storage | 1.2.840.10008.5.1.4.1.1.88.67 | No | Yes | Yes | No |
| Encapsulated PDF Storage SOP Class | 1.2.840.10008.5.1.4.1.1.104.1 | No | Yes | Yes | No |
| Positron Emission Tomography Image Storage | 1.2.840.10008.5.1.4.1.1.128 | No | Yes | Yes | Yes |

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|---|-------------------------------|-----------------------|-----|---------------------------|-----|
| RT Image Storage | 1.2.840.10008.5.1.4.1.1.481.1 | No | Yes | Yes | Yes |
| RT Dose Storage | 1.2.840.10008.5.1.4.1.1.481.2 | No | Yes | Yes | No |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | No | Yes | Yes | No |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | No | Yes | Yes | No |
| RT Plan Storage | 1.2.840.10008.5.1.4.1.1.481.5 | No | Yes | Yes | No |
| RT Treatment Summary Record Storage | 1.2.840.10008.5.1.4.1.1.481.7 | No | Yes | Yes | No |
| Transfer (Private SOP Class) | | | | | |
| Syngo Non-Image Storage | | Yes | | Yes | |
| Storage Commitment | | | | | |
| Storage Commitment Push Model SOP Class | 1.2.840.10008.1.20.1 | Yes | | Yes | |
| Worklist Management | | | | | |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Yes | | No | |
| Modality Performed Procedure Step SOP Class | 1.2.840.10008.3.1.2.3.3 | Yes | | No | |
| Query/Retrieve | | | | | |
| Patient Root Q/R Information Model - FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Yes | | Yes | |
| Patient Root Q/R - Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Yes | | Yes | |
| Study Root Q/R - Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | | Yes | |
| Study Root Q/R - Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | | Yes | |
| Patient/Study Only Q/R - Information Model FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | | Yes | |
| Patient/Study Only Q/R - Information Model MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | Yes | | Yes | |
| Print Management | | | | | |
| Basic Grayscale Print Management Meta | | Yes | | No | |
| Basic Color Print Management Meta | | Yes | | No | |
| Print Job | | Yes | | No | |
| Presentation LUT | | Yes | | No | |

Table 2 - Media Services

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| Compact Disk - Recordable | | |
| STD-GEN-CD | Yes | Yes |
| AUG-GEN-CD | Yes | Yes |
| | | |
| DVD | | |
| AUG-GEN-DVD | Yes | Yes |
| AUG- GEN-DVD-J2K | Yes | Yes |

| | | |
|------------------|-----|-----|
| STD-GEN-DVD | Yes | Yes |
| STD-GEN-DVD-J2K | Yes | Yes |
| | | |
| USB | | |
| AUG- GEN-USB-J2K | Yes | Yes |
| STD-GEN-USB-J2K | Yes | Yes |
| | | |

Table 3 - Implementation Identifying Information

| Name | Value |
|-----------------------------|------------------------|
| Application Context Name | 1.2.840.100008.3.1.1.1 |
| Implementation Class UID | 1.3.12.2.1107.5.2 |
| Implementation Version Name | "SYNGO_MR_XA10A" |

2 Table of Contents

| | | |
|------------|---|-----------|
| 1 | Conformance Statement Overview | 2 |
| 2 | Table of Contents | 6 |
| 3 | Introduction..... | 8 |
| 3.1 | Scope and Field | 8 |
| 3.2 | Audience | 8 |
| 3.3 | Remarks | 8 |
| 3.4 | Definitions, Terms and Abbreviations | 8 |
| 3.5 | References | 9 |
| 4 | Networking | 10 |
| 4.1 | Implementation Model..... | 10 |
| 4.1.1 | Application Data Flow | 10 |
| 4.1.2 | Functional Definitions of Application Entities | 12 |
| 4.1.3 | Sequencing of Activities..... | 13 |
| 4.2 | Application Entity Specification | 16 |
| 4.2.1 | Verification AE Specification | 16 |
| 4.2.2 | Storage AE Specification | 18 |
| 4.2.3 | Storage Commitment AE Specification | 23 |
| 4.2.4 | Query/Retrieve AE Specification | 26 |
| 4.2.5 | Modality Worklist AE Specification | 31 |
| 4.2.6 | Modality Performed Procedure Step AE Specification | 37 |
| 4.2.7 | Print AE Specification | 40 |
| 4.2.8 | Print SCU Specification | 40 |
| 4.3 | Network Interfaces | 49 |
| 4.3.1 | Physical Network Interface | 49 |
| 4.3.2 | Additional Protocols | 49 |
| 4.3.3 | IPv4 and IPv6 Support..... | 49 |
| 4.4 | Configuration..... | 49 |
| 4.4.1 | AE Title/Presentation Address Mapping | 49 |
| 4.4.2 | Parameters | 50 |
| 5 | Media Interchange | 51 |
| 5.1 | Implementation Model..... | 51 |
| 5.1.1 | Application Data Flow Diagram | 51 |
| 5.1.2 | Functional definitions of AEs | 52 |
| 5.1.3 | Sequencing of Real-World Activities | 52 |
| 5.1.4 | File Meta Information for Implementation Class and Version..... | 53 |
| 5.2 | AE SPECIFICATIONS..... | 53 |
| 5.2.1 | Media Storage AE – Specification | 53 |
| 5.3 | AUGMENTED AND PRIVATE APPLICATION PROFILES..... | 55 |
| 5.3.1 | Augmented Application Profiles..... | 55 |
| 5.4 | MEDIA CONFIGURATION..... | 55 |

| | | |
|---|--|------------------|
| 6 | <i>Support of Extended Character Sets</i> | <i>56</i> |
| 7 | <i>Attribute confidentiality profiles.....</i> | <i>59</i> |
| 7.1 | De-identification | 59 |
| 8 | <i>Security</i> | <i>65</i> |
| 8.1 | Security Profiles | 65 |
| 8.2 | Association Level Security | 65 |
| 8.3 | Application Level Security | 65 |
| 9 | <i>Annexes.....</i> | <i>66</i> |
| 9.1 | IOD Contents | 66 |
| 9.1.1 | Enhanced MR Image IOD..... | 66 |
| 9.1.2 | Enhanced MR Color Image IOD | 75 |
| 9.1.3 | MR Spectroscopy IOD | 75 |
| 9.1.4 | Raw Data | 83 |
| 9.1.5 | Evidence Documents..... | 83 |
| 9.2 | Data Dictionary of Private Attributes | 83 |
| 9.3 | Grayscale Image Consistency | 84 |
| <i>Annex A: Index of Tables.....</i> | | <i>85</i> |
| <i>Annex B: Table of Figures</i> | | <i>87</i> |

3 Introduction

3.1 Scope and Field

This DICOM Conformance Statement refers to SIEMENS MR Products using software *syngo*® MR XA10A. The following table relates *syngo*® MR XA10A software versions to SIEMENS *syngo*® MR products.

| Software Name | SIEMENS MR Product |
|----------------|--------------------|
| syngo MR XA10A | MAGNETOM Vida |

The syngo® MR product is a “*syngo*®-speaking^a” Imaging Modality or workstation. The syngo® MR product is designed to be integrated into an environment of medical DICOM-based devices.

3.2 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

3.3 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between syngo® MR and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [1]. DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of conformance statements is the first step towards assessing interconnectivity and interoperability between syngo® MR and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens representative for the most recent product information.

3.4 Definitions, Terms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

| | |
|-------|--|
| AE | DICOM Application Entity |
| AET | Application Entity Title |
| ASCII | American Standard Code for Information Interchange |
| CSE | Customer Service Engineer |
| DCS | DICOM Conformance Statement |
| DICOM | Digital Imaging and Communications in Medicine |
| FSC | File Set Creator |

^a *syngo* is a registered trademark of Siemens Healthcare GmbH

| | |
|-------|---|
| FSR | File Set Reader |
| FSU | File Set Updater |
| GSDF | Grayscale Standard Display Function |
| IOD | DICOM Information Object Definition |
| ISO | International Standard Organization |
| MOD | Magneto-optical Disk |
| n. a. | not applicable |
| NEMA | National Electrical Manufacturers Association |
| O | Optional Key Attribute |
| PDU | DICOM Protocol Data Unit |
| R | Required Key Attribute |
| SCU | DICOM Service Class User (DICOM client) |
| SCP | DICOM Service Class Provider (DICOM Server) |
| SOP | DICOM Service-Object Pair |
| SPS | Scheduled Procedure Step |
| SR | Structured Report |
| TFT | Thin Film Transistor (Display) |
| TID | Template ID |
| U | Unique Key Attribute |
| UID | Unique Identifier |
| UTF-8 | Unicode Transformation Format-8 |
| VR | Value Representation |

3.5 References

[NEMA PS3] Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>^b

[IHE] Integrating the Healthcare Enterprise – IHE Radiology Technical Framework – www.ihe.net

• ^b The DICOM Standard is under continuous maintenance, the current official version is available at <http://dicom.nema.org>

4 Networking

4.1 Implementation Model

syngo® MR supports storing DICOM images to remote nodes like workstations or Archiving Systems. Using the Storage Commitment Service it can request safe keeping of previously stored instances from an Archiving system. Additionally the syngo® MR can query remote notes, retrieve and store selected instances from that node. Using the Modality Worklist service the syngo® MR can query a HIS/RIS for scheduled procedures. Performed procedure status and other procedure data can be returned to the HIS/RIS using the Modality Performed Procedure Step (MPPS) Service.

4.1.1 Application Data Flow

The following figure provides a functional overview of the syngo® MR Application Entities (AE). Relationships are shown between user-invoked activities (in the circles at the left of the AEs) and the associated real-world activities provided by DICOM service providers (in the circles at the right of the AEs)




Figure 1: syngo® MR DICOM Data Flow diagram

4.1.2 Functional Definitions of Application Entities

The SCP components of the Application Entities of the syngo® MR operate as background server processes. They exist as soon as the system is powered up and wait for association requests. Upon accepting an association with a negotiated Presentation Context they start to receive and process the requests described in the following sections.

The SCU components of the Application Entity are invoked upon requests from the user interface or indirect by trigger from internal processes.

4.1.2.1 Functional Definition of Verification AE

Syngo® MR supports the Verification service as a SCP and SCU. As a SCU, Verification can be activated from the Admin Portal during system configuration. As a SCP the product processes and responds to incoming

As a SCP of the Verification Service the syngo® MR processes and responds to incoming verification requests.

4.1.2.2 Functional Definition of Storage AE

The syngo® MR Storage SCU is invoked either directly by the user, by an auto-archive trigger or internally by the Query/Retrieve Application Entity that is responsible for processing retrieve requests. The job consists of data describing the composite objects selected for storage and the destination Application Entity Title. An association is negotiated with the destination Application Entity and the image data is transferred using the DIMSE C-STORE -Request. The transfer status is reported to the initiator of the Storage request.

The Storage SCP component of syngo® MR starts to receive the Composite Image Objects and import them into the database after accepting an association with a negotiated Presentation Context. The system can be configured in such a way, that Responses to the Storage Request are sent immediately after reception of the data, after persistent storage on the hard disc or after storage and indexing in the local database.

4.1.2.3 Functional Definition of the Storage Commitment AE

If configured, syngo® MR can serve as a SCU for the DICOM Storage Commitment service. Upon successful completion of a storage job, the system uses the DIMSE N-ACTION Request to request storage commitment from a remote DICOM storage commitment SCP. This can either be the same as the storage destination or a different system depending on the system configuration. Storage Commitment Requests are sent after a configurable delay after storing the objects.

4.1.2.4 Functional Definition of Query/Retrieve AE

Syngo® MR supports DICOM Query/Retrieve as a SCU: The user can initiate a query to a remote node using the C-Find Request. After matching the specified keys, the remote Query /Retrieve SCP uses the C-Find Response to return the results of its search, which will be displayed to the user. Depending on user action the syngo® MR Query/Retrieve DICOM SCU sends a C-MOVE request to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query/Retrieve SCP) to the system's Storage SCP.

Syngo® MR supports the following query models:

- Study Root Query Model.
- Patient Root Query Model
- Patient/Study Only Query Model

Furthermore the SCU services may issue relational queries, if supported by the remote Query/Retrieve SCP node and required by the querying Application.

The syngo® MR DICOM Query/Retrieve SCP accepts C-Find Request, queries the local database based on the provided matching keys and returns the matches using the C-FIND Response. Depending on further request from the remote Query/Retrieve SCU, syngo® MR responds to C-MOVE requests by initiating a C-STORE sub-operation to send image objects to the Storage SCP of the querying system.

4.1.2.5 Functional Definition of Modality Worklist AE

The syngo® MR Modality Worklist SCU issues DICOM Modality Worklist requests using C-FIND requests. The results in the C-FIND response are stored in internal database. The provided Patient and Procedure information is used for patient registration prior to starting an exam.

4.1.2.6 Functional Definition of Modality Performed Procedure Step SCU AE

The syngo® MR MPPS SCU uses the N-CREATE request to inform an Information System that a procedure step is IN PROGRESS.

Syngo® MR MPPS SCU uses the N-SET request to inform the Information System about the finalization of the Procedure Step, using either a status of COMPLETED or DISCONTINUED.

4.1.2.7 Functional Definition of Print AE

The Print SCU of syngo® MR is invoked by the user interface to setup film-sheet layout and whenever an image is ready to be printed on film. The Print SCU will hold and maintain all data needed to compile a complete film-sheet from the data (images, layout, configuration) received. Whenever a film-sheet is ready to print the related data is used to supply the Information to the SOP Classes of the Print Management Service Class. A queue is maintained, in order to immediately store several film-sheets in case of resource problems on printer. The SCU will only supply and require the mandatory SOP Classes of the Print Management Service Class.

4.1.3 Sequencing of Activities

This section describes the sequencing of Real-World Activities performed by syngo® MR Entities using a UML sequence diagram. Real-World Activities are depicted as vertical bars and arrows show the events exchanged between them.

4.1.3.1 System Configuration




Figure 2: Sequence Diagram for Real World Activities - System Configuration

4.1.3.2 Acquisition Workflow




Figure 3: Sequence Diagram for Real World Activities -Acquisition workflow

4.2 Application Entity Specification

This section outlines the specifications for each of the Application Entities that are part of syngo® MR.

4.2.1 Verification AE Specification

4.2.1.1 SOP Classes

The Verification AE of syngo® MR provides standard conformance to the Verification SOP Class listed in "Table 1: Network Services" section "Verification" in the ["Conformance Statement Overview"](#).

4.2.1.2 Association Policy

The syngo® MR Admin Portal attempts to open an association for verification request whenever the "Echo" function is activated.

Table 4: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

4.2.1.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any size is supported.

Table 5: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.1.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to "Table 3 - Implementation Identifying Information" in the ["Conformance Statement Overview"](#).

4.2.1.3 Association Initiation Policy

4.2.1.3.1 Activity – "Send Verification" Request

4.2.1.3.1.1 Description and Sequencing of Activity

Syngo® MR serves as a SCU of the Verification Service Class. A C-ECHO request is initiated by the Admin Portal whenever "Verification" is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than "Success" this will be indicated to the user and the association is closed.

¹ Default, the value is configurable

4.2.1.3.1.2 Proposed Presentation Contexts

Table 6 - Presentation Context Table "Verification" below lists the supported presentation contexts for verification requests.

| Table 6 - Presentation Context Table "Verification" | | | | | |
|---|-------------------|---------------------------|-------------------|------|----------------------|
| Presentation Context Table – "Verification" | | | | | |
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

4.2.1.3.1.3 SOP Specific Conformance – Verification SCU

The application conforms to the definition of the Verification SCU in accordance to the DICOM Standard.

4.2.1.4 Association Acceptance Policy

4.2.1.4.1 Activity – "Receive Verification Request"

4.2.1.4.1.1 Description and Sequencing of Activity

Syngo® MR serves as a SCP of the Verification Service Class. If the Verification SCP accepts an association, it will respond to C-ECHOC requests. If the Called AE Title does not match any pre-configured AE Title shared by SCP, the association will be rejected.

4.2.1.4.1.2 Accepted Presentation Contexts

The syngo® MR DICOM application will accept Presentation Contexts as shown in the following table:

| Table 7 - Presentation Context Table "Verification" | | | | | |
|---|-------------------|---------------------------|-------------------|------|----------------------|
| Presentation Context Table – "Verification" | | | | | |
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |

4.2.1.4.1.3 SOP Specific Conformance – Verification SCP

The ECHO-SCP provides standard conformance to the Verification Service Class.

4.2.2 Storage AE Specification

4.2.2.1 SOP Classes

The Storage AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “SOP Classes Created by syngo® MR” and “SOP Classes Managed by syngo® MR” in the [“Conformance Statement Overview”](#).

4.2.2.2 Association Policy

Table 8: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

Syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

4.2.2.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 9: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.2.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.2.3 Association Initiation Policy

4.2.2.3.1 Activity – “Send Storage Request”

4.2.2.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU of the Storage Service Class. The Storage SCU is triggered by the transfer job queue or by an external retrieve request. An association request is sent to the destination AE. Upon successful negotiation of a Presentation Context, the transfer is started. Objects will be transferred sequentially on the same open association

¹ Default, the value is configurable

Syngo® MR does not provide any automated retry mechanism.

4.2.2.3.1.2 Proposed Presentation Contexts

For all Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with “yes” in the Image Objects Column of the table below are supported.

For all Non-Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with “yes” in the Non-Image Objects Column of the table below are supported.

For a distinction between Image and Non-Image Objects please refer to the DICOM Standard PS3.3 Section A.1.4 “Overview of the Composite IOD Module Content”.

Table 10: Proposed Presentation Contexts for Storage

| UID value | Transfer Syntax | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2 | Implicit Value Representation Little Endian native | yes | yes |
| 1.2.840.10008.1.2.1 | Explicit Value Representation Little Endian native | yes | yes |
| 1.2.840.10008.1.2.2 | Explicit Value Representation Big Endian | yes | yes |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed | yes | no |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed | yes | no |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | yes | no |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed | yes | no |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed | yes | no |
| 1.2.840.10008.1.2.5 | RLE Lossless compressed | yes | no |

Depending on the configuration, the Storage SCU will choose a compressed or uncompressed Transfer Syntax among those accepted by the SCP. The Transfer Syntax chosen is the preferred one among the compressed and uncompressed ones. The preference order is the order of occurrence in the configuration. It is possible to configure for a specific node, which Transfer Syntax shall be used and which one shall be excluded. The configuration can even be extended, based on the combination of SOP Classes and supported Transfer Syntaxes.

An instance will be JPEG lossless (Process 1 and Process 2+4) compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME or RGB or YBR_FULL or YBR_FULL_422
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) is >8
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D

An instance will be JPEG lossy compressed only if it fulfills the following criteria:

- is an image and not already compressed
- photometric interpretation (0028,0004) is MONOCHROME or RGB
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) equal to 12'D or 8'D
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D

An instance will be JPEG 2000 lossless compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric interpretation (0028,0004) not MONOCHROME or RGB or YBR_FULL or YBR_FULL_422
- Bits Allocated (0028,0100) not 16'D or 8'D

An instance will be JPEG 2000 lossy compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric interpretation (0028,0004) is MONOCHROME or RGB
- Bits Stored (0028,0101) equal to 12'D or 8'D

There is no extended negotiation as an SCU.

4.2.2.3.1.3 SOP specific Conformance for SOP classes

Syngo® MR will not add or change private attributes by default, even in case objects are compressed or the image header is updated according to the IHE Patient Information Reconciliation Profile.

The behavior of syngo® MR when encountering status codes in a C-STORE response is summarized in Table 11:

Table 11: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|------------------------------|------------|--|
| Error | Any other DIMSE Error Status | 0xFFFF | Send is continued till the end. Log message is created. |
| Success | Image is successfully stored | 0000 | If configured, Storage Commitment is requested for successfully stored instances |

Table 12 below indicates the behavior if exceptions occur:

Table 12: DICOM Command Communication Failure Behavior

| Exception | Behavior |
|---------------------|--|
| Timeout | Log message is created (Timeout configurable; default 30s) |
| Association Aborted | Send is failed. Log message is created. |

4.2.2.4 Association Acceptance Policy

4.2.2.4.1 Activity – “Receive Storage Request”

4.2.2.4.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCP of the Storage Service Class. The storage SCP accepts incoming C-Store Request from any configured AE Title, receives supported objects transmitted on that association and stores them in the local database.

4.2.2.4.1.2 Accepted Presentation Contexts

For all supported Transfer Objects (see “Table 1: Network Services” section “SOP Classes Created by syngo® MR” and “SOP Classes Managed by syngo® MR” in the [“Conformance Statement Overview”](#).) the appropriate Transfer Syntaxes are supported.

Generally all Presentation Contexts are accepted as long as they contain at least one suitable Transfer Syntax. All other Presentation Contexts are rejected.

There is no Extended Negotiation as an SCP

4.2.2.4.1.3 SOP-specific Conformance Statement for Storage SOP classes

Syngo® MR conforms to the Full Storage Class at Level 2.

In case of a successful C-STORE operation, the image has successfully been written on disk either in Explicit Little Endian format or in the compression format received.

The Storage AE of syngo® MR returns the status “success” when the data is stored to disk and a minimal image header validation has been performed.

The following header attributes must be available and filled:

- SOP Class UID,
- Study Instance UID,
- Series Instance UID and
- SOP Instance UID.

Table 13 below list the status codes that syngo® MR can return:

Table 13: Storage C-STORE Response Status

| Service Status | Further Meaning | Error Code | Reason |
|----------------|-----------------|------------|--------|
|----------------|-----------------|------------|--------|

| | | | |
|---------|-----------------------------------|------|--|
| Success | success | 0000 | Image received correctly (success notification is done after receiving, before indexing and storing) |
| Failure | Out-of-resource | A700 | No resource left in the Short Term Storage |
| Failure | Unable to Process | Cxxx | Error during instance reception |
| Failure | Data set does not match SOP Class | A9xx | The data set is not conform to the SOP Class contained in the resource. |

Restriction: successful operation does not guarantee storage of header data in the database.

4.2.2.4.1.4 Other SOP specific behavior

If an image is received that is already stored in the database - identified by the SOP Instance UID
- the new image will be ignored. The existing instance is not superseded.

4.2.3 Storage Commitment AE Specification

4.2.3.1 SOP Classes

The Storage Commitment AE of syngo® MR provides standard conformance to the SOP Class listed in “Table 1: Network Services” section “Storage Commitment” in the [“Conformance Statement Overview”](#).

4.2.3.2 Association Policy

Table 14: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

Syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

4.2.3.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 15: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.3.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.3.3 Association Initiation Policy

4.2.3.3.1 Activity “Send Initial Storage Commitment”

4.2.3.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU of the Storage Commitment Service Class. After successful transfer of Imaging Objects to a configured Archive, the Storage Commitment SCU initiates an N-Action Request, if Storage Commitment is configured. This request will be sent on a different association than the storage request.

¹ Default, the value is configurable

The Storage Commitment Request will be sent out with a delay, in order to ensure that the remote node properly indexes received instances. The delay time is configurable with a default delay of 10 minutes.

The system may issue one N-ACTION Request for a complete set (bundle) of instances or issue one N-ACTION-Request per instance. This behavior is configurable; the default value is “bundled”.

Syngo® MR will accept the N-Event-Report-Request on the same association if sent immediately after the N-ACTION-Response. However it will not wait for it. The association is closed after three seconds.

4.2.3.3.1.2 Proposed Presentation Contexts

Table 16 below lists the supported presentation contexts for storage commitment.

Table 16: Proposed Presentation Contexts for Storage Commitment

| Presentation Context Table | | | | | |
|-------------------------------|----------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.3.3.1.3 SOP specific Conformance for SOP classes

The behavior of syngo® MR when encountering status codes in an N-ACTION response is summarized in Table 17:

Table 17: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|-------------------|---|
| Error | Any failure that occurs | Any non null Code | Failure reported to user; corresponding object(s) will be marked as “Archived failed” |
| Success | All Instances are available on the remote node | 0000 | Success reported to user; in case failures exist, the corresponding instances will be marked as “Archived failed” |

Table 18: DICOM Command Communication Failure Behavior

| Exception | Behavior |
|---------------------|---|
| Timeout | Failure reported to user (Timeout configurable; default 30s); the request will be retried |
| Association Aborted | Failure reported to user; the request will be retried |

4.2.3.4 Association Acceptance Policy

4.2.3.4.1 Activity “Receive Reply to Initial Storage Commitment”

4.2.3.4.1.1 Description and Sequencing of Activities

Syngo® MR supports the reverse role negotiation of the Storage Commitment Service Class as SCU. If the N-ACTION request was not answered in 3 seconds, it closes the association. It accepts incoming N-EVENT-REPORT Request in a new association.

4.2.3.4.1.2 Accepted Presentation Contexts

The syngo® MR DICOM application supports the presentation contexts listed in the following table for the Storage Commitment Service Class.

Table 19 - Presentation Context Table "Update Flag Information"

| Presentation Context Table – “Update Flag Information” | | | | |
|--|--|---|------|-----------|
| Abstract Syntax | Transfer Syntax | | Role | Ext. Neg. |
| Description | Name List | UID List | | |
| 1.2.840.10008.1.20.1 Storage Commitment Push Model | Explicit VR Little Endian Explicit VR Big Endian Implicit VR Little Endian | 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2 1.2.840.10008.1.2 | SCP | None |

4.2.3.4.1.3 SOP-specific Conformance Statement for Storage SOP classes

The behavior of syngo® MR when encountering status codes in an N-EVENT-REPORT response is summarized in the following table:

Table 20: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------------------------|--------------------|--------------------------------------|
| Error | Storage Commitment Reply ignored. | Any none null Code | Storage Commitment will be repeated. |
| Success | Storage Commitment Reply noticed. | 0000 | Success reported to user. |

4.2.4 Query/Retrieve AE Specification

4.2.4.1 SOP Classes

The Query/Retrieve AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Query/Retrieve” in the [“Conformance Statement Overview”](#).

4.2.4.2 Association Policy

Table 21: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

Syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

4.2.4.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 22: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.4.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.4.3 Association Initiation Policy

4.2.4.3.1 Activity “Querying a Remote Node” for Instances

4.2.4.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU for the following SOP Classes

¹ Default, the value is configurable

- Patient Root Q/R Information Model - FIND SOP Class
- Study Root Q/R Information Model –FIND SOP Class
- Patient/Study only Q/R Information Model – FIND SOP Class.

Using the attributes specified by the user as Query Keys (in accordance with the query model) the Query SCU initiates a C-FIND Request and display the responses to the user.

4.2.4.3.1.2 Proposed Presentation Contexts

Syngo® MR will propose Presentation Contexts as shown in the following table:

Table 23: Proposed Presentation Contexts for Query

| Presentation Context Table | | | | | |
|---|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | Yes |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Study Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | Yes |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Patient/Study Only Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

Table 24: Extended Negotiation as an SCU

| Name | UID | Extended Negotiation |
|--|-----------------------------|---|
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Relational Query will be negotiated if necessary as defined in DICOM PS3.4. |
| Study Root Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Relational Query will be negotiated if necessary as defined in DICOM PS3.4. |

4.2.4.3.1.3 SOP Specific Conformance Statement to Query SOP classes

Syngo® MR checks for the following status codes in the Query SCP's C-FIND-Response:

Table 25: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|--------------------|------------------------------------|
| Failure | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process | Any none null Code | Failure reported to user |
| Pending | All optional keys are supported the same manner as Required Keys. | FF00 | Pending state is indicated to user |
| | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01 | Pending state is indicated to user |
| Success | Query has been performed successfully. | 0000 | Success reported to user |

Table 26: DICOM Command Communication Failure Behavior

| exception | Behavior |
|---------------------|---|
| Timeout | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user |

Syngo® MR supports the following query levels:

- Study
- Series

Matching Keys on Instance Level is not supported by the syngo® MR SCU.

The following table lists the various attributes at Study and Series levels, which can be used for hierarchical queries as well as return values for display. The display capabilities are highly configurable and “yes” indicates that it is possible to configure display of the data:

A “yes” in the UI column will indicate that the attribute may be visualized when browsing the Query results with the Browser. The Browser display is additionally influenced by the related Browser configuration

Table 27: Attributes supported for instance Query - SCU

| Attribute Name | Tag | Type | User input | UI |
|----------------------------|-------------|------|-------------|-----|
| Study Level | | | | |
| Patient's Name | (0010,0010) | O | enter value | yes |
| Patient ID | (0010,0020) | O | enter value | yes |
| Patient's Birth Date | (0010,0030) | O | enter value | yes |
| Patient's Birth Time | (0010,0032) | O | enter value | yes |
| Patient's Sex | (0010,0040) | O | enter value | yes |
| Accession Number | (0008,0050) | O | enter value | yes |
| Study ID | (0020,0010) | O | enter value | yes |
| Study Instance UID | (0020,000D) | U | enter value | yes |
| Study Date | (0008,0020) | O | enter value | yes |
| Study Time | (0008,0030) | O | enter value | yes |
| Referring Physician's Name | (0008,0090) | O | enter value | yes |

| Attribute Name | Tag | Type | User input | UI |
|---|---------------------------|------|-------------|-----|
| Study Description | (0008,1030) | O | enter value | yes |
| Number of Study related Instances | (0020,1208) | O | - | yes |
| Modalities in Study | (0008,0061) | O | enter value | yes |
| Number of Study Related Series | (0020,1206) | O | - | yes |
| Series Level | | | | |
| Modality | (0008,0060) | O | enter value | yes |
| Series Date | (0008,0021) | O | enter value | yes |
| Series Time | (0008,0031) | O | enter value | yes |
| Number of Series related Instances | (0020,1209) | O | - | yes |
| Series Number | (0020,0011) | O | enter value | yes |
| Series Description | (0008,103E) | O | enter value | yes |
| Request Attributes Sequence \ Requested Procedure ID | (0040,0275) \ (0040,1001) | O | enter value | yes |
| Request Attributes Sequence \ Scheduled Procedure Step ID | (0040,0275) \ (0040,0009) | O | enter value | yes |
| Performed Procedure Step Start Date | (0040,0244) | O | enter value | yes |
| Performed Procedure Step Start Time | (0040,0245) | O | enter value | yes |
| Series Instance UID | (0020,000E) | U | - | yes |

4.2.4.3.1 Activity “Retrieve Instances from a remote node”

4.2.4.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU for the following SOP Classes

- Patient Root Q/R Information Model - MOVE SOP Class
- Study Root Q/R Information Model – MOVE SOP Class
- Patient/Study only Q/R Information Model – MOVE SOP Class.

The C-MOVE-Request is used to retrieve the selected imaging objects. The Retrieve AE supports the query model Study Root.

4.2.4.3.1.2 Proposed Presentation Contexts

Syngo® MR will propose Presentation Contexts as shown in the following table:

Table 28: Proposed Presentation Contexts for Retrieve and Activity “MOVE SCU”

| Presentation Context Table | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Study Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.4.3.1.3 SOP Specific Conformance Statement for Move SCU Classes

The presentation context is negotiated at association establishment time. When the C-MOVE-Request is processed, the Move Destination attribute (receiver of images) is ignored. However the Move Destination AE must conform to the DICOM conventions (value representation AE).

The behavior of syngo® MR when encountering status codes in a C-MOVE response is summarized in Table 29

Table 29: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|-------------------|-----------------------------------|
| Error | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process; Move destination unknown | Any non null Code | Failure reported to user |
| Pending | Move Operation continues | FF00 | Operation continues in background |
| Success | Move has been performed successfully. | 0000 | Success reported to user |

Table 30: DICOM Command Communication Failure Behavior

| exception | Behavior |
|---------------------|---|
| Timeout | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user |

4.2.4.4 Association Acceptance Policy

Syngo® MR does not provide SCP functionality.

4.2.5 Modality Worklist AE Specication

4.2.5.1 SOP Classes

The Modality Worklist AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Worklist Management” in the [“Conformance Statement Overview”](#).

4.2.5.2 Association Policy

Table 31: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

Syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

4.2.5.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 32: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.5.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.5.3 Association Initiation Policy

4.2.5.3.1 Activity “Querying a Remote Node” for Modality Worklist

4.2.5.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU of the Modality Worklist service. It performs worklist queries by issuing a C-FIND request at regular intervals. In addition a worklist request can be triggered manually.

¹ Default, the value is configurable

4.2.5.3.1.2 Proposed Presentation Contexts

Syngo® MR will propose Presentation Contexts as shown in the following table:

Table 33: Proposed Presentation Contexts for Worklist

| Presentation Context Table | | | | | |
|----------------------------|------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Modality Worklist-FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.5.3.1.3 SOP Specific Conformance for SOP Classes

Search Key Attributes of the Worklist C-FIND

Syngo® MR Modality Worklist SCU supports “broad worklist queries” with all required search keys. The following tables describe the “broad query” search keys that the SCU supports. The list is configurable in ‘DICOM Modality Worklist Query’.

Table 34: Broad Query search keys

| Attribute Name | Tag | Matching Key Type | Query Value |
|--|-------------|-------------------|-------------------------------|
| Scheduled Procedure Step | | | |
| Scheduled Procedure Step Sequence | (0040,0100) | R | |
| >Modality | (0008,0060) | R | “*” or <configured Modality> |
| >Scheduled Station AE Title | (0040,0001) | R | <own AET> or “*” ^c |
| >Scheduled Procedure Step Start Date | (0040,0002) | R | Range from UI ^d |
| >Scheduled Procedure Step Description | (0040,0007) | O | |
| >Scheduled Station Name | (0040,0010) | O | |
| >Scheduled Procedure Step Location | (0040,0011) | O | |
| >Scheduled Procedure Step Status | (0040,0020) | O | |
| >Scheduled Performing Physician's Name | (0040,0006) | O | |
| >Scheduled Protocol Code Sequence | (0040,0008) | O | |

^c This depends on user configuration (Administration Portal->Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist) if the “own AET” is provided or not.

^d A time window can be configured by defining how many days to look into the past and into the future (Administration Portal-> Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist)

| Attribute Name | Tag | Matching Key Type | Query Value |
|-----------------------------------|--------------|-------------------|-------------|
| >>Code Value | (0008,0100) | O | |
| Requested Procedure Description | (0032,1060) | O | |
| Requested Procedure Priority | (0040,1003) | O | |
| Patient Transport Arrangements | (0040,1004) | O | |
| Requested Procedure Comments | (0040,1400) | O | |
| Requested Procedure Code Sequence | (0032,1064) | O | |
| >Code Value | (0008,0100) | O | |
| Requesting Physician | (0032,1032) | O | |
| Referring Physicians Name | (0008,0090) | O | |
| Current Patient Location | (0038,0300) | O | |
| Pregnancy Status | (0010, 21C0) | O | |
| Medical Alerts | (0010,2000) | O | |
| Allergies | (0010,2110) | O | |

Return Key Attributes of the Modality Worklist C-FIND

The syngo® MR Modality Worklist SCU supports worklist queries with return key attributes of all types. The following tables describe the return keys that the SCU supports.

An “x” in the **UI** column indicates that the attribute may be visualized when browsing the Worklist results with the Browser. The Browser display is additionally influenced by the related Browser configuration.

Table 35: Modality Worklist C-Find Return keys

| Attribute Name | Tag | Return Key Type | UI | Notes |
|--|-------------|-----------------|----|---|
| SOP Common | | | | |
| Specific Character Set | (0008,0005) | 1C | - | |
| Scheduled Procedure Step | | | | |
| Scheduled Procedure Step Sequence | (0040,0100) | 1 | | |
| >Modality | (0008,0060) | 1 | x | |
| >Scheduled Station AE Title | (0040,0001) | 1 | | “Scheduled Station AE Title” is taken as default for “Performed Station AE Title” |
| >Scheduled Procedure Step Start Date | (0040,0002) | 1 | - | |
| >Scheduled Procedure Step Start Time | (0040,0003) | 1 | - | |
| >Scheduled Procedure Step End Date | (0040,0004) | 3 | - | |
| >Scheduled Procedure Step End Time | (0040,0005) | 3 | - | |
| >Scheduled Performing Physician's Name | (0040,0006) | 1 | x | “Scheduled Performing Physician's Name” is taken as default for “Performing Physician's Name” |
| >Scheduled Procedure Step Description | (0040,0007) | 1C | x | “Scheduled Procedure Step Description” is taken as |

| Attribute Name | Tag | Return Key Type | UI | Notes |
|---|-------------|-----------------|----|--|
| | | | | default for "Performed Procedure Step Description" |
| >Scheduled Protocol Code Sequence ** | (0040,0008) | 1C | - | Uses universal sequence match "Scheduled Protocol Code Sequence" is taken as default for "Performed Protocol Code Sequence" |
| >>Code Value | (0008,0100) | 1C | - | |
| >>Coding Scheme Designator | (0008,0102) | 1C | - | |
| >>Coding Scheme Version | (0008,0103) | 3 | - | |
| >>Code Meaning | (0008,0104) | 3 | - | |
| >Scheduled Procedure Step ID | (0040,0009) | 1 | x | "Scheduled Procedure Step ID" is taken as default for "Performed Procedure Step ID" |
| >Scheduled Station Name | (0040,0010) | 2 | x | |
| >Scheduled Procedure Step Location | (0040,0011) | 2 | - | "Scheduled Procedure Step Location" is taken as default for "Performed Location" |
| >Scheduled Procedure Step Status | (0040,0020) | 3 | - | |
| >Comments on the Scheduled Procedure Step | (0040,0400) | 3 | - | |
| Requested Procedure | | | | |
| Study Date | (0008,0020) | 3 | x | |
| Study Time | (0008,0030) | 3 | x | |
| Referenced Study Sequence ** | (0008,1110) | 2 | - | Uses universal sequence match |
| >Referenced SOP Class UID | (0008,1150) | 1C | - | |
| >Referenced SOP Instance UID | (0008,1155) | 1C | - | |
| Study Instance UID | (0020,000D) | 1 | - | |
| Requested Procedure Description | (0032,1060) | 1C | x | |
| Requested Procedure Code Sequence ** | (0032,1064) | 1C | - | Uses universal sequence match "Requested Procedure Code Sequence" is taken as default for "Procedure Code Sequence" |
| >Code Value | (0008,0100) | 1C | - | |
| >Coding Scheme Designator | (0008,0102) | 1C | - | |
| >Coding Scheme Version | (0008,0103) | 3 | - | |
| >Code Meaning | (0008,0104) | 3 | - | |
| Requested Procedure ID | (0040,1001) | 1 | x | "Requested Procedure ID" is taken as default for "Study ID" |
| Reason for the Requested Procedure | (0040,1002) | 3 | - | |
| Requested Procedure Priority | (0040,1003) | 2 | x | |
| Patient Transport Arrangements | (0040,1004) | 2 | - | |
| Confidentiality Code | (0040,1008) | 3 | - | |
| Reporting Priority | (0040,1009) | 3 | x | |
| Names of intended Recipients of Results | (0040,1010) | 3 | - | |
| Requested Procedure Comments | (0040,1400) | 3 | - | |
| Imaging Service Request | | | | |
| Accession Number | (0008,0050) | 2 | x | |

| Attribute Name | Tag | Return Key Type | UI | Notes |
|---|-------------|-----------------|----|--|
| Referring Physician's Name | (0008,0090) | 2 | x | |
| Requesting Physician | (0032,1032) | 2 | x | |
| Requesting Service | (0032,1033) | 3 | - | |
| Issuing Date of Imaging Service Request | (0040,2004) | 3 | - | |
| Issuing Time of Imaging Service Request | (0040,2005) | 3 | - | |
| Placer Order Number / Imaging Service Request * | (0040,2016) | 3 | - | Old tag (0040,2006) is retired and not used. |
| Filler Order Number / Imaging Service Request * | (0040,2017) | 3 | - | Old tag (0040,2007) is retired and not used. |
| Order entered by ... | (0040,2008) | 3 | - | |
| Order Enterer's location | (0040,2009) | 3 | - | |
| Order Callback Phone Number | (0040,2010) | 3 | - | |
| Imaging Service Request Comments | (0040,2400) | 3 | - | |
| Visit Identification | | | | |
| Admission ID | (0038,0010) | 2 | x | |
| Issuer of Admission ID | (0038,0011) | 3 | - | |
| Visit Status | | | | |
| Current Patient Location | (0038,0300) | 2 | x | |
| Visit Admission | | | | |
| Admitting Diagnosis Description | (0008,1080) | 3 | x | |
| Admitting Date | (0038,0020) | 3 | - | |
| Patient Identification | | | | |
| Patient's Name | (0010,0010) | 1 | x | |
| Patient ID | (0010,0020) | 1 | x | |
| Issuer of Patient ID | (0010,0021) | 3 | x | |
| Other Patient IDs | (0010,1000) | 3 | x | |
| Other Patient Names | (0010,1001) | 3 | x | |
| Patient's Birth Name | (0010,1005) | 3 | - | |
| Patient Demographic | | | | |
| Patient's Birth Date | (0010,0030) | 2 | x | |
| Patient's Birth Time | (0010,0032) | 3 | x | |
| Patient's Sex | (0010,0040) | 2 | x | |
| Patient's Insurance Plan Code Sequence ** | (0010,0050) | 3 | - | Uses universal sequence match |
| >Code Value | (0008,0100) | 1C | - | |
| >Coding Scheme Designator | (0008,0102) | 1C | - | |
| >Coding Scheme Version | (0008,0103) | 3 | - | |
| >Code Meaning | (0008,0104) | 3 | - | |
| Patient's Age | (0010,1010) | 3 | - | |
| Patient's Size | (0010,1020) | 3 | x | |
| Patient's Weight | (0010,1030) | 2 | x | |
| Patient's Address | (0010,1040) | 3 | x | |
| Military Rank | (0010,1080) | 3 | x | |
| Branch of Service | (0010,1081) | 3 | - | |
| Ethnic Group | (0010,2160) | 3 | x | |
| Patient Comments | (0010,4000) | 3 | x | |
| Patient Medical | | | | |
| Medical Alerts | (0010,2000) | 2 | x | |
| Allergies | (0010,2110) | 2 | x | |
| Pregnancy Status | (0010,21C0) | 2 | x | |
| Smoking Status | (0010,21A0) | 3 | | |
| Last Menstrual Date | (0010,21D0) | 3 | | |
| Additional Patient History | (0010,21B0) | 3 | | |

| Attribute Name | Tag | Return Key Type | UI | Notes |
|----------------|-------------|-----------------|----|-------|
| Special Needs | (0038,0050) | 2 | x | |

The behavior of syngo® MR when encountering status codes in a C-FIND response is summarized in Table 36:

Table 36: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|--------------------|------------------------------------|
| Error | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process | Any none null Code | Failure reported to user |
| Pending | All optional keys are supported the same manner as Required Keys. | FF00 | Pending state is indicated to user |
| | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01 | Pending state is indicated to user |
| Success | Query has been performed successfully. | 0000 | Success reported to user |

Table 37: DICOM Command Communication Failure Behavior

| exception | Behavior |
|---------------------|---|
| Timeout | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user |

4.2.5.4 Association Acceptance Policy

Syngo® MR does not provide the functionality of a SCP of the Modality Worklist – Find SOP Class.

4.2.6 Modality Performed Procedure Step AE Specification

4.2.6.1 SOP Classes

The Modality Performed Procedure Step AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Worklist Management” in the [“Conformance Statement Overview”](#).

4.2.6.2 Association Policy

Table 38: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

Syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

4.2.6.2.1 Asynchronous Nature

Syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 39: Asynchronous Nature as an Association Initiator

| | |
|--|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|--|----|

4.2.6.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.6.3 Association Initiation Policy

4.2.6.3.1 Activity “Create Modality Performed Procedure Step”

4.2.6.3.1.1 Description and Sequencing of Activities

Syngo® MR serves as a SCU of the Modality Performed Procedure Step SOP Class. It sends N-CREATE request to inform the Information System that a Procedure Step has been started.

4.2.6.3.1.2 Accepted Presentation Contexts

Syngo® MR proposes Presentation Contexts as shown in the following table:

¹ Default, the value is configurable

Table 40: Acceptable Presentation Contexts Activity “Create MPPS”

| Presentation Context Table | | | | | |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.6.3.1.3 SOP specific Conformance for MPPS SOP class

The behavior of syngo® MR when encountering status codes in an N-CREATE-RSP response is summarized in Table 41:

Table 41: MPPS N-CREATE Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|--------------------|----------------------|
| Error | MPPS creation request could not be processed. | Any none null Code | MPPS is not created. |
| Success | MPPS creation request processed successfully. | 0000 | MPPS is created. |

4.2.6.3.2 Activity “Update Modality Performed Procedure Step”

4.2.6.3.2.1 Description and Sequencing of Activities

When the procedure step has been finished, syngo® MR sends N-SET request to inform the Information System about the finalization of the procedure step (completed or discontinued).

4.2.6.3.2.2 Proposed Presentation Contexts

Syngo® MR proposes Presentation Contexts as shown in the following table:

Table 42: Acceptable Presentation Contexts Activity “Update MPPS”

| Presentation Context Table | | | | | |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.6.3.2.3 SOP specific Conformance for MPPS SOP class

The behavior of syngo® MR when encountering status codes in an N-SET-RSP response is summarized in Table 43:

Table 43: MPPS N-SET Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|--------------------|----------------------|
| Error | MPPS update request could not be processed. | Any none null Code | MPPS is not updated. |
| Success | MPPS update request could processed successfully. | 0000 | MPPS is updated. |

4.2.6.4 Association Acceptance Policy

Syngo® MR does not provide the functionality of a SCP of the Modality Performed Procedure Step SOP Class.

4.2.7 Print AE Specification

4.2.7.1 SOP Classes

4.2.7.2 Association Policy

4.2.7.3 Association Initiation Policy

4.2.7.4 Association Acceptance Policy

4.2.8 Print SCU Specification

4.2.8.1 SOP Classes

For SOP Classes supported, please refer to “Table 1: Network Services” section “Print Management” in the [“Conformance Statement Overview”](#).

4.2.8.2 Association Policies

4.2.8.2.1 General

Whenever a film-sheet is completely set up and printed by command or automated rule, the job is prepared for processing. As soon as the queue is ready to process the job, it is activated and worked according to the processing data. The Print application will initiate an association to the print destination and process the printing.

The default PDU size used will be used.

4.2.8.2.2 Number of Associations

The syngo® MR DICOM application initiates one association at a time for each different print device configured.

4.2.8.2.3 Asynchronous Nature

The syngo® MR DICOM print application does not support asynchronous communication (multiple outstanding transactions over a single association).

4.2.8.2.4 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.8.3 Association Initiation Policy

Triggered by the Print job queue the Print Management SCU establishes an association by using the DICOM association services. An N-GET request determines the printer status prior to printing. If the printer status is “normal”, the print job is started.

After the last film is printed from queue, the Print application will leave open the association for another 60 seconds. If a new film job is ready for printing within this time-limit, the job will be immediately processed over the still open association. If there is no new job, the association is closed.

During the “idle-time” (no open association to printer) the Print application will issue a cyclic camera status request (using N-GET of the Printer SOP Class) every 5 minutes.

4.2.8.3.1 Activity - Print Film

4.2.8.3.1.1 Description and Sequencing of Activity

The film sheet is internally processed, converted to a Standard/1-1 page and then the page image is sent. Status is controlled by awaiting any N-EVENT message through the transfer until the last image or film-sheet is sent.

If the response from the remote application contains a status other than Success or Warning the association is aborted.

4.2.8.3.1.2 Proposed Presentation Context

The syngo® MR DICOM application will propose Presentation Contexts as shown in the following table:

Table 44 - Presentation Context Table "Print Film"

| Presentation Context Table | | | | | |
|---|------------------------|--|--|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | Implicit VR Little Endian Explicit VR Little Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 | SCU | None |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 | Implicit VR Little Endian Explicit VR Little Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 | SCU | None |
| Print Job SOP Class | 1.2.840.10008.5.1.1.14 | Implicit VR Little Endian Explicit VR Little Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 | SCU | None |
| Presentation LUT SOP Class | 1.2.840.10008.5.1.1.23 | Implicit VR Little Endian Explicit VR Little Endian | 1.2.840.10008.1.2 1.2.840.10008.1.2.1 | SCU | None |

4.2.8.3.1.3 SOP Specific Conformance

The syngo® MR DICOM print management SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class and to Basic Color Print Management Meta SOP Class.

The application uses a setting platform to define the properties of the connected DICOM SCP, e.g.:

- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP
- LUT type to be attached.

The printing is suspended in the case of a failure return status of the SCP or when the user cancels the job.

Basic Film Session SOP Class

The Basic Film Session information object definition describes all the user-defined parameters, which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes that are printed on one hardcopy printer.

The syngo® MR DICOM print management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the following attributes:

Table 45 - Basic Film Session N-CREATE attributes

| Attribute Name | Tag | Usage SCU | Supported Values |
|------------------|-------------|-----------|----------------------------------|
| Number of Copies | (2000,0010) | U | "1" |
| Medium Type | (2000,0030) | U | BLUE FILM CLEAR FILM PAPER |

U = User Option

The number of copies sent to the DICOM Printer is always 1, the job is sent n times for n copies.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session – see below:

| Attribute Name | Tag | Source of Information |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000) →(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Basic Film Session |

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

The Basic Film Session SOP Class interprets the following status codes (from N-CREATE-RSP, N-DELETE-RSP messages):

Table 46 - Basic Film Session Status Codes

| Service Status | Meaning | Error Codes |
|----------------|--|-------------|
| Failure | Film session SOP instances hierarchy does not contain film box SOP instances | C600 |
| | Unable to create print job, print queue is full | C601 |
| | Image size is larger than images box size | C603 |
| Warning | Memory allocation not supported | B600 |
| | Film session printing is not supported | B601 |
| | Film box does not contain image box (empty page) | B602 |
| Success | Film belonging to the film session are accepted for printing | 0000 |

Basic Film Box SOP Class

The Basic Film Box information object definition describes all the user-defined parameter of one film of the film session. The Basic Film Box information description defines the presentation parameters, which are common for all images on a given sheet of film.

The Basic Film Box refers to one or more Image Boxes.

Supported Service Elements as SCU are:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP Class N-CREATE-RQ message uses the following attributes (the actual values for each attribute depend on DICOM printer configuration within the syngo® MR DICOM print management SCU):

Table 47 - Basic Film Box N-CREATE attributes

| Attribute Name | Tag | Usage SCU | Supported Values |
|--------------------------------------|-------------|-----------|--|
| Image Display Format | (2010,0010) | M | STANDARD\1,1 |
| Referenced Film Session Sequence | (2010,0500) | M | n. a. |
| > Referenced SOP Class UID | (0008,1150) | M | 1.2.840.10008.5.1.1.1 |
| > Referenced SOP Instance UID | (0008,1155) | M | |
| Film Orientation | (2010,0040) | M | PORTRAIT |
| Film Size ID | (2010,0050) | M | 8INX10IN, 10INX12IN, 10INX14IN, 11INX14IN,, 14INX14IN, 14INX17IN, 24CMX24CM, 24CMX30CM |
| Magnification Type | (2010,0060) | M | BILINEAR, CUBIC, NONE, REPLICATE |
| Border Density | (2010,0100) | U | BLACK, WHITE |
| Max Density | (2010,0130) | U | > 0 |
| Min Density | (2010,0120) | U | 50 > value > 0 |
| Illumination | (2010,015E) | U | > 0 Required if Presentation LUT is present. |
| Reflective Ambient Light | (2010,0160) | U | > 0 Required if Presentation LUT is present. |
| Referenced Presentation LUT Sequence | (2050,0500) | U | |

M = Mandatory, **U** = User Option

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally and used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the DICOM print manager will issue a N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (N-ACTION-RQ and N-DELETE-RQ) on the Basic Film Box - see below:

| Attribute Name | Tag | Source of Information |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000) →(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Basic Film Box |

The Basic Film Box SOP Class interprets the following status codes:

Table 48 - Basic Film Box Status Codes

| Service Status | Meaning | Error Codes |
|----------------|---|-------------|
| Failure | Unable to create print job, print queue is full | C601 |
| | Image size is larger than images box size | C603 |

| Service Status | Meaning | Error Codes |
|----------------|---|-------------|
| Warning | Film box does not contain image box (empty page) | B603 |
| | Requested MinDensity or MaxDensity outside of Printer's operating range | B605 |
| Success | Film accepted for printing | 0000 |

Basic Grayscale Image Box SOP Class

The Basic Grayscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Grayscale Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the following attributes:

Table 49 - Basic Grayscale Image Box N-SET attributes

| Attribute Name | Tag | Usage SCU | Supported Values |
|--------------------------------|-------------|-----------|-----------------------|
| Image Position | (2020,0010) | M | 1 |
| Basic Grayscale Image Sequence | (2020,0110) | M | n.a. |
| > Samples per Pixel | (0028,0002) | M | 1 |
| > Photometric Interpretation | (0028,0004) | M | MONOCHROME2 |
| > Rows | (0028,0010) | M | <Printer/Film config> |
| > Columns | (0028,0011) | M | <Printer/Film config> |
| > Pixel Aspect Ratio | (0028,0034) | M | (1:1) |
| > Bits Allocated | (0028,0100) | M | 8, 16 |
| > Bits Stored | (0028,0101) | M | 8, 12 |
| > High Bit | (0028,0102) | M | 7, 11 |
| > Pixel Representation | (0028,0103) | M | 0 |
| > Pixel Data | (7FE0,0010) | M | |

M = Mandatory

The Grayscale Image Box SOP Class interprets the following status codes:

Table 50 - Basic Grayscale Image Box Status Codes

| Service Status | Meaning | Error Codes |
|----------------|---|-------------|
| Failure | Image contains more pixel than printer can print in Image Box | C603 |
| | Insufficient memory in printer to store the image | C605 |
| Warning | Requested MinDensity or MaxDensity outside of Printer's operating range | B605 |
| Success | | 0000 |

Basic Color Image Box SOP Class

The Basic Color Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Color Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Color Image Box SOP Class uses only the N-SET-RQ with the following attributes:

Table 51 - Basic Color Image Box N-SET attributes

| Attribute Name | Tag | Usage SCU | Supported Values |
|----------------------------|-------------|-----------|------------------|
| Image Position | (2020,0010) | M | 1 |
| Basic Color Image Sequence | (2020,0111) | M | n.a. |
| > Samples per Pixel | (0028,0002) | M | 3 |

| Attribute Name | Tag | Usage SCU | Supported Values |
|------------------------------|-------------|-----------|-----------------------|
| > Photometric Interpretation | (0028,0004) | M | RGB |
| >Planar Configuration | (0028,0006) | M | 1 |
| > Rows | (0028,0010) | M | <Printer/Film config> |
| > Columns | (0028,0011) | M | <Printer/Film config> |
| > Pixel Aspect Ratio | (0028,0034) | M | (1:1) |
| > Bits Allocated | (0028,0100) | M | 8 |
| > Bits Stored | (0028,0101) | M | 8, |
| > High Bit | (0028,0102) | M | 7 |
| > Pixel Representation | (0028,0103) | M | 0 |
| > Pixel Data | (7FE0,0010) | M | |

M = Mandatory

The Color Image Box SOP Class interprets the following status codes:

Table 52 - Basic Color Image Box Status Codes

| Service Status | Meaning | Error Codes |
|----------------|---|-------------|
| Failure | Image contains more pixel than printer can print in Image Box | C603 |
| | Insufficient memory in printer to store the image | C605 |
| Warning | Requested MinDensity or MaxDensity outside of Printer's operating range | B605 |
| Success | | 0000 |

Presentation LUT SOP Class

The Presentation LUT tailors image hardcopy printing for specific modalities, applications and user preferences.

The output of the Presentation LUT is Presentation Values (P-Values). P-Values are approximately related to human perceptual response. They are intended to facilitate common input for hardcopy. P-Values are intended to be independent of the specific class or characteristics of the hardcopy device.

The Presentation LUT SOP Class uses only the N-CREATE-RQ with the following attributes:

Table 53 - Presentation LUT N-CREATE attribute

| Attribute Name | Tag | Usage SCU | Supported Values |
|------------------------|-------------|-----------|------------------|
| Presentation LUT Shape | (2050,0020) | U | IDENTITY |

U = User Option

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and is used for later requests on the Basic Film Box (N-CREATE-RQ) and on the Presentation LUT (N-DELETE-RQ) - see below:

| Attribute Name | Tag | Source of Information |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000) →(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Presentation LUT |

The Presentation LUT SOP Class interprets the following status codes:

Table 54 - Presentation LUT Status Codes

| Service Status | Meaning | Codes |
|----------------|---------|-------|
|----------------|---------|-------|

| Service Status | Meaning | Codes |
|----------------|--|-------|
| Success | Presentation LUT successfully created | 0000 |
| Warning | Requested MinDensity or MaxDensity outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | B605 |

Printer SOP Class

The Printer SOP Class allows to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The SCU uses the mandatory N-EVENT Report DIMSE service to monitor the changes of the printer status in an asynchronous way.

The following returned information is supported:

Table 55 - Used Printer N-EVENT Report attributes

| Event-type Name | Event | Attributes | Tag | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal | 1 | | | |
| Warning | 2 | Printer Status Info | (2110,0020) | U |
| Failure | 3 | Printer Status Info | (2110,0020) | U |

U = User Option

Table 56 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes

| Attribute Name | Tag | Usage SCP | Supported Values |
|---------------------|-------------|-----------|----------------------------------|
| Printer Status | (2110,0010) | M | NORMAL, FAILURE, WARNING |
| Printer Status Info | (2110,0020) | M | See tables in Annex for details. |

M = Mandatory

Note: For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the appropriate Annex section.

Printer Job SOP Class

The Print Job SOP Class allows to monitor the execution of the print process.

The syngo® MR DICOM Print Management application supports the optional N-EVENT-REPORT DIMSE Service to receive the changes of the Print Job Status in an asynchronous way.

It can receive Events from the Print SCP asynchronously:

- N-EVENT-REPORT

Note: The underlying syngo DICOM Print AE does not support receiving of N-EVENT-REPORT messages from camera during open print sessions. This is typically configurable in the camera setup.

The following information is supported:

Table 57 - Used Print Job N-EVENT Report attributes

| Event-type Name | Event | Attributes | Tag | Usage SCU |
|-----------------|-------|-----------------------|-------------|--|
| Normal | 1 | Execution Status Info | (2100,0030) | U |
| | | Print Job ID | (2100,0010) | -- (Print Queue Management SOP Class not supported) |
| | | Film Session Label | (2000,0050) | U |

| Event-type Name | Event | Attributes | Tag | Usage SCU |
|-----------------|-------|-----------------------|-------------|--|
| | | Printer Name | (2110,0030) | U |
| Printing | 2 | Execution Status Info | (2100,0030) | U |
| | | Print Job ID | (2100,0010) | -- (Print Queue Management SOP Class not supported) |
| | | Film Session Label | (2000,0050) | U |
| | | Printer Name | (2110,0030) | U |
| | | | | |
| Done | 3 | Execution Status Info | (2100,0030) | U |
| | | Print Job ID | (2100,0010) | -- (Print Queue Management SOP Class not supported) |
| | | Film Session Label | (2000,0050) | U |
| | | Printer Name | (2110,0030) | U |
| | | | | |
| Failure | 4 | Execution Status Info | (2100,0030) | U |
| | | Print Job ID | (2100,0010) | -- (Print Queue Management SOP Class not supported) |
| | | Film Session Label | (2000,0050) | U |
| | | Printer Name | (2110,0030) | U |
| | | | | |

U = User Option

Note: For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the appropriate Annex section".

4.2.8.3.2 Activity - Show Device Status

4.2.8.3.2.1 Description and Sequencing of Activity

With no printing activity ongoing ("idle time"), the syngo® MR DICOM Print SCU application will cyclically request the printer status to update the related printer state in the Printing UI.

4.2.8.3.2.2 Proposed Presentation Context

The syngo® MR DICOM application will propose Presentation Contexts as shown in the following table:

Table 58 - Presentation Context Table "Show Device Status"

| Presentation Context Table | | | | | |
|----------------------------|------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.8.3.2.3 SOP Specific Conformance

The Printer SOP Class allows to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The Print SCU AE application will cyclically "ask" the Printer (SCP) for its status synchronously:

- N-GET as SCU

The following information is supported:

Table 59 - Used Printer N-EVENT Report attributes

| Event-type Name | Event | Attributes | Tag | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal | 1 | | | |
| Warning | 2 | Printer Status Info | (2110,0020) | U |
| Failure | 3 | Printer Status Info | (2110,0020) | U |

U = User Option

<modify If detailed status displays of the DICOM Print SCU are described in the Annex, add a reference to the Annex section below

Table 60 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes

| Attribute Name | Tag | Usage SCP | Supported Values |
|---------------------|-------------|-----------|----------------------------------|
| Printer Status | (2110,0010) | M | NORMAL, FAILURE, WARNING |
| Printer Status Info | (2110,0020) | M | See tables in Annex for details. |

M = Mandatory

Note: For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the Annex section....

>

4.2.8.4 Association Acceptance Policy

The syngo® MR DICOM application does not support Print Management Services as an SCP.

4.3 Network Interfaces

4.3.1 Physical Network Interface

Syngo® MR provides DICOM 3.0 TCP/IP network communication support as defined in Part 8 of the DICOM Standard. The network communication is independent from the physical medium over which TCP/IP executes; it inherits this from the Windows OS system upon which it executes.

4.3.2 Additional Protocols

none

4.3.3 IPv4 and IPv6 Support

IPv4 and IPv6 are supported. Regarding IPv6 please note, that the complete networking infrastructure in the hospital (firewalls, DNS-Servers ...) must support IPv6 in order to get a functioning communication.

4.4 Configuration

4.4.1 AE Title/Presentation Address Mapping

AE Titles shall be unique within the hospital. A common way to achieve that is to use the host-name as part of the AE Titles. The string can be up to 16 characters and must not contain any extended characters. Only 7-bit ASCII characters (excluding Control Characters) are allowed according to the DICOM Standard.

4.4.1.1 Local AE Titles

Syngo® MR allows configuring AETitles, Ports and Services in any wished way. Default delivery is that all services are using the same AE title and only one port number. In case the connected systems cannot handle this default, the customer service engineer is able to configure for each service its own AE title and Port number.

| Parameter | Configurable | Default Value |
|------------------|--------------|--|
| Default AE title | Yes | hostname in uppercase characters; limited to 16 characters |
| Default Port | Yes | 104 |

4.4.1.2 Remote AE Title/Presentation Address Mapping

4.4.1.2.1 Remote Association Initiators

All relevant remote applications that may setup DICOM associations towards syngo® MR need to be configured in syngo® MR, before the association can be established. This behavior is configurable but it is recommended, not to change this behavior.

The mapping of external AE Titles to TCP/IP addresses and ports is configurable and initially set at the time of installation by Installation Personnel. Changes can later on also be performed by the local system administrator. The Application Entity Titles and supported transfer syntaxes need to be known for configuration.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

4.4.1.2.2 Remote SCP's

For remote applications that shall be able to accept DICOM associations from syngo® MR, the following information needs to be available:

- Application Entity Title
- Host Name / IP address on which the remote application service runs
- Port number on which the remote application accepts association requests.

The remote system will be indicated in the UI of syngo® MR with a logical name, that is also entered when configuring the node in the administration UI.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

4.4.2 Parameters

The next table lists configuration parameters, which are true for all Application Entities.

Table 61: Parameter List

| Parameter | Configurable | Default Value |
|--|--------------|---------------|
| max PDU size | Yes | 32768 Bytes |
| time-out for accepting/rejecting an association request | Yes | 30 s |
| time-out for responding to an association open/close request | Yes | 30 s |
| time-out for accepting a message over network | Yes | 30 s |
| time-out for waiting for data between TCP/IP-packets | Yes | 5 s |
| time-outs for waiting for a Service Request/Response message from the remote node (Storage SCP/SCU) | Yes | 30 s |
| time-outs for waiting for a Service Request/Response message from the remote node (Query/Retrieve SCP/SCU) | Yes | 30 s |
| time-out for waiting for a C-MOVE-RSP | No | 1200 s |
| number of image collection before saving to database | Yes | 20 |
| max matches query limit | Yes | 100 |
| max number of parallel receiving associations | Yes | 12 |

5 Media Interchange

5.1 Implementation Model

5.1.1 Application Data Flow Diagram




Figure 4: Media Interchange Application Data Flow Diagram

Syngo® MR provides the functionality to Import or Export DICOM Instances from and to the File System. During export, a DicomDIR may also be generated (user selection). A complete ISO Image ready-to-burn can also be generated. All SOP Classes defined in Table 1 are supported for the Import/Export functionality.

5.1.2 Functional definitions of AEs

The syngo® MR application is capable of

- creating a new File-set in the File System (Export to ...)
- importing SOP Instances from the medium onto local storage
- writing the File-sets DICOMDIR information into the file system and joining it to an ISO image.

5.1.3 Sequencing of Real-World Activities

Whenever data is written to an external media, syngo® MR creates a DICOMDIR from the selected data and creates an ISO image of the selected data on the local hard disk. Depending on the selected data and options (selected media size, with or without compression) either General Purpose CD profile or DVD-J2K profile are used.




Figure 5: Sequence diagram – Media creation

5.1.4 File Meta Information for Implementation Class and Version

This section describes the values assigned to the File Meta Information attributes that pertain to the Implementation Class and Version. The implementation Class UID and the Implementation Version name in the File Meta Header are the same as the values specified for networking.

Table 62: Implementation Class/Version Name - Media Interchange

| | |
|--------------------------------------|----------------------------------|
| File Meta Information Version | 0001 |
| Implementation Class UID | 1.3.12.2.1107.5.8.15.10.20090701 |
| Implementation Version Name | |

5.2 AE SPECIFICATIONS

5.2.1 Media Storage AE – Specification

Syngo® MR provides conformance to the following Application Profiles as an FSC as well as an FSR. The FSU role is only supported only on a non-optical storage device (e.g. USB stick).

In addition augmented conformance is provided to store extra data attributes important for the full feature support of the syngo®-based products. Details are listed below:

Table 63: Media - Application Profiles and Real-World Activities

| Application Profiles Supported | Real-World Activity | Role | Service Class Option |
|---------------------------------------|--|--------------------|-----------------------------|
| AUG-GEN-CD | Browse Directory Information Import into Application Export to local Archive Media | FSR , FSC , FSU | Interchange |
| AUG-GEN-DVD | | | |
| AUG- GEN-DVD-J2K | | | |
| AUG- GEN-USB-J2K | | | |
| STD-GEN-CD | Browse Directory Information Import into Application Export to local Archive Media | FSR , FSC , FSU | Interchange |
| STD-GEN-DVD | | | |
| STD-GEN-DVD-J2K | | | |
| STD-GEN-USB-J2K | | | |

5.2.1.1 Real-World Activities

5.2.1.1.1 Activity “Browse Directory Information”

Syngo® MR acts as FSR using the interchange option when requested to read the media directory.

Syngo® MR will read the DICOMDIR and insert those directory entries that are valid for the application profiles supported, into a local database. The database then is used for browsing media contents.

Note: The "Icon Image Sequence" is also supported in DICOMDIR. But only those Icon Images with "Bits Allocated" (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the Browser.

5.2.1.1.2 Real World Activity "Import into Application"

Syngo® MR application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the application.

The SOP Instance selected from the media directory will be copied into the running Application. Only SOP Instances, that are valid for the application profile supported and supported by syngo® MR can be retrieved from media.

5.2.1.1.3 Real-World Activity "Export to local Archive Media"

The syngo® MR application acts as FSU (for media with existing DICOM file-set) or FSC (media not initialized) using the interchange option when requested to copy SOP Instances from the local storage to local Archive Medium. The activity as FSU is only possible as long as the local burning SW of syngo® MR has not already processed the generated ISO file.

The syngo® MR application will receive a list of SOP Instances to be copied to the local archive medium. Depending on the profile selected (Standard: uncompressed, with DICOMDIR; Patient: compressed with DICOMDIR) the SOP Instances will be taken and an ISO file is being generated that includes the DICOMDIR and the corresponding objects.

It is then up to the syngo® MR local configuration (if equipped with a local media burner) to burn the ISO file to the appropriate media.

5.2.1.2 SOP Classes and Transfer Syntaxes

These Application Profiles are based on the Media Storage Service Class with the Interchange Option. Syngo® MR provides Standard Conformance to the the SOP Classes listed in "Table 1: Network Services" section "SOP Classes Created by syngo® MR" and "SOP Classes Managed by syngo® MR" in the ["Conformance Statement Overview"](#).

Using the Application Profiles supporting compression (AUG- GEN-DVD-J2K, AUG- GEN-USB-J2K, STD-GEN-DVD-J2K, STD-GEN-USB-J2K) the following Transfer Syntaxes are supported:

Table 64: Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K

| UID value | Transfer Syntax | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2.1 | Explicit Value Representation Little Endian native | yes | yes |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed | yes | no |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed | yes | no |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | yes | no |

| UID value | Transfer Syntax | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed | yes | no |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed | yes | no |
| 1.2.840.10008.1.2.5 | RLE Lossless compressed | yes | no |

Using the Application Profiles that do not support compression (AUG- GEN-DVD, AUG- GEN-USB, STD-GEN-DVD, STD-GEN-USB) only Explicit Value Representation Little Endian (1.2.840.10008.1.2.1) is supported.

5.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

5.3.1 Augmented Application Profiles

The standard application profiles are augmented with private object Siemens CSA Non-Image.

Table 65: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles

| Information Object Definition | SOP Class UID | Transfer Syntax UID | FSC | FSR |
|-------------------------------|---------------------|---|-----|-----|
| CSA Non-Image Storage | 1.3.12.2.1107.5.9.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 | O | M |

The Siemens non-image is typically used for raw data and 3D private data.

5.4 MEDIA CONFIGURATION

none

6 Support of Extended Character Sets

The syngo® MR DICOM application supports the following character sets as defined in the four tables below:

Table 66: Single-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|----------------------|
| Default repertoire | none | ISO_IR 6 | ISO 646: |
| Latin alphabet No. 1 | ISO_IR 100 | ISO_IR 100 | Supplementary set |
| | | ISO_IR 6 | ISO 646: |
| Latin alphabet No. 2 | ISO_IR 101 | ISO_IR 101 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 3 | ISO_IR 109 | ISO_IR 109 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 4 | ISO_IR 110 | ISO_IR 110 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Cyrillic | ISO_IR 144 | ISO_IR 144 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Arabic | ISO_IR 127 | ISO_IR 127 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Greek | ISO_IR 126 | ISO_IR 126 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Hebrew | ISO_IR 138 | ISO_IR 138 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 5 | ISO_IR 148 | ISO_IR 148 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Japanese | ISO_IR 13 | ISO_IR 13 | JIS X 0201: Katakana |
| | | ISO_IR 14 | JIS X 0201: Romaji |
| Thai | ISO_IR166 | ISO_IR166 | TIS 620-253 (1990) |
| | | ISO_IR 6 | ISO 646 |

Table 67: Single-Byte Characters Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------|-------------------------|-------------------|
| Default repertoire | ISO 2022 IR 6 | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.1 | ISO 2022 IR 100 | ISO 2022 | ESC 02/13 04/01 | ISO-IR 100 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.2 | ISO 2022 IR 101 | ISO 2022 | ESC 02/13 04/02 | ISO-IR 101 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.3 | ISO 2022 IR 109 | ISO 2022 | ESC 02/13 04/03 | ISO-IR 109 | Supplementary set |

| | | | | | |
|---------------------|-----------------|----------|-----------------|------------|-------------------|
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.4 | ISO 2022 IR 110 | ISO 2022 | ESC 02/13 04/04 | ISO-IR 110 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |

Multi-Byte Character Sets without Code Extension

Table 68: Multi-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|---|
| Unicode | ISO_IR 192 | ISO 10646 | Unicode in UTF-8 |
| Chinese | GB18030 | GB18030 | GB 18030-2000 (China Association for Standardization) |

Table 69: Multi-Byte Character Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------------|-------------------------|-------------------------------------|
| Japanese | ISO 2022 IR 159 | ISO 2022 | ESC 02/04 02/08 04/04 | ISO-IR 159 | JIS X 0212: Supplementary Kanji set |
| Korean | ISO 2022 IR 149 | ISO 2022 | ESC 02/04 02/09 04/03 | ISO-IR 149 | KS X 1001: Hangul and Hanja |

All Special Character Sets (SCS) listed above are supported for incoming Data. When creating new Instances, the system will use the default SCS (or SCS List) configured on the machine.

When there is a mismatch between the given character set in attribute (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM conform:

- Convert each illegal character to a '?'.

There are three categories of character sets which have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO_IR 6, ISO 2022 IR 6, ISO_IR 100, etc.
→ encoded in ISO 2022
- ISO_IR 192 → encoded in UTF-8
- GB18030 → encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO_IR 192 \leftrightarrow (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030 \leftrightarrow (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022 \leftrightarrow (0008,0005) contains ISO_IR 192
- An attribute value is encoded in ISO 2022 \leftrightarrow (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies is not DICOM conform. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

Syngo® MR supports Kanji characters in the byte zones after 74 (79, 7A, 7B and 7C).

7 Attribute confidentiality profiles

7.1 De-identification

The syngo® MR application can de-identify attributes, when exporting to Media. Three different levels of de-identification are supported. The user needs to select the appropriate de-identification level during export.

For full and reduced anonymization private attributes are not included in anonymized Studies. For service anonymization all private attributes are included in anonymized Studies.

Note: reduced anonymization applies the following options: Retain UIDs , Patient Chars, Long Full Dates

In the following table for attributes marked with:

- 'Yes' - data are anonymized
- 'No' - data are kept

Table 70: Application Level Confidentiality Profile Attributes (standard tags)

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|---|------|---------|---------|
| (0000,1000) | Affected SOP Instance UID | Yes | No | No |
| (0000,1001) | Requested SOP Instance UID | Yes | No | No |
| (0002,0003) | Media Storage SOP Instance UID | Yes | No | No |
| (0004,1511) | Referenced SOP Instance UID in File | Yes | No | No |
| (0008,0014) | Instance Creator UID | Yes | No | No |
| (0008,0015) | Instance Coercion DateTime | Yes | No | No |
| (0008,0018) | SOP Instance UID | Yes | No | Yes |
| (0008,0020) | Study Date | Yes | No | No |
| (0008,0021) | Series Date | Yes | No | No |
| (0008,0022) | Acquisition Date | Yes | No | No |
| (0008,0023) | Content Date | Yes | No | No |
| (0008,0024) | Overlay Date | Yes | No | No |
| (0008,0025) | Curve Date | Yes | No | No |
| (0008,002A) | Acquisition DateTime | Yes | No | No |
| (0008,0030) | Study Time | Yes | No | No |
| (0008,0031) | Series Time | Yes | No | No |
| (0008,0032) | Acquisition Time | Yes | No | No |
| (0008,0033) | Content Time | Yes | No | No |
| (0008,0034) | Overlay Time | Yes | No | No |
| (0008,0035) | Curve Time | Yes | No | No |
| (0008,0050) | Accession Number | Yes | Yes | No |
| (0008,0058) | Failed SOP Instance UID List | Yes | No | No |
| (0008,0080) | Institution Name | Yes | Yes | No |
| (0008,0081) | Institution Address | Yes | Yes | No |
| (0008,0082) | Institution Code Sequence | Yes | Yes | No |
| (0008,0090) | Referring Physician's Name | Yes | Yes | Yes |
| (0008,0092) | Referring Physician's Address | Yes | Yes | Yes |
| (0008,0094) | Referring Physician's Telephone Numbers | Yes | Yes | Yes |
| (0008,0096) | Referring Physician's Identification Sequence | Yes | Yes | No |
| (0008,010D) | Context Group Extension Creator UID | Yes | No | No |
| (0008,0201) | Timezone Offset From UTC | Yes | No | No |
| (0008,1010) | Station Name | Yes | Yes | Yes |
| (0008,1030) | Study Description | Yes | Yes | No |
| (0008,103E) | Series Description | Yes | Yes | No |
| (0008,1040) | Institutional Department Name | Yes | Yes | No |
| (0008,1048) | Physician(s) of Record | Yes | Yes | Yes |

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|---|------|---------|---------|
| (0008,1049) | Physician(s) of Record Identification Sequence | Yes | Yes | No |
| (0008,1050) | Performing Physicians' Name | Yes | Yes | Yes |
| (0008,1052) | Performing Physicians' Identification Sequence | Yes | Yes | No |
| (0008,1060) | Name of Physician(s) Reading Study | Yes | Yes | Yes |
| (0008,1062) | Physician Reading Study Identification Sequence | Yes | Yes | No |
| (0008,1070) | Operators' Name | Yes | Yes | Yes |
| (0008,1072) | Operators' Identification Sequence | Yes | Yes | No |
| (0008,1080) | Admitting Diagnoses Description | Yes | Yes | No |
| (0008,1084) | Admitting Diagnoses Code Sequence | Yes | Yes | No |
| (0008,1110) | Referenced Study Sequence | Yes | No | No |
| (0008,1111) | Referenced Performed Procedure Step Sequence | Yes | No | No |
| (0008,1120) | Referenced Patient Sequence | Yes | Yes | No |
| (0008,1140) | Referenced Image Sequence | Yes | No | No |
| (0008,1155) | Referenced SOP Instance UID | Yes | No | No |
| (0008,1195) | Transaction UID | Yes | No | No |
| (0008,2111) | Derivation Description | Yes | No | No |
| (0008,2112) | Source Image Sequence | Yes | No | No |
| (0008,3010) | Irradiation Event UID | Yes | No | No |
| (0008,4000) | Identifying Comments | Yes | Yes | No |
| (0008,9123) | Creator Version UID | Yes | No | No |
| (0010,0010) | Patient's Name | Yes | Yes | Yes |
| (0010,0020) | Patient ID | Yes | Yes | Yes |
| (0010,0021) | Issuer of Patient ID | Yes | Yes | No |
| (0010,0030) | Patient's Birth Date | Yes | Yes | No |
| (0010,0032) | Patient's Birth Time | Yes | Yes | No |
| (0010,0040) | Patient's Sex | Yes | No | No |
| (0010,0050) | Patient's Insurance Plan Code Sequence | Yes | Yes | No |
| (0010,0101) | Patient's Primary Language Code Sequence | Yes | Yes | No |
| (0010,0102) | Patient's Primary Language Modifier Code Sequence | Yes | Yes | No |
| (0010,1000) | Other Patient IDs | Yes | Yes | Yes |
| (0010,1001) | Other Patient Names | Yes | Yes | Yes |
| (0010,1002) | Other Patient IDs Sequence | Yes | Yes | No |
| (0010,1005) | Patient's Birth Name | Yes | Yes | Yes |
| (0010,1010) | Patient's Age | Yes | No | No |
| (0010,1020) | Patient's Size | Yes | No | No |
| (0010,1030) | Patient's Weight | Yes | No | No |
| (0010,1040) | Patient Address | Yes | Yes | Yes |
| (0010,1050) | Insurance Plan Identification | Yes | Yes | No |
| (0010,1060) | Patient's Mother's Birth Name | Yes | Yes | Yes |
| (0010,1080) | Military Rank | Yes | Yes | No |
| (0010,1081) | Branch of Service | Yes | Yes | No |
| (0010,1090) | Medical Record Locator | Yes | Yes | No |
| (0010,1100) | Referenced Patient Photo Sequence | Yes | Yes | No |
| (0010,2000) | Medical Alerts | Yes | Yes | No |
| (0010,2110) | Allergies | Yes | Yes | No |
| (0010,2150) | Country of Residence | Yes | Yes | No |
| (0010,2152) | Region of Residence | Yes | Yes | No |
| (0010,2154) | Patient's Telephone Number | Yes | Yes | Yes |
| (0010,2160) | Ethnic Group | Yes | No | No |
| (0010,2180) | Occupation | Yes | Yes | No |
| (0010,21A0) | Smoking Status | Yes | No | No |
| (0010,21B0) | Additional Patient's History | Yes | Yes | Yes |
| (0010,21C0) | Pregnancy Status | Yes | No | No |
| (0010,21D0) | Last Menstrual Date | Yes | No | No |
| (0010,21F0) | Patient's Religious Preference | Yes | Yes | No |
| (0010,2203) | Patient Sex Neutered | Yes | No | No |
| (0010,2297) | Responsible Person | Yes | Yes | No |

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|---|------|---------|---------|
| (0010,2299) | Responsible Organization | Yes | Yes | No |
| (0010,4000) | Patient Comments | Yes | Yes | Yes |
| (0018,0010) | Contrast Bolus Agent | Yes | Yes | No |
| (0018,1000) | Device Serial Number | Yes | Yes | No |
| (0018,1002) | Device UID | Yes | No | No |
| (0018,1004) | Plate ID | Yes | Yes | No |
| (0018,1005) | Generator ID | Yes | Yes | No |
| (0018,1007) | Cassette ID | Yes | Yes | No |
| (0018,1008) | Gantry ID | Yes | Yes | No |
| (0018,1030) | Protocol Name | Yes | Yes | No |
| (0018,1400) | Acquisition Device Processing Description | Yes | Yes | No |
| (0018,2042) | Target UID | Yes | No | No |
| (0018,4000) | Acquisition Comments | Yes | Yes | No |
| (0018,700A) | Detector ID | Yes | Yes | No |
| (0018,9424) | Acquisition Protocol Description | Yes | Yes | No |
| (0018,9516) | Start Acquisition DateTime | Yes | No | No |
| (0018,9517) | End Acquisition DateTime | Yes | No | No |
| (0018,A003) | Contribution Description | Yes | Yes | No |
| (0020,000D) | Study Instance UID | Yes | No | Yes |
| (0020,000E) | Series Instance UID | Yes | No | Yes |
| (0020,0010) | Study ID | Yes | Yes | No |
| (0020,0052) | Frame of Reference UID | Yes | No | No |
| (0020,0200) | Synchronization Frame of Reference UID | Yes | No | No |
| (0020,3401) | Modifying Device ID | Yes | Yes | No |
| (0020,3404) | Modifying Device Manufacturer | Yes | Yes | No |
| (0020,3406) | Modified Image Description | Yes | Yes | No |
| (0020,4000) | Image Comments | Yes | Yes | No |
| (0020,9158) | Frame Comments | Yes | Yes | No |
| (0020,9161) | Concatenation UID | Yes | No | No |
| (0020,9164) | Dimension Organization UID | Yes | No | No |
| (0028,1199) | Palette Color Lookup Table UID | Yes | No | No |
| (0028,1214) | Large Palette Color Lookup Table UID | Yes | No | No |
| (0028,4000) | Image Presentation Comments | Yes | Yes | No |
| (0032,0012) | Study ID Issuer | Yes | Yes | No |
| (0032,1020) | Scheduled Study Location | Yes | Yes | No |
| (0032,1021) | Scheduled Study Location AE Title | Yes | Yes | No |
| (0032,1030) | Reason for Study | Yes | Yes | No |
| (0032,1032) | Requesting Physician | Yes | Yes | No |
| (0032,1033) | Requesting Service | Yes | Yes | No |
| (0032,1060) | Requested Procedure Description | Yes | Yes | No |
| (0032,1070) | Requested Contrast Agent | Yes | Yes | No |
| (0032,4000) | Study Comments | Yes | Yes | No |
| (0038,0004) | Referenced Patient Alias Sequence | Yes | Yes | No |
| (0038,0010) | Admission ID | Yes | Yes | No |
| (0038,0011) | Issuer of Admission ID | Yes | Yes | No |
| (0038,001E) | Scheduled Patient Institution Residence | Yes | Yes | No |
| (0038,0020) | Admitting Date | Yes | No | No |
| (0038,0021) | Admitting Time | Yes | No | No |
| (0038,0040) | Discharge Diagnosis Description | Yes | Yes | No |
| (0038,0050) | Special Needs | Yes | Yes | No |
| (0038,0060) | Service Episode ID | Yes | Yes | No |
| (0038,0061) | Issuer of Service Episode ID | Yes | Yes | No |
| (0038,0062) | Service Episode Description | Yes | Yes | No |
| (0038,0300) | Current Patient Location | Yes | Yes | No |
| (0038,0400) | Patient's Institution Residence | Yes | Yes | No |
| (0038,0500) | Patient State | Yes | Yes | No |
| (0038,4000) | Visit Comments | Yes | Yes | No |

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|---|------|---------|---------|
| (0040,0001) | Scheduled Station AE Title | Yes | Yes | No |
| (0040,0002) | Scheduled Procedure Step Start Date | Yes | No | No |
| (0040,0003) | Scheduled Procedure Step Start Time | Yes | No | No |
| (0040,0004) | Scheduled Procedure Step End Date | Yes | No | No |
| (0040,0005) | Scheduled Procedure Step End Time | Yes | No | No |
| (0040,0006) | Scheduled Performing Physician Name | Yes | Yes | No |
| (0040,0007) | Scheduled Procedure Step Description | Yes | Yes | No |
| (0040,000B) | Scheduled Performing Physician Identification Sequence | Yes | Yes | No |
| (0040,0010) | Scheduled Station Name | Yes | Yes | No |
| (0040,0011) | Scheduled Procedure Step Location | Yes | Yes | No |
| (0040,0012) | Pre-Medication | Yes | Yes | No |
| (0040,0241) | Performed Station AE Title | Yes | Yes | No |
| (0040,0242) | Performed Station Name | Yes | Yes | No |
| (0040,0243) | Performed Location | Yes | Yes | No |
| (0040,0244) | Performed Procedure Step Start Date | Yes | No | No |
| (0040,0245) | Performed Procedure Step Start Time | Yes | No | No |
| (0040,0248) | Performed Station Name Code Sequence | Yes | Yes | No |
| (0040,0250) | Performed Procedure Step End Date | Yes | No | No |
| (0040,0251) | Performed Procedure Step End Time | Yes | No | No |
| (0040,0253) | Performed Procedure Step ID | Yes | Yes | No |
| (0040,0254) | Performed Procedure Step Description | Yes | Yes | No |
| (0040,0275) | Request Attributes Sequence | Yes | Yes | No |
| (0040,0280) | Comments on Performed Procedure Step | Yes | Yes | No |
| (0040,0555) | Acquisition Context Sequence | Yes | Yes | No |
| (0040,1001) | Requested Procedure ID | Yes | Yes | No |
| (0040,1004) | Patient Transport Arrangements | Yes | Yes | No |
| (0040,1005) | Requested Procedure Location | Yes | Yes | No |
| (0040,1010) | Names of Intended Recipient of Results | Yes | Yes | No |
| (0040,1011) | Intended Recipients of Results Identification Sequence | Yes | Yes | No |
| (0040,1101) | Person Identification Code Sequence | Yes | Yes | No |
| (0040,1102) | Person Address | Yes | Yes | No |
| (0040,1103) | Person Telephone Numbers | Yes | Yes | No |
| (0040,1400) | Requested Procedure Comments | Yes | Yes | No |
| (0040,2001) | Reason for Imaging Service Request | Yes | Yes | No |
| (0040,2008) | Order Entered By | Yes | Yes | No |
| (0040,2009) | Order Enterer Location | Yes | Yes | No |
| (0040,2010) | Order Callback Phone Number | Yes | Yes | No |
| (0040,2016) | Placer Order Number of Imaging Service Request | Yes | Yes | No |
| (0040,2017) | Filler Order Number of Imaging Service Request | Yes | Yes | No |
| (0040,2400) | Imaging Service Request Comments | Yes | Yes | No |
| (0040,3001) | Confidentiality Constraint on Patient Data Description | Yes | Yes | No |
| (0040,4005) | Scheduled Procedure Step Start DateTime | Yes | No | No |
| (0040,4010) | Scheduled Procedure Step Modification DateTime | Yes | No | No |
| (0040,4011) | Expected Completion Date Time | Yes | No | No |
| (0040,4023) | Referenced General Purpose Scheduled Procedure Step Transaction UID | Yes | No | No |
| (0040,4025) | Scheduled Station Name Code Sequence | Yes | Yes | No |
| (0040,4027) | Scheduled Station Geographic Location Code Sequence | Yes | Yes | No |
| (0040,4030) | Performed Station Geographic Location Code Sequence | Yes | Yes | No |
| (0040,4034) | Scheduled Human Performers Sequence | Yes | Yes | No |
| (0040,4035) | Actual Human Performers Sequence | Yes | Yes | No |
| (0040,4036) | Human Performers Organization | Yes | Yes | No |
| (0040,4037) | Human Performers Name | Yes | Yes | No |
| (0040,4050) | Performed Procedure Step Start DateTime | Yes | No | No |
| (0040,4051) | Performed Procedure Step End DateTime | Yes | No | No |
| (0040,4052) | Procedure Step Cancellation DateTime | Yes | No | No |
| (0040,A027) | Verifying Organization | Yes | Yes | No |

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|---|------|---------|---------|
| (0040,A073) | Verifying Observer Sequence | Yes | Yes | No |
| (0040,A075) | Verifying Observer Name | Yes | Yes | No |
| (0040,A078) | Author Observer Sequence | Yes | Yes | No |
| (0040,A07A) | Participant Sequence | Yes | Yes | No |
| (0040,A07C) | Custodial Organization Sequence | Yes | Yes | No |
| (0040,A088) | Verifying Observer Identification Code Sequence | Yes | Yes | No |
| (0040,A123) | Person Name | Yes | Yes | No |
| (0040,A124) | UID | Yes | Yes | No |
| (0040,A171) | Observation UID | Yes | No | No |
| (0040,A172) | Referenced Observation UID (Trial) | Yes | No | No |
| (0040,A192) | Observation Date (Trial) | Yes | No | No |
| (0040,A193) | Observation Time (Trial) | Yes | No | No |
| (0040,A307) | Current Observer (Trial) | Yes | Yes | No |
| (0040,A352) | Verbal Source (Trial) | Yes | Yes | No |
| (0040,A353) | Address (Trial) | Yes | Yes | No |
| (0040,A354) | Telephone Number (Trial) | Yes | Yes | No |
| (0040,A358) | Verbal Source Identifier Code Sequence (Trial) | Yes | Yes | No |
| (0040,A402) | Observation Subject UID (Trial) | Yes | No | No |
| (0040,A730) | Content Sequence | Yes | Yes | No |
| (0040,DB0C) | Template Extension Organization UID | Yes | No | No |
| (0040,DB0D) | Template Extension Creator UID | Yes | No | No |
| (0070,0001) | Graphic Annotation Sequence | Yes | Yes | No |
| (0070,0084) | Content Creator's Name | Yes | Yes | No |
| (0070,0086) | Content Creator's Identification Code Sequence | Yes | Yes | No |
| (0070,031A) | Fiducial UID | Yes | No | No |
| (0088,0140) | Storage Media Fileset UID | Yes | No | No |
| (0088,0200) | Icon Image Sequence | Yes | Yes | No |
| (0088,0904) | Topic Title | Yes | Yes | No |
| (0088,0906) | Topic Subject | Yes | Yes | No |
| (0088,0910) | Topic Author | Yes | Yes | No |
| (0088,0912) | Topic Keywords | Yes | Yes | No |
| (0400,0100) | Digital Signature UID | Yes | Yes | No |
| (0400,0402) | Referenced Digital Signature Sequence | Yes | Yes | No |
| (0400,0403) | Referenced SOP Instance MAC Sequence | Yes | Yes | No |
| (0400,0404) | MAC | Yes | Yes | No |
| (0400,0550) | Modified Attributes Sequence | Yes | Yes | No |
| (0400,0561) | Original Attributes Sequence | Yes | Yes | No |
| (2030,0020) | Text String | Yes | Yes | No |
| (3006,0024) | Referenced Frame of Reference UID | Yes | No | No |
| (3006,00C2) | Related Frame of Reference UID | Yes | No | No |
| (3008,0105) | Source Serial Number | No | No | No |
| (300A,0013) | Dose Reference UID | Yes | No | No |
| (300E,0008) | Reviewer Name | Yes | Yes | No |
| (4000,0010) | Arbitrary | Yes | Yes | No |
| (4000,4000) | Text Comments | Yes | Yes | No |
| (4008,0042) | Results ID Issuer | Yes | Yes | No |
| (4008,0102) | Interpretation Recorder | Yes | Yes | No |
| (4008,010A) | Interpretation Transcriber | Yes | Yes | No |
| (4008,010B) | Interpretation Text | Yes | Yes | No |
| (4008,010C) | Interpretation Author | Yes | Yes | No |
| (4008,0111) | Interpretation Approver Sequence | Yes | Yes | No |
| (4008,0114) | Physician Approving Interpretation | Yes | Yes | No |
| (4008,0115) | Interpretation Diagnosis Description | Yes | Yes | No |
| (4008,0118) | Results Distribution List Sequence | Yes | Yes | No |
| (4008,0119) | Distribution Name | Yes | Yes | No |
| (4008,011A) | Distribution Address | Yes | Yes | No |
| (4008,0202) | Interpretation ID Issuer | Yes | Yes | No |

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|-------------|-----------------------------|------|---------|---------|
| (4008,0300) | Impressions | Yes | Yes | No |
| (4008,4000) | Results Comments | Yes | Yes | No |
| (50xx,xxxx) | Curve Data | Yes | Yes | No |
| (60xx,0100) | Overlay Bits Allocated | Yes | Yes | No |
| (60xx,0102) | Overlay Bit Position | Yes | Yes | No |
| (60xx,3000) | Overlay Data | Yes | Yes | No |
| (60xx,4000) | Overlay Comments | Yes | Yes | No |
| (FFFA,FFFA) | Digital Signatures Sequence | Yes | Yes | No |
| (FFFC,FFFC) | Data Set Trailing Padding | Yes | Yes | No |

Table 71: Application Level Confidentiality Profile Attributes (private tags)

| DICOM Tag | Attribute Name | Full | Reduced | Service |
|--|----------------------------------|------|---------|---------|
| (0019, SIEMENS CT VAO COAD, 90) | Osteo offset | Yes | No | No |
| (0019, SIEMENS CT VAO COAD, 92) | Osteo Regression Line Slope | Yes | No | No |
| (0019, SIEMENS CT VAO COAD, 93) | Osteo Regression Line Intercept | Yes | No | No |
| (0019, SIEMENS CT VAO COAD, 96) | Osteo Phantom Number | Yes | No | No |
| (0043, GEMS_PARM_01, 1E) | GE Delta Start Time | Yes | No | No |
| (0029, SIEMENS CSA ENVELOPE, 10) | Syngo Report Data | Yes | No | No |
| (0029, SIEMENS CSA ENVELOPE, 11) | Syngo Report Presentation | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 08) | Modality Image Header Type | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 09) | Modality Image Header Version | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 10) | Modality Image Header Info | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 18) | Modality Series Header Type | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 19) | Modality Series Header Version | Yes | No | No |
| (0029, SIEMENS CSA HEADER, 20) | Modality Series Header Info | Yes | No | No |
| (0029, SIEMENS MEDCOM HEADER, 40) | Application Header Sequence | Yes | No | No |
| (0029, SIEMENS MEDCOM HEADER, 41) | Application Header Type | Yes | No | No |
| (0029, SIEMENS MEDCOM HEADER, 42) | Application Header ID | Yes | No | No |
| (0029, SIEMENS MEDCOM HEADER, 43) | Application Header Version | Yes | No | No |
| (0029, SIEMENS MEDCOM HEADER, 44) | Application Header Info | Yes | No | No |
| (0029, SIEMENS CT APPL DATASET, 00) | Dual Energy Algorithm Parameters | Yes | No | No |
| (0029, SIEMENS CT APPL ALG PARAMS, 20) | Perfusion Result Set Id | Yes | No | No |

8 Security

8.1 Security Profiles

Time Synchronization Profiles: syngo® MR acts as an NTP Client in the Maintain Time Transaction.

8.2 Association Level Security

It is possible to configure whether the SCP will only answer to known AETs or to any AET.

8.3 Application Level Security

- User must login with own password
- For configuration and Maintenance, Service Technician must login with a separate password.

9 Annexes

9.1 IOD Contents

9.1.1 Enhanced MR Image IOD

This chapter describes the DICOM attributes of Enhanced MR Image Instances performed by the MR acquisition.

9.1.1.1 Patient Module

| Attribute Name | Tag | Supported Values |
|--------------------------|-------------|-------------------------------|
| Patient's Name | (0010,0010) | RIS defined or set by creator |
| Patient ID | (0010,0020) | RIS defined or set by creator |
| Issuer of Patient ID | (0010,0021) | RIS defined |
| Patient's Birth Date | (0010,0030) | RIS defined or set by creator |
| Patient's Birth Time | (0010,0032) | RIS defined or set by creator |
| Patient's Sex | (0010,0040) | RIS defined or set by creator |
| Other Patient IDs | (0010,1000) | RIS defined or set by creator |
| Other Patient Names | (0010,1001) | RIS defined or set by creator |
| Ethnic Group | (0010,2160) | RIS defined or set by creator |
| Patient Comments | (0010,4000) | RIS defined or set by creator |
| Other Patient Names | (0010,1001) | RIS defined or set by creator |
| Patient Identity Removed | (0012,0062) | set by creator |

9.1.1.2 General Study Module

| Attribute Name | Tag | Supported Values |
|----------------------------|-------------|-------------------------------|
| Study Instance UID | (0020,000D) | RIS defined or set by creator |
| Study Date | (0008,0020) | set by creator |
| Study Time | (0008,0030) | set by creator |
| Referring Physician's Name | (0008,0090) | RIS defined or set by creator |
| Study ID | (0020,0010) | set by creator |
| Accession Number | (0008,0050) | RIS defined or set by creator |
| Study Description | (0008,1030) | set by creator |
| Procedure Code Sequence | (0008,1032) | RIS defined |
| >Code Value | (0008,0100) | RIS defined |
| >Code Scheme Designator | (0008,0102) | RIS defined |
| >Code Meaning | (0008,0104) | RIS defined |

9.1.1.3 Patient Study Module

| Attribute Name | Tag | Supported Values |
|------------------|-------------|-------------------------------|
| Patient's Age | (0010,1010) | RIS defined or set by creator |
| Patient's Size | (0010,1020) | RIS defined or set by creator |
| Patient's Weight | (0010,1030) | RIS defined or set by creator |

9.1.1.4 General Series Module

| Attribute Name | Tag | Supported Values |
|--|-------------|-------------------------------|
| Modality | (0008,0060) | MR |
| Series Instance UID | (0020,000E) | set by creator |
| Series Number | (0020,0011) | set by creator |
| Series Date | (0008,0021) | set by creator |
| Series Time | (0008,0031) | set by creator |
| Performing Physicians' Name | (0008,1050) | RIS defined or set by creator |
| Protocol Name | (0018,1030) | set by creator |
| Series Description | (0008,103E) | set by creator |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator |
| >Referenced SOP Class UID | (0008,1150) | set by creator |
| >Referenced SOP Instance UID | (0008,1155) | set by creator |
| Body Part Examined | (0018,0015) | set by creator |
| Patient Position | (0018,5100) | set by creator |
| Request Attributes Sequence | (0040,0275) | RIS defined |
| >Requested Procedure ID | (0040,1001) | RIS defined |
| >Accession Number | (0008,0050) | RIS defined |
| >Study Instance UID | (0020,000D) | RIS defined |
| >Requested Procedure Description | (0032,1060) | RIS defined |
| >Scheduled Procedure Step ID | (0040,0009) | RIS defined |
| >Scheduled Procedure Step Description | (0040,0007) | RIS defined |
| > Referenced Study Sequence | (0008,1110) | RIS defined |
| >> Referenced SOP Class UID | (0008,1150) | RIS defined |
| >> Referenced SOP Instance UID | (0008,1155) | RIS defined |
| > Requested Procedure Code Sequence | (0032,1064) | RIS defined |
| >> Code Value | (0008,0100) | RIS defined |
| >> Code Scheme Designator | (0008,0102) | RIS defined |
| >> Code Meaning | (0008,0104) | RIS defined |
| > Scheduled Protocol Code Sequence | (0040,0008) | RIS defined |
| >> Code Value | (0008,0100) | RIS defined |
| >> Code Scheme Designator | (0008,0102) | RIS defined |
| >> Code Meaning | (0008,0104) | RIS defined |
| Performed Procedure Step ID | (0040,0253) | set by creator |
| Performed Procedure Step Start Date | (0040,0244) | set by creator |
| Performed Procedure Step Start Time | (0040,0245) | set by creator |
| Performed Procedure Step Description | (0040,0254) | set by creator |

9.1.1.5 MR Series Module

| Attribute Name | Tag | Supported Values |
|----------------|-------------|------------------|
| Modality | (0008,0060) | MR |

9.1.1.6 Frame of Reference Module

| Attribute Name | Tag | Supported Values |
|------------------------------|-------------|------------------|
| Frame of Reference UID | (0020,0052) | set by creator |
| Position Reference Indicator | (0020,1040) | empty |

9.1.1.7 General Equipment Module

| Attribute Name | Tag | Supported Values |
|---------------------------|-------------|-------------------------------|
| Manufacturer | (0008,0070) | Siemens HealthCare GmbH |
| Institution Name | (0008,0080) | RIS defined or set by creator |
| Institution Address | (0008,0081) | RIS defined set by creator |
| Manufacturer's Model Name | (0008,1090) | set by creator |
| Device Serial Number | (0018,1000) | set by creator |
| Software Versions | (0018,1020) | set by creator |

9.1.1.8 Enhanced General Equipment Module

| Attribute Name | Tag | |
|----------------|-------------|-------------------------|
| Manufacturer | (0008,0070) | Siemens HealthCare GmbH |

9.1.1.9 Image Pixel Module

| Attribute Name | Tag | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel | (0028,0002) | 1 |
| Photometric Interpretation | (0028,0004) | MONOCHROME2 |
| Rows | (0028,0010) | set by creator |
| Columns | (0028,0011) | set by creator |
| Bits Allocated | (0028,0100) | 16 |
| Bits Stored | (0028,0101) | 12 |
| High Bit | (0028,0102) | 11 |
| Pixel Representation | (0028,0103) | 0 |
| Pixel Data | (7FE0,0010) | set by creator |
| Smallest Image Pixel Value | (0028,0106) | set by creator |
| Largest Image Pixel Value | (0028,0107) | set by creator |

9.1.1.10 Enhanced Contrast/Bolus Module

| Attribute Name | Tag | Supported Values |
|---|-------------|------------------|
| Contrast/Bolus Agent Sequence | (0018,0012) | set by creator |
| >Code Value | (0008,0100) | set by creator |
| >Code Scheme Designator | (0008,0102) | set by creator |
| >Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Agent Number | (0018,9337) | 1 |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator |
| >>Code Value | (0008,0100) | set by creator |
| >>Code Scheme Designator | (0008,0102) | set by creator |
| >>Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Ingredient Code Sequence | (0018,9338) | set by creator |
| >>Code Value | (0008,0100) | set by creator |
| >>Code Scheme Designator | (0008,0102) | set by creator |
| >>Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Volume | (0018,1041) | set by creator |
| >Contrast/Bolus Ingredient Concentration | (0018,1049) | set by creator |

9.1.1.11 Multi-frame Functional Groups Module

| Attribute Name | Tag | Supported Values |
|------------------|-------------|------------------|
| Instance Number | (0020,0013) | set by creator |
| Content Date | (0008,0023) | set by creator |
| Content Time | (0008,0033) | set by creator |
| Number of Frames | (0028,0008) | set by creator |

9.1.1.12 Multi-frame Dimension Module

| Attribute Name | Tag | Supported Values |
|---------------------------------|-------------|------------------|
| Dimension Organization Sequence | (0020,9221) | set by creator |
| Dimension Index Sequence | (0020,9222) | set by creator |

9.1.1.13 Enhanced MR Image Functional Groups

| Functional Group Macro | Attribute Name | Tag | Supported Values |
|------------------------|--------------------------------------|-------------|------------------|
| Pixel Measures | Pixel Measures Sequence | (0028,9110) | set by creator |
| | >Pixel Spacing | (0018,0030) | set by creator |
| | >Slice Thickness | (0028,0050) | set by creator |
| Frame Content | Frame Content Sequence | (0020,9111) | set by creator |
| | >Frame Acquisition Number | (0020,9156) | set by creator |
| | >Frame Reference Date Time | (0018,9151) | set by creator |
| | >Frame Acquisition Date Time | (0018,9074) | set by creator |
| | >Frame Acquisition Duration | (0018,9220) | set by creator |
| | >Dimension Index Values | (0020,9157) | set by creator |
| | >Temporal Position Index | (0020,9128) | set by creator |
| | >Stack ID | (0020,9056) | set by creator |
| | >In-Stack Position Number | (0020,9057) | set by creator |
| | >Frame Comments | (0020,9158) | set by creator |
| Plane Position | Plane Position Sequence | (0020,9113) | set by creator |
| | >Image Position (Patient) | (0020,0032) | set by creator |
| Plane Orientation | Plane Orientation Sequence | (0020,9116) | set by creator |
| | >Image Orientation (Patient) | (0020,0037) | set by creator |
| Referenced Image | Referenced Image Sequence | (0008,1140) | set by creator |
| | >Referenced SOP Class UID | (0008,1150) | set by creator |
| | >Referenced SOP Instance UID | (0008,1155) | set by creator |
| | >Referenced Frame Number | (0008,1160) | set by creator |
| | >Purpose of Referenced Code Sequence | (0040,A170) | set by creator |
| | >>Code Value | (0008,0100) | e.g. 121311 |
| | >>Coding Scheme Designator | (0008,0102) | e.g. DCM |
| | >>Code Meaning | (0008,0104) | e.g. Localizer |
| Derivation Image | Derivation Image Sequence | (0008,9124) | set by creator |
| | >Derivation Code Sequence | (0008,9215) | set by creator |
| | >>Code Value | (0008,0100) | set by creator |
| | >>Coding Scheme Designator | (0008,0102) | set by creator |
| | >>Code Meaning | (0008,0104) | set by creator |
| | >Source Image Sequence | (0008,2112) | set by creator |
| | >>Referenced SOP Class UID | (0008,1150) | set by creator |
| | >>Referenced SOP Instance UID | (0008,1155) | set by creator |

| | | | |
|----------------------------------|---|-------------|----------------|
| | >>Referenced Frame Number | (0008,1160) | set by creator |
| | >>Purpose of Referenced Code Sequence | (0040,A170) | set by creator |
| | >>>Code Value | (0008,0100) | set by creator |
| | >>>Coding Scheme Designator | (0008,0102) | set by creator |
| | >>>Code Meaning | (0008,0104) | set by creator |
| Cardiac Synchronization | Cardiac Synchronization Sequence | (0018,9118) | set by creator |
| | >Nominal Percentage of Cardiac Phase | (0020,9241) | set by creator |
| | >Nominal Cardiac Trigger Delay Time | (0020,9153) | set by creator |
| | >Actual Cardiac Trigger Delay Time | (0020,9252) | set by creator |
| | >Nominal Cardiac Trigger Time Prior to R-peak | (0020,9154) | set by creator |
| | >Actual Cardiac Trigger Time Prior to R-peak | (0020,9155) | set by creator |
| | >Intervals Acquired | (0018,1083) | set by creator |
| | >Intervals Rejected | (0018,1084) | set by creator |
| | >Heart Rate | (0018,1088) | set by creator |
| | >R-R Interval Time Nominal | (0020,9251) | set by creator |
| | >Low R-R Value | (0018,1081) | set by creator |
| | >High R-R Value | (0018,1082) | set by creator |
| Frame Anatomy | Frame Anatomy Sequence | (0020,9071) | set by creator |
| | >Frame Laterality | (0020,9072) | set by creator |
| | >Anatomic Region Sequence | (0008,2218) | set by creator |
| | >>Code Value | (0008,0100) | set by creator |
| | >>Coding Scheme Designator | (0008,0102) | set by creator |
| | >>Code Meaning | (0008,0104) | set by creator |
| Pixel Value Transformation | Pixel Value Transformation Sequence | (0028,9145) | set by creator |
| | >Rescale Intercept | (0028,1052) | set by creator |
| | >Rescale Slope | (0028,1053) | set by creator |
| | >Rescale Type | (0028,1054) | set by creator |
| Frame VOI LUT | Frame VOI LUT Sequence | (0028,9132) | set by creator |
| | >Window Center | (0028,1050) | set by creator |
| | >Window Width | (0028,1051) | set by creator |
| Real World Value Mapping | Real World Value Mapping Sequence | (0040,9096) | set by creator |
| | >Real World Value Intercept | (0040,9224) | set by creator |
| | >Real World Value Slope | (0040,9225) | set by creator |
| | >Measurement Units Code Sequence | (0040,08EA) | set by creator |
| | >>Code Value | (0008,0100) | set by creator |
| | >>Coding Scheme Designator | (0008,0102) | set by creator |
| | >>Code Meaning | (0008,0104) | set by creator |
| Contrast/Bolus Usage | Contrast/Bolus Usage Sequence | (0018,9341) | set by creator |
| | >Contrast/Bolus Agent Number | (0018,9337) | 1 |
| | >Contrast/Bolus Agent Administered | (0018,9342) | set by creator |
| | >Contrast/Bolus Agent Detected | (0018,9343) | set by creator |
| | >Contrast/Bolus Agent Phase | (0018,9344) | set by creator |
| MR Image Frame Type | MR Image Frame Type Sequence | (0018,9226) | set by creator |
| | >Frame Type | (0008,9007) | set by creator |
| | >Pixel Presentation | (0008,9205) | set by creator |
| | >Volumetric Properties | (0008,9206) | set by creator |
| | >Volume Based Calculation Technique | (0008,9207) | set by creator |
| | >Complex Image Component | (0008,9208) | set by creator |
| | >Acquisition Contrast | (0008,9209) | set by creator |
| | >Functional Settling Phase Frames Present | (0018,9622) | set by creator |
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence | (0018,9112) | set by creator |
| | >Repetition Time | (0018,0080) | set by creator |
| | >Flip Angle | (0018,1314) | set by creator |

| | | | |
|---------------------|---|---|--|
| | >Echo Train Length >RF Echo Train Length >Gradient Echo Train Length >Specific Absorption Rate Sequence >>Specific Absorption Rate Definition >>Specific Absorption Rate Value >Gradient Output Type >Gradient Output >Operation Mode Sequence >>Operating Mode Type >>Operating Mode | (0018,0091) (0018,9240) (0018,9241) (0018,9239) (0018,9179) (0018,9181) (0018,9180) (0018,9182) (0018,9176) (0018,9177) (0018,9178) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR FOV/Geometry | MR FOV/Geometry Sequence >In-plane Phase Encoding Direction >MR Acquisition Frequency Encoding Steps >MR Acquisition Phase Encoding Steps in-plane >MR Acquisition Phase Encoding Steps out-of-plane >Percent Sampling >Percent Phase Field of View | (0018,9125) (0018,1312) (0018,9058) (0018,9231) (0018,9232) (0018,0093) (0018,0094) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Echo | MR Echo Sequence >Effective Echo Time | (0018,9114) (0018,9082) | set by creator set by creator |
| MR Modifier | MR Modifier Sequence >Inversion Recovery >Inversion Times >Flow Compensation >Flow Compensation Direction >Spoiling >T2 Preparation >Spectrally Selected Excitation >Spatial Pre-saturation >Partial Fourier >Partial Fourier Direction >Parallel Acquisition >Parallel Acquisition Technique >Parallel Reduction Factor In-plane >Parallel Reduction Factor out-of-plane | (0018,9115) (0018,9009) (0018,9079) (0018,9010) (0018,9183) (0018,9016) (0018,9021) (0018,9026) (0018,9027) (0018,9081) (0018,9036) (0018,9077) (0018,9078) (0018,9069) (0018,9155) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Imaging Modifier | MR Imaging Modifier Sequence >Magnetization Transfer >Blood Signal Nulling >Tagging >Tag Spacing First Dimension >Tag Spacing Second Dimension >Tag Angle First Axis >Tag Angle Second Axis >Tag Thickness >Tagging Delay >Transmitter Frequency >Pixel Bandwidth | (0018,9006) (0018,9020) (0018,9022) (0018,9028) (0018,9030) (0018,9218) (0018,9019) (0018,9219) (0018,9035) (0018,9184) (0018,9098) (0018,0095) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Receive Coil | MR Receive Coil Sequence >Receive Coil Name >Receive Coil Manufacturer Name >Receive Coil Type >Quadrature Receive Coil >Multi-Coil Definition Sequence | (0018,9042) (0018,1250) (0018,9041) (0018,9043) (0018,9044) (0018,9045) | set by creator set by creator set by creator set by creator set by creator set by creator |

| | | | |
|---------------------------|--|-------------|----------------|
| | >>Multi-Coil Element Name | (0018,9047) | set by creator |
| | >>Multi-Coil Element Used | (0018,9048) | set by creator |
| MR Transmit Coil | MR Transmit Coil Sequence | (0018,9049) | set by creator |
| | >Transmit Coil Name | (0018,1251) | set by creator |
| | >Transmit Coil Manufacturer Name | (0018,9050) | set by creator |
| | >Transmit Coil Type | (0018,9051) | set by creator |
| MR Diffusion | MR Diffusion Sequence | (0018,9117) | set by creator |
| | >Diffusion b-value | (0018,9087) | set by creator |
| | >Diffusion Directionality | (0018,9075) | set by creator |
| | >Diffusion Gradient Direction Sequence | (0018,9076) | set by creator |
| | >>Diffusion Gradient Orientation | (0018,9089) | set by creator |
| | >Diffusion b-matrix Sequence | (0018,9601) | set by creator |
| | >>Diffusion b-value XX | (0018,9602) | set by creator |
| | >>Diffusion b-value XY | (0018,9603) | set by creator |
| | >>Diffusion b-value XZ | (0018,9604) | set by creator |
| | >>Diffusion b-value YY | (0018,9605) | set by creator |
| | >>Diffusion b-value YZ | (0018,9606) | set by creator |
| | >>Diffusion b-value ZZ | (0018,9607) | set by creator |
| MR Averages | MR Averages Sequence | (0018,9119) | set by creator |
| | >Number of Averages | (0018,0083) | set by creator |
| MR Velocity Encoding | MR Velocity Encoding Sequence | (0018,9197) | set by creator |
| | >Velocity Encoding Direction | (0018,9090) | set by creator |
| | >Velocity Encoding Minimum Value | (0018,9091) | set by creator |
| | >Velocity Encoding Maximum Value | (0018,9217) | set by creator |
| MR Arterial Spin Labeling | MR Arterial Spin Labeling Sequence | (0018,9251) | set by Creator |
| | >ASL Technique Description | (0018,9252) | set by Creator |
| | >ASL Context | (0018,9257) | set by Creator |
| | >ASL Slab Sequence | (0018,9260) | set by Creator |
| | >>ASL Slab Number | (0018,9253) | set by Creator |
| | >>ASL Slab Thickness | (0018,9254) | set by Creator |
| | >>ASL Slab Orientation | (0018,9255) | set by Creator |
| | >>ASL Mid Slab Position | (0018,9256) | set by Creator |
| | >>ASL Pulse Train Duration | (0018,9258) | set by Creator |
| | >ASL Crusher Flag | (0018,9259) | set by Creator |
| | >ASL Crusher Flow Limit | (0018,925A) | set by Creator |
| | >ASL Crusher Description | (0018,925B) | set by Creator |
| | >ASLBolus Cut-off Flag | (0018,925C) | set by Creator |
| | >ASLBolus Cut-off Timing Sequence | (0018,925D) | set by Creator |
| | >>ASL Bolus Cut-off Delay Time | (0018,925F) | set by Creator |
| | >>ASL Bolus Cut-off Technique | (0018,925E) | set by Creator |

9.1.1.14 Cardiac Synchronization Module

| Attribute Name | Tag | Supported Values |
|-----------------------------------|-------------|------------------|
| Cardiac Synchronization Technique | (0018,9037) | set by creator |
| Cardiac Signal Source | (0018,9085) | set by creator |
| Cardiac RR Interval Specified | (0018,9070) | set by creator |
| Low R-R Value | (0018,1081) | set by creator |
| High R-R Value | (0018,1082) | set by creator |
| Intervals Acquired | (0018,1083) | set by creator |
| Intervals Rejected | (0018,1084) | set by creator |

9.1.1.15 Respiratory Synchronization Nodule

| Attribute Name | Tag | Supported Values |
|---|-------------|------------------|
| Respiratory Motion Compensation Technique | (0018,9170) | set by creator |

9.1.1.16 Bulk Motion Synchronization Module

| Annex A: Attribute Name | Annex B: Tag | Annex C: Supported Values |
|------------------------------------|--------------|---------------------------|
| Bulk Motion Compensation Technique | (0018,9172) | set by creator |
| Bulk Motion Signal Source | (0018,9173) | set by creator |

9.1.1.17 Supplemental Palette Color Lookup Table Module

| Annex D: Attribute Name | Annex E: Tag | Annex F: Supported Values |
|---|--------------|---------------------------|
| Red Palette Color Lookup Table Descriptor | (0028,1101) | set by creator |
| Green Palette Color Lookup Table Descriptor | (0028,1102) | set by creator |
| Blue Palette Color Lookup Table Descriptor | (0028,1103) | set by creator |
| Red Palette Color Lookup Table Data | (0028,1201) | set by creator |
| Green Palette Color Lookup Table Data | (0028,1202) | set by creator |
| Blue Palette Color Lookup Table Data | (0028,1203) | set by creator |

9.1.1.18 Acquisition Context Module

| Attribute Name | Tag | Supported Values |
|------------------------------|-------------|------------------|
| Acquisition Context Sequence | (0040,0555) | TID 99_3100 |

9.1.1.18.1 MR Acquisition Context TID 99_3100

The MR Acquisition Context TID 99_3100 is used in the Acquisition Context Sequence (0040,0555) of Enhanced MR Images.

Type: Extensible

| | NL | VT | Concept Name | VM | Req Ty- pe | Condition | Value Set Constraint |
|---|----|------|---------------------------------------|----|---------------|-----------|----------------------|
| 1 | | CODE | EV (A-52, 99SMS_CTMR, "MR Technique") | 1 | M | | DCID(A-200) |

9.1.1.18.2 MR Technique Context Group A-200

Type: Extensible

| Coding Scheme Designator (0008,0102) | Coding Scheme Version (0008,0103) | Code Value (0008,0100) | Code Meaning (0008,0104) |
|--------------------------------------|-----------------------------------|------------------------|--------------------------|
| 99SMS_CTMR | 1.0 | A-200 | GRASP |
| 99SMS_CTMR | 1.0 | A-202 | FastView |
| 99SMS_CTMR | 1.0 | A-204 | Angio |

9.1.1.19 Enhanced MR Image Module

| Attribute Name | Tag | Supported Values |
|------------------------------------|-------------|------------------|
| Acquisition Number | (0020,0012) | set by creator |
| Acquisition Date Time | (0008,002A) | set by creator |
| Acquisition Duration | (0018,9073) | set by creator |
| Content Qualification | (0018,9004) | set by creator |
| Resonant Nucleus | (0018,9100) | set by creator |
| k-space Filtering | (0018,9064) | set by creator |
| Magnetic Field Strength | (0018,0087) | set by creator |
| Applicable Safety Standard Agency | (0018,9174) | set by creator |
| Image Comments | (0020,4000) | set by creator |
| Image Type | (0008,0008) | set by creator |
| Pixel Presentation | (0008,9205) | set by creator |
| Volumetric Properties | (0008,9206) | set by creator |
| Volume Based Calculation Technique | (0008,9207) | set by creator |
| Complex Image Component | (0008,9208) | set by creator |
| Acquisition Contrast | (0008,9209) | set by creator |
| Samples per Pixel | (0028,0002) | 1 |
| Photometric Interpretation | (0028,0004) | MONOCHROME2 |
| Bits Allocated | (0028,0100) | 16 |
| Bits Stored | (0028,0101) | 12 |
| High Bit | (0028,0102) | 11 |
| Pixel Representation | (0028,0103) | 0 |
| Burned In Annotation | (0028,0301) | NO |
| Lossy Image Compression | (0028,2110) | 00 |
| Presentation LUT Shape | (2050,0020) | IDENTITY |

9.1.1.20 MR Pulse Sequence Module

| Attribute Name | Tag | Supported Values |
|-------------------------------------|-------------|------------------|
| Pulse Sequence Name | (0018,9005) | set by creator |
| MR Acquisition Type | (0018,0023) | set by creator |
| Echo Pulse Sequence | (0018,9008) | set by creator |
| Multiple Spin Echo | (0018,9011) | set by creator |
| Multi-planar Excitation | (0018,9012) | set by creator |
| Phase Contrast | (0018,9014) | set by creator |
| Time of Flight Contrast | (0018,9015) | set by creator |
| Steady State Pulse Sequence | (0018,9017) | set by creator |
| Echo Planar Pulse Sequence | (0018,9018) | set by creator |
| Saturation Recovery | (0018,9024) | set by creator |
| Spectral Selected Suppression | (0018,9025) | set by creator |
| Oversampling Phase | (0018,9029) | set by creator |
| Geometry of k-Space | (0018,9032) | set by creator |
| Rectilinear Phase Encode Reordering | (0018,9034) | set by creator |
| Segmented k-Space | (0018,9033) | set by creator |
| Coverage of k-Space | (0018,9094) | set by creator |
| Number of k-Space Trajectories | (0018,9093) | set by creator |

9.1.1.21 SOP Common Module

| Attribute Name | Tag | Supported Values |
|------------------------|-------------|-----------------------------------|
| SOP Class UID | (0008,0016) | Enhanced MR Storage SOP Class UID |
| SOP Instance UID | (0008,0018) | set by creator |
| Specific Character Set | (0008,0005) | set by creator |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

9.1.2 Enhanced MR Color Image IOD

The Enhanced MR Color Image IOD supports the same DICOM attributes as the Enhanced MR Image IOD with the exception of the Photometric Interpretation which is RGB.

9.1.2.1 Image Pixel Module

| Attribute Name | Tag | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel | (0028,0002) | 3 |
| Photometric Interpretation | (0028,0004) | RGB |
| Rows | (0028,0010) | set by creator |
| Columns | (0028,0011) | set by creator |
| Bits Allocated | (0028,0100) | 8 |
| Bits Stored | (0028,0101) | 8 |
| High Bit | (0028,0102) | 7 |
| Pixel Representation | (0028,0103) | 0 |

9.1.3 MR Spectroscopy IOD

This chapter describes the DICOM attributes of MR Spectroscopy Instances performed by the MR acquisition.

9.1.3.1 Patient Module

| Attribute Name | Tag | Supported Values |
|--------------------------|-------------|-------------------------------|
| Patient's Name | (0010,0010) | RIS defined or set by creator |
| Patient ID | (0010,0020) | RIS defined or set by creator |
| Issuer of Patient ID | (0010,0021) | RIS defined |
| Patient's Birth Date | (0010,0030) | RIS defined or set by creator |
| Patient's Birth Time | (0010,0032) | RIS defined or set by creator |
| Patient's Sex | (0010,0040) | RIS defined or set by creator |
| Other Patient IDs | (0010,1000) | RIS defined or set by creator |
| Other Patient Names | (0010,1001) | RIS defined or set by creator |
| Ethnic Group | (0010,2160) | RIS defined or set by creator |
| Patient Comments | (0010,4000) | RIS defined or set by creator |
| Other Patient Names | (0010,1001) | RIS defined or set by creator |
| Patient Identity Removed | (0012,0062) | set by creator |

9.1.3.2 General Study Module

| Attribute Name | Tag | Supported Values |
|----------------------------|-------------|-------------------------------|
| Study Instance UID | (0020,000D) | RIS defined or set by creator |
| Study Date | (0008,0020) | set by creator |
| Study Time | (0008,0030) | set by creator |
| Referring Physician's Name | (0008,0090) | RIS defined or set by creator |
| Study ID | (0020,0010) | set by creator |
| Accession Number | (0008,0050) | RIS defined or set by creator |
| Study Description | (0008,1030) | set by creator |
| Procedure Code Sequence | (0008,1032) | RIS defined |
| >Code Value | (0008,0100) | RIS defined |
| >Code Scheme Designator | (0008,0102) | RIS defined |
| >Code Meaning | (0008,0104) | RIS defined |

9.1.3.3 Patient Study Module

| Attribute Name | Tag | Supported Values |
|------------------|-------------|-------------------------------|
| Patient's Age | (0010,1010) | RIS defined or set by creator |
| Patient's Size | (0010,1020) | RIS defined or set by creator |
| Patient's Weight | (0010,1030) | RIS defined or set by creator |

9.1.3.4 General Series Module

| Attribute Name | Tag | Supported Values |
|--|-------------|-------------------------------|
| Modality | (0008,0060) | MR |
| Series Instance UID | (0020,000E) | set by creator |
| Series Number | (0020,0011) | set by creator |
| Series Date | (0008,0021) | set by creator |
| Series Time | (0008,0031) | set by creator |
| Performing Physicians' Name | (0008,1050) | RIS defined or set by creator |
| Protocol Name | (0018,1030) | set by creator |
| Series Description | (0008,103E) | set by creator |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator |
| >Referenced SOP Class UID | (0008,1150) | set by creator |
| >Referenced SOP Instance UID | (0008,1155) | set by creator |
| Body Part Examined | (0018,0015) | set by creator |
| Patient Position | (0018,5100) | set by creator |
| Request Attributes Sequence | (0040,0275) | RIS defined |
| >Requested Procedure ID | (0040,1001) | RIS defined |
| >Accession Number | (0008,0050) | RIS defined |
| >Study Instance UID | (0020,000D) | RIS defined |
| >Requested Procedure Description | (0032,1060) | RIS defined |
| >Scheduled Procedure Step ID | (0040,0009) | RIS defined |
| >Scheduled Procedure Step Description | (0040,0007) | RIS defined |
| > Referenced Study Sequence | (0008,1110) | RIS defined |
| >> Referenced SOP Class UID | (0008,1150) | RIS defined |
| >> Referenced SOP Instance UID | (0008,1155) | RIS defined |
| > Requested Procedure Code Sequence | (0032,1064) | RIS defined |
| >> Code Value | (0008,0100) | RIS defined |

| | | |
|--------------------------------------|-------------|----------------|
| >> Code Scheme Designator | (0008,0102) | RIS defined |
| >> Code Meaning | (0008,0104) | RIS defined |
| > Scheduled Protocol Code Sequence | (0040,0008) | RIS defined |
| >> Code Value | (0008,0100) | RIS defined |
| >> Code Scheme Designator | (0008,0102) | RIS defined |
| >> Code Meaning | (0008,0104) | RIS defined |
| Performed Procedure Step ID | (0040,0253) | set by creator |
| Performed Procedure Step Start Date | (0040,0244) | set by creator |
| Performed Procedure Step Start Time | (0040,0245) | set by creator |
| Performed Procedure Step Description | (0040,0254) | set by creator |

9.1.3.5 MR Series Module

| Attribute Name | Tag | Supported Values |
|----------------|-------------|------------------|
| Modality | (0008,0060) | MR |

9.1.3.6 Frame of Reference Module

| Attribute Name | Tag | Supported Values |
|------------------------------|-------------|------------------|
| Frame of Reference UID | (0020,0052) | set by creator |
| Position Reference Indicator | (0020,1040) | empty |

9.1.3.7 General Equipment Module

| Attribute Name | Tag | Supported Values |
|---------------------------|-------------|-------------------------|
| Manufacturer | (0008,0070) | Siemens HealthCare GmbH |
| Institution Name | (0008,0080) | set by creator |
| Institution Address | (0008,0081) | set by creator |
| Manufacturer's Model Name | (0008,1090) | set by creator |
| Device Serial Number | (0018,1000) | set by creator |
| Software Versions | (0018,1020) | set by creator |

9.1.3.8 Enhanced General Equipment Module

| Attribute Name | Tag | Supported Values |
|----------------|-------------|-------------------------|
| Manufacturer | (0008,0070) | Siemens HealthCare GmbH |

9.1.3.9 Enhanced Contrast/Bolus Module

| Attribute Name | Tag | Supported Values |
|---|-------------|------------------|
| Contrast/Bolus Agent Sequence | (0018,0012) | set by creator |
| >Code Value | (0008,0100) | set by creator |
| >Code Scheme Designator | (0008,0102) | set by creator |
| >Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Agent Number | (0018,9337) | 1 |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator |
| >>Code Value | (0008,0100) | set by creator |
| >>Code Scheme Designator | (0008,0102) | set by creator |
| >>Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Ingredient Code Sequence | (0018,9338) | set by creator |

| | | |
|--|-------------|----------------|
| >>Code Value | (0008,0100) | set by creator |
| >>Code Scheme Designator | (0008,0102) | set by creator |
| >>Code Meaning | (0008,0104) | set by creator |
| >Contrast/Bolus Volume | (0018,1041) | set by creator |
| >Contrast/Bolus Ingredient Concentration | (0018,1049) | set by creator |

9.1.3.10 Multi-frame Functional Groups Module

| Attribute Name | Tag | Supported Values |
|-----------------|-------------|------------------|
| Instance Number | (0020,0013) | set by creator |
| Content Date | (0008,0023) | set by creator |
| Content Time | (0008,0033) | set by creator |

9.1.3.11 Multi-frame Dimension Module

| Attribute Name | Tag | Supported Values |
|---------------------------------|-------------|------------------|
| Dimension Organization Sequence | (0020,9221) | set by creator |
| Dimension Index Sequence | (0020,9222) | set by creator |

9.1.3.12 MR Spectroscopy Functional Groups

| Functional Group Macro | Attribute Name | Tag | Supported Values |
|---------------------------|--------------------------------------|-------------|------------------|
| Pixel Measures | Pixel Measures Sequence | (0028,9110) | set by creator |
| | >Pixel Spacing | (0018,0030) | set by creator |
| | >Slice Thickness | (0028,0050) | set by creator |
| Frame Content | Frame Content Sequence | (0020,9111) | set by creator |
| | >Frame Acquisition Number | (0020,9156) | set by creator |
| | >Frame Reference Date Time | (0018,9151) | set by creator |
| | >Frame Acquisition Date Time | (0018,9074) | set by creator |
| | >Frame Acquisition Duration | (0018,9220) | set by creator |
| | >Dimension Index Values | (0020,9157) | set by creator |
| | >Temporal Position Index | (0020,9128) | set by creator |
| | >Stack ID | (0020,9056) | set by creator |
| | >In-Stack Position Number | (0020,9057) | set by creator |
| Plane Position | >Frame Comments | (0020,9158) | set by creator |
| | Plane Position Sequence | (0020,9113) | set by creator |
| Plane Orientation | >Image Position (Patient) | (0020,0032) | set by creator |
| | Plane Orientation Sequence | (0020,9116) | set by creator |
| Referenced Image | >Image Orientation (Patient) | (0020,0037) | set by creator |
| | Referenced Image Sequence | (0008,1140) | set by creator |
| | >Referenced SOP Class UID | (0008,1150) | set by creator |
| | >Referenced SOP Instance UID | (0008,1155) | set by creator |
| | >Referenced Frame Number | (0008,1160) | set by creator |
| | >Purpose of Referenced Code Sequence | (0040,A170) | set by creator |
| | >>Code Value | (0008,0100) | set by creator |
| | >>Coding Scheme Designator | (0008,0102) | set by creator |
| Derivation Image | >>Code Meaning | (0008,0104) | set by creator |
| | Derivation Image Sequence | (0008,9124) | set by creator |
| | >Derivation Code Sequence | (0008,9215) | set by creator |
| | >>Code Value | (0008,0100) | set by creator |
| | >>Coding Scheme Designator | (0008,0102) | set by creator |

| | | | |
|----------------------------------|---|--|--|
| | >>Code Meaning >Source Image Sequence >>Referenced SOP Class UID >>Referenced SOP Instance UID >>Referenced Frame Number >>Purpose of Referenced Code Sequence >>>Code Value >>>Coding Scheme Designator >>>Code Meaning | (0008,0104) (0008,2112) (0008,1150) (0008,1155) (0008,1160) (0040,A170) (0008,0100) (0008,0102) (0008,0104) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| Frame Anatomy | Frame Anatomy Sequence >Frame Laterality >Anatomic Region Sequence >>Code Value >>Coding Scheme Designator >>Code Meaning | (0020,9071) (0020,9072) (0008,2218) (0008,0100) (0008,0102) (0008,0104) | set by creator set by creator set by creator set by creator set by creator set by creator |
| Contrast/Bolus Usage | Contrast/Bolus Usage Sequence >Contrast/Bolus Agent Number >Contrast/Bolus Agent Administered >Contrast/Bolus Agent Detected >Contrast/Bolus Agent Phase | (0018,9341) (0018,9337) (0018,9342) (0018,9343) (0018,9344) | set by creator 1 set by creator set by creator set by creator |
| MR Spectroscopy Frame Type | MR Spectroscopy Frame Type Sequence >Frame Type >Volumetric Properties >Volume Based Calculation Technique >Complex Image Component >Acquisition Contrast | (0018,9227) (0008,9007) (0008,9206) (0008,9207) (0008,9208) (0008,9209) | set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence >Repetition Time >Flip Angle >Echo Train Length >RF Echo Train Length >Gradient Echo Train Length >Specific Absorption Rate Sequence >>Specific Absorption Rate Definition >>Specific Absorption Rate Value >Gradient Output Type >Gradient Output >Operation Mode Sequence >>Operating Mode Type >>Operating Mode | (0018,9112) (0018,0080) (0018,1314) (0018,0091) (0018,9240) (0018,9241) (0018,9239) (0018,9179) (0018,9181) (0018,9180) (0018,9182) (0018,9176) (0018,9177) (0018,9178) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Spectroscopy FOV/Geometry | MR Spectroscopy FOV/Geometry Sequence >Spectroscopy Acquisition Data Columns >Spectroscopy Acquisition Phase Rows >Spectroscopy Acquisition Phase Columns >Spectroscopy Acquisition Out-of-Plane Phase Steps >Percent Sampling >Percent Phase Field of View | (0018,9103) (0018,9127) (0018,9095) (0018,9234) (0018,9159) (0018,0093) (0018,0094) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Echo | MR Echo Sequence >Effective Echo Time | (0018,9114) (0018,9082) | set by creator set by creator |
| MR Modifier | MR Modifier Sequence >Inversion Recovery >Inversion Times >Flow Compensation | (0018,9115) (0018,9009) (0018,9079) (0018,9010) | set by creator set by creator set by creator set by creator |

| | | | |
|----------------------|--|--|--|
| | >Flow Compensation Direction >Spoiling >T2 Preparation >Spectrally Selected Excitation >Spatial Pre-saturation >Partial Fourier >Partial Fourier Direction >Parallel Acquisition >Parallel Acquisition Technique >Parallel Reduction Factor In-plane >Parallel Reduction Factor out-of-plane | (0018,9183) (0018,9016) (0018,9021) (0018,9026) (0018,9027) (0018,9081) (0018,9036) (0018,9077) (0018,9078) (0018,9069) (0018,9155) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Receive Coil | MR Receive Coil Sequence >Receive Coil Name >Receive Coil Manufacturer Name >Receive Coil Type >Quadrature Receive Coil >Multi-Coil Definition Sequence >>Multi-Coil Element Name >>Multi-Coil Element Used | (0018,9042) (0018,1250) (0018,9041) (0018,9043) (0018,9044) (0018,9045) (0018,9047) (0018,9048) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Transmit Coil | MR Transmit Coil Sequence >Transmit Coil Name >Transmit Coil Manufacturer Name >Transmit Coil Type | (0018,9049) (0018,1251) (0018,9050) (0018,9051) | set by creator set by creator set by creator set by creator |
| MR Diffusion | MR Diffusion Sequence >Diffusion b-value >Diffusion Directionality >Diffusion Gradient Direction Sequence >>Diffusion Gradient Orientation >Diffusion b-matrix Sequence >>Diffusion b-value XX >>Diffusion b-value XY >>Diffusion b-value XZ >>Diffusion b-value YY >>Diffusion b-value YZ >>Diffusion b-value ZZ | (0018,9117) (0018,9087) (0018,9075) (0018,9076) (0018,9089) (0018,9601) (0018,9602) (0018,9603) (0018,9604) (0018,9605) (0018,9606) (0018,9607) | set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator set by creator |
| MR Averages | MR Averages Sequence >Number of Averages | (0018,9119) (0018,0083) | set by creator set by creator |
| MR Velocity Encoding | MR Velocity Encoding Sequence >Velocity Encoding Direction >Velocity Encoding Minimum Value >Velocity Encoding Maximum Value | (0018,9197) (0018,9090) (0018,9091) (0018,9217) | set by creator set by creator set by creator set by creator |

9.1.3.13 Cardiac Synchronization Module

| Attribute Name | Tag | Supported Values |
|-----------------------------------|-------------|------------------|
| Cardiac Synchronization Technique | (0018,9037) | set by creator |
| Cardiac Signal Source | (0018,9085) | set by creator |
| Cardiac RR Interval Specified | (0018,9070) | set by creator |
| Low R-R Value | (0018,1081) | set by creator |
| High R-R Value | (0018,1082) | set by creator |
| Intervals Acquired | (0018,1083) | set by creator |
| Intervals Rejected | (0018,1084) | set by creator |

9.1.3.14 Respiratory Synchronization Module

| Attribute Name | Tag | Supported Values |
|---|-------------|------------------|
| Respiratory Motion Compensation Technique | (0018,9170) | set by creator |

9.1.3.15 Bulk Motion Synchronization Module

| Attribute Name | Tag | Supported Values |
|------------------------------------|-------------|------------------|
| Bulk Motion Compensation Technique | (0018,9172) | set by creator |
| Bulk Motion Signal Source | (0018,9173) | set by creator |

9.1.3.16 Acquisition Context Module

| Attribute Name | Tag | Supported Values |
|------------------------------|-------------|------------------|
| Acquisition Context Sequence | (0040,0555) | empty |

9.1.3.17 MR Spectroscopy Module

| Attribute Name | Tag | Supported Values |
|--------------------------------------|-------------|------------------|
| Acquisition Number | (0020,0012) | set by creator |
| Acquisition Date Time | (0018,002A) | set by creator |
| Acquisition Duration | (0018,9073) | set by creator |
| Referenced Raw Data Sequence | (0008,9121) | set by creator |
| Content Qualification | (0018,9004) | PRODUCT |
| Resonant Nucleus | (0018,9100) | set by creator |
| k-space Filtering | (0018,9064) | set by creator |
| Magnetic Field Strength | (0018,0087) | set by creator |
| Applicable Safety Standard Agency | (0018,9174) | set by creator |
| Image Comments | (0020,4000) | set by creator |
| Image Type | (0008,0008) | set by creator |
| Volumetric Properties | (0008,9206) | set by creator |
| Volume Based Calculation Technique | (0008,9207) | set by creator |
| Complex Image Component | (0008,9208) | set by creator |
| Acquisition Contrast | (0008,9209) | set by creator |
| Transmitter Frequency | (0018,9098) | set by creator |
| Spectral Width | (0018,9052) | set by creator |
| Chemical Shift Reference | (0018,9053) | set by creator |
| Volume Localisation Technique | (0018,9054) | set by creator |
| Volume Localization Sequence | (0018,9126) | set by creator |
| >Slab Thickness | (0018,9104) | set by creator |
| >Slab Orientation | (0018,9105) | set by creator |
| >Mid SlabPosition | (0018,9106) | set by creator |
| De-coupling | (0018,9059) | set by creator |
| De-coupling Nucleus | (0018,9060) | set by creator |
| De-coupling Frequency | (0018,9061) | set by creator |
| De-coupling Method | (0018,9062) | set by creator |
| De-coupling Chemical Shift Reference | (0018,9063) | set by creator |
| Time Domain Filtering | (0018,9065) | set by creator |
| Number Of Zero Fills | (0018,9066) | set by creator |
| Baseline Correction | (0018,9067) | set by creator |

| | | |
|-------------------------------------|-------------|----------------|
| Frequency Correction | (0018,9101) | set by creator |
| First Order Phase Correction | (0018,9198) | set by creator |
| Water Referenced Phase Correction | (0018,9199) | set by creator |
| Water Reference Acquisition | (0018,9297) | set by creator |
| Referenced Instance Sequence | (0008,114A) | set by creator |
| >Purpose of Reference Code Sequence | (0040,A170) | set by creator |

9.1.3.18 MR Spectroscopy Pulse Sequence Module

| Attribute Name | Tag | Supported Values |
|-------------------------------------|-------------|------------------|
| Pulse Sequence Name | (0018,9005) | set by creator |
| MR Spectroscopy Acquisition Type | (0018,9200) | set by creator |
| Echo Pulse Sequence | (0018,9008) | set by creator |
| Multi Spin Echo | (0018,9011) | set by creator |
| Multi-planar Excitation | (0018,9012) | set by creator |
| Steady State Pulse Sequence | (0018,9017) | set by creator |
| Echo Planar Pulse Sequence | (0018,9018) | set by creator |
| Spectrally Selected Suppression | (0018,9025) | set by creator |
| Geometry of k-Space Traversal | (0018,9032) | set by creator |
| Rectilinear Phase Encode Reordering | (0018,9034) | set by creator |
| Segmented k-Space Traversal | (0018,9033) | set by creator |
| Coverage of k-Space | (0018,9094) | set by creator |
| Number of k-Space Trajectories | (0018,9093) | set by creator |
| Echo Peak Position | (0018,9298) | set by creator |

9.1.3.19 MR Spectroscopy Data Module

| Attribute Name | Tag | Supported Values |
|------------------------------------|-------------|------------------|
| Rows | (0028,0010) | set by creator |
| Columns | (0028,0011) | set by creator |
| Data Point Rows | (0028,9001) | set by creator |
| Data Point Columns | (0028,9002) | set by creator |
| Data Representation | (0028,9108) | set by creator |
| Signal Domain Columns | (0028,9003) | set by creator |
| Signal Domains Rows | (0028,9235) | set by creator |
| First Order Phase Correction Angle | (5600,0010) | set by creator |
| Spectroscopy Data | (5600,0020) | set by creator |

9.1.3.20 SOP Common Module

| Attribute Name | Tag | Supported Values |
|------------------------|-------------|-----------------------------------|
| SOP Class UID | (0008,0016) | MR Spectroscopy SOP Class UID |
| SOP Instance UID | (0008,0018) | set by creator |
| Specific Character Set | (0008,0005) | set by creator |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

9.1.4 Raw Data

Neuro Diffusion Tensor Imaging data are encoded in the Raw Data SOP Class.

9.1.5 Evidence Documents

Evidence Documents will be created by applications e.g. cardiac analysis to store evaluated results. The Evidence Documents are encoded in the SR Enhanced SOP Class.

9.1.5.1 Evidence Document Templates

The finding related results of applications are stored into the Content Sequence of a Structured Evidence Document. The Content Items of a Content Sequence are specified in Structured Reporting Templates.

Examples of Structured Reporting Templates are:

- BI-RADS reporting
- Cardiac reporting
- PI-RADS 2 reporting
- PhoenixZip documentation.

9.2 Data Dictionary of Private Attributes

The following table Table 72: Private Data Element Dictionary lists all private attributes created by syngo® MR which may be included in the generated instances. These private attributes may be deprecated or replaced with standard attributes in the future.

Table 72: Private Data Element Dictionary

| DICOM Tag | Name | VR | VM |
|---|---------------------------------|----|----|
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,01) | Business Unit Code | CS | 1 |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,02) | Application Type | LO | 1 |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,03) | Application Attributes Sequence | SQ | 1 |
| (0029,SIEMENS SYNGO FUNCTION ASSIGNMENT,01) | Data Reference | LO | 1 |
| (0009,SIEMENS SYNGO INDEX SERVICE,20) | Object Insertion Date | DA | 1 |
| (0009,SIEMENS SYNGO INDEX SERVICE,A0) | Sender System Device Name | LO | 1 |
| (0029,SIEMENS SYNGO VOLUME,12) | Slices | US | 1 |
| (0029,SIEMENS SYNGO VOLUME,14) | Volume Histogram | OB | 1 |
| (0029,SIEMENS SYNGO VOLUME,18) | Volume Level | IS | 1 |
| (0029,SIEMENS SYNGO VOLUME,30) | Voxel Spacing | DS | 3 |
| (0029,SIEMENS SYNGO VOLUME,32) | Volume Position (Patient) | DS | 3 |
| (0029,SIEMENS SYNGO VOLUME,37) | Volume Orientation (Patient) | DS | 9 |

| DICOM Tag | Name | VR | VM |
|--------------------------------------|--------------------------------|----|-----|
| (0029,SIEMENS SYNGO VOLUME,40) | Resampling Flag | CS | 1 |
| (0029,SIEMENS SYNGO VOLUME,42) | Normalization Flag | CS | 1 |
| (0029,SIEMENS SYNGO VOLUME,44) | SubVolume Sequence | SQ | 1-n |
| (0071,SIEMENS SYNGO REGISTRATION,20) | Registered Image Sequence | SQ | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,21) | Registration Is Validated Flag | CS | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,20) | Registered Image Sequence | SQ | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,21) | Registration Is Validated Flag | CS | 1 |
| (7FDF,SIEMENS SYNGO DATA PADDING,FC) | Pixel Data Leading Padding | OB | 1 |

Interpretation of the DICOM Tags from the above table:

(gggg, pp,ee) -> (gggg, ppee)

gggg - odd group number

pp - private creator identification code

ee - private element

9.3 Grayscale Image Consistency

The high resolution TFT display monitor option of syngo® MR comes with a DICOM Grayscale Standard Display Function (GSDF) compliant factory pre-setting. A typical working environment setup is assumed for ambient light.

Annex A: Index of Tables

| | |
|---|----|
| Table 1: Network Services..... | 2 |
| Table 2 - Media Services..... | 4 |
| Table 3 - Implementation Identifying Information | 5 |
| Table 4: Association Policies..... | 16 |
| Table 5: Asynchronous Nature as an Association Initiator..... | 16 |
| Table 6 - Presentation Context Table "Verification" | 17 |
| Table 7 - Presentation Context Table "Verification" | 17 |
| Table 8: Association Policies..... | 18 |
| Table 9: Asynchronous Nature as an Association Initiator..... | 18 |
| Table 10: Proposed Presentation Contexts for Storage..... | 19 |
| Table 11: DICOM Command Response Status Handling Behavior..... | 20 |
| Table 12: DICOM Command Communication Failure Behavior | 21 |
| Table 13: Storage C-STORE Response Status | 21 |
| Table 14: Association Policies..... | 23 |
| Table 15: Asynchronous Nature as an Association Initiator..... | 23 |
| Table 16: Proposed Presentation Contexts for Storage Commitment | 24 |
| Table 17: DICOM Command Response Status Handling Behavior..... | 24 |
| Table 18: DICOM Command Communication Failure Behavior | 24 |
| Table 19 - Presentation Context Table "Update Flag Information" | 25 |
| Table 20: DICOM Command Response Status Handling Behavior..... | 25 |
| Table 21: Association Policies..... | 26 |
| Table 22: Asynchronous Nature as an Association Initiator..... | 26 |
| Table 23: Proposed Presentation Contexts for Query | 27 |
| Table 24: Extended Negotiation as an SCU..... | 27 |
| Table 25: DICOM Command Response Status Handling Behavior..... | 28 |
| Table 26: DICOM Command Communication Failure Behavior | 28 |
| Table 27: Attributes supported for instance Query - SCU..... | 28 |
| Table 28: Proposed Presentation Contexts for Retrieve and Activity "MOVE SCU" | 29 |
| Table 29: DICOM Command Response Status Handling Behavior..... | 30 |
| Table 30: DICOM Command Communication Failure Behavior | 30 |
| Table 31: Association Policies..... | 31 |
| Table 32: Asynchronous Nature as an Association Initiator..... | 31 |
| Table 33: Proposed Presentation Contexts for Worklist | 32 |
| Table 34: Broad Query search keys..... | 32 |
| Table 35: Modality Worklist C-Find Return keys | 33 |
| Table 36: DICOM Command Response Status Handling Behavior..... | 36 |
| Table 37: DICOM Command Communication Failure Behavior | 36 |
| Table 38: Association Policies..... | 37 |
| Table 39: Asynchronous Nature as an Association Initiator..... | 37 |
| Table 40: Acceptable Presentation Contexts Activity "Create MPPS" | 38 |
| Table 41: MPPS N-CREATE Response Status Handling Behavior..... | 38 |
| Table 42: Acceptable Presentation Contexts Activity "Update MPPS" | 38 |
| Table 43: MPPS N-SET Response Status Handling Behavior | 39 |
| Table 44 - Presentation Context Table "Print Film" | 41 |
| Table 45 - Basic Film Session N-CREATE attributes | 42 |
| Table 46 - Basic Film Session Status Codes | 42 |
| Table 47 - Basic Film Box N-CREATE attributes | 43 |
| Table 48 - Basic Film Box Status Codes..... | 43 |
| Table 49 - Basic Grayscale Image Box N-SET attributes | 44 |
| Table 50 - Basic Grayscale Image Box Status Codes | 44 |
| Table 51 - Basic Color Image Box N-SET attributes..... | 44 |
| Table 52 - Basic Color Image Box Status Codes | 45 |
| Table 53 - Presentation LUT N-CREATE attribute..... | 45 |
| Table 54 - Presentation LUT Status Codes..... | 45 |
| Table 55 - Used Printer N-EVENT Report attributes..... | 46 |
| Table 56 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes | 46 |
| Table 57 - Used Print Job N-EVENT Report attributes | 46 |
| Table 58 - Presentation Context Table "Show Device Status" | 47 |

| | |
|--|----|
| Table 59 - Used Printer N-EVENT Report attributes..... | 48 |
| Table 60 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes | 48 |
| Table 61: Parameter List | 50 |
| Table 62: Implementation Class/Version Name - Media Interchange | 53 |
| Table 63: Media - Application Profiles and Real-World Activities | 53 |
| Table 64: Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K | 54 |
| Table 65: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles | 55 |
| Table 66: Single-Byte Character Sets without Code Extension | 56 |
| Table 67: Single-Byte Characters Sets with Code Extension | 56 |
| Table 68: Multi-Byte Character Sets without Code Extension | 57 |
| Table 69: Multi-Byte Character Sets with Code Extension | 57 |
| Table 70: Application Level Confidentiality Profile Attributes (standard tags)..... | 59 |
| Table 71: Application Level Confidentiality Profile Attributes (private tags) | 64 |
| Table 72: Private Data Element Dictionary | 83 |

Annex B: Table of Figures

| | |
|---|----|
| Figure 1: syngo® MR DICOM Data Flow diagram | 11 |
| Figure 2: Sequence Diagram for Real World Activities - System Configuration | 14 |
| Figure 3: Sequence Diagram for Real World Activities -Acquisition workflow | 15 |
| Figure 4: Media Interchange Application Data Flow Diagram | 51 |
| Figure 5: Sequence diagram – Media creation | 52 |

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 84-0
siemens.com/healthcare

HRT H00D05162002530837 | © Siemens Healthcare GmbH, 07. 2017

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.