

The Siemens logo is displayed in a bold, teal, sans-serif font within a white rectangular box in the top left corner of the advertisement.A Siemens MAGNETOM RT Pro edition MRI scanner is shown in a clinical setting. The machine is white and blue, with a large circular bore. A patient bed is partially visible, extending from the bore. The background is a light blue wall.

MAGNETOM RT Pro edition

The time is now for MR in RT.

usa.siemens.com/mri-in-rt

From 2006 to 2013, the utilization of MR images in radiotherapy treatment planning has **increased from 6% to 21% of all cases**¹. This is due to the fact that MRI's excellent soft-tissue contrast, exact information on the tumor extent, and information about functional parameters can add substantial clinical value to the therapy process.

The MAGNETOM RT Pro edition is a comprehensive package including hardware and software tailored to the needs of radiation oncology.

Easily scan patients in the treatment position



MRI Scanner

The top-of-the-line 70 cm Open Bore systems MAGNETOM Aera 1.5T and MAGNETOM Skyra 3T are the sound basis for MR in RT. Excellent patient access, a large Field of View, and trendsetting imaging applications enabled by the core technologies Tim 4G and Dot make the difference. The optional Tim Dockable Table improves your process efficiency by enabling the preparation of a patient outside the MR exam room.



Coils

The MAGNETOM RT Pro edition is equipped with a set of additional Flex coils and a Body coil with an extended cable for more flexible positioning of the coils. This ensures excellent signal-to-noise ratio for brain, head & neck, and body (pelvic) imaging when immobilization devices are in place.



RT Positioning Packages

To help ensure accurate and reproducible patient positioning, Siemens has partnered with CIVCO Medical Solutions and Orfit Industries, leading vendors for RT patient positioning solutions. Patients can be immobilized with thermoplastic masks to assure that imaging is performed in exactly the same position for simulation and treatment.




MR Compatible Laser Bridge

The optional LAP DORADOnova MR3T Laser Bridge is an external laser system supporting virtual simulation at the MR with

red or green laser lights. For high accuracy the marked position is shifted to the center of the scanner when using the RT Dot Engine, so there is no need to use the laser of the MR scanner in addition.

Consistently rely on intuitive workflows and spatial integrity




RT Dot Engine

Designed to enable consistent data acquisition and geometric integrity, the RT Dot Engine supports the user in the acquisition of suitable RT planning

images, which can be further processed in external RT applications. It contains features like automatic distortion correction, laser quality assurance, and automatic axial image reconstruction.

Instructions to perform quality assurance with the commercially available ACR MRI accreditation phantom can be downloaded at usa.siemens.com/mr-rt-qa.




Sequences & Applications

Siemens offers a variety of optional imaging applications which are tailored to address some of the main challenges in

radiation oncology. Techniques like BLADE, Advanced WARP³, Siemens-unique FREEZEit and RESOLVE help to maintain geometric integrity of the acquired data and to deal with typical sources of image artifacts.

Effectively establish multi-modality image viewing and precise contouring



Post-processing

syngo.via RT Image Suite is a dedicated RT software solution enabling advanced multi-modality image viewing, registration and target contouring. For

monitoring and therapy control, *syngo.MR OncoCare* is the tool to visualize functional changes in order to evaluate treatment response.

¹ IMV 2014 Radiation Therapy Market Summary Report

² Data on file. Results may vary.

³ The MRI restrictions (if any) of the metal implant must be considered prior to patient undergoing MRI exam. MR imaging of patients with metallic implants brings specific risks. However, certain implants are approved by the governing regulatory bodies to be MR conditionally safe. For such implants, the previously mentioned warning may not be applicable. Please contact the implant manufacturer for the specific conditional information. The conditions for MR safety are the responsibility of the implant manufacturer, not of Siemens.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Local Contact Information

Siemens Medical Solutions USA, Inc.
40 Liberty Boulevard
Malvern, PA 19355-9998
USA
Phone: +1-888-826-9702
usa.siemens.com/healthcare