



MR-integrated workflows **in Radiation Therapy** **for MAGNETOM Systems**

Thank you!

We would like to express our sincere gratitude to the following senior experts and their respective institutions who made this book possible.



Yue Cao, Ph.D.
University of Michigan,
Ann Arbor, USA




Gary Liney, Ph.D. and Robba Rai, MHLthSc
Liverpool and Macarthur Cancer Therapy Centre,
Ingham Institute for Applied Medical Research,
Sydney, Australia



Eric Paulson, Ph.D.
Medical College of
Wisconsin, Milwaukee,
USA



**Cynthia Ménard, M.D., FRCPC and
David Roberge, M.D., FRCPC**
Centre hospitalier de l'Université de Montréal,
Canada



Maja Sohlín, Ph.D.
Sahlgrenska University
Hospital, Gothenburg,
Sweden

Foreword

Radiation Oncology is experiencing growth in the use of magnetic resonance imaging in treatment planning and response assessment. While it can seem daunting to start integrating this technology in Radiation Therapy departments, this trend is inevitable and has been embraced by physicists and physicians alike. Clinical teams of several institutions have created their own solutions to optimize images, in an effort to add valuable clinical insights in treatment planning. This booklet is a compilation of study protocols* and practical tips for several body regions provided by these experts.

This compilation is being shared to foster clinical knowledge among medical professionals. You can download the .edx and .exar1 files at usa.siemens.com/magnetom-world-mr-in-rt

We would like to thank all the contributors who have worked tirelessly and generously shared their knowledge and expertise.

We look forward to hearing your feedback and suggestions so that we at Siemens Healthineers can continuously improve and do our part in the care of our patients.

Note: MAGNETOM systems are indicated for use as magnetic resonance diagnostic devices (MRDD) that produce transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that display the internal structure and/or function of the head, body or extremities. These images and/or spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician, yield information that may assist in diagnosis. The study protocols and another information contained in this brochure were developed entirely by the medical professionals identified at the beginning of each section. These medical professionals did not receive financial support or any other kind of assistance from Siemens Healthineers.

Patient results may vary depending on many factors including system, software version, options, coils, technique and patient condition. Siemens Healthineers is not responsible or liable for results.

Discover ways to optimize the use of your MR simulation applications.

Explore the latest clinical cases—from experts for experts!

*Results may vary depending on multiple factors like system, software version, options, coils and technique besides patient condition. Siemens Healthineers will not be liable for any damages."

Contents

Brain

- 04 Glioma 1**
Liverpool and Macarthur Cancer Therapy Centre,
Ingham Institute for Applied Medical Research,
Sydney, Australia
- 06 Glioma 2**
Sahlgrenska University Hospital,
Gothenburg, Sweden
- 08 Small Targets for Stereotactic Treatment**
Sahlgrenska University Hospital,
Gothenburg, Sweden

Head & Neck

- 10 Head & Neck 1**
Liverpool and Macarthur Cancer Therapy Centre,
Ingham Institute for Applied Medical Research,
Sydney, Australia
- 12 Head & Neck 2**
Medical College of Wisconsin,
Milwaukee, WI, USA
- 14 Head & Neck 3**
Sahlgrenska University Hospital,
Gothenburg, Sweden
- 16 Head & Neck 4**
University of Michigan,
Ann Arbor, MI, USA

Thorax

- 18 Esophagus**
Medical College of Wisconsin,
Milwaukee, WI, USA
- 20 Lung**
Liverpool and Macarthur Cancer Therapy Centre,
Ingham Institute for Applied Medical Research,
Sydney, Australia

Abdomen

- 22 Pancreas**
Medical College of Wisconsin,
Milwaukee, WI, USA
- 24 Cholangiocarcinoma**
Medical College of Wisconsin,
Milwaukee, WI, USA
- 26 Liver**
Centre hospitalier de l'Université de Montréal,
Canada

Pelvis

- 30 Prostate 1**
Centre hospitalier de l'Université de Montréal,
Canada
- 34 Prostate 2**
Sahlgrenska University Hospital,
Gothenburg, Sweden
- 38 Prostate 3**
Medical College of Wisconsin,
Milwaukee, WI, USA
- 40 Cervix – Brachytherapy**
Medical College of Wisconsin,
Milwaukee, WI, USA

Extremities

- 42 Musculoskeletal**
Liverpool and Macarthur Cancer Therapy Centre,
Ingham Institute for Applied Medical Research,
Sydney, Australia

Glioma 1

Robba Rai; Gary Liney
Liverpool and Macarthur Cancer Therapy Centre, Ingham Institute for Applied Medical Research,
Sydney, Australia

Overview

This protocol uses a 20-channel head & neck coil without thermoplastic mask immobilisation to maintain homogenous and high image quality acquired with a site specific volume coil. Care should be taken when positioning the patient to ensure they are straight and head tilt matches CT.

Patient

Patient Presentation	<p>High Grade Glioma: Considered in all patients, especially if interval from Sim to diagnostic imaging >3–4 weeks</p> <p>Low Grade Glioma: Considered if no post-operative MRI available, poor quality or non-contemporaneous diagnostic scan</p>
Treatment Prescribed	<p>Varies on grade and location of tumor:</p> <ul style="list-style-type: none">• 3D Conformal• IMRT (intensity-modulated radiation therapy)• VMAT (volumetric modulated arc therapy)• Tomotherapy

Imaging Study

Scanner	MAGNETOM Skyra 3T
Equipment Required	<ul style="list-style-type: none">• Head/Neck 20 coil• Immobilisation sponges
Patient Prep Needed	22-gauge intravenous cannula required prior to imaging
Setup and Landmark	<ul style="list-style-type: none">• Patient positioned head first supine.• Straighten patient through mid forehead, mid nose and mid chin using external lasers.• Use immobilisation sponges on the sides of patient’s ears and top of head to prevent movement during imaging.• Landmark over the lower orbital margin (OM) using the bore lasers.



1 Brain setup using 20-channel head and neck coil.

Sequences


Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:15 min		
T1w MPRAGE 3D sagittal	5:21 min	Anatomical delineation	
T2w FLAIR 2D	4:26 min	Post-operative changes Target delineation	Use 2 concatenations to minimize cross talk artifacts and maintain grey-white matter contrast
T1w MPRAGE 3D sagittal + Gad	5:21 min	Tumor delineation Residual enhancement post-op	<p>Injection Rate Hand bolus</p> <p>Gad Vol 0.1 mmol/kg (Gadobutrol)</p> <p>Saline Flush 10 mL before and after Gad</p>

Special considerations

Tips and Tricks	<ul style="list-style-type: none">• When positioning patient copy chin to chest measurements from CT Sim to match head tilt position to minimize variation between CT and MRI.• Use the same slice position from reconstructed transverse T1w MPRAGE and transverse T2w FLAIR to minimize errors between CT-MRI registration.
Preparing Images for Planning	<ul style="list-style-type: none">• Reconstruct the T1w MPRAGE pre and post-contrast in the transverse plane only for planning.
RTP Requirements	<ul style="list-style-type: none">• Use 3D distortion correction where applicable. If unavailable use 2D.• Maintain a receiver bandwidth with a ≤ 1 pixel fat-water shift (≥ 400 Hz/Px at 3T).• 0 mm interslice gap to eliminate interpolation in radiotherapy treatment planning (RTP) system.• 2 mm slices should be used for post-operative and high/low grade Glioma cases.• 1 mm slices should be used for stereotactic patients.

Example Images



2 T2w FLAIR showing gross tumor volume (GTV), clinical target volume (CTV), and planning target volume (PTV) margins (2A), gadolinium-enhanced T1w MPRAGE demonstrating peripheral enhancement of post surgical margins with clear definition of brain structures (2B), co-registration of non-contrast planning CT and post-contrast T1w MPRAGE, showing accurate registration of bone and soft tissue anatomy (2C).

Glioma 2

Maja Sohlín, Ph.D.; Christian Gustafsson, Ph.D.
Sahlgrenska University Hospital, Gothenburg, Sweden
Reproduced with permission from the “Method book for the use of MRI in radiotherapy” (Version 3) of the Swedish Vinova project.

Patient

Treatment Prescribed	Varies on grade and location of tumor.
----------------------	--

Imaging Study

Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">Two Flex 4 Large coilsVelcro fastenersSand bags

Setup and Landmark	<ul style="list-style-type: none">Patient positioned head first supine.Line things up based on the markings previously made on the immobilizing mask, if such exist. To ensure that the patient is positioned as straight as possible in the mask, it is particularly important to check that the sagittal laser follows the previously made marking, or is centered on the patient.The coils are inserted and clamped in place under the head area of the flat table top so they also cover the rear of the head. They are held together in the front with Velcro fasteners. The spine coil elements are not used with this setup.The coils are supported with sand bags on the sides to get closer to the head and prevent them from pressing on the patient’s nose and face.
--------------------	--




Sequences

Sequence	Scan time	Utility for RT	Characteristics
Axial T1w	5:25 min	Matching to CT plus anatomy before contrast	
Axial T2w	5:47 min	Delineating of target and risk organs	
Axial T2w FLAIR	5:02 min	Outlining of vasogenic edema/ infiltrative glioma (bright)	Potential artifacts from CSF pulsation
Axial T1w with contrast	5:25 min	Outlining of areas with defective blood brain barrier and neovascularization (bright)	Post-operative blood products – compare with pre-contrast T1w

Special Considerations

Tips and Tricks	<ul style="list-style-type: none">If possible, center the target volume in the isocenter to maximize field homogeneity and efficiency of any fat suppression and to minimize distortion caused by non-linear gradients.The T1-weighted image is used for matching to the CT and must therefore cover as large an area of the head as possible while maintaining correct registration.Other images must cover the target and any risk organs to be delineated.
-----------------	---

Example Images



2A T1 TSE tra (4:54 min),
2B T1 TSE Gd tra, (4:54 min),
2C T2 TSE tra (4:43 min),
2D T2 TIRM FLAIR tra (5:45 min).

Small targets for stereotactic treatment

Maja Sohlin, Ph.D.; Christian Gustafsson, Ph.D.
Sahlgrenska University Hospital, Gothenburg, Sweden

Reproduced with permission from the "Method book for the use of MRI in radiotherapy" (Version 3) of the Swedish Vinova project.

Patient

Treatment Prescribed	Stereotactic treatment is used for small target areas that often require high-resolution MR images with thin slices for target delineating. This means that it may be advantageous to use 3D sequences for stereotactic patients.
----------------------	---

Imaging Study

Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">Two Flex 4 Large coilsVelcro fastenersSand bags
Setup and Landmark	<ul style="list-style-type: none">Patient positioned head first supine.Line things up based on the markings previously made on the immobilizing mask, if such exist. To ensure that the patient is positioned as straight as possible in the mask, it is particularly important to check that the sagittal laser follows the previously made marking, or is centered on the patient.A holder was specially fabricated for the stereotactic frame (CIVCO trUpoint ARCH™ SRS/SRT System, CIVCO Medical Solutions, Kalona, IA, USA). To get the holder in a stable position, it was made to be fastened below the protruding head section of the flat table top. Because of this, the patient must be positioned well down on the examination table, which could potentially be problematic for unusually tall people. The advantage of this setup is that the spine coil elements can be used to cover the rear part of the cranium.Two Flex 4 Large coils are wound around the stereotactic frame and are held in place with sand bags on the sides.



1 Coil positioning for MAGNETOM Aera 1.5T system in Gothenburg with two Flex 4 Large coils covering the front part of the head. The rear part of the head is covered by the spine elements in the table.


Sequences

Sequence	Scan time	Utility for RT	Characteristics
Axial T1w	5:25 min	Matching to CT plus anatomy before contrast	
Axial T2w	5:47 min	Delineating of target and risk organs	
Axial T2w FLAIR	5:02 min	Outlining of vasogenic edema/ infiltrative glioma (bright)	Potential artifacts from CSF pulsation
Axial T1w with contrast	5:25 min	Outlining of areas with defective blood brain barrier and neovascularization (bright)	Post-operative blood products – compare with pre-contrast T1w

Special Considerations

Tips and Tricks	<ul style="list-style-type: none">If possible, center the target volume in the isocenter to maximize field homogeneity and efficiency of any fat suppression and to minimize distortion caused by non-linear gradients.The T1-weighted image is used for matching to the CT and must therefore cover as large an area of the head as possible while maintaining correct registration.Other images must cover the target and any risk organs to be delineated.
-----------------	---

Example Images



2 (2A) T1 TSE tra (5:18 min), (2B) T1 Gd tra (5:18 min), (2C) T2 TSE tra (6:20 min), (2D) T2 TIRM FLAIR tra (5:22 min).

Head & Neck 1

Robba Rai; Gary Liney
Liverpool and Macarthur Cancer Therapy Centre, Ingham Institute for Applied Medical Research,
Sydney, Australia

Overview

Protocol is designed for imaging in thermoplastic immobilisation mask. Coil arrangement has been developed to ensure homogenous signal throughout imaging volume whilst minimizing geometric distortions.

Patient

Patient Presentation	<p>Mucosal Primary:</p> <ul style="list-style-type: none">• All nasopharyngeal primaries• All other mucosal primaries where extent of disease uncertain after endoscopy/PET <p>Non-mucosal Primary:</p> <ul style="list-style-type: none">• Consider if primary or nodal disease is close to critical structures.
Treatment Prescribed	IMRT or Tomotherapy
Imaging Study	
Scanner	MAGNETOM Skyra 3T
Equipment Required	<ul style="list-style-type: none">• Two Body 18 surface coils• RF coil bridges• Flat table overlay• Velcro straps• Patient-specific immobilization device/s
Patient prep needed	20-gauge intravenous cannula required prior to imaging
Setup and Landmark	<ul style="list-style-type: none">• Adjust the flat table overlay by moving it down the MRI bed to ensure the top of the head of the thermoplastic mask is covered by the superior coil element of the integrated spine coil.• Index patient-specific immobilization devices to flat table overlay.• Patient to be positioned as per CT Simulation. Head first supine.• Straighten the patient according to the anatomical landmarks with the external laser system.• Secure coil bridges to table and use Velcro straps to stabilize the Body 18 surface coil to the bridges. The coil should cover the vertex of the skull to the suprasternal notch.• Lower the bridges so that it is close to the patient’s mask but not touching.• Place the second Body 18 surface coil over the patient’s torso.• Use longer Velcro straps to attach the two coils together, ensuring that there is minimal separation between the coils to minimize signal drop off in the images.• Landmark over the chin.



1 Head & Neck planning MRI setup using two 18-channel surface coils.

Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Sequence	Utility for RT	Characteristics
Localizer	0:15 min		
T2 TSE Fast Dixon	3:17 min	Primary and nodal delineation (particularly where contrast not given)	Plan with in-phase and water images
T1 TSE Fast Dixon	3:05 min	Primary and nodal delineation (particularly where contrast not given)	In-phase images used for bone registration in RTP system
T1 TSE Fast Dixon + Gad	3:05 min	Primary and nodal GTV delineation	<p>Injection Rate 2 mL/s</p> <p>Gad Vol 0.1 mmol/kg (Gadobutrol)</p> <p>Saline Flush 20 mL before and after Gad</p>

Special Considerations

Tips and Tricks	<ul style="list-style-type: none">• Dixon provides more robust and homogenous fat suppression as well as minimizing geometric distortion in the head and neck and should be considered for Head & Neck protocols.• Adjust the coil bridges so the patient head is equidistant to the spine and surface coil. This will help create a more homogenous image.• Landmark over the patient’s chin using the in bore lasers before securing any coils in place as the coils will obstruct the view of the patient’s external anatomy.
RTP Requirements	<ul style="list-style-type: none">• Use 3D distortion correction where applicable. If unavailable use 2D.• Maintain a receiver bandwidth with a ≤ 1 pixel fat-water shift ($\geq 400\text{Hz/Px}$ at 3T).• 0 mm interslice gap to eliminate interpolation in radiotherapy treatment planning (RTP) system.• All sequences should have the same centre slice to maintain consistency with CT-MRI rigid registration in RTP system.• 2–3 mm slice thickness to encompass entire neck in an ideal acquisition time to maintain image quality.

Example Images

2A

2B

2

T2 Dixon in-phase (2A), T2 Dixon water demonstrating robust fat suppression and delineation of tumor volume (2B), post-contrast T1 Dixon water showing enhancement of tumor and surrounding vasculature (2C), co-registration of planning CT and planning MRI (2D).

2C

2D

Head & Neck 2

Jared Robbins, M.D.; Eric Paulson, Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Squamous cell carcinoma
Treatment Prescribed	Volumetric Modulated Arc Therapy

Imaging Study

Scanner	MAGNETOM Verio 3T
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coils4-channel extremity coilSpine array coilThree Nylon strapsFlat table overlayHead restPatient-specific immobilization device
Patient Prep Needed	Change into gown and robe.
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine array.Index patient-specific immobilization devices to flat table overlay.Setup patient head-first supine, using external lasers to align and straighten.Wrap two 6-channel Body Matrix coils around patient’s head and secure with two Nylon straps.Confirm that the coils overlap anteriorly, are positioned as inferiorly as possible, and are symmetric.Position 4-channel extremity coil over neck and secure with Nylon strap Landmark over patient chin.



1 Immobilization device and RF coil positioning.

Sequences


Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT
Localizer	0:19 min	
Ax STIR TSE	7:00 min	Delineation of gross tumor volume and associated edema
Ax T1 TSE	4:20 min	Organs at risk delineation
Ax DWI RESOLVE	5:30 min	Localization of gross tumor volume
Ax fat-suppressed T1 TSE + Gad	5:30 min	Delineation of gross tumor volume Injection Rate: Hand bolus Gad Vol: 0.1 mmol/kg (Gadobenic acid)

Special Considerations

Acquisition	<ul style="list-style-type: none">Contiguous 3 mm slices with 1 mm in-plane resolution.Readout bandwidths adjusted to ensure WFS <1 pixel.
Preparing Images for Treatment Planning	<ul style="list-style-type: none">Apply 3D distortion correction to all images.Apply image standardization to all non-diffusion images.

Example Images




2 Multiparametric images used for target and OAR delineation: STIR (2A), T1w (2B), fat-suppressed T1+Gad (2C), RESOLVE DWI b=800 s/mm² (2D), ADC, apparent diffusion coefficient map (2E).

Head & Neck 3

Maja Sohlin, Ph.D.; Christian Gustafsson, Ph.D.
Sahlgrenska University Hospital, Gothenburg, Sweden
Reproduced with permission from the "Method book for the use of MRI in radiotherapy" (Version 3) of the Swedish Vinova project.

Imaging Study

Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">• Body 18 Long coil• Flex 4 Small coil• Coil holders• Flat table overlay• Patient-specific immobilization device
Setup and Landmark	<ul style="list-style-type: none">• Position flat table overlay over spine array.• Line things up based on the markings previously made on the immobilizing mask, if such exist. To ensure that the patient is positioned as straight as possible in the mask, it is particularly important to check that the sagittal laser follows the previously made marking, or is centered on the patient.• Coverage area see Figure 1.



1 Coverage areas (sag, cor) for head and neck. The large FOV used gives the skin contour of the patient and is used for better control of the matching.



2 Coil positioning for Head & Neck with Body 18 Long. There is also a Flex 4 Small coil positioned posteriorly (under the table top) to increase the SNR from the rear skull parts.


Sequences

Sequence	Scan time	Utility for RT	Characteristics
Axial T2 STIR	6:41 min	Differential edema (bright)	Potential swallowing artifacts, flow artifacts
Axial T1	6:36 min	Delineating of nerves and teeth	
Axial ADC		Plotting of hypercellularity (dark)	
Axial fat-suppressed post-contrast T1	6:36 min	Plotting of faulty, leaky tissue (bright)	Potential swallowing artifacts, flow artifacts

Special Considerations

If the patient has dental fillings¹, these will appear as streak artifacts in CT. This problem is often avoided in MRI as the artifacts instead become small, local, signal-poor parts. The size of the artifacts in CT depend on what type of material the dental filling is made out of. The same applies to MRI, but the artifact there can instead take the form of a geometric distortion, signal loss and/or signal shift. With such severe geometric distortion, great care should be taken when defining the target and risk organs in or near the artifact.

Example Images



3 (3A) T1 TSE tra 3 mm (6:36 min), (3B) T1 TSE tra 3 mm FatSat (6:36 min), (3C) T1 TSE tra 3 mm GD (6:36 min), (3D) T2 TSE tra 3 mm (6:41 min). ALL run with Neck shim and WARP.

¹The MRI restrictions (if any) of the metal implant must be considered prior to patient undergoing MRI exam. MR imaging of patients with metallic implants brings specific risks. However, certain implants are approved by the governing regulatory bodies to be MR conditionally safe. For such implants, the previously mentioned warning may not be applicable. Please contact the implant manufacturer for the specific conditional information. The conditions for MR safety are the responsibility of the implant manufacturer, not of Siemens Healthineers.

Head & Neck 4

Yue Cao, Ph.D.; James Balter, Ph.D.
University of Michigan, Ann Arbor, MI, USA

Patient

Patient Presentation A patient is positioned in the treatment configuration using five-point face mask and bit bar as physician indication.

Treatment Prescribed Standard 35 fractions of 70 Gy. Protocol scans include local boosting.

Imaging Study

Scanner MAGNETOM Skyra 3T

- Equipment Required**
- Flat table top
 - Body 18 long and Spine 32 matrix coils
 - Patient-specific mask, neck rest and bite block/dental guards when indicated
 - Sand bags
 - Head/neck board
 - Arm pulls if/as needed

- Setup and Landmark**
- Flat table top is usually already in place, but place if needed.
 - Mask can be made in CT or in Zone 3 (sag laser on wall).
 - Patient zeroed at External Auditory Meatus, mask is positioned using indexed head/neck board using the exact bar.
 - Body 18 placed over mask, sand bags used near head region.
 - Scan range varies by disease site and diagnosis.




Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT
Localizer	0:13 min	
2D T1 TSE FS pre Gd	2:08 min	Assisting GTV delineation
2D DWI RESOLVE	5:00 min	Boost target delineation
2D T2 TSE fs	3:14 min	Assisting GTV delineation
2D T1 TSE post Gd	5:35 min	GTV definition

Example Images



Esophagus

Candice Johnstone, M.D.; Paul Knechtges, M.D.; Eric Paulson, Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Esophageal cancer
Treatment Prescribed	3D conformal radiation therapy

Imaging Study

Scanner	MAGNETOM Verio 3T
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coilsSpine array coilTwo adjustable coil bridgesTwo Nylon strapsFlat table overlayHeadphonesECG leads and Bluetooth transceiverPatient-specific immobilization device/s
Patient Prep Needed	<ul style="list-style-type: none">Change into gown and robe.Attach ECG leads and Bluetooth transceiver.Place 22-gauge intravenous cannula into antecubital vein.
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine array.Index patient-specific immobilization devices to flat table overlay.Setup patient head-first supine, using external lasers to align and straighten, and confirm pelvic crest is above S9 on spine array coil.Place RF coil bridges over thorax and conform to patient surface anatomy.Wrap two 6-channel Body Matrix coils circumferentially over RF coil bridges and secure with two Nylon straps.Place headphones on patient.Landmark over thorax.



1 Adjustable RF coil bridges and RF coil positioning.

Sequences


Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer – free breathing	0:30 min		
Localizer – breath-hold	0:17 min		
Ax Restore T2 TSE	6:00 min	Delineation of gross tumor volume	Cardiac or respiratory gating
Ax BH T1 Dixon CAIPI pre	0:12 min	Reference for gross tumor volume delineation	
Care Bolus			Injection Rate: 2 cc/sec Gad Vol: 0.1 mmol/kg (Gadobenic acid)
Ax BH T1 Dixon CAIPI + Gad arterial	0:12 min	Reference for gross tumor volume	
Ax BH T1 Dixon CAIPI + Gad venous	0:12 min	Reference for gross tumor volume delineation	
Ax BH T1 Dixon CAIPI + Gad equilibrium	0:12 min	Reference for gross tumor volume delineation	

Special Considerations

Acquisition	<ul style="list-style-type: none">Respiratory gating performed for tumors close to gastroesophageal junction (GEJ).Cardiac gating performed for tumors close to heart.Contrast injection performed using bolus tracking technique.Readout bandwidths adjusted to ensure WFS <1 pixel.
Preparing Images for Treatment Planning	<ul style="list-style-type: none">Apply 3D distortion correction to all images.Apply image standardization to non-DWI images.Perform local rigid registration over gross tumor volume to align reference MR image with planning CT images.

Example Images



2 Planning CT (2A) and cardiac and respiratory gated Restore T2 TSE (2B).

Lung

Gary Liney; Robba Rai
Liverpool and Macarthur Cancer Therapy Centre, Ingham Institute for Applied Medical Research,
Sydney, Australia

Overview


This is a free breathing procedure incorporating morphological and functional imaging protocols. Ultrashort echo time (UTE) imaging is incorporating for the purposes of MRI only planning and lung parenchyma assessment.

Patient


Patient Presentation	Treatment response research – free-breathing protocol
Treatment Prescribed	IMRT or 3D Conformal (Multifield)

Imaging Study

Scanner	MAGNETOM Skyra 3T
Equipment Required	<ul style="list-style-type: none">• Body 18, 18-channel surface coil• Spine 32, 32-channel spine coil• RF coil bridges• Flat table overlay• Velcro straps• Patient-specific immobilization device/s
Patient Prep Needed	20-gauge intravenous cannula required prior to imaging
Setup and Landmark	<ul style="list-style-type: none">• Index patient-specific immobilization devices to flat table overlay.• Patient to be positioned as per CT Simulation. Head first supine.• Straighten the patient to AP and rotation tattoos with the external laser system.• Secure coil bridges to table and strap 18-channel surface coil to the bridges.• Lower the bridges so that it is close to the patient’s torso but not deforming the external contours.• Ensure that the coil covers the top of lung apices to lower costophrenic angles.• Landmark over the cross hair of the coil (center of the thorax).



2 Free-breathing T2 HASTE using phase navigation for respiratory triggering (2A), free-breathing T1 StarVIBE (2B), free-breathing coronal T1 StarVIBE perfusion used for tumor perfusion analysis (2C), free-breathing diffusion-weighted imaging, b-value 750 s/mm² (2D) and ADC map (2E).



Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:09 min		
T2_HASTE_transverse	5:00 min	Target delineation	Respiratory Navigation
T1 StarVIBE_transverse	3:40 min	Node and bone delineation	
T2_TRUFI_coronal	0:19 min	Tumor motion	
T2_TRUFI_sagittal	0:13 min	Tumor motion	
ZOOMit_EPI_transverse	3:18 min	Target delineation	
Radial VIBE coronal	9:04 min	Lung parenchyma delineation and MR-only planning	
T1 StarVIBE coronal 2 degree flip angle	0:14 min	T1 Mapping	
T1 StarVIBE coronal 15 degree flip angle	0:14 min	T1 Mapping	
T1 StarVIBE coronal dynamic + Gad	5:42 min	Tumor perfusion	Injection Rate: 4 mL/s Gad Vol: 0.05 mmol/kg (Gadobutrol) Saline Flush: 20 mL inject after 3 measurements
T1 StarVIBE SPAIR post Gad	4:13 min	Target delineation and nodal enhancement	Injection Rate: 2 mL/s Gad Vol: 0.05 mmol/kg (Gadobutrol) Saline Flush: 20 mL scan after 20 sec. post Gad injection

Special Considerations

Tips & Tricks	<ul style="list-style-type: none">• Use phase navigation for T2 HASTE for optimal image quality.• Place cardiac shim box over the heart for TrueFISP imaging to minimize off-resonance and flow related artifacts.
RTP Requirements	<ul style="list-style-type: none">• Use 3D distortion correction where applicable. If unavailable use 2D.• Maintain a receiver bandwidth with a ≤1 pixel fat-water shift (≥400 Hz/Px at 3T).• 0 mm interslice gap to eliminate interpolation in radiotherapy treatment planning (RTP) system.• All sequences should have the same centre slice to maintain consistency with CT-MRI rigid registration in RTP system.• 3–4 mm slice thickness to encompass entire thorax in an ideal acquisition time to maintain image quality.

Pancreas

Beth Erickson, M.D.; Paul Knechtges, M.D.; Eric Paulson, Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Locally advanced unresectable pancreatic adenocarcinoma
Treatment Prescribed	Respiratory-gated intensity modulated radiation therapy

Imaging Study

Scanner	MAGNETOM Verio 3T
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coilsSpine array coilTwo adjustable coil bridgesTwo Nylon strapsFlat table overlayHeadphonesPatient-specific immobilization device/s
Patient Prep Needed	<ul style="list-style-type: none">Change into gown and robe.Place 22-gauge intravenous cannula into antecubital vein.
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine array.Index patient-specific immobilization devices to flat table overlay.Setup patient head-first supine, using external lasers to align and straighten, and confirm pelvic crest is above S9 on spine array coil.Place RF coil bridges over abdomen and conform to patient surface anatomy.Wrap two 6-channel Body Matrix coils circumferentially over RF coil bridges and secure with two Nylon straps.Place headphones on patient.Landmark over diaphragm.




1 Adjustable RF coil bridges and RF coil positioning.

Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer – free breathing	0:30 min		
Localizer – breath-hold	0:17 min		
Ax PACE T2 HASTE	2:30 min	Organs at risk delineation	
Ax PACE fat-suppressed T2 HASTE	2:30 min	Localization of gross tumor volume	
Ax PACE DWI	5:00 min	Localization of gross tumor volume	
Ax BH T1 Dixon CAIPI pre	0:12 min	Reference for gross tumor volume delineation	
Care Bolus			Injection Rate: 2 cc/sec Gad Vol: 0.1 mmol/kg (Gadobenic acid)
Ax BH T1 Dixon CAIPI + Gad arterial	0:12 min	Delineation of gross tumor volume	
Ax BH T1 Dixon CAIPI + Gad venous	0:12 min	Reference for gross tumor volume delineation	
Ax BH T1 Dixon CAIPI + Gad equilibrium	0:12 min	Reference for gross tumor volume delineation	

2 Multiparametric images used for target and OAR localization and delineation. Planning CT (2A), IV contrast CT (2B), T2 HASTE (2C), fat-suppressed T2 HASTE (2D), pre-contrast Dixon T1 (2E), arterial phase Dixon T1 + Gad (2F), venous phase Dixon T1 + Gad (2G), apparent diffusion coefficient (ADC) map (2H).



Cholangiocarcinoma

Jared Robbins, M.D.; Paul Knechtges, M.D.; Eric Paulson, Ph.D.; X. Allen Li, Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Intrahepatic cholangiocarcinoma
Treatment Prescribed	Liver stereotactic body radiation therapy

Imaging Study

Scanner	MAGNETOM Verio 3T
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coilsSpine array coilTwo adjustable coil bridgesTwo Nylon strapsFlat table overlayHeadphonesPatient-specific immobilization device/s
Patient Prep Needed	<ul style="list-style-type: none">Change into gown and robe.Place 22-gauge intravenous cannula into antecubital vein.
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine array.Index patient-specific immobilization devices to flat table overlay.Setup patient head-first supine, using external lasers to align and straighten, and confirm pelvic crest is above S9 on spine array coil.Place RF coil bridges over abdomen and conform to patient surface anatomy.Wrap two 6-channel Body Matrix coils circumferentially over RF coil bridges and secure with two Nylon straps.Place headphones on patient.Landmark over diaphragm.



1 Adjustable RF coil bridges and RF coil positioning.

Sequences


Sequence	Scan time	Utility for RT	Characteristics
Localizer – free breathing	0:30 min		
Localizer – breath-hold	0:17 min		
Ax BH T1 Dixon CAIPI pre	0:12 min	Reference for gross tumor volume delineation; organs at risk delineation	
Care Bolus			Injection Rate: 2 cc/sec Agent Vol: 10 mL (Gadoxetate disodium)
Ax PACE fat-suppressed T2 TSE	4:00 min	Localization of gross tumor volume	
Ax PACE DWI	5:00 min	Localization of gross tumor volume	
Ax BH T1 Dixon CAIPI 20 min delay	0:12 min	Delineation of gross tumor volume	
Ax BH T1 Dixon CAIPI 30 min delay	0:12 min	Delineation of gross tumor volume	

Special Considerations

Acquisition	<ul style="list-style-type: none">Contrast injection performed using bolus tracking technique.Administer 1 mg Glucagon IV (1/2 at start of exam, 1/2 midway through exam) to suppress peristalsis.Respiratory gating and breath-holds performed at end expiration to match phase of gating window used for treatment delivery.Readout bandwidths adjusted to ensure WFS <1 pixel.
Preparing Images for Treatment Planning	<ul style="list-style-type: none">Apply 3D distortion correction to all images.Apply image standardization to non-DWI images.Identify a reference MR image and perform local MR-MR rigid registration over gross tumor volume to correct for inter-sequence motion.Perform local rigid registration over gross tumor volume to align reference MR image with planning CT images.

Example Images

2 Multiparametric images used for target and OAR localization and delineation. Planning CT (2A), IV contrast CT (2B), fat-suppressed T2 TSE (2C), 30-minute delayed Dixon T1 + contrast medium (2D), DWI b=500 s/mm² (2E), ADC map (2F).



Liver

Karim Boudam; Jean-Charles Côté; Marie-Pierre Campeau; Guila Delouya;
David Donath; David Roberge
Centre hospitalier de l'Université de Montréal, Canada

Patient

Patient Presentation	58-year-old man presenting with metastatic liver lesions from a primary colon cancer.
Treatment Prescribed	Stereotactic Body Radiation Therapy (SBRT) 45 Gy in 3 fractions

Imaging Study

Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">• Spine 32 surface coil• Body 18 surface coil• Wireless respiratory monitoring system
Patient Prep Needed	Immobilization device: Full-body vacuum cushion (BodyFIX, Elekta Instruments, Crawley, UK) + abdominal compression belt (Aktina Medical, Congers, NY, USA)
Setup and Landmark	The patients are provided the compression belt and immobilized in treatment position using the vacuum cushion. The positioning on the couch is carried out using the Siemens Healthineers MRI built in lasers.



1 Patient immobilization setup.

Imaging Protocol	All the sequences, except the 2D cine, are corrected for geometry distortion using the built in Siemens Healthineers 3D correction algorithm.
------------------	---

Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:48 min		
T1 VIBE tra 3D	0:15 min	Positioning of the 2D cine slab	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
2D cine without compression	1:00 min	Tumor motion assessment without compression	1.0 x 1.0 x 10 mm
2D cine with compression	1:00 min	Tumor motion assessment with compression	1.0 x 1.0 x 10 mm
T1 VIBE tra 3D pre-contrast sequence (Breath-held expiration)	0:15 min	Differentiation of hepatocytes from tumor (dark)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
T1 VIBE tra 3D pre-contrast sequence (Breath-held inspiration)	0:15 min	Differentiation of hepatocytes from tumor (dark)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
T2 FS TSE tra 2D (Triggered expiration)	5:00 min	Delineation of tumor (bright)	1.0 x 1.0 x 3.0 mm
T2 FS TSE tra 2D mBH	0:45 min	Optional sequence Delineation of tumor (bright)	2.0 x 2.0 x 3.0 mm
T1 VIBE tra 3D post-contrast: Gadobutrol – 30 seconds delay arterial phase (Breath-held expiration)	0:15 min	Delineation of hypervascular tumors (bright)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
T1 VIBE tra 3D post-contrast: Gadobutrol – venous phase (Breath-held inspiration)	0:15 min	Delineation of hypervascular tumors (bright)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
T1 VIBE tra 3D post-contrast: Gadobutrol – 5 min delay (Breath-held expiration)	0:15 min	Delineation of hypovascular tumors (dark)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm
T1 VIBE tra 3D post-contrast: Gadobutrol – 5 min delay (Breath-held inspiration)	0:15 min	Delineation of hypovascular tumors (dark)	1.7 x 1.7 x 4.0 mm Resample: 1.7 x 1.7 x 2.0 mm

Sequences


The MR imaging workflow starts with a crude localization sequence followed by a 3D T1 VIBE sequence to identify the target area and thus fine tune the positioning of the 2D cine sequences. Tumor motion assessment is first done without compression then, multiple 2D cine sequences are acquired with increasing abdominal pressure. The belt pressure is chosen to obtain a clinically satisfying reduction in tumor motion within the limits of patient comfort.

Following this iterative evaluation of belt pressure, T1 and T2 sequences are acquired with the selected pressure. The T1 sequences are acquired in breath-hold at expiration and inspiration with and without contrast (30 seconds and 5 minutes post injection). The T2 sequence is obtained at expiration using the triggering module. A T2 multiple breath-hold sequence is also acquired if the tumor is not visible on the T1 inspiration images.


Special Considerations

Preparing Images for Planning


Figure 2 describes the treatment planning registration workflow. A 4D CT scan is acquired in treatment position and the maximum expiration and inspiration phases are identified amongst the 10 binned phases. The MRI sequences are then anatomically registered using a manual 6D rigid transformation. The acquisition of sequences in expiration and inspiration allows us to minimize the differences between the CT and MR images.



Example Images



3 T1w axial fs breath-hold sequence acquired in inspiration (3A). 2D cine (2B).



4 T1w axial fs breath-hold sequences acquired in expiration and inspiration at 30 seconds (4A, B) and 5 min post-contrast injection (4C, D). T2w axial fat sat triggered at expiration (4E).

Prostate 1

.....

Cynthia Ménard¹; Jean-Charles Côté¹; Guila Delouya¹; Maroie Barkati¹; Karim Boudam¹; Carole Lambert¹; Gary Liney²
¹Centre hospitalier de l'Université de Montréal, Canada
²Liverpool and Macarthur Cancer Therapy Centre, Ingham Institute for Applied Medical Research, Sydney, Australia

Patient

Patient Presentation	67-year-old man presenting with intermediate-risk prostate cancer, clinical stage T1c, Gleason VII (3+4) with 7/18 biopsies positive, PSA 5.46.
Treatment Prescribed	SBRT 36.25 Gy in 5 fractions, one fraction per week.
Imaging Study	
Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">• Spine 32 surface coil• Body 18 long surface coil• Foot rest
Patient Prep Needed	<ul style="list-style-type: none">• Fiducial markers (4 in-house 1 x 3mm platinum cylinders) are implanted under transrectal ultrasound (TRUS) guidance at least 7 days prior to treatment planning.• Patients are instructed to empty their bladder before simulation. Bowel prep includes a fleet enema before imaging.• Buscopan 20 mg IM at 1ml/min, 10 to 15 minutes before MRI (unless contraindicated).
Setup and Landmark	Patients are positioned supine directly on the spine coil with foot rest.



Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:18 min		
T2 SPACE tra 3D	6:22 min	Prostate segmentation	1.0 x 1.0 x 1.0 mm
PD T2 TSE tra 2D	5:41 min	Marker-based registration and corresponding diagnostic T2w TSE image for tumor identification	0.8 x 0.8 x 3.0 mm
DWI tra b=50, 400, 800 s/mm², 2D	5:14 min	Tumor segmentation	1.8 x 1.8 x 4.0 mm ADC and calculated b-value of 1400
T1 SPACE cor 3D	5:12 min	Diagnostic image highlighting intraprostatic blood (assist tumor segmentation) as well as pelvic bones and lymph nodes	1.1 x 1.1 x 1.2 mm

- Table in ISO positioning mode
- 3D distortion correction applied for anatomical images (2D for diffusion)
- Bandwidth of at least 220 Hz/mm (1 mm fat-water separation at 1.5 Tesla)
- Include prostate, seminal vesicles and a part of bladder, rectum
- Save separately PD and T2 of PD T2 TSE TRA 2D
- Save an axial MPR of the T1 SPACE COR 3D

This imaging protocol strives to balance treatment planning requirements (segmentation and registration), while complying with PIRADS v2.0 to permit valid diagnostic reporting. For the purpose of PIRADS reporting, the T2 SPACE image can be reconstructed in sagittal and coronal planes.


The examination is designed to keep each individual acquisition under 6 minutes to minimize motion related artifacts, and the overall examination under 30 minutes. The imaging range is the prostate gland, with the exception of the final T1-weighted image, where the range is the entire pelvis.

A dual-echo (PD/T2) 2D acquisition is at the center of a rigid marker-based MRI-CT registration workflow. A volumetric T2w acquisition can be anatomically registered to the 2D T2w dataset if motion is observed between adjacent acquisitions. In that case, anatomic registration (T2 to T2) can be performed. The last two image sets can be used by clinicians to assist in tumor segmentation.

Special Considerations


Preparing Images for Planning

- Send images from MR to PACS.
- Load images (MR and CT) from PACS on fusion console (MIM Vista, MIM Software Inc., Cleveland, OH, USA).
- Register MR to CT (workflow below).
- Send registered images to contouring server.
- Contour volumes on contouring client (Somavision Varian Medical Systems, Palo Alto, CA, USA).
- Send CT and contour from contouring server to planning console.




1 Manual point-based registration (PD to CT) using 4 implanted markers (sup and inf) for a total of 8 points.
2 If required (due to motion) anatomic automated registration prostate-to-prostate using mutual information algorithm.

Example Images



3 Slice corresponding to implanted marker right mid-gland best visualized on PD image (3A).




4 Slice corresponding to tumor at right anterior apex best visualized on calculated b=1400 DWI (4C).

Prostate 2

Maja Sohlin, Ph.D.; Christian Gustafsson, Ph.D.
Sahlgrenska University Hospital, Gothenburg, Sweden
Reproduced with permission from the "Method book for the use of MRI in radiotherapy" (Version 3) of the Swedish Vinova project.


Imaging Study

Scanner	MAGNETOM Aera 1.5T
Equipment Required	<ul style="list-style-type: none">• Body 18 Long coil• Coil holders• Flat table overlay (Civco Medical Solutions, Kalona, IA, USA)• Knee support on an in-house moveable rail
Setup and Landmark	<ul style="list-style-type: none">• If the patient comes in for the MRI after undergoing the CT scan and creating a CT treatment plan, then there will already be body markings that were applied to the patient's hip during the CT visit. The position of these body markings are determined by a laser with preset definitions that are the same as the preset definitions for CT and the treatment unit. By matching this laser system with these markings, the patient can be positioned the same way in each radiotherapy fraction and you get a reproducible position for the patient. This avoids rotation of the pelvis, thereby minimizing the prostate's position deviation between different radiotherapy fractions.• When the patient comes in to MRI it is advisable for the same positioning procedure to be used. There is therefore usually a laser system in the MRI room. If the laser system does not align with the body markings, you can rotate the patient's hip by grasping it and rotating.• Once this is done, it is advisable to ask the patient to lift their pelvis straight up in the air slightly and then lower it again. Then check that the laser aligns with the body markings as shown below:



1 (1A) Patient's position requires correction, (1B) patient's position is correct.

- Because of the manual rotation performed with the hip, there is a risk of the buttocks will end up slightly asymmetrical. By asking the patient to lift their hips slightly, you avoid having the buttocks end up in different positions between the CT image and the MR image, which could make registration between these images more difficult. The separation between the buttocks can be used as a landmark to check whether any asymmetry has occurred in the position:



2 Separation of the buttocks on (2A) CT image, and (2B) MR image.

- Coverage areas (cor, sag, tra) for a transverse T2-weighted series: cover the entire prostate, approx. 3–4 slices, from above the bladder base down to the penis root.
- Coverage areas (cor, sag, tra) for a transverse T1-weighted series: use a large FOV to visualize the skin contour for better quality control of the registration.



3 Coil positioning with Body 18 Long coil. The coil holders prevent the coil from pressing on the skin contour. The knee support sits on a moveable rail created at the workshop. The knee support should be the same type as in CT to obtain the best reproducible patient positioning.


Sequences

Sequence	Scan time	Utility for RT	Characteristics
Sag T2		Delineating of rectum and bladder	
Axial T2	6:17 min	Delineating of prostate and extracapsular disease (dark)	Potential post-biopsy hemorrhage, difficult to see the marker
Axial fat suppressed T2		Delineating of extracapsular disease (dark), lymph nodes (bright)	Potential post-biopsy hemorrhage
Axial T1	3:38 min	Detection of post-biopsy haemorrhaging (bright), visualization of markers	
Axial Diffusion ADC		Delineating of tumor (dark)	

Special Considerations

- The use of contrast medium is governed by local procedures. The procedures for when the patient should drink and/or urinate are governed by local regulations, but the recommendation is to use the same method for CT, MRI and the radiotherapy treatment unit for each radiotherapy fraction.
- It is important to pay attention to what type of ink is used for the body markings as certain ink types have been shown to produce artifacts in the form of signal loss on the MR image. This is not acceptable when MRI is used for treatment planning, but in connection with matching the MR image to a CT image it is not as critical since only a small proportion of the skin has disappeared. After automatic matching, always check that the results are correct.
- For patients with MRI-compatible hip prostheses¹, the protocols can be adjusted so that metal artifacts are as small as possible. This is best done by applying any metal artefact reduction functions, such as WARP. If this is not available, a higher bandwidth and reduced echo time could limit the spread of metal artifacts. The artifacts are generally smaller on 1.5T compared to 3T. On CT, streak artifacts penetrate the prostate, while the prostate is intact on the MRI. Bear in mind that automatic matching with CT, if used, may be incorrect. Carefully check the matching manually.
- In order for a T2-weighted image series to be matched against CT data, you often rely on an image series where the gold markers are clearly visible and instead register these to the CT. This image series is often taken with some type of gradient echo-based imaging technique to show the position of the gold seeds in the prostate more clearly. We call this image series the differentiation series. The registration between the T2-weighted image series and the CT data then follows automatically since the T2-weighted image series and the differentiation series are taken in the same frame of reference. A disadvantage of this approach is that there is a risk that the patient and/or the patient’s anatomy has moved between the image capture of the T2-weighted image series and the image capture of the differentiation series. This can be addressed by directly visualizing the gold markers in the T2-weighted image series. It may be challenging to visualize the gold markers directly on the T2-weighted image series because the gold markers themselves are small, and the T2-weighted image series is often based on the turbo spine echo technique (TSE), which automatically minimizes and to some degree compensates the susceptibility artifacts from the gold markers. In TSE sequences, there is also smearing in the image as a result of the repeated 180° pulses it uses. On a MAGNETOM Aera 1.5T system, the following adjustments have been shown to visualize the gold seeds directly on the T2-weighted image series:

- Increase the spatial resolution of the image. Approx. 0.7 x 0.7 mm, depending on the size of the gold markers used.
- Reduce the slice thickness to reduce partial volume artifacts. Approx. 2.5 mm.
- Reduce echo train smearing by reducing the echo train length in the TSE sequence.
- Restrict the number of NEX to just 1 to reduce the smearing caused by averaging.
- Eliminate imaging acceleration to increase SNR.
- Restrict imaging to 1 package.
- Use movement correction, if available.
- These measures often lead to a reduction in SNR and thereby an increase in scan time to obtain acceptable image quality.



4 (4A) T1 VIBE tra 2 mm (3:38 min), (4B) T2 TSE tra 2 mm (6:17 min).

¹The MRI restrictions (if any) of the metal implant must be considered prior to patient undergoing MRI exam. MR imaging of patients with metallic implants brings specific risks. However, certain implants are approved by the governing regulatory bodies to be MR conditionally safe. For such implants, the previously mentioned warning may not be applicable. Please contact the implant manufacturer for the specific conditional information. The conditions for MR safety are the responsibility of the implant manufacturer, not of Siemens Healthineers.

Prostate 3

Colleen Lawton, M.D.; Eric Paulson, Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Prostate adenocarcinoma
Treatment Prescribed	Volumetric modulated arc therapy

Imaging Study

Scanner	MAGNETOM 3T Verio
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coilsSpine array coilTwo adjustable coil bridgesTwo Nylon strapsFlat table overlayPatient-specific immobilization device(s)
Patient Prep Needed	<ul style="list-style-type: none">Change into gown and robeBladder filling protocol
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine arrayIndex patient-specific immobilization devices to flat table overlaySetup patient head-first supine, using external lasers to align and straighten, and confirm pelvic crest is above S9 on spine array coilPlace RF coil bridges over pelvis and conform to patient surface anatomyWrap two 6-channel Body Matrix coils circumferentially over RF coil bridges and secure with two Nylon strapsLandmark over diaphragm



1 Adjustable RF coil bridges and RF coil positioning.

Sequences


Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:30		
Sag Restore T2 TSE	3:00	Delineation of gland and rectum / gland interface	
Ax Restore T2 TSE	5:00	Delineation of gland	
Cor Restore T2 TSE	3:00	Delineation of gland and localization of urethra	
Ax DWI	5:00	Localization of gross tumor volume	Acquired b-values: 200, 600, 1000 s/mm ² ; calculated b-value: 1400 s/mm ²

Special Considerations

Acquisition	<ul style="list-style-type: none">Readout bandwidths adjusted to ensure WFS < 1 pixel
Preparing Images for Treatment Planning	<ul style="list-style-type: none">Apply 3D distortion correction to all imagesApply image standardization to non-DWI imagesPerform local rigid registration over prostate gland to align reference MR image with planning CT images

Example Images



2 Multiparametric images used for target and OAR localization and delineation. Planning CT (2A), sagittal T2 TSE (2B), axial T2 TSE (2C), and coronal T2 TSE (2D).

Cervix—Brachytherapy

Beth Erickson M.D.; Jason Rownd MS; Eric Paulson Ph.D.
Medical College of Wisconsin, Milwaukee, WI, USA

Patient

Patient Presentation	Cervix cancer
Treatment Prescribed	High dose rate brachytherapy boost

Imaging Study

Scanner	MAGNETOM 3T Verio
Equipment Required	<ul style="list-style-type: none">Two 6-channel Body Matrix coilsSpine array coilTwo adjustable coil bridgesTwo Nylon strapsFlat table overlayMR-compatible physiological monitoring equipment
Patient Prep Needed	<ul style="list-style-type: none">Brachytherapy preparation (including anesthesia)
Setup and Landmark	<ul style="list-style-type: none">Position flat table overlay over spine arrayTransfer patient from brachytherapy suite to MR table, setting up patient foot-first supine and confirm pelvic crest is below S9 on spine array coilAttached physiological monitoring devicesPlace RF coil bridges over pelvis and conform to patient surface anatomyWrap two 6-channel Body Matrix coils circumferentially over RF coil bridges and secure with two Nylon strapsLandmark over diaphragm



1 Adjustable RF coil bridges and RF coil positioning.

Sequences

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:30		
3PL T2 HASTE	0:59	Confirmation of applicator positioning and absence of perforation	Aligned to applicator in three planes
Ax 3D T2 SPACE	12:00	Planning MR image, delineation of targets and OAR	1 x 1 x 1.5 mm ³ voxels
Ax T2 TSE	4:00	Delineation of gross tumor volume and high risk clinical target volume	Para-axial, parallel to ring
Sag T2 TSE	3:00	Delineation of gross tumor volume and high risk clinical target volume	Para-sagittal, parallel to tandem
Cor T2 TSE	3:00	Delineation of gross tumor volume and high risk clinical target volume	Para-coronal, parallel to tandem

Special Considerations

Acquisition	<ul style="list-style-type: none">Administer 1 mg Glucagon IV (1/2 at start of exam, 1/2 midway through exam) to suppress peristalsisReadout bandwidths adjusted to ensure WFS < 1 pixel
Preparing Images for Treatment Planning	<ul style="list-style-type: none">Apply 3D distortion correction to all imagesApply image standardization to all images

Example Images



2 Multiparametric images used for target and OAR localization and delineation. **2A:** Axial 3D T2 SPACE planning MR images displayed with sagittal and coronal reformats. A saline filled dummy catheter, inserted into the tandem lumen, is visible on the sagittal image. Registered 3D T2 SPACE (**2B**), para-sagittal 2D T2 TSE (**2C**), para-axial 2D T2 TSE (**2D**), and para-coronal 2D T2 TSE (**2E**) images. A saline filled dummy catheter, inserted into the ring lumen, is visible on the para-axial images. The saline filled dummy catheters assist in aligning applicator templates during brachytherapy treatment planning.

Musculoskeletal

Robba Rai; Gary Liney
Liverpool and Macarthur Cancer Therapy Centre, Ingham Institute for Applied Medical Research,
Sydney, Australia

Overview

Coil and protocol selection will vary depending on the anatomical site of interest and patient positioning. For areas greater than 30 cm in scan length, TimCT should be considered to maintain geometrical integrity of scans for accurate CT-MRI rigid registration in radiotherapy treatment planning system.

Patient


Patient Presentation	Any peripheral primary or secondary lesions <ul style="list-style-type: none">SarcomaOsteosarcomaMetastatic disease
Treatment Prescribed	Dependent on location and disease extent

Imaging Study

Scanner	MAGNETOM Skyra 3T
Equipment Required	<ul style="list-style-type: none">Flex 4 Large and Small (Site dependant)Body 18 surface coil (Site dependant)Flat table overlay (CIVCO)Velcro strapsPatient-specific immobilization device/s
Patient Prep Needed	22-gauge intravenous cannula required prior to imaging if gadolinium is requested
Setup and Landmark	<ul style="list-style-type: none">Index patient-specific immobilization devices to flat table overlayPatient to be positioned as per CT SimulationStraighten the patient's anatomy of interest using the rotation tattoos with the external laser systemDepending on the region of interest



1 MR Sim setup using vacuum bag for immobilisation over tricep Sarcoma. 18-channel surface coil used over centre over affected arm. Patient was offset to position arm in the isocentre of the MRI scanner.



2 MR Sim setup in prone position for gluteal sarcoma. 18-channel surface coil used over centre of affected area.



3 MR Sim setup for TimCT of femur using a combination of coils including 18-channel surface coils and two 4-channel flex coils.

Sequences

Small region of interest

Download the .exar1 file at usa.siemens.com/magnetom-world-mr-in-rt

Sequence	Scan time	Utility for RT	Characteristics
Localizer	0:09 min		
T1 TSE	2:04 min	Bone registration to CT	
T2 TSE Dixon	2:40 min	Delineation of gross tumor volume, nodes and associated edema	
T1 TSE Dixon + Gad	2:40 min	Delineation of gross tumor volume	Injection Rate: Hand bolus Gad Vol: 0.1 mmol/kg (Gadobutrol) Saline Flush: 10 ml before and after Gad

Long bones and larger regions of interest (TimCT)

Sequence	Scan time	Utility for RT
Localizer	0:09 min	
T1 FLASH TimCT	4:48 min	Bone registration to CT
T2 HASTE FS TimCT	6:43 min	Delineation of gross tumor volume and nodes
T2 HASTE TimCT	6:43 min	Delineation of gross tumor volume and associated edema

Special Considerations


Tips & Tricks

- Use MRI markers over reference tattoos to help with CT-MRI registration in radiotherapy treatment planning system.
- Use immobilisation methods such as sponges and tape to ensure patient does not move during image acquisition.
- Run two localizers for upper extremities e.g. arm with the second localizer offset to the affected side. Acquire second localizer in ISO positioning mode.
- If possible, MR staff to be present at CT Sim to assist with positioning and provide advice on immobilisation to ensure compatibility with MRI.
- Use Dixon where applicable for robust fat suppression particularly for upper extremities.
- For areas longer than 30 cm, such as long bones, consider TimCT to maintain geometrical integrity of anatomy.

RTP Requirements


- Use 3D distortion correction where applicable. If unavailable use 2D.
- Maintain a receiver bandwidth with a ≤ 1 pixel fat-water shift (≥ 400 Hz/Px at 3T).
- 0 mm interslice gap to eliminate interpolation in radiotherapy treatment planning (RTP) system.
- All sequences should have the same centre slice to maintain consistency with CT-MRI rigid registration in RTP system.
- 2–3 mm slice thickness will vary on anatomical site of interest. Ensure that the slices encompass entire extremity in an ideal acquisition time to maintain image quality.

Example Images



4 Tricep Sarcoma

Figure 4: T2 Dixon in-phase (4A), T1 Dixon in-phase (4B), T2 Dixon water (4C), co-registration of post-contrast T1 Dixon water (thermal window) with planning CT, showing agreement between two scans and good registration of tumor volume.



5 Gluteal Sarcoma

Figure 5: T2 Dixon in-phase (5A), co-registration of planning CT and T2 Dixon water sensitive planning MRI, showing robust fat suppression over a large field of view and excellent registration of bony and soft tissue anatomy (5B).




6 Metastatic disease in femur using TimCT for a 60 cm scan length

Figure 6: Reconstructed T2 HASTE with SPAIR showing hyperintense metastatic disease in bone marrow of femoral shaft (6A), coronal reconstructed T1 FLASH used for bone registration to planning CT (6B), accurate co-registration of planning CT with T2 HASTE planning MRI.

MAGNETOM World Community

Discover ways to optimize the use of your MR simulation applications.

Explore the latest clinical cases—from experts for experts!



The MAGNETOM World is the community of Siemens Healthineers MR users worldwide, providing you with relevant clinical information. Here you will find application tips and protocols developed by medical professionals that may optimize your daily work. Lectures and presentations from experts in the field will allow you to be exposed to new ideas and alternative clinical approaches.

Put the advantages of the MAGNETOM World
to work for you!

MAGNETOM systems are indicated for use as magnetic resonance diagnostic devices (MRDD) that produce transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that display the internal structure and/or function of the head, body or extremities. These images and/or

spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician, yield information that may assist in diagnosis.. The study protocols, application tips, lectures and presentations contained in this document were developed entirely by the medical professionals identified in the respective articles. These medical professionals did not receive financial support or any other kind of assistance from Siemens Healthineers.

Patient results may vary depending on many factors including system, software version, options, coils, technique and patient condition. Siemens Healthineers is not responsible or liable for results.

usa.siemens.com/magnetom-world-mr-in-rt

Notes

[illegible]

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

For USA Distribution Only.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens.com/healthineers

Local Contact Information

Siemens Medical Solutions USA, Inc.
40 Liberty Boulevard
Malvern, PA 19355-9998, USA
Phone: +1-888-826-9702
usa.siemens.com/healthineers