

# Left atrial appendage closure

*syngo DynaCT Cardiac has the potential to facilitate and accelerate the procedure, and make the treatment safer. syngo DynaCT Cardiac provides additional value when measuring the ostium and length of the LAA, and helps in selecting the correct device size.*

## Structural Heart Disease

### Study Protocol

#### Courtesy of

Christian Schlundt, MD,  
Vice Chairman of Cardiology,  
University Erlangen, Germany

#### Supported by

- syngo DynaCT Cardiac
- syngo iPilot
- syngo InSpace3D

#### System & Software

Artis Q.zen  
syngo X Workplace

### Case Description

#### Patient history

The patient was a 74-year-old male weighing 91 kg, with a height of 1.67 m.

#### Diagnosis

Permanent atrial fibrillation (AFib), contraindication for oral anti-coagulation due to severe bleeding events. Pre-procedural echo showed no thrombus in the left atrial appendage (LAA).

#### Treatment

LAA closure with occluder (Boston Scientific Watchman™)


#### Tips and tricks

- TEE probe was pulled back during acquisition.
- Arms can be kept in normal position beside the body.
- Use breathhold to reduce respiratory motion artifacts.
- No rapid pacing (if patient is in permanent AFib, rapid pacing can be skipped).
- Pigtail catheter with side holes (close to the center of the LAA ostium).


## Left atrial appendage closure

| Acquisition protocol  | 5s DCT Card                          |
|-----------------------|--------------------------------------|
| Number of projections | 248                                  |
| System dose           | 0.36 µGy/f                           |
| Increment in degrees  | 0.8 °/f                              |
| Injection protocol    |                                      |
| Catheter position     | Left atrial appendage                |
| Contrast medium (CM)  | 350 mg iodine/mL                     |
| Test bolus            | w/o                                  |
| Dilution              | 50 %                                 |
| Injection volume      | 65 mL<br>(32.5 mL CM/32.5 mL saline) |
| Injection rate        | 15 mL/s over 4.33 s                  |
| Duration of injection | 5 s                                  |
| X-ray delay           | 2 s                                  |
| Power injector used   | Yes                                  |


| Reconstructions       |           |
|-----------------------|-----------|
| VOI size              | Full      |
| Slice matrix          | 512 × 512 |
| Kernel type           | HU        |
| Image characteristics | Normal    |
| Reconstruction mode   | Nat fill  |
| Viewing preset        | Golden    |



AP projection from rotational angiography with direct injection into the LAA



syngo iPilot with overlay of 3D structure and landing zone onto the live fluoro image



syngo InSpace reconstruction of the LAA, with measurements of the ostium and the length of the LAA in MPRs

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.

All rights reserved.

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[siemens.com/healthineers](https://www.siemens.com/healthineers)