




Study Protocol

Visualization of clipped aneurysm

Neuro Interventions



DynaCT with intravenous contrast injection is a fast and easy procedure for non-invasive arterial imaging.

Courtesy of

Andrey Sergeev, MD,
Department of Neurosurgery,
The Federal State Budgetary
Institution "Almazov National
Medical Research Centre"
of the Ministry of Health,
Russia

Supported by

syngo DualVolume
syngo DynaCT DSA

System & Software

Artis zeego, VC21
syngo X-Workplace, VB21

Case Description

Patient history

Patient hemorrhaged from a giant ACom aneurysm. Urgent admission and surgery was performed the next day.

Diagnosis

Giant aneurysm of the ACom artery. Acute subarachnoid hemorrhage from the aneurysm, Hunt-Hess II.

Treatment

Acute craniotomy with intra-intracranial anastomosis A2-A2 and aneurysm clipping.

General comments

DynaCT with intravenous contrast injection is a fast and easy procedure for non-invasive arterial imaging. Intraoperative imaging of cerebral vessels allows neurosurgeons to visualize the effectiveness of clipped aneurysms and patency of bypasses immediately after the surgery.

Tips and tricks

The injection of the iodine contrast takes place directly after the mask run. We observe the contrast bolus using DSA. As soon as the carotid artery is filled with contrast, we start the robotic C-arm rotation around the patient's head. The head is held in position with fixing straps to reduce patient movement during 3D acquisition. syngo DynaCT acquired in Zoom 1 (42 cm) delivers better resolution due to the smaller voxel size compared to syngo DynaCT in Zoom 0 (48 cm), and results in better and more detailed visualization of the clips and vessels.

Visualization of clipped aneurysm

Acquisition protocol	10s DSA DCT Head, zoom 42 cm
----------------------	------------------------------

Injection protocol	
Catheter position	Cubital vein 18 g
Contrast medium (CM)	370 mg iodine/mL
Dilution (CM/Saline):	No
Injection volume	100 mL
Injection rate	5 mL/s
Duration of injection	20 s
X-ray delay	DSA bolus tracking
Power injector used	Yes

Reconstructions	Primary	Secondary
Name	DynaCT Head NatFill HU Normal HU	DynaCT Head Sub MoCo HU Very smooth
VOI size	Full	Full
Slice matrix	512 × 512	512 × 512
Kernel type	HU	HU
Image characteristics	Auto	Very smooth
Reconstruction mode	NatFill	Sub with MoCo (motion correction algorithm)
Viewing preset	Dual Volume 2	Dual Volume 2

Clinical Images



Figure 1: Clipped aneurysm. Intra-intracranial anastomosis A2-A2 with retrograde filling of right A2 vessel (secondary reconstruction)




Figure 2: 3D visualization of the clipped aneurysm (secondary reconstruction)



Figure 3: Transversal MIP with vessels and clips (primary reconstruction)

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.

All rights reserved.