

ARFI Imaging

Managing Patients with Liver Disease

Paul S. Sidhu, FRCR, Eleni Konstantatou, MD, Dylan Lewis, FRCR
Department of Radiology, King's College Hospital, Denmark Hill, London, United Kingdom




Figure 1. A focal liver lesion (arrow) on a background of mildly increased echogenicity. The patient had a previous sonographic examination suggesting the presence of a hemangioma in the liver; the appearances are not typical for a hemangioma, but in the presence of a fatty liver a hemangioma can be of low echogenicity.

Acoustic Radiation Force Impulse (ARFI) Imaging

A 62-year-old man was referred to the Ultrasound Department from his primary-care physician for a right upper quadrant and liver examination as part of investigations for continuing upper abdominal pain following an injury. His past medical history included diet-controlled diabetes and hypertension.

The sonographic examination demonstrated mild hepatosteatorosis, and a 36 x 36 mm mixed reflective lesion in the central aspect of the liver (Figure 1).

On questioning, the patient had a previous sonographic examination of the liver approximately 3 years previously and was told that he had a “benign” liver hemangioma.




Figure 2a. An arterialized lesion (arrow) is demonstrated on the CECT examination, suggestive but not conclusive of a HCC.




Figure 2b. Minor central washout (arrow) was seen on the portal venous image.

Further hematological work-up by the primary-care physician indicated a positive result for hepatitis C and a normal feto-protein (AFP). Further imaging with contrast-enhanced computed tomography (CECT)(Figure 2) and magnetic resonance (MR)(Figure 3) imaging indicated that the abnormality in the right liver lobe was vascular, but retained hepatospecific contrast at 20 minutes, atypical of a hepato-cellular carcinoma (HCC).

The indeterminate focal lesion needed to be subject to confirmation by histology for further appropriate clinical management; an HCC was suspected. The clinical management pathway needed to ascertain the presence of any underlying chronic liver fibrosis or cirrhosis; advanced underlying liver disease would alter the patient management pathway, as advanced fibrosis will not regress on treatment. The grey-scale sonographic evaluation alone demonstrated a fatty liver but no other evidence of chronic liver disease.

The need for containing the size of the HCC with a trans-arterial chemoembolization (TACE), and resection or then consideration of eventual liver transplant, is dependent on the extent of underlying diffuse liver disease. The underlying liver was subject to an ultrasound scan using Virtual Touch quantification (acoustic radiation force impulse [ARFI] imaging), with sampling from the right liver lobe demonstrating a median of 1.48 m/sec, (Figure 4). This result indicated that a biopsy of the underlying liver was also required. When the biopsy of the liver was performed, it confirmed the degree of fibrosis estimated from the ARFI measurement.

The HCC was subject to a TACE procedure with reduction in the size of the HCC (Figure 5). Liver transplantation is awaited.

This case indicates the importance of the role of ARFI in encouraging the need to stage the patient adequately for correct management, and in this case, biopsy of both the focal liver lesion (confirmed as HCC) and the background liver (METAVIR score F3) was performed to optimize clinical management.






Figure 3a. Gadolinium-enhanced MR imaging in the arterial phase demonstrates vascular enhancement of the lesion (arrow).

Figure 3b. Portal venous phase images indicate the possibility of an HCC.

Figure 3c. Hepatic-specific contrast imaging at 20 minutes shows wash-out of contrast (arrow), which is not a typical feature of HCC.




Figure 4. The ARFI measurement from deep within the right liver lobe indicates a reading that is higher than normal (1.47 m/sec). Ten separate readings are obtained and the median calculated.




Figure 5. The arterial phase of the TACE procedure demonstrates the lesion vascularity (arrow).

Bibliography

Jaffer OS, Lung PFC, Bosanac D, Patel V, Ryan SM, Heneghan MA, Quaglia A, Sidhu PS. acoustic radiation force impulse (ARFI) quantification: repeatability of measurements in selected liver segments and influence of age, body mass index (BMI) and liver capsule-to-box (CB) distance. *Br J Radiol.* 2012;85:e858-e863.

Friedrich-Rust M, Nierhoff J, Lupsor M, Sporea I, Fierbinteanu-Braticevici C, Strobel D, Takahashi H, Yoneda M, Suda T, Zeuzem S, Herrmann E. Performance of acoustic radiation force impulse imaging for the staging of liver fibrosis: a pooled meta-analysis. *J Viral Hepat.* 2012;19:e212-9.

Karlas T, Pfrepper C, Wiegand J, Wittekind C, Neuschulz M, Mössner J, Berg T, Tröltzsch M, Keim V. Acoustic radiation force impulse imaging (ARFI) for non-invasive detection of liver fibrosis: examination standards and evaluation of interlobe differences in healthy subjects and chronic liver disease. *Scand J Gastroenterol.* 2011;46:1458-67.

Ellis EL, Mann DA. Clinical evidence for the regression of liver fibrosis. *J Hepatol.* 2012;56:1171-1180.

Some/all of the features and products described here may not be available in the United States on all systems. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Please contact your local Siemens Healthineers sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Virtual Touch is a trademark of Siemens Medical Solutions USA, Inc.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 913184-0
siemens.com/healthineers

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
685 East Middlefield Road
Mountain View, CA 94043
USA
Phone: 1-888-826-9702
siemens.com/ultrasound