

Shear Wave Detection Imaging in the Liver

The Palpation Factor

New technologies in ultrasound elastography imaging, including Acoustic Radiation Force Impulse (ARFI) and shear wave detection imaging, are providing a new dimension of diagnostic information for physicians. At the Interdisciplinary Ultrasound Center of Germany's Munich University Hospital, Professor Dirk-André Clevert, MD, explores the possibilities of these new technologies and what they mean for the detection, diagnosis, and treatment of liver disease.

By Andrea McMurray

"An entirely new dimension of information from ultrasound is provided through elastography, or tissue strain imaging," says Professor Dirk Clevert. "In addition to 'seeing' and 'hearing' with ultrasound using B-mode and Doppler ultrasound imaging, ultrasound elastography adds the tactile part, allowing us to obtain information about the mechanical stiffness of tissue. In the past, we were only able to get this information via physical palpation – the age-old proven medical technique that has been used by physicians for over a thousand years."

Palpation requires a lot of experience, even when used in superficial organs and structures, as for example the breast or the thyroid. It is very difficult, if not impossible, to be applied with confidence in deeper regions of the body. The advent of new ultrasound elastography technologies is about to change that: Siemens Healthineers suite of Virtual Touch tissue strain analysis applications provides both qualitative and quantitative information about the stiffness of tissue, including deep abdominal tissue such as the liver.

This marks an important step in the development of elastography that has the potential to redefine how we use ultrasound in the diagnosis, treatment, and therapy of liver disease today.

How It Works

Unlike traditional, manual-compression ultrasound elastography, which has been around for several years now, Virtual Touch tissue strain analysis technologies use sound waves to compress tissue, a method known as Acoustic Radiation Force Impulse (ARFI) imaging. ARFI applies an ultrasound push pulse to a defined region of interest that causes displacement of the respective tissue. The extent of displacement and/or compression varies with the specific stiffness properties of different tissues: Stiff tissue will not be displaced/compressed as much as soft tissue.



“In addition to visual and audio information, I can now ‘virtually touch’ the tissue, allowing me to determine its texture.”

Professor Dirk-André Clevert, MD,
Section Head of the Interdisciplinary Ultrasound Center
at the University Hospital of Munich, Germany, considers
tissue strain analytics as the most important development
in ultrasound technology since the advent of Doppler imaging.




Figure 1. Tissue strain analysis provides a new dimension of diagnostic information from ultrasound technologies.





Figure 2. ARFI works almost like physical compression, but can penetrate deeper tissue. A push pulse generates lateral shear waves. The speed of shear waves is closely correlated to stiffness and can be accurately measured with new technologies. The ability to exactly measure the shear waves improves the consistency of exams.

i Virtual Touch imaging


Virtual Touch imaging is indicated for assessment of tissue stiffness in abdominal and small parts (breast, thyroid, testes, etc.).

Virtual Touch imaging uses a push pulse and creates a more detailed and reproducible qualitative elastogram than is possible with manual compression elastography.



q Virtual Touch quantification

Virtual Touch quantification is indicated for assessment of tissue stiffness in abdominal (including liver) and small parts.



The compression of tissue using ARFI push pulses generates shear waves that propagate perpendicularly to the push wave. An advanced, proprietary algorithm measures the velocity of these shear waves, which correlates with the stiffness of tissue at a defined region of interest.

Virtual Touch tissue strain analysis offers both ARFI and shear wave imaging capabilities that allow the assessment of tissue stiffness.

There are three different Virtual Touch technologies:

- Virtual Touch™ imaging produces a qualitative elastogram that is more detailed and reproducible than a conventional manual-compression elastogram.
- Virtual Touch™ quantification offers quantitative results by measuring the velocity of shear waves that are produced through push pulse compression.
- Virtual Touch™ IQ combines compression imaging and shear wave measurements on a single screen, allowing both the qualitative and quantitative assessment of tissue.

“The combination of all these technologies allows us to gather valuable additional parameters of diagnostic information that we didn’t have before,” says Clevert.


The applications provide stiffness measurements at a depth of up to 8 centimeters and stiffness tissue maps at up to 10 centimeters below the skin.

Integrating the Workflow

Clevert outlines the typical examination workflow when using tissue strain analysis. A thorough exam of the liver using B-mode and Doppler imaging provides an overall status of the liver tissue as well as blood flow in the hepatic and portal veins and other vessels. “This will tell me if there are anomalies that are signs of liver disease: fibrotic and cirrhotic changes, fluid, open vessels, general flow,” he says. That accomplished, the application of Virtual Touch quantification reveals a numerical value that allows a classification of tissue stiffness. “If I see a lesion in the standard B-mode, for example, a hemangioma, with Virtual Touch imaging I can immediately determine whether it is stiff or soft.”



Virtual Touch IQ is indicated for assessment of tissue stiffness in **small parts (breasts, thyroid, etc.)**.



Unlike conventional elastography, which provides the information on the basis of manual compression, Virtual Touch technologies use intelligent algorithms that automatically calculate the values. At the same time, they address the challenges related to user dependence and variability, as well as consistency and reproducibility. "That's a great advantage and an absolute prerequisite for integrating these applications into the clinical routine."

The Third Dimension of Information

According to Clevert, tissue strain analysis is the most important development in ultrasound technology since the advent of Doppler imaging. "In addition to visual and audio information, I can now 'virtually touch' the tissue, allowing me to determine its texture. In addition, Virtual Touch quantification lets me obtain a numerical value to confirm and exactly determine what I see." For Clevert, this is a decisive advantage over conventional elastography imaging, which provides visual information only. "It takes us beyond such statements as 'a little red or a little blue,' as provided by a conventional elastogram. It gives us a numerical value related to tissue stiffness which allows us to assess the state and development of the disease."

Visual information provided by B-mode and Doppler ultrasound imaging is an indispensable and valuable tool in liver diagnosis today; however, it has its limitations. Often, there are no visible anomalies, but the laboratory data continue to show that something is wrong. "Virtual Touch quantification allows us to measure whether there are abnormalities in the liver parenchyma indicating an increased stiffness of the liver before it becomes visible in the B-mode image. It is the increased stiffness as compared to the stiffness of a normal liver that adds the information to help with the diagnosis of abnormality,"


says Clevert. "With Virtual Touch applications, we have a very reliable and fast tool to assess tissue stiffness with the additional advantage of being able to see what we are doing. By combining numerical values with the visual ultrasound information of the tissue and flow, we can make sure that we have really assessed it all." In addition, Virtual Touch technologies can be used for a wide variety of patients, including those with advanced ascites or distinct portal hypertension.

Therapy Control

Determining the stiffness of the liver will be increasingly requested by referring physicians. "Virtual Touch analysis helps us detect anomalies in the liver at a very early stage, when we haven't seen anything on standard ultrasound yet. It tells us whether the liver is fibrotic or cirrhotic."

At Clevert's lab, the technology is also used to determine the success of interferon therapy used to treat chronic and acute hepatitis C. "The therapy with interferon is expensive; however, not all patients respond and it can have adverse effects. Virtual Touch imaging lets us monitor the success of the treatment, and change gear if necessary."



Clevert thinks that tissue strain analysis will be implemented in the clinical routine in the near future. He says that an experienced user can make the measurements in less than 3 minutes. "Who wouldn't want to make use of the additional information if it can be obtained so fast and in such an uncomplicated manner?" As the areas of use for ARFI and shear wave detection imaging expand, they will be integrated into a broader range of transducers. "These are great applications, and expanding their use will benefit all, patients and physicians alike."



A: Virtual Touch imaging delivers a qualitative assessment of tissue

B: Virtual Touch quantification shows a shear wave velocity of 0.70 m/s in the tissue surrounding the hemangioma

C: Inside the hemangioma, the shear wave velocity is clearly increased at 1.10 m/s



Andrea McMurray is a marketing communications manager at Siemens Healthineers Ultrasound Business Area in Mountain View, California.

Article Reprint: Siemens Medical Solutions Magazine, August 2009.

Because of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/all of the features and products described here may not be available in the United States on all systems. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described here without prior notice.

Please contact your local Siemens Healthineers sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Virtual Touch is a trademark of Siemens Medical Solutions USA, Inc.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens.com/healthineers

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
685 East Middlefield Road
Mountain View, CA 94043
USA
Phone: 1-888-826-9702
siemens.com/ultrasound