

MAGNETOM Flash

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vol.

14

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David Hillier

MR Application Specialist, Siemens Healthineers Great Britain and Ireland (GB&I)

10 Whole-body Diffusion-weighted MR Image Analysis with *syngo.via* Frontier MR Total Tumor Load

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Robert Grimm¹, Anwar R. Padhani²

¹ Siemens Healthineers, Erlangen, Germany

² Paul Strickland Scanner Centre, Mount Vernon Hospital, Northwood, Middlesex, UK

13 Metastatic Prostate Cancer in Practice – the MET-RADS-P Imaging Response System Using Whole-body MRI

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Anwar R. Padhani¹; Nina Tunariu²

¹ Paul Strickland Scanner Centre, Mount Vernon Hospital, Northwood, Middlesex, UK

² Royal Marsden NHS Foundation Trust & Institute of Cancer Research and Cancer Research UK Cancer Imaging Centre, Sutton, Surrey, UK

22 Observing Endocrine Therapy Resistance in Metastatic Breast Cancer with Whole-body MRI

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Anwar R. Padhani

Paul Strickland Scanner Centre, Mount Vernon Hospital, Northwood, Middlesex, UK

The information presented in MAGNETOM Flash is for illustration only and is not intended to be relied upon by the reader for instruction as to the practice of medicine. Any health care practitioner reading this information is reminded that they must use their own learning, training and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens Healthcare to be used for any purpose in that regard. The treating physician bears the sole responsibility for the diagnosis and treatment of patients, including drugs and doses prescribed in connection with such use. The Operating Instructions must always be strictly followed when operating the MR System. The source for the technical data is the corresponding data sheets.

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References

- 1 <https://www.nice.org.uk/guidance/ng35/chapter/recommendations>
- 2 Padhani, A.R., Douis, H. and Gall, P. (2012), Whole-Body Diffusion-Weighted MRI of the Bone Marrow in Health and Disease, MAGNETOM Flash, Vol 50, Issue 2 RSNA Edition, pp60-65.
- 3 Koh, D-M., Blackledge, M., Padhani, A.R., Takahara, T., Kwee, T.C., Leach, M.D. and Collins, S.J. (2012), Whole-Body Diffusion-Weighted MRI: Tips, Tricks and Pitfalls, American Journal of Roentgenology, Vol 199, pp252-262.
- 4 Smith, G. (2015), Using Remote Assist to Expand Your Expertise, MAGNETOM Flash, Vol 63, Issue 3 RSNA Edition, pp 87-89.
- 5 https://w5.siemens.com/web/ua/ru/medecine/detection_diagnosis/magnetic_resonans/035-15-MRI-scanners/MAGNETOM-Aera/Documents/mri_magnetom-aera_productbrochure_en-00044418.pdf
- 6 Furuta, A., Isoda, H., Yamashita, R., Ohno, T., Kawahara, S., Shimizu, H., Shibata, T. and Togashi, K. (2014), Comparison of Monopolar and Bipolar Diffusion Weighted Imaging Sequences for Detection of Small Hepatic Metastases, European Journal of Radiology, Vol 83 Issue 9, pp1626-1630.
- 7 Giles, S.L., deSouza, N.M., Collins, D.J., Morgan, V.A., West, S., Davies, F.E., Morgan, G.J. and Messiou, C. (2015), Assessing Myeloma Bone Disease With Whole-Body Diffusion-Weighted Imaging: Comparison With X-Ray Skeletal Survey by Region and Relationship With Laboratory Estimates of Disease Burden, Clinical Radiology, Vol 70, Issue 6, pp614-621.

Contact

David Hillier
Siemens Healthineers
Great Britain and Ireland (GB&I)
HC WEA GBR SV CS BD AM MR

Mobile: +44 7808 828554
david.hillier@siemens.com

Whole-body Diffusion-weighted MR Image Analysis with syngo.via Frontier MR Total Tumor Load

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References

- 1 Pawlyn C, Fowkes L, Otero S, et al. Whole-body diffusion-weighted MRI: a new gold standard for assessing disease burden in patients with multiple myeloma? *Leukemia*. 2016;30(6):1446-1448. doi:10.1038/Leu.2015.338.
- 2 Padhani AR, Lecouvet FE, Tunariu N, Koh DM, De Keyzer F, et al. (2016) Rationale for Modernising Imaging in Advanced Prostate Cancer. *European Urology Focus*, Available online. doi:10.1016/j.euf.2016.06.018.
- 3 Padhani AR and Tunariu N, Metastatic response system in using whole body MRI for prostate cancer in practice. See also page 76 in this issue of MAGNETOM Flash.
- 4 Blackledge MD, Collins DJ, Tunariu N, Orton MR, Padhani AR, et al. (2014) Assessment of Treatment Response by Total Tumor Volume and Global Apparent Diffusion Coefficient Using Diffusion-Weighted MRI in Patients with Metastatic Bone Disease: A Feasibility Study. *PLOS ONE* 9(4): e91779. doi: 10.1371/journal.pone.0091779.
- 5 Blackledge MD, Tunariu N, Orton MR, Padhani AR, Collins DJ, et al. (2016) Inter- and Intra-Observer Repeatability of Quantitative Whole-Body, Diffusion-Weighted Imaging (WBDWI) in Metastatic Bone Disease. *PLOS ONE* 11(4): e0153840. doi: 10.1371/journal.pone.0153840
- 6 Amlani A, Ghosh-Ray S, van Ree K, Makris A, D'Sa S, Ostler P, Anyamene N, Padhani AR. Relationships Between Diffusion Weighted Signal Intensity, ADC and Water/fat Content of Malignant Bone Marrow. Program Number: 3447. Proceedings of the International Society of Magnetic Resonance in Medicine. 20th Scientific Meeting and Exhibition, Salt Lake City, May 2013.
- 7 Blackledge MD, Leach M, Rata M, Tunariu N, Koh DM, et al. (2016) Visualizing whole-body treatment response heterogeneity using multi-parametric magnetic resonance imaging. *Journal of Algorithms & Computation Technology*, in press. doi: 10.1177/1748301816668024.

Contact	
	<p>Robert Grimm, Ph.D. Siemens Healthcare GmbH HC DI MR PI TIO ONCO</p> <p>Postbox 32 60 91050 Erlangen Germany Phone: +49 (0) 9131 84-2717 robertgrimm@siemens.com</p>

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Metastatic Prostate Cancer in Practice – the MET-RADS-P Imaging Response System Using Whole-body MRI

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1A Physician: Exam date current: Exam date comparison:

1A MET-RADS Prostate Report

Primary Involved RAC:
 Pelvic nodes Involved RAC:
 Retroperitoneal Involved RAC:
 Other nodes Involved RAC:
 Liver Involved RAC:
 Lungs Involved RAC:
 Other sites Involved RAC:

Bones MET-RADS:
 Skull Involved RAC:
 Cervical spine Involved RAC:
 Dorsal spine Involved RAC:
 Lumbar spine Involved RAC:
 Pelvis Involved RAC:
 Thorax Involved RAC:
 Limbs Involved RAC:

OVERALL ASSESSMENT

	No dis	CR	PR	SD	PD	DISCORD
Primary						
Nodes						
Viscera						
Bones						

Comments:
 Radiologist:
 Date:

1B RAC

Region	MET-RADS-P Description
1	Consistent with RECIST v1.1/PCWG criteria for unequivocal response (partial/complete). Return of normal marrow in areas previously infiltrated by focal/diffuse metastatic infiltration. Decrease in number/size of focal lesions. Evolution diffuse neoplastic pattern to focal lesions. Decreasing soft tissue associated with bone disease. Dense lesion sclerosis (edge to edge), sharply defined, very thin/disappearance of hyperintense rim on T2W-FS images. The emergence of intra/peritumoural fat within/around lesions (fat dot/halo sign). Previously evident lesion shows increase in ADC from $\leq 1400 \mu\text{m}^2/\text{s}$ to $>1400 \mu\text{m}^2/\text{s}$. Previously evident lesion shows increase in ADC from $\leq 1000 \mu\text{m}^2/\text{s}$ to $>1000 \mu\text{m}^2/\text{s}$. Morphological findings consistent with stable or responding disease.
2	Changes depicting tumour response that do not meet RECIST v1.1/PCWG criteria for partial or complete response (see below). Evidence of improvement, but not enough to fulfill criteria for RAC 1. For example: Previously evident lesions showing increases in ADC from $\leq 1000 \mu\text{m}^2/\text{s}$ to $>1000 \mu\text{m}^2/\text{s}$ but $<40\%$ increase in ADC from baseline with corresponding decrease in high b-value SI and morphological findings consistent with stable or responding disease.
3	No observable change. Changes depicting tumour progression that do not meet RECIST v1.1/PCWG criteria for progression. Evidence of worsening disease, but not enough to fulfill criteria for RAC 5. Equivocal appearance of new lesion(s). No change in size but increasing SI on high b-value images (with ADC values $<1400 \mu\text{m}^2/\text{s}$) consistent with possible disease progression. Relapse disease: re-emergence of lesion(s) that previously disappeared or enlargement of lesion(s) lesions that had partially regressed/stabilized with prior treatments. Imaging depicted bone lesions that might be clinically significant (therefore excludes asymptomatic fractures in non-critical bones). Soft tissue in spinal canal causing narrowing not associated with neurological findings and not requiring radiotherapy.
4	Tumour progression that meet RECIST v1.1/PCWG criteria for unequivocal progression. New critical fracture(s)/cord compression requiring radiotherapy/surgical intervention \rightarrow only if confirmed as malignant by MRI signal characteristics. Unequivocal new focal/diffuse area(s) of metastatic infiltration in regions of prior normal marrow. Unequivocal increase in number/size of focal lesions. Evolution of focal lesions to diffuse neoplastic pattern. Appearance/increasing soft tissue associated with bone disease. New lesion/regions of high signal intensity on high b-value images with ADC value between $800-1000 \mu\text{m}^2/\text{s}$.
5	Highly likely to be progressing.

Response Assessment Category (RAC) allocation rules – compare to relevant baseline scan

The primary RAC value is based on the dominant response of more than half of the disease within the region; The secondary RAC value is for the second most frequent pattern of response.
For a single lesion in a region only the primary number category is assessed. Regions with multiple lesions/diffuse disease, all with the same RAC, both the primary and secondary have the same values.
When equal numbers of lesions are of higher and lower RACs then the primary pattern allocation is reserved for the higher RAC.

RECIST v1.1 categories

- Complete Response (CR): Disappearance of all target lesions.
- Partial Response (PR): At least a 30% decrease in the sum of the longest diameter (LD) of target lesions, taking as reference the baseline sum LD.
- Stable Disease (SD): Neither sufficient shrinkage to qualify for PR nor sufficient increase to qualify for PD, taking as reference the smallest sum LD since the treatment started.
- Progressive Disease (PD): At least a 30% increase in the sum of the LD of target lesions, taking as reference the smallest sum LD recorded since the treatment started or the appearance of one or more new lesions.

Progression of focal prostate disease: Use RECIST v1.1 for progression criteria above applied to local disease.
Progression of nodes (short axis)

- <1 cm nodes have to have grown by at least 5mm in from baseline or otherwise radii and be ≥ 1 cm to be considered to have progressed.
- For nodes that are 1.0-1.5 cm that have grown by at least 5 mm from baseline or treatment radii and are ≥ 1.0 cm in short axis can be considered to have progressed.
- For nodes ≥ 1.5 cm short axis use RECIST v1.1 progression criteria.

Progression of visceral disease: Use RECIST v1.1 progression criteria above applied to visceral disease.

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2B 30/08/2016 MRI Whole body

Clinical details: mCRPC. Restaging post urinary diversion. On Abiraterone and zoledronic acid.

Technique: A whole-body MRI scan with whole body diffusion sequences. Comparison is made to the previous whole-body MRI scan dated 29/04/2016.

Findings:

Cervical and dorsal spine:
The intervertebral bony alignment is normal. Regrettably, there is marked disease progression. New metastases are seen throughout the cervical and dorsal spine with multifocal lesions. No interval loss of vertebral height. The craniocervical junction is normal. The cervical and dorsal cord outline normally.

Lumbosacral spine:
The intervertebral bony alignment remains normal with no interval loss of vertebral height. Degenerative spinal stenosis at L3/L4 as noted on the previous occasion. Regrettably, there is marked disease progression in the lumbosacral spine since the previous study.

Body scan:
No skull vault deposits have emerged. Normal sinonasal airways.
No supraclavicular fossa lymphadenopathy. Progressive left axillary lymphadenopathy also.
There is marked disease progression left scapula bone with extra osseous soft tissue disease now visible. Marked disease progression the ribs bilaterally also. There are artefacts are sternotomy wires.
The central mediastinal and hilar regions are normal. No lung abnormalities are detected.
The liver and spleen are homogeneous. Normal pancreas and adrenal glands. Both kidneys are unobstructed with bilateral renal stents in situ. Small retroperitoneal lymph nodes are also detected.
Extensive metastatic bone disease is present in the sacrum predominantly on the left side, right and left hemipelvis bone disease also.
Nodal disease in the common iliac regions bilaterally. Right obturator region with extra nodal tumour spread.
There is large locally advanced prostate carcinoma with bladder and ureteric involvement. No tumour involvement of the rectum or rectosigmoid junction.
Metastatic disease in the right proximal femur also.

Impression:
There is marked disease progression since 29/04/2016. Disease progression is seen locally the prostate gland with extensive bladder invasion. There is disease progression within pelvic and retroperitoneal lymph nodes with extra nodal tumour spread. There is disease progression in the left axilla also. Bone disease progression additionally throughout the spine with extra-osseous soft tissue disease also visible. No new visceral relapse of disease.
Please see graphical MET-RADS-P report also.

2C Physician: PO, Exam 30/8/16, Exam 29/4/16

2C MET-RADS Prostate Report

Primary Involved RAC: 5/5
Rectal and bladder involvement
Pelvic nodes Involved RAC: 5/5
Retroperitoneal Involved RAC: 5/5
Other nodes Involved RAC: 5/5
Liver Involved RAC: 5/5
Lungs Involved RAC: 5/5
Other sites Involved RAC: 5/5

Bones MET-RADS:
 Skull Involved RAC: N
 Cervical spine Involved RAC: 5/5
 Dorsal spine Involved RAC: 5/5
 Lumbar spine Involved RAC: 5/5
 Pelvis Involved RAC: 5/5
 Thorax Involved RAC: 5/5
 Limbs Involved RAC: 5/5

OVERALL ASSESSMENT

	No dis	CR	PR	SD	PD	DISCORD
Primary						
Nodes						
Viscera						
Bones						

Comments: Bladder and rectal involvement by primary tumor.

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4D Physician: Po, Exam: 10/3/15

Soft Tissues RECIST

- Primary Involved RAC: 3/4
- Pelvic nodes Involved RAC: 5/2
- Retroperitoneal Involved RAC: 5/5
- Other nodes Involved RAC: 5/5 (axilla)
- Liver Involved RAC: N
- Lungs Involved RAC: N
- Other sites Involved RAC: N

Bones MET-RADS

- Skull Involved RAC: N
- Cervical spine Involved RAC: 2/3
- Dorsal spine Involved RAC: 2/3
- Lumbosacral spine Involved RAC: 5/5
- Pelvis Involved RAC: 5/5
- Thorax Involved RAC: 3/5
- Limbs Involved RAC: 5/2

OVERALL ASSESSMENT

	No dis	CR	PR	SD	PD	DISCORD
Primary						
Nodes						Not in RECIST category
Viscera	✓					TRUS/US/US/PS/US/PS/US/PS
Bones				2		major

Comments: Progression in nodes. Response (pseudoprogression) in dorsal spine but true progression in lumbar spine.

Response assessment categories (RAC): 1 Highly likely to be responding; 2 Likely to be responding; 3 No change; 4 Likely to be progressing; 5 Highly likely to be progressing. Single lesion? RAC only; 22 lesions/diffuse disease use both RACs.

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- Lecouvet FE. Whole-Body MR Imaging: Musculoskeletal Applications. *Radiology*. 2016;279:345–365.
- Padhani AR, Lecouvet FE, Tunariu N, et al. Rationale for Modernising Imaging in Advanced Prostate Cancer. *Eur Urol Focus*. European Association of Urology; 2016;44:198–205.
- Dimopoulos M a., Hillengass J, Usmani S, et al. Role of magnetic resonance imaging in the management of patients with multiple myeloma: a consensus statement. *J Clin Oncol*. 2015;33:657–664.
- Padhani AR, Lecouvet FE, Tunariu N, et al. METastasis Reporting and Data System for Prostate Cancer: Practical Guidelines for Acquisition, Interpretation, and Reporting of Whole-body Magnetic Resonance Imaging-based Evaluations of Multiorgan Involvement in Advanced Prostate Cancer. *Eur Urol*. European Association of Urology; 2017;71:81–92.
- Mottet N, Bellmunt J, Bolla M, et al. EAU-ESTRO-SIOG Guidelines on Prostate Cancer. Part 1: Screening, Diagnosis, and Local Treatment with Curative Intent. *Eur Urol*. 2016;1–12.
- Cornford P, Bellmunt J, Bolla M, et al. EAU-ESTRO-SIOG Guidelines on Prostate Cancer. Part II: Treatment of Relapsing, Metastatic, and Castration-Resistant Prostate Cancer. *Eur Urol*. 2016;1–13.
- Eisenhauer EA, Therasse P, Bogaerts J, et al. New response evaluation criteria in solid tumours: revised RECIST guideline (version 1.1). *Eur J Cancer*. 2008/12/23. 2009;45:228–247.
- Scher HI, Morris MJ, Stadler WM, et al. Trial Design and Objectives for Castration-Resistant Prostate Cancer: Updated Recommendations From the Prostate Cancer Clinical Trials Working Group 3. *J Clin Oncol*. 2016;34:1402–1418.

Contact

Professor Anwar R. Padhani
MB, BS, FRCP, FRCR
Paul Strickland Scanner Centre
Mount Vernon Hospital

Rickmansworth Road
Northwood, Middlesex HA6 2RN
United Kingdom
anwar.padhani@stricklandscanner.org.uk

Contact

Nina Tunariu
MD MRCP FRCR MDR
Royal Marsden NHS Foundation Trust &
Institute of Cancer Research
Cancer Research UK Cancer Imaging Centre


Downs Road, Sutton, Surrey, SM2 5PT
United Kingdom
Nina.Tunariu@icr.ac.uk

Visit us at

www.siemens.com/magnetom-world

to download the MET-RADS-P template form.

Observing Endocrine Therapy Resistance in Metastatic Breast Cancer with Whole-body MRI

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References

- Jeselson R, Buchwalter G, De Angelis C, Brown M, Schiff R. ESR1 mutations – a mechanism for acquired endocrine resistance in breast cancer. *Nat Rev Clin Oncol*. 2015;12:573–583.
- Padhani AR, Lecouvet FE, Tunariu N, et al. METastasis Reporting and Data System for Prostate Cancer: Practical Guidelines for Acquisition, Interpretation, and Reporting of Whole-body Magnetic Resonance Imaging-based Evaluations of Multiorgan Involvement in Advanced Prostate Cancer. *Eur Urol*. European Association of Urology; 2017;71:81–92.

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Rickmansworth Road
Northwood, Middlesex HA6 2RN
United Kingdom
anwar.padhani@
stricklandscanner.org.uk

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