

# **Case Study**

# Robust quantification of neuroendocrine tumor imaging: xSPECT Quant with <sup>111</sup>In tracers enables new possibilities in SPECT/CT imaging

By Mario Jreige, Gilles Allenbach and J. Prior Lausanne University Hospital, Lausanne, Switzerland

## History

A patient with a well-differentiated neuroendocrine tumor stage IV [HEP, LYM] of the pancreas G2, MIB1 ~5% in March 2014, was treated with different chemotherapy regimens with progressive disease until January 2016. After two peptide receptor radionuclide therapy (PRRT) cycles, the patient was referred for restaging in August 2016.

Two- and 24-hour post-injection scans were performed after <sup>111</sup>In-octreotide using a Symbia Intevo<sup>™1</sup> 16 SPECT/CT scanner with xSPECT Quant<sup>™1</sup> technology.

### Diagnosis

This scan demonstrated the metastatic disease with focal mass uptake in the pancreatic region as well as the disseminated spread, as shown in figures

1 and 2. Compared to a prior examination (not shown), the final results were concluded as stable disease after PRRT.

#### Comments

of several different diseases, but are mostly used in the early detection, staging and potential treatment of neuroendocrine tumors of the upper and mid-gut. With the advent of new therapy approaches in neuroendocrine tumors such as PRRT, the necessity for accurate uptake measurements has emerged. Quantitative response evaluation in neuroendocrine tumors is primarily in the domain of PET imaging, but 111 In-labelled octreotide is still one of the most commonly used agents for imaging of neuroendocrine tumors in

diagnosis and staging, in particular, where PET tracers are not available.2 Although it is known that there is a correlation of somatostatin receptor expression and the visibility in somatostatin receptor imaging, due to the lack of robust and repeatable quantification, visual scores were recommended.3 The feasibility of 99mTc-MDP quantification and its potential clinical impact was already demonstrated.4 Now, with the expansion of accurate and robust quantification using xSPECT Quant<sup>™1</sup> for additional tracers, clinicians have started to define the meaning of these uptake values.

For this case, considering the tracers' distribution and uptake after 2 hours and 24 hours and, in particular, the uptake values calculated by SUV, one






Figure 1: Coronal view of a fused xSPECT/CT image with measurements of one of the liver lesions as well as the hilar lymph node metastasis. Note the automated SUV measurements using syngo®.via. The tracer distribution, retention and elimination can be identified by comparing the two time points.

Data courtesy of Department of Nuclear Medicine and Molecular Imaging


Cervical


2-hour post-injection

24-hour post-injection

Figure 2: Coronal view of the maximum intensity projection (MIP) of the xSPECT data from both time points, showing four different lesions. Selection of the lesions is for demonstration purposes only.

Data courtesy of Department of Nuclear Medicine and Molecular Imaging





2-hour post-injection

24-hour post-injection




Figure 3: Plotted SUV from the four selected lesions over time [2-hr., 24-hr.], demonstrate the quantified changes of the lesions.

can identify and quantify the biology of the tumor's behavior in terms of octreotide receptor binding, tracer target density, as well as the accumulation in potential organs at risk. As an example, Figure 3 illustrates the changes of the SUV between time point 1 at 2 hours, and time point 2 at 24 hours post-injection of all lesions measured in the overview in Figure 2. Using 111In, as presented here, physicians worldwide can overcome the initial and inherent limitations of regular SPECT/CT and use accurate and reproducible quantification for uptake calculations and precise follow-up evaluation.

# Conclusion

xSPECT Quant-enabled quantification for 111 In helps to measure precisely the initial uptake and changes over time in primary diagnosis, treatment follow-up and treatment planning.

<sup>1</sup> Symbia Intevo and xSPECT Quant are not commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

#### References:

- <sup>2</sup> Kwekkeboom D, Kam B, van Essen M, et al. Somatostatin receptor-based imaging and therapy of gastroenteropancreatic neuroendocrine tumors. Endocrine-Related Cancer 2010;17(1):R53-73.
- <sup>3</sup> Diakatou E, Alexandraki K, Tsolakis A, et al. Somatostatin and dopamine receptor expression in neuroendocrine neoplasms: correlation of immunohistochemical findings with somatostatin receptor scintigraphy visual scores. Clin Endocrinol 2015;83(3):420-8.
- Salaun P-Y, Abgral R, Laroche RD. xSPECT Quant in Treatment Monitoring. 2017. siemens.com/mi | xSPECT Case Study.

# **Examination Protocol**

Scanner: Symbia Intevo 16

SPECT	
Injected dose (representative scan)	180 MBq (4.9 mCi) 111In-octreotide
Scan delay	2-hour, 24-hour
Acquisition	xSPECT Quant

CT		2-hour	24-hour
	KVP	130 kV	130 kV
	Current	28 mAs	24 mAs
	CTDIvol	1.8816	1.6128

The statements by Siemens customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prrior notice. Some/all of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the mose current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

"Siemens Healthineers" is considered a brand name. Its use is not intended to represent the legal entity to which this product is registered. Please contact your local Siemens organization for further details.

### Siemens Healthcare Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen Germany Phone: +49 9131 84-0 siemens.com/healthcare

# **Global Business Line**

Siemens Medical Solutions USA, Inc. Molecular Imaging 2501 North Barrington Road Hoffman Estates, IL 60192 USA

Phone: +1 847 304-7700 siemens.com/mi