

A Pancreatic Squamous Cell Carcinoma

By Silva, C.F.G, BS; Bertolazzi, P, BS; Mendes, R.E.V, BS; Oliveira, F.F.; Viana, P.C.C, MD; Horvat, N, MD
Radiology Department, Sírío Libanês Hospital, Sao Paulo, Brazil

History

A 67-year-old male patient, suffering from loss of appetite and weight within the last 2 months, came to the hospital for a checkup. His medical history was unremarkable. A Dual Energy (DE) CT was performed to investigate a suspected pancreatic lesion.

Diagnosis


DECT images showed an irregular, heterogeneously enhanced lesion in the head of the pancreas, measuring 145.9 HU in density and 4.7 mg/mL in iodine uptake. The image contrast was significantly enhanced using Monoenergetic Plus (Mono+) at 45 keV.

The 3D image perspective could be better demonstrated using cinematic VRT (cVRT) for improved depth and shape perceptions. A pancreas tumor was suspected and subsequently confirmed by a PET-CT. A biopsy was performed which revealed a squamous cell carcinoma.


The patient underwent a pancreatectomy, followed by chemotherapy. In a 3-month follow-up, he was in good general condition with improved symptoms.

Comments


Pancreatic squamous cell carcinoma is an extremely rare subtype of pancreatic cancer of ductal origin. It is characterized by enhancement on contrast CT due to its hypervascularity. The dataset acquired by DECT can be processed to improve the contrast enhancement using *syngo*.CT DE Monoenergetic Plus, to both present the iodine map as well as to quantify the iodine uptake using *syngo*.CT DE Virtual Unenhanced. Using cVRT allows a more photo-realistic representation of the CT image data, with an enhanced and more natural depth and shape perception. ●



1 Mixed image (Fig. 1a), Mono+ image (Fig. 1b) and iodine map (Fig. 1c) show an irregular, heterogeneously enhanced lesion (arrows) in the head of the pancreas, with an iodine uptake of 4.7 mg/mL (Fig. 1c), which is almost double that in the adjacent normal pancreatic tissues (2.4 mg/mL). The measurement is normalized to the enhancement of the aorta. The enhancement is significantly improved in the Mono+ image display.



2 Quantitative measurements and graphical presentation of the CT attenuation in accordance with energy (keV) levels. Comparing 45 keV with 70 keV, the contrast difference between the lesion and the normal tissue is almost threefold.



3 cVRT image (Fig. 3a) demonstrates a better 3D perspective with improved depth and shape perceptions. A PET-CT image (Fig. 3b) confirmed the lesion in the head of the pancreas.

Examination Protocol

Scanner	SOMATOM Definition Flash		
Scan area	Upper abdomen	Slice collimation	32 × 0.6 mm
Scan mode	Dual Source Dual Energy	Slice width	3 mm
Scan length	302.4 mm	Reconstruction increment	2.4 mm
Scan direction	Cranio-caudal	Reconstruction kernel	D30f
Scan time	13 s	Contrast	350 mg/mL
Tube voltage	80 / Sn140 kV	Volume	120 mL + 50 mL saline
Effective mAs	145 / 67 mAs	Flow rate	5 mL/s
CTDI _{vol}	5.96 mGy	Start delay	Bolus tracking in the descending aorta with a threshold of 150 HU and an additional delay of 10 s
DLP	191 mGy cm		
Rotation time	0.5 s		
Pitch	0.6		

The outcomes by Siemens Healthineers customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.