MAGNETOM Flash

The Magazine of MRI 翻訳版

Issue Number 1/2016 | SCMR Edition

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
Page 15

The MAGNETOM World

Your portal to talks, articles and case studies

siemens.com/magnetom-world

Siemens' global MRI community offers peer-to-peer support and information. Radiologists, physicists, cardiologists, and technologists, have all contributed with publications, presentations, training documents, case studies, and more – all freely available to you via this unique network.


Don't miss the >300 lectures and presentations by international and renowned experts on all aspects of MRI that will allow you to be exposed to new ideas and alternative approaches.

Put the advantages of the MAGNETOM World to work for you!

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The centerpiece of the MAGNETOM World Internet platform consists of MAGNETOM users' results. Here you will find case reports, articles and application tips allowing you to optimize your daily work.

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University of Kentucky Medical Center, Gill Heart Institute, Lexington, KY, USA

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References

- White SK, Sado DM, Flett AS, Moon JC. Characterising the myocardial interstitial space: the clinical relevance of non-invasive imaging. Heart 2012;98:773-9.
- Dazzi H, Kaufmann, Follath. Anthracycline-induced acute cardiotoxicity in adults treated for leukaemia. Analysis of the clinico-pathological aspects of documented acute anthracycline induced cardiotoxicity in patients treated for acute leukaemia at the University Hospital of Zurich, Switzerland, between 1990 and 1996. Ann Oncol 2001; 12(7): 963-6.
- Bristow MR, Thompson PD, Martin RP, et al. Early anthracycline cardiotoxicity. Am J Med 1978; 65(5): 823-32.

Contact

Steve W. Leung, M.D., FACC
University of Kentucky Medical Center
Division of Cardiovascular Medicine
Gill Heart Institute
326 Charles T. Wethington Building
Lexington, KY 40536
USA
Phone: +1 859-323-8040
steve.leung@uky.edu

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Jeremy D. Collins, M.D.¹; Bruce Spottiswoode, Ph.D.²; James C. Carr, M.D.¹

¹ Feinberg School of Medicine, Department of Radiology, Sections of Cardiovascular Imaging and Interventional Radiology, Northwestern University, Chicago, IL, USA

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References

1. Lundin M, Ugander M. Clinical Utility of Cardiac T1- and Extracellular Volume (ECV) Mapping. A Brief Review. MAGNETOM Flash 1/2015: 18-20.
2. Maestrini V, Abdel-Gadir A, Herrey AS, Moon JC. New Generation Cardiac Parametric Mapping: the Clinical Role of T1 and T2 Mapping. MAGNETOM Flash 5/2013: 104-107.
3. Schelbert EB, Wong TC. Clinical Benefits of T1 and ECV Mapping. MAGNETOM Flash 1/2015: 12-17.
4. Xue H, Greiser A, Zuehlsdorff S, Jolly MP, Guehring J, Arai AE et al. Phase-Sensitive Inversion Recovery for Myocardial T1 Mapping with Motion Correction and Parametric Fitting. Magn Reson Med. 2013 May; 69(5): 1408-1420.
5. Gahl S, Shah S, Xui H, Chung YC, Pennell ML, Guehring J et al. Myocardial T2 Mapping With Respiratory Navigator and Automatic Nonrigid Motion Correction. Magn Reson Med. 2012 November; 68(5): 1570-1578.

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
References

- 1 Simonetti OP, Kim RJ, Fieno DS, Hillenbrand HB, Wu E, Bundy JM, Finn JP, Judd RM. An improved MR imaging technique for the visualization of myocardial infarction. *Radiology* 2001;218:215-223. doi: 10.1148/radiology.218.1.r01ja50215.
- 2 Kim RJ, Fieno DS, Parrish TB, Harris K, Chen E-L, Simonetti O, Bundy J, Finn JP, Klocke FJ, Judd RM. Relationship of MRI Delayed Contrast Enhancement to Irreversible Injury, Infarct Age, and Contractile Function. *Circulation* 1999;100:1992-2002. doi: 10.1161/01.CIR.100.19.1992.
- 3 Viallon M, Jacquier A, Rotaru C, Delattre BMA, Mewton N, Vincent F, Croisille P. Head-to-head comparison of eight late gadolinium-enhanced cardiac MR (LGE CMR) sequences at 1.5 tesla: from bench to bedside. *J Magn Reson Imaging* 2011;34:1374-1387. doi: 10.1002/jmri.22783.
- 4 Kellman P, Arai AE, McVeigh ER, Aletras AH. Phase-sensitive inversion recovery for detecting myocardial infarction using gadolinium-delayed hyperenhancement. *Magn Reson Med* 2002;47:372-383. doi: 10.1002/mrm.10051.
- 5 Varga-Szemes A, van der Geest RJ, Spottiswoode BS, Muscogiuri G, De Cecco CN, Suranyi P, Rehwald WG, Schoepf JU. Quantification of myocardial late gadolinium enhancement using synthetic inversion recovery imaging. *Journal of Cardiovascular Magnetic Resonance* 2015;17(Suppl 1):O8.
- 6 Kim HW, Farzaneh-Far A, Kim RJ. Cardiovascular Magnetic Resonance in Patients with Myocardial Infarction. *Journal of the American College of Cardiology* 2010;55. doi:10.1016/j.jacc.2009.06.059.

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