# MAGNETOM Flash

The Magazine of MRI 翻訳版

Issue Number 2/2016 | ISMRM Edition

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
Page 18

# The MAGNETOM World


Your portal to talks, articles and case studies

[siemens.com/magnetom-world](https://www.siemens.com/magnetom-world)

Siemens' global MRI community offers peer-to-peer support and information. Radiologists, physicists, cardiologists, and technologists, have all contributed with publications, presentations, training documents, case studies, and more – all freely available to you via this unique network.

Don't miss the >300 lectures and presentations by international and renowned experts on all aspects of MRI that will allow you to be exposed to new ideas and alternative approaches.

The centerpiece of the MAGNETOM World Internet platform consists of MAGNETOM users' results. Here you will find case reports, articles and application tips allowing you to optimize your daily work.

Put the advantages of the MAGNETOM World to work for you!

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Cover image courtesy of Vikas Gulani, M.D.

( Department of Radiology Case Western Reserve University University Hospitals Case Medical Center 11100 Euclid Ave- Bolwell Building, Room B120 Cleveland, OH 44106 USA)

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Departments of Radiology (1) and Biomedical Engineering (2), Case Western Reserve University, University Hospitals Case Medical Center, Cleveland, OH, USA

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## References

- 1 Larsson H, Frederiksen J. Assessment of demyelination, edema, and gliosis by in vivo determination of T1 and T2 in the brain of patients with acute attack of multiple sclerosis. *Magn Reson Med* [Internet]. 1989;348:337–48. Available from: <http://onlinelibrary.wiley.com/doi/10.1002/mrm.1910110308/abstract>
- 2 Usman AA, Taimen K, Wasielewski M, et al. Cardiac magnetic resonance T2 mapping in the monitoring and follow-up of acute cardiac transplant rejection: a pilot study. *Circ Cardiovasc Imaging* [Internet]. 2012;5(6):782–90. Available

- from: <http://www.ncbi.nlm.nih.gov/pubmed/23071145>
- 3 Payne AR, Berry C, Kellman P, et al. Bright-blood T 2-weighted MRI has high diagnostic accuracy for myocardial hemorrhage in myocardial infarction a preclinical validation study in swine. *Circ Cardiovasc Imaging*. 2011;4(6):738–45.
  - 4 Van Heeswijk RB, Feliciano H, Bongard C, et al. Free-breathing 3 T magnetic resonance T2-mapping of the heart. *JACC Cardiovasc Imaging* [Internet]. 2012 Dec;5(12):1231–9. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/23236973>
  - 5 lles L, Pfluger H, Phrommintikul A, et al. Evaluation of Diffuse Myocardial Fibrosis in Heart Failure With Cardiac Magnetic Resonance Contrast-Enhanced T1 Mapping. *J Am Coll Cardiol*. 2008;52(19):1574–80.
  - 6 Hernando D, Levin YS, Sirlin CB, Reeder SB. Quantification of liver iron with MRI: state of the art and remaining challenges. *J Magn Reson Imaging* [Internet]. 2014 Nov;40(5):1003–21. Available from: <http://www.pubmed-central.nih.gov/articlerender.fcgi?artid=4308740&tool=pmcentrez&rendertype=abstract>
  - 7 Reeder SB, Cruite I, Hamilton G, Sirlin CB. Quantitative assessment of liver fat with magnetic resonance imaging and spectroscopy. *J Magn Reson Imaging* [Internet]. 2011 Oct;34(4):729–49. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/21928307>
  - 8 Look DC, Locker DR. Time saving in measurement of NMR and EPR relaxation times. *Rev Sci Instrum*. 1970;41(2):250–1.
  - 9 Huang TY, Liu YJ, Stemmer A, Poncelet BP. T2 measurement of the human myocardium using a T 2-prepared transient-state trueFISP sequence. *Magn Reson Med*. 2007;57:960–6.
  - 10 Mehta BB, Chen X, Bilchick KC, Salerno M, Epstein FH. Accelerated and navigator-gated look-locker imaging for cardiac T1 estimation (ANGIE): Development and application to T1 mapping of the right ventricle. *Magn Reson Med*. 2015;73(1):150–60.
  - 11 Zhu DC, Penn RD. Full-brain T1 mapping through inversion recovery fast spin echo imaging with time-efficient slice ordering. *Magn Reson Med*. 2005;54(3):725–31.
  - 12 Cheng H-LM, Wright GA. Rapid High-Resolution T1 Mapping by Variable Flip Angles: Accurate and Precise Measurements in the Presence of Radiofrequency Field Inhomogeneity. *Magn Reson Med*. 2006;55:566–76.
  - 13 Doneva M, Börnert P, Eggers H, Stehning C, S  n  gas J, Mertins A. Compressed sensing reconstruction for magnetic resonance parameter mapping. *Magn Reson Med*. 2010;64(4):1114–20.
  - 14 Deoni SCL, Rutt BK, Peters TM. Rapid combined T1 and T2 mapping using

- gradient recalled acquisition in the steady state. *Magn Reson Med*. 2003;49(3):515–26.
- 15 Blume U, Lockie T, Stehning C, et al. Interleaved T1 and T2 relaxation time mapping for cardiac applications. *J Magn Reson Imaging*. 2009;29:480–7.
  - 16 Warntjes JBM, Dahlqvist O, Lundberg P. Novel method for rapid, simultaneous T1, T2\*, and proton density quantification. *Magn Reson Med*. 2007;57(3):528–37.
  - 17 Schmitt P, Griswold MA, Jakob PM, et al. Inversion recovery TrueFISP: quantification of T(1), T(2), and spin density. *Magn Reson Med*. 2004;51(4):661–7.
  - 18 Ehses P, Seiberlich N, Ma D, et al. IR TrueFISP with a golden-ratio-based radial readout: Fast quantification of T1, T2, and proton density. *Magn Reson Med*. 2013;69(1):71–81.
  - 19 Ma D, Gulani V, Liu K, et al. Magnetic resonance fingerprinting. *Nature* [Internet]. 2013;495(7440):187–92.
  - 20 Jiang Y, Ma D, Seiberlich N, Gulani V, Griswold M a. MR fingerprinting using fast imaging with steady state precession (FISP) with spiral readout. *Magn Reson Med* [Internet]. 2014;. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/25491018>
  - 21 Chen Y, Jiang Y, Pahva Shivani, et al. MR Fingerprinting for Rapid Quantitative Abdominal Imaging. *Radiology*. 2016;000(0):1–9.
  - 22 Schmitt P, Griswold MA, Gulani V, Haase A, Flentje M, Jakob PM. A simple geometrical description of the TrueFISP ideal transient and steady-state signal. *Magn Reson Med*. 2006;55(1):177–86.
  - 23 Perlin K. An image synthesizer. *ACM SIGGRAPH Comput Graph*. 1985;19(3):287–96.
  - 24 Ma D, Pierre EY, Jiang Y, et al. Music-based magnetic resonance fingerprinting to improve patient comfort during MRI examinations. *Magn Reson Med* [Internet]. 2015; Available from: <http://doi.wiley.com/10.1002/mrm.25818>
  - 25 Lee JH, Hargreaves B a., Hu BS, Nishimura DG. Fast 3D Imaging Using Variable-Density Spiral Trajectories with Applications to Limb Perfusion. *Magn Reson Med*. 2003;50(6):1276–85.
  - 26 Bloch F. Nuclear induction. *Phys Rev*. 1946;70:460–85.
  - 27 Jiang Y, Wright KL, Seiberlich N, Gulani V, Griswold MA. Simultaneous T1, T2, diffusion and proton density quantification with MR fingerprinting. In: In proceedings of the 22nd annual meeting of ISMRM meeting & exhibition in Milan, Italy. 2014. p. 28.
  - 28 Wright KL, Ma D, Jiang Y, Gulani V, Griswold MA, Luis H-G. Theoretical framework for MR fingerprinting with ASL: simultaneous quantification of CBF, transit time, and T1. In: In proceedings of the 22nd annual meeting of ISMRM meeting & exhibition in Milan, Italy. 2014. p. 417.
  - 29 Weigel M, Schwenk S, Kiselev VG, Scheffler K, Hennig J. Extended phase graphs with anisotropic diffusion. *J Magn Reson*. 2010;205(2):276–85.

- 30 Deshmane AV, Ma D, Jiang Y, et al. Validation of Tissue Characterization in Mixed Voxels Using MR Fingerprinting. In: In proceedings of the 22nd annual meeting of ISMRM meeting & exhibition in Milan, Italy. 2014. p. 0094.
- 31 McGivney D, Pierre E, Ma D, et al. SVD Compression for Magnetic Resonance Fingerprinting in the Time Domain. *IEEE Trans Med Imaging*. 2014;0062(12):1–13.
- 32 Cauley SF, Setsonpomp K, Ma D, et al. Fast group matching for MR fingerprinting reconstruction. *Magn Reson Med* [Internet]. 2014;00. Available from: <http://dx.doi.org/10.1002/mrm.25439>
- 33 Hamilton JI, Jiang Y, Chen Y, et al. MRF for Rapid Quantification of Myocardial T1, T2, and Proton Spin Density. *Magn Reson Med*. 2016;In press.
- 34 Deoni SCL, Peters TM, Rutt BK. High-resolution T1 and T2 mapping of the brain in a clinically acceptable time with DESPOT1 and DESPOT2. *Magn Reson Med*. 2005;53(1):237–41.
- 35 Jamin Y, Tucker ER, Poon ES, et al. Evaluation of Clinically Translatable MR Imaging Biomarkers of Therapeutic Response in the TH- MYCN Transgenic Mouse Model of Neuroblastoma. *Radiology*. 2012.
- 36 Weidensteiner C, Allegrini PR, Sticker-Jantschkeff M, Romanet V, Ferretti S, McSheehy PM. Tumour T1 changes in vivo are highly predictive of response to chemotherapy and reflect the number of viable tumour cells - a preclinical MR study in mice. *BMC Cancer*. 2014;14(1):88.
- 37 Christen T, Pannetier NA, Ni WW, et al. MR vascular fingerprinting: A new approach to compute cerebral blood volume, mean vessel radius, and oxygenation maps in the human brain. *Neuroimage* [Internet]. 2014 Apr 1 [cited 2015 Nov 19];89:262–70. Available from: <http://www.pubmedcentral.nih.gov/articlerender.fcgi?artid=3940168&tool=pmcentrez&rendertype=abstract>
- 38 Pan S, Mao D, Peiyang L, Yang L, Babu W, Hanzhang L. Arterial Spin Labeling without control/label pairing and post-labeling delay: an MR fingerprinting implementation. In: In proceedings of the 23rd annual meeting of ISMRM meeting & exhibition in Toronto, Canada. 2015. p. 0276.

- 39 Hamilton JI, Deshmane AV, Stephanie H, Griswold MA, Seiberlich N. Magnetic Resonance Fingerprinting with Chemical Exchange (MRF-X) for Quantification of Subvoxel T1, T2, Volume Fraction, and Exchange Rate. In: In proceedings of the 23rd annual meeting of ISMRM meeting & exhibition in Toronto, Canada. 2015. p. 0329.
- 40 Vymazal J, Righini A, Brooks RA, et al. T1 and T2 in the Brain of Healthy Subjects, Patients with Parkinson Disease, and Patients with Multiple System Atrophy: Relation to Iron Content. *Radiology*. 1999;211(2):489–95.
- 41 Whittall KP, MacKay AL, Graeb DA, Nugent RA, Li DK, Paty DW. In vivo measurement of T2 distributions and water contents in normal human brain. *Magn Reson Med*. 1997;37(1):34–43.
- 42 Poon CS, Henkelman RM. Practical T2 quantitation for clinical applications. *J Magn Reson Imaging*. 1992;2(5):541–53.
- 43 Wansapura JP, Holland SK, Dunn RS, Ball WS. NMR relaxation times in the human brain at 3.0 tesla. *J Magn Reson Imaging*. 1999;9(4):531–8.
- 44 de Bazelaire CM, Duhamel GD, Rofsky NM, Alsop DC. MR imaging relaxation times of abdominal and pelvic tissues measured in vivo at 3.0 T: preliminary results. *Radiology*. 2004;230(3):652–9.
- 45 Chen Y, Lee GR, Aandal G, et al. Rapid volumetric t1 mapping of the abdomen using three-dimensional through-time spiral GRAPPA. *Magnetic Resonance in Medicine*. 2016 Apr;75(4):1457-65.
- 46 Stanisz GJ, Odobina EE, Pun J, et al. T1, T2 relaxation and magnetization transfer in tissue at 3T. *Magn Reson Med* [Internet]. 2005/08/09 ed. 2005;54(3):507–12. Available from: <http://www.ncbi.nlm.nih.gov/pubmed/16086319>.

## Contact

Vikas Gulani, M.D.  
Department of Radiology  
Case Western Reserve University  
University Hospitals Case Medical Center  
11100 Euclid Ave  
Bolwell Building, Room B120  
Cleveland, OH 44106  
USA  
[vxg46@case.edu](mailto:vxg46@case.edu)

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Otto Rapalino, M.D.<sup>1</sup>; Keith Heberlein, Ph.D.<sup>2</sup>

<sup>1</sup> Neuroradiology Division, Department of Radiology, Massachusetts General Hospital, Boston, MA, USA


<sup>2</sup> Director of Neuro Applications Development, Siemens Healthcare USA

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References

- 1
- Mehan WA, Gonzales RG, Buchbinder BR, Chen JW, Copen WA, Gupta R, et al. Optimal brain MRI protocol for new neurological complaint. PLOS One 2014; 9(10): e110803. doi:10.1371/journal.pone.0110803.
- 2
- Prakkamakul S, Witzel T, Huang S, Boulter D, Borja MJ, Schaefer P, et al. Comparative analysis of image quality and diagnostic concordance between 5-minute and conventional magnetic resonance protocols for brain imaging. Manuscript submitted for publication 2016.
- 3
- Griswold MA, Jakob PM, Heidemann RM, Nittka M, Jellus V, Wang J, et al. GeneRalized Autocalibrating Partially Parallel Acquisitions (GRAPPA). Magn Reson Med 2002; 47:1202-1210.
- 4
- Keil B, Wald LL. Massively parallel MRI detector arrays. J Magn Reson 2013; 229: 75-89.

MAGNETOM Skyra, 3T, Head/Neck 20											Total exam time: 5:11		
Plane	TR (ms)	TE (ms)	FOV (mm)	Phase FOV (%)	Slices	Slice (mm)	Gap (%)	Matrix	Phase Direc-tions	iPAT factor	b-values	Direc-tions (no)	TA (mins)
AutoAlign Head Scout													
3D	3.15	1.37	260	100	128	1.6	20	160	A-P	3	NA	NA	0:14
T1 GRE FLASH													
Sag	240	2.46	220	100	35	4.0	20	256	A-P	2	NA	NA	0:41
T2 TSE													
Axial	6200	78	220	87.5	25	5.0	20	256	R-L	3	NA	NA	1:02
T2 TSE DarkFluid													
Axial	8000	119	220	87.5	25	5.0	20	256	R-L	2	NA	NA	1:52
ep2d Diffusion													
Axial	4200	72	240	100	31	5.0	12	160	A-P	2	0,800	12	1:16
ep2d T2*-weighted													
Axial	6120	30	220	100	25	5.0	20	128	A-P	1	NA	NA	0:06

MAGNETOM Aera, 1.5T, Head/Neck 20 (XJ and XQ gradients)											Total exam time: 5:56		
Plane	TR (ms)	TE (ms)	FOV (mm)	Phase FOV (%)	Slices	Slice (mm)	Gap (%)	Matrix	Phase Direc-tions	iPAT factor	b-values	Direc-tions (no)	TA (mins)
AutoAlign Head Scout													
3D	4.52	2.38	260	100	128	1.6	20	160	A-P	3	NA	NA	0:19
T1 SE													
Sag	595	11	230	100	27	5.0	20	256	A-P	2	NA	NA	1:11
T2 TSE													
Axial	4700	101	220	87.5	25	5.0	20	256	R-L	2	NA	NA	0:56
T2 TSE DarkFluid													
Axial	5700	80	220	87.5	25	5.0	20	256	R-L	1	NA	NA	2:03
ep2d Diffusion													
Axial	4500	89	240	100	31	5.0	12	128	A-P	2	0,800	12	1:21
ep2d T2*-weighted													
Axial	6120	75	220	100	25	5.0	20	128	A-P	1	NA	NA	0:06

- 5
- Wiggins GC, Triantafyllou C, Potthast A, Reykowski A, Nittka M, Wald LL. 32-channel 3 Tesla receive-only phased-array head coil with soccer-ball element geometry. Magn Reson Med 2006; 56(1): 216-223.
- 6
- Setsompop K, Cohen-Adad J, Gagoski BA, Raij T, Yendiki A, Keil B, et al. Improving diffusion MRI using simultaneous multi-slice echo planar imaging. Neuroimage 2012; 63:569-580.
- 7
- Setsompop K, Cauley S, Wald L. (2015). Advancing diffusion MRI using SMS EPI. MAGNETOM Flash (Special SMS Supplement) 2015; 63:16-22.
- 8
- Smith SM, Beckmann CF, Andersson J, Auerbach EJ, Bijsterbosch J, Douaud G, et al. Resting state fMRI in the Human Connectome Project. Neuroimage 2013; 80:144-168.
- 9
- Ugurbil K, Auerbach EJ, Modeller S, Xu J, Vu A, Glasser MF, et al. Slice acceleration in the 3 Tesla component of the human connectome project. MAGNETOM Flash (Special SMS Supplement) 2015; 63:49-56.

Contact

Otto Rapalino, M.D.  
Neuroradiology Division  
Department of Radiology  
Massachusetts General Hospital  
55 Fruit Street  
Boston MA, 02114-2696  
USA  
Phone: +1 617-726-8320

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<sup>1</sup>Department of Radiology, China-Japan Friendship Hospital, Beijing, China

<sup>2</sup>Department of Radiology, Peking Union Medical College Hospital, Peking Union Medical College and Chinese Academy of Medical Sciences, Beijing, China

<sup>3</sup>MR Application-Predevelopment, Siemens Healthcare, Erlangen, Germany

<sup>4</sup>MR Collaboration, Siemens Healthcare, Shanghai, China

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

### References

- 1 Padhani AR, Liu G, Koh DM, et al. Diffusion-weighted magnetic resonance imaging as a cancer biomarker: consensus and recommendations. Neoplasia 2009;11(2):102-125.
- 2 Padhani AR, Koh DM, Collins DJ. Whole-body diffusion-weighted MR imaging in cancer: current status and research directions. Radiology 2011;261(3):700-718.
- 3 Attariwala R, Picker W. Whole body MRI: improved lesion detection and characterization with diffusion weighted techniques. J MagnReson Imaging 2013;38(2):253-268.
- 4 Petralia G, Padhani A, Summers P, et al. Whole-body diffusion-weighted imaging: is it all we need for detecting metastases in melanoma patients? EurRadiol. 2013;23(12):3466-3476.
- 5 Brioli A, Morgan GJ, Durie B, Zamagni E. The utility of newer imaging techniques as predictors of clinical outcomes in multiple myeloma. Expert RevHematol 2014;7(1):13-16.
- 6 Mayerhoefer ME, Karanikas G, Kletter K, et al. Evaluation of diffusion-weighted MRI for pretherapeutic assessment and staging of lymphoma: results of a prospective study in 140 patients. Clin Cancer Res 2014;20(11):2984-2993.
- 7 Littooj AS, Kwee TC, Barber I, et al. Whole-body MRI for initial staging of paediatric lymphoma: prospective comparison to an FDG-PET/CT-based reference standard. EurRadiol 2014;24(5):1153-1165.

- 8 Klenk C, Gawande R, Uslu L, et al. Ionising radiation-free whole-body MRI versus (18)F-fluorodeoxyglucose PET/CT scans for children and young adults with cancer: a prospective, non-randomised, single-centre study. *Lancet Oncol* 2014;15(3):275-285.
- 9 Zamagni E, Nanni C, Patriarca F, et al. A prospective comparison of 18F-fluoro-deoxyglucose positron emission tomography-computed tomography, magnetic resonance imaging and whole-body planar radiographs in the assessment of bone disease in newly diagnosed multiple myeloma. *Haematologica*.2007; 92(1):50-55.
- 10 Schmidt H, Brendle C, Schraml C, et al. Correlation of simultaneously acquired diffusion-weighted imaging and 2-deoxy-[18F] fluoro-2-D-glucose positron emission tomography of pulmonary lesions in a dedicated whole-body magnetic resonance/positron emission tomography system. *Invest Radiol*. 2013;48(5):247-55.
- 11 Bogner W, Pinker-Domenig K, Bickel H, et al. Readout-segmented echo-planar imaging improves the diagnostic performance of diffusion-weighted MR breast examinations at 3.0 T. *Radiology* 2012;263(1):64-76.
- 12 Porter DA, Heidemann RM. High resolution diffusion-weighted imaging using readout-segmented echo-planar imaging, parallel imaging and a two-dimensional navigator-based reacquisition. *Magn Reson Med*. 2009;62(2):468-475.
- 13 Chen NK, Guidon A, Chang HC, et al. A robust multi-shot scan strategy for high-resolution diffusion weighted MRI enabled by multiplexed sensitivity-encoding (MUSE). *Neuroimage* 2013;72:41-7.
- 14 Riffel P, Michaely HJ, Morelli JN, et al. Zoomed EPI-DWI of the head and neck with two-dimensional, spatially-selective radiofrequency excitation pulses, *EurRadiol* 2014; 24(10):2507-12.
- 15 Thierfelder KM, Scherr MK, Notohami-prodjo M, et al. Diffusion-weighted MRI of the prostate: advantages of Zoomed EPI with parallel-transmit-accelerated 2D-selective excitation imaging. *EurRadiol* 2014; 24(12):3233-41.
- 16 Morrell G, Spielman D. Dynamic shimming for multi-slice magnetic resonance imaging. *MagnResonMed* 1997;38(3):477-83.
- 17 Lee SK, Tan ET, Govenkar A, et al. Dynamic slice-dependent shim and center frequency update in 3 T breast diffusion weighted imaging. *MagnReson Med* 2014;71(5):1813-1818.
- 18 Zhang H, Xue H, Stemmer A, et al. Integrated Shimming Improves Lesion Detection in Whole-Body Diffusion-Weighted Examinations of Patients With Plasma Disorder at 3 T. *Investigative Radiology*, 2015 [Epub ahead of print]
- 19 Stemmer A and Kiefer B. Combination of integrated slice-specific dynamic shimming and pixel-wise unwarping of residual EPI distortions. *Proc. Intl. Soc. Mag. Reson. Med*. 2015; 23:3729.
- 20 Jezzard P and Balaban RS. Correction for geometric distortion in echo planar images from B<sub>0</sub> field variations. *Magnetic Resonance in Medicine* 1995; 34:65-73.
- 21 Takahara T, Imai Y, Yamashita T, et al. Diffusion weighted whole body imaging with background body signal suppression (DWIBS): technical improvement using free breathing, STIR and high resolution 3D display. *Radiat Med*. 2004;22:275–282.
- 22 Thomas C. Kwee, Taro Takahara, Reiji Ochiai, et al. Diffusion-weighted whole-body imaging with background body signal suppression (DWIBS): features and potential applications in oncology". *EurRadiol*. 2008; 18:1937–1952
- 23 Messiou C, Giles S, Collins D J, et al. Assessing response of myeloma bone disease with diffusion-weighted MRI [J]. *The British journal of radiology*, 2012, 85(1020): e1198-203.
- 24 Giles S L, Messiou C, Collins D J, et al. Whole-body diffusion-weighted MR imaging for assessment of treatment response in myeloma [J]. *Radiology*, 2014, 271(3): 785-94.

## Contact

Professor Xue Huadan, M.D.  
Department of Radiology  
Peking Union Medical College  
Hospital, Peking Union  
Medical College and Chinese  
Academy of Medical Sciences  
Shuaifuyuan 1#, Wangfujing  
Street, Dongcheng District  
Beijing, China, 100730  
bjdanna95@hotmail.com

## Contact

Zhang Haibo, M.D.  
Department of Radiology  
Chinese-Japan Friendship  
Hospital  
Yinghuayuan East Street #2  
Chaoyang District  
Beijing, China, 100029  
zhh\_hello@163.com