

Workflow syngo EVAR Guidance

Preparation of Preprocedural CT Dataset

- 1. Use preprocedural CT dataset for EVAR planning
 - ◆ Load respective CT dataset into 4D task card



♦ Select application *syngo* EVAR Guidance

2. Creation of vessel tree

- ♦ The software detects the vessel wall of the aorta and all main branching vessels. Centerlines are calculated and all identified vessels are automatically labelled, shown in the **Create Vessel Tree** blind.
- 3. Complete the vessel tree by adding additional vessels



- ◆ Activate Add vessel
- ♦ Use 1-click segmentation by clicking on the most distal end of the branching vessel in one of the MPR segments.

 The vessel segments will be automatically detected.
- ♦ Rename the vessels added in the input field








4. Verification of landmarks

♦ Select the **Define Landmarks** tab card



- Switch along landmarks to verify the perpendicular views
- ♦ If needed, adjust the landmark position in the MPR segment directly by dragging to the correct position




2D/3D Registration

- 1. Start syngo 2D/3D fusion
 - ♦ Select the **Overlay** tab card
 - ♦ Perform two fluoro projections in the frontal and lateral positions

2. Alignment of the fluoro projections with the CTA volume

- ♦ Align the volume with two 2D images by moving or rotating the volume in the side-by-side layout
- ♦ Accept 2D/3D alignment




Fusion Imaging

1. Overlay of 3D contours




◆ To display contour information from the 3D volume on the live fluoroscopy image, select the **Overlay** tab card and activate the **Live** button




2. Choose the target vessel from the heads-up display

- Use the heads-up display menu and select the VRT onscreen menu and the submenu Adjust to landmark
- Select the required vessel. The corresponding projection angulation is sent automatically to the C-arm.
- ♦ Deflect the joystick to drive the C-arm to the defined optimal viewing angle without any radiation



- 3. Choose the target vessel from the syngo X-Workplace
 - ♦ Select the **Overlay** tab card
 - ◆ Choose the required vessel/landmark from the list
 - ♦ Deflect the joystick to drive the C-arm to the defined optimal viewing angle without any radiation



3D Assessment

- 1. Perform syngo DynaCT for immediate 3D assessment of the stent position
- 2. In case of endoleak type II detection, use the DynaCT data for treatment planning with syngo Needle Guidance