
## **Study Protocol**

# Visualization of prostate arteries

supported by Artis Q and syngo Embolization Guidance

Vascular



### **Case Description**

#### Patient history

A 50-year-old man with benign prostate hyperplasia and acute urinary retention was referred for prostate artery embolization. He had been experiencing lower urinary tract symptoms for nearly three years and had been treated medically with no clinical improvement.

#### Diagnosis

Benign prostate hyperplasia and acute urinary retention.

#### Treatment

With the catheter positioned in the abdominal aorta above the bifurcation, we got perfect visualization of the prostate and supplying vessels from both the left and right side. With the help of syngo Embolization Guidance we obtained a map of supplying vessel that was used during embolization. Successful vascular occlusion was then accomplished with microspheres (300–500 mm in size) for embolization.

#### **General comments**

syngo Embolization Guidance with arterial contrast injection facilitates selective prostate artery catheterization, which has the potential to reduce dose and contrast medium. We marked the distal point of the prostatic artery and the proximal internal artery. syngo Embolization Guidance automatically creates the track between the two points, which can be overlaid to the live fluoro images.

#### Tips & Tricks:

First, the power injector was readied. An injection protocol with an X-ray delay was chosen to get good visualization of the iliac arteries and prostate. 6s DynaCT Body acquisition protocol was used. The injector is triggered by the system and, according to the injection protocol, the C-arm starts DynaCT rotation with a 7 s delay.

#### Courtesy of

Bulat Sharafutdinov, MD, Edgar Gaziev, MD, Clinics of Kazan University, Kazan city, Russia

#### Supported by

- syngo Embolization Guidance
- syngo DynaCT
- syngo 3D Roadmap

#### System & Software

Artis Q floor VD11 syngo X Workplace VD10




# Protocol


Acquisition protocol	6sDCT Body
Injection protocol	
Catheter position	Pigtail 5F; abdominal aorta above bifurcation
Contrast medium (CM)	350 mg iodine/mL
Dilution	No
Injection volume	65 mL
Injection rate	5 mL/s
Duration of injection	13 s
X-ray delay	7 s
Power injector used	Yes


Reconstructions	1. Reconstruction
Name	DynaCT Body NatFill HU Normal HU
VOI size	Full
Slice matrix	512x512
Kernel type	HU
Image characteristics	Auto
Reconstruction mode	NatFill
Viewing preset	DynaCT Body

# Clinical Images












Fig. 2 a, b. 75 mm MIPs (left and right side) with vessel trajectory identified by *syngo* Embolization Guidance.

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.

All rights reserved.

#### Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany

Phone: +49 9131 84-0 siemens-healthineers.com