SOMATOM On.site

Bringing critical care imaging to your patient

Focus: Intensive Care Units

siemens-healthineers.com/somatom-onsite



International version.

Not for distribution or use in the U.S.

Product does not yet fulfill neither MDD (93/42/EEC) nor MDR 2017/745 and is therefore not yet available in the EU.





Providing CT scans of the head for ICU patients is crucial but often challenging

ICU (Intenesive Care Unit) patients suffering from acute and critical head conditions face a high risk of sudden health status deterioration. When this happens, fast evaluation is essential for immediate treatment decision-making. CT scans are ideally suited for this task – but providing imaging for ICU patients is often cumbersome and not without risk.



CT scans of the head are in high demand

Head scans are among the most common procedures for ICU patients and account for 30% of all CT scans.¹



ICU patient transport involves logistical and staffing challenges

Transporting patients from the ICU to the radiology department is a delicate and complex process that requires up to five staff members² to leave their other duties behind.



ICU staff shortage limits patient care

Frequent bed closures affect 76% of German ICU floors, mostly due to staff shortage.³



Transport may negatively affect the patient's health status

Providing CT imaging requires transporting ICU patients to the radiology department. However, for 1 in 5 patients, transport makes existing health conditions worse.⁴



ICU patient scans complicate radiology department schedules

Patients with acute and critical conditions tie up CT scanners for at least twice as long as most other patients. ^{5,6} Transports are also frequently delayed, ² making it necessary to reschedule other patients.





Bringing critical care imaging to your patient

Our mobile head CT scanner SOMATOM On.site provides access to consistent and reliable SOMATOM image quality, right at your patient's bedside. Thanks to the latest advances in digital healthcare, the intelligent software myExam Companion together with our GO technologies guides users through the scanning process. With this workflow guidance the system is easy to use and it enables greater consistency in image acquisition across your team.

SOMATOM On.site will change the way you care for ICU patients suffering from acute and critical head conditions. Eliminating transport to the radiology department may help prevent related complications and allows your staff to concentrate on what matters most: providing optimal care for every patient at the ward. With this easy access to CT imaging of the head directly at the point of care, you can get images to help diagnose and treat sudden changes in the patient's condition or arising complications faster. Additionally, offering CT imaging directly in the ICU allows the radiology department to maintain their scheduling.

With SOMATOM On.site, you can transform care delivery.



See for yourself

Our workflow video demonstrates how point-of-care CT scanning works with SOMATOM On.site.



Transform care delivery with SOMATOM On.site

Minimize patient transports

SOMATOM On.site permits CT imaging of the head directly in the ICU ward, so patients need not encounter any transport-associated risks. Acquiring images directly at the patient's bedside may speed up diagnostics and following decisions and steps. This may have an impact on hospital length of stay and infection control.

Optimize staff deployment

SOMATOM On.site reduces the need for transporting patients to the radiology department for CT scans of the head. ICU staff and doctors can stay in the ICU. This means they can gain additional time for providing bedside patient care and focus more on their core tasks.

Achieve high image quality at the point of care

SOMATOM On.site brings SOMATOM image quality to the patient's bedside. It is combining proven technologies from existing SOMATOM scanners with newly designed components.

Point-of-care imaging reduces the number of ICU patient scans that need to be performed in the radiology department, potentially resulting in more streamlined utilization of stationary systems and staff, along with more predictable radiology appointment schedules.





Experience convenient and reliable bedside imaging

SOMATOM On.site simplifies your entire imaging workflow, enabling you to provide access to SOMATOM image quality for critically ill patients directly in the ICU ward.

Real mobility for convenient driving and maneuvering

Thanks to the intelligent driving concept, SOMATOM On.site helps you provide imaging to patients at virtually any location. The motorized trolley features four driving wheels along with an ergonomic drive handle that enable easy maneuverability even in small patient rooms.

For safe handling in busy environments, SOMATOM On.site is equipped with an integrated drive camera that guides the user while driving. The integrated computers lender further flexibility for patient scans in small rooms where additional trolleys would be a hindrance.

Eliminating common patient positioning challenges

SOMATOM On.site features a telescopic gantry that remains retracted until you are ready to start scanning. This gives you ample space for patient positioning with tubes and lines in place. The scanner also has an integrated shoulder board for bridging the gap between the bed and scanner. Sliding the patient's upper body onto the shoulder board permits convenient patient positioning in preparation for scans. Furthermore, the integrated head holder of SOMATOM On.site helps for correct isocentric positioning of the patient's head inside the gantry bore.



Simple steps suitable for cross-trained staff



Select your patient via the intuitive user interface of SOMATOM On.site.



Choose your preferred scan protocol for image acquisition.



Scan&GO allows you to supervise the scan right at your patient's bedside.



Recon&GO handles all postprocessing for you and automatically uploads the images to the PACS.

Fast and intuitive scan preparation

SOMATOM On.site leverages key digital health-care innovations to unleash the full potential of modern imaging technology – independent of the operator's experience. The scanner offers an intuitive user interface and, thanks to myExam Companion, is equipped with builtin expertise to guide you through the complete imaging examination.

In addition SOMATOM On.site features the proven GO technologies. Scan&GO enables techs to control the entire scan process via the integrated Touch UI and stay right next to the patient throughout the image acquisition process. A dedicated Recon&GO functionality performs automated postprocessing and uploads the images to PACS without any user interaction. As a result, radiologists receive comprehensive results and consistent image quality. Together, these features make scanning intuitive and easy.



Safe scanning in ICU patient rooms

Traditional stationary CT scanners require techs to return to the control room during scans, because a safety distance is necessary for radiation protection. SOMATOM On.site allows you to stay close to patients and remain at their side.

The self-shielded scanner design, along with the detachable radiation shields* covering the front and back of the gantry bore, protect staff and neighboring patients from scattered radiation.





Scan start position

Scan end position

Telescopic gantry for hassle-free acquisitions

SOMATOM On.site features a specially designed telescopic gantry with two gantry compartments and an internally moving tube and detector. Thanks to this technology, the systems' trolley position remains fixed throughout the scan, eliminating common problems associated with scanners that may move during image acquisition. The many lines and drains typically in place for ICU patients remain unaffected while scanning with SOMATOM On.site, as the gantry cover close to the patient head stays fixed throughout the acquisition. Taking scanner, table, and patient movement out of the equation means that you do not have to make any compromises when it comes to image quality.

Technical specifications

Key data

Slices	32
Z-coverage	24 mm (32 x 0.75 mm)
Detector material	Stellar detector UFC (Ultra Fast Ceramics)
Rotation time	1.0 s
Max. scan speed	36 mm/s
Spatial resolution	17 lp/cm

kV settings	80–120 kV	
Current	up to 25 mA	
Scanner dimensions (height / width / depth)	1.55 / 1.60 / 0.74 m	
Iterative reconstruction	SAFIRE, iMAR	
Scanner weight	~ 890 kg	
Motorized trolley with integrated drive camera		

^{*} Radiation shields are currently under development and their future availability cannot be guaranteed.

CT imaging performance

SOMATOM On.site delivers reliable and consistent image quality with the advantages of mobile CT scanning thanks to a combination of existing innovations from SOMATOM scanners (e.g., low-noise detectors) and newly designed components (e.g., the tube).



Stellar detector

SOMATOM On.site is the first mobile CT system outfitted with our SOMATOM Stellar detector. Stellar detector's technology helps to keep electronic noise low, increases dose efficiency, and improves spatial resolution. It's, high channel density and integrated design of the detector results in the grey-white matter differentiation that is crucial for evaluating CT scans of the head.



Dedicated X-ray tube

Our X-ray tube was designed specifically for neuro CT imaging. It features the resistance to vibrations and shock that is needed in mobile CT scanners. It is built to reduce maintenance downtime and improve performance in demanding environments. The cooling system is a special design combining water and air.


State-of-the-art image reconstruction algorithms

SOMATOM On.site incorporates iterative reconstruction algorithms to help you get the most out of your point-of-care CT scan.



Metal artifact reduction with iMAR


iMAR (iterative metal artifact reduction) allows for reconstructed artifact reduction without losing valuable information. iMAR is designed to yield images with a reduced level of metal artifacts compared to conventional reconstruction if the underlying CT data is distorted by the presence of metal in the scanned object.




Noise removal with SAFIRE

This iterative reconstruction algorithm delivers natural looking images at the lowest possible dose. Multiple iterations are performed when processing the raw data. During these iterations the noise is reduced without compromising image quality.

Imaging at the patient's bedside



Apply known tools and applications



Maximum intensity projection (MIP), e.g., for intervention planning



Cinematic Rendering to enhance education and communication

This product does not yet fulfill all the essential requirements according to the European Medical Device Directive (93/42/EEC) and its national implementations or all applicable General Safety and Performance Requirements according to the European Medical Device Regulation 2017/745. It is not yet commercially available in the European Union.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide.

Availability and packaging may vary by country and is subject to change without prior notice. Some or all of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

- 1 Data on file
- 2 Barnes E. Portable CT brings lower costs, reduce risk to ICU. AuntMinnie.com; 2008.
- 3 Karagiannidis C, Kluge S, Riessen R, Krakau M, Bein T, Janssens U. Auswirkungen des Pflegepersonalmangels auf die intensivmedizinische Versorgungskapazität in Deutschland. Med Klin Intensivmed Notfmed. 2018 Jul 9;114(4):327–33. Available from: https://www.springermedizin.de/auswirkungen-des-pflegepersonalmangels-auf-die-intensivmedizinis/15927204.
- 4 Anticipating unique disruptions. PwC's 21st CEO Survey: key findings from the healthcare industry [homepage on the internet]. PwC Turkey; 2018 [cited 2019 Aug 6]. Available from: https://www.pwc.com.tr/en/yayinlar/arastirmalar/ceo-survey/21-ceo-arastirmasi/saglik-sektoru-sonuclari.html.
- 5 Agrawal, D et al. (2016). Bedside computed tomography in traumatic brain injury: Experience of 10,000 consecutive cases in neurosurgery at a level 1 trauma center in India. Neurol India. 2016 Jan–Feb; 64(1):62– 365. Available from: https://www.ncbi.nlm.nih.gov/pubmed/26754994.
- 6 Masaryk et al. (2008). The economic and clinical benefits of portable head/neck CT imaging in the intensive care unit. Radiol Manage. 2008 March–April; 30(2):50–54. Available from: https://www.ncbi.nlm.nih. qov/pubmed/18431942.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone: +49 9131 84-0

Phone: +49 9131 84-0 siemens-healthineers.com