



Transducers

ACUSON Sequoia ultrasound system

Version VA11 SW

siemens-healthineers.com/ultrasound



Taking ultrasound to new heights

ACUSON Sequoia offers 11 transducers leveraging Siemens Healthineers unique new InTune transducer design and architecture for optimal performance.



Contents

| | |
|--------------------------------------|----|
| Curved | 3 |
| Linear | 4 |
| Endocavity | 6 |
| Vector | 7 |
| Pencil | 9 |
| Selectable Frequencies | 10 |
| Cable Length | 11 |
| Connector Type | 11 |
| Needle Guide | 12 |
| Fusion – eTRAX Needle Tracking | 13 |
| Advanced Applications | 13 |

Curved



9C3 Transducer

| | |
|------------------------------|--------------------------|
| Form factor | Curved |
| Design | 1D, Hanafy, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 2.2–9.2 MHz |
| Axial and lateral resolution | 0.56 and 0.96 mm |
| Maximum depth | Up to 30 cm |
| Field of view | 79 deg |
| Physical footprint | 69.6 x 20.5 mm |
| Total weight | 780 g |



5C1 Transducer

| | |
|------------------------------|--------------------|
| Form factor | Curved |
| Design | 1D, Single Crystal |
| Gesture detection | Yes |
| Bandwidth | 1.0–5.7 MHz |
| Axial and lateral resolution | 0.67 and 1.2 mm |
| Maximum depth | Up to 30 cm |
| Field of view | 72 deg |
| Physical footprint | 63.3 x 18.2 mm |
| Total weight | 743 g |



DAX Transducer

| | |
|------------------------------|-----------------------|
| Form factor | Curved |
| Design | Multi-D, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 1.0–3.5 MHz |
| Axial and lateral resolution | 0.8 and 2.3 mm |
| Maximum depth | Up to 40 cm |
| Field of view | 50 deg |
| Physical footprint | 57.7 x 30.2 mm |
| Total weight | 848 g |

Linear



18H6 Transducer

| | |
|------------------------------|------------------|
| Form factor | Linear |
| Design | 1D, Piezoceramic |
| Gesture detection | No |
| Bandwidth | 5.5–21.10 MHz |
| Axial and lateral resolution | 0.2 and 0.23 mm |
| Maximum depth | Up to 6 cm |
| Field of view | 28 mm |
| Physical footprint | 13.6 x 40.4 mm |
| Total weight | 630 g |



18L6 Transducer

| | |
|------------------------------|--------------------------|
| Form factor | Linear |
| Design | 1D, Hanafy, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 4.6–17.8 MHz |
| Axial and lateral resolution | 0.3 and 0.43 mm |
| Maximum depth | Up to 8 cm |
| Field of view | 58 mm |
| Physical footprint | 69.2 x 16.5 mm |
| Total weight | 762 g |



14L5 Transducer

| | |
|------------------------------|-----------------------|
| Form factor | Linear |
| Design | Multi-D, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 4.8–13.6 MHz |
| Axial and lateral resolution | 0.3 and 0.38 mm |
| Maximum depth | Up to 8 cm |
| Field of view | 38 mm |
| Physical footprint | 49.6 x 12.9 mm |
| Total weight | 727 g |



10L4 Transducer

| | |
|------------------------------|-----------------------|
| Form factor | Linear |
| Design | Multi-D, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 2.9–9.9 MHz |
| Axial and lateral resolution | 0.3 and 0.52 mm |
| Maximum depth | Up to 14 cm |
| Field of view | 38 mm |
| Physical footprint | 49.3 x 18.9 mm |
| Total weight | 723 g |

Endocavity



9EC4 Transducer

| | |
|------------------------------|------------------|
| Form factor | Curved |
| Design | 1D, Piezoceramic |
| Gesture detection | Yes |
| Bandwidth | 2.9–8.1 MHz |
| Axial and lateral resolution | 0.46 and 0.8 mm |
| Maximum depth | Up to 14 cm |
| Field of view | 176 deg |
| Physical footprint | 17.0 x 22.0 mm |
| Total weight | 700 g |

Vector



10V4 Transducer

| | |
|------------------------------|--------------------------|
| Form factor | Vector |
| Design | 1D, Hanafy, Piezoceramic |
| Gesture detection | N/A |
| Bandwidth | 3.4–10.4 MHz |
| Axial and lateral resolution | 0.34 and 0.62 mm |
| Maximum depth | Up to 14 cm |
| Field of view | 90 deg |
| Physical footprint | 22.6 x 14.3 mm |
| Total weight | 585 g |



8V3 Transducer

| | |
|------------------------------|--------------------------|
| Form factor | Vector |
| Design | 1D, Hanafy, Piezoceramic |
| Gesture detection | N/A |
| Bandwidth | 2.1–8.3 MHz |
| Axial and lateral resolution | 0.59 and 0.79 mm |
| Maximum depth | Up to 24 cm |
| Field of view | 90 deg |
| Physical footprint | 26.9 x 16.6 mm |
| Total weight | 644 g |



5V1 Transducer

| | |
|------------------------------|--------------------|
| Form factor | Vector |
| Design | 1D, Single Crystal |
| Gesture detection | Yes |
| Bandwidth | 1.1–4.9 MHz |
| Axial and lateral resolution | 1.06 and 3.72 mm |
| Maximum depth | Up to 30 cm |
| Field of view | 90 deg |
| Physical footprint | 27.2 x 18.7 mm |
| Total weight | 640 g |



4V1 Transducer

| | |
|------------------------------|--------------------------|
| Form factor | Vector |
| Design | 1D, Hanafy, Piezoceramic |
| Gesture detection | N/A |
| Bandwidth | 1.4–5.1 MHz |
| Axial and lateral resolution | 0.9 and 1.1 mm |
| Maximum depth | Up to 30 cm |
| Field of view | 90 deg |
| Physical footprint | 35.5 x 20.2 mm |
| Total weight | 639 g |

Pencil



CW2 Transducer

| | |
|------------------------------|---------|
| Form factor | Pencil |
| Design | N/A |
| Gesture detection | N/A |
| Bandwidth | N/A |
| Axial and lateral resolution | N/A |
| Maximum depth | N/A |
| Field of view | N/A |
| Physical footprint | 17.0 mm |
| Total weight | N/A |

Table 1: Selectable Frequencies¹

| Transducer | 2D | THI | Color Doppler | PW Doppler | Contrast |
|------------|---------------------|---------------------|--------------------------|----------------|----------------|
| 9C3 | Low, Mid, High | Low, Mid, High | Pen, Mid, Res | Low, Mid, High | Low, Mid |
| 5C1 | Pen, Low, Mid, High | Pen, Low, Mid, High | Pen, Low, Mid, High, Res | Low, Mid | Low, Mid, High |
| DAX | Pen, Low, Mid | Pen, Low, Mid | Pen, Mid, Res | Pen, Low | Pen, Mid, High |
| 18H6 | Low, Mid, High | High | Pen, Mid, Res | Low, Mid, High | Low, Mid |
| 18L6 | Low, Mid, High | Low, Mid, High | Pen, Mid, Res | Low, Mid, High | Low, Mid |
| 14L5 | Low, Mid, High | Low, Mid, High | Pen, Mid, Res | Low, Mid | Low, Mid |
| 10L4 | Low, Mid, High | Low, Mid, High | Pen, Mid, High, Res | Low, Mid | Low, Mid |
| 9EC4 | Low, Mid, High | Low, Mid, High | Pen, Mid, Res | Low, Mid, High | Low, Mid |
| 10V4 | Low, Mid, High, Res | Mid, High | Pen, Mid, High, Res | Low, Mid, High | |
| 8V3 | Low, Mid, High, Res | Mid, High | Pen, Mid, High, Res | Low, Mid, High | |
| 5V1 | Pen | Low, Mid, High | Pen, Mid, Res | Low | |
| 4V1 | Low, Mid, High | Low, Mid, High | Pen, Mid, Res | Low, Mid | Low, Mid, High |

¹ System specific

Table 2: Cable Length

| Transducer | Cable Length |
|------------|--------------|
| 9C3 | 2.1 m |
| 5C1 | 2.1 m |
| DAX | 2.7 m |
| 18L6 | 2.1 m |
| 14L5 | 2.1 m |
| 10L4 | 2.1 m |
| 9EC4 | 2.2 m |
| 8V3 | 2.2 m |
| 5V1 | 2.1 m |
| 4V1 | 1.9 m |

Table 3: Connector Type

| Transducer | Connector |
|------------|---------------------------|
| 9C3 | Compact pinless connector |
| 5C1 | Compact pinless connector |
| DAX | Compact pinless connector |
| 18L6 | Compact pinless connector |
| 14L5 | Compact pinless connector |
| 10L4 | Compact pinless connector |
| 9EC4 | Compact pinless connector |
| 8V3 | Compact pinless connector |
| 5V1 | Compact pinless connector |
| 4V1 | Compact pinless connector |
| CW2 | Hirose |

Table 4: Needle Guide

| Transducer | Product Description | Guidance Angle Selection – Depth |
|------------|---|---|
| 9C3 | Ultra-Pro II Bracket Stater Kit | A – 5 cm B – 10 cm |
| 5C1 | Verza Tracking Bracket Starter Kit | 1 – 2.2 cm 2 – 3.8 cm 3 – 6.1 cm 4 – 9.9 cm 5 – 15.0 cm |
| DAX | Verza Tracking Bracket Starter Kit | 1 – 2.4 cm 2 – 4.1 cm 3 – 6.4 cm 4 – 9.9 cm 5 – 15 cm |
| 18L6 | Ultra-Pro II Bracket Starter Kit | A – 2.1 cm B – 5.4 cm |
| 14L5 | Verza Bracket Starter Kit | 1 – 1.8 cm 2 – 3.0 cm 3 – 4.3 cm 4 – 6.4 cm 5 – 8.9 cm |
| 10L4 | Verza Tracking Bracket Starter Kit | 1 – 2.2 cm 2 – 3.6 cm 3 – 5.6 cm 4 – 8.6 cm 5 – 13 cm |
| 9EC4 | Disposable Endocavity Guide Kit – 24 pack | 1° Needle Path angle |
| 9EC4 | Reusable Endocavity Guide | 1° Needle Path angle |
| 4V1 | Ultra-Pro II Tracking Bracket Starter Kit | A – 5 cm B – 10 cm |

Table 5: Fusion – eTRAX Needle Tracking

Product Description

| |
|------------------------|
| 12GA eTRAX Starter Kit |
| 14GA eTRAX Starter Kit |
| 16GA eTRAX Starter Kit |
| 18GA eTRAX Starter Kit |

Table 6: Advanced Applications

| Transducer | Strain Elastography | Point Shear Wave Elastography | 2D Shear Wave Elastography | Contrast Imaging | Fusion Imaging |
|------------|---------------------|-------------------------------|----------------------------|------------------|----------------|
| 9C3 | N/A | N/A | N/A | Yes | N/A |
| 5C1 | N/A | Yes | Yes | Yes | Yes |
| DAX | N/A | Yes | Yes | Yes | Yes |
| 18L6 | Yes | N/A | N/A | Yes | N/A |
| 14L5 | Yes | N/A | N/A | Yes | N/A |
| 10L4 | Yes | Yes | Yes | Yes | Yes |
| 9EC4 | Yes | N/A | N/A | Yes | N/A |
| 8V3 | N/A | N/A | N/A | N/A | N/A |
| 5V1 | N/A | N/A | N/A | N/A | N/A |
| 4V1 | N/A | Yes | N/A | Yes | Yes |

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons, their future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

ACUSON Sequoia and InTune are trademarks of Siemens Medical Solutions USA, Inc.

eTRAX is a trademark of CIVCO. CIVCO is a registered trademark of CIVCO Medical Solutions.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey towards expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated 5 million patients globally everyday benefit from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 170 years of experience and 18,000 patents globally. With more than 48,000 dedicated colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029, USA
Phone: 1-888-826-9702
siemens-healthineers.com/ultrasound