



DICOM Conformance Statement

syngo.via View&GO
VA15A

CONFORMANCE STATEMENT OVERVIEW

The **syngo.via View&GO** is lightweight reading software solution with connectivity to DICOM modalities and healthcare information systems. By default one **syngo.via View&GO (AE)** is used. It is possible to configure usage of multiple different AEs for the individual DICOM services.

The **syngo.via View&GO**:

- requests storage of objects (images, reports, encapsulated PDF)
- support query and retrieval of objects from a remote node
- displays images to a user
- imports objects from portable interchange media
- exports objects to non-optical storage device (e.g. USB stick)

The **syngo.via View&GO** conforms to the DICOM Standard and supports the network services as described in Table 0-1 and the media services as described in Table 0-2.

Table 0-1 Network Services

| SOP Classes | Service Class User (SCU) | Service Class Provider (SCP) |
|--|--------------------------|------------------------------|
| Verification | | |
| Verification Service | Yes | Yes |
| Transfer(Image SOP Class) | | |
| Breast Tomosynthesis Image Storage | Yes | Yes |
| Computed Radiography Image Storage | Yes | Yes |
| CT Image Storage | Yes | Yes |
| Digital Mammography X-Ray Image Storage - For Presentation | Yes | Yes |
| Digital Mammography X-Ray Image Storage - For Processing | Yes | Yes |
| Digital X-Ray Image Storage - For Presentation | Yes | Yes |
| Digital X-Ray Image Storage - For Processing | Yes | Yes |
| Enhanced CT Image Storage | Yes | Yes |
| Enhanced MR Image Storage | Yes | Yes |
| Enhanced MR Color Image Storage | Yes | Yes |
| Enhanced XA Image Storage | Yes | Yes |
| Enhanced XRF Image Storage | Yes | Yes |
| MR Image Storage | Yes | Yes |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | Yes | Yes |
| Multi-frame Grayscale Word Secondary Capture Image Storage | Yes | Yes |
| Multi-frame Single Bit Secondary Capture Image Storage | Yes | Yes |
| Multi-frame True Color Secondary Capture Image Storage | Yes | Yes |

| SOP Classes | Service Class User (SCU) | Service Class Provider (SCP) |
|--|--------------------------|------------------------------|
| Nuclear Medicine Image Storage | Yes | Yes |
| PET Image Storage | Yes | Yes |
| RT Image Storage | Yes | Yes |
| Secondary Capture Image Storage | Yes | Yes |
| Segmentation Storage | Yes | Yes |
| Ultrasound Image Storage | Yes | Yes |
| Ultrasound Image Storage (Retired) | Yes | Yes |
| Ultrasound Multi-Frame Image Storage | Yes | Yes |
| Ultrasound Multi-Frame Image Storage (Retired) | Yes | Yes |
| X-Ray 3D Angiographic Image Storage | Yes | Yes |
| X-Ray Angiographic Image Storage | Yes | Yes |
| X-Ray Radio-Fluoroscopic Image Storage | Yes | Yes |
| Transfer(Non-Image SOP Class) | | |
| 12-lead ECG Waveform Storage | Yes | Yes |
| Ambulatory ECG Waveform Storage | Yes | Yes |
| Basic Text SR Storage | Yes | Yes |
| Blending Softcopy Presentation State Storage | Yes | Yes |
| Cardiac Electrophysiology Waveform Storage | Yes | Yes |
| Color Softcopy Presentation State Storage (store & forward only) | Yes | Yes |
| Comprehensive SR Storage | Yes | Yes |
| Deformable Spatial Registration Storage | Yes | Yes |
| Encapsulated PDF Storage | Yes | Yes |
| Enhanced SR Storage | Yes | Yes |
| General ECG Waveform Storage | Yes | Yes |
| Grayscale Softcopy Presentation State Storage (store & forward only) | Yes | Yes |
| Hemodynamic Waveform Storage | Yes | Yes |
| Key Object Selection Document Storage | Yes | Yes |
| Mammography CAD SR Storage | Yes | Yes |
| MR Spectroscopy Storage | Yes | Yes |
| Pseudo-Color Softcopy Presentation State Storage | Yes | Yes |
| Procedure Log Storage | Yes | Yes |
| Raw DataStorage | Yes | Yes |
| Real World Value Mapping Storage | Yes | Yes |
| RT Beams Treatment Record Storage | Yes | Yes |
| RT Dose Storage | Yes | Yes |
| RT Ion Beams Treatment Record Storage | Yes | Yes |
| RT Ion Plan Storage | Yes | Yes |
| RT Plan Storage | Yes | Yes |
| RT Structure Set Storage | Yes | Yes |
| RT Treatment Summary Record Storage | Yes | Yes |
| Spatial Fiducials Storage | Yes | Yes |
| Spatial Registration Storage | Yes | Yes |
| Surface Segmentation Storage | Yes | Yes |

| SOP Classes | Service Class User (SCU) | Service Class Provider (SCP) |
|--|--------------------------|------------------------------|
| X-Ray Radiation Dose SR Storage | Yes | Yes |
| Transfer(Private SOP Class) | | |
| CSA Non-Image Storage | Yes | Yes |
| Query / Retrieve | | |
| Patient Root – Query/Retrieve Information Model – FIND | Yes | No |
| Patient Root – Query/Retrieve Information Model – MOVE | Yes | No |
| Study Root – Query/Retrieve Information Model – FIND | Yes | No |
| Study Root – Query/Retrieve Information Model – MOVE | Yes | No |
| Patient/Study Only – Query/Retrieve Information Model – FIND | Yes | No |
| Patient/Study Only – Query/Retrieve Information Model – MOVE | Yes | No |
| Workflow Management | | |
| Storage Commitment Push Model | No | No |
| Modality Worklist Information Model – FIND | No | No |
| Print Management | | |
| Basic Grayscale Print Management Meta SOP Class | No | No |
| Basic Film Session SOP Class | No | No |
| Basic Film Box SOP Class | No | No |
| Basic Grayscale Image Box SOP Class | No | No |
| Printer SOP Class | No | No |
| Print Job SOP Class | No | No |
| Presentation LUT SOP Class | No | No |
| Basic Color Print Management Meta SOP Class | No | No |
| Basic Color Image Box SOP Class | No | No |

Table 0-2 Media Services

| Media Storage Application Profile | Write Files (FSC) | Read Files (FSR) |
|--|-------------------|------------------|
| Compact Disk – Recordable | | |
| STD-GEN-CD (augmented, see 3.2.1) | No | Yes |
| DVD – Recordable | | |
| STD-GEN-DVD (augmented, see 3.2.1) | No | Yes |
| STD-GEN-DVD-J2K (augmented, see 3.2.1) | No | Yes |
| USB | | |
| STD-GEN-USB-J2K (augmented, see 3.2.1) | Yes | Yes |

TABLE OF CONTENTS

| | |
|--|-----------|
| CONFORMANCE STATEMENT OVERVIEW | 3 |
| TABLE OF CONTENTS | 7 |
| List of Tables | 9 |
| List of Figures | 10 |
| 1 INTRODUCTION | 11 |
| 1.1 REVISION HISTORY | 11 |
| 1.2 GENERAL | 11 |
| 1.3 AUDIENCE | 11 |
| 1.4 REMARKS | 11 |
| 1.5 TERMS AND DEFINITIONS | 11 |
| 1.6 ABBREVIATIONS | 12 |
| 1.7 REFERENCES | 12 |
| 1.8 SCOPE AND FIELD OF APPLICATION | 13 |
| 2 NETWORKING | 14 |
| 2.1 IMPLEMENTATION MODEL | 14 |
| 2.1.1 Application Data Flow | 14 |
| 2.1.2 Functional Definition of AE's | 14 |
| 2.1.3 Sequencing of Real-World Activities | 16 |
| 2.2 AE SPECIFICATIONS | 18 |
| 2.2.1 syngo.via View&GO AE | 18 |
| 2.3 NETWORK INTERFACES | 31 |
| 2.3.1 Physical Network Interface | 31 |
| 2.3.2 Additional Protocols | 31 |
| 2.3.3 IPv4 and IPv6 Support | 31 |
| 2.4 CONFIGURATION | 31 |
| 2.4.1 AE Title/Presentation Address Mapping | 31 |
| 3 MEDIA INTERCHANGE | 33 |
| 3.1 IMPLEMENTATION MODELS | 33 |
| 3.1.1 Application Data Flow Diagram | 33 |
| 3.1.2 Functional definitions of AEs | 34 |
| 3.1.3 Sequencing of Real-World Activities | 34 |
| 3.1.4 File Meta Information for Implementation Class and Version | 35 |
| 3.2 AE SPECIFICATIONS | 35 |
| 3.2.1 Media Storage AE – Specification | 35 |
| 3.3 AUGMENTED AND PRIVATE APPLICATION PROFILES | 45 |
| 3.3.1 Augmented Application Profiles | 45 |
| 3.4 MEDIA CONFIGURATION | 46 |
| 3.5 Attribute confidentiality profiles | 46 |
| 3.5.1 De-identification | 46 |
| 4 SUPPORT OF CHARACTER SETS | 51 |
| 4.1 CHARACTER SETS FOR syngo.via View&GO | 51 |

| | | |
|----------|---|-----------|
| 5 | SECURITY | 54 |
| 5.1 | SECURITY PROFILES | 54 |
| 5.2 | ASSOCIATION LEVEL SECURITY | 54 |
| 6 | Annexes | 55 |
| 6.1 | SOP Classes supported | 55 |
| 6.2 | IOD CONTENTS | 58 |
| 6.2.1 | Created SOP Instance(s) | 58 |
| 6.2.2 | Data Dictionary of Private Attributes | 58 |
| 6.2.3 | Usage of Attributes from received IODs | 59 |
| 6.2.4 | Coerced / Modified fields | 59 |
| 6.3 | CODED TERMINOLOGY AND TEMPLATES | 59 |
| 6.3.1 | Context Groups | 59 |
| 6.3.2 | Template Specifications | 59 |
| 6.3.3 | Private Code definitions | 60 |
| 6.4 | GRAYSCALE IMAGE CONSISTENCY | 60 |
| 6.5 | STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES | 60 |

List of Tables

| | |
|--|----|
| Table 0-1 Network Services | 3 |
| Table 0-2 Media Services | 6 |
| Table 2-1: Association Policies | 19 |
| Table 2-2: Asynchronous Nature as an Association Initiator | 19 |
| Table 2-3: DICOM Implementation Class and Version | 19 |
| Table 2-4: Association initiation policies | 19 |
| Table 2-5: Acceptable Presentation Contexts for Verification Activity | 20 |
| Table 2-6: Proposed Presentation Contexts for Storage | 21 |
| Table 2-7: DICOM Command Response Status Handling Behavior | 22 |
| Table 2-8: DICOM Command Communication Failure Behavior | 22 |
| Table 2-9: Proposed Presentation Contexts for Query | 24 |
| Table 2-10: Extended Negotiation as an SCU | 25 |
| Table 2-11: DICOM Command Response Status Handling Behavior | 25 |
| Table 2-12: DICOM Command Communication Failure Behavior | 25 |
| Table 2-13: Attributes supported for instance Query - SCU | 26 |
| Table 2-14: Proposed Presentation Contexts for Retrieve and Activity "MOVE SCU" | 27 |
| Table 2-15: DICOM Command Response Status Handling Behavior | 27 |
| Table 2-16: DICOM Command Communication Failure Behavior | 29 |
| Table 2-17: Acceptable Presentation Contexts for Verification Activity | 29 |
| Table 2-18: Storage C-STORE Response Status | 30 |
| Table 3-1: Implementation Class/Version Name - Media Interchange | 35 |
| Table 3-2: Media - Application Profiles and Real-World Activities | 35 |
| Table 3-3: SOP Classes and Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K | 37 |
| Table 3-4: SOP Classes and Transfer Syntaxes for STD-GEN-CD and STD-GEN-DVD Profile | 43 |
| Table 3-5: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles | 46 |
| Table 3-6: Application Level Confidentiality Profile Attributes (standard tags) | 46 |
| Table 4-1: Single-Byte Character Sets without Code Extension | 51 |
| Table 4-2: Single-Byte Characters Sets with Code Extension | 52 |
| Table 4-3: Multi-Byte Character Sets without Code Extension | 52 |
| Table 4-4: Multi-Byte Character Sets with Code Extension | 52 |
| Table 6-1 SOP CLASSES for Storage | 55 |
| Table 6-2: Supported Non-Storage SOP Classes | 57 |
| Table 6-3: List of created SOP Classes | 58 |
| Table 6-4: Private Data Element Dictionary | 58 |

List of Figures

| | |
|--|----|
| Figure 1.8-1: Overview about DICOM capabilities of <i>syngo.via</i> View&GO VA15A..... | 13 |
| Figure 2.1-1: Application Data Flow Diagram..... | 14 |
| Figure 2.1-2: Sequence diagram – Verify..... | 16 |
| Figure 2.1-3: Sequence diagram – Storage | 17 |
| Figure 2.1-4: Sequence diagram – Query/Retrieve | 18 |
| Figure 3.1-1: Media Interchange Application Data Flow Diagram | 33 |
| Figure 3.1-2: Sequence diagram – Media creation on USB..... | 34 |

1 INTRODUCTION

1.1 REVISION HISTORY

| Version/ Status | Date of Issue | Product / Version | Author | Change & Reason of Change |
|--------------------|---------------|-------------------------------|-------------------------|--|
| 1.0 | 2018-12-01 | syngo.via View&GO VA15A | CT DD DS EU HU O 4 3 | Release version. |
| 0.1 | 2018-11-15 | syngo.via View&GO VA15A | CT DD DS EU HU O 4 3 | Draft version for syngo.via View&GO VA15A based on 2.0 version of VA10A. Anonymization profile update. C-Store SCP behavior update. |

1.2 GENERAL

The Conformance Statement describes the DICOM interface for the Siemens **syngo.via View&GO** in terms of DICOM Part PS3.2 2016a [1].

1.3 AUDIENCE

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

1.4 REMARKS

DICOM, by itself, does not guarantee interoperability. However, the Conformance Statement facilitates a first-level validation for interoperability between different applications supporting the same DICOM functionality as SCU and SCP, respectively.

This Conformance Statement is not intended to replace validation with other DICOM equipment to ensure proper exchange of information intended.

The scope of this Conformance Statement is to facilitate communication with Siemens and other vendors' medical equipment. The Conformance Statement should be read and understood in conjunction with the DICOM PS3.1-3.20 2016a Standard [1]. However, by itself it is not guaranteed to ensure the desired interoperability and a successful interconnectivity.

The user should be aware of the following important issues:

- The comparison of different conformance statements is the first step towards assessing interconnectivity.
- Test procedures should be defined and tests should be performed by the user to validate the connectivity desired. DICOM itself and the conformance parts do not specify this.
- The standard will evolve to meet the users' future requirements. Siemens is actively involved in developing the standard further and therefore reserves the right to make changes to its products or to discontinue its delivery.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens representative for the most recent product information.

1.5 TERMS AND DEFINITIONS

Terms used in this document shall be interpreted as defined in the DICOM Standard.

1.6 ABBREVIATIONS

| | |
|-------|--|
| ACR | American College of Radiology |
| AE | DICOM Application Entity |
| ASCII | American Standard Code for Information Interchange |
| DB | Database |
| DCS | DICOM Conformance Statement |
| DSA | Digital Subtraction Angiography |
| IIDC | Image-Intensifier Distortion Correction |
| IOD | DICOM Information Object Definition |
| ISO | International Standard Organization |
| FSC | File-set Creator |
| FSR | File-set Reader |
| FSU | File-set Updater |
| MWL | Modality Worklist |
| NEMA | National Electrical Manufacturers Association |
| O | Optional Key Attribute |
| PDU | DICOM Protocol Data Unit |
| R | Required Key Attribute |
| RIS | Radiology Information System |
| SC | Storage Commitment |
| SCU | DICOM Service Class User |
| SCP | DICOM Service Class Provider |
| SOP | DICOM Service-Object Pair |
| SCS | Specific Character Set |
| U | Unique Key Attribute |

1.7 REFERENCES

- [1] Digital Imaging and Communications in Medicine (DICOM PS3.1-3.20 2016a), National Electrical Manufacturers Association (NEMA), <http://medical.nema.org/>
- [2] IHE Radiology Technical Framework, Vol. I – IV, http://www.ihe.net/Technical_Frameworks

1.8 SCOPE AND FIELD OF APPLICATION




Figure 1.8-1: Overview about DICOM capabilities of **syngo.via View&GO VA15A**

2 NETWORKING

2.1 IMPLEMENTATION MODEL

2.1.1 Application Data Flow

The Application Data Flow diagram in Figure 2.1-1 depicts the DICOM data flow to and from the individual applications within **syngo.via View&GO**.




Figure 2.1-1: Application Data Flow Diagram

2.1.2 Functional Definition of AE's

The SCP components of the Application Entities of the **syngo.via View&GO** operate as background server processes. They exist as soon as the application starts and wait for association requests. Upon accepting an association with a negotiated Presentation Context they start to receive and process the requests described in the following sections.

The SCU components of the Application Entity are invoked upon requests from the user interface or indirectly triggered from internal processes.

2.1.2.1 Verification

Verification requests (SCP) will be processed and responded by the **syngo.via View&GO** AE. The **syngo.via View&GO** AE can also initiate an association and request verification to a remote AE (SCU).

2.1.2.2 Storage

The **syngo.via View&GO** Storage SCU is invoked directly by the user, or by the Query/Retrieve Application Entity that is responsible for processing retrieve requests. The request consists of data describing the composite objects selected for storage and the destination AET. An association is negotiated with the destination AE and the image data is transferred using the DIMSE C-STORE-Service.

The **syngo.via View&GO** Storage SCP starts to receive the Composite Image Objects and import them into the database after accepting an association with a negotiated Presentation Context. The system responses to the Storage Request immediately after reception of the Data.

2.1.2.3 Storage Commitment

The **syngo.via View&GO** does not support Storage Commitment.

2.1.2.4 Query

The C-FIND request to the remote SCP is invoked directly by the user. The remote SCP returns a list of responses with matching the request which are displayed to the user. The user can decide to start retrieving any of the responses or to issue another query.

As an SCU **syngo.via View&GO** supports

- Study Root Query/Retrieve Information Model (FIND)
- Patient Root Query/Retrieve Information Model (FIND)
- Patient/Study Only Query/Retrieve Information Model (FIND)
- Furthermore the SCU services may issue relational queries, if supported by the SCP node and required by the querying Application.

2.1.2.5 Retrieve

The **syngo.via View&GO** initiates a C-MOVE request to the remote Retrieve SCP. The remote Retrieve SCP in turn starts C-STORE sub operations to the **syngo.via View&GO** Storage SCP.

As an SCP **syngo.via View&GO** supports

- Study Root Query/Retrieve Information Model (MOVE,GET)
- Patient Root Query/Retrieve Information Model (MOVE,GET)
- Patient/Study Only Query/Retrieve Information Model (MOVE,GET)

2.1.3 Sequencing of Real-World Activities

Verify:

The communication between **syngo.via View&GO** and an external DICOM node in case of Verify is depicted in Figure 2.1-2 in more detail.




Figure 2.1-2: Sequence diagram – Verify

Storage:

The communication between **syngo.via View&GO** and an external DICOM node in case of triggering the transfer or accepting storage requests is depicted in more detail.




Figure 2.1-3: Sequence diagram – Storage

Query and Retrieval:

The communication between **syngo.via View&GO** and an external DICOM node in case of querying of objects from a remote DICOM node and retrieval to **syngo.via View&GO** is depicted in Figure 2.1-4 in more detail.




Figure 2.1-4: Sequence diagram – Query/Retrieve

2.2 AE SPECIFICATIONS

This section outlines the specifications for each of the Application Entities that are part of the **syngo.via View&GO** solution.

2.2.1 syngo.via View&GO AE

2.2.1.1 SOP Classes supported

This Application Entity provides Standard Conformance to the SOP Classes listed in Table 6-1 SOP CLASSES and

Table 6-2: Supported Non-Storage SOP Classes.

2.2.1.2 Association Establishment Policies

Table 2-1: Association Policies

| | |
|---|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 32 kB ¹ |
| Maximum number of simultaneous associations as an association acceptor | 12 ¹ |
| Maximum number of simultaneous associations as an association initiator | unlimited |

The **syngo.via View&GO** AE contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system. Nevertheless, transfer jobs to one distinct remote system (Send, Retrieve) will be run sequentially one after the other.

2.2.1.2.1 Asynchronous Nature

The **syngo.via View&GO** supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

Table 2-2: Asynchronous Nature as an Association Initiator

| | |
|---|----------|
| Maximum number of outstanding asynchronous transactions | Infinite |
|---|----------|

2.2.1.2.2 Implementation Identifying Information

Table 2-3: DICOM Implementation Class and Version

| | |
|-----------------------------|----------------------------------|
| Implementation Class UID | 1.3.12.2.1107.5.8.15.10.20090701 |
| Implementation Version Name | syngo.ViewAndGO |

2.2.1.3 Association Initiation Policy

syngo.via View&GO initiates associations while processing the service operations and internal messages as shown below:

Table 2-4: Association initiation policies

| Operation or Real-World Activity | Association for |
|----------------------------------|-----------------|
| Verification | C-ECHO |
| Send / Receive Instance | C-STORE |
| Querying a remote node | C-FIND |
| Retrieval of Instances | C-MOVE |

¹ Default, the value is configurable

2.2.1.3.1

2.2.1.3.2 Activity "Verification"

2.2.1.3.2.1 Description and Sequencing of Activities

The **syngo.via View&GO** verification process will accept an association, send ECHO request and receive reply.

2.2.1.3.2.2 Accepted Presentation Contexts

Table 2-5: Acceptable Presentation Contexts for Verification Activity

| Presentation Context Table | | | | | |
|----------------------------|-------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.2.1.3.2.3 SOP specific Conformance for SOP classes

syngo.via View&GO provides standard conformance to the Verification service SCU.

2.2.1.3.3 Activity "Send To"

2.2.1.3.3.1 Description and Sequencing of Activities

Storage of DICOM object is triggered internally in the **syngo.via View&GO**.

If an association to a remote Application Entity could successfully be established, each image will be transferred one after another via the same open association.

In case the transfer fails for a permanent reason (rejection permanent reported by SCP, all Presentation Contexts are refused, ...) the transfer will not be retried.

2.2.1.3.3.2 Proposed Presentation Contexts

For all supported Transfer objects (see SOP Classes in Table 6-1) the following Transfer Syntaxes are supported:

Table 2-6: Proposed Presentation Contexts for Storage

| UID value | Transfer Syntax | Type |
|------------------------|--|-------------------|
| 1.2.840.10008.1.2.1 | Explicit Value Representation Little Endian native | Image / non-image |
| 1.2.840.10008.1.2 | Implicit Value Representation Little Endian native | Image / non-image |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | Image |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed | Image |
| 1.2.840.10008.1.2.5 | RLE Lossless compressed | Image |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed | Image |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed | Image |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed | Image |

Storage SCU Service will choose a compressed or uncompressed Transfer Syntax among those accepted by the SCP for images.

An instance will be JPEG lossless compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME or RGB or YBR_FULL or YBR_FULL_422
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) equal to 12'D or 8'D
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D

An instance will be JPEG lossy compressed only if it fulfills the following criteria:

- is an image and not already compressed
- photometric interpretation (0028,0004) is MONOCHROME or RGB
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) equal to 12'D or 8'D
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D

An instance will be JPEG 2000 lossless compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric interpretation (0028,0004) not MONOCHROME or RGB or YBR_FULL or YBR_FULL_422
- Bits Allocated (0028,0100) either 16'D nor 8'D

An instance will be JPEG 2000 lossy compressed only if it fulfills the following criteria:

- is an image and not already compressed
- Photometric interpretation (0028,0004) is MONOCHROME or RGB
- Bits Stored (0028,0101) equal to 12'D or 8'D

There is no extended negotiation as an SCU.

2.2.1.3.3.3 SOP specific Conformance for SOP classes

The **syngo.via View&GO** will not add or change private attributes by default, even in case objects are compressed or image header is updated according to IHE [2] Patient Information Reconciliation. The behavior of **syngo.via View&GO** when encountering status codes in a C-STORE response is summarized in Table 2-7:

Table 2-7: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|------------|--|
| Error | Duplicate SOP Instance UID: some of the instances sent to the SCP were already available there. | 0111 | Job is continued till the end. |
| Error | Out-Of-Resources: The remote node has run out of resources (storage resources for example) | A7XX | Job is continued till the end. The result can be checked by checking the availability of the data on the target node. |
| Error | Any other DIMSE Error Status | XXXX | User does not get notified. Job is continued till the end. Error is logged in the system log. The result can be checked by checking the availability of the data on the target node. |
| Success | Image is successfully stored on file system. | 0000 | User does not get notified. The result can be checked by checking the availability of the data on the target node. |

Table 2-8: DICOM Command Communication Failure Behavior

| Exception | Behavior |
|---------------------|--|
| Timeout | User does not get notified. The result can be checked by checking the availability of the data on the target node. |
| Association Aborted | User does not get notified. The result can be checked by checking the availability of the data on the target node. |

2.2.1.3.3.4 CT Derived object for LungCAD findings generated by **syngo.via View&GO**

| Attribute | Tag | Type | Description |
|------------------------|--------------|------|--|
| Specific Character Set | (0008, 0005) | 1C | Copied from Input Image |
| Image Type | (0008,0008) | 1 | Value 1: DERIVED Value 2: SECONDARY Value 3: AXIAL Value 4: AlgorithmName_AlgorithmVersion_DO |
| Instance Creation Date | (0008,0012) | 3 | DO ¹ instance creation date |
| Instance Creation Time | (0008,0013) | 3 | DO ¹ instance creation time |

| Attribute | Tag | Type | Description |
|------------------------------------|-------------|------|--|
| SOP Class UID | (0008,0016) | 1 | Copied from Input Image |
| SOP Instance UID | (0008,0018) | 1 | 1.3.12.2.1107.5.99.3.99.UID ² |
| Study Date | (0008,0020) | 2 | Copied from Input Image |
| Series Date | (0008,0021) | 3 | Creation Date of the DO ¹ in <YYYYMMDD> format |
| Acquisition Date | (0008,0022) | 3 | Copied from Input Image |
| Study Time | (0008,0030) | 2 | Copied from Input Image |
| Series Time | (0008,0031) | 3 | Creation Time of the DO ¹ in <HHMMSS> format |
| Acquisition Time | (0008,0032) | 3 | Copied from Input Image |
| Accession Number | (0008,0050) | 2 | Copied from Input Image |
| Modality | (0008,0060) | 1 | Copied from Input Image |
| Manufacturer | (0008,0070) | 2 | Copied from Input Image |
| Institution Name | (0008,0080) | 3 | Copied from Input Image |
| Referring Physician's Name | (0008,0090) | 2 | Copied from Input Image |
| Study Description | (0008,1030) | 3 | Copied from Input Image |
| Series Description | (0008,103E) | 3 | Series description as configured by the application |
| Name Of Physician(s) Reading Study | (0008,1060) | 3 | Copied from Input Image |
| Manufacturer's Model Name | (0008,1090) | 3 | Copied from Input Image |
| Derivation description | (0008,2111) | 3 | Set AlgorithmName_AlgorithmVersion_DO value. For example: LUNGCAD_VD10C_DO) |
| Source Image Sequence | (0008,2112) | 3 | Set input image SOP class UID and input image SOP instance UID from original image |
| >Reference SOP Class UID | (0008,1150) | | From original image's SOP Class UID |
| >Reference SOP Instance UID | (0008,1155) | | From original image's SOP Instance UID |
| Patient Name | (0010,0010) | 2 | Copied from Input Image |
| Patient ID | (0010,0020) | 2 | Copied from Input Image |
| Patient's Birth Date | (0010,0030) | 2 | Copied from Input Image |
| Patient's Sex | (0010,0040) | 2 | Copied from Input Image |
| Patient's Size | (0010,1020) | 3 | Copied from Input Image |
| Patient's Weight | (0010,1030) | 3 | Copied from Input Image |
| Pregnancy Status | (0010,21C0) | 3 | Copied from Input Image |
| Imager Pixel Spacing | (0018,1164) | 3 | Copied from Input Image |
| Patient Position | (0018,5100) | 2C | Copied from Input Image |
| View Position | (0018,5101) | 3 | Copied from Input Image |
| Detector Element Physical Size | (0018,7020) | 3 | Copied from Input Image |
| Detector Element Spacing | (0018,7022) | 3 | Copied from Input Image |
| Study Instance UID | (0020,000D) | 1 | Copied from Input Image |
| Series Instance UID | (0020,000E) | 1 | 1.3.12.2.1107.5.99.3.99.UID ² |
| Study ID | (0020,0010) | 2 | Copied from Input Image |
| Series Number | (0020,0011) | 2 | Series Number as configured by the user |
| Patient Orientation | (0020,0020) | 2 | Copied from Input Image |
| Image Laterality | (0020,0062) | 1 | Copied from Input Image |

| Attribute | Tag | Type | Description |
|-----------------------------|-------------|------|--|
| Samples PerPixel | (0028,0002) | 1 | Copied from input image |
| Photometric Interpretation | (0028,0004) | 1 | Copied from input image |
| Rows | (0028,0010) | 1 | Copied from input image |
| Columns | (0028,0011) | 1 | Copied from input image |
| Pixel Spacing | (0028,0030) | 1 | Copied from input image |
| Bits Allocated | (0028,0100) | 1 | Copied from input image |
| Bits Stored | (0028,0101) | 1 | Copied from input image |
| High Bit | (0028,0102) | 1 | Copied from input image |
| Pixel Representation | (0028,0103) | 1 | Copied from input image |
| Window Center | (0028,1050) | 1 | Copied from Input Image |
| Window Width | (0028,1051) | 1 | Copied from Input Image |
| Rescale Intercept | (0028,1052) | 1 | Copied from input image |
| Rescale Slope | (0028,1053) | 1 | Copied from input image |
| Rescale Type | (0028,1054) | 1 | Copied from input image |
| Presentation Creator's Name | (0070,0084) | 2 | Algorithm Name_Version_DO ¹ |
| Study Comments | (0032,4000) | 3 | Copied from input image |
| Pixel Data | (7FE0,0010) | 1 | Sets the Derived Pixel data |

1 – Derived Object

2 – UID generated by syngo.via
View@GO

2.2.1.3.4 Activity “Querying a Remote Node” for Instances

2.2.1.3.4.1 Description and Sequencing of Activities

The associated Real-World activity is a C-Find request initiated by the user (see also Figure 2.1-4). The user specifies some attributes and will send a C-Find request (according to the query model) and will then return the results to the initiating application.

2.2.1.3.4.2 Proposed Presentation Contexts

The **syngo.via View&GO** will propose Presentation Contexts as shown in the following table:

Table 2-9: Proposed Presentation Contexts for Query

| Presentation Context Table | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | Yes |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Study Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | Yes |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |

| | | | | | |
|---|-----------------------------|---------------------------|---------------------|-----|----|
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Patient/Study Only Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

Table 2-10: Extended Negotiation as an SCU

| Name | UID | Extended Negotiation |
|--|-----------------------------|--|
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Relational Query will be negotiated as defined in DICOM PS3.4 2016a. |
| Study Root Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Relational Query will be negotiated as defined in DICOM PS3.4 2016a. |

2.2.1.3.4.3 SOP Specific Conformance Statement to Query SOP classes

The **syngo.via View&GO** checks for the following status codes in the Query SCP's C-FIND-Response:

Table 2-11: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|--------------------|------------------------------------|
| Failure | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process | Any none null Code | Failure reported to user |
| Pending | All optional keys are supported the same manner as Required Keys. | FF00 | Pending state is indicated to user |
| | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01 | Pending state is indicated to user |
| Success | Query has been performed successfully. | 0000 | Success reported to user |

Table 2-12: DICOM Command Communication Failure Behavior

| exception | Behavior |
|---------------------|--|
| Timeout | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user |

The **syngo.via View&GO** supports the following query levels:

- Study

- Series

Matching Keys on Instance Level is not supported by **syngo.via View&GO** as SCU.

The following table lists the various attributes at Study and Series levels, which can be used for **hierarchical** queries as well as return values for display. The display capabilities are highly configurable and “yes” indicates that it is possible to configure display of the data:
A “yes” in the **UI** column will indicate that the attribute may be visualized when browsing the Query results with the Browser. The Browser display is additionally influenced by the related Browser configuration

Table 2-13: Attributes supported for instance Query - SCU

| Attribute Name | Tag | Type | User input | UI |
|---|---------------------------|------|-------------|-----|
| Study Level | | | | |
| Patient's Name | (0010,0010) | O | enter value | yes |
| Patient ID | (0010,0020) | O | enter value | yes |
| Issuer of Patient ID | (0010,0021) | O | enter value | yes |
| Patient's Birth Date | (0010,0030) | O | enter value | yes |
| Patient's Birth Time | (0010,0032) | O | enter value | yes |
| Patient's Sex | (0010,0040) | O | enter value | yes |
| Accession Number | (0008,0050) | O | enter value | yes |
| Study ID | (0020,0010) | O | enter value | yes |
| Study Instance UID | (0020,000D) | U | enter value | yes |
| Study Date | (0008,0020) | O | enter value | yes |
| Study Time | (0008,0030) | O | enter value | yes |
| Referring Physician' s Name | (0008,0090) | O | enter value | yes |
| Study Description | (0008,1030) | O | enter value | yes |
| Number of Study related Instances | (0020,1208) | O | - | yes |
| Modalities in Study | (0008,0061) | O | enter value | yes |
| Number of Study Related Series | (0020,1206) | O | - | yes |
| Series Level | | | | |
| Modality | (0008,0060) | O | enter value | yes |
| Series Date | (0008,0021) | O | enter value | yes |
| Series Time | (0008,0031) | O | enter value | yes |
| Number of Series related Instances | (0020,1209) | O | - | yes |
| Series Number | (0020,0011) | O | enter value | yes |
| Series Description | (0008,103E) | O | enter value | yes |
| Request Attributes Sequence \ Requested Procedure ID | (0040,0275) \ (0040,1001) | O | enter value | yes |
| Request Attributes Sequence \ Scheduled Procedure Step ID | (0040,0275) \ (0040,0009) | O | enter value | yes |
| Performed Procedure Step Start Date | (0040,0244) | O | enter value | yes |
| Performed Procedure Step Start Time | (0040,0245) | O | enter value | yes |
| Series Instance UID | (0020,000E) | U | - | yes |

2.2.1.3.5 Activity Retrieving Instances from a remote node

2.2.1.3.5.1 Description and Sequencing of Activities

The C-MOVE-RQs are used to retrieve the referenced images. The Retrieve AE supports the query model Study Root.

2.2.1.3.5.2 Accepted Presentation Contexts

Table 2-14: Proposed Presentation Contexts for Retrieve and Activity “MOVE SCU”

| Presentation Context Table | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Patient Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Study Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Patient Study Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.2.1.3.5.3 SOP Specific Conformance Statement for Move SCU Classes

At association establishment time the C-MOVE presentation context shall be negotiated. When the C-MOVE-RQ is processed, the Move Destination attribute (receiver of images) is ignored. However the Move Destination AE must conform to the DICOM conventions (value representation AE).

The behavior of **syngo.via View&GO** when encountering status codes in a C-MOVE response is summarized in Table 2-15:

Table 2-15: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|-------------------|--------------------------|
| Error | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process; Move destination unknown | Any non null Code | Failure reported to user |
| Pending | Move Operation continues | FF00 | Operation continues in |

| | | | |
|---------|---------------------------------------|------|--------------------------|
| | | | background |
| Success | Move has been performed successfully. | 0000 | Success reported to user |

Table 2-16: DICOM Command Communication Failure Behavior

| exception | Behavior |
|---------------------|---|
| Timeout | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user |

2.2.1.4 Association Acceptance Policy

The **syngo.via View&GO** attempts to accept a new association for

- DIMSE C-STORE
service operations.

Generally associations are accepted if all of the following conditions are true:

- The "called AET" matches one of the configured Application Entity Titles of the **syngo.via View&GO**.
- The "calling AET" is known (configured) at **syngo.via View&GO**.
- The maximum number of incoming associations is not reached.
- At least one Presentation Context with a minimum of one suitable transfer syntax has been proposed as defined by the "Presentation Context Tables" in the following subsections.
- The system has enough available resources to perform the service requested (e.g. enough free disk space, less than the max. number of associations are already in use)

2.2.1.4.1 Activity "Verification"

2.2.1.4.1.1 Description and Sequencing of Activities

The **syngo.via View&GO** verification process will accept an association, receive ECHO request and send reply.

2.2.1.4.1.2 Accepted Presentation Contexts

Table 2-17: Acceptable Presentation Contexts for Verification Activity

| Presentation Context Table | | | | | |
|----------------------------|-------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

2.2.1.4.1.3 SOP specific Conformance for SOP classes

syngo.via View&GO provides standard conformance to the Verification service SCP.

2.2.1.4.2 Activity “Receive Instances”

2.2.1.4.2.1 Description and Sequencing of Activities

The **syngo.via View&GO** receiving process will accept C-STORE association requests, receive any objects transmitted on that association and store the objects on disk.

2.2.1.4.2.2 Accepted Presentation Contexts

For all supported Transfer objects (see SOP Classes in Table 6-1) the Transfer Syntaxes described in Table 2-6 are supported.

Generally all Presentation Contexts are accepted as long as they contain at least one suitable Transfer Syntax. All other Presentation Contexts are rejected.

There is no Extended Negotiation as an SCP

2.2.1.4.2.3 SOP-specific Conformance Statement for Storage SOP classes

The **syngo.via View&GO** conforms to the Full Storage Class at Level 2.

In case of a successful C-STORE operation, the image has successfully been received in either Explicit Little Endian format or in the compression format received.

syngo.via View&GO AE returns the status “success” when the data is received and a minimal image header validation has been performed.

The following header attributes must be available and filled:

- SOP Class UID,
- Study Instance UID,
- Series Instance UID and
- SOP Instance UID.

Table 2-18: Storage C-STORE Response Status

| Service Status | Further Meaning | Error Code | Reason |
|----------------|-------------------|------------|---|
| Success | success | 0000 | Image received correctly (success notification is done after receiving, before indexing and storing) |
| Failure | Out-of-resource | A700 | No resource left in the Short Term Storage |
| Failure | Unable to Process | Cxxx | Error during instance reception |

| | | | |
|---------|-----------------------------------|------|---|
| Failure | Data set does not match SOP Class | A9xx | The data set is not conform to the SOP Class contained in the resource. |
|---------|-----------------------------------|------|---|

Restriction: successful operation does not guarantee storage on disk and storage of header data in the database.

2.2.1.4.2.4 Other SOP specific behavior

- If an image is received that is already stored in the database - identified by the SOP Instance UID - the new image will be ignored. The existing instance is not superseded.

2.3 NETWORK INTERFACES

2.3.1 Physical Network Interface

The **syngo.via View&GO** provides DICOM 3.0 TCP/IP network communication support as defined in Part 8 of DICOM [1]. The network communication is independent from the physical medium over which TCP/IP executes; it inherits this from the Windows OS system upon which it executes.

2.3.2 Additional Protocols

none

2.3.3 IPv4 and IPv6 Support

IPv4 and IPv6 are supported. Regarding IPv6 please note, that the complete networking infrastructure in the hospital (firewalls, DNS-Servers, ...) must support IPv6 in order to get a functioning communication.

2.4 CONFIGURATION

2.4.1 AE Title/Presentation Address Mapping

AE Titles shall be unique within the hospital. A common way to achieve that is to use the hostname as part of the AE Titles. The string can be up to 16 characters and must not contain any extended characters. Only 7-bit ASCII characters (excluding Control Characters) are allowed according to DICOM [1].

2.4.1.1 Local AE Titles

The **syngo.via View&GO** allows to configure AETitles, Ports any way. Default delivery is that all services are using the same AE title and only one port number.

| Parameter | Configurable | Default Value |
|------------------|--------------|--|
| Default AE title | Yes | hostname in uppercase characters; limited to 16 characters |
| Default Port | Yes | 104 |

2.4.1.2 Remote AE Title/Presentation Address Mapping

2.4.1.2.1 Remote Association Initiators

All relevant remote applications that may setup DICOM associations towards **syngo.via View&GO** need to be configured in **syngo.via View&GO**, before the association can be established. This behavior is configurable but it is recommended, not to change this behavior.

The mapping of external AE Titles to TCP/IP addresses and ports is configurable. The Application Entity Titles and supported transfer syntaxes need to be known for configuration.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

2.4.1.2.2 Remote SCP's

For remote applications that shall be able to accept DICOM associations from **syngo.via View&GO**, the following information needs to be available:

- Application Entity Title
- Host Name / IP address on which the remote application service runs
- Port number on which the remote application accepts association requests.

The remote system will be indicated in the UI of **syngo.via View&GO** with a logical name that is also entered when configuring the node in the configuration UI.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

3 MEDIA INTERCHANGE

3.1 IMPLEMENTATION MODELS

3.1.1 Application Data Flow Diagram




Figure 3.1-1: Media Interchange Application Data Flow Diagram

The **syngo.via View&GO** provides the functionality to Import or Export DICOM Instances from and to the File System. During export, a DICOMDIR may also be generated (user selection). All SOP Classes defined in Table 3-3 and Table 3-4 are supported for the Import/Export functionality.

3.1.2 Functional definitions of AEs

The **syngo.via View&GO** application is capable of

- creating a new File-set in the File System (Export to ...)
- importing SOP Instances from the medium onto local storage
- writing the File-sets DICOMDIR information into the file system

3.1.3 Sequencing of Real-World Activities

Whenever data is written to an external media, **syngo.via View&GO** a DICOMDIR may be created for the selected data.




Figure 3.1-2: Sequence diagram – Media creation on USB




Figure 5.1-3: Sequence diagram – Media import

3.1.4 File Meta Information for Implementation Class and Version

This section describes the values assigned to the File Meta Information attributes (see DICOM [1] part PS 3.10) that pertain to the Implementation Class and Version. The implementation Class UID and the Implementation Version name in the File Meta Header are the same as the values specified for networking.

Table 3-1: Implementation Class/Version Name - Media Interchange

| | |
|-------------------------------|----------------------------------|
| File Meta Information Version | 0001 |
| Implementation Class UID | 1.3.12.2.1107.5.8.15.10.20090701 |
| Implementation Version Name | syngo.ViewAndGO |

3.2 AE SPECIFICATIONS

3.2.1 Media Storage AE – Specification

The **syngo.via View&GO** provides conformance to the following Application Profiles as an FSR. FSC and FSU are supported only on a non-optical storage device (e.g. USB stick).

In addition augmented conformance is provided to store extra data attributes important for the full feature support of the **syngo®**-based products. Details are listed below:

Table 3-2: Media - Application Profiles and Real-World Activities

| Application Profiles Supported | Real-World Activity | Role | Service Class Option |
|--------------------------------|---------------------|------|----------------------|
|--------------------------------|---------------------|------|----------------------|

| Application Profiles Supported | Real-World Activity | Role | Service Class Option |
|--------------------------------|--|--------------------|----------------------|
| AUG-GEN-CD | Browse Directory Information Import into Application | FSR | Interchange |
| AUG-GEN-DVD | | | |
| AUG- GEN-DVD-J2K | | | |
| AUG- GEN-USB-J2K | Browse Directory Information Import into Application Export to local Archive Media | FSR , FSC , FSU | |
| STD-GEN-CD | Browse Directory Information Import into Application | FSR | Interchange |
| STD-GEN-DVD | | | |
| STD-GEN-DVD-J2K | | | |
| STD-GEN-USB-J2K | Browse Directory Information Import into Application Export to local Archive Media | FSR , FSC , FSU | |

3.2.1.1 Real-World Activities

3.2.1.1.1 Activity “Browse Directory Information”

The **syngo.via View&GO** acts as FSR using the interchange option when requested to read the media directory.

The **syngo.via View&GO** will read the DICOMDIR and insert those directory entries that are valid for the application profiles supported, into a local database. The database then is used for browsing media contents.

Note: The “Icon Image Sequence” is also supported in DICOMDIR. But only those Icon Images with “Bits Allocated” (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the Browser.

3.2.1.1.1.1 Media Storage Application Profiles

See Table 3-2 for the Application Profiles listed that invoke this Application Entity for the Browse Directory Information

3.2.1.1.2 Activity “Import into Application”

The **syngo.via View&GO** application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the application.

The SOP Instance selected from the media directory will be copied into the running Application. Only SOP Instances, that are valid for the application profile supported and supported by **syngo.via View&GO** (see Table 6-1), can be retrieved from media.

3.2.1.1.3 Media Storage Application Profiles

See Table 3-2 for the Application Profiles listed that invoke this Application Entity for the local Archive Media Real-World Activity.

3.2.1.2 SOP Classes and Transfer Syntaxes

These Application Profiles are based on the Media Storage Service Class with the Interchange Option. In the table below (Table 3-3) the Transfer Syntax UID “RLE Lossless” only applies for decompression.

Table 3-3: SOP Classes and Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|--|--------------------------------|---|
| Basic Directory | 1.2.840.10008.1.3.10 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| 12-lead ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ambulatory ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Blending Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Breast Tomosynthesis Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Cardiac Electrophysiology Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.3.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Color Softcopy Presentation State Storage (store & forward only) | 1.2.840.10008.5.1.4.1.1.11.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Comprehensive SR | 1.2.840.10008.5.1.4.1.1.88.33 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| CR Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Deformable Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|---|-------------------------------|---|
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| Digital Mammography X-Ray Image Storage- For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | Explicit VR Big Endian Uncompressed 1.2.840.10008.1.2.2 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| Digital Mammography X-Ray Image Storage- For Presentation | 1.2.840.10008.5.1.4.1.1.1.2 | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| Encapsulated PDF | 1.2.840.10008.5.1.4.1.1.104.1 | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Enhanced CT Image Storage | 1.2.840.10008.5.1.4.1.1.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Enhanced MR Color Image Storage | 1.2.840.10008.5.1.4.1.1.4.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|---|--------------------------------|--|
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Enhanced Structured Report | 1.2.840.10008.5.1.4.1.1.88.22 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced XA Image Storage | 1.2.840.10008.5.1.4.1.1.12.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 RLE Lossless 1.2.840.10008.1.2.5 JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Enhanced XRF Image Storage | 1.2.840.10008.5.1.4.1.1.12.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 RLE Lossless 1.2.840.10008.1.2.5 JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| General ECG Waveform | 1.2.840.10008.5.1.4.1.1.9.1.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Grayscale Softcopy Presentation State | 1.2.840.10008.5.1.4.1.1.11.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Hemodynamic Waveform | 1.2.840.10008.5.1.4.1.1.9.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Key Object Selection Document | 1.2.840.10008.5.1.4.1.1.88.59 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Mammography CAD SR Storage | 1.2.840.10008.5.1.4.1.1.88.50 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 RLE Lossless 1.2.840.10008.1.2.5 JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 RLE Lossless 1.2.840.10008.1.2.5 JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Multi-frame Grayscale Byte SC Image Storage | 1.2.840.10008.5.1.4.1.1.7.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 RLE Lossless 1.2.840.10008.1.2.5 JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Multi-frame Grayscale Word SC Image Storage | 1.2.840.10008.5.1.4.1.1.7.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 JPEG Lossless Process 14 (selection value 1) |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|--|-------------------------------|--|
| | | 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) |
| | | 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless |
| | | 1.2.840.10008.1.2.5 |
| Multi-frame Single Bit SC Image Storage | 1.2.840.10008.5.1.4.1.1.7.1 | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) |
| | | 1.2.840.10008.1.2.4.70 |
| Multi-frame True Color SC Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | JPEG Lossy (baseline or extended) |
| | | 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless |
| | | 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) |
| | | 1.2.840.10008.1.2.4.70 |
| | | RLE Lossless |
| Positron Emission Tomography Image Storage | 1.2.840.10008.5.1.4.1.1.128 | 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) |
| Procedure Log Storage | 1.2.840.10008.5.1.4.1.1.88.40 | 1.2.840.10008.1.2.4.70 |
| | | RLE Lossless |
| | | 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| Pseudo-Color Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.3 | Explicit VR Little Endian Uncompressed |
| Raw DataStorage | 1.2.840.10008.5.1.4.1.1.66 | 1.2.840.10008.1.2.1 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) |
| | | 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless |
| | | 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed |
| | | 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) |
| | | 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) |
| | | 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless |
| | | 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|---------------------------------------|-------------------------------|---|
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Dose Storage | 1.2.840.10008.5.1.4.1.1.481.2 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Image Storage | 1.2.840.10008.5.1.4.1.1.481.1 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Ion Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.9 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Ion Plan Storage | 1.2.840.10008.5.1.4.1.1.481.8 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Plan Storage | 1.2.840.10008.5.1.4.1.1.481.5 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 |
| | | 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|-------------------------------------|-------------------------------|---|
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Treatment Summary Record Storage | 1.2.840.10008.5.1.4.1.1.481.7 | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| Spatial Fiducials Storage | 1.2.840.10008.5.1.4.1.1.66.2 | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | |
| | | |
| | | |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|--|--------------------------------|---|
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| Ultrasound Multi-frame Storage (retired) | 1.2.840.10008.5.1.4.1.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | |
| Ultrasound Multi-frame Image storage | 1.2.840.10008.5.1.4.1.1.3.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | |
| X-Ray 3D Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.1 | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| X-Ray Radiation Dose Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.67 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | JPEG Lossless Process 14 (selection value 1) 1.2.840.10008.1.2.4.70 |
| | | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| | | JPEG Lossy (baseline or extended) 1.2.840.10008.1.2.4.50 1.2.840.10008.1.2.4.51 |
| | | RLE Lossless 1.2.840.10008.1.2.5 |
| | | JPEG 2000 Lossless 1.2.840.10008.1.2.4.90 |
| | | JPEG 2000 Lossy 1.2.840.10008.1.2.4.91 |
| | | |

Table 3-4: SOP Classes and Transfer Syntaxes for STD-GEN-CD and STD-GEN-DVD Profile

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|-------------------------------|----------------------|---|
| Basic Directory | 1.2.840.10008.1.3.10 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|--|--------------------------------|--|
| 12-lead ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ambulatory ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Blending Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Breast Tomosynthesis Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Cardiac Electrophysiology Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.3.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Deformable Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Digital Mammography X-Ray Image Storage- For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Digital Mammography X-Ray Image Storage- For Presentation | 1.2.840.10008.5.1.4.1.1.1.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Encapsulated PDF Storage | 1.2.840.10008.5.1.4.1.1.104.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced CT Image Storage | 1.2.840.10008.5.1.4.1.1.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced MR Color Image Storage | 1.2.840.10008.5.1.4.1.1.4.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced XA Image Storage | 1.2.840.10008.5.1.4.1.1.12.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Enhanced XRF Image Storage | 1.2.840.10008.5.1.4.1.1.12.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| General ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Grayscale Softcopy Presentation State StorageSOP Class | 1.2.840.10008.5.1.4.1.1.11.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Hemodynamic Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.2.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Key Object Selection Document Storage | 1.2.840.10008.5.1.4.1.1.88.59 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Mammography CAD SR Storage | 1.2.840.10008.5.1.4.1.1.88.50 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Multi-frame Single Bit SC Image | 1.2.840.10008.5.1.4.1.1.7.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Multi-frame True Color SC Image | 1.2.840.10008.5.1.4.1.1.7.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |

| Information Object Definition | SOP Class UID | Transfer Syntax UID |
|--|--------------------------------|---|
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Pseudo-Color Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Positron Emission Tomography Image Storage | 1.2.840.10008.5.1.4.1.1.128 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Procedure Log Storage | 1.2.840.10008.5.1.4.1.1.88.40 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Raw DataStorage | 1.2.840.10008.5.1.4.1.1.66 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Dose Storage | 1.2.840.10008.5.1.4.1.1.481.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Image Storage | 1.2.840.10008.5.1.4.1.1.481.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Ion Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.9 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Ion Plan Storage | 1.2.840.10008.5.1.4.1.1.481.8 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Plan Storage | 1.2.840.10008.5.1.4.1.1.481.5 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| RT Treatment Summary Record Storage | 1.2.840.10008.5.1.4.1.1.481.7 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Spatial Fiducials Storage | 1.2.840.10008.5.1.4.1.1.66.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ultrasound Multi-frame (retired) | 1.2.840.10008.5.1.4.1.1.3 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| Ultrasound Multi-frame Storage Image | 1.2.840.10008.5.1.4.1.1.3.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| X-Ray 3D Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| X-Ray Radiation Dose Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.67 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 |

3.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

3.3.1 Augmented Application Profiles

The standard application profiles are augmented with private object Siemens CSA Non-Image.

Table 3-5: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles

| Information Object Definition | SOP Class UID | Transfer Syntax UID | FSC | FSR |
|-------------------------------|---------------------|---|-----|-----|
| CSA Non-Image Storage | 1.3.12.2.1107.5.9.1 | Explicit VR Little Endian Uncompressed 1.2.840.10008.1.2.1 | O | M |

The Siemens non-image is typically used for raw data and 3D private data.

3.4 MEDIA CONFIGURATION

none

3.5 Attribute confidentiality profiles

3.5.1 De-identification

The **syngo.via View&GO** application can de-identify attributes. During export to file system it is the user responsibility to select anonymization.

For anonymization private attributes are not included in anonymized Studies.

The following table contains the attributes that are anonymized:

Table 3-6: Application Level Confidentiality Profile Attributes (standard tags)

| DICOM Tag | Attribute Name |
|-------------|--|
| (0000,1001) | Requested SOP Instance UID |
| (0002,0003) | Media Storage SOP Instance UID |
| (0004,1511) | Referenced SOP Instance UID in File |
| (0008,0014) | Instance Creator UID |
| (0008,0015) | Instance Coercion DateTime |
| (0008,0018) | SOP Instance UID |
| (0008,0020) | Study Date |
| (0008,0021) | Series Date |
| (0008,0022) | Acquisition Date |
| (0008,0023) | Content Date |
| (0008,0024) | Overlay Date |
| (0008,0025) | Curve Date |
| (0008,002A) | Acquisition DateTime |
| (0008,0030) | Study Time |
| (0008,0031) | Series Time |
| (0008,0032) | Acquisition Time |
| (0008,0033) | Content Time |
| (0008,0034) | Overlay Time |
| (0008,0035) | Curve Time |
| (0008,0050) | Accession Number |
| (0008,0058) | Failed SOP Instance UID List |
| (0008,0080) | Institution Name |
| (0008,0081) | Institution Address |
| (0008,0082) | Institution Code Sequence |
| (0008,0090) | Referring Physician's Name |
| (0008,0092) | Referring Physician's Address |
| (0008,0094) | Referring Physician's Telephone Numbers |
| (0008,0096) | Referring Physician's Identification Sequence |
| (0008,010D) | Context Group Extension Creator UID |
| (0008,0201) | Timezone Offset From UTC |
| (0008,1010) | Station Name |
| (0008,1030) | Study Description |
| (0008,103E) | Series Description |
| (0008,1040) | Institutional Department Name |
| (0008,1048) | Physician(s) of Record |
| (0008,1049) | Physician(s) of Record Identification Sequence |

| DICOM Tag | Attribute Name |
|-------------|---|
| (0008,1050) | Performing Physicians' Name |
| (0008,1052) | Performing Physicians' Identification Sequence |
| (0008,1060) | Name of Physician(s) Reading Study |
| (0008,1062) | Physician Reading Study Identification Sequence |
| (0008,1070) | Operators' Name |
| (0008,1072) | Operators' Identification Sequence |
| (0008,1080) | Admitting Diagnoses Description |
| (0008,1084) | Admitting Diagnoses Code Sequence |
| (0008,1110) | Referenced Study Sequence |
| (0008,1111) | Referenced Performed Procedure Step Sequence |
| (0008,1120) | Referenced Patient Sequence |
| (0008,1140) | Referenced Image Sequence |
| (0008,1155) | Referenced SOP Instance UID |
| (0008,1195) | Transaction UID |
| (0008,2111) | Derivation Description |
| (0008,2112) | Source Image Sequence |
| (0008,3010) | Irradiation Event UID |
| (0008,4000) | Identifying Comments |
| (0008,9123) | Creator Version UID |
| (0010,0010) | Patient's Name |
| (0010,0020) | Patient ID |
| (0010,0021) | Issuer of Patient ID |
| (0010,0030) | Patient's Birth Date |
| (0010,0032) | Patient's Birth Time |
| (0010,0040) | Patient's Sex |
| (0010,0050) | Patient's Insurance Plan Code Sequence |
| (0010,0101) | Patient's Primary Language Code Sequence |
| (0010,0102) | Patient's Primary Language Modifier Code Sequence |
| (0010,1000) | Other Patient IDs |
| (0010,1001) | Other Patient Names |
| (0010,1002) | Other Patient IDs Sequence |
| (0010,1005) | Patient's Birth Name |
| (0010,1010) | Patient's Age |
| (0010,1020) | Patient's Size |
| (0010,1030) | Patient's Weight |
| (0010,1040) | Patient Address |
| (0010,1050) | Insurance Plan Identification |
| (0010,1060) | Patient's Mother's Birth Name |
| (0010,1080) | Military Rank |
| (0010,1081) | Branch of Service |
| (0010,1090) | Medical Record Locator |
| (0010,1100) | Referenced Patient Photo Sequence |
| (0010,2000) | Medical Alerts |
| (0010,2110) | Allergies |
| (0010,2150) | Country of Residence |
| (0010,2152) | Region of Residence |
| (0010,2154) | Patient's Telephone Number |
| (0010,2160) | Ethnic Group |
| (0010,2180) | Occupation |
| (0010,21A0) | Smoking Status |
| (0010,21B0) | Additional Patient's History |
| (0010,21C0) | Pregnancy Status |
| (0010,21D0) | Last Menstrual Date |
| (0010,21F0) | Patient's Religious Preference |
| (0010,2203) | Patient Sex Neutered |
| (0010,2297) | Responsible Person |
| (0010,2299) | Responsible Organization |
| (0010,4000) | Patient Comments |
| (0018,0010) | Contrast Bolus Agent |
| (0018,1000) | Device Serial Number |
| (0018,1002) | Device UID |
| (0018,1004) | Plate ID |
| (0018,1005) | Generator ID |
| (0018,1007) | Cassette ID |
| (0018,1008) | Gantry ID |
| (0018,1030) | Protocol Name |
| (0018,1400) | Acquisition Device Processing Description |
| (0018,2042) | Target UID |
| (0018,4000) | Acquisition Comments |
| (0018,700A) | Detector ID |

| DICOM Tag | Attribute Name |
|-------------|--|
| (0018,9424) | Acquisition Protocol Description |
| (0018,9516) | Start Acquisition DateTime |
| (0018,9517) | End Acquisition DateTime |
| (0018,A003) | Contribution Description |
| (0020,000D) | Study Instance UID |
| (0020,000E) | Series Instance UID |
| (0020,0010) | Study ID |
| (0020,0052) | Frame of Reference UID |
| (0020,0200) | Synchronization Frame of Reference UID |
| (0020,3401) | Modifying Device ID |
| (0020,3404) | Modifying Device Manufacturer |
| (0020,3406) | Modified Image Description |
| (0020,4000) | Image Comments |
| (0020,9158) | Frame Comments |
| (0020,9161) | Concatenation UID |
| (0020,9164) | Dimension Organization UID |
| (0028,1199) | Palette Color Lookup Table UID |
| (0028,1214) | Large Palette Color Lookup Table UID |
| (0028,4000) | Image Presentation Comments |
| (0032,0012) | Study ID Issuer |
| (0032,1020) | Scheduled Study Location |
| (0032,1021) | Scheduled Study Location AE Title |
| (0032,1030) | Reason for Study |
| (0032,1032) | Requesting Physician |
| (0032,1033) | Requesting Service |
| (0032,1060) | Requested Procedure Description |
| (0032,1070) | Requested Contrast Agent |
| (0032,4000) | Study Comments |
| (0038,0004) | Referenced Patient Alias Sequence |
| (0038,0010) | Admission ID |
| (0038,0011) | Issuer of Admission ID |
| (0038,001E) | Scheduled Patient Institution Residence |
| (0038,0020) | Admitting Date |
| (0038,0021) | Admitting Time |
| (0038,0040) | Discharge Diagnosis Description |
| (0038,0050) | Special Needs |
| (0038,0060) | Service Episode ID |
| (0038,0061) | Issuer of Service Episode ID |
| (0038,0062) | Service Episode Description |
| (0038,0300) | Current Patient Location |
| (0038,0400) | Patient's Institution Residence |
| (0038,0500) | Patient State |
| (0038,4000) | Visit Comments |
| (0040,0001) | Scheduled Station AE Title |
| (0040,0002) | Scheduled Procedure Step Start Date |
| (0040,0003) | Scheduled Procedure Step Start Time |
| (0040,0004) | Scheduled Procedure Step End Date |
| (0040,0005) | Scheduled Procedure Step End Time |
| (0040,0006) | Scheduled Performing Physician Name |
| (0040,0007) | Scheduled Procedure Step Description |
| (0040,000B) | Scheduled Performing Physician Identification Sequence |
| (0040,0010) | Scheduled Station Name |
| (0040,0011) | Scheduled Procedure Step Location |
| (0040,0012) | Pre-Medication |
| (0040,0241) | Performed Station AE Title |
| (0040,0242) | Performed Station Name |
| (0040,0243) | Performed Location |
| (0040,0244) | Performed Procedure Step Start Date |
| (0040,0245) | Performed Procedure Step Start Time |
| (0040,0248) | Performed Station Name Code Sequence |
| (0040,0250) | Performed Procedure Step End Date |
| (0040,0251) | Performed Procedure Step End Time |
| (0040,0253) | Performed Procedure Step ID |
| (0040,0254) | Performed Procedure Step Description |
| (0040,0275) | Request Attributes Sequence |
| (0040,0280) | Comments on Performed Procedure Step |
| (0040,0555) | Acquisition Context Sequence |
| (0040,1001) | Requested Procedure ID |
| (0040,1004) | Patient Transport Arrangements |
| (0040,1005) | Requested Procedure Location |

| DICOM Tag | Attribute Name |
|-------------|---|
| (0040,1010) | Names of Intended Recipient of Results |
| (0040,1011) | Intended Recipients of Results Identification Sequence |
| (0040,1101) | Person Identification Code Sequence |
| (0040,1102) | Person Address |
| (0040,1103) | Person Telephone Numbers |
| (0040,1400) | Requested Procedure Comments |
| (0040,2001) | Reason for Imaging Service Request |
| (0040,2008) | Order Entered By |
| (0040,2009) | Order Enterer Location |
| (0040,2010) | Order Callback Phone Number |
| (0040,2016) | Placer Order Number of Imaging Service Request |
| (0040,2017) | Filler Order Number of Imaging Service Request |
| (0040,2400) | Imaging Service Request Comments |
| (0040,3001) | Confidentiality Constraint on Patient Data Description |
| (0040,4005) | Scheduled Procedure Step Start DateTime |
| (0040,4010) | Scheduled Procedure Step Modification DateTime |
| (0040,4011) | Expected Completion Date Time |
| (0040,4023) | Referenced General Purpose Scheduled Procedure Step Transaction UID |
| (0040,4025) | Scheduled Station Name Code Sequence |
| (0040,4027) | Scheduled Station Geographic Location Code Sequence |
| (0040,4030) | Performed Station Geographic Location Code Sequence |
| (0040,4034) | Scheduled Human Performers Sequence |
| (0040,4035) | Actual Human Performers Sequence |
| (0040,4036) | Human Performers Organization |
| (0040,4037) | Human Performers Name |
| (0040,4050) | Performed Procedure Step Start DateTime |
| (0040,4051) | Performed Procedure Step End DateTime |
| (0040,4052) | Procedure Step Cancellation DateTime |
| (0040,A027) | Verifying Organization |
| (0040,A073) | Verifying Observer Sequence |
| (0040,A075) | Verifying Observer Name |
| (0040,A078) | Author Observer Sequence |
| (0040,A07A) | Participant Sequence |
| (0040,A07C) | Custodial Organization Sequence |
| (0040,A088) | Verifying Observer Identification Code Sequence |
| (0040,A123) | Person Name |
| (0040,A124) | UID |
| (0040,A171) | Observation UID |
| (0040,A172) | Referenced Observation UID (Trial) |
| (0040,A192) | Observation Date (Trial) |
| (0040,A193) | Observation Time (Trial) |
| (0040,A307) | Current Observer (Trial) |
| (0040,A352) | Verbal Source (Trial) |
| (0040,A353) | Address (Trial) |
| (0040,A354) | Telephone Number (Trial) |
| (0040,A358) | Verbal Source Identifier Code Sequence (Trial) |
| (0040,A402) | Observation Subject UID (Trial) |
| (0040,A730) | Content Sequence |
| (0040,DB0C) | Template Extension Organization UID |
| (0040,DB0D) | Template Extension Creator UID |
| (0070,0001) | Graphic Annotation Sequence |
| (0070,0084) | Content Creator's Name |
| (0070,0086) | Content Creator's Identification Code Sequence |
| (0070,031A) | Fiducial UID |
| (0088,0140) | Storage Media Fileset UID |
| (0088,0200) | Icon Image Sequence |
| (0088,0904) | Topic Title |
| (0088,0906) | Topic Subject |
| (0088,0910) | Topic Author |
| (0088,0912) | Topic Keywords |
| (0400,0100) | Digital Signature UID |
| (0400,0402) | Referenced Digital Signature Sequence |
| (0400,0403) | Referenced SOP Instance MAC Sequence |
| (0400,0404) | MAC |
| (0400,0550) | Modified Attributes Sequence |
| (0400,0561) | Original Attributes Sequence |
| (2030,0020) | Text String |
| (3006,0024) | Referenced Frame of Reference UID |
| (3006,00C2) | Related Frame of Reference UID |

| DICOM Tag | Attribute Name |
|-------------|--------------------------------------|
| (3008,0105) | Source Serial Number |
| (300A,0013) | Dose Reference UID |
| (300E,0008) | Reviewer Name |
| (4000,0010) | Arbitrary |
| (4000,4000) | Text Comments |
| (4008,0042) | Results ID Issuer |
| (4008,0102) | Interpretation Recorder |
| (4008,010A) | Interpretation Transcriber |
| (4008,010B) | Interpretation Text |
| (4008,010C) | Interpretation Author |
| (4008,0111) | Interpretation Approver Sequence |
| (4008,0114) | Physician Approving Interpretation |
| (4008,0115) | Interpretation Diagnosis Description |
| (4008,0118) | Results Distribution List Sequence |
| (4008,0119) | Distribution Name |
| (4008,011A) | Distribution Address |
| (4008,0202) | Interpretation ID Issuer |
| (4008,0300) | Impressions |
| (4008,4000) | Results Comments |
| (50xx,xxxx) | Curve Data |
| (60xx,0100) | Overlay Bits Allocated |
| (60xx,0102) | Overlay Bit Position |
| (60xx,3000) | Overlay Data |
| (60xx,4000) | Overlay Comments |
| (FFFA,FFFA) | Digital Signatures Sequence |
| (FFFC,FFFC) | Data Set Trailing Padding |

Table 3-7: Application Level Confidentiality Profile Attributes (private tags)

| DICOM Tag | Attribute Name |
|--|----------------------------------|
| (0019,0005) | Multiphase UID |
| (0019, SIEMENS CT VA0 COAD, 90) | Osteo offset |
| (0019, SIEMENS CT VA0 COAD, 92) | Osteo Regression Line Slope |
| (0019, SIEMENS CT VA0 COAD, 93) | Osteo Regression Line Intercept |
| (0019, SIEMENS CT VA0 COAD, 96) | Osteo Phantom Number |
| (0021, SIEMENS MR SDS 01, 19) | MR Phoenix Protocol |
| (0029, SIEMENS CSA ENVELOPE, 10) | Syngo Report Data |
| (0029, SIEMENS CSA ENVELOPE, 11) | Syngo Report Presentation |
| (0029, SIEMENS CSA HEADER, 08) | Modality Image Header Type |
| (0029, SIEMENS CSA HEADER, 09) | Modality Image Header Version |
| (0029, SIEMENS CSA HEADER, 10) | Modality Image Header Info |
| (0029, SIEMENS CSA HEADER, 18) | Modality Series Header Type |
| (0029, SIEMENS CSA HEADER, 19) | Modality Series Header Version |
| (0029, SIEMENS CSA HEADER, 20) | Modality Series Header Info |
| (0029, SIEMENS MEDCOM HEADER, 40) | Application Header Sequence |
| (0029, SIEMENS MEDCOM HEADER, 41) | Application Header Type |
| (0029, SIEMENS MEDCOM HEADER, 42) | Application Header ID |
| (0029, SIEMENS MEDCOM HEADER, 43) | Application Header Version |
| (0029, SIEMENS MEDCOM HEADER, 44) | Application Header Info |
| (0029, SIEMENS CT APPL DATASET, 00) | Dual Energy Algorithm Parameters |
| (0029, SIEMENS CT APPL ALG PARAMS, 20) | Perfusion Result Set Id |
| (0029, SIEMENS CT EXAM IMAGE, 49) | Metal Artifact Reduction Type |
| (0043, GEMS_PARM_01, 1E) | GE Delta Start Time |

4 SUPPORT OF CHARACTER SETS

4.1 CHARACTER SETS FOR *syngo.via* View&GO

The ***syngo.via* View&GO** DICOM application supports the following character sets as defined in the three tables below.

Table 4-1: Single-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|----------------------|
| Default repertoire | none | ISO_IR 6 | ISO 646: |
| Latin alphabet No. 1 | ISO_IR 100 | ISO_IR 100 | Supplementary set |
| | | ISO_IR 6 | ISO 646: |
| Latin alphabet No. 2 | ISO_IR 101 | ISO_IR 101 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 3 | ISO_IR 109 | ISO_IR 109 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 4 | ISO_IR 110 | ISO_IR 110 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Cyrillic | ISO_IR 144 | ISO_IR 144 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Arabic | ISO_IR 127 | ISO_IR 127 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Greek | ISO_IR 126 | ISO_IR 126 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Hebrew | ISO_IR 138 | ISO_IR 138 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 5 | ISO_IR 148 | ISO_IR 148 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Japanese | ISO_IR 13 | ISO_IR 13 | JIS X 0201: Katakana |
| | | ISO_IR 14 | JIS X 0201: Romaji |
| Thai | ISO_IR166 | ISO_IR166 | TIS 620-253 (1990) |
| | | ISO_IR 6 | ISO 646 |

Table 4-2: Single-Byte Characters Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------|-------------------------|-------------------|
| Default repertoire | ISO 2022 IR 6 | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.1 | ISO 2022 IR 100 | ISO 2022 | ESC 02/13 04/01 | ISO-IR 100 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.2 | ISO 2022 IR 101 | ISO 2022 | ESC 02/13 04/02 | ISO-IR 101 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.3 | ISO 2022 IR 109 | ISO 2022 | ESC 02/13 04/03 | ISO-IR 109 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.4 | ISO 2022 IR 110 | ISO 2022 | ESC 02/13 04/04 | ISO-IR 110 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |

Multi-Byte Character Sets without Code Extension

Table 4-3: Multi-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|---|
| Unicode | ISO_IR 192 | ISO 10646 | Unicode in UTF-8 |
| Chinese | GB18030 | GB18030 | GB 18030-2000 (China Association for Standardization) |

Table 4-4: Multi-Byte Character Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------------|-------------------------|-------------------------------------|
| Japanese | ISO 2022 IR 159 | ISO 2022 | ESC 02/04 02/08 04/04 | ISO-IR 159 ISO-IR 87 | JIS X 0212: Supplementary Kanji set |
| Korean | ISO 2022 IR 149 | ISO 2022 | ESC 02/04 02/09 04/03 | ISO-IR 149 | KS X 1001: Hangul and Hanja |

All Specific Character Set listed above are supported for incoming Data. When creating new Instances, the system will use the default SCS (or SCS List) configured on the machine.

When there is a mismatch between the given character set in attribute (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM conform:

- Convert each illegal character to a '?'.

There are now three categories of character sets which have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO_IR 6, ISO 2022 IR 6, ISO_IR 100, etc.
→ encoded in ISO 2022
- ISO_IR 192 → encoded in UTF-8
- GB18030 → encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO_IR 192 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains ISO_IR 192
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies is not DICOM conform. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

The **syngo.via View&GO** supports Kanji characters in the byte zones after 74 (79, 7A, 7B and 7C).

5 SECURITY

5.1 SECURITY PROFILES

Time Synchronization Profiles: **syngo.via View&GO** acts as an NTP Client in the Maintain Time Transaction.

5.2 ASSOCIATION LEVEL SECURITY

The Echo SCP will answer to AETs.

6 Annexes

6.1 SOP Classes supported

Table 6-1 SOP CLASSES for Storage

| SOP Class Name | SOP Class UID | User of Service (SCU) | Provider of Service (SCP) | Visualization |
|--|--------------------------------|-----------------------|---------------------------|---------------|
| Supported Storage SOP Classes | | | | |
| 12-lead ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.1 | Yes | Yes | No |
| Ambulatory ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.3 | Yes | Yes | No |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | Yes | Yes | No |
| Blending Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.4 | Yes | Yes | No |
| Breast Tomosynthesis Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.3 | Yes | Yes | Yes |
| Cardiac Electrophysiology Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.3.1 | Yes | Yes | No |
| Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.2 | Yes | Yes | No |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | Yes | Yes | No |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Yes | Yes | Yes |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Yes | Yes | Yes |
| Deformable Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.3 | Yes | Yes | No |
| Digital Mammography X-Ray Image Storage- For Presentation | 1.2.840.10008.5.1.4.1.1.1.2 | Yes | Yes | Yes |
| Digital Mammography X-Ray Image Storage- For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | Yes | Yes | No |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Yes | Yes | Yes |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | Yes | Yes | No |
| Encapsulated PDF Storage | 1.2.840.10008.5.1.4.1.1.104.1 | Yes | Yes | No |
| Enhanced CT Image Storage | 1.2.840.10008.5.1.4.1.1.2.1 | Yes | Yes | Yes |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | Yes | Yes | Yes |
| Enhanced MR Color Image Storage | 1.2.840.10008.5.1.4.1.1.4.3 | Yes | Yes | Yes |
| Enhanced Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.22 | Yes | Yes | No |
| Enhanced XA Image Storage | 1.2.840.10008.5.1.4.1.1.12.1.1 | Yes | Yes | Yes |
| Enhanced XRF Image Storage | 1.2.840.10008.5.1.4.1.1.12.2.1 | Yes | Yes | No |
| General ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.2 | Yes | Yes | No |
| Grayscale Softcopy Presentation State Storage (store & forward only) | 1.2.840.10008.5.1.4.1.1.11.1 | Yes | Yes | Yes |
| Hemodynamic Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.2.1 | Yes | Yes | No |
| Key Object Selection Document Storage | 1.2.840.10008.5.1.4.1.1.88.59 | Yes | Yes | Yes |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | Yes | Yes | Yes |

| SOP Class Name | SOP Class UID | User of Service (SCU) | Provider of Service (SCP) | Visualization |
|--|--------------------------------|-----------------------|---------------------------|---------------|
| Mammography CAD SR Storage | 1.2.840.10008.5.1.4.1.1.88.50 | Yes | Yes | No |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | Yes | Yes | Yes |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2 | Yes | Yes | Yes |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3 | Yes | Yes | Yes |
| Multi-frame Single Bit Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.1 | Yes | Yes | Yes |
| Multi-frame True Color Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | Yes | Yes | Yes |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | Yes | Yes | Yes |
| Pseudo-Color Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.3 | Yes | Yes | No |
| PET Image Storage | 1.2.840.10008.5.1.4.1.1.128 | Yes | Yes | Yes |
| Procedure Log Storage | 1.2.840.10008.5.1.4.1.1.88.40 | Yes | Yes | No |
| Raw DataStorage | 1.2.840.10008.5.1.4.1.1.66 | Yes | Yes | No |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | Yes | Yes | No |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | Yes | Yes | No |
| RT Dose Storage | 1.2.840.10008.5.1.4.1.1.481.2 | Yes | Yes | Yes |
| RT Image Storage | 1.2.840.10008.5.1.4.1.1.481.1 | Yes | Yes | Yes |
| RT Ion Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.9 | Yes | Yes | No |
| RT Ion Plan Storage | 1.2.840.10008.5.1.4.1.1.481.8 | Yes | Yes | No |
| RT Plan Storage | 1.2.840.10008.5.1.4.1.1.481.5 | Yes | Yes | Yes |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | Yes | Yes | Yes |
| RT Treatment Summary Record Storage | 1.2.840.10008.5.1.4.1.1.481.7 | Yes | Yes | No |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Yes | Yes | Yes |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | Yes | Yes | No |
| Spatial Fiducials Storage | 1.2.840.10008.5.1.4.1.1.66.2 | Yes | Yes | No |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | Yes | Yes | No |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | Yes | Yes | No |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | Yes | Yes | Yes |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | Yes | Yes | Yes |
| Ultrasound Multi-frame Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.3 | Yes | Yes | Yes |
| X-Ray 3D Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.1 | Yes | Yes | Yes |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | Yes | Yes | Yes |
| X-Ray Radiation Dose Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.67 | Yes | Yes | No |
| X-Ray Radio-Fluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | Yes | Yes | Yes |
| Supported private Storage SOP Classes | | | | |
| CSA Non-Image Storage | 1.3.12.2.1107.5.9.1 | No | Yes | No |

Table 6-2: Supported Non-Storage SOP Classes

| SOP Class Name | SOP Class UID | User of Service (SCU) | Provider of Service (SCP) |
|---|-----------------------------|-----------------------|---------------------------|
| Supported Verification SOP Classes | | | |
| Verification | 1.2.840.10008.1.1 | Yes | Yes |
| Supported Storage Commitment SOP Classes | | | |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | No | No |
| Storage Commitment Push Model well known SOP Instance | 1.2.840.10008.1.20.1.1 | No | No |
| Supported Query/Retrieve-FIND SOP Classes | | | |
| Patient Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Yes | Yes |
| Study Root Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | Yes |
| Patient /Study Only Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Yes | Yes |
| Supported Query/Retrieve-MOVE SOP Classes | | | |
| Study Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Yes | Yes |
| Patient Root Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | Yes | Yes |
| Patient/Study Only Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | Yes | Yes |
| Modality Worklist Information SOP Class | | | |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Yes | No |
| Grayscale Print Management META SOP classes | | | |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | Yes | No |
| - Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Yes | No |
| - Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Yes | No |
| - Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | Yes | No |
| - Printer SOP Class | 1.2.840.10008.5.1.1.16 | Yes | No |
| Print Job SOP Class | 1.2.840.10008.5.1.1.14 | Yes | No |
| Presentation LUT SOP Class | 1.2.840.10008.5.1.1.23 | Yes | No |
| Color Print Management META SOP classes | | | |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 | Yes | No |
| - Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Yes | No |
| - Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Yes | No |
| - Basic Color Image Box SOP Class | 1.2.840.10008.5.1.1.4.1 | Yes | No |
| - Printer SOP Class | 1.2.840.10008.5.1.1.16 | Yes | No |
| Print Job SOP Class | 1.2.840.10008.5.1.1.14 | Yes | No |

6.2 IOD CONTENTS

6.2.1 Created SOP Instance(s)

The applications from **syngo.via View&GO** create objects of the following SOP Classes during Transferring, Post-Processing and Reading:

Table 6-3: List of created SOP Classes

| SOP Class Name | SOP Class UID | Internally used (neither SCU nor SCP is applicable) |
|---|-------------------------------|---|
| Basic Text Structured Report | 1.2.840.10008.5.1.4.1.1.88.11 | No |
| Comprehensive Structured Report Storage | 1.2.840.10008.5.1.4.1.1.88.33 | No |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | No |
| Encapsulated PDF Storage | 1.2.840.10008.5.1.4.1.1.104.1 | No |
| Enhanced Structured Report | 1.2.840.10008.5.1.4.1.1.88.22 | No |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | No |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | No |
| Raw Data Storage | 1.2.840.10008.5.1.4.1.1.66 | No |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | No |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | No |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | No |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | No |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | No |

6.2.2 Data Dictionary of Private Attributes

The following table Table 6-4: Private Data Element Dictionary lists all private attributes created by **syngo.via View&GO**, which may be included in the generated instances. These private attributes may be deprecated or replaced with standard attributes in the future.

Table 6-4: Private Data Element Dictionary

| DICOM Tag | Name | VR | VM |
|---|---------------------------------|----|----|
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,01) | Business Unit Code | CS | 1 |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,02) | Application Type | LO | 1 |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,03) | Application Attributes Sequence | SQ | 1 |

| DICOM Tag | Name | VR | VM |
|---|--------------------------------|----|-----|
| (0029,SIEMENS SYNGO FUNCTION ASSIGNMENT,01) | Data Reference | LO | 1 |
| (0009,SIEMENS SYNGO INDEX SERVICE,20) | Object Insertion Date | DA | 1 |
| (0009,SIEMENS SYNGO INDEX SERVICE,A0) | Sender System Device Name | LO | 1 |
| (0029,SIEMENS SYNGO VOLUME,12) | Slices | US | 1 |
| (0029,SIEMENS SYNGO VOLUME,14) | Volume Histogram | OB | 1 |
| (0029,SIEMENS SYNGO VOLUME,18) | Volume Level | IS | 1 |
| (0029,SIEMENS SYNGO VOLUME,30) | Voxel Spacing | DS | 3 |
| (0029,SIEMENS SYNGO VOLUME,32) | Volume Position (Patient) | DS | 3 |
| (0029,SIEMENS SYNGO VOLUME,37) | Volume Orientation (Patient) | DS | 9 |
| (0029,SIEMENS SYNGO VOLUME,40) | Resampling Flag | CS | 1 |
| (0029,SIEMENS SYNGO VOLUME,42) | Normalization Flag | CS | 1 |
| (0029,SIEMENS SYNGO VOLUME,44) | SubVolume Sequence | SQ | 1-n |
| (0071,SIEMENS SYNGO REGISTRATION,20) | Registered Image Sequence | SQ | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,21) | Registration Is Validated Flag | CS | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,20) | Registered Image Sequence | SQ | 1 |
| (0071,SIEMENS SYNGO REGISTRATION,21) | Registration Is Validated Flag | CS | 1 |
| (7FDF,SIEMENS SYNGO DATA PADDING,FC) | Pixel Data Leading Padding | OB | 1 |

Interpretation of the Dicom Tags from the above table:

(gggg, pp,ee) -> (gggg, ppee)

gggg - odd group number

pp - private creator identification code

ee - private element

6.2.3 Usage of Attributes from received IODs

N/A

6.2.4 Coerced / Modified fields

N/A

6.3 CODED TERMINOLOGY AND TEMPLATES

See application specific annexes.

6.3.1 Context Groups

See application specific annexes.

6.3.2 Template Specifications

See application specific annexes.

6.3.3 Private Code definitions

See application specific annexes.

6.4 GRAYSCALE IMAGE CONSISTENCY

N/A

6.5 STANDARD EXTENDED / SPECIALIZED / PRIVATE SOP CLASSES

N/A

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Legal Manufacturer
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

Siemens Healthcare Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Telephone +49 9131 84-0
[siemens.com/healthcare](https://www.siemens.com/healthcare)