

# COVID-19 Laboratory Testing Guide

siemens-healthineers.com/covid-19









#### What is COVID-19?

COVID-19 (coronavirus disease 2019) is the disease resulting from infection with a newly emerged coronavirus named SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2).<sup>1</sup> Coronaviruses are a family of RNA viruses usually found in animals. Mutations in the virus can result in human infection and subsequent spread.<sup>2</sup>

SARS-CoV-2 is closely related to the SARS virus identified in 2003 (SARS-CoV) and (to a lesser degree) the MERS-CoV virus from 2012. All three can produce a severe respiratory syndrome and associated mortality.<sup>2-4</sup> While both SARS-CoV and the MERS-CoV viruses seem to have a higher comparative rate of mortality (especially MERS-CoV), the newly emerged SARS-CoV-2 appears much more infectious, with significant human-to-human transmission.<sup>5-7</sup> Asymptomatic individuals may transmit the virus, challenging infection control.<sup>8,9</sup> Its rapid spread has produced a true global pandemic.

#### Diagnosis of COVID-19:

Diagnosis cannot be made solely on signs or symptoms as these overlap with other respiratory illness, so confirmation of the presence of the virus is essential. Table 1 describes the range and percent of symptoms seen in confirmed COVID-19 infections.

Sign or symptom	% of patients
Fever	83–99
Cough	59–82
Fatigue	44–70
Anorexia	40-84
Sputum production	28-33
Shortness of breath	31–40
Myalgia (muscle aches)	11–35

**Table 1.** Range and percent of symptoms seen in confirmed COVID-19 infections.<sup>13</sup>

In some populations, a loss of taste or smell is also a widely reported symptom.<sup>10</sup> Molecular testing specific for the SARS-CoV-2 RNA is used to confirm presence of the virus, and serology testing can identify antibodies to the virus.<sup>11,12</sup> In addition, antibody testing is a promising approach to assess prevalence of infection and potentially identify immunity. It remains to be confirmed whether antibodies to SARS-CoV-2 offer protection (immunity) from subsequent exposure.

#### Conclusion:

Testing is critical to differentiate COVID-19 from other respiratory diseases such as influenza or RSV (respiratory syncytial virus). Widespread testing will both inform a COVID-19 diagnosis and aid a greater understanding of disease prevalence, especially in infections that are asymptomatic.



### Clinical classification and transmission of COVID-1914-27

# **Clinical classification**



### Human-to-human transmission



#### **Disease spectrum**



Incubation period median 5 days, range 1–27 days Patients with mild symptoms in week 1...

...may progress in week 2

Progression may be rapid and sudden.

Possible transmission during asymptomatic phase 4–6 days before the onset of symptoms.

Possible transmission

Highest risk of transmission during symptomatic phase

Possible transmission after remission of the symptoms

Possible transmission

Remission

#### **Asymptomatic**

No symptoms\* **1.2–17.6%** 

Infectious virus can be shed, viral loads may be comparable to symptomatic patients

### Mild to Moderate

80%

No signs (mild) or some signs (moderate) of pneumonia on imaging

#### **Symptomatic**

Severe **10–15**%

#### Critical

2-5%

- Respiratory distress of ≥30 breath per minute
- Oxygen saturation ≤93% at rest
- PaO2/FiO2 ≤300 mmHg
- Lung lesion progression >50% in 24–48h
- Respiratory failure requiring mechanical ventilation
- Shock
- Any organ failure requiring ICU care
- Case/fatality rate current estimation 0.3 to 4%

<sup>\*</sup>Recent SARS-CoV-2 testing studies demonstrate a wide variety of asymptomatic patients between 1% and >80% partially due to different stages of outbreak, testing availability and population studied

### Course of patients at risk and/or infected by COVID-1914-27



# Identification of high risk patients



- Age >65 years
- Any age with severe obesity BMI >40
- Pre-existing pulmonary disease (chronic lung disease or moderate to severe asthma, COPD, lung cancer, pulmonary hypertension, emphysema (smoking, A1AT deficiency))
- Chronic kidney disease
- Diabetes
- History of:
- hypertension (treated and untreated)
- cardiovascular disease
- liver disease
- transplants or other immunosuppression (ex. cancer treatment)
- Patients with HIV or other immunocompromised conditions
- Patients with endocrine pathologies
- Use of biologic drugs



# Population triage<sup>2,3,4</sup>

## **Epidemiologic surveillance Laboratory testing:**

- RT-PCR
- Anti-SARS-CoV-2 antibodies



# Diagnosis and Determination of Disease Status

#### **Disease detection:**

• Molecular testing (RT-PCR)

#### **Immune Response:**

• Anti-SARS-CoV-2 antibodies (IgG, Total)

#### **Imaging:**

- Chest CT
- Ultrasound
- X-ray



#### Recovery

#### **Laboratory testing:**

 Anti-SARS-CoV-2 antibodies (IgG, Total)

### Confirm immunization

(under investigation)



#### Monitoring

#### **Laboratory testing:**

- Arterial blood gas
- Complete blood count
- Acute phase proteins, inflammation & coagulation
- Liver, kidney, cardiac function biomarkers
- Additional testing related to comorbidities

#### **Imaging:**

- Chest CT (ground glass opacification)
- Ultrasound

# The essential role of laboratory diagnostics in SARS-CoV-2 infection<sup>14-27</sup>



#### **Essential lab testing**

Viral serologies HBV, HCV, HIV Ab/Ag

Daily labs

CBC with differential (trend total lymphocyte count)
Comprehensive metabolic panel:
Electrolytes: Na, K, Total CO<sub>2</sub>, Chloride
Total protein and Albumin
Creatinine
Bilirubin, ALT, AST
CPK (total creatine kinase)
Lactate

D-dimer, Ferritin, CRP, ESR, LDH, Cardiac troponin

# Frequent laboratory abnormalities in patients with COVID-19<sup>†</sup>



#### **Decreased**

Blood lymphocyte count (35–75%) Albumin (50–98%) Hemoglobin (41–50%)



#### Increased

Neutrophil count
Erythrocyte sedimentation rate (ESR; up to 85%)
C-reactive protein (CRP; 75–93%)

Lactate dehydrogenase (LDH; 27–92%)

Alanine aminotransferase (ALT)

Aspartate

aminotransferase (AST) Total bilirubin

Cardiac troponin

Procalcitonin (6–25%) Prothrombin time (PT)

D-dimer (36-43%)

†Approximate percentage of patients

Additional essential lab testing	Test	Potential clinical significance
Arterial blood gas	pH, PaCO2, PaO2, and aHCO3	For ventilator adjustments
Hematology	Lymphopenia with atypical lymphocytes Leukocytosis, Neutrophilia, low eosinophils Thrombocytopenia: Platelet count	Decreased immunological response to the virus Bacterial (super) infection Consumption (disseminated) coagulopathy
Hemostasis	Prothrombin time, D-Dimer	Activation of blood coagulation and/or disseminated coagulopathy PT and D-dimer are significant predictors of disease severity
Inflammation/Infection	CRP, Ferritin, IL6, TNFα, SAA Procalcitonin	Severe viral infection/viremia Bacterial (super) infection
Cardiac	High-sensitivity troponin, CK-MB, BNP/NT-proBNP	Increased levels may be associated with higher mortality <sup>‡</sup>
Liver	Albumin, ALT, AST, Bilirubin	Impairment of liver function, Liver injury
Renal	Creatinine, Cystatin C	Kidney injury

## **Impact of Comorbidities on COVID-19 Patients**

Patients with comorbidities (like diabetes mellitus, hypertension, cardiovascular, chronic lung and chronic kidney disease) are particularly susceptible to COVID-19 infection and are likely to have more severe illness<sup>14-26</sup>

	Comorbidities	Additional testing	Impact of COVID-19
6	Cardiovascular disease	Troponin, Natriuretic peptides, CKMB	Precipitates cardiac complications like: acute heart failure, myocardial infarction, myocardial injury, cardiac arrest.
6,3	Chronic kidney disease	Blood: Creatinine, Cystatin C, eGFR Urine: Albumin	Challenges for patients on dialysis, in particular, in-center hemodialysis; uremic patients are particularly vulnerable to infection and may exhibit greater variations in clinical symptoms and infectivity.
	Heart/liver/kidney transplant	Immunosuppressant Drugs: Mycophenolate, Cyclosporine, Tacrolimus, Sirolimus, Everolimus	Patients may be more vulnerable due to immunocompromised status.
楽	Viral co-infection	Hepatitis B serologies (anti-HBs, anti-HBc, and HBsAg) Hepatitis C serology (anti-HCV), unless positive in past HIV 1/2/O, CD4 count	Viral serologies assist in interpretation of ALT elevations, present in ~25% of COVID-19 patients. HIV patients may get severe side effects when taking Tocilizumab (drug being used for COVID-19 pneumonia).§ Hepatitis patients are at higher risk for liver complications.
	Diabetes	Blood glucose	Patients with diabetes who are infected with COVID-19 may see their glycemic control deteriorate during the illness.
	Chronic lung disease	LDH	Patients may be more vulnerable due to lung function insufficiency.
	Chronic liver disease	Albumin, AST, ALT, Total Protein, Bilirubin, PT INR	Patients may be more vulnerable due to liver function insufficiency.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Siemens Healthineers' SARS-CoV-2 molecular and serology tests have not been FDA cleared or approved. These tests have been authorized by FDA under an EUA for use by authorized laboratories. The molecular test has been authorized only for detecting the presence of antibodies against SARS-CoV-2, not for any other viruses or pathogens. The serology test has been authorized only for detecting the presence of antibodies against SARS-CoV-2, not for any other viruses or pathogens. The serology test has been authorized only for detecting the presence of antibodies against SARS-CoV-2, not for any other viruses or pathogens. These tests are only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

#### Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone: +49 9131 84-0 siemens-healthineers.com

#### Published by

Siemens Healthcare Diagnostics Inc. Laboratory Diagnostics 511 Benedict Avenue Tarrytown, NY 10591-5005 USA

Phone: +1 914-631-8000

#### References:

- 1. Lai CC, Shih TP, Ko WC, et al. Int J Antimicrob Agents. 2020 Mar;55(3):105924.
- 2. Chen, Y. et al. J Med Virol. 2020 Apr;92(4):418-423
- 3. Guarner J. Am J Clin Pathol. 2020 Mar 9;153(4):420-421.
- 4. Liu J. et al. J Med Virol. 2020;92:491-494
- Riou J, Althous CL. Euro Surveill. 2020;25(4):2000058. https://www.eurosurveillance.org/content/10.2807/1560-7917.ES.2020.25.4.2000058
- 6. Liu Y, Gayle AA, Wilder-Smith A, et al. J Travel Med. 2020;27(2). https://doi.org/10.1093/jtm/taaa021
- 7. Guo YR et al. Mil Med Res. 2020;7(1):11. https://mmrjournal.biomedcentral.com/articles/10.1186/s40779-020-00240-0
- 8. Bai Y. et al. JAMA. Published online February 21, 2020. doi:10.1001/jama.2020.2565
- https://www.cdc.gov/coronavirus/2019-ncov/hcp/faq.html#transmission. Accessed April 6, 2020.
- https://www.the-hospitalist.org/hospitalist/article/219556/coronavirus-updates/reportsincreasingly-suggest-anosmia/hyposmia-can Accessed April 6, 2020.
- 11. https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/testing.html Accessed April 6, 2020.
- 12. https://www.cdc.gov/coronavirus/2019-ncov/about/testing.html Accessed April6, 2020.
- 13. https://www.cdc.gov/coronavirus/2019-ncov/hcp/clinical-guidance-management-patients.html
- 14. World Health Organization. Coronavirus disease 2019 (COVID-19) Situation Report 44. https://www.who.int/emergencies/diseases/novel-coronavirus-2019/situation-reports
- 15. https://www.ifcc.org/ifcc-news/2020-03-26-ifcc-information-guide-on-covid-19/
- 16. Report of the WHO-China Joint Mission on Coronavirus Disease 2019 (COVID-19)
- 17. March 22, 2020 Content source: National Cent Immunization and Respiratory Diseases (NCIRD) Division of Viral Diseases
- 18. JAMA.2020 Feb7. doi:10,1001/jama,20201585
- 19. Ann Transl Med 2020;8(3):48 | http://dx.doi.org/10.21037/atm.2020.02.06
- 20. Lippi G, PlebaniM. Laboratory abnormalities in patients with COVID-2019 infection. Clin Chem Lab Med. 2020 Mar 3. doi: 10.1515/cclm-2020-0198.
- 21. https://www.massgeneral.org/news/coronavirus/coronavirus-latest-updates
- 22. Jingyuan Liu, https://doi.org/10.1101/2020.02.10.20021584
- 23. Hematologic parameters in patients with COVID-19 infection 10.1002/ajh.25774
- 24. Ruan et al: https://doi.org/10.1007/s00134-020-05991-x
- 25. Boettler T et al. Care of patients with liver disease during the COVID-19 pandemic: EASL-ESCMID Position Paper. https://easl.eu/wp-content/uploads/2020/04/EASL-ESCMID-Position-Paper-on-COVID-19-and-the-liver-2-April-2020.pdf April 2, 2020.
- 26. ifcc.org/ifcc-news/2020-03-26-ifcc-information-guide-on-covid-19/
- 27. https://www.cebm.net/covid-19/covid-19-what-proportion-are-asymptomatic/