


Angiogenic Factors

Key Biomarkers for Pregnancy Management in Women at Risk of Preeclampsia



Background

Bringing a child safely into the world is one of life's most profound moments, yet preeclampsia (PE) remains a serious threat to mothers and babies worldwide. PE is a complex, life-threatening pregnancy disorder that typically develops after 20 weeks of gestation, marked by new-onset hypertension and organ dysfunction (See table 1 for detailed diagnostic criteria for PE). For mothers, complications can include eclampsia, stroke, HELLP syndrome, organ failure, and even death. For the fetus, risks include growth restriction, prematurity, and loss. Women with a history of PE face a 15–20% chance of recurrence¹ and an increased long-term risk of cardiovascular disease² and diabetes.³

PE has also a significant public health impact globally. It has been estimated that it affects approximately 2–8% of pregnancies worldwide. Globally, there are around 46,000 maternal deaths per year and 0.5 million fetal or newborn deaths per year due to PE.⁴ The incidence varies by region with higher rates in low and middle-income countries at around 10% of maternal deaths in Asia and Africa and 25% in Latin America due to preeclampsia and eclampsia.⁴ In the United States, Canada and Western Europe, the incidence of PE is reported to be between 2% and 5%.⁵ In France, Olié et al. reported that PE complicates around 2% of pregnancies.⁶

Table 1. Definition of preeclampsia (de novo) by the International Society for the Study of Hypertension in Pregnancy¹

Gestational hypertension (sBP \geq 140 mm Hg and/or dBP \geq 90 mm Hg) accompanied by one or more of the following new-onset conditions at \geq 20 weeks' gestation:

1. Proteinuria: 24-h urine protein \geq 300 mg/d or spot urine protein/creatinine ratio \geq 0.30 mg/mg or urine dipstick testing \geq 2+
2. Other maternal end-organ dysfunction, including:
 - Neurological complications (e.g., eclampsia, altered mental status, blindness, stroke, clonus, severe headaches, or persistent visual scotomata)
 - Pulmonary oedema
 - Hematological complications (e.g., platelet count $<$ 150,000/ μ L, Disseminated Intravascular Coagulation, hemolysis)
 - Acute kidney injury (creatinine \geq 90 μ mol/L; 1 mg/dL)
 - Liver involvement (e.g., elevated transaminases such as ALT or AST $>$ 40 IU/L) with or without right upper quadrant or epigastric abdominal pain)
3. Uteroplacental dysfunction (e.g., placental abruption, angiogenic imbalance, fetal growth restriction, abnormal umbilical artery Doppler waveform analysis, or intrauterine fetal death)

Not for distribution in the United States

Furthermore, PE has a significant impact on healthcare systems. In the United States, where the incidence of PE has increased in the past years,^{7,33} Stevens et al. estimated at \$2.18 billion in the United States the economic burden of PE during the first 12 months after delivery, comprising \$1.03 billion in excess maternal health care expenditures and \$1.15 billion in additional infant-related costs.⁸




Currently, starting preventive treatment with low-dose aspirin before 16 weeks of gestation is recommended to reduce the risk of preeclampsia in high-risk pregnant women (such as those with a history of PE).¹ Delivery, however, remains the only definitive treatment. Management prior to delivery depends on disease severity and gestational age, which is why PE is often classified based on the timing of gestation: early-onset (before 34 weeks of gestation) and late-onset (at or after 34 weeks).⁹ Identifying those at risk for early-onset PE is critical not only to inform expectant mothers, but also to personalize antenatal care and reduce the likelihood of severe complications.

The key role of angiogenic factors

The physiopathology of PE is not completely known but involves placental dysfunction and angiogenic imbalance.¹⁰⁻¹² This imbalance is characterized by abnormal concentrations of angiogenic factors in maternal circulation with a decrease of the pro-angiogenic placental growth factor (PlGF) and an increase in the anti-angiogenic soluble fms-like tyrosine kinase-1 (sFlt-1).

Placental Growth Factor (PlGF)


PlGF is a key member of the vascular endothelial growth factor (VEGF) family, predominantly expressed by placental cells. It plays a central role in placental development by promoting cytotrophoblast invasion, spiral artery remodeling, and angiogenesis, which are essential for fetal vascularization.^{10,11} PlGF concentrations increase progressively from around week 10, peak at week 30, and decline toward term.¹²



Its dynamic expression is critical for maintaining proper placental function and supporting fetal growth. It has been reported that from 13 to 16 weeks onward, PlGF is significantly lower in women who developed PE compared to those with normal and healthy pregnancies.¹² (See figures 1 and 2)

Soluble fms-like tyrosine kinase receptor-1 (sFlt-1)

sFlt-1 is a splice variant of the VEGF receptor Flt-1 and lacks membrane-binding domains, allowing it to act as a decoy receptor. It binds circulating VEGF and PlGF, preventing their interaction with cell surface receptors and thereby inhibiting angiogenesis.



Produced by the placenta, sFlt-1 regulates vascular growth during pregnancy.¹¹ In PE, elevated sFlt-1 concentrations—especially early in gestation—are associated with impaired placental vascularization. Concentrations of sFlt-1 increased earlier in gestation and to higher concentrations in individuals who later developed PE. (See figures 1 and 2)

Figure 1. PIGF and sFlt-1 production during normal pregnancy. Adapted from reference 12.

Free PIGF concentrations begin rising around week 11, peak between weeks 29–32 and supports placental development. sFlt-1 concentrations remain stable until weeks 33–36, then increase by approximately 145 pg/mL per week until delivery.





Figure 2. PIGF and sFlt-1 production in PE. Adapted from reference 12.

PIGF concentrations follow a similar pattern in all three groups but are significantly lower in pregnant women who develop PE (orange/red lines) compared to controls (green line). sFlt-1 concentrations are elevated in PE groups, with the highest concentrations observed in women with established PE (red line).



Elevated sFlt-1 and lower PIGF concentrations leading to an elevated sFlt-1/PIGF ratio is observed in cases of placental dysfunction and is associated with pregnancies complicated by PE. (Figures 3 and 4)

Figure 3. sFlt-1/PIGF ratio in women who develop PE later in pregnancy. Adapted from reference 12.

Representation of both PIGF and sFlt-1 concentrations over time reveals a significantly higher sFlt-1/PIGF ratio in women who develop PE compared to controls.





Figure 4. sFlt-1/PIGF ratio in women with established PE. Adapted from reference 12.

The sFlt-1/PIGF ratio is markedly elevated in women with established PE compared to controls, supporting its potential use in predicting disease onset



Clinical and cost-effectiveness evidence supporting the value of using the sFlt-1/PIGF ratio in clinical practice

Once the role of angiogenic factors in the pathophysiology of PE was uncovered, researchers sought to better define the clinical utility of the measurements of those biomarkers and how they could support clinicians in accurately diagnosing the disease and improving patient care. Extensive research has explored the interest of measuring the sFlt-1/PIGF ratio in three key areas: to improve diagnostic accuracy,¹³ to predict short-term development of PE,¹⁴⁻¹⁶ and to assess the prognosis of adverse maternal and fetal outcomes.¹⁷ This section highlights a selection of studies supporting the clinical and economic value of measuring the sFlt-1/PIGF ratio.¹⁸⁻²⁰ While not exhaustive, these examples illustrate the growing evidence supporting the use of the sFlt-1/PIGF ratio.

Diagnostic accuracy across gestational ages

Verlohren et al.¹³ conducted a multicenter case–control study involving 1149 patients to define gestational age–specific sFlt-1/PIGF cutoffs for the diagnosis of PE. Their results showed that using two cutoffs, one focusing on high sensitivity and the other on high specificity, improved the diagnostic precision across gestational ages:

Early gestation (20+0 – 33+6 weeks)

Low-risk cutoff ≤33	Sensitivity 95%	Specificity 94%
High-risk cutoff ≥85	Specificity 99.5%	Sensitivity 88%

Late gestation (≥34 weeks)

Low-risk cutoff ≤33	Sensitivity 89.6%	Specificity 73.1%
High-risk cutoff ≥110	Specificity 95.5%	Sensitivity 58.2%

The authors also indicate that sFlt-1/PIGF ratio values falling between these cutoffs for the early and late gestational phases are considered equivocal and should help identifying patients at intermediate risk.¹³

Short-term prediction of PE

The PROGNOSIS study¹⁴ aimed to determine whether the sFlt-1/PIGF ratio could predict the absence of PE within 1 week and the presence of PE within 4 weeks. The study was conducted as an international prospective, multicenter observational trial involving 1273 pregnant women (aged 18 years and older) at 24 to 36+6 weeks of gestation with singleton pregnancies in whom PE was suspected. The results demonstrated that an sFlt-1/PIGF ratio ≤38 helped to rule out PE within one week (NPV 99.3%), while a ratio >38 helped to predict PE within four weeks (PPV 36.7%). These findings were confirmed by PROGNOSIS Asia¹⁵ (2019; n = 764), reporting similar performance (NPV 98.6%, PPV 30.3%). The INSPIRE trial¹⁶ included 370 women (18 years old or older) with singleton pregnancy between 24+0 and 37+0 weeks of gestation and a clinical suspicion of PE. This study further validated that using the sFlt-1/PIGF ratio improved clinical precision showing 100% sensitivity and NPV for PE within one week when the ratio was ≤38. However, the use of the sFlt-1/PIGF ratio didn't contribute to a significant reduction of hospitalizations (26% vs. 32%; p = 0.192).

Prediction of adverse outcomes

Rana et al.¹⁷ measured the sFlt-1/PLGF ratio at presentation in 616 women with singleton pregnancies (28.6% before gestational week 34) with signs and symptoms of preeclampsia. They reported that the individuals who developed adverse events had a significantly higher sFlt-1/PLGF ratio, with a difference being especially pronounced before 34 weeks of gestation. They concluded that adding the sFlt-1/PIGF ratio to hypertension and proteinuria for women presenting before 34 weeks gestation could predict adverse outcomes, with an area under the curve (AUC) of 0.93 (vs. 0.84 for hypertension and proteinuria alone; P=0.001).

A sFlt-1/PIGF ratio ≥85 before 34 weeks was associated with an 86% of delivery within two weeks, compared to only 15.8% when the ratio was <85.

In addition to clinical utility and benefits, several other studies have employed economic models to evaluate the impact of implementing the sFlt-1/PIGF ratio in terms of reduction of unnecessary hospitalizations and overall healthcare expenditures. Such studies concluded that the implementation of the sFlt-1/PIGF ratio could result in cost savings per patient, estimated at £344 (UK),¹⁸ €361 (Germany)¹⁹ and €670 (Italy).²⁰ Those potential cost-savings would be driven primarily by enabling earlier intervention and avoiding unnecessary hospitalizations.

sFlt-1/PIGF Ratio in Practice: Insights from International Societies and Clinical Guidelines

The NICE guidelines,²¹ the European Society of Cardiology (ESC) guidelines,^{22,23} the combined guidelines on Hypertensive Disorders in Pregnancy²⁴ (German Society of Gynecology and Obstetrics (DGGG), Austrian Society of Gynecology and Obstetrics (OEGGG), Swiss Society of Gynecology and Obstetrics (SGGG)), the Spanish Society of Gynecology and Obstetrics (SEGO)²⁵ recommend measuring the sFlt-1/PIGF ratio to help diagnose and/or predict development of PE. In contrast, the French College of Obstetricians and Gynecologists²⁶ does not recommend routine use of the sFlt-1/PIGF ratio solely to reduce maternal or perinatal morbidity in suspected or confirmed PE, citing moderate to low evidence and limited test availability—as STAT tests available 24/7—in maternity centers across France. Similarly, the Danish Society for Obstetrics and Gynecology²⁷ excludes angiogenic assays from its guidelines due to limited clinical implementation.

Conclusion

While limited availability and access to angiogenic assays might remain a challenge in some countries, numerous studies provide evidence supporting the value of measuring the sFlt-1/PlGF ratio for managing pregnant women with signs and symptoms of preeclampsia.

About Atellica IM sFlt-1 and Atellica IM PlGF assays*

The Atellica IM PlGF and Atellica IM sFlt 1 assays are for in vitro diagnostic use in the quantitative determination of the placental growth factor and of the soluble fms-like tyrosine kinase 1 in human serum and plasma (dipotassium EDTA).²⁸⁻³¹ These two assays are used in combination on Atellica IM and CI Analyzers to determine the sFlt-1/PlGF ratio. The Atellica IM sFlt-1/PlGF ratio is intended for use:


- in the prognosis of preterm delivery and adverse outcomes in women presenting with signs and symptoms of preeclampsia.
- as an aid in the diagnosis of preeclampsia, in conjunction with signs and symptoms of the disease.
- as an aid in short-term prediction of preeclampsia (rule-out and rule-in) in pregnant women with suspicion of preeclampsia in conjunction with other diagnostic and clinical information.

The results of the sFlt-1/PlGF ratio should always be interpreted in conjunction with the patient’s medical history, clinical presentation, and other findings.

See Instructions for Use²⁸⁻³¹ and below for detailed information on each intended use of the sFlt-1/PlGF ratio determined on Atellica IM and CI Analyzers.

Prognosis of preterm delivery and adverse outcomes in women presenting with signs and symptoms of preeclampsia†

Gestational Age 20+0 – 35+6 (weeks + days)




Clinical Performance		
PE Ratio (sFlt-1/PlGF) cutoff ≥ 38	Estimate	95% CI
Sensitivity	88.42% (84/95)	81.99–94.86%
Specificity	87.65% (213/243)	83.52–91.79%
Positive Predictive Value (PPV)	73.68% (84/114)	65.60–81.77%
Negative Predictive Value (NPV)	95.09% (213/224)	92.26–97.92%

A total of 338 specimens collected from 12 sites across the United States were tested at 1 clinical testing site. Specimens were from individuals with singleton pregnancies between gestational weeks 20+0 days and 35+6 days with signs and symptoms of preeclampsia.³² The clinical performance of the PE ratio was evaluated for the prognosis of preterm delivery and adverse outcomes within 2 weeks using a cut-off value of 38. Of the 338 subjects, 95 subjects experienced preterm delivery or an adverse outcome within 2 weeks of presentation.²⁸⁻³¹

Aid in short-term prediction of preeclampsia in pregnant women with suspicion of preeclampsia in conjunction with other diagnostic and clinical information†

Gestational Age 24+0 – 36+6/7 (weeks + days)




Clinical Performance			
Parameter	N ^a	Estimate	95% CI
Positive Percent Agreement	27	100% (27/27)	(87.5%, 100%)
Negative Percent Agreement	105	96.2% (101/105)	(90.6%, 98.5%)

^aNumber of measurements.

A study enrolled 132 pregnant women aged 18 years or older with singleton pregnancies between gestational weeks 24+0 days and 36+6 days, all presenting with signs and symptoms of preeclampsia.³² Specimens were used to evaluate the diagnostic concordance between the sFlt 1/PIGF ratio obtained from the Atellica IM sFlt 1 and Atellica IM PIGF assays and a commercially available sFlt 1/PIGF ratio, using a PE ratio cut-off value of 38 for short-term prediction of preeclampsia. Results were established using the Atellica IM Analyzer.²⁸⁻³¹


Aid in the diagnosis of preeclampsia, in conjunction with signs and symptoms of the disease†

Early Gestation 20+0 – 33+6 (weeks + days)



Clinical Performance				
PE Ratio (sFlt-1/PIGF) cutoff	Sensitivity	95% CI	Specificity	95% CI
≤33 (Rule-Out)	91.6% (120/131)	85.6–95.2%	97.6% (240/246)	94.8–98.9%
≥85 (Rule-In)	81.7% (107/131)	74.2–87.4%	99.2% (244/246)	97.1–99.8%

Late Gestation 34+0 – delivery (weeks + days)



Clinical Performance				
PE Ratio (sFlt-1/PIGF) cutoff	Sensitivity	95% CI	Specificity	95% CI
≤33 (Rule-Out)	85.3% (81/95)	76.8–91.0%	84.5% (212/251)	79.5–88.4%
≥110 (Rule-In)	55.8% (53/95)	45.8–65.4%	96.8% (243/251)	93.8–98.4%

Testing was performed on 723 prospectively collected serum specimens from 18 sites across the United States at 1 site using 1 reagent lot for each assay. The 2:1 case control study included women with singleton pregnancy diagnosed with preeclampsia at the time of enrollment, and women with singleton pregnancy who did not develop preeclampsia/eclampsia/hemolysis, elevated liver enzyme levels, and low platelet count or intrauterine growth restriction during pregnancy. Preeclampsia diagnoses were confirmed by adjudication using the International Society for the Study of Hypertension in Pregnancy (ISSHP) definition for preeclampsia.¹ Separate PE ratio (sFlt 1/PlGF) cut-off values were verified for both early gestational phase (week 20+0 days and week 33+6 days) and late gestational phase (≥ 34 weeks) using the ruleout and rule-in PE ratio (sFlt 1/PlGF) cut-off values of 33 and 85 for subjects in early gestation and 33 and 110 for subjects in late gestation: Results were established using the Atellica IM Analyzer.²⁸⁻³¹

Talk to your Siemens Healthineers representative to learn more about the sFlt-1 and PlGF tests available on Atellica Analyzers. Do not miss the opportunity to enhance your lab's clinical value. By adding angiogenic markers to your laboratory's test menu, you can help improve access for obstetricians and gynecologists, supporting them in the timely and effective care of preeclampsia. Click here to explore Siemens Healthineers comprehensive test portfolio for reproductive endocrinology.

* All trademarks are the property of their respective owners. Atellica IM sFlt-1 and Atellica IM PlGF are not available for sale in the U.S.A. The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed.

† See products' Instructions for Use for detailed information on study design and clinical performance supporting each intended use.²⁸⁻³¹

References

- Magee LA, Brown MA, Hall DR, Gupte S, Hennessy A, Karumanchi SA, Kenny LC, McCarthy F, Myers J, Poon LC, Rana S, Saito S, Staff AC, Tsigas E, von Dadelszen P. The 2021 International Society for the Study of Hypertension in Pregnancy classification, diagnosis & management recommendations for international practice. *Pregnancy Hypertens* 2022; 27: 148–169
- Wu P, Haththotuwa R, Kwok CS, Babu A, Kotronias RA, Rushton C, Zaman A, Fryer AA, Kadam U, Chew-Graham CA, Mamas MA. Preeclampsia and Future Cardiovascular Health: A Systematic Review and Meta-Analysis. *Circ Cardiovasc Qual Outcomes*. 2017 Feb;10(2):e003497. doi: 10.1161/CIRCOUTCOMES.116.003497. Epub 2017 Feb 22. PMID: 28228456.
- Wu P, Kwok CS, Haththotuwa R, Kotronias RA, Babu A, Fryer AA, Myint PK, Chew-Graham CA, Mamas MA. Pre-eclampsia is associated with a twofold increase in diabetes: a systematic review and meta-analysis. *Diabetologia*. 2016 Dec;59(12):2518-2526. doi: 10.1007/s00125-016-4098-x. Epub 2016 Sep 19. PMID: 27646865; PMCID: PMC6518071.
- World Health Organization, 2025 <https://www.who.int/news-room/fact-sheets/detail/pre-eclampsia> (Consulted 11/05/2025)
- <https://www.preeclampsia.org/the-news/legislative-advocacy/preeclampsia-and-maternal-mortality-a-global-burden> (Consulted 11/05/2025)
- Olié V, Moutengou E, Grave C, Deneux-Tharaux C, Regnault N, Kretz S, Gabet A, Mounier-Vehier C, Tsatsaris V, Plu-Bureau G, Blacher J. Prevalence of hypertensive disorders during pregnancy in France (2010-2018): The Nationwide CONCEPTION Study. *J Clin Hypertens (Greenwich)*. 2021 Jul;23(7):1344-1353. doi: 10.1111/jch.14254. Epub 2021 May 27. PMID: 34042277; PMCID: PMC8678732.
- Stegers EA, von Dadelszen P, Duvekot JJ, Pijnenborg R Pre-eclampsia. *Lancet* 2010; 376: 631-44 [https://doi.org/10.1016/S0140-6736\(10\)60279-6](https://doi.org/10.1016/S0140-6736(10)60279-6)
- Stevens, Warren et al. *American Journal of Obstetrics & Gynecology*, Volume 217, Issue 3, 237 - 248. e16 Am J Obstet Gynecol. 2017 September ; 217(3): 235–236. doi:10.1016/j.ajog.2017.06.011.
- Dimitriadis, E., Rolnik, D.L., Zhou, W. et al. Pre-eclampsia. *Nat Rev Dis Primers* 9, 8 (2023). <https://doi.org/10.1038/s41572-023-00417-6>
- Chau K, Hennessy A, Makris A. Placental growth factor and pre-eclampsia. *J Hum Hypertens*. 2017 Dec;31(12):782-786. doi: 10.1038/jhh.2017.61. Epub 2017 Aug 24. PMID: 29115294; PMCID: PMC5680413.
- Phipps EA, Thadhani R, Benzing T, Karumanchi SA. Pre-eclampsia: pathogenesis, novel diagnostics and therapies. *Nat Rev Nephrol*. 2019 May;15(5):275-289. doi: 10.1038/s41581-019-0119-6. Erratum in: *Nat Rev Nephrol*. 2019 Jun;15(6):386. doi: 10.1038/s41581-019-0156-1. PMID: 30792480; PMCID: PMC6472952.
- Levine RJ, Maynard SE, Qian C, Lim KH, England LJ, Yu KF, Schisterman EF, Thadhani R, Sachs BP, Epstein FH, Sibai BM, Sukhatme VP, Karumanchi SA. Circulating angiogenic factors and the risk of preeclampsia. *N Engl J Med*. 2004 Feb 12;350(7):672-83. doi: 10.1056/NEJMoa031884. Epub 2004 Feb 5. PMID: 14764923.
- Verlohren S, Herraiz I, Lapaire O, et al. New gestational phase-specific cutoff values for the use of the soluble fms-like tyrosine kinase-1/placental growth factor ratio as a diagnostic test for preeclampsia. *Hypertension* 2014;63:346-52
- Harald Zeisler, M.D., Elisa Llurba, M.D., Ph.D., Frederic Chantraine, M.D., Ph.D., Manu Vatish, M.B., Ch.B., D.Phil., Anne Cathrine Staff, M.D., Ph.D., Maria Sennström, M.D., Ph.D., Matts Olovsson, M.D., Ph.D. and Stefan Verlohren, M.D., Ph.D. Predictive Value of the sFlt-1:PlGF Ratio in Women with Suspected Preeclampsia; *N Engl J Med* 2016;374:13-22; DOI: 10.1056/NEJMoa1414838
- Bian X, Biswas A, Huang X, Lee KJ, Li TK, Masuyama H, Ohkuchi A, Park JS, Saito S, Tan KH, Yamamoto T, Dietl A, Grill S, Verhagen-Kamerbeek WDJ, Shim JY, Hund M. Short-Term Prediction of Adverse Outcomes Using the sFlt-1 (Soluble fms-Like Tyrosine Kinase 1)/PlGF (Placental Growth Factor) Ratio in Asian Women With Suspected Preeclampsia. *Hypertension*. 2019 Jul;74(1):164-172. doi: 10.1161/HYPERTENSIONAHA.119.12760. Epub 2019 Jun 12. PMID: 31188674; PMCID: PMC6587370.
- Cerdeira AS, O'Sullivan J, Ohuma EO, Harrington D, Szafranski P, Black R, Mackillop L, Impey L, Greenwood C, James T, Smith I, Papageorgiou AT, Knight M, Vatish M. Randomized Interventional Study on Prediction of Preeclampsia/Eclampsia in Women With Suspected Preeclampsia: INSPIRE. *Hypertension*. 2019 Oct;74(4):983-990. doi: 10.1161/HYPERTENSIONAHA.119.12739. Epub 2019 Aug 12. PMID: 31401877; PMCID: PMC6756298.
- Rana S, Powe CE, Salahuddin S, Verlohren S, Perschel FH, Levine RJ, Lim KH, Wenger JB, Thadhani R, Karumanchi SA. Angiogenic factors and the risk of adverse outcomes in women with suspected preeclampsia. *Circulation*. 2012 Feb 21;125(7):911-9. doi: 10.1161/CIRCULATIONAHA.111.054361. Epub 2012 Jan 18. PMID: 22261192; PMCID: PMC3319742.
- Vatish M, Strunz-McKendry T, Hund M, Allegranza D, Wolf C, Smare C. sFlt-1/PlGF ratio test for pre-eclampsia: an economic assessment for the UK. *Ultrasound Obstet Gynecol*. 2016 Dec;48(6):765-771. doi: 10.1002/uog.15997. Epub 2016 Nov 8. PMID: 27300726; PMCID: PMC5215695.
- Schlembach D, Hund M, Schroer A, Wolf C. Economic assessment of the use of the sFlt-1/PlGF ratio test to predict preeclampsia in Germany. *BMC Health Serv Res*. 2018 Aug 6;18(1):603. doi: 10.1186/s12913-018-3406-1. PMID: 30081904; PMCID: PMC6080558.
- Frusca, T., Gervasi, M. T., Paolini, D., Dionisi, M., Ferre, F., & Cetin, I. (2017). Budget impact analysis of sFlt-1/PlGF ratio as prediction test in Italian women with suspected preeclampsia. *The Journal of Maternal-Fetal & Neonatal Medicine*, 30(18), 2166–2173. <https://doi.org/10.1080/14767058.2016.1242122>

References continued

21. PLGF-based testing to help diagnose suspected preterm pre-eclampsia. Diagnostics guidance [DG49] National Institute for Health and Care Excellence (NICE) Published: 27 July 2022 www.nice.org.uk/guidance/dg49
22. Regitz-Zagrosek V, Roos-Hesselink JW, Bauersachs J, Blomström-Lundqvist C, Cifková R, De Bonis M, et al. 2018 ESC Guidelines for the management of cardiovascular diseases during pregnancy. *Eur Heart J* 2018;39:3165–241. <https://doi.org/10.1093/eurheartj/ehy340>
23. Julie De Backer, Kristina H Haugaa, Nina Eide Hasselberg, Michèle de Hosson, Margarita Brida, Silvia Castelletti, Matthew Cauldwell, Elisabetta Cerbai, Lia Crotti, Natasja M S de Groot, Mette-Elise Estensen, Eva S Goossens, Bernhard Haring, Donata Kurpas, Carmel M McEniery, Sanne A E Peters, Amina Rakisheva, Antonia Sambola, Oliver Schlager, Florian S Schoenhoff, Tommaso Simoncini, Françoise Steinbach, Isabella Sudano, Lorna Swan, Anne Marie Valente, ESC Scientific Document Group, 2025 ESC Guidelines for the management of cardiovascular disease and pregnancy: Developed by the task force on the management of cardiovascular disease and pregnancy of the European Society of Cardiology (ESC). Endorsed by the European Society of Gynecology (ESG), *European Heart Journal*, 2025;, ehaf193, <https://doi.org/10.1093/eurheartj/ehaf193>
24. Pecks U, Baumann M, Binder J, Contini C, Dathan-Stumpf A, Dechend R, Enna-Kirchmair B, Fischer T, Girard T, Greve S, Groten T, Hartung A, Kehl S, Koch M, Köbke A, Kranke P, Lapaire O, Mader S, Rump LC, Sperling A, Stepan H, Stracke S, Verlohren S, von Versen-Höyneck F, Winkler K, Zemlin M, Schlembach D. Hypertensive Disorders in Pregnancy (HDP): Diagnostics and Therapy. Guideline of the DGGG, OEGGG and SGGG (S2k-Level, AWMF Registry No. 015/018, June 2024). *Geburtshilfe Frauenheilkd.* 2025 Jun 2;85(8):810-850. doi: 10.1055/a-2522-2347. PMID: 40771833; PMCID: PMC12324857.
25. Trastornos hipertensivos en la gestación. Sociedad Española de Ginecología y Obstetricia. *Trastornos hipertensivos en la gestación. Prog Obstet Ginecol* 2020;63:244-272
26. La pré-éclampsie : recommandations pour la pratique clinique du Collège national des gynécologues obstétriciens français. Loïc Sentilhes, Thomas Schmitz, Chloé Arthuis, Tiphaine Barjat, Paul Berveiller, Céline Camilleri, Alizée Froeliger, Charles Garabedian, Paul Guerby, Diane Korb, Edouard Lecarpentier, Aurélien Mattuizzi, Jeanne Sibuide, Marie-Victoire Sénat, Vassilis Tsatsaris. *Gynécologie Obstétrique Fertilité & Sénologie*. Volume 52, Issue 1, January 2024, Pages 3-44.
27. Danish Society for Obstetrics and Gynecology: Danish 2018 national clinical guideline for hypertensive disorders in pregnancy and preeclampsia. *Obstetrik — DSOG* <https://static1.squarespace.com/static/5467abcce4b056d72594db79/t/5bac84e7652dea0a1b5fb489/1538032877105>
28. Placental Growth Factor (PIGF). Instructions for Use 11205121_EN Rev. 02, 2025-08 (Atellica IM Analyzer)
29. Placental Growth Factor (PIGF). Instructions for Use 11205028_EN Rev.02, 2025-08 (Atellica CI Analyzer)
30. Soluble fms-like Tyrosine Kinase-1 (sFlt 1) Instructions for Use 11205117_EN Rev. 02, 2025-08 (Atellica IM Analyzer)
31. Soluble fms-like Tyrosine Kinase-1 (sFlt1) Instructions for Use 111205055_EN Rev.02, 2025-08 (Atellica CI Analyzer)
32. Hypertension in pregnancy. Report of the American College of Obstetricians and Gynecologists' Task doi:10.1097/01.AOG.0000437382.03963.88 *Force on Hypertension in Pregnancy. Obstet Gynecol.* 2013;122(5):1122-1131.
33. Ayyash MK, McLaren R Jr, Shaman M, Al-Kouatly HB. Trends in preeclampsia risk factors in the US from 2010 to 2021. *JAMA*. Published online June 10, 2024. doi:10.1001/jama.2024.8931

All trademarks are the property of their respective owners. Not available for sale in the U.S.A. The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed.

Not for distribution in the United States

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Core Lab Solutions
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000