

## Customer Case Study



# Analytical performance evaluation of a new integrated clinical chemistry analyzer

Like all clinical laboratories, MVZ Labor Dr. Limbach is challenged with reporting accurate results under operational constraints

## MVZ Labor Dr. Limbach & Kollegen GbR



Heidelberg, Germany



This laboratory is part of Limbach group SE, and according to its website, the largest network of laboratories in Germany. Limbach group SE is one of the **largest networks of laboratories in Europe**.



### Challenge

Deliver precise, clinically relevant, and timely results while managing staff shortage and budgets limitations is the challenge that most of clinical laboratories are facing.



### Risks

- Patient misdiagnosis
- Delayed treatment
- Regulatory / Legal
- Reputation damage
- Increased costs



### Solution: Atellica CI Analyzer

- Standardization for increased operational efficiency
- Intelligent sample aspiration
- Maintained throughput for rack-based processing
- Fast turnaround time



**Up to 1000 TPH\***

(600 PMs, 400 ISEs\*)  
with 70 reagent positions



**Up to 120 TPH\***

with 20 primary and  
20 ancillary reagent positions



1

Instill confidence in quality and accuracy of results

2

Increase testing efficiency

3

Address staff shortages and budget constraints effectively



**Study objectives:** validation of analytical performances of Atellica CI key assays in real world condition

### 17 commonly used assays

- ✓ **13 chemistry tests (CC):** Albumin, Alkaline Phosphatase, ALAT, AST, Amylase, Calcium, Cholesterol, Glucose, HDL Cholesterol, LDL Cholesterol, Triglycerides, Uric Acid, Urea Nitrogen
- ✓ **4 immunoassay tests (IA):** TSH, FT4, PSA, HCG

### Assessment of precision

- ✓ CLSI - EP15-A3
- ✓ QC materials or patient pools
- ✓ Repeatability and within-laboratory

### Method comparisons vs. Roche COBAS 6000

- ✓ CLSI EP09-A3
- ✓ Remnant anonymized serum patient samples
- ✓ Correlation coefficients of  $\geq 0.95$  were considered acceptable



## Outcomes



### Precision

CC	Repeatability CVs	from 0 % to 1.9 %		from 0.5 % to 5.1 %	
	Within-laboratory CVs	from 0.5 % to 3.7 %		from 0.9 % to 5.2 %	
IA	Repeatability CVs	from 1.1 % to 2.5 %		from 0.6 % to 2.3 %	
	Within-laboratory CVs	from 1.6 % to 3.8 %		from 1.6 % to 5.4 %	
In total	Within-laboratory CVs*	< 3%	< 4%	< 3%	< 4%
	% of CC and IA tests	76%	100%	65%	85%

- ✓ Results of the Atellica CI Analyzer precision studies met manufacturer's specifications.
- ✓ Evaluations using pooled patient samples and QC materials showed comparable results.

### Method comparisons

Passing Bablok	Correlation coef. (r) values from 0.98 to 1 (CC: 0.989 to 1.0 – IA: 0.98 to 0.999).
	Slopes: CC tests from 0.946 to 1.143. - IA tests from 0.798 to 1.276.
Bland Altman	CC tests : mean bias from -0.4% (smallest) to 12,4% (largest).
	IA tests: mean bias from -3.9% (smallest for TSH) to 26.2% (largest for HCG).

- ✓ Overall, results of method comparisons were acceptable.

1

Analytical performances meet specifications

2

Atellica CI Analyzer is suitable for routine clinical use

## References

Ruffing, S. Mickeler, M. Kraft, P. Findeisen, Analytical performance evaluation of a new integrated clinical chemistry and immunoassay analyzer, *Practical Laboratory Medicine*,  
<https://doi.org/10.1016/j.plabm.2024.e00427>

The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Because there is no "typical" hospital or laboratory and many variables exist (e.g., hospital size, samples mix, case mix, level of IT and/or automation adoption) there can be no guarantee that other customers will achieve the same results).

Product availability varies by country and is subject to varying regulatory requirements. Their future availability cannot be guaranteed. Please contact your local representative for availability.

This study was funded by Siemens Healthineers.

---

### Siemens Healthineers Headquarters

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 180  
[siemens-healthineers.com](https://www.siemens-healthineers.com)

### Published by

Siemens Healthcare Diagnostics Inc.  
Core Lab Solutions  
511 Benedict Avenue  
Tarrytown, NY 10591-5005  
USA  
Phone: +1 914-631-8000