



Study Protocol

Control *syngo* DynaCT directly following endovascular aneurysm sealing procedure

Interventional Surgery

Immediate post-procedural DynaCT improves technical success of EVAR procedures.

Courtesy of

Prof. Dietmar Böckler, MD,
Department for Vascular and
Endovascular Surgery,
University Hospital Heidelberg,
Germany

Supported by

syngo DynaCT DSA

System & Software

Artis zeego VC14
syngo X Workplace VB21

Case Description

Patient history

73-year-old male

Diagnosis

Asymptomatic patient with progressive, penetrating atherosclerotic ulcer (PAU) in the infrarenal aorta (max. diameter 40 mm).

Treatment

Endovascular aneurysm sealing using Endologix Nellix™.

Tips & Tricks

syngo DynaCT in portrait orientation (24 cm vertical coverage) has a large anatomical coverage that allows visualization of the renal arteries, partially suprarenal aorta, infrarenal aorta and its branches, and the iliac arteries during EVAR procedures. Immediate post-procedural control DynaCT allows the physician to evaluate his treatment so that the physician can see the complications early and change the treatment

strategy. *syngo* DynaCT removes the need for further postprocedural CT examinations.

General comments

Control *syngo* DynaCT directly following aneurysm sealing enables assesment of correct deployment of stent graft, patency of renal, visceral, and iliac arteries. Complications such as endoleaks and stent stenosis also excluded thanks to *syngo* DynaCT.

Control *syngo* DynaCT directly following endovascular aneurysm sealing procedure

Acquisition protocol	5s DCT Body CARE Portrait (prototype calibration)
----------------------	---------------------------------------------------

Injection protocol

Catheter position	Above renal arteries
Contrast medium (CM)	300 mg iodine/mL
Dilution	66 %
Injection volume	70 mL
Injection rate	10 mL/s
Duration of injection	7 s
X-ray delay	2 s
Power injector used	Yes

Reconstructions

Primary

Name	DynaCT Body Nat Fill HU
VOI size	Full
Slice matrix	512 × 512
Kernel type	HU
Image characteristics	Normal
Reconstruction mode	Nat fill
Viewing preset	Auto

Clinical Images



Figure 1: Coronal MIP 50 mm shows optimal deployment of stent and patency of renal arteries




Figure 2: Coronal VRT shows correct deployment of stent and patency of vessels

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.

All rights reserved.