

# Quantitative Determination of Apixaban: Method Comparison Studies Between Different Assays and Platforms

C. KIK<sup>1</sup>, M. WILKENS<sup>1</sup>, S. PILGRIM<sup>1</sup>, M. BOEHM-WEIGERT<sup>1</sup>, R.H. CHRISTENSON<sup>2</sup>, R. FRANCIS<sup>3</sup>, A. MOHAMMAD<sup>4</sup>, A. SIEGEMUND<sup>5</sup>, E. WALTERS<sup>6</sup>

<sup>1</sup>Siemens Healthcare Diagnostics Products GmbH, Marburg, Germany; <sup>2</sup>University of Maryland School of Medicine, Baltimore, MD, USA; <sup>3</sup>Columbia University Medical Center, New York, NY, USA;

<sup>4</sup>Baylor Scott & White Health Center, Temple, TX, USA ; <sup>5</sup>MVZ Limbach Magdeburg GmbH, Magdeburg, Germany; <sup>6</sup>Loma Linda University Medical Center, Loma Linda, CA, USA

## Introduction and Aim

Apixaban is a direct oral anticoagulant (DOAC) that specifically inhibits Factor Xa (FXa). Although laboratory monitoring of DOACs is not required per se, there is a growing recognition of the need to test in specific clinical situations and patient populations. The quantitative determination of DOACs is possible with chromogenic anti-Xa assays originally developed for quantitative determination of heparin enoxaparin concentration using drug-specific calibrators.

The INNOVANCE Anti-Xa assay from Siemens Healthineers is an FXa-based chromogenic assay used with drug-specific standards and controls for the quantitative determination of heparin activity (unfractionated heparin [UFH] and low molecular weight heparin [LMWH]), as well as rivaroxaban\*, edoxaban\*, and apixaban concentration in human citrated plasma. Respective applications are available on various hemostasis systems.

This study evaluates the comparability of the INNOVANCE Anti-Xa assay to measure apixaban by employing apixaban-specific standards and controls with liquid chromatography–tandem mass spectrometry (LC-MS/MS) and other commercial methods. The INNOVANCE Anti-Xa assay comparability on different systems from Siemens Healthineers was also investigated.

## Methods

### Assays

Samples were obtained from patients receiving apixaban therapy or suspected to be under therapy with apixaban and tested with the following assays/devices:

- INNOVANCE Anti-Xa assay from Siemens Healthineers for the quantitative determination of apixaban (performed on the CS-2500, CS-5100, CN-3000\*, and CN-6000\* Systems)
- Liquid chromatography–tandem mass spectrometry (LC-MS/MS; Swiss BioQuant, Reinach, Switzerland)
- STA-Liquid Anti-Xa assay (Stago; Asnières-sur-Seine, France) for the quantitative determination of apixaban (performed on the STA Compact Max System, Stago)
- HemosIL Liquid Anti-Xa assay (Werfen; Barcelona, Spain) for the quantitative determination of apixaban (performed on the ACL TOP Family 50 Series analyzer, Werfen).

### Measurements

#### Method comparison of INNOVANCE Anti-Xa assay versus LC-MS/MS

A method comparison study was performed using the INNOVANCE Anti-Xa assay on the CS-2500 System versus the LC-MS/MS method for measuring apixaban concentrations of 101 frozen samples (n = 47 native samples, n = 52 spiked samples, and n = 2 controls). The sample measurements with the INNOVANCE Anti-Xa assay were done at Siemens Healthineers, and with the LC-MS/MS method at Swiss BioQuant.

#### Method comparison of INNOVANCE Anti-Xa assay versus STA Liquid Anti-Xa assay

Using a different set of 107 frozen samples (n = 102 native samples from the intended use population and n = 5 spiked samples), a method comparison study was performed comparing the INNOVANCE Anti-Xa assay on the CS-2500 System to the STA-Liquid Anti-Xa assay on the STA Compact Max system at Siemens Healthineers.

#### Method comparison of INNOVANCE Anti-Xa assay versus HemosIL Liquid Anti-Xa assay

Method comparisons using samples from patients under therapy or suspected to be under therapy with apixaban (n = 306) were enrolled at three clinical sites. N = 18 samples were diluted or spiked to cover the measuring interval of the INNOVANCE Anti-Xa assay (20–700 ng/mL). Samples were measured fresh or after frozen storage at the enrollment sites using the INNOVANCE Anti-Xa assay on the CS-2500 System and the HemosIL Liquid Anti-Xa assay on the ACL TOP Family 50 Series analyzer.

#### Method comparison of INNOVANCE Anti-Xa assay on different Siemens Healthineers systems

The samples described above (n = 306) were also measured at Siemens Healthineers with the INNOVANCE Anti-Xa assay on CS-5100, CN-3000\*, and CN-6000\* systems. All samples were measured following frozen storage.

### Statistical analysis

The method comparison studies were performed and evaluated by Passing-Bablok regression using webMC (method comparison) in accordance with CLSI EP09-A3, “Measurement Procedure Comparison and Bias Estimation.” Correlation coefficient (r), intercept, slope and bias at apixaban concentrations of 30 ng/mL and 50 ng/mL were calculated for the entire measuring interval (20–700 ng/mL) as well as for the lower concentration range (20–100 ng/mL). We also performed the regression in the lower concentration range to better evaluate the data in that range and achieve a more precise and correct estimate for the predicted bias at 30 ng/mL and 50 ng/mL.

## Results

The summarized statistical results from the regression analysis are listed in Table 1 and are graphically displayed in Figures 1–4.

Table 1. Method comparison results.

Method 2 (y-axis)		Method 1 (x-axis)		Evaluated interval (ng/mL)	n	Correlation Coefficient (r)	Slope	Intercept (ng/mL)	Deviation at:	
Assay	System	Assay	System						30 ng/mL	50 ng/mL
INNOVANCE Anti-Xa	CS-2500	n/a	LC-MS/MS	20–350 (Figure 1)	101	0.987	1.034	0.759	1.8	2.5
				20–100	65	0.969	0.982	3.651	3.1	2.7
INNOVANCE Anti-Xa	CS-2500	STA Liquid-Anti-Xa	STA Compact Max	20–400 (Figure 2)	107	0.975	1.081	-2.329	0.1	1.7
				20–100	41	0.943	1.029	-0.178	0.7	1.3
INNOVANCE Anti-Xa	CS-2500	HemosIL Liquid Anti-Xa	ACL TOP Family 50 Series	20–600 (Figure 3)	301	0.989	1.026	-8.590	-7.8	-7.3
				20–100	90	0.935	0.801	8.889	2.9	-1.1
INNOVANCE Anti-Xa	CS-5100	INNOVANCE Anti-Xa	CS-2500	20–600 (Figure 4)	305	0.999	1.007	0.169	0.4	0.5
				20–100	103	0.998	1.018	-0.554	0.0	0.3
INNOVANCE Anti-Xa	CN-3000*	INNOVANCE Anti-Xa	CS-2500	20–600 (Figure 4)	306	0.998	1.023	1.023	1.7	2.2
				20–100	104	0.997	1.053	-0.211	1.4	2.4
INNOVANCE Anti-Xa	CN-6000*	INNOVANCE Anti-Xa	CS-2500	20–600 (Figure 4)	305	0.998	1.022	-1.532	-0.9	-0.5
				20–100	103	0.998	1.018	-1.036	-0.5	-0.1




Figure 1. Passing-Bablok regression analysis for the INNOVANCE Anti-Xa assay on the CS-2500 System versus results obtained using the liquid chromatography–tandem mass spectrometry (LC-MS/MS).




Figure 2. Passing-Bablok regression analysis for the INNOVANCE Anti-Xa assay on the CS-2500 System versus results obtained using the STA-Liquid Anti-Xa assay on the STA Compact Max System.




Figure 3. Passing-Bablok regression analysis for the INNOVANCE Anti-Xa assay on the CS-2500 System versus the HemosIL Liquid Anti-Xa assay on the ACL TOP Family 50 Series analyzer.




Figure 4. Passing-Bablok regression analysis for the INNOVANCE Anti-Xa assay performed on: A) the CS-2500 System versus the CS-5100 System, B) the CS-2500 System versus the CN-3000 System\*, and C) the CS-2500 System versus the CN-6000 System\*.

## Conclusion

The method comparison studies demonstrate that the INNOVANCE Anti-Xa assay, in combination with the INNOVANCE Apixaban Standards and INNOVANCE Apixaban Controls on the CS-2500 System, compares very well with the LC-MS/MS method and with assays from Stago and Werfen, with correlation coefficients  $\geq 0.975$ , intercepts  $\leq 8.590$  ng/mL, and slopes between 1.026 and 1.081.

The INNOVANCE Anti-Xa assay for Apixaban is available on various hemostasis systems from Siemens Healthineers. The method comparison study of the INNOVANCE Anti-Xa assay on the CS-2500 System versus the INNOVANCE Anti-Xa assay on the CS-5100 System, the CN-3000 System\*, and the CN-6000 System\* from Siemens Healthineers demonstrates an excellent agreement with correlation coefficients  $\geq 0.998$ , intercepts  $\leq 1.532$  ng/mL, and slopes between 1.007 and 1.023.

In summary, the data demonstrate the accuracy and the comparability among systems for the INNOVANCE Anti-Xa assay in combination with INNOVANCE Apixaban Standards and INNOVANCE Apixaban Controls for the quantitative determination of apixaban in human-citrated plasma.

All trademarks are the property of their respective owners.

(\*) Not available for sale in the U.S. The availability of products/features/applications varies by countries. Their future availability cannot be guaranteed.

LC-MS/MS refers to liquid chromatography–tandem mass spectrometry.

CS-2500 and CS-5100 Systems refer to Automated Blood Coagulation Analyzer CS-2500 and Automated Blood Coagulation Analyzer CS-5100, respectively.

CN-3000 and CN-6000 Systems refer to Automated Blood Coagulation Analyzer CN-3000 and Automated Blood Coagulation Analyzer CN-6000, respectively.

Abstract Presentation/Board Number: PB0189

Contact information

Corinna Kik: corinna.kik@siemens-healthineers.com

Dr. Matthias Wilkens: matthias.wilkens@siemens-healthineers.com