MAGNETOM Flash

The Magazine of MR 翻訳版

SIEMENS

Clinical Cardiovascular MRI

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Contact

Igor Klem, M.D.
 Duke Cardiovascular Magnetic Resonance Center
 Duke University Medical Center 3934
 Durham, NC 27710
 USA
 igorklem@duke.edu

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- 1 Kim RJ, Shah DJ, Judd RM. How we perform delayed enhancement imaging. JCMR 2003;5:505–514.
- 2 Kim RJ, Fieno DS, Parrish TB, et al. Relationship of MRI delayed contrast enhancement to irreversible injury, infarct age, and contractile function. Circulation 1999;100:1992–2002.
- 3 Kim RJ, Wu E, Rafael A, et al. The use of contrast-enhanced magnetic resonance imaging to identify reversible myocardial dysfunction. N Engl J Med 2000;343:1445–53.
- 4 Wagner A, Mahrholdt H, Holly TA, et al. Contrast-enhanced MRI and routine single photon emission computed tomography (SPECT) perfusion imaging for detection of subendocardial myocardial infarcts: an imaging study. Lancet 2003;361:374–79.
- 5 Choi KM, Kim RJ, Gubernikoff G, et al. Transmural extent of acute myocardial infarction predicts long-term improvement in contractile function. Circulation 2001;104:1101–7.
- 6 Klem I, Heitner JF, Shah DJ, et al. Improved detection of coronary artery disease by stress perfusion cardiovascular magnetic resonance with the use of delayed enhancement infarction imaging. J Am Coll Cardiol 2006;47:1630–8.

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- Likoff MJ, Chandler SL, Kay HR. Clinical determinants of mortality in chronic congestive heart failure secondary to idiopathic dilated or to ischemic cardiomyopathy. *Am J Cardiol* 1987;59:634-8.
- Bart BA, Shaw LK, McCants Jr CB, et al. Clinical determinants of mortality in patients with angiographically diagnosed ischemic or nonischemic cardiomyopathy. *J Am Coll Cardiol* 1997;30: 1002-8.
- Allman KC, Shaw LJ, Hachamovitch R, et al. Myocardial viability testing and impact of revascularization on prognosis in patients with coronary artery disease and left ventricular dysfunction: a meta-analysis. *J Am Coll Cardiol* 2002;39:1151-8.
- Casolo G, Minneci S, Manta R, et al. Identification of the ischaemic etiology of heart failure by cardiovascular magnetic resonance imaging: diagnostic accuracy of late gadolinium enhancement. *Am Heart J* 2006;151:101-108.
- McCrohon JA, Moon JC, Prasad SK, et al. Differentiation of heart failure related to dilated cardiomyopathy and coronary artery disease using gadolinium-enhanced cardiovascular magnetic resonance. *Circulation* 2003;108:54-9.
- Bello D, Shah DJ, Farah GM, et al. Gadolinium cardiovascular magnetic resonance predicts reversible myocardial dysfunction and remodeling in patients with heart failure undergoing beta-blocker therapy. *Circulation* 2003;108:1945-53.
- Soriano CJ, Ridocci F, Estornell J, et al. Late gadolinium-enhanced cardiovascular magnetic resonance identifies patients with standardized definition of ischemic cardiomyopathy: a single centre experience. *Int J Cardiol* 2007;116:167-173.
- Bulkley BH, Hutchins GM, Bailey I, et al. Thallium 201 imaging and gated cardiac blood pool scans in patients with ischemic and idiopathic congestive cardiomyopathy: a clinical and pathological study. *Circulation* 1977;55:753-760.
- Wagner A, Mahrholdt H, Holly TA, Elliott M.D., Regenfus M, Parker M, Klocke FJ, Bonow RO, Kim RJ, Judd RM. Contrast-enhanced MRI and routine single photon emission computed tomography (SPECT) perfusion imaging for detection of sub-endocardial myocardial infarcts: an imaging study. *Lancet* 2003; 361: 374-9.
- Barkhausen J, Hunold P, Jochims M, et al. Imaging of myocardial perfusion with magnetic resonance. *J Magn Reson Imaging* 2004;19:750-7.
- Assomull RG, Prasad SK, Lyne J, et al. Cardiovascular magnetic resonance, fibrosis, and prognosis in dilated cardiomyopathy. *J Am Coll Cardiol* 2006;48:1977-1985.
- Keenan NG, Pennell DJ. CMR of ventricular function. *Echocardiography* 2007;24:185-193.
- Bottini PB, Carr AA, Prisant M, et al: Magnetic resonance imaging compared to echocardiography to assess left ventricular mass in the hypertensive patient. *Am J Hypertens* 1995;8:221-228.
- Olivieri NF, Nathan DG, MacMillan JH, Wayne AS, Liu PP, McGee A, Martin M, Koren G, Cohen AR. Survival in medically treated patients with homozygous beta-thalassemia. *N Engl J Med*. 1994;331: 574 -578.
- Borgna-Pignatti C, Rugolotto S, De Stefano P, Zhao H, Cappellini M.D., Del Vecchio GC, Romeo MA, Forni GL, Gamberini MR, Ghilardi R, Piga A, Cnaan A. Survival and complications in patients with thalassemia major treated with transfusion and deferoxamine. *Haematologica*. 2004;89:1187-1193.

- 16 Modell B, Khan M, Darlison M. Survival in beta thalassaemia major in the UK: data from the UK Thalassaemia Register. *Lancet*. 2000; 355:2051–2052.
- 17 Wacker Ph.D., Balmer-Ruedin D, Oberhansli I, Wyss M. Regression of cardiac insufficiency after ambulatory intravenous deferoxamine in thalassaemia major. *Chest* 1993; 103: 1276–8.
- 18 Anderson LJ, Holden S, Davis B, et al. Cardiovascular T2-star (*) magnetic resonance for the early diagnosis of myocardial iron overload. *Eur Heart J* 2001;21:2171–2179.
- 19 Tanner MA, Galanello R, Dessi C, et al. Myocardial iron loading in patients with thalassaemia major on deferoxamine chelation. *J Cardiovasc Magn Reson* 2006;8:543–547.
- 20 Anderson LJ, Westwood MA, Holden S, et al. Myocardial iron clearance during reversal of siderotic cardiomyopathy with intravenous desferrioxamine: a prospective study using T2* cardiovascular magnetic resonance. *Br J Haematol* 2004;127:348–355.
- 21 Anderson LJ, Wonke B, Prescott E, et al. Comparison of effects of oral deferiprone and subcutaneous desferrioxamine on myocardial iron concentrations and ventricular function in beta-thalassaemia. *Lancet* 2002;360:516–520.
- 22 Pennell DJ, Berdoukas V, Karagiorga M, et al. Randomised controlled trial of deferiprone or deferoxamine in beta-thalassaemia major in patients with asymptomatic myocardial siderosis. *Blood* 2006;107:3738–3744.
- 23 Tanner MA, Galanello R, Dessi C, et al. A randomized, placebo-controlled, double-blind trial of the effect of combined therapy with deferoxamine and deferiprone on myocardial iron in thalassaemia major using cardiovascular magnetic resonance. *Circulation* 2007;115:1876–1884.
- 24 Modell B, Khan M, Darlison M, Westwood MA, Pennell DJ, David Ingram D. Falling mortality and changing causes of death in thalassaemia major: data from the UK thalassaemia register and possible contribution from T2* cardiovascular magnetic resonance. *J Cardiovasc Magn Reson* In review.
- 25 Telfer P, Coen PG, Christou S, Hadjigavriel M, Kolnakou A, Pangalou E, Pavlides N, Psiloines M, Simamonian K, Skordos G, Sitarou M, Angastiniotis M. Survival of medically treated thalassaemia patients in Cyprus. Trends and risk factors over the period 1980–2004. *Haematologica* 2006; 91: 1187–92.
- 26 Falk RH, Skinner M. The systemic amyloidoses: an overview. *Adv Intern Med*. 2000;45:107–137.
- 27 Maceira AM, Joshi J, Prasad SK, et al. Cardiovascular magnetic resonance in cardiac amyloidosis. *Circulation* 2005;111: 186–193.
- 28 Iwai K, Tachibana T, Takemura T, et al. Pathological studies on sarcoidosis autopsy. I. Epidemiological features of 320 cases in Japan. *Acta Pathol Jpn*. Jul–Aug 1993;43(7–8):372–376.
- 29 Feldman AM, McNamara D. Myocarditis. *N Engl J Med*. Nov 9 2000;343(19):1388–1398.
- 30 Kawai C. From myocarditis to cardiomyopathy: mechanisms of inflammation and cell death: learning from the past for the future. *Circulation*. Mar 2 1999;99(8):1091–1100.
- 31 Magnani JW, Dec GW. Myocarditis: current trends in diagnosis and treatment. *Circulation*. Feb 14 2006;113(6):876–890.
- 32 Friedrich MG, Stroh O, Schulz-Menger J, et al. Contrast media-enhanced magnetic resonance imaging visualizes myocardial changes in the course of viral myocarditis. *Circulation*. May 12 1998;97(18):1802–1809.
- 33 Mahrholdt H, Goedecke C, Wagner A, et al. Cardiovascular magnetic resonance assessment of human myocarditis: a comparison to histology and molecular pathology. *Circulation*. Mar 16 2004;109(10):1250–1258.
- 34 Mahrholdt H, Wagner A, Deluigi CC, et al. Presentation, patterns of myocardial damage, and clinical course of viral myocarditis. *Circulation*. Oct 10 2006;114(15):1581–1590.
- 35 Abdel-Aty H, Boye P, Zagrosek A, et al. Diagnostic performance of cardiovascular magnetic resonance in patients with suspected acute myocarditis: comparison of different approaches. *J Am Coll Cardiol*. Jun 7 2005;45(11): 1815–1822.
- 36 Schulz-Menger J, Wassmuth R, Abdel-Aty H, et al. Patterns of myocardial inflammation and scarring in sarcoidosis as assessed by cardiovascular magnetic resonance. *Heart*. Mar 2006;92(3):399–400.
- 37 Smedema JP, Snoep G, van Kroonenburgh MP, et al. Cardiac involvement in patients with pulmonary sarcoidosis assessed at two university medical centers in the Netherlands. *Chest*. Jul 2005;128(1):30–35.
- 38 Vignaux O, Dhote R, Duboc D, et al. Detection of myocardial involvement in patients with sarcoidosis applying T2-weighted, contrast-enhanced, and cine magnetic resonance imaging: initial results of a prospective study. *J Comput Assist Tomogr*. Sep-Oct 2002;26(5):762–767.
- 39 Singh JA, Woodard PK, Davila-Roman VG, et al. Cardiac magnetic resonance imaging abnormalities in systemic lupus erythematosus: a preliminary report. *Lupus*. 2005;14(2):137–144.
- 40 Rickers C, Wilke NM, Jerosch-Herold M, et al. Utility of cardiac magnetic resonance imaging in the diagnosis of hypertrophic cardiomyopathy. *Circulation*. Aug 9 2005;112(6):855–861.
- 41 Schulz-Menger J, Abdel-Aty H, Busjahn A, et al. Left ventricular outflow tract planimetry by cardiovascular magnetic resonance differentiates obstructive from non-obstructive hypertrophic cardiomyopathy. *J Cardiovasc Magn Reson*. 2006;8(5):741–746.
- 42 Schulz-Menger J, Gross M, Messroghli D, et al. Cardiovascular magnetic resonance of acute myocardial infarction at a very early stage. *J Am Coll Cardiol*. Aug 6 2003;42(3):513–518.
- 43 van Dockum WG, ten Cate FJ, ten Berg JM, et al. Myocardial infarction after percutaneous transluminal septal myocardial ablation in hypertrophic obstructive cardiomyopathy: evaluation by contrast-enhanced magnetic resonance imaging. *J Am Coll Cardiol*. Jan 7 2004;43(1):27–34.
- 44 van Dockum WG, Beek AM, ten Cate FJ, et al. Early onset and progression of left ventricular remodeling after alcohol septal ablation in hypertrophic obstructive cardiomyopathy. *Circulation*. May 17 2005;111(19):2503–2508.
- 45 Moon JC, McKenna WJ, McCrohon JA, et al. Toward clinical risk assessment in hypertrophic cardiomyopathy with gadolinium cardiovascular magnetic resonance. *J Am Coll Cardiol*. May 7 2003;41(9):1561–1567.
- 46 Tandri H, Friedrich MG, Calkins H, et al. MRI of arrhythmogenic right ventricular cardiomyopathy/dysplasia. *J Cardiovasc Magn Reson*. 2004;6(2):557–563.
- 47 Bluemke DA, Krupinski EA, Ovitt T, et al. MR Imaging of arrhythmogenic right ventricular cardiomyopathy: morphologic findings and interobserver reliability. *Cardiology*. 2003;99(3):153–162.
- 48 Tandri H, Saranathan M, Rodriguez ER, et al. Noninvasive detection of myocardial fibrosis in arrhythmogenic right ventricular cardiomyopathy using delayed-enhancement magnetic resonance imaging. *J Am Coll Cardiol*. Jan 4 2005;45(1):98–103.

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- Fleischmann, K.E., et al., Exercise echocardiography or exercise SPECT imaging? A meta-analysis of diagnostic test performance. *J Nucl Cardiol*, 2002. 9(1): p. 133–4.
- Klocke, F.J., et al., ACC/AHA/ASNC guidelines for the clinical use of cardiac radionuclide imaging-executive summary: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (ACC/AHA/ASNC Committee to Revise the 1995 Guidelines for the Clinical Use of Cardiac Radionuclide Imaging). *Circulation*, 2003. 108(11): p. 1404–18.
- Klem, I., et al., Improved detection of coronary artery disease by stress perfusion cardiovascular magnetic resonance with the use of delayed enhancement infarction imaging. *J Am Coll Cardiol*, 2006. 47(8): p. 1630–8.
- Ingkanisorn, W.P., et al., Prognosis of negative adenosine stress magnetic resonance in patients presenting to an emergency department with chest pain. *J Am Coll Cardiol*, 2006. 47(7): p. 1427–32.
- Okuda, S., et al., Evaluation of ischemic heart disease on a 1.5 Tesla scanner: combined first-pass perfusion and viability study. *Radiat Med*, 2005. 23(4): p. 230–5.
- Sakuma, H., et al., Diagnostic accuracy of stress first-pass contrast-enhanced myocardial perfusion MRI compared with stress myocardial perfusion scintigraphy. *AJR Am J Roentgenol*, 2005. 185(1): p. 95–102.
- Plein, S., et al., Coronary artery disease: myocardial perfusion MR imaging with sensitivity encoding versus conventional angiography. *Radiology*, 2005. 235(2): p. 423–30.
- Takase, B., et al., Whole-heart dipyridamole stress first-pass myocardial perfusion MRI for the detection of coronary artery disease. *Jpn Heart J*, 2004. 45(3): p. 475–86.
- Paetsch, I., et al., Myocardial perfusion imaging using OMNISCAN: a dose finding study for visual assessment of stress-induced regional perfusion abnormalities. *J Cardiovasc Magn Reson*, 2004. 6(4): p. 803–9.
- Paetsch, I., et al., Comparison of dobutamine stress magnetic resonance, adenosine stress magnetic resonance, and adenosine stress magnetic resonance perfusion. *Circulation*, 2004. 110(7): p. 835–42.
- Wolff, S.D., et al., Myocardial first-pass perfusion magnetic resonance imaging: a multicenter dose-ranging study. *Circulation*, 2004. 110(6): p. 732–7.
- Thiele, H., et al., Color-encoded semiautomatic analysis of multi-slice first-pass magnetic resonance perfusion: comparison to tetrofosmin single photon emission computed tomography perfusion and X-ray angiography. *Int J Cardiovasc Imaging*, 2004. 20(5): p. 371–84; discussion 385–7.
- Plein, S., et al., Assessment of non-ST-segment elevation acute coronary syndromes with cardiac magnetic resonance imaging. *J Am Coll Cardiol*, 2004. 44(11): p. 2173–81.
- Bunce, N.H., et al., Combined coronary and perfusion cardiovascular magnetic resonance for the assessment of coronary artery stenosis. *J Cardiovasc Magn Reson*, 2004. 6(2): p. 527–39.
- Nagel, E., et al., Magnetic resonance perfusion measurements for the noninvasive detection of coronary artery disease. *Circulation*, 2003. 108(4): p. 432–437.
- Ishida, N., et al., Noninfarcted myocardium: correlation between dynamic first-pass contrast-enhanced myocardial MR imaging and quantitative coronary angiography. *Radiology*, 2003. 229(1): p. 209–16.
- Doyle, M., et al., The impact of myocardial flow reserve on the detection of coronary artery disease by perfusion imaging methods: an NHLBI WISE study. *J Cardiovasc Magn Reson*, 2003. 5(3): p. 475–85.
- Kinoshita, M., et al., Myocardial perfusion magnetic resonance imaging for diagnosing coronary arterial stenosis. *Jpn Heart J*, 2003. 44(3): p. 323–34.
- Ibrahim, T., et al., Assessment of coronary flow reserve: comparison between contrast-enhanced magnetic resonance imaging and positron emission tomography. *J Am Coll Cardiol*, 2002. 39(5): p. 864–70.
- Schwittler, J., et al., Assessment of myocardial perfusion in coronary artery disease by magnetic resonance: a comparison with positron emission tomography and coronary angiography. *Circulation*, 2001. 103(18): p. 2230–2235.
- Panting, J.R., et al., Echo-planar magnetic resonance myocardial perfusion imaging: parametric map analysis and comparison with thallium SPECT. *J Magn Reson Imaging*, 2001. 13(2): p. 192–200.
- Al-Saadi, N., et al., Noninvasive detection of myocardial ischemia from perfusion reserve based on cardiovascular magnetic resonance. *Circulation*, 2000. 101(12): p. 1379–83.
- Christian, T.F., et al., Absolute myocardial perfusion in canines measured by using dual-bolus first-pass MR imaging. *Radiology*, 2004. 232(3): p. 677–84.
- Gatehouse, P.D., et al., Accurate assessment of the arterial input function during high-dose myocardial perfusion cardiovascular magnetic resonance. *J Magn Reson Imaging*, 2004. 20(1): p. 39–45.
- Panting, J.R., et al., Abnormal subendocardial perfusion in cardiac syndrome X detected by cardiovascular magnetic resonance imaging. *N Engl J Med*, 2002. 346(25): p. 1948–53.
- Jahnke, C., et al., Prognostic value of cardiac magnetic resonance stress tests: adenosine stress perfusion and dobutamine stress wall motion imaging. *Circulation*, 2007. 115(13): p. 1769–76.

Contact

Andrew E. Arai, MD
Senior Investigator
National Heart,
Lung and Blood Institute
National Institutes
of Health
Bethesda, MD, 20892
USA

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- 1 Schwitter J, Nanz D, Kneifel S, et al. Assessment of myocardial perfusion in coronary artery disease by magnetic resonance: a comparison with positron emission tomography and coronary angiography. *Circulation*. 2001;103(18):2230–5.
- 2 Klem I, Heitner JF, Shah DJ, et al. Improved detection of coronary artery disease by stress perfusion cardiovascular magnetic resonance with the use of delayed enhancement infarction imaging. *J Am Coll Cardiol*. 2006;47(8):1630–8.
- 3 Schwitter J. MR IMPACT: First-pass MR perfusion test superior to SPECT for detection of CAD. Oral presentation. Annual meeting of European Society of Cardiology 2005, Stockholm, Sweden.
- 4 Pilz G, Bernhardt P, Klos M, et al. Clinical implication of adenosine-stress cardiac magnetic resonance imaging as potential gatekeeper prior to invasive examination in patients with AHA/ACC class II indication for coronary angiography. *Clin Res Cardiol*. 2006;95(10):531–8.
- 5 Ingkanisorn WP, Kwong RY, Bohme NS, et al. Prognosis of negative adenosine stress magnetic resonance in patients presenting to an emergency department with chest pain. *J Am Coll Cardiol*. 2006;47(7):1427–32.
- 6 Jahnke C, Nagel E, Gebker R, et al. Prognostic value of cardiac magnetic resonance stress tests: adenosine stress perfusion and dobutamine stress wall motion imaging. *Circulation*. 2007;115(13):1769–76.
- 7 Hardung D, v. der Recke G. Event free survival in patients with suspected coronary artery disease and a negative test result in adenosine stress perfusion magnetic resonance imaging. Oral presentation. SCMR Ninth Annual Scientific Sessions 2006, Miami, USA.
- 8 Miszalski K, v. der Recke G. A prospective study in 644 consecutive patients to assess safety and feasibility of adenosine stress magnetic resonance perfusion for diagnosis of myocardial ischemia. Oral presentation. SCMR Ninth Annual Scientific Sessions 2006, Miami, USA.

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