

siemens.com/somatom-drive

Absolute Dose Values in Computed Tomography

SOMATOM Drive

| Reference values | | Switzerland ¹ | Germany ² | European Union ³ | USA ⁴ |
|------------------|---------------------------|--------------------------|----------------------|-----------------------------|------------------|
| Head Routine | CTDI _{vol} [mGy] | 65 | 60 | 60 | 75 |
| Thorax Routine | CTDI _{vol} [mGy] | 15 | 10 | 30 | 21 |
| Abdomen Routine | CTDI _{vol} [mGy] | 15 | 15 | 35 | 25 |

| Default Siemens Protocol | | Standard IR values*/**/**** | Put your current values here |
|--------------------------|---------------------------|-----------------------------|------------------------------|
| Head Routine | CTDI _{vol} [mGy] | 37.9 | |
| Thorax Routine | CTDI _{vol} [mGy] | 4.22 | |
| Lung Low Dose SPS**** | CTDI _{vol} [mGy] | 0.61 | |
| Abdomen Routine | CTDI _{vol} [mGy] | 9.39 | |

1 Bundesamt für Gesundheit (Merkblatt R-06-06, Diagnostische Referenzwerte in der Computertomographie, 01.04.2010).

2 Bundesamt für Strahlenschutz (Diagnostische Referenzwerte für diagnostische und interventionelle Röntgenanwendungen vom 22. Juni 2016).

3 European Guidelines on Quality Criteria for Computed Tomography (<http://www.drs.dk/guidelines/ct/quality/htmlindex.htm>).

4 American College of Radiology (CT Accreditation Program Requirements, Clinical Image Quality Guide, Amended 2014).

*Values are based on the default protocols of the SOMATOM Drive with syngo CT VA62A and an average sized patient of 1.75 m and 75 kg

**Iterative Reconstruction ADMIRE is used

****In clinical practice, the use of ADMIRE may reduce CT patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task. Image quality as defined by low contrast detectability using a model observer method for evaluation. Equivalent low contrast detectability can be achieved with 80% to 85% less dose using ADMIRE at highest strength level for thin (0.6 mm) reconstruction slices in measured and simulated body and head phantoms for low contrast objects with different contrasts. See ADMIRE data sheet for further information.

****SPS: Selective Photon Shield/Tin Filter is used

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

Note: Any technical data contained in this document may vary within defined tolerances.

International version. Not for distribution in the US.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen

Germany

Phone: +49 9131 84 0

siemens.com/healthineers