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Ileana Jelescu, Ph.D.

Microstructure Mapping Lab, Department of Radiology,
Lausanne University Hospital (CHUV) and University of Lausanne
(UNIL), Lausanne, Switzerland

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Simona Schiavi, Ph.D.¹; Matilde Inglese, M.D., Ph.D.^{1,2}


¹ Department of Neuroscience, Rehabilitation, Ophthalmology,
Genetics, Maternal and Child Health (DINOEMI),
University of Genoa, Italy

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Josef Pfeuffer, Thorsten Feiweier, Heiko Meyer

Application Development, Siemens Healthineers AG, Erlangen,
Germany

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Wenjun Fan^{1,2}; Wen Shen¹

¹ Department of Radiology, Tianjin First Center Hospital,
Tianjin, China

² Armed Police Corps Hospital of Henan, Zhengzhou,
Henan, China

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¹ Department of Radiological, Oncological, and Pathological
Sciences, Sapienza University of Rome, Italy

² CNR-ISC Roma Sapienza, Physics Department, Rome, Italy

³ CNR-ISC Roma Sapienza, Physics Department and SAIMLAL
Department, Morphogenesis & Tissue Engineering, Sapienza
University of Rome, Italy

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¹ Department of Neuroscience, Rehabilitation, Ophthalmology,
Genetics, Maternal and Child Health (DINOEMI), University of
Genoa, Italy

¹ Institute of Diagnostic and Interventional Radiology, University
Hospital Zurich, University of Zurich, Switzerland

² Siemens Shenzhen Magnetic Resonance Ltd., Shenzhen, China

³ Siemens Healthineers, MR Application Development,
Erlangen, Germany

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Ileana Jelescu, Ph.D.

Microstructure Mapping Lab, Department of Radiology, Lausanne University Hospital (CHUV) and University of Lausanne (UNIL), Lausanne, Switzerland

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References

- Moseley ME, de Crespigny AJ, Roberts TP, Kozniowska E, Kucharczyk J. Early detection of regional cerebral ischemia using high-speed MRI. *Stroke*. 1993;24(12 Suppl):160–5.
- Lansberg MG, Thijs VN, O'Brien MW, Ali JO, de Crespigny AJ, Tong DC, et al. Evolution of apparent diffusion coefficient, diffusion-weighted, and T2-weighted signal intensity of acute stroke. *AJNR Am J Neuroradiol*. 2001;22(4):637–644.
- Panesar SS, Abhinav K, Yeh FC, Jacquesson T, Collins M, Fernandez-Miranda J. Tractography for Surgical Neuro-Oncology Planning: Towards a Gold Standard. *Neurotherapeutics*. 2019;16(1):36–51.
- Shahbandi A, Sattari SA, Haghshomar M, Shab-Bidar S, Lawton MT. Application of diffusion tensor-based tractography in treatment of brain arteriovenous malformations: a systematic review. *Neurosurg Rev*. 2023;46(1):115.
- Basser PJ, Mattiello J, LeBihan D. Estimation of the effective self-diffusion tensor from the NMR spin echo. *J Magn Reson B*. 1994;103(3):247–254.
- Chandarana H, Taouli B. Diffusion and perfusion imaging of the liver. *Eur J Radiol*. 2010;76(3):348–358.
- Girardet R, Dubois M, Manasseh G, Jreige M, Du Pasquier C, Canniff E, et al. The combination of non-contrast abbreviated MRI and alpha foetoprotein has high performance for hepatocellular carcinoma screening. *Eur Radiol*. 2023;33(10):6929–6938.
- Sobeh T, Inbar Y, Apter S, Soffer S, Anteby R, Kraus M, et al. Diffusion-weighted MRI for predicting and assessing treatment response of liver metastases from CRC – A systematic review and meta-analysis. *Eur J Radiol*. 2023;163:110810.
- Padhani AR, Schoots IG. Prostate cancer screening-stepping forward with MRI. *Eur Radiol*. 2023;33(10):6670–6676.
- lima M, Honda M, Sigmund EE, Ohno Kishimoto A, Kataoka M, Togashi K. Diffusion MRI of the breast: Current status and future directions. *J Magn Reson Imaging*. 2020;52(1):70–90.
- Partridge SC, Nissan N, Rahbar H, Kitsch AE, Sigmund EE. Diffusion-weighted breast MRI: Clinical applications and emerging techniques. *J Magn Reson Imaging*. 2017;45(2):337-355.
- Jelescu IO, Palombo M, Bagnato F, Schilling KG. Challenges for biophysical modeling of microstructure. *J Neurosci Methods*. 2020;344:108861.
- Griswold MA, Jakob PM, Heidemann RM, Nittka M, Jellus V, Wang J, et al. Generalized autocalibrating partially parallel acquisitions (GRAPPA). *Magn Reson Med*. 2002;47(6):1202–1210.
- Feinberg DA, Setsompop K. Ultra-fast MRI of the human brain with simultaneous multi-slice imaging. *J Magn Reson*. 2013;229:90–100.
- Moeller S, Pisharady PK, Ramanna S, Lenglet C, Wu X, Dowdle L, et al. NOise reduction with DIstribution Corrected (NORDIC) PCA in dMRI with complex-valued parameter-free locally low-rank processing. *Neuroimage*. 2021;226:117539.
- Veraart J, Novikov DS, Christiaens D, Ades-Aron B, Sijbers J, Fieremans E. Denoising of diffusion MRI using random matrix theory. *Neuroimage*. 2016;142:394–406.
- Alexander DC, Dyrby TB, Nilsson M, Zhang H. Imaging brain microstructure with diffusion MRI: practicality and applications. *NMR Biomed*. 2019;32(4):e3841.
- Novikov DS, Fieremans E, Jespersen SN, Kiselev VG. Quantifying brain microstructure with diffusion MRI: Theory and parameter estimation. *NMR Biomed*. 2019;32(4):e3998.
- Novikov DS, Kiselev VG, Jespersen SN. On modeling. *Magn Reson Med*. 2018;79(6):3172–3193.
- Jelescu IO, Fieremans E. Chapter 2 - Sensitivity and specificity of diffusion MRI to neuroinflammatory processes. In: Laule C, Port JD (Eds.). *Advances in Magnetic Resonance Technology and Applications. Imaging Neuroinflammation*. Academic Press. 2023:31–50. <https://doi.org/10.1016/B978-0-323-91771-1.00010-1>
- Jensen JH, Helpert JA, Ramani A, Lu H, Kaczynski K. Diffusional kurtosis imaging: the quantification of non-gaussian water

diffusion by means of magnetic resonance imaging.

Magn Reson Med. 2005;53(6):1432–1440.

- Jelescu IO, Veraart J, Adisetiyo V, Milla SS, Novikov DS, Fieremans E. One diffusion acquisition and different white matter models: how does microstructure change in human early development based on WMTI and NODDI? *Neuroimage*. 2015;107:242–256.
- Paydar A, Fieremans E, Nwankwo JI, Lazar M, Sheth HD, Adisetiyo V, et al. Diffusional kurtosis imaging of the developing brain. *AJNR Am J Neuroradiol*. 2014;35(4):808–814.
- Benitez A, Jensen JH, Falangola MF, Nietert PJ, Helpert JA. Modeling white matter tract integrity in aging with diffusional kurtosis imaging. *Neurobiol Aging*. 2018;70:265–275.
- Cho KIK, Kwak YB, Hwang WJ, Lee J, Kim M, Lee TY, et al. Microstructural Changes in Higher-Order Nuclei of the Thalamus in Patients With First-Episode Psychosis. *Biol Psychiatry*. 2019;85(1):70–78.
- Kochunov P, Rowland LM, Fieremans E, Veraart J, Jahanshad N, Eskandar G, et al. Diffusion-weighted imaging uncovers likely sources of processing-speed deficits in schizophrenia. *Proc Natl Acad Sci U S A*. 2016;113(47):13504–13509.
- Pavan T, Alemán-Gómez Y, Jenni R, Steullet P, Schilliger Z, Dwir D, et al. White Matter Microstructure Alterations and Their Link to Symptomatology in Early Psychosis and Schizophrenia. 2024. <https://doi.org/10.1101/2024.02.01.24301979>
- Dong JW, Jelescu IO, Ades-Aron B, Novikov DS, Friedman K, Babb JS, et al. Diffusion MRI biomarkers of white matter microstructure vary nonmonotonically with increasing cerebral amyloid deposition. *Neurobiol Aging*. 2020;89:118–128.
- de Kouchkovsky I, Fieremans E, Fleyshter L, Herbert J, Grossman RI, Inglese M. Quantification of normal-appearing white matter tract integrity in multiple sclerosis: a diffusion kurtosis imaging study. *J Neurol*. 2016;263(6):1146–1155.
- Bai Y, Lin Y, Tian J, Shi D, Cheng J, Haacke EM, et al. Grading of Gliomas by Using Monoexponential, Biexponential, and Stretched Exponential Diffusion-weighted MR Imaging and Diffusion Kurtosis MR Imaging. *Radiology*. 2016;278(2):496–504.
- Hui ES, Fieremans E, Jensen JH, Tabesh A, Feng W, Bonilha L, et al. Stroke assessment with diffusional kurtosis imaging. *Stroke*. 2012;43(11):2968–2273.
- Umesh Rudrapatna S, Wieloch T, Beirup K, Ruscher K, Mol W, Yanev P, et al. Can diffusion kurtosis imaging improve the sensitivity and specificity of detecting microstructural alterations in brain tissue chronically after experimental stroke? Comparisons with diffusion tensor imaging and histology. *Neuroimage*. 2014;97:363–373.
- Chen J, Chung S, Li T, Fieremans E, Novikov DS, Wang Y, et al. Identifying relevant diffusion MRI microstructure biomarkers relating to exposure to repeated head impacts in contact sport athletes. *Neuroradiol J*. 2023;36(6):693–701.
- Karlsen RH, Einarsen C, Moe HK, Håberg AK, Vik A, Skandsen T, et al. Diffusion kurtosis imaging in mild traumatic brain injury and postconcussional syndrome. *J Neurosci Res*. 2019;97(5):568–581.
- Rosenkrantz AB, Padhani AR, Chenevert TL, Koh DM, De Keyser F, Taouli B, et al. Body diffusion kurtosis imaging: Basic principles, applications, and considerations for clinical practice. *J Magn Reson Imaging*. 2015;42(5):1190–1202.
- Guo F, Leemans A, Viergever MA, Dell'Acqua F, De Luca A. Generalized Richardson-Lucy (GRL) for analyzing multi-shell diffusion MRI data. *Neuroimage*. 2020;218:116948.
- Jeurissen B, Tournier JD, Dhollander T, Connelly A, Sijbers J. Multi-tissue constrained spherical deconvolution for improved analysis of multi-shell diffusion MRI data. *Neuroimage*. 2014;103:411–426.
- Jelescu IO, Budde MD. Design and validation of diffusion MRI models of white matter. *Front Phys*. 2017;28:61.
- Guglielmetti C, Veraart J, Roelant E, Mai Z, Daans J, Van Audekerke J, et al. Diffusion kurtosis imaging probes cortical alterations and white matter pathology following cuprizone

induced demyelination and spontaneous remyelination.

Neuroimage. 2016;125:363–377.

- Jelescu IO, Zurek M, Winters KV, Veraart J, Rajaratnam A, Kim NS, et al. In vivo quantification of demyelination and recovery using compartment-specific diffusion MRI metrics validated by electron microscopy. *Neuroimage*. 2016;132:104–114.
- Mollink J, Kleinnijenhuis M, Cappellen van Walsum AV, Sotiropoulos SN, Cottaar M, Mirfin C, et al. Evaluating fibre orientation dispersion in white matter: Comparison of diffusion MRI, histology and polarized light imaging. *Neuroimage*. 2017;157:561–574.
- Salo RA, Belevich I, Jokitalo E, Gröhn O, Sierra A. Assessment of the structural complexity of diffusion MRI voxels using 3D electron microscopy in the rat brain. *Neuroimage*. 2021;225:117529.
- Sepehrband F, Clark KA, Ullmann JF, Kurniawan ND, Leverage G, Reutens DC, et al. Brain tissue compartment density estimated using diffusion-weighted MRI yields tissue parameters consistent with histology. *Hum Brain Mapp*. 2015;36(9):3687–3702.
- Veraart J, Fieremans E, Novikov DS. On the scaling behavior of water diffusion in human brain white matter. *Neuroimage*. 2019;185:379–387.
- Coelho S, Baete SH, Lemberskiy G, Ades-Aron B, Barrol G, Veraart J, et al. Reproducibility of the Standard Model of diffusion in white matter on clinical MRI systems. *Neuroimage*. 2022;257:119290.
- Zhang H, Schneider T, Wheeler-Kingshott CA, Alexander DC. NODDI: practical in vivo neurite orientation dispersion and density imaging of the human brain. *Neuroimage*. 2012;61(4):1000–1016.
- Fieremans E, Jensen JH, Helpert JA. White matter characterization with diffusional kurtosis imaging. *Neuroimage*. 2011;58(1):177–188.
- Diao Y, Jelescu I. Parameter estimation for WMTI-Watson model of white matter using encoder-decoder recurrent neural network. *Magn Reson Med*. 2023;89(3):1193–1206.
- Jespersen SN, Olesen JL, Hansen B, Shemesh N. Diffusion time dependence of microstructural parameters in fixed spinal cord. *Neuroimage*. 2018;182:329–342.
- Assaf Y, Freidlin RZ, Rohde GK, Basser PJ. New modeling and experimental framework to characterize hindered and restricted water diffusion in brain white matter. *Magn Reson Med*. 2004;52(5):965–978.
- Novikov DS, Veraart J, Jelescu IO, Fieremans E. Rotationally-invariant mapping of scalar and orientational metrics of neuronal microstructure with diffusion MRI. *Neuroimage*. 2018;174:518–538.
- Tristão Pereira C, Diao Y, Yin T, da Silva AR, Lanz B, Pierzchala K, et al. Synchronous nonmonotonic changes in functional connectivity and white matter integrity in a rat model of sporadic Alzheimer's disease. *Neuroimage*. 2021;225:117498.
- Jelescu IO, Veraart J, Fieremans E, Novikov DS. Degeneracy in model parameter estimation for multi-compartmental diffusion in neuronal tissue. *NMR Biomed*. 2016;29(1):33–47.
- Liao Y, Coelho S, Chen J, Ades-Aron B, Pang M, Osorio R, et al. Mapping tissue microstructure of brain white matter in vivo in health and disease using diffusion MRI. *ArXiv [Preprint]*. 2023:arXiv:2307.16386v2.
- Edwards LJ, Pine KJ, Ellerbrock I, Weiskopf N, Mohammadi S. NODDI-DTI: Estimating Neurite Orientation and Dispersion Parameters from a Diffusion Tensor in Healthy White Matter. *Front Neurosci*. 2017;11:720.
- Jelescu IO, de Skowronski A, Geffroy F, Palombo M, Novikov DS. Neurite Exchange Imaging (NEXI): A minimal model of diffusion in gray matter with inter-compartment water exchange. *Neuroimage*. 2022;256:119277.
- Olesen JL, Østergaard L, Shemesh N, Jespersen SN. Diffusion time dependence, power-law scaling, and exchange in gray matter. *Neuroimage*. 2022;251:118976.
- Lee HH, Olesen JL, Tian Q, Ramos Llorden G, Jespersen SN, Huang SY. Revealing diffusion time-dependence and exchange

effect in the in vivo human brain gray matter by using high gradient diffusion MRI. Proc. Intl. Soc. Mag. Reson. Med. 2022. London, UK, p. 254.

59 Uhl Q, Pavan T, Molendowska M, Jones DK, Palombo M, Jelescu I. Quantifying human gray matter microstructure using Neurite Exchange Imaging (NEXI) and 300 mT/m gradients. Imaging Neuroscience. 2024;2:1–19.

60 Uhl Q, Pavan T, Feiweier T, Piredda GF, Jespersen SN, Jelescu I. NEXI for the quantification of human gray matter microstructure on a clinical MRI scanner. Proc. Intl. Soc. Mag. Reson. Med. 2024. Singapore.

61 Uhl Q, Pavan T, Feiweier T, Canales-Rodríguez EJ, Jelescu I. Optimizing The NEXI Acquisition Protocol For Quantifying Human Gray Matter Microstructure On A Clinical MRI Scanner Using Explainable AI. Proc. Intl. Soc. Mag. Reson. Med. 2023. Presented at the Annual Meeting of the ISMRM, Toronto, Canada, p. 948.

62 Williamson NH, Ravin R, Cai TX, Falgairolle M, O'Donovan MJ, Basser PJ. Water exchange rates measure active transport and homeostasis in neural tissue. PNAS Nexus. 2023;2(3):pgad056.

63 Hertanu A, Uhl Q, Pavan T, Lamy CM, Jelescu IO. Quantifying features of human gray matter microstructure postmortem using NeuriteExchange Imaging (NEXI) at ultra-high field. Proc. Intl. Soc. Mag. Reson. Med. 2023. Toronto, Canada.

64 Xu J, Does MD, Gore JC. Quantitative characterization of tissue microstructure with temporal diffusion spectroscopy. J Magn Reson. 2009;200(2):189–197.

65 Merboldt KD, Hänicke W, Bruhn H, Gyngell ML, Frahm J. Diffusion imaging of the human brain in vivo using high-speed STEAM MRI. Magn Reson Med. 1992;23(1):179–192.

66 Szczepankiewicz F, Sjölund J, Ståhlberg F, Lätt J, Nilsson M. Tensor-valued diffusion encoding for diffusional variance decomposition (DIVIDE): Technical feasibility in clinical MRI systems. PLoS One. 2019;14(3):e0214238.

67 Jiang X, Li H, Xie J, McKinley ET, Zhao P, Gore JC, et al. In vivo imaging of cancer cell size and cellularity using temporal diffusion spectroscopy. Magn Reson Med. 2017;78(1):156–164.

68 Reynaud O, Winters KV, Hoang DM, Wadghiri YZ, Novikov DS, Kim SG. Pulsed and oscillating gradient MRI for assessment of cell size and extracellular space (POMACE) in mouse gliomas. NMR Biomed. 2016 Oct;29(10):1350–1363.

69 Zhang H, Liu K, Ba R, Zhang Z, Zhang Y, Chen Y, et al. Histological and molecular classifications of pediatric glioma with time-dependent diffusion MRI-based microstructural mapping. Neuro Oncol. 2023;25(6):1146–1156.

70 Zhu A, Shih R, Huang RY, DeMarco JK, Bhushan C, Morris HD, et al. Revealing tumor microstructure with oscillating diffusion encoding MRI in pre-surgical and post-treatment glioma patients. Magn Reson Med. 2023;90(5):1789–1801.

71 Wu D, Jiang K, Li H, Zhang Z, Ba R, Zhang Y, et al. Time-Dependent Diffusion MRI for Quantitative Microstructural Mapping of Prostate Cancer. Radiology. 2022;303(3):578–587.

72 Szczepankiewicz F, van Westen D, Englund E, Westin CF, Ståhlberg F, Lätt J, et al. The link between diffusion MRI and tumor heterogeneity: Mapping cell eccentricity and density by diffusional variance decomposition (DIVIDE). Neuroimage. 2016;142:522–532.

73 Topgaard D. Multidimensional diffusion MRI.

J Magn Reson. 2017;275:98–113.

74 Westin CF, Knutsson H, Pasternak O, Szczepankiewicz F, Özsarslan E, van Westen D, et al. Q-space trajectory imaging for multidimensional diffusion MRI of the human brain. Neuroimage. 2016;135:345–362.

75 Szczepankiewicz F, Lasič S, van Westen D, Sundgren PC, Englund E, Westin CF, et al. Quantification of microscopic diffusion anisotropy disentangles effects of orientation dispersion from microstructure: applications in healthy volunteers and in brain tumors. Neuroimage. 2015;104:241–252.

76 Nilsson M, Englund E, Szczepankiewicz F, van Westen D, Sundgren PC. Imaging brain tumour microstructure. Neuroimage. 2018;182:232–250.

77 Nery F, Szczepankiewicz F, Kerkelä L, Hall MG, Kaden E, Gordon I, et al. In vivo demonstration of microscopic anisotropy in the human kidney using multidimensional diffusion MRI. Magn Reson Med. 2019;82(6):2160–2168.

78 Fokkinga E, Hernandez-Tamames JA, Ianus A, Nilsson M, Tax CMW, Perez-Lopez R, et al. Advanced Diffusion-Weighted MRI for Cancer Microstructure Assessment in Body Imaging, and Its Relationship With Histology. J Magn Reson Imaging. 2023. doi: 10.1002/jmri.29144. Epub ahead of print.

79 Reynaud O. Time-Dependent Diffusion MRI in Cancer: Tissue Modeling and Applications. Front. Phys. 2017;5. https://doi.org/10.3389/fphy.2017.00058

80 Lemberskiy G, Fieremans E, Veraart J, Deng FM, Rosenkrantz AB, Novikov DS. Characterization of prostate microstructure using water diffusion and NMR relaxation. Front Phys. 2018;6:91.

81 Lemberskiy G, Rosenkrantz AB, Veraart J, Taneja SS, Novikov DS, Fieremans E. Time-Dependent Diffusion in Prostate Cancer. Invest Radiol. 2017;52(7):405–411.

82 Lemberskiy G, Feiweier T, Gyftopoulos S, Axel L, Novikov DS, Fieremans E. Assessment of myofiber microstructure changes due to atrophy and recovery with time-dependent diffusion MRI. NMR Biomed. 2021;34(7):e4534.

83 McDowell AR, Feiweier T, Muntoni F, Hall MG, Clark CA. Clinically feasible diffusion MRI in muscle: Time dependence and initial findings in Duchenne muscular dystrophy. Magn Reson Med. 2021;86(6):3192–3200.

84 Piredda GF, Hilbert T, Granziera C, Bonnier G, Meuli R, Molinari F, et al. Quantitative brain relaxation atlases for personalized detection and characterization of brain pathology. Magn Reson Med. 2020;83(1):337–351.

85 Pavan T, Alemán-Gómez Y, Jenni R, Cleusix M, Alameda L, Do KQ, et al. Normative Modelling of White Matter Microstructure in Early Psychosis and Schizophrenia. Proc. Intl. Soc. Mag. Reson. Med. 2023. Toronto, Canada, p. 651.

86 Szczepankiewicz F, Lasič S, van Westen D, Sundgren PC, Englund E, Westin CF, et al. Quantification of microscopic diffusion anisotropy disentangles effects of orientation dispersion from microstructure: applications in healthy volunteers and in brain tumors. Neuroimage. 2015;104:241-252.

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Simona Schiavi, Ph.D.¹; Matilde Inglese, M.D., Ph.D.^{1,2}

¹Department of Neuroscience, Rehabilitation, Ophthalmology, Genetics, Maternal and Child Health (DINOgMI), University of Genoa, Italy

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References

- Inglese M, Petracca M. MRI in multiple sclerosis: Clinical and research update. Curr Opin Neurol. 2018;31(3):249–255.
- Basser PJ, Mattiello J, LeBihan D. MR diffusion tensor spectroscopy and imaging. Biophys J. 1994;66(1):259–267.
- Novikov DS, Kiselev VG, Jespersen SN. On modeling. Magn Reson Med. 2018;79(6):3172–3193.
- Alexander DC, Dyrby TB, Nilsson M, Zhang H. Imaging brain microstructure with diffusion MRI: practicality and applications. NMR Biomed. 2019;32(4):e3841.
- Jelescu IO, Budde MD. Design and Validation of Diffusion MRI Models of White Matter. Front Phys. 2017;28:61. Epub 2017 Nov 28.
- Basser PJ, Pajevic S, Pierpaoli C, Duda J, Aldroubi A. In vivo fiber tractography using DT-MRI data. Magn Reson Med. 2000;44(4):625–632.
- Zhang F, Daducci A, He Y, Schiavi S, Seguin C, Smith RE, et al. Quantitative mapping of the brain’s structural connectivity using diffusion MRI tractography: A review. Neuroimage. 2022;249:118870.
- Cercignani M, Lannucci G, Filippi M. Diffusion-weighted imaging in multiple sclerosis. Ital J Neurol Sci. 1999;20(5 Suppl):S246–249.
- Inglese M, Bester M. Diffusion imaging in multiple sclerosis: Research and clinical implications. NMR Biomed. 2010;23(7):865–872.
- Filippi M, Cercignani M, Inglese M, Horsfield MA, Comi G. Diffusion tensor magnetic resonance imaging in multiple sclerosis. Neurology. 2001;56(3):304–311.
- Hygino da Cruz LC, Batista RR, Domingues RC, Barkhof F. Diffusion magnetic resonance imaging in multiple sclerosis. Neuroimaging Clin N Am. 2011;21(1)71–88, vii-viii.
- Cercignani M, Gandini Wheeler-Kingshott C. From micro- to macro-structures in multiple sclerosis: what is the added value of diffusion imaging. NMR Biomed. 2019;32(4):e3888.
- Pagani E, Bammer R, Horsfield MA, Rovaris M, Gass A, Ciccarelli O, et al. Diffusion MR imaging in multiple sclerosis: technical aspects and challenges. AJNR Am J Neuroradiol. 2007;28(3):411–420.
- Zhang H, Schneider T, Wheeler-Kingshott CA, Alexander DC. NODDI: Practical in vivo neurite orientation dispersion and density imaging of the human brain. Neuroimage. 2012;61(4):1000–1016.
- Kaden E, Kelm ND, Carson RP, Does MD, Alexander DC. Multi-compartment microscopic diffusion imaging. Neuroimage. 2016;139:346–359.
- Lakhani DA, Schilling KG, Xu J, Bagnato F. Advanced Multicompartment Diffusion MRI Models and Their Application in Multiple Sclerosis. AJNR Am J Neuroradiol. 2020;41(5):751–757.
- Rahmanzadeh R, Lu PJ, Barakovic M, Weigel M, Maggi P, Nguyen TD, et al. Myelin and axon pathology in multiple sclerosis assessed by myelin water and multi-shell diffusion imaging. Brain. 2021;144(6):1684–1696.
- Wang X, Cusick MF, Wang Y, Sun P, Libbey JE, Trinkaus K, et al. Diffusion basis spectrum imaging detects and distinguishes coexisting subclinical inflammation, demyelination and axonal injury in experimental autoimmune encephalomyelitis mice. NMR Biomed. 2014;27(7):843–852.
- Chiang CW, Wang Y, Sun P, Lin TH, Trinkaus K, Cross AH, et al. Quantifying white matter tract diffusion parameters in the presence of increased extra-fiber cellularity and vasogenic edema.

Neuroimage. 2014;101:310–319.

- Schiavi S, Petracca M, Sun P, Fleysher L, Cocozza S, El Mendili MM, et al. Non-invasive quantification of inflammation, axonal and myelin injury in multiple sclerosis. Brain. 2021;144(1):213–223.
- Wang Y, Wang Q, Haldar JP, Yeh FC, Xie M, Sun P, et al. Quantification of increased cellularity during inflammatory demyelination. Brain. 2011;134(Pt 12):3590–3601.
- Jensen JH, Helpert JA, Ramani A, Lu H, Kaczynski K. Diffusional kurtosis imaging: The quantification of non-gaussian water diffusion by means of magnetic resonance imaging. Magn Reson Med. 2005;53(6):1432–1440.
- de Kouchkovsky I, Fieremans E, Fleysher L, Herbert J, Grossman RI, Inglese M. Quantification of normal-appearing white matter tract integrity in multiple sclerosis: a diffusion kurtosis imaging study. J Neurol. 2016;263(6):1146–1155.
- Todea RA, Lu PJ, Fartaria MJ, Bonnier G, Du Pasquier R, Krueger G, et al. Evolution of Cortical and White Matter Lesion Load in Early-Stage Multiple Sclerosis: Correlation With Neuroaxonal Damage and Clinical Changes. Front Neurol. 2020;11:973.
- Palombo M, Ianus A, Guerrieri M, Nunes D, Alexander DC, Shemesh N, et al. SANDI: A compartment-based model for non-invasive apparent soma and neurite imaging by diffusion MRI. Neuroimage. 2020;215:116835.
- Genc S, Chamberland M, Koller K, Tax CMW, Zhang H, Palombo M, et al. Repeatability of soma and neurite metrics in cortical and subcortical grey matter. bioRxiv. 2020.
- Xu J, Moeller S, Auerbach EJ, Strupp J, Smith SM, Feinberg DA, et al. Evaluation of slice accelerations using multiband echo planar imaging at 3 T. Neuroimage. 2013; 83:991–1001.
- Feinberg DA, Moeller S, Smith SM, Auerbach E, Ramanna S, Gunther M, et al. Multiplexed echo planar imaging for sub-second whole brain fMRI and fast diffusion imaging. PLoS One. 2010;5(12):e15710.
- Moeller S, Yacoub E, Olman CA, Auerbach E, Strupp J, Harel N, et al. Multiband multislice GE-EPI at 7 tesla, with 16-fold acceleration using partial parallel imaging with application to high spatial and temporal whole-brain fMRI. Magn Reson Med. 2010;63(5):1144–1153.
- Lipp I, Parker GD, Tallantyre EC, Goodall A, Grama S, Patitucci E, et al. Tractography in the presence of multiple sclerosis lesions. Neuroimage. 2019;116471.
- Petracca M, Schiavi S, Battocchio M, El Mendili MM, Fleysher L, Daducci A, et al. Streamline density and lesion volume reveal a postero–anterior gradient of corpus callosum damage in multiple sclerosis. Eur J Neurol. 2020;27(6):1076–1082.
- Schiavi S, Petracca M, Battocchio M, El Mendili MM, Paduri S, Fleysher L, et al. Sensory-motor network topology in multiple sclerosis: Structural connectivity analysis accounting for intrinsic density discrepancy. Hum Brain Mapp. 2020 ;41(11):2951–2963.
- Schiavi S, Azzari A, Mensi A, et al. Classification of multiple sclerosis patients based on structural disconnection: a robust feature selection approach. J Neuroimaging. 2022.
- Bosticardo S, Schiavi S, Schaedelin S, Lu PJ, Barakovic M, Weigel M, et al. Microstructure-Weighted Connectomics in Multiple Sclerosis. Brain Connect. 2021;12(1):6–17.

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Josef Pfeuffer, Thorsten Feiweier, Heiko Meyer

Application Development, Siemens Healthineers AG, Erlangen, Germany

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References

- 1 Jezzard P, Balaban RS. Correction for geometric distortion in echo planar images from B0 field variations. Magn Reson Med. 1995;34(1):65–73.
- 2 Pfeuffer J, Van de Moortele PF, Ugurbil K, Hu X, Glover GH. Correction of physiologically induced global off-resonance effects in dynamic echo-planar and spiral functional imaging. Magn Reson Med. 2002;47(2):344–353.

- 3 Andersson JL, Skare S, Ashburner J. How to correct susceptibility distortions in spin-echo echo-planar images: application to diffusion tensor imaging. Neuroimage. 2003;20(2):870–888.
- 4 Raj D, Paley DP, Anderson AW, Kennan RP, Gore JC. A model for susceptibility artefacts from respiration in functional echo-planar magnetic resonance imaging. Phys Med Biol. 2000;45(12):3809–3820.
- 5 Glover GH, Lai S. Self-navigated spiral fMRI: interleaved versus single-shot. Magn Reson Med. 1998;39(3):361–368.
- 6 Noll DC, Pauly JM, Meyer CH, Nishimura DG, Macovski A. Deblurring for non-2D Fourier transform magnetic resonance imaging. Magn Reson Med. 1992;25(2):319–333.
- 7 Irarrazabal P, Meyer CH, Nishimura DG, Macovski A. Inhomogeneity correction using an estimated linear fieldmap. Magn Reson Med 1996;35(2):278–282.
- 8 Andersson JL, Hutton C, Ashburner J, Turner R, Friston K. Modeling geometric deformations in EPI time series. Neuroimage. 2001;13(5):903–919.
- 9 Durand E, Van De Moortele PF, Pachot-Clouard M, Le Bihan D. Artifact due to B(0) fluctuations in fMRI: correction using the k-space central line. Magn Reson Med. 2001;46(1):198–201.
- 10 Liu HL, Kochunov P, Lancaster JL, Fox PT, Gao JH. Comparison of navigator echo and centroid corrections of image displacement induced by static magnetic field drift on echo planar functional MRI. J Magn Reson Imaging. 2001;13(2):308–312.
- 11 Holland D, Kuperman JM, Dale AM. Efficient correction of inhomogeneous static magnetic field-induced distortion in Echo Planar Imaging. Neuroimage. 2010;50(1):175–183.
- 12 Hutton C, Bork A, Josephs O, Deichmann R, Ashburner J, Turner R. Image distortion correction in fMRI: A quantitative evaluation. Neuroimage. 2002;16(1):217–240.
- 13 Tong A, Lemberskiy G, Huang C, Shanbhogue K, Feiweier T, Rosenkrantz AB. Exploratory study of geometric distortion correction of prostate diffusion-weighted imaging using B0 map acquisition. J Magn Reson Imaging. 2019;50(5):1614–1619.
- 14 Graunhan NF, Grünebach N, Brockstedt L, Sanner A, Feiweier T, Schöffling V, et al. Reduction of Distortion Artifacts in Brain MRI Using a Field Map-based Correction Technique in Diffusion-weighted Imaging: A Prospective Study. Clin Neuroradiol. 2023. doi: 10.1007/s00062-023-01338-3. Epub ahead of print.
- 15 Dagher J, Reese T, Bilgin A. High-resolution, large dynamic range field map estimation. Magn Reson Med. 2014;71(1):105–117.
- 16 Fautz HP, Gross P, Gumbrecht R. Method and computer to determine a B0 field map with a magnetic resonance apparatus. United States patent US 20150077115A1. 2015.
- 17 Gumbrecht R, Koehler M, Schneider R. Method and magnetic resonance apparatus for correction of a B0 map for chemical shifts. United States patent US 9830711B2. 2015.
- 18 Bredfeldt JS, Miao X, Kaza E, Schneider M, Requardt M, Feiweier T, et al. Patient specific distortion detection and mitigation in MR images used for stereotactic radiosurgery. Phys Med Biol. 2022;67(6):10.1088/1361-6560/ac508e.
- 19 Jin J, Tourell M, Sati P, Patil S, Liu K, Derbyshire J, et al. Segmented 3D EPI with CAIPIRINHA for Fast, High-Resolution T2*-weighted Imaging. Proceedings ISMRM 2021.
- 20 Jin J, O'Brien K, et al. Fast, High-resolution Whole Brain SWI and QSM with CAIPIRINHA 3D-EPI and Deep Learning Reconstruction. Proceedings ISMRM 2024.

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


References

- 1Levey AS, Bosch JP, Lewis JB, Greene T, Rogers N, Roth D: A more accurate method to estimate glomerular filtration rate from serum creatinine: a new prediction equation. Modification of Diet in Renal Disease Study Group. Ann Intern Med 1999, 130(6):461-470.
- 2aimes C, Darge K, Khrichenko D, Carson RH, Berman JI: Diffusion tensor imaging and tractography of the kidney in children: feasibility and preliminary experience. Pediatr Radiol 2013.
- 3Lanzman RS, Ljmani A, Pentang G, Zgoura P, Zenginli H, Kropil P, Heusch P, Schek J, Miese FR, Blondin D et al: Kidney transplant: functional assessment with diffusion-tensor MR imaging at 3T. Radiology 2013, 266(1):218-225.

Contact

Shen Wen, M.D.
Chief Radiologist
Dept. of Radiology
Tianjin First Center Hospital
24 Fukang road, Nankai district,
Tianjin, China, 300192
Phone: +86-11-22-23626583
shenwen66happy@163.com

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¹ Department of Radiological, Oncological, and Pathological Sciences, Sapienza University of Rome, Italy
² CNR-ISC Roma Sapienza, Physics Department, Rome, Italy
³ CNR-ISC Roma Sapienza, Physics Department and SAIMLAL Department, Morphogenesis & Tissue Engineering, Sapienza University of Rome, Italy

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¹ Institute of Diagnostic and Interventional Radiology, University Hospital Zurich, University of Zurich, Switzerland
² Siemens Shenzhen Magnetic Resonance Ltd., Shenzhen, China
³ Siemens Healthineers, MR Application Development, Erlangen, Germany
⁴ Kantonsspital Münsterlingen, Switzerland
⁵ Inselspital Bern, Switzerland

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
References

- Regnault TR, Galan HL, Parker TA, Anthony RV. Placental development in normal and compromised pregnancies—a review. Placenta. 2002; 23 Suppl A:S119-29.
- Benirschke K, Burton GJ, Baergen RN. Pathology of the Human Placenta. 6th ed. Springer; 2012.
- Serov AS, Salafia CM, Filoche M, Grebenkov DS. Analytical theory of oxygen transport in the human placenta. J Theor Biol. 2015; 368:133-44.
- Broses I, Pijnenborg R, Vercruyssen L, Romero R. The “Great Obstetrical Syndromes” are associated with disorders of deep placentation. Am J Obstet Gynecol. 2011; 204(3):193-201.
- Baschat AA, Hecher K. Fetal growth restriction due to placental disease. Semin Perinatol. 2004; 28(1):67-80.
- Serov AS, Salafia C, Grebenkov DS, Filoche M. The role of morphology in mathematical models of placental gas exchange. J Appl Physiol. 2016; 120:17-28.
- Siauve N, Chalouhi GE, Deloison B, Alison M, Clement O, Ville Y, Salomon LJ. Functional imaging of the human placenta with magnetic resonance. Am J Obstet Gynecol. 2015; 213 Suppl 4:S103-14
- Blaicher W, Brugger PC, Mittermayer C, Schwindt J, Deutinger J, Bernaschek G, Prayer D. Magnetic resonance imaging of the normal placenta. Eur J Radiol. 2006; 57(2):256-60.
- Basser JP. Inferring microstructural features and the physiological state of tissues from diffusion-weighted images. NMR Biomed. 1995; 8(7-8):333-44.
- Moore R, Issa B, Tokarczuk P, Duncan KR, Boulby P, Baker PN, Bowtell RW, Worthington BS, Johnson IR, Gowland PA. In vivo intravoxel incoherent motion measurements in the human placenta using echo-planar imaging at 0.5T. Magn Reson Med. 2000; 43(2):295-302.
- Moore RJ, Strachan BK, Tyler DJ, Duncan KR, Boulby P, Baker PN, Worthington BS, Johnson IR, Gowland PA. In utero Perfusing Fraction Maps in Normal and Growth Restricted Pregnancy Measured Using IVIM Echo-Planar MRI. Placenta. 2000; 21:726-32.
- Linduska N, Dekan S, Messerschmidt A, Kasprian G, Brugger PC, Chalubinski K, Weber M, Prayer D. Placental pathologies in fetal MRI with pathohistological correlation. Placenta. 2009; 30(6):555-9.
- Manganaro L, Fierro F, Tomei A, La Barbera L, Savelli S, Sollazzo P, Sergi ME, Vinci V, Ballesio L, Marini M. MRI and DWI: feasibility of DWI and ADC maps in the evaluation of placental changes during gestation. Prenat Diagn. 2010; 30(12-13):1178-84.
- Sivrioğlu AK, Özcan ÜA, Türk A, Ulus S, Yıldız ME, Sönmez G, Mutlu. Evaluation of the placenta with relative apparent diffusion coefficient and T2 signal intensity analysis. Diagn Interv Radiol. 2013; 19:495-500.
- Sohlberg S, Mulic-Lutvica A, Lindgren P, Ortiz-Nieto F, Wikström AK, Wikström J. Placental perfusion in normal pregnancy and early and late preeclampsia: a magnetic resonance imaging study. Placenta. 2014; 35(3):202-6.

- Villar J, Ismail LC, Victora CG, Ohuma EO, Bertino E, Altman DG, Lambert A, Papageorgiou AT, Carvalho M, Jaffer YA, Gravett MG, Purwar M, Frederick IO, Noble AJ, Pang R, Barros FC, Chumlea C, Bhutta ZA, Kennedy SH; International Fetal and Newborn Growth Consortium for the 21st Century (INTERGROWTH-21st). International standards for newborn weight, length, and head circumference by gestational age and sex: the Newborn Cross-Sectional Study of the INTERGROWTH-21st Project. Lancet. 2014; 384(9946):857-68.
- Niendorf T, Dijkhuizen RM, Norris DG, Van Lookeren Campagne M, Nicolay K. Biexponential diffusion attenuation in various states of brain tissue: implications for diffusion-weighted imaging. Magn Reson Med. 1996; 36(6):847-57.
- Mulkern RV, Gudbjartsson H, Westin CF, Zengingonul HP, Gartner W, Guttman CR, Holtzman D, Kyriakos W, Jolesz FA, Maier SE. Multi-component apparent diffusion coefficients in human brain. NMR Biomed. 1999; 12:51-62.
- De Santis S, Gabrielli A, Palombo M, Maraviglia B, Capuani S. Non-Gaussian diffusion imaging: a brief practical review. Magn Reson Imaging. 2011; 29(10):1410-6.
- Nair G, Hu XP. Manifestation and Post-hoc Correction of Gradient Cross-term Artifacts in DTI. Magn Reson Imaging. 2012; 30(6):764-73.
- Rossi C, Boss A, Martirosian P, Steidle G, Capuani S, Claussen CD, Maraviglia B, Schick F. Influence of steady background gradients on the accuracy of molecular diffusion anisotropy measurements. Magn Reson Imaging 2008; 26:1250-8.
- De Santis S, Rebuzzi M, Di Pietro G, Fasano F, Maraviglia B, Capuani S. In vitro and in vivo MR evaluation of internal gradient to assess trabecular bone density. Phys Med Biol. 2010; 55(19):5767.
- Capuani S, Guerreri M, Antonelli A, Bernardo S, Porpora MG, Giacotti A, Catalano C, Manganaro L. Diffusion and perfusion quantified by Magnetic Resonance Imaging are markers of human placenta development in normal pregnancy. Placenta. 2017; 58:33-9.
- Baergen RN. Manual of Pathology of the Human Placenta. 2nd ed. Springer; 2011.
- Avni R, Neeman M, Garbow JR. Functional MRI of the placenta – From rodents to humans. Placenta. 2015; 36(6):615-22.
- Solomon E, Avni R, Hadas R, Raz T, Garbow JR, Bendel P, Frydman L. Major mouse placental compartments revealed by diffusion-weighted MRI, contrast-enhanced MRI, and fluorescence imaging. Proc Natl Acad Sci U S A. 2014; 111(28):10353-8.
- Capuani S, Manenti G, Iundusi R, Tarantino U. Focus on Diffusion MR Investigations of Musculoskeletal Tissue to Improve Osteoporosis Diagnosis: A Brief Practical Review. BioMed Res Int. 2015; 2015:948610.
- Le Bihan D, Breton E, Lallemand D, Aubin ML, Vignaud J, Laval-Jeantet M. Separation of diffusion and perfusion in intravoxel incoherent motion MR imaging. Radiology. 1988; 168(2):497-505.

Contact

Lucia Manganaro, M.D.
Department of Radiological,
Oncological and Pathological
Sciences
Sapienza University of Rome

V.le Regina Elena 324
00161 Roma
Italy
Phone +39 06 49972639
lucia.manganaro@uniroma1.it

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
Parameter	RESOLVE (no slice acceleration) 4 min 21 sec	SMS RESOLVE (2-fold acceleration) 2 min 35 sec	SMS RESOLVE (3-fold acceleration) 1 min 44 sec
Echo time [ms]	56	55	55
Distance factor	20%		
Readout segments	5		
Field-of-view [mm²]	340 x 136		
Echo spacing [ms]	0.38		
Averages	b=50 (1), b=800 (4)		
Table 1: Imaging parameters of conventional RESOLVE and SMS RESOLVE at our institution.			

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
References

- 1 Bickel H, Pinker-Domenig K, Bogner W, Spick C, Bago-Horvath Z, Weber M, et al. Quantitative apparent diffusion coefficient as a noninvasive imaging biomarker for the differentiation of invasive breast cancer and ductal carcinoma in situ. Invest Radiol. 2015;50(2):95-100.
- 2 Bogner W, Gruber S, Pinker K, Grabner G, Stadlbauer A, Weber M, et al. Diffusion-weighted MR for differentiation of breast lesions at 3.0 T: how does selection of diffusion protocols affect diagnosis? Radiology. 2009;253(2):341-51.
- 3 Pinker K, Bogner W, Baltzer P, Karanikas G, Magometschnigg H, Brader P, et al. Improved differentiation of benign and malignant breast tumors with multiparametric 18fluorodeoxyglucose positron emission tomography magnetic resonance imaging: a feasibility study. Clin Cancer Res. 2014;20(13):3540-9.
- 4 Iacconi C, Giannelli M, Marini C, Cilotti A, Moretti M, Viacava P, et al. The role of mean diffusivity (MD) as a predictive index of the response to chemotherapy in locally advanced breast cancer: a preliminary study. Eur Radiol. 2010;20(2):303-8.
- 5 Liu S, Ren R, Chen Z, Wang Y, Fan T, Li C, et al. Diffusion-weighted imaging in assessing pathological response of tumor in breast cancer subtype to neoadjuvant chemotherapy. Journal of magnetic resonance imaging : JMRI. 2015;42(3):779-87.
- 6 Porter DA, Heidemann RM. High resolution diffusion-weighted imaging using readout-segmented echo-planar imaging, parallel imaging and a two-dimensional navigator-based reacquisition. Magnetic resonance in medicine. 2009;62(2):468-75.
- 7 Bogner W, Pinker-Domenig K, Bickel H, Chmelik M, Weber M, Helbich TH, et al. Readout-segmented echo-planar imaging improves the diagnostic performance of diffusion-weighted MR breast examinations at 3.0 T. Radiology. 2012;263(1):64-76.
- 8 Frost R, Jezard P, Douaud G, Clare S, Porter DA, Miller KL. Scan time reduction for readout-segmented EPI using simultaneous multislice acceleration: Diffusion-weighted imaging at 3 and 7 Tesla. Magnetic resonance in medicine. 2014.
- 9 Setsompop K, Gagoski BA, Polimeni JR, Witzel T, Wedeen VJ, Wald LL. Blipped-controlled aliasing in parallel imaging for simultaneous multislice echo planar imaging with reduced g-factor penalty. Magnetic resonance in medicine. 2012;67(5):1210-24.
- 10 Filli L, Ghafoor S, Kenkel D, Liu W, Weiland E, Andreisek G, et al. Simultaneous multi-slice readout-segmented echo planar imaging for accelerated diffusion-weighted imaging of the breast. Eur J Radiol. 2016;85(1):274-8.

Contact

Lukas Filli, M.D.
Institute of Diagnostic and Interventional Radiology
University Hospital Zurich, University of Zurich
Raemistrasse 100
CH-8091 Zurich
Switzerland
Phone: +41 44 225 11 11
lukas.filli@usz.ch

MAGNETOM Cima.X with Gemini Gradients

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