

**DICOM Conformance Statement**  
***AI-Rad Companion Chest CT***



This Statement is also valid for higher software versions - until a document for that higher version is published.



# 1 Overview

AI-Rad Companion Chest CT is a cloud-based/on-premises(edge) application, which communicates indirectly with other DICOM nodes since it makes use of the network services provided by teamplay DICOM Hub.

AI-Rad Companion Chest CT:

- Receives the input DICOM data from teamplay DICOM Hub after a storage request to teamplay Receiver using the configurable AET (E.g., "AIRC").
- Displays images to a user (browser-based viewer application).
- Generate result objects CT Image, Comprehensive SR and Secondary Capture Image.
- Stores result DICOM data via teamplay DICOM Hub and teamplay Receiver to one or several target DICOM nodes configured in teamplay DICOM Hub. AI-Rad Companion Chest CT conforms to the DICOM Standard [2] and supports a subset of the storage SOP classes supported by teamplay DICOM Hub, as described in Table 1: Network Services.

Please refer to the DICOM Conformance Statements of teamplay DICOM Hub [1] for further information on the provided network services.

NOTE: This DICOM Conformance Statement is applicable for AI-Rad Companion Chest CT (Pulmonary) VA40A, Chest CT (Cardiovascular) VA22A and Chest CT (Musculoskeletal) VA22A and later until superseded by a more recent document applicable to a more recent version.

Table 1: Network Services

SOP Classes	SOP Class UID	User of Service (SCU)		Provider of Service (SCP)	
		Create	Send	Store	Display
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Yes	No <sup>1</sup>	No <sup>1</sup>	Yes
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	Yes	No <sup>1</sup>	No <sup>1</sup>	No
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Yes	No <sup>1</sup>	No <sup>1</sup>	No

<sup>1</sup> Network communication is performed by teamplay DICOM hub. Please refer to the DICOM Conformance Statements of teamplay DICOM Hub [1] for further information.

Table 2: Media Services

Media Storage Application Profile	Write Files (FSC or FSU)	Read Files (FSR)
NOT APPLICABLE		

Table 3: Implementation Identifying Information

Name	Value
Implementation Class UID	1.3.12.2.1107.5.8.21
Implementation Version Name	AI-Rad Companion

## 2 Table of Contents

1	OVERVIEW .....	2
2	TABLE OF CONTENTS .....	3
3	INTRODUCTION .....	5
3.1	Revision History.....	5
3.2	Audience .....	5
3.3	Remarks.....	5
3.4	Definitions, Terms and Abbreviations.....	5
3.5	References.....	6
4	NETWORKING .....	7
5	MEDIA INTERCHANGE.....	8
6	TRANSFORMATIONS OF DICOM TO CDA .....	9
7	SUPPORT OF EXTENDED CHARACTER SETS.....	10
7.1	Character sets for AI-Rad Companion Chest CT .....	10
8	ATTRIBUTE CONFIDENTIALITY PROFILES .....	13
8.1	De-identification .....	13
9	SECURITY .....	14
9.1	Security Profiles .....	14
9.2	Association Level Security.....	14
9.3	Application Level Security.....	14
10	ANNEXES.....	15
10.1	IOD Contents .....	15
10.1.1	Created SOP Instances.....	15
10.1.1.1	Comprehensive SR.....	15
10.1.1.2	Secondary Capture Report .....	21
10.1.1.3	VRT Results & MPR Series (Secondary Capture) .....	25
10.1.1.4	CT Image Storage .....	30
10.1.2	Usage of Attributes from Received IODs .....	32
10.1.3	Attribute Mapping .....	32
10.1.4	Coerced/Modified Fields .....	32
10.2	Data Dictionary of Private Attributes .....	32
10.3	Coded Terminology and Templates .....	32
10.3.1	Context Groups .....	32



10.3.2	Template Specifications.....	35
10.3.3	Private Code definitions.....	45
<b>10.4</b>	<b>Grayscale Image Consistency.....</b>	<b>48</b>
<b>10.5</b>	<b>Standard Extended / Specialized / Private SOP Classes .....</b>	<b>48</b>
<b>10.6</b>	<b>Private Transfer Syntaxes.....</b>	<b>48</b>
	<b>INDEX OF TABLES.....</b>	<b>49</b>

# 3 Introduction

## 3.1 Revision History

Version	Date	Change
1.0	10/12/2023	First version of document valid from AI-Rad Companion Chest CT Pulmonary VA40A, Cardiovascular VA22A, Musculoskeletal VA22A
2.0	18/01/2024	Updated Lesion Review Status in Table 51 and adapted review comments.

## 3.2 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

## 3.3 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between AI-Rad Companion Chest CT and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [2]. DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of conformance statements is the first step towards assessing interconnectivity and interoperability between AI-Rad Companion Chest CT and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility. Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens Healthineers representative for the most recent product information.

## 3.4 Definitions, Terms and Abbreviations

Definitions, terms, and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

AE	DICOM Application Entity
AET	Application Entity Title
ASCII	American Standard Code for Information Interchange
DCS	DICOM Conformance Statement
DICOM	Digital Imaging and Communications in Medicine
FSC	File Set Creator
FSR	File Set Reader
FSU	File Set Updater

IOD	DICOM Information Object Definition
ISO	International Standard Organization
NEMA	National Electrical Manufacturers Association
O	Optional Key Attribute
PDU	DICOM Protocol Data Unit
R	Required Key Attribute
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM Server)
SOP	DICOM Service-Object Pair
SR	Structured Report
TID	Template ID
U	Unique Key Attribute
UID	Unique Identifier
UTF-8	Unicode Transformation Format-8
VR	Value Representation

### 3.5 References

- [1] DICOM Conformance Statements of teamplay DICOM Hub – <https://www.siemens-healthineers.com/en-in/services/it-standards/dicom-conformance-statements-digital-and-automation/teamplay>
- [2] NEMA PS3 / ISO 12052, Digital Imaging and Communications in Medicine (DICOM) Standard, National Electrical Manufacturers Association, Rosslyn, VA, USA (available free at <https://www.dicomstandard.org/>)
- [3] AI-Rad Companion Data Privacy and Security White Paper – Based on on-demand request from the end users.
- [4] Teamplay Data Privacy and Security White Paper – Based on on-demand request from the end users.
- [5] Integrating the Healthcare Enterprise – IHE Radiology Technical Framework – <http://www.ihe.net>

# 4 Networking

Please refer to the latest version of DICOM Conformance Statement of teamplay DICOM Hub [1] for further information on the provided networking capabilities for AI-Rad Companion Chest CT.

AI-Rad Companion Chest CT only supports subset of transfer syntaxes supported by teamplay and the following table lists the supported transfer syntaxes.

Table 4: Supported Transfer Syntaxes

Abstract Syntax		Transfer Syntax		Read	Create	Role	Extension Negotiation
SOP Classes	SOP Class UID	Name List	UID List				
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Explicit VR Little Endian	1.2.840.10008.1.2.1	No	Yes	SCP/SCU	None
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Explicit VR Little Endian	1.2.840.10008.1.2.1	No	Yes	SCP	None
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.8.33	Explicit VR Little Endian	1.2.840.10008.1.2.1	No	Yes	SCP	None

# 5 Media Interchange

AI-Rad Companion Chest CT is not providing any means for media interchange.



## **6 Transformations of DICOM to CDA**

NOT APPLICABLE

# 7 Support of Extended Character Sets

## 7.1 Character sets for AI-Rad Companion Chest CT

AI-Rad Companion Chest CT supports the following character sets as defined in the tables in this section.

Table 5: Single-Byte Character Sets without Code Extension

Character Set Description	Defined Term	ISO Registration Number	Character Set
Default repertoire	None	ISO_IR 6	ISO 646
Latin alphabet No. 1	ISO_IR 100	ISO_IR 100	Supplementary set
		ISO_IR 6	ISO 646
Latin alphabet No. 2	ISO_IR 101	ISO_IR 101	Supplementary set
		ISO_IR 6	ISO 646
Latin alphabet No. 3	ISO_IR 109	ISO_IR 109	Supplementary set
		ISO_IR 6	ISO 646
Latin alphabet No. 4	ISO_IR 110	ISO_IR 110	Supplementary set
		ISO_IR 6	ISO 646
Latin alphabet No. 5	ISO_IR 148	ISO_IR 148	Supplementary set
		ISO_IR 6	ISO 646
Cyrillic	ISO_IR 144	ISO_IR 6	Supplementary set
		ISO_IR 6	ISO 646
Arabic	ISO_IR 127	ISO_IR 127	Supplementary set
		ISO_IR 6	ISO 646
Greek	ISO_IR 126	ISO_IR 126	Supplementary set
		ISO_IR 6	ISO 646
Hebrew	ISO_IR 138	ISO_IR 138	Supplementary set
		ISO_IR 6	ISO 646
Japanese	ISO_IR 13	ISO_IR 13	JIS X 0201: Katakana
		ISO_IR 14	JIS X 0201: Romaji
Thai	ISO_IR 166	ISO_IR 166	TIS 620-253 (1990)
		ISO_IR 6	ISO 646

Table 6: Single-Byte Characters Sets with Code Extension

Character Set Description	Defined Term	Standard for Code Extension	ESC sequence	ISO Registration Number	Character Set
Default repertoire	ISO 2022 IR 6	ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646
Latin alphabet No.1	ISO 2022 IR 100	ISO 2022	ESC 02/13 04/01	ISO-IR 100	Supplementary set
		ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646
Latin alphabet No.2	ISO 2022 IR 101	ISO 2022	ESC 02/13 04/02	ISO-IR 101	Supplementary set
		ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646
Latin alphabet No.3	ISO 2022 IR 109	ISO 2022	ESC 02/13 04/03	ISO-IR 109	Supplementary set
		ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646
Latin alphabet No. 4	ISO 2022 IR 110	ISO 2022	ESC 02/13 04/04	ISO-IR 110	Supplementary set
		ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646
Latin alphabet No. 5	ISO 2022 IR 148	ISO 2022	ESC 02/13 04/13	ISO-IR 148	Supplementary set
		ISO 2022	ESC 02/08 04/02	ISO-IR 6	ISO 646

Table 7: Multi-Byte Character Sets without Code Extension

Character Set Description	Defined Term	ISO Registration Number	Character Set
Unicode	ISO_IR 192	ISO 10646	Unicode in UTF-8
Chinese	GB18030	GB18030	GB 18030-2023 (China Association for Standardization)

Table 8: Multi-Byte Character Sets with Code Extension

Character Set Description	Defined Term	Standard for Code Extension	ESC sequence	ISO Registration Number	Character Set
Japanese	ISO 2022 IR 159	ISO 2022	ESC 02/04 02/08 04/04	ISO-IR 159	JIS X 0212: Supplementary Kanji set
Korean	ISO 2022 IR 149	ISO 2022	ESC 02/04 02/09 04/03	ISO-IR 149	KS X 1001: Hangul and Hanja

All SCS (Special Character Sets) listed above are supported for incoming Data.

Three categories of character sets have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO\_IR 6, ISO 2022 IR 6, ISO\_IR 100, etc. → Encoded in ISO 2022
- ISO\_IR 192 → Encoded in UTF-8
- GB18030 → Encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO\_IR 192  $\leftarrow \rightarrow$  (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030  $\leftarrow \rightarrow$  (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022  $\leftarrow \rightarrow$  (0008,0005) contains ISO\_IR 192
- An attribute value is encoded in ISO 2022  $\leftarrow \rightarrow$  (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies does not conform to DICOM. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

AI-Rad Companion Chest CT supports Kanji characters in the byte zone after 74 (79, 7A, 7B and 7C).

# 8 Attribute confidentiality profiles

## 8.1 De-identification

Please refer to the AI-Rad Companion Data Privacy and Security White Paper [3] for further information on the support of de-identification of attributes natively for AI-Rad Companion Cloud deployment. In case of edge (on-premises) deployment, no de-identification happens from teamplay.

# 9 Security

## 9.1 Security Profiles

AI-Rad Companion Chest CT does not support any specific security measures.

## 9.2 Association Level Security

NOT APPLICABLE

## 9.3 Application Level Security

NOT APPLICABLE

# 10 Annexes

The tables in these section uses a number of abbreviations.

The abbreviations used in the “Presence” column are:

VNAP	Value is Not Always Present. Attribute is sent zero length if no value is present.
ANAP	Attribute Not Always Present.
ALWAYS	Attribute and Value are always present.
EMPTY	Attribute is sent zero length.

The abbreviations used in the “Source” Column are:

MWL	The attribute value is copied from Modality Worklist.
USER	The attribute value is entered by the user.
AUTO	The attribute value is generated by the system.
CONFIG	The attribute value is obtained by configuration

## 10.1 IOD Contents

### 10.1.1 Created SOP Instances

#### 10.1.1.1 Comprehensive SR

Table 9: Comprehensive SR IOD Modules

Information Entity	Module	Reference	Presence of Module
Patient	Patient	Table 10	ALWAYS
Study	General Study	Table 11	ALWAYS
	Patient Study	Table 12	ALWAYS
Series	SR Document Series	Table 13	ALWAYS
Equipment	General Equipment	Table 14	ALWAYS
SR Document	SR Document General	Table 15	ALWAYS
	SR Document Content	Table 16	ALWAYS
	SOP Common	Table 17	ALWAYS

Table 10 lists all Attributes that are supported in the Patient Module.

Table 10: Patient Module

Attribute	Tag	Source	Value	Presence	Comments
Patient`s Name	(0010,0010)	AUTO	Copied from source image	ALWAYS	
Patient ID	(0010,0020)	AUTO	Copied from source image	ALWAYS	
Patient`s Birth Date	(0010,0030)	AUTO	Copied from source image	VNAP	
Patient`s Sex	(0010,0040)	AUTO	Copied from source image	VNAP	
Other Patient IDs	(0010,1000)	AUTO	Copied from source image	ANAP	

Table 11 lists all Attributes that are supported in the General Study Module

Table 11: General Study Module

Attribute	Tag	Source	Value	Presence	Comments
Study Instance UID	(0020,000D)	AUTO	Copied from source image	ALWAYS	
Study Date	(0008,0020)	AUTO	Copied from source image	ALWAYS	
Study Time	(0008,0030)	AUTO	Copied from source image	VNAP	
Referring Physician's Name	(0008,0090)	AUTO	Copied from source image	VNAP	
Study ID	(0020,0010)	AUTO	Copied from source image	ALWAYS	
Accession Number	(0008,0050)	AUTO	Copied from source image	VNAP	
Study Description	(0008,1030)	AUTO	Copied from source image	ANAP	

Table 12 lists all Attributes that are supported in the Patient Study Module

Table 12: Patient Study Module

Attribute	Tag	Source	Value	Presence	Comments
Patient`s Age	(0010,1010)	AUTO	Copied from source image	VNAP	

Table 13 lists all Attributes that are supported in the SR Document Series Module

Table 13: SR Document Series Module

Attribute	Tag	Source	Value	Presence	Comments
Modality	(0008, 0060)	AUTO	SR	ALWAYS	
Series Instance UID	(0020, 000E)	AUTO	Unique identifier of the Series	ALWAYS	
Series Number	(0020, 0011)	AUTO	1999	ALWAYS	



Attribute	Tag	Source	Value	Presence	Comments
Series Date	(0008,0021)	AUTO	Date the Series started.	ALWAYS	
Series Time	(0008,0031)	AUTO	Time the Series started.	ALWAYS	
Protocol Name	(0018, 1030)	AUTO	Copied from source image	ALWAYS	
Series Description	(0018, 103E)	AUTO	Description of the SR series Example values AI-Rad Companion Structured Report	ALWAYS	
Referenced Performed Procedure Step Sequence	(0008,1111)	AUTO		EMPTY	

Table 14 lists all Attributes that are supported in the General Equipment Module

Table 14: General Equipment Module

Attribute	Tag	Source	Value	Presence	Comments
Manufacturer	(0008,0070)	AUTO	Siemens Healthineers	ALWAYS	
Manufacturer's Model Name	(0008,1090)	AUTO	AI-Rad Companion Chest CT	ALWAYS	
Device Serial Number	(0018,1000)	AUTO	9999	ANAP	
Software Versions	(0018,1020)	AUTO	Value can be one of the following: PULMO VAXXXISR-1.0 MSK VAXXXISR-1.0 CV VAXXXISR-1.0  where XXX indicates the used version of AI-Rad Companion Chest CT. For Example: - For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40AISR-1.0 For AI Rad Companion Chest CT CV version VA22A the value is CV VA22AISR-1.0 For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22AISR-1.0	ALWAYS	

Table 15 lists all Attributes that are supported in the SR Documents General Module

Table 15: SR Document General Module

Attribute	Tag	Source	Value	Presence	Comments
Instance Number	(0020,0013)	AUTO	Unique identifier	ALWAYS	
Verification Flag	(0040, A493)	AUTO	UNVERIFIED	ALWAYS	

Attribute	Tag	Source	Value	Presence	Comments
Completion Flag	(0040, A491)	AUTO	COMPLETE	ALWAYS	
Content Date	(0008,0023)	AUTO	The date the document content creation started.	ALWAYS	
Content Time	(0008,0033)	AUTO	The time the document content creation started.	ALWAYS	
Performed Procedure Code Sequence	(0040,A372)	AUTO		ALWAYS	
> Code Value	(0008,0100)	AUTO	24627-2	ALWAYS	
> Coding Scheme Designator	(0008,0102)	AUTO	LN	ALWAYS	LOINC code
> Coding Scheme Version	(0008,0103)	AUTO	CT Chest	ALWAYS	
Current Requested Procedure Evidence Sequence	(0040,A375)	AUTO		ALWAYS	
> Study Instance UID	(0020,000D)	AUTO	Study instance UID of the referenced study	ALWAYS	
> Referenced Series Sequence	(0008,1115)	AUTO		ALWAYS	
>> Series Instance UID	(0020,000E)	AUTO	Series instance UID of the referenced series	ALWAYS	
>> Referenced SOP Sequence	(0008,1199)	AUTO		ALWAYS	
>>> Referenced SOP Class UID	(0008,1150)	AUTO	SOP Class UID of the referenced Instance	ALWAYS	
>>> Referenced SOP Instance UID	(0008,1155)	AUTO	SOP Instance UID of the referenced Instance	ALWAYS	

Table 16 lists all Attributes that are supported in the SR Document Content Module

Table 16: SR Document Content Module

Attribute	Tag	Source	Value	Presence	Comments
Value Type	(0040,A040)	AUTO	CONTAINER	ALWAYS	
Concept Name Code Sequence	(0040,A043)	AUTO		ALWAYS	
> Code Value	(0008,0100)	AUTO	126000	ALWAYS	
> Coding Scheme Designator	(0008,0102)	AUTO	DCM	ALWAYS	
> Code Meaning	(0008,0104)	AUTO	Imaging Measurement Report	ALWAYS	
Content Sequence	(0040,A730)	AUTO		ALWAYS	

Refer to Table 51 for details regarding SR document content module.

Table 17 lists all Attributes that are supported in the SOP Common Module

Table 17: SOP Common Module

Attribute	Tag	Source	Value	Presence	Comments
Specific Character Set	(0008,0005)	AUTO	ISO_IR 192	ALWAYS	
Instance Creation Date	(0008,0012)	AUTO	Date when series is created	ALWAYS	
Instance Creation Time	(0008,0013)	AUTO	Time when series is created	ALWAYS	
SOP Class UID	(0008,0016)	AUTO	SOP Class UID of the respective DICOM Object	ALWAYS	
SOP Instance UID	(0008,0018)	AUTO	Unique Identifier of the instance	ALWAYS	
Original Specialized SOP Class UID	(0008,001B)	AUTO	1.3.12.2.1107.5.99.3.10	ANAP	
Timezone Offset from UTC	(0008,0201)	AUTO	Copied from the original image. If not present, then it is +0000	ALWAYS	
Contributing Equipment Sequence	(0018,A001)	AUTO		ALWAYS	
> Manufacturer	(0008,0070)	AUTO	The value can be either of the following: Siemens Healthineers SIEMENS	ALWAYS	
> Manufacturer Model Name	(0008,1090)	AUTO	AI-Rad Companion Chest CT	ALWAYS	
> Software Versions	(0018,1020)	AUTO	The value can be either of the following: MSK VAXXX CV VAXXX PULMO VAXXX	ALWAYS	

Attribute	Tag	Source	Value	Presence	Comments
			<p>Where XXX is the version number.</p> <p>For Example: -</p> <p>For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40A.</p> <p>For AI Rad Companion Chest CT CV version VA22A the value is CV VA22A.</p> <p>For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22A.</p>		
> Purpose Of Reference Code Sequence	(0040,A170)	AUTO		ALWAYS	
>> Code Value	(0008,0100)	AUTO	Newcode1	ALWAYS	
>> Coding Scheme Designator	(0008,0102)	AUTO	99SHSAIRC	ALWAYS	
>> Code Meaning	(0008,0104)	AUTO	Processing Algorithm	ALWAYS	
Instance Number	(0020,0013)	AUTO	A number that identifies this Composite object instance.	ALWAYS	

### 10.1.1.2 Secondary Capture Report

Table 18: Secondary Capture IOD Modules

Information Entity	Module	Reference	Presence of Module
Patient	Patient	Table 10	ALWAYS
Study	General Study	Table 11	ALWAYS
	Patient Study	Table 12	ALWAYS
Series	General Series	Table 19	ALWAYS
Equipment	General Equipment	Table 14	ALWAYS
	SC Equipment	Table 20	ALWAYS
Acquisition	General Acquisition	Table 24	ALWAYS
Image	General Image	Table 21	ALWAYS
	General Reference	Table 25	ALWAYS
	Image Pixel	Table 22	ALWAYS
	Device	Table 26	ALWAYS
	SC Image	Table 23	ALWAYS
	SOP Common	Table 17	ALWAYS

Table 19 lists all Attributes that are supported in the General Series Module

Table 19: General Series Modules

Attribute	Tag	Source	Value	Presence	Comments
Modality	(0008,0060)	AUTO	CT	ALWAYS	
Series Instance UID	(0020,000E)	AUTO	Unique identifier of the series.	ALWAYS	
Series Number	(0020,0011)	AUTO	AI-Rad Companion Results - 2000 AI-Rad Companion Pulmonary Lesions - 2010 AI-Rad Companion Pulmonary Lesions Current Time Point - 2010 AI-Rad Companion Pulmonary Lesions Prior Time Point - 2011 AI-Rad Companion Pulmonary Lesions Overview - 2019 AI-Rad Companion Pulmonary Parenchyma - 2020 AI-Rad Companion Pulmonary Parenchyma Overview - 2021 AI-Rad Companion Pulmonary Density - 2022 AI-Rad Companion Pulmonary Density Overview - 2023 AI-Rad Companion Musculoskeletal - 2030 AI-Rad Companion Musculoskeletal Overview - 2031 AI-Rad Companion Cardiovascular Aorta - 2040 AI-Rad Companion Cardiovascular Aorta Overview - 2041 AI-Rad Companion Cardiovascular Heart - 2042 AI-Rad Companion Cardiovascular Heart Overview - 2043	ALWAYS	

Attribute	Tag	Source	Value	Presence	Comments
			AI-Rad Companion Structured Report - 1999		
Series Date	(0008,0021)	AUTO	Date when series is created	ALWAYS	
Series Time	(0008,0031)	AUTO	Time when series is created	ALWAYS	
Series Description	(0008,103E)	AUTO	AI-Rad Companion Results	ALWAYS	
Patient Position	(0018,5100)	AUTO	Copied from input image	VNAP	

Table 20 lists all Attributes that are supported in the SC Equipment Module

Table 20: SC Equipment Modules

Attribute	Tag	Source	Value	Presence	Comments
Conversion Type	(0008,0064)	AUTO	WSD	ALWAYS	Synthetic Image
Secondary Capture Device ID	(0018,1010)	AUTO	AI-Rad Companion	VNAP	
Secondary Capture Device Manufacturer	(0018,1016)	AUTO	Siemens Healthineers	VNAP	
Secondary Capture Device Software Versions	(0018,1019)	AUTO	PULMO VAXXX, CV VAXXX, MSK VAXXX  VA XXX where XXX indicates the used version of AI-Rad Companion Chest CT.  For Example: - VA40A. For Example: -  For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40A.  For AI Rad Companion Chest CT CV version VA22A the value is CV VA22A.  For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22A.	VNAP	

Table 21 lists all Attributes that are supported in the General Image Module

Table 21: General Image Modules

Attribute	Tag	Source	Value	Presence	Comments
Instance Number	(0020,0013)	AUTO	A number that identifies this image. Which is a running number of the instance.	ALWAYS	
Patient Orientation	(0020,0020)	AUTO	Copied from input image	ANAP	

Attribute	Tag	Source	Value	Presence	Comments
Content Date	(0008, 0023)	AUTO	Date when series is created	ALWAYS	
Content Time	(0008,0033)	AUTO	Time when series is created	ALWAYS	
Image Type	(0008,0008)	AUTO	Value can be any of the following: In case of secondary capture series value will be: DERIVED\SECONDARY\0\RESULT In case of CT image storage series value will be: DERIVED\SECONDARY\AXIAL\MPR THICK In case of VRT Results and MPR series value will be: DERIVED\SECONDARY\0\CVRT	ALWAYS	
Image Comments	(0020,4000)	AUTO	auto-created	ALWAYS	

Table 22 lists all Attributes that are supported in the Image Pixel Module

Table 22: Image Pixel Modules

Attribute	Tag	Source	Value	Presence	Comments
Samples per Pixel	(0028,0002)	AUTO	3	ALWAYS	
Photometric Interpretation	(0028,0004)	AUTO	Copied from input image	ALWAYS	
Planar Configuration	(0028,0006)	AUTO	0	ALWAYS	
Rows	(0028,0010)	AUTO	1024	ALWAYS	
Columns	(0028,0011)	AUTO	1024	ALWAYS	
Bits Allocated	(0028,0100)	AUTO	8	ALWAYS	
Bits Stored	(0028,0101)	AUTO	8	ALWAYS	
High Bit	(0028,0102)	AUTO	7	ALWAYS	
Pixel Representation	(0028,0103)	AUTO	0	ALWAYS	
Pixel Data	(7FE0,0010)	AUTO	A data stream of the pixel samples that comprise the Image	ALWAYS	

Table 23 lists all Attributes that are supported in the SC Image Module

Table 23: SC Image Modules

Attribute	Tag	Source	Value	Presence	Comments
Date of Secondary Capture	(0018,1012)	AUTO	The date the Secondary Capture Image was captured.	ALWAYS	
Time of Secondary Capture	(0018,1014)	AUTO	The time the Secondary Capture Image was captured.	ALWAYS	
Pixel Spacing	(0028,0030)	AUTO	Copied from input image	ALWAYS	

Table 24 lists all Attributes that are supported in the General Acquisition Module

Table 24: General Acquisition Modules

Attribute	Tag	Source	Value	Presence	Comments
Acquisition Number	(0020,0012)	AUTO	Copied from input image	ALWAYS	
Acquisition Date	(0008,0022)	AUTO	Copied from input image	VNAP	
Acquisition Time	(0008,0032)	AUTO	Copied from input image	VNAP	

Table 25 lists all Attributes that are supported in the General Reference Module

Table 25: General Reference Modules

Attribute	Tag	Source	Value	Presence	Comments
Derivation Description	(0008,2111)	AUTO	Secondary Capture	ALWAYS	
Derivation Code Sequence	(0008,9215)	AUTO		ANAP	Root node
>> Code Value	(0008,0100)	AUTO	Identifier of the coded entry	ANAP	
>> Coding Scheme Designator	(0008,0102)	AUTO	DCM	ANAP	
>> Coding Scheme Version	(0008,0103)	AUTO	1.0	ANAP	
>> Code Meaning	(0008,0104)	AUTO	Volume rendering	ANAP	

Table 26 lists all Attributes that are supported in the Device Module

Table 26: Device Modules

Attribute	Tag	Source	Value	Presence	Comments
Manufacturer	(0008,0070)	AUTO	Siemens Healthineers		
Manufacturer's Model Name	(0008,1090)	AUTO	AI-Rad Companion Chest CT		
Device Serial Number	(0018,1000)	AUTO	9999		



Table 27 lists all Attributes that are supported in the General Equipment Module

Table 27: General Equipment Module

Attribute	Tag	Source	Value	Presence	Comments
Manufacturer	(0008,0070)	AUTO	Siemens Healthineers	ALWAYS	
Manufacturer's Model Name	(0008,1090)	AUTO	AI-Rad Companion Chest CT	ALWAYS	
Device Serial Number	(0018,1000)	AUTO	9999	ANAP	
Software Versions	(0018,1020)	AUTO	PULMO VAXXX, CV VAXXX, MSK VAXXX  where XXX indicates the used version of AI-Rad Companion Chest CT. For Example: -  For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40A.  For AI Rad Companion Chest CT CV version VA22A the value is CV VA22A.  For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22A.	ALWAYS	

### 10.1.1.3 Secondary Capture IOD (VRT Results & MPR Series (Secondary Capture)

Table 28: Secondary Capture IOD Modules (VRT Results & MPR Series)

Information Entity	Module	Reference	Presence of Module
Patient	Patient	Table 10	ALWAYS
Study	General Study	Table 11	ALWAYS
	Patient Study	Table 12	ALWAYS
Series	General Series	Table 29	ALWAYS
Equipment	General Equipment	Table 30	ALWAYS
	SC Equipment	Table 31	ALWAYS
Acquisition	General Acquisition	Table 24	ALWAYS
Image	General Image	Table 32	ALWAYS
	General Reference	Table 25	ALWAYS
	Image Pixel	Table 33	ALWAYS
	Device	Table 26	ALWAYS
	SC Image	Table 23	ALWAYS
	SOP Common	Table 17	ALWAYS

Table 29 lists all attributes that are supported in the General Series Module

Table 29: General Series Modules

Attribute	Tag	Source	Value	Presence	Comments
Modality	(0008,0060)	AUTO	CT	ALWAYS	
Series Instance UID	(0020,000E)	AUTO	Unique identifier of the series.	ALWAYS	
Series Number	(0020,0011)	AUTO	AI-Rad Companion Pulmonary Lesions - 2010 AI-Rad Companion Pulmonary Lesions Current Time Point - 2010 AI-Rad Companion Pulmonary Lesions Prior Time Point - 2011 AI-Rad Companion Pulmonary Lesions Overview - 2019 AI-Rad Companion Pulmonary Parenchyma - 2020 AI-Rad Companion Pulmonary Parenchyma Overview - 2021 AI-Rad Companion Pulmonary Density - 2022 AI-Rad Companion Pulmonary Density Overview - 2023 AI-Rad Companion Musculoskeletal - 2030 AI-Rad Companion Musculoskeletal Overview - 2031 AI-Rad Companion Cardiovascular Aorta - 2040 AI-Rad Companion Cardiovascular Aorta Overview - 2041 AI-Rad Companion Cardiovascular Heart - 2042 AI-Rad Companion Cardiovascular Heart Overview - 2043	ALWAYS	
Series Date	(0008,0021)	AUTO	Date when series is created	ALWAYS	
Series Time	(0008,0031)	AUTO	Time when series is created	ALWAYS	
Series Description	(0008,103E)	AUTO	Values can be one of the below: AI-Rad Companion Pulmonary Lesions Overview AI-Rad Companion Pulmonary Parenchyma Overview AI-Rad Companion Pulmonary Density Overview AI-Rad Companion Musculoskeletal Overview AI-Rad Companion Cardiovascular Aorta Overview AI-Rad Companion Cardiovascular Heart Overview AI-Rad Companion Pulmonary Lesions Current Time Point AI-Rad Companion Pulmonary Lesions Prior Time Point	ALWAYS	
Body Part Examined	(0018,0015)	AUTO	CHEST	ALWAYS	
Patient Position	(0018,5100)	AUTO	Copied from input image	VNAP	

Table 30 lists all Attributes that are supported in the General Equipment Module

Table 30: General Equipment Module

Attribute	Tag	Source	Value	Presence	Comments
Manufacturer	(0008,0070)	AUTO	SIEMENS	ALWAYS	

Attribute	Tag	Source	Value	Presence	Comments
Manufacturer's Model Name	(0008,1090)	AUTO	SOMATOM Definition AS	ALWAYS	
Device Serial Number	(0018,1000)	AUTO	9999	ALWAYS	
Software Versions	(0018,1020)	AUTO	<p>Value can be one of the following: PULMO VAXXX, MSK VAXXX, CV VAXXX</p> <p>where XXX indicates the used version of AI-Rad Companion Chest CT. For Example: -</p> <p>For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40A.</p> <p>For AI Rad Companion Chest CT CV version VA22A the value is CV VA22A.</p> <p>For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22A.</p>	ALWAYS	

Table 31 lists all Attributes that are supported in the SC Equipment Module

Table 31: SC Equipment Modules

Attribute	Tag	Source	Value	Presence	Comments
Conversion Type	(0008,0064)	AUTO	WSD	ALWAYS	Synthetic Image
Modality	(0008,0060)	AUTO	CT		
Secondary Capture Device ID	(0018,1010)	AUTO	AI-Rad Companion		
Secondary Capture Device Manufacturer	(0018,1016)	AUTO	Siemens Healthineers		
Secondary Capture Device Software Versions	(0018,1019)	AUTO	<p>Value can be one of the following: PULMO VAXXX, MSK VAXXX, CV VAXXX</p> <p>where XXX indicates the used version of AI-Rad Companion Chest CT. For Example: -</p> <p>For AI Rad Companion Chest CT Pulmo version VA40A the value is PULMO VA40A.</p>		

Attribute	Tag	Source	Value	Presence	Comments
			For AI Rad Companion Chest CT CV version VA22A the value is CV VA22A. For AI Rad Companion Chest CT MSK version VA22A the value is MSK VA22A.		

Table 32 lists all Attributes that are supported in the General Image Module

Table 32: General Image Modules

Attribute	Tag	Source	Value	Presence	Comments
Instance Number	(0020,0013)	AUTO	A number that identifies this image. Which is a running number of the instance.	ALWAYS	
Patient Orientation	(0020,0020)	AUTO	Copied from input image	ANAP	
Content Date	(0008, 0023)	AUTO	Date when series is created	ALWAYS	
Content Time	(0008,0033)	AUTO	Time when series is created	ALWAYS	
Image Type	(0008,0008)	AUTO	The below values is applicable for Pulmonary Lesions, Pulmonary Parenchyma, Pulmonary Density, Musculoskeletal, Cardiovascular Aorta and Cardiovascular Heart. DERIVED\SECONDARY\0\CVRT	ALWAYS	
Image Comments	(0020,4000)	AUTO	auto-created	ALWAYS	

Table 33 lists all Attributes that are supported in the Image Pixel Module

Table 33: Image pixel Modules

Attribute	Tag	Source	Value	Presence	Comments
Samples per Pixel	(0028,0002)	AUTO	Value can be: <ul style="list-style-type: none"> <li>3 – in case of burnt-in graphics</li> <li>1 – in case of overlay graphics</li> </ul>	ALWAYS	
Photometric Interpretation	(0028,0004)	AUTO	Copied from input image	ALWAYS	
Planar Configuration	(0028,0006)	AUTO	0	ALWAYS	
Rows	(0028,0010)	AUTO	512	ALWAYS	
Columns	(0028,0011)	AUTO	512	ALWAYS	
Bits Allocated	(0028,0100)	AUTO	Value can be:	ALWAYS	

Attribute	Tag	Source	Value	Presence	Comments
			<ul style="list-style-type: none"> <li>8 – in case of burnt-in graphics</li> <li>16 – in case of overlay graphics</li> </ul>		
Bits Stored	(0028,0101)	AUTO	Value can be: <ul style="list-style-type: none"> <li>8 – in case of burnt-in graphics</li> <li>16 – in case of overlay graphics</li> </ul>	ALWAYS	
High Bit	(0028,0102)	AUTO	Value can be: <ul style="list-style-type: none"> <li>7 – in case of burnt-in graphics</li> <li>15 – in case of overlay graphics</li> </ul>	ALWAYS	
Pixel Representation	(0028,0103)	AUTO	0	ALWAYS	
Pixel Data	(7FE0,0010)	AUTO	A data stream of the pixel samples that comprise the Image	ALWAYS	

#### 10.1.1.4 CT Image Storage

Table 34: CT Image Storage IOD Modules

Information Entity	Module	Reference	Presence of Module
Patient	Patient	Table 10	ALWAYS
Study	General Study	Table 11	ALWAYS
	Patient Study	Table 12	ALWAYS
Series	General Series	Table 29	ALWAYS
Frame of Reference	Frame of Reference	Table 35	ALWAYS
Equipment	General Equipment	Table 30	ALWAYS
Acquisition	General Acquisition	Table 24	ALWAYS
Image	General Image	Table 21	ALWAYS
	General Reference	Table 25	ALWAYS
	Image Plane	Table 36	ALWAYS
	Image Pixel	Table 33	ALWAYS
	CT Image	Table 37	ALWAYS
	SOP Common	Table 17	ALWAYS

Table 35 lists all the attributes of Frame of Reference Module

Table 35: Frame of Reference Module Attributes

Attribute	Tag	Source	Value	Presence	Comments
Frame of Reference UID	(0020,0052)	AUTO	Copied from source image	ALWAYS	
Position Reference Indicator	(0020,1040)	AUTO	Copied from source image	ALWAYS	

Table 36 lists all the attributes of Image Plane module

Table 36: Image Plane Module Attributes

Attribute	Tag	Source	Value	Presence	Comments
Pixel Spacing	(0028,0030)	AUTO	Copied from source image	ALWAYS	
Image Orientation (Patient)	(0020,0037)	AUTO	Copied from source image	ALWAYS	
Image Position (Patient)	(0020,0032)	AUTO	Copied from source image	ALWAYS	

Table 37 lists all the attributes of CT Image module

Table 37: CT Image Module Attributes

Attribute	Tag	Source	Value	Presence	Comments
Image Type	(0008,0008)	AUTO	The value can be any one of the below based on the output image type: DERIVED SECONDARY AXIAL MPR THICK	ALWAYS	
Samples per Pixel	(0028,0002)	AUTO	1	ALWAYS	
Photometric Interpretation	(0028,0004)	AUTO	MONOCHROME2	ALWAYS	
Bits Allocated	(0028,0100)	AUTO	16	ALWAYS	
Rescale Intercept	(0028,1052)	AUTO	-1024	ALWAYS	
Rescale Slope	(0028,1053)	AUTO	Copied from source image	ALWAYS	
Rescale Type	(0028,1054)	AUTO	Copied from source image	ANAP	
KVP	(0018,0060)	AUTO	120	VNAP	
Acquisition Number	(0020,0012)	AUTO	Copied from source image	ALWAYS	
Reconstruction Diameter	(0018,1100)	AUTO	Copied from source image	ANAP	

Attribute	Tag	Source	Value	Presence	Comments
Distance Source To Detector	(0018,1110)	AUTO	Copied from source image	ANAP	
Distance Source To Patient	(0018,1111)	AUTO	Copied from source image	ANAP	
Gantry Detector Tilt	(0018,1120)	AUTO	Copied from source image	ANAP	
Table Height	(0018,1130)	AUTO	Copied from source image	ANAP	
Rotation Direction	(0018,1140)	AUTO	Copied from source image	ANAP	
Exposure Time	(0018,1150)	AUTO	Copied from source image	ANAP	
Filter Type	(0018,1160)	AUTO	Copied from source image	ANAP	
Generator Power	(0018,1170)	AUTO	Copied from source image	ANAP	
Focal Spots	(0018,1190)	AUTO	Copied from source image	ANAP	
Convolution Kernel	(0018,1210)	AUTO	Copied from source image	ANAP	

### 10.1.2 Usage of Attributes from Received IODs

NOT APPLICABLE

### 10.1.3 Attribute Mapping

NOT APPLICABLE

### 10.1.4 Coerced/Modified Fields

NOT APPLICABLE

## 10.2 Data Dictionary of Private Attributes

NOT APPLICABLE

## 10.3 Coded Terminology and Templates

### 10.3.1 Context Groups

Table 38: CID 43. Numeric Value Failure Qualifier

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
DCM	114006	Measurement failure



Table 39: CID 270. Observer Type

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
DCM	121007	Device

Table 40: CID 7000. Diagnostic Imaging Report Document Title

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
LN	24627-2	CT Chest

Table 41: CID 6110. Lung Anatomy Finding or Feature

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	31094006	Lung lobes
DCM	121401	Derivation

Table 42: CID 9000. Physical Quantity Descriptor

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	363698007	Finding site

Table 43: CID 210. Qualitative Evaluation Modifier Type

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	106233006	Topographical modifier

Table 44: CID 6126. Location in Lung

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	72481006	Middle lobe of right lung

Table 45: CID 6129. Chest Site Involvement

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	39607008	Lung

Table 46: CID 7470. Linear Measurement

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	121207	Height

Table 47: CID 6141. Attenuation Coefficient Measurement

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
DCM	112031	Attenuation coefficient

Table 48: CID 219. Geometry Graphical Representation

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
DCM	111041	Contour

Table 49: CID 7470. Linear Measurements

Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)
SCT	103339001	Long Axis
SCT	103340004	Short Axis

### 10.3.2 Template Specifications

AI-Rad Companion Chest CT will generate the results in the form of Comprehensive DICOM SR represented in TID 1500 Measurement Report format. Please see the below tables for an overview of DICOM attributes and their values corresponding to this Measurement Report template.

Table 50: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
		CONTAINER	EV (126000, DCM, "Imaging Measurement Report")	Imaging Measurement Report	1500
>	HAS CONCEPT MOD	CODE	EV (113011, DCM, "Document Title Modifier")	("CHESTCT0999", "99SHSAIRC", "AI-Rad CT Lung Lesion"/"AI-Rad CT Lung Parenchyma"/"AI-Rad CT Cardio"/"AI-Rad CT Vascular Aorta"/"AI-RAD CT Spine"/"AI-Rad CT Pulmonary Density")	2010
>	HAS CONCEPT MOD	CODE	EV ("121049", "DCM", "Language of Content Item and Descendants")	("eng", "RFC5646", "English")	1204
>>	HAS CONCEPT MOD	CODE	EV ("121046", "DCM", "Country of Language")	("US", "ISO3166_1", "United States")	1204
>	HAS CONCEPT MOD	CODE	EV ("121058", "DCM", "Procedure Reported")	("24627-2", "LN", "CT Chest")	1500
>	CONTAINS	CONTAINER	EV ("111028", "DCM", "Image Library")		1600
>>	CONTAINS	CONTAINER	EV ("126200", "DCM", "Image Library Group")		1600
>>>	HAS ACQ CONTEXT	CODE	EV ("121139", "DCM", "Modality")	("CT", "DCM", "Computed Tomography")	1602
>>>	HAS ACQ CONTEXT	DATE	EV ("111060", "DCM", "Study Date")	Copied from input image	1602
>>>	HAS ACQ CONTEXT	TIME	EV ("111061", "DCM", "Study Time")	Copied from input image	1602
>	CONTAINS	CONTAINER	EV ("126010", "DCM", "Image Measurements")		1500
>>	CONTAINS	CONTAINER	EV ("125007", "DCM", "Image Measurement Group")	Refer to subsequent tables for details of values	1411

Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Lesion

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	<ol style="list-style-type: none"> <li>1. Lesion, If more than one lesions are found then identifier is named as L1, L2 etc.</li> <li>2. "No findings" if there are no lesions detected</li> </ol>	1411
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	Unique Identifiers	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	Different Codes and Meaning of findings: <ol style="list-style-type: none"> <li>1. ("RID50149", RADLEX, "Pulmonary nodule").</li> <li>2. ("C35552", NCIt, "Cardiovascular System Finding")</li> <li>3. ("C110937", NCIt, "Musculoskeletal Finding")</li> <li>4. ("RID28530", DCM, "Opacities")</li> <li>5. ("RID50149", RADLEX, "Pulmonary nodule")</li> <li>6. ("RID6039", RADLEX, "Low attenuation values in lung")</li> </ol>	1411
>>>	HAS CONCEPT MOD	CODE	EV ("363698007", "SCT", "Finding site")	Different Codes and Meaning of finding sites: <ol style="list-style-type: none"> <li>1. ("RID1302", RADLEX, "Right lung")</li> <li>2. ("RID1326", RADLEX, "Left lung")</li> <li>3. ("RID1327", RADLEX, "Upper lobe of left lung")</li> <li>4. ("RID1338", RADLEX, "Lower lobe of left lung")</li> <li>5. ("RID1303", RADLEX, "Upper lobe of right lung")</li> <li>6. ("RID1310", RADLEX, "Middle lobe of lung")</li> <li>7. ("RID1315", RADLEX, "Lower lobe of right lung")</li> <li>8. ("RID1301", RADLEX, "Lung")</li> <li>9. ("39607008", SCT, "Both lungs")</li> </ol>	1419
>>>	HAS OBS CONTEXT	TEXT	EV (CHESTCT0102,99SHSAIRC,"Lesion Review Status")	Measurement auto-confirmed/Measurement accepted/ To be measured.	1411

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers to result image SOP Instance UID	1411
>>>	CONTAINS	NUM	EV (103339001, SCT, "Long Axis")	UNITS = EV (mm,UCUM,milimeter)	1400
>>>>	INFERRED FROM	SCoord		Graphic data (0070, 0022) and Graphic type (0070,0023) attributes are added.	320
>>>	CONTAINS	NUM	EV (LOJK, IBSI, "Maximum 3D Diameter of a Mesh")	UNITS = EV (mm,UCUM,milimeter)	1400
>>>	CONTAINS	NUM	EV (103340004, SCT, "Short Axis")	UNITS = EV (mm,UCUM,milimeter)	1400
>>>>	INFERRED FROM	SCoord		Graphic data (0070, 0022) and Graphic type (0070,0023) attributes are added.	320
>>>	CONTAINS	NUM	EV (RID50155, RADLEX, "Mean 2D diameter")	UNITS = EV (mm,UCUM,milimeter)	1411
>>>	CONTAINS	NUM	EV (RID28668, RADLEX, "Volume")	UNITS = EV (mm3,UCUM,cubic milimeter)	1411
>>>	CONTAINS	Text	EV (121106, DCM, "Comment")	A warning indicating that slice thickness is outside optimum range would be added in case slice thickness of input data is >2.5 mm for US Institutes.	1410
>	CONTAINS	CONTAINER	EV (126011, DCM, Derived Imaging Measurements)		1500
>>	CONTAINS	NUM	EV (CHESTCT0103, 99SHSAIRC, Maximum 2D Diameter Change)	UNITS = EV (% , UCUM, Percent)	1411
>>	CONTAINS	NUM	EV (CHESTCT0104, 99SHSAIRC, Maximum 3D diameter Change)	UNITS = EV (% , UCUM, Percent)	1411
>>	CONTAINS	NUM	EV (CHESTCT0105, 99SHSAIRC, Maximum perpendicular 2D diameter Change)	UNITS = EV (% , UCUM, Percent)	1411
>>	CONTAINS	NUM	EV (CHESTCT0106, 99SHSAIRC, Mean 2D diameter Change)	UNITS = EV (% , UCUM, Percent)	1411
>>	HAS OBS CONTEXT	NUM	EV (CHESTCT0108, 99SHSAIRC, Volume Change)	UNITS = EV (d, UCUM, Day);  If volume change exceeds 999 days, then UNITS = EV(% , UCUM, Percent)	1411

Table 52: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Parenchyma

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	<ol style="list-style-type: none"> <li>LeftUpperLobe / LeftLowerLobe / RightUpperLobe / RightMiddleLobe / RightLowerLobe / LeftLung / RightLung / BothLungs / Undefined / Lung Applied Range</li> <li>"No findings" if there is no parenchyma results detected</li> </ol>	1411
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	<ol style="list-style-type: none"> <li>Unique identifier</li> </ol>	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	<p>Different Codes and Meaning of findings:</p> <ol style="list-style-type: none"> <li>("RID6039", RADLEX, "Low attenuation in lung") when parenchyma results were found</li> <li>("CHECTCT0006", 99SHSAIRC, "No parenchyma results available or all results rejected")</li> <li>("C98451", NCI, "Chronic Lung Disorder")</li> <li></li> </ol>	1411
>>>	HAS CONCEPT MOD	CODE	EV ("363698007", "SCT", "Finding site")	<p>Different Codes and Meaning of finding sites:</p> <ol style="list-style-type: none"> <li>("RID1302", RADLEX, "Right lung")</li> <li>("RID1326", RADLEX, "Left lung")</li> <li>("RID1327", RADLEX, "Upper lobe of left lung")</li> <li>("RID1338", RADLEX, "Lower lobe of left lung")</li> <li>("RID1303", RADLEX, "Upper lobe of right lung")</li> <li>("RID1310", RADLEX, "Middle lobe of lung")</li> <li>("RID1315", RADLEX, "Lower lobe of right lung")</li> <li>("RID1301", RADLEX, "Lung")</li> <li>("39607008", SCT, "Both lungs")</li> </ol>	1419
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers result image SOP Instance UID	1411
>>>	CONTAINS	CODE	EV (130400, DCM, "Geometric purpose of region")	(111041, DCM, "Contour")	1411

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	CONTAINS	NUM	EV (CHESTCT0201, 99SHSAIRC, "LAV950")	UNITS = EV (% ,UCUM,Percent)	1411
>>>	CONTAINS	CODE	EV (CHECTCT0001, 99SHSAIRC, "Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX , "Green") 2. ("RID39037", RADLEX , "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	
>>>	CONTAINS	CODE	EV (CHESTCT0202, 99SHSAIRC, "Lung Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX , "Green") 2. ("RID39037", RADLEX , "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411
>>>	CONTAINS	IMAGE	EV (121232, DCM, "Source series for segmentation")	Original image series	1411

Table 53: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Cardio

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	1. Heart/Calcium Score/Calcium Applied Range 2. "No findings" if there are no cardiac results detected	1411
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	Unique identifier	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	Different Codes and Meaning of findings: 1. ("C35552", NCI, "Cardio Vascular System Finding") 2. ("CHECTCT0006", 99SHSAIRC, "No cardiac results available or all results rejected")	1411
>>>	HAS CONCEPT MOD	CODE	EV ("363698007", "SCT", "Finding site")	Different Codes and Meaning of finding sites: 1. ("RID1385", "RADLEX", "Heart") 2. ("C12843", "NCIt", "Coronary Artery")	1419

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers result image SOP Instance UID	1411
>>>	CONTAINS	NUM	EV (CHESTCT0301, 99SHSAIRC, "Heart Volume")	UNITS = EV (mL,UCUM,milimeter )	1411
>>>	CONTAINS	NUM	EV (CHESTCT0302, 99SHSAIRC, "Coronary Calcium")	UNITS = EV (mm3, UCUM,cubic milimeter )	1411
>>>	CONTAINS	CODE	EV (CHECTCT0001, 99SHSAIRC, "Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411
>>>	CONTAINS	CODE	EV (CHESTCT0303, 99SHSAIRC, "Coronary Calcium Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411

Table 54: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Vascular

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	1. AorticSinus / Sinoturbularjunction / MidAscending / ProximalArch / MidArch / ProximalDescending / MidDescending / Diaphragm / Abdomina / Aorta Applied Range 2. "No Findings", If vascular results are not found	1411
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	Unique Identifiers	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	Different Codes and Meaning of findings: 1. ("C35552", NCI, "Cardio Vascular System Finding") 2. ("CHECTCT0006", 99SHSAIRC, "No aorta	1411



Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
				results available or all results rejected")	
>>>	HAS CONCEPT MOD	CODE	EV ("363698007", "SCT", "Finding site")	Different Codes and Meaning of finding sites: 1. ("C33557", NCIt, "Sinus of Valsalva") 2. ("RID579", RADLEX, "Sinotubular Junction") 3. ("CHESTCT0401", 99SHSAIRC, "Mid Ascending Aorta") 4. ("CHESTCT0402", 99SHSAIRC, "Proximal Aortic Arch") 5. ("CHESTCT0403", 99SHSAIRC, "Mid Aortic Arch") 6. ("CHESTCT0404", 99SHSAIRC, "Proximal Descending Thoracic Aorta") 7. ("CHESTCT0405", 99SHSAIRC, "Mid Descending Thoracic Aorta") 8. ("CHESTCT0406", 99SHSAIRC, "Aorta at Diaphragm") 9. ("RID905", RADLEX, "Abdominal Aorta") 10. ("RID480", RADLEX, "Aorta")	1419
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers to result image SOP Instance UID	1411
>>>	CONTAINS	NUM	EV (RID13432, RADLEX, "Diameter")	UNITS = EV (mm, UCUM, milimeter)	1411
>>>	CONTAINS	CODE	EV (CHECTCT0001, 99SHSAIRC, "Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411
>>>	CONTAINS	CODE	EV (CHESTCT0407, 99SHSAIRC, "Aorta Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411

Table 55: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Spine

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	<ol style="list-style-type: none"> <li>1. Thoracic vertebrae labels / Spine Applied Range</li> <li>2. "No Findings", if algorithm failed to determine spine labels</li> </ol>	1411
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	Unique identifiers	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	<p>Different Codes and Meaning of findings:</p> <ol style="list-style-type: none"> <li>1. ("C110937", NCIt, "Musculoskeletal Finding")</li> <li>2. ("CHECTCT0006", 99SHSAIRC, "No spine results available or all results rejected")</li> </ol>	1411
>>>	HAS CONCEPT MOD	CODE	EV ("363698007", "SCT", "Finding site")	<p>Different Codes and Meaning of finding sites:</p> <ol style="list-style-type: none"> <li>1. ("RID29198", "RADLEX", "First thoracic vertebra")</li> <li>2. ("RID29199", "RADLEX", "Second thoracic vertebra")</li> <li>3. ("RID29200", "RADLEX", "Third thoracic vertebra")</li> <li>4. ("RID29201", "RADLEX", "Fourth thoracic vertebra")</li> <li>5. ("RID29202", "RADLEX", "Fifth thoracic vertebra")</li> <li>6. ("RID29203", "RADLEX", "Sixth thoracic vertebra")</li> <li>7. ("RID29204", "RADLEX", "Seventh thoracic vertebra")</li> <li>8. ("RID31704", "RADLEX", "Eighth thoracic vertebra")</li> <li>9. ("RID29206", "RADLEX", "Ninth thoracic vertebra")</li> <li>10. ("RID29207", "RADLEX", "Tenth thoracic vertebra")</li> <li>11. ("RID29208", "RADLEX", "Eleventh thoracic vertebra")</li> <li>12. ("RID29209", "RADLEX", "Twelfth thoracic vertebra")</li> <li>13. ("RID29154", "RADLEX", "Vertebra")</li> </ol>	1419
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers result image SOP Instance UID	1411

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	CONTAINS	NUM	EV (121207, DCM, "Height")	UNITS = EV (mm,UCUM,milimeter)	1411
>>>>	HAS CONCEPT MOD	CODE	DT (106233006, SCT, "Topographical Modifier")	Different codes and meanings of modifiers in RADLEX: 1. ("RID5818", RADLEX, "Anterior") 2. ("RID5820", RADLEX, "Medial") 3. ("RID5819", RADLEX, "Posterior")	1419
>>>>	HAS CONCEPT MOD	CODE	EV (CHECTCT0001, 99SHSAIRC, "Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411
>>>	CONTAINS	NUM	EV (112031, DCM, "Attenuation coefficient")	UNITS = EV (hnsfU,UCUM,Hounsfield unit)	1411
>>>>	HAS CONCEPT MOD	CODE	EV (121401, DCM, "Derivation")	( C53319,NClt,"Mean") )	1411
>>>	CONTAINS	CODE	EV (CHECTCT0501, 99SHSAIRC, "Spine Range")	Different Codes and Meaning of ranges: 1. ("RID39089", RADLEX, "Green") 2. ("RID39037", RADLEX, "Yellow") 3. ("CHECTCT0002", 99SHSAIRC, "Red") 4. ("CHECTCT0003", 99SHSAIRC, "Orange")	1411

Table 56: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Pulmonary density

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
Refer Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT					
>>>	HAS OBS CONTEXT	TEXT	DT ("112039", "DCM", "Tracking Identifier")	1. LeftUpperLobe, LeftLowerLobe, RightUpperLobe, RightMiddleLobe, RightLowerLobe, LeftLung, RightLung, BothLungs, Lung 2. "Not found or Not confirmed", if lung opacities are not found 3. "No Findings", if algorithm fails to detect opacities in lung	1411

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	HAS OBS CONTEXT	UIDREF	EV ("112040", "DCM", "Tracking Unique Identifier")	Unique identifiers	1411
>>>	CONTAINS	CODE	EV ("121071", "DCM", "Finding")	Different Codes and Meaning of findings: 1. ("RAD28530", DCM, "Opacities") 2. ("CHECTCT0006", 99SHSAIRC, "Algorithm failed: no pulmonary density results available") 3. ("CHECTCT0006", 99SHSAIRC, "Lung opacities not found or not confirmed")	1411
>>>	HAS CONCEPT MOD	CODE	EV (31094006, SCT, "Lung lobes")	Different Codes and Meaning of lung lobes: 1. ("RID1302", RADLEX, "Right lung") 2. ("RID1326", RADLEX, "Left lung") 3. ("RID1327", RADLEX, "Upper lobe of left lung") 4. ("RID1338", RADLEX, "Lower lobe of left lung") 5. ("RID1303", RADLEX, "Upper lobe of right lung") 6. ("RID1310", RADLEX, "Middle lobe of lung") 7. ("RID1315", RADLEX, "Lower lobe of right lung") 8. ("RID1301", RADLEX, "Lung") 9. ("39607008", SCT, "Both lungs")	1411
>>>	CONTAINS	CODE	EV (130400, DCM, "Geometric purpose of region")	(111041, DCM, "Outline")	1411
>>>	CONTAINS	IMAGE	EV ("130401", "DCM", "Visual representation")	Refers to result image SOP Instance UID	1411
>>>	CONTAINS	IMAGE	EV (121232, DCM, "Source series for segmentation")	Original image series	1411
>>>	CONTAINS	NUM	EV (CHESTCT0601, 99SHSAIRC, "Opacity score")	UNITS = EV ({Number}, UCUM, Number)	1411
>>>	CONTAINS	NUM	EV (CHESTCT0602, 99SHSAIRC, "Total Volume")	UNITS = EV (mL, UCUM, milliliter)	1411
>>>	CONTAINS	NUM	EV (CHESTCT0603, 99SHSAIRC, "Opacity volume")	UNITS = EV (mL, UCUM, milliliter)	1411
>>>	CONTAINS	NUM	EV (CHESTCT0604, 99SHSAIRC, "Opacity percentage")	UNITS = EV (%{vol}, UCUM, VolumePercent)	1411
>>>	CONTAINS	NUM	EV (CHESTCT0605, 99SHSAIRC, "High opacity volume")	UNITS = EV (mL, UCUM, milliliter)	1411

Nesting Level	Relationship Type	Value Type	Concept Name	Values	TID
>>>	CONTAINS	NUM	EV (CHESTCT0606, 99SHSAIRC, "High opacity percentage")	UNITS = EV (%{vol},UCUM,VolumePercent)	1411
>>>	CONTAINS	NUM	EV (CHESTCT0607, 99SHSAIRC, "Mean HU total")	UNITS = EV (hnsfU,UCUM,Hoursfield unit)	1411

### 10.3.3 Private Code definitions

The following tables list all private attributes created by AI-Rad Companion Chest CT which may be included in the generated instances.

Table 57: Private Coded Entry Attributes

Tag	Attribute	Source	Value Type	Presence	Comments
(0029, 0010)	Private Creator	SIEMENS MEDCOM HEADER	LO	1	
(0029, 0011)	Private Creator	SIEMENS SYNGO ADVANCED PRESENTATION	LO	1	
(0029, 1041)	Application Header Type	VIA_NO_VOLUME	CS	1	
(0029, 1042)	Application Header ID	NOT FOR VOLUME WORKFLOW	LO	1	
(0029, 1043)	Application Header Version	V1 20120620	LO	1	

Table 58: Private Code definitions

Code Value	Code Meaning	Definition	Notes
CHECTCT0003	Orange		
CHECTCT0002	Red		
CHECTCT0001	Range		
Newcode1	Processing Algorithm		
CHECTCT0007	Undefined		
CHECTCT0004	Light Yellow		
CHECTCT0005	Unknown		
CHECTCT0006	Message description		
CHECTCT0501	Spine Range		
CHECTCT0502	AI-Rad CT Spine		
CHECTCT0611	AI-Rad CT Pulmonary Density		

Code Value	Code Meaning	Definition	Notes
CHESTCT0601	Opacity score		
CHESTCT0602	Total volume		
CHESTCT0603	Opacity volume		
CHESTCT0604	Opacity percentage		
CHESTCT0605	High opacity volume		
CHESTCT0606	High opacity percentage		
CHESTCT0607	Mean HU total		
CHESTCT0608	Mean HU of opacity		
CHESTCT0609	Standard deviation total		
CHESTCT0610	Standard deviation of opacity		
CHESTCT0201	LAV950		
CHESTCT0202	Lung Range		
CHESTCT0203	AI-Rad CT Lung Parenchyma		
CHESTCT0304	AI-Rad CT Cardio		
CHESTCT0301	Heart Volume		
CHESTCT0302	Coronary Calcium		
CHESTCT0303	Coronary Calcium Range		
CHESTCT0103	Maximum 2D diameter Change		
CHESTCT0104	Maximum 3D diameter Change		
CHESTCT0105	Maximum perpendicular 2D diameter Change		
CHESTCT0106	Mean 2D diameter Change		
CHESTCT0108	Volume Change		
CHESTCT0999	AI-Rad CT Lung Lesion		
CHESTCT0102	Lesion Review Status		
CHESTCT0401	Mid Ascending Aorta		
CHESTCT0402	Proximal Aortic Arch		
CHESTCT0403	Mid Aortic Arch		
CHESTCT0404	Proximal Descending Thoracic Aorta		
CHESTCT0405	Mid Descending Thoracic Aorta		
CHESTCT0406	Aorta at Diaphragm		
CHESTCT0408	Maximum of Ascending Aorta		
CHESTCT0409	Maximum of Descending Aorta		

Code Value	Code Meaning	Definition	Notes
CHESTCT0407	Aorta Range		
CHESTCT0410	AI-Rad CT Vascular Aorta		

The Coding Scheme Designator for all the above private codes is "99SHSAIRC".

10.4 Grayscale Image Consistency

NOT APPLICABLE

10.5 Standard Extended / Specialized / Private SOP Classes

NOT APPLICABLE

10.6 Private Transfer Syntaxes

NOT APPLICABLE



# Index of Tables

Table 1: Network Services.....	2
Table 2: Media Services .....	2
Table 3: Implementation Identifying Information .....	2
Table 4: Supported Transfer Syntaxes .....	7
Table 5: Single-Byte Character Sets without Code Extension.....	10
Table 6: Single-Byte Characters Sets with Code Extension .....	11
Table 7: Multi-Byte Character Sets without Code Extension .....	11
Table 8: Multi-Byte Character Sets with Code Extension .....	11
Table 9: Comprehensive SR IOD Modules .....	15
Table 10: Patient Module.....	16
Table 11: General Study Module .....	16
Table 12: Patient Study Module .....	16
Table 13: SR Document Series Module .....	16
Table 14: General Equipment Module .....	17
Table 15: SR Document General Module .....	17
Table 16: SR Document Content Module.....	19
Table 17: SOP Common Module .....	19
Table 18: Secondary Capture IOD Modules.....	21
Table 19: General Series Modules .....	21
Table 20: SC Equipment Modules.....	22
Table 21: General Image Modules.....	22
Table 22: Image Pixel Modules.....	23
Table 23: SC Image Modules.....	24
Table 24: General Acquisition Modules .....	24
Table 25: General Reference Modules .....	24
Table 26: Device Modules.....	24
Table 27: General Equipment Module .....	25
Table 28: Secondary Capture IOD Modules (VRT Results & MPR Series) .....	25
Table 29: General Series Modules .....	26
Table 30: General Equipment Module .....	26
Table 31: SC Equipment Modules.....	27
Table 32: General Image Modules.....	28
Table 33: Image pixel Modules .....	28
Table 34: CT Image Storage IOD Modules.....	30
Table 35: Frame of Reference Module Attributes.....	31
Table 36: Image Plane Module Attributes.....	31
Table 37: CT Image Module Attributes .....	31
Table 39: CID 43. Numeric Value Failure Qualifier .....	32
Table 40: CID 270. Observer Type .....	33
Table 41: CID 7000. Diagnostic Imaging Report Document Title .....	33
Table 42: CID 6110. Lung Anatomy Finding or Feature.....	33
Table 43: CID 9000. Physical Quantity Descriptor .....	33
Table 44: CID 210. Qualitative Evaluation Modifier Type .....	33
Table 45: CID 6126. Location in Lung.....	33
Table 46: CID 6129. Chest Site Involvement.....	34
Table 47: CID 7470. Linear Measurement.....	34
Table 48: CID 6141. Attenuation Coefficient Measurement .....	34
Table 49: CID 219. Geometry Graphical Representation .....	34
Table 50: CID 7470. Linear Measurements .....	34
Table 51: TID 1500 Measurement Report for Comprehensive DICOM SR – AI-Rad Companion Chest CT.....	35
Table 52: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Lesion .....	36
Table 53: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Parenchyma .....	37
Table 54: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Cardio .....	39
Table 55: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Vascular .....	40
Table 56: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Spine.....	42
Table 57: TID 1500 Measurement Report for Comprehensive DICOM SR – Chest CT Pulmonary density.....	43
Table 58: Private Coded Entry Attributes .....	45

Table 59: Private Code definitions .....45

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

---

Siemens Healthineers Headquarters  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com



Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany