



ACUSON Sequoia ultrasound system VA25 SW release

Virtual Touch 2D SWE

Quick Reference Card


siemens-healthineers.com/ultrasound



Using 2D shear wave elastography (SWE) to measure shear velocity (Vs) in meters per second (m/s) and elasticity (E) in kilopascal (kPa) | Small Parts Application

Step 1

- Select compatible transducer and exam preset¹



Step 2

- Press **VT** (Virtual Touch) located on control panel



Step 3

- **SWE** tab is default application on touch screen



Step 4

- Roll trackball to reposition the region of interest (ROI)
- Press right or left **Set** key and roll trackball to resize ROI
- Use minimal scanning pressure to prevent pre-compression effect
- Press **Update** on control panel to begin acquisition



¹ 2D SWE imaging is compatible with 5C1, DAX: General, Abdomen;
10L4: General, Thyroid, Breast, MSK; 18L6: General, Thyroid, Breast.

Step 5

- When acquisition is complete, system automatically freezes, emits one audible “beep” and activates Velocity or **Elasticity** display mode
- Adjust **Min** and/or **Max Velocity** or **Elasticity**
 - Increase range if measured value exceeds range by displaying “High” or,
 - Decrease range if better visualization of velocity of stiffness differences is needed
- Rotate **Shear Wave** to view **Quality** and **Displacement** displays
 - Measurements should be taken in regions with high quality (green to yellow-green) in the **Quality** display
 - Correlate the areas in the Velocity or Elasticity image with shear wave signal-to-noise ratio and **Displacement** amplitude

Note: During freeze, a cooling timer indicates remaining seconds until transducer is available for image acquisition. During a short cooling time (e.g., 1–2 seconds), the timer appears and then disappears very quickly.

Shear Wave	Min Velocity	Max Velocity
Velocity	0.5	6.5

Step 6

- If desired, select **Unit m/s** to display the 2D SWE image as an elasticity image in kPa
- If desired, select **Live Dual** to display the SWE image in full-screen format
- If desired, select **SWE** to hide / display 2D SWE or,
- Rotate **VT** on control panel to change transparency of 2D SWE

Unit
m/s

SWE

Live Dual

VT

Measurements

Shear Velocity and Elasticity

Measurements can be performed or viewed in the Velocity (or Elasticity), Quality or Displacement display modes.

- Press **Caliper** on the control panel
- **Shear Velocity** is default measurement tool
- Roll trackball to position measurement ROI. In **Live Dual**, the measurement ROI can be placed on either image.
- If needed, rotate **ROI Diameter** to resize
- Shear wave measurements and measurement ROI depth and diameter are displayed near image (Fig. 1)
- To display measurement in report, select measurement label
- If desired, press **Set** key and roll trackball to perform additional measurements
- Select **Delete** on the control panel to delete the measurement
- Select **Report** on left side of touch screen to view, store, print or transfer shear wave measurements

Note: If measurement ROI value exceeds the current maximum velocity or elasticity setting, the measured results display **HIGH** for the shear wave value. If measurement ROI value is lower than the current minimum velocity and elasticity setting, the measured results display **LOW**. The measured results display **X.XX** (for velocity) or **XX.X** (for elasticity) if the shear wave was not detected within the ROI or if a measurement is performed outside the ROI.





Fig. 1: Velocity display mode in Live Dual format and shear wave measurements of a breast lesion.

Shadow

Duplicates and displays a measurement on the adjacent image for a comparison of lesion size and/or location.

- In **Live Dual**, press **Caliper** on the control panel
- Select measurement tool
- **Shadow** is activated as a default on the touch screen
- Measure the right or left image
- The Shadow measurement is displayed in the measured results

Tip: Shadow may be used with the Shear Velocity tool when **Live Dual** is activated.



For the proper use of the software or hardware, please always use the Operator Manual or Instructions for Use (hereinafter collectively "Operator Manual") issued by Siemens Healthineers. This material is to be used as training material only and shall by no means substitute the Operator Manual. Any material used in this training will not be updated on a regular basis and does not necessarily reflect the latest version of the software and hardware available at the time of the training. The Operator Manual shall be used as your main reference, in particular for relevant safety information like warnings and cautions.

Note: Some functions shown in this material are optional and might not be part of your system.

Certain products, product related claims or functionalities described in the material (hereinafter collectively "Functionality") may not (yet) be commercially available in your country. Due to regulatory requirements, the future availability of said functionalities in any specific country is not guaranteed. Please contact your local Siemens Healthineers sales representative for the most current information. The reproduction, transmission or distribution of this training or its contents is not permitted without express written authority. Offenders will be liable for damages.

All names and data of patients, parameters and configuration dependent designations are fictional and examples only. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

ACUSON Sequoia and Virtual Touch are trademarks of Siemens Medical Solutions USA, Inc.

Siemens Healthineers Ultrasound owns the rights to all images.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey towards expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated 5 million patients globally everyday benefit from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 170 years of experience and 18,000 patents globally. With more than 48,000 dedicated colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029, USA
Phone: 1-888-826-9702
siemens-healthineers.com/ultrasound