

ADVIA Centaur Immunoassay System

High-Sensitivity Troponin I Assay



siemens-healthineers.com/TNIH











Intended Use

The ADVIA Centaur® High-Sensitivity Troponin I (TNIH) Assay is for in vitro diagnostic use in the quantitative measurement of cardiac troponin I in human serum or plasma using the ADVIA Centaur XP/XPT Immunoassay Systems. The assay can be used to aid in the diagnosis of acute myocardial infarction (AMI).

Assay Description

The Siemens Healthineers ADVIA Centaur TNIH Assay is a three-site sandwich immunoassay using direct chemiluminometric technology. The Solid Phase reagent is magnetic latex particles conjugated with streptavidin, with two bound biotinylated capture monoclonal antibodies each recognizing a unique cTnI epitope.

The Lite reagent comprises a conjugate whose architecture consists of a proprietary acridinium ester and a recombinant antihuman cTnI sheep Fab covalently attached to bovine serum albumin (BSA) for chemiluminescent detection. The accumulated light signal is directly related to the sample cTnI concentration.

Siemens Healthineers offers a proven, true high-sensitivity cardiac troponin I assay with fast, accurate, and actionable test results:

- Demonstrates clinical proof that the ADVIA Centaur High-Sensitivity Troponin I Assay is a true highsensitivity troponin I assay according to the IFCC guidelines.¹
- Offers fast, accurate, and actionable results that enable use of faster protocols.
- Provides confidence in results due to high tolerance to common interferents:

ADVIA Centaur TNIH

Hemolysis	Biotin (ng/mL)				
visible at naked eye	supplementation level				
<10% bias	<10% bias				
500 mg/dL Hb	3500 ng/mL				




ADVIA Centaur High-Sensitivity Troponin I Assay

Performance Summary

System	Sample Type	Sample Volume	Assay Range	LoB	LoD	LoQ (20% CV)	LoQ (10% CV)	Onboard Stability		99th percentile (n = 2010)
ADVIA Centaur XP/XPT	Human serum, plasma (lithium heparin)	100 μL	2.50– 25,000.00 pg/mL (ng/L)	0.50 pg/mL (ng/L)	1.60 pg/mL (ng/L)	2.50 pg/mL (ng/L)	6.0 pg/mL (ng/L)	28 days	18 min	Combined: 47.34 pg/mL (ng/L)* Male: 57.27 pg/mL (ng/L) Female: 36.99 pg/mL (ng/L)

^{*99}th percentile value determined using combined gender data and lithium heparin sample type.

ADVIA Centaur XP and XPT TNIH Assay Precision Curve



Ordering Information

System	SMN No.	Tests per Kit	Contents				
	10994774	100	1 ReadyPack® 1 vial High/Low Calibrator ADVIA Centaur TNIH Master Curve Card ADVIA Centaur TNIH Calibrator Assigned Value Card and bar-code labe				
ADVIA Centaur XP/XPT	10994775	500	 5 ReadyPacks 2 vials High/Low Calibrator ADVIA Centaur TNIH Master Curve Card ADVIA Centaur TNIH Calibrator Assigned Value Card and bar-code labels 				
	10994776		ADVIA Centaur TNIH Master Curve Material, 5 x 1.0 mL ADVIA Centaur MCM lot-specific value sheet				

ADVIA Centaur, ReadyPack, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone: +49 9131 84-0 siemens-healthineers.com

Legal Manufacturer

Siemens Healthcare Diagnostics Inc. Laboratory Diagnostics 511 Benedict Avenue Tarrytown, NY 10591-5005 USA

Phone: +1 914-631-8000

Reference:

1. Apple F, Sandoval Y, Jaffe A, Ordonez-Llanos J. Cardiac troponin assays: guide to understanding analytical characteristics and their impact on clinical care. Report of the IFCC Task Force on Clinical Applications of Cardiac Bio-Markers. Clinical Chemistry. 2016 Oct 10;63(1):73-81.