

## Study Protocol

# Transarterial thoracic chemotherapy using *syngo* Dyna3D HighSpeed

Interventional Oncology

The very short acquisition time of less than 3 s with *syngo* Dyna3D HighSpeed allows acquiring 3D datasets without breathing motion artifacts even in very sick patients.

---

### Courtesy of

Prof. Thomas J. Vogl,  
PD Dr. Stefan Zangos,  
Department of Radiology,  
University of Frankfurt,  
Germany

### Supported by

*syngo* Dyna3D HighSpeed

### System & Software

Artis zeego VC21  
*syngo* X Workplace VB21  
(Prototype software)

---

## Case Description

### Patient history

61-year-old female.  
Adenocarcinoma of the lung,  
recurrent tumor after surgery,  
radiation therapy and systemic  
chemotherapy. Actual third line  
therapy protocol.

### Diagnosis

No systemic metastases, local  
intrapulmonary infiltration.

### Treatment

Transarterial thoracic chemotherapy  
with a mix of Mitomycin,  
Gemcitabine and Cisplatin injected  
in the ascending aorta performed  
after using the *syngo* Dyna3D  
HighSpeed protocol.

### General comments

The very short acquisition time  
of less than 3 s makes it possible  
to acquire 3D datasets without  
breathing motion artifacts even  
in very sick patients.

*syngo* Dyna3D HighSpeed allowed  
performing the 3D acquisition with  
only 25 mL of contrast and helped  
to save nearly 30% contrast media  
compared to a regular 5 s protocol.  
(A regular 5 s protocol would have  
required a 7 s injection protocol,  
resulting in a total volume of 105 mL  
with 35 mL of contrast).

## Transarterial thoracic chemotherapy

Acquisition protocol	3sDR HighSpeed
<b>Injection protocol</b>	
Catheter position	Ascending aorta
Contrast medium (CM)	350 mg iodine/mL
Dilution (CM/Saline):	33 %
Injection volume	75 mL
Injection rate	15 mL/s
Duration of injection	5 s
X-ray delay	2 s
Power injector used	Yes
<b>Reconstructions</b>	
<b>Primary</b>	
Name	DynaCT Body Nat Fill HU
VOI size	Full
Slice matrix	512 × 512
Kernel type	HU
Image characteristics	Normal
Reconstruction mode	Nat fill
Viewing preset	DynaCT Body

## Clinical Images




Figure 1: MIP 6 mm

The presented images show very good and sharp delineation of pulmonary arteries and their small branches



Figure 2: Coronar MIP 29 mm

Visualization of intercostal arteries (white arrow)




Figure 3: Thick MIP

Visualization of pulmonary arteries with the tumor

---

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen, Germany

Phone: +49 9131 84-0

siemens-healthineers.com

*The statements by Siemens' customers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.*

*On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens sales organization worldwide.*

*All rights reserved.*