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Siemens Healthineers Hämostase-Workshop

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

The Key Players

How does the anti-PF4/heparin antibody cause HIT?

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


	Number of patients with TEC / number of evaluable patients (%)	Odds Ratio	95% con- fidence in- terval	p-value
Gender (406 patients)				
Male	80/167 (47.9%)	1.0		
Female	146/239 (61.1%)	1.08	0.63–1.87	0.77
Age (405 patients)				
0–39 years	16/34 (47.1%)	1.0		
40–59 years	57/103 (55.3%)	2.71	0.96–7.61	0.06
60–79 years	134/232 (57.8%)	1.94	0.76–4.96	0.16
≥ 80 years	19/36 (52.8%)	0.98	0.29–3.29	0.98
Underlying disease (400 patients)				
Internal medicine	62/143 (43.4%)	1.0		
General surgery	20/53 (37.7%)	0.72	0.34–1.51	0.38
Cardio-vascular Surgery	12/31 (38.7%)	0.62	0.23–1.68	0.35
Orthopedic and trauma surgery	98/125 (78.4%)	5.34	2.67–10.68	<0.001
Others	30/48 (62.5%)	1.70	0.76–3.84	0.20
Platelet count decrease (319 patients)				
≤ 30%	5/19 (26.3%)	1.0		
> 30–50%	9/29 (31.0%)	1.72	0.40–7.47	0.47
> 50–70%	27/57 (47.4%)	3.63	0.95–13.86	0.06
> 70–90%	101/158 (63.9%)	8.11	2.29–28.66	0.001
> 90%	38/56 (67.9%)	8.79	2.26–34.17	0.002
§ All predictors were simultaneously entered into the model. Hosmer-Lemeshow Chi-square=11.5 (8 df, p=0.18); Nagelkerke R Square=0.27.				

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1-3 = niedrig; 4-5 = mittel, 6-8 = hoch (*Warkentin, Curr Hematol Rep 2003*)

Diagnosis, Step 1: The 4Ts Score

4Ts Score²¹ = Pre-test scoring system to assess the probability of HIT
Calculated as a sum of points from 4 components:

① Thrombocytopenia

Points:

- 2** — Platelet fall > 50% *and* platelet nadir $\geq 20 \times 10^9/L$
- 1** — Platelet fall 30%-50% *or* Platelet nadir $10-19 \times 10^9/L$
- 0** — Platelet fall < 30% *or* Platelet nadir < $10 \times 10^9/L$


② Timing of platelet count fall

Points:

- 2** — Fall starts days 5-10 (timing is clear) *or* Fall starts ≤ 1 day + heparin exposure within 30 days
- 1** — Fall starts days 5-10 (timing is unclear) *or* Fall starts after day 10 *or* Fall starts ≤ 1 day + heparin exposure 30-100 days ago
- 0** — Fall starts ≤ 4 days without recent heparin exposure

Visualization of points:

Line color correlates with # points listed above

1 Calculate using peak platelet count *after* heparin started (i.e., not necessarily the day heparin started)

2 % decline is used, not an absolute platelet threshold, as HIT can occur without thrombocytopenia if the peak count is high (e.g., after surgery, myeloproliferative neoplasm).²¹

Total 4Ts Score
= probability of HIT²²

- ≥6 High (64%)**
- 4-5 Intermediate (14%)**
- ≤3 Low (0.2%)**


④ Other causes of thrombocytopenia

Points:

- 2** None apparent
 - 1** Possible
 - 0** Definite
- 1** This variable can be challenging to determine and leads to variation between providers in score determination.²³

3 **HEP score**²⁶ An alternate score with similar criteria but better performance with less experienced clinicians, likely due to inclusion of defined 'other causes of thrombocytopenia'.²⁷

1. Chronic thrombocytopenia
2. Newly initiated non-heparin med known to cause thrombocytopenia
3. Severe infection
4. Severe disseminated intravascular coagulation (DIC)
5. Indwelling intra-arterial device
6. Cardiopulmonary bypass within 96 hours

1 Calculated based on the first day of consistent platelet fall.

2 In many patients, platelet fall can be *biphasic*, with an initial decline due to a non-HIT insult (e.g. surgery), followed by a rebound and a second, HIT-induced decline.

3 It takes 5-10 days to produce new IgG antibodies. Immediate HIT onset can occur in those with recent exposure due to presence of previously formed antibodies.^{23,24}

③ Thrombosis or other sequelae

Points:

- 2** New thrombosis; skin necrosis; acute systemic reaction post-intravenous UFH bolus
- 1** Progressive or recurrent thrombosis; non-necrotizing (erythematous) skin lesions at heparin injection sites; suspected thrombosis (not proven)
- 0** None

1 Includes thrombosis that occurs or progresses in the period of platelet count decline *after* heparin exposure (i.e., not thrombosis that occurred before heparin exposure)

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Diagnosis, Step 2a: Laboratory diagnosis, *Immunoassays*

The laboratory diagnosis of HIT is classically divided into 2 steps:
(a) an **immunoassay** and (b) a **functional assay**.

Enzyme Immunoassay (EIA)

The EIA is the prototypical immunoassay, and is used as the standard immunoassay throughout this manuscript

Purpose: Determine if anti-PF4/heparin antibodies are *present*, and if so, their abundance

Performance: High sensitivity (0.97, 95% CI 0.95-0.99), limited specificity (0.82, 95% CI 0.90-0.84)²⁸

💡 Specificity varies based on assay used, increasing with IgG-specific assays^{28,29} and higher OD thresholds.³⁰

💡 Higher OD correlates with higher HIT likelihood³⁰ and thrombotic risk.³¹

How it works: (A) The plate well is coated with multimolecular complexes of PF4/heparin (or PF4/heparin-like polyanions). Patient serum or plasma is added and anti-PF4/heparin antibodies bind to the complex. (B) Tagged anti-human antibody is added and binds anti-PF4/heparin antibodies. A substrate is added to elicit color change. (C) Degree of color change (optical density, OD) increases with increasing concentration of anti-PF4/heparin antibodies.

*If EIA is positive, proceed with functional assay (see next page, "Laboratory diagnosis, Functional assays").
If EIA is negative, HIT is essentially excluded.*

Additional immunoassay options: **Rapid Immunoassays**

Multiple types of immunoassays exist, all with different performance characteristics and different international availability. Selected examples of rapid immunoassays are included below.

Purpose: Similar to EIA but provide results in ≤30 minutes (EIA is often batched, takes hours to perform)

Performance: Varies by assay

		Sensitivity (95% CI)	Specificity (95% CI)	How it works: Method of anti-PF4/heparin detection	Ab class detected
Particle gel immunoassay ³²	PaGIA	0.98 (0.94-1.00)	0.88 (0.91-0.96)	Polymers coated with PF4/heparin complexes	IgG1
IgG-specific chemiluminescence assay ³²	IgG CIA	0.95 (0.90-1.00)	0.94 (0.89-0.99)	Magnetic particles coated with PF4/polyvinyl sulfonate	IgG
Lateral flow immunoassay ³²	LFIA	0.97 (0.92-1.00)	0.91 (0.86-0.96)	PF4-polyanion complexes	IgG
Latex immunoturbidimetric assay ³³	LIA	0.974	0.94	Competitive inhibition of nanoparticles coated with HIT-like antibodies	IgG, IgA, IgM

There is limited guidance on how rapid immunoassays should be incorporated into diagnostic algorithms, and whether additional confirmatory laboratory testing is needed.³²⁻³³

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Diagnosis, Step 2b: Laboratory diagnosis, *Functional assays*


¹⁴C-serotonin release assay (SRA)

The SRA is the "gold standard" functional assay^{34,35} and is used as the standard functional assay throughout this manuscript

Purpose: Determine if anti-PF4/heparin antibodies are *pathogenic* (i.e. activate platelets)

Performance: High sensitivity (~0.95), high specificity (~0.95)³⁵

Due to variation in assay performance and reporting practices,³⁶ real-world performance may vary from literature reports.


How it works:^{34,35} (A) Patient serum is combined with washed donor platelets radiolabeled with ¹⁴C-serotonin. (B) At least two different doses of heparin are added. At therapeutic dose, the HIT antigen forms and is recognized by the anti-PF4/heparin antibody. At high dose, heparin disrupts the multimolecular HIT antigen^{34,39} so the anti-PF4/heparin antibody cannot bind. (C) When the antigen/antibody complex binds to the platelet Fc receptor, platelets activate and release granules containing ¹⁴C-serotonin. Radioactivity is measured and converted to a % ¹⁴C-serotonin release. (D) A diagnosis of HIT (i.e. positive SRA) is generally characterized by strong activation (high % serotonin release) with therapeutic dose heparin *plus* inhibited activation (low % serotonin release) with high dose heparin. A negative SRA shows no activation at either concentration. An indeterminate SRA shows activation at both heparin concentrations.

Additional functional assays

Multiple functional assays exist that differ by the method used for measuring platelet activation:


Heparin-induced platelet activation assay (HIPA)⁴⁰

Measurement of activation:
Visual assessment of decreased turbidity with platelet activation

P-selectin expression assay (PEA)⁴¹⁻⁴³

Measurement of activation:
Platelet P-selectin (CD62p) expression by flow cytometry

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LIMITATIONS

LEGAL NOTICE

The TORADI-HIT algorithm is a validated diagnostic instrument to calculate the probability of HIT in individual patients.

Platelet count (nadir) [$10^9/L$]

66

CRP (C-reactive protein; if available) [mg/L]

22

Type of immunoassay

Acustar HIT-IgG (CLIA) [U/mL]

Numeric immunoassay test result

3.6

For the Diamed ID-H/PF4 use Titer⁻¹. E.g. use 64 instead of 1:64

Did the patient receive unfractionated heparin during the current hospital stay?


Yes

Timing of the thrombocytopenia

Day 5-10, or day 1 if recent heparin exposure

Other causes of thrombocytopenia

Results

Cutoff

HIT less likely

HIT more likely

HIT highly probable

HIT highly improbable

96 %

Interpretation

The classification probability of HIT is 96 %. The probability of HIT is very high with values in the red area (90% of HIT-positive patients have values in the red area).

Personalized Clinical Diagnostics

von: www.toradi-hit.dbmr.unibe.ch

Diagnosis: Putting it all together

Why use risk scores before laboratory testing?

- Thrombocytopenia is common in hospitalized patients, and it is rarely due to HIT⁴⁴
- The 4Ts Score has excellent *negative* predictive value (can exclude HIT)²²
- EIA is sensitive, not specific. False positives are common.

★ Choosing Wisely recommendation to avoid laboratory testing in patients with low probability 4Ts Score^{45,46}

Why not wait for the laboratory results to change anticoagulation?

Without appropriate therapy in a patient with HIT, the daily risk of complications (combined risk of thromboembolism, amputation, and death) is 6.1%.⁴⁷


Risk of HIT can be further stratified based on EIA optical density (OD) and 4Ts Score:

EIA OD	4Ts Score		
	4-5	6-8	
	14%	64%	without EIA OD
<0.6	0%	0%	
0.6-1.49	5-25%	16%	68%
1.5-1.99	50%	53%	93%
2+	90%	92%	99%
			without 4Ts Score

Lighter shades indicate probability of the individual component (4Ts Score without EIA OD,²² EIA OD without 4Ts Score⁴⁸). Darker shades indicate probability accounting for combination of 4Ts Score and EIA OD.^{49,50} Exact probabilities will vary with different EIA assays.

★ In patients with a high probability 4Ts score and a strongly positive EIA, a functional assay may not be necessary.⁵¹

★ Rarely, functional assays are falsely negative.⁵² In patients with a high 4Ts score and high EIA OD, treatment for HIT can be considered despite a negative functional assay.⁵¹

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Selection of non-heparin anticoagulation & other therapies

Anticoagulation selection in Acute HIT

Dependent on 4 clinical criteria

		1) Clinically unstable	2) Life- or limb-threatening ischemia	3) Renal dysfunction CrCl <30mL/min	4) Liver dysfunction Bilirubin >1.5mg/dL
Oral	Oral Xa inhibitors*	Red	Red	Yellow	Yellow
	Vitamin K antagonist (VKA)	Red	Red	Red	Red
IV	Argatroban	Green	Green	Green	Red
	Bivalirudin	Green	Green	Red ^	Green
SQ	Danaparoid	Green	Green	Green	Green
	Fondaparinux*	Yellow	Green	Red ~	Green


Bili, bilirubin; CrCl, creatinine clearance; IV, intravenous; SQ, subcutaneous

- Green circle: Agent is preferred
- Yellow circle: Agent is not preferred, but can be considered based on availability, risk/benefit
- Red circle: Agent is not recommended


Laboratory monitoring required

* Existing data with rivaroxaban, apixaban
 ^ If argatroban not available, can use with close monitoring due to accumulation risk
 ~ Trivial risk of reported HIT,⁵⁴ but has been demonstrated to be safe for use in acute HIT^{55,56}
 ~ Use in renal dysfunction has been reported^{57,58}

Oral anticoagulation pearls^{51,59}

Additional treatment considerations



May JM.
 Res Pract Thromb Haemost 2023;
 7: e100283

Autoimmune HIT (aHIT)

What is autoimmune HIT?


- Rarely, the clinical manifestations of HIT can develop via **heparin-independent platelet activation**.
- Autoimmune HIT (aHIT) refers to any HIT-related platelet count fall that begins or persists in the absence of heparin.^{70,71}
- Six subtypes of aHIT have been described.
- Notably, spontaneous HIT has been reported to occur most commonly after knee replacement or infection.⁷¹

Pathophysiology

(1) Positively charged PF4 tetramers repel one another. (2) In aHIT, a non-heparin molecule overcomes this electrostatic repulsion, bringing individual PF4 tetramers together to form large multimolecular PF4 complexes. The exact non-heparin molecule able to perform this action is debated and may vary based on aHIT subtype, but proposed examples include HIT "super-antibodies" that bridge two PF4 tetramers and/or highly anionic molecules which could include chondroitin sulfate (released from joint cartilage in joint surgery), polyphosphate, and/or nucleotides such as DNA and RNA.^{70,71} (3) The PF4 complexes then serve as antigens bound by anti-PF4/polyanion IgG, which (4) trigger activation of platelets and other cell types.

Principles of diagnosis & management of aHIT syndromes are **the same as classical HIT**, with some unique considerations:

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