

## DICOM Conformance Statement

# *syngo*® MR XA50A

### Product Name

|                     |                   |
|---------------------|-------------------|
| MAGNETOM Avanto Fit | MAGNETOM Free.Max |
| MAGNETOM Free.Star  | MAGNETOM Lumina   |
| MAGNETOM Skyra Fit  | MAGNETOM Vida     |
| MAGNETOM Vida Fit   |                   |



# DICOM Conformance Statement Overview

The *syngo*® MR XA50A conforms to the DICOM Standard and supports the network services as described in Table 1: Network Services and the media services as described in Table 2: Media Services

**Table 1: Network Services**

| SOP Classes                                                | SOP Class UID                 | User of Service (SCU) |      | Provider of Service (SCP) |         |
|------------------------------------------------------------|-------------------------------|-----------------------|------|---------------------------|---------|
| Verification                                               |                               |                       |      |                           |         |
| Verification                                               | 1.2.840.10008.1.1             | Yes                   |      | Yes                       |         |
| SOP Classes created by syngo® MR XA50A                     |                               |                       |      |                           |         |
|                                                            |                               | Create                | Send | Store                     | Display |
| MR Image Storage                                           | 1.2.840.10008.5.1.4.1.1.4     | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Image Storage                                  | 1.2.840.10008.5.1.4.1.1.4.1   | Yes                   | Yes  | Yes                       | Yes     |
| MR Spectroscopy Storage                                    | 1.2.840.10008.5.1.4.1.1.4.2   | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Color Image Storage                            | 1.2.840.10008.5.1.4.1.1.4.3   | Yes                   | Yes  | Yes                       | Yes     |
| Secondary Capture Image Storage                            | 1.2.840.10008.5.1.4.1.1.7     | Yes                   | Yes  | Yes                       | Yes     |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3   | Yes                   | Yes  | Yes                       | Yes     |
| Grayscale Softcopy Presentation State Storage              | 1.2.840.10008.5.1.4.1.1.11.1  | Yes                   | Yes  | Yes                       | Yes     |
| Raw Data Storage                                           | 1.2.840.10008.5.1.4.1.1.66    | Yes                   | Yes  | Yes                       | No      |
| Spatial Registration Storage                               | 1.2.840.10008.5.1.4.1.1.66.1  | Yes                   | Yes  | Yes                       | Yes     |
| Segmentation Storage                                       | 1.2.840.10008.5.1.4.1.1.66.4  | Yes                   | Yes  | Yes                       | Yes     |
| Real World Value Mapping Storage                           | 1.2.840.10008.5.1.4.1.1.67    | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced SR Storage                                        | 1.2.840.10008.5.1.4.1.1.88.22 | Yes                   | Yes  | Yes                       | Yes     |
| Comprehensive SR Storage                                   | 1.2.840.10008.5.1.4.1.1.88.33 | Yes                   | Yes  | Yes                       | Yes     |
| Encapsulated PDF Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.104.1 | Yes                   | Yes  | Yes                       | Yes     |
| RT Image Storage                                           | 1.2.840.10008.5.1.4.1.1.481.1 | Yes                   | Yes  | Yes                       | Yes     |
| RT Dose Storage                                            | 1.2.840.10008.5.1.4.1.1.481.2 | Yes                   | Yes  | Yes                       | Yes     |
| RT Structure Set Storage                                   | 1.2.840.10008.5.1.4.1.1.481.3 | Yes                   | Yes  | Yes                       | Yes     |
| RT Plan Storage                                            | 1.2.840.10008.5.1.4.1.1.481.5 | Yes                   | Yes  | Yes                       | Yes     |

| SOP Classes                                                | SOP Class UID                 | User of Service (SCU) |      | Provider of Service (SCP) |         |
|------------------------------------------------------------|-------------------------------|-----------------------|------|---------------------------|---------|
| SOP Classes managed by <i>syngo</i> ® MR XA50A             |                               |                       |      |                           |         |
|                                                            |                               | Create                | Send | Store                     | Display |
| Computed Radiography Image Storage                         | 1.2.840.10008.5.1.4.1.1.1     | No                    | Yes  | Yes                       | Yes     |
| Digital X-Ray Image Storage<br>- For Presentation          | 1.2.840.10008.5.1.4.1.1.1.1   | No                    | Yes  | Yes                       | Yes     |
| Digital X-Ray Image Storage – For Processing               | 1.2.840.10008.5.1.4.1.1.1.1.1 | No                    | Yes  | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | No                    | Yes  | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage – For Processing   | 1.2.840.10008.5.1.4.1.1.1.2.1 | No                    | Yes  | Yes                       | Yes     |
| Digital Intra-Oral X-Ray Image - for Presentation - IMAGE  | 1.2.840.10008.5.1.4.1.1.1.3   | No                    | Yes  | Yes                       | Yes     |
| Digital Intra-Oral X-Ray Image - for Processing - IMAGE    | 1.2.840.10008.5.1.4.1.1.1.3.1 | No                    | Yes  | Yes                       | Yes     |
| CT Image Storage                                           | 1.2.840.10008.5.1.4.1.1.2     | No                    | Yes  | Yes                       | Yes     |
| Enhanced CT Image Storage                                  | 1.2.840.10008.5.1.4.1.1.2.1   | No                    | Yes  | Yes                       | Yes     |
| Ultrasound Multi-frame Image Storage                       | 1.2.840.10008.5.1.4.1.1.3.1   | No                    | Yes  | Yes                       | Yes     |
| MR Image Storage                                           | 1.2.840.10008.5.1.4.1.1.4     | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Image Storage                                  | 1.2.840.10008.5.1.4.1.1.4.1   | Yes                   | Yes  | Yes                       | Yes     |
| MR Spectroscopy Storage                                    | 1.2.840.10008.5.1.4.1.1.4.2   | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Color Image Storage                            | 1.2.840.10008.5.1.4.1.1.4.3   | Yes                   | Yes  | Yes                       | Yes     |
| Ultrasound Image Storage                                   | 1.2.840.10008.5.1.4.1.1.6.1   | No                    | Yes  | Yes                       | Yes     |
| Enhanced US Volume Storage                                 | 1.2.840.10008.5.1.4.1.1.6.2   | No                    | Yes  | Yes                       | Yes     |
| Secondary Capture Image Storage                            | 1.2.840.10008.5.1.4.1.1.7     | Yes                   | Yes  | Yes                       | Yes     |
| Multi-frame Single Bit Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.1   | No                    | Yes  | Yes                       | Yes     |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2   | No                    | Yes  | Yes                       | Yes     |

| SOP Classes                                                  | SOP Class UID                 | User of Service (SCU) |     | Provider of Service (SCP) |     |
|--------------------------------------------------------------|-------------------------------|-----------------------|-----|---------------------------|-----|
| Multi-frame Grayscale Word Secondary Capture Image Storage   | 1.2.840.10008.5.1.4.1.1.7.3   | Yes                   | Yes | Yes                       | Yes |
| Multi-frame True Color Secondary Capture Image Storage       | 1.2.840.10008.5.1.4.1.1.7.4   | No                    | Yes | Yes                       | Yes |
| 12-lead ECG Waveform Storage                                 | 1.2.840.10008.5.1.4.1.1.9.1.1 | No                    | Yes | Yes                       | No  |
| General ECG Waveform Storage                                 | 1.2.840.10008.5.1.4.1.1.9.1.2 | No                    | Yes | Yes                       | No  |
| Ambulatory ECG Waveform Storage                              | 1.2.840.10008.5.1.4.1.1.9.1.3 | No                    | Yes | Yes                       | No  |
| Hemodynamic Waveform Storage                                 | 1.2.840.10008.5.1.4.1.1.9.2.1 | No                    | Yes | Yes                       | No  |
| Cardiac Electrophysiology Waveform Storage                   | 1.2.840.10008.5.1.4.1.1.9.3.1 | No                    | Yes | Yes                       | No  |
| Basic Voice Audio Waveform Storage                           | 1.2.840.10008.5.1.4.1.1.9.4.1 | No                    | Yes | Yes                       | No  |
| General Audio Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.4.2 | No                    | Yes | Yes                       | No  |
| Arterial Pulse Waveform Storage                              | 1.2.840.10008.5.1.4.1.1.9.5.1 | No                    | Yes | Yes                       | No  |
| Respiratory Waveform Storage                                 | 1.2.840.10008.5.1.4.1.1.9.6.1 | No                    | Yes | Yes                       | No  |
| Grayscale Softcopy Presentation State Storage                | 1.2.840.10008.5.1.4.1.1.11.1  | Yes                   | Yes | Yes                       | Yes |
| Color Softcopy Presentation State Storage SOP Class          | 1.2.840.10008.5.1.4.1.1.11.2  | No                    | Yes | Yes                       | Yes |
| Pseudo-Color Softcopy Presentation State Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.11.3  | No                    | Yes | Yes                       | Yes |
| Blending Softcopy Presentation State Storage                 | 1.2.840.10008.5.1.4.1.1.11.4  | No                    | Yes | Yes                       | Yes |
| Grayscale Planar MPR Volumetric Presentation State Storage   | 1.2.840.10008.5.1.4.1.1.11.6  | No                    | Yes | Yes                       | No  |
| Compositing Planar MPR Volumetric Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.7  | No                    | Yes | Yes                       | No  |
| Advanced Blending Presentation State Storage                 | 1.2.840.10008.5.1.4.1.1.11.8  | No                    | Yes | Yes                       | No  |
| Volume Rendering Volumetric Presentation State Storage       | 1.2.840.10008.5.1.4.1.1.11.9  | No                    | Yes | Yes                       | No  |

| SOP Classes                                                      | SOP Class UID                  | User of Service (SCU) |     | Provider of Service (SCP) |     |
|------------------------------------------------------------------|--------------------------------|-----------------------|-----|---------------------------|-----|
| Segmented Volume Rendering Volumetric Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.10  | No                    | Yes | Yes                       | No  |
| Multiple Volume Rendering Volumetric Presentation State Storage  | 1.2.840.10008.5.1.4.1.1.11.11  | No                    | Yes | Yes                       | No  |
| X-Ray Angiographic Image Storage                                 | 1.2.840.10008.5.1.4.1.1.12.1   | No                    | Yes | Yes                       | Yes |
| Enhanced XA Image Storage                                        | 1.2.840.10008.5.1.4.1.1.12.1.1 | No                    | Yes | Yes                       | Yes |
| X-Ray Radiofluoroscopic Image Storage                            | 1.2.840.10008.5.1.4.1.1.12.2   | No                    | Yes | Yes                       | Yes |
| Enhanced XRF Image Storage                                       | 1.2.840.10008.5.1.4.1.1.12.2.1 | No                    | Yes | Yes                       | Yes |
| X-Ray 3D Angiographic Image Storage                              | 1.2.840.10008.5.1.4.1.1.13.1.1 | No                    | Yes | Yes                       | Yes |
| Breast Tomosynthesis Image Storage                               | 1.2.840.10008.5.1.4.1.1.13.1.3 | No                    | Yes | Yes                       | Yes |
| Nuclear Medicine Image Storage                                   | 1.2.840.10008.5.1.4.1.1.20     | No                    | Yes | Yes                       | Yes |
| Raw Data Storage                                                 | 1.2.840.10008.5.1.4.1.1.66     | Yes                   | Yes | Yes                       | No  |
| Spatial Registration Storage                                     | 1.2.840.10008.5.1.4.1.1.66.1   | Yes                   | Yes | Yes                       | Yes |
| Spatial Fiducials Storage                                        | 1.2.840.10008.5.1.4.1.1.66.2   | No                    | Yes | Yes                       | No  |
| Deformable Spatial Registration SOP Class                        | 1.2.840.10008.5.1.4.1.1.66.3   | No                    | Yes | Yes                       | No  |
| Segmentation Storage                                             | 1.2.840.10008.5.1.4.1.1.66.4   | Yes                   | Yes | Yes                       | Yes |
| Surface Segmentation Storage                                     | 1.2.840.10008.5.1.4.1.1.66.5   | No                    | Yes | Yes                       | Yes |
| Real World Value Mapping Storage                                 | 1.2.840.10008.5.1.4.1.1.67     | Yes                   | Yes | Yes                       | Yes |
| Basic Text SR Storage                                            | 1.2.840.10008.5.1.4.1.1.88.11  | No                    | Yes | Yes                       | Yes |
| Enhanced SR Storage                                              | 1.2.840.10008.5.1.4.1.1.88.22  | Yes                   | Yes | Yes                       | Yes |
| Comprehensive SR Storage                                         | 1.2.840.10008.5.1.4.1.1.88.33  | Yes                   | Yes | Yes                       | Yes |
| Procedure Log Storage                                            | 1.2.840.10008.5.1.4.1.1.88.40  | No                    | Yes | Yes                       | No  |
| Mammography CAD SR Storage                                       | 1.2.840.10008.5.1.4.1.1.88.50  | No                    | Yes | Yes                       | No  |
| Key Object Selection Document Storage                            | 1.2.840.10008.5.1.4.1.1.88.59  | No                    | Yes | Yes                       | No  |
| X-Ray Radiation Dose SR Storage                                  | 1.2.840.10008.5.1.4.1.1.88.67  | No                    | Yes | Yes                       | No  |
| Radiopharmaceutical Radiation Dose SR                            | 1.2.840.10008.5.1.4.1.1.88.68  | No                    | Yes | Yes                       | No  |

| SOP Classes                                     | SOP Class UID                 | User of Service (SCU) |     | Provider of Service (SCP) |     |
|-------------------------------------------------|-------------------------------|-----------------------|-----|---------------------------|-----|
| Encapsulated PDF Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.104.1 | Yes                   | Yes | Yes                       | Yes |
| Positron Emission Tomography Image Storage      | 1.2.840.10008.5.1.4.1.1.128   | No                    | Yes | Yes                       | Yes |
| RT Image Storage                                | 1.2.840.10008.5.1.4.1.1.481.1 | Yes                   | Yes | Yes                       | Yes |
| RT Dose Storage                                 | 1.2.840.10008.5.1.4.1.1.481.2 | Yes                   | Yes | Yes                       | Yes |
| RT Structure Set Storage                        | 1.2.840.10008.5.1.4.1.1.481.3 | Yes                   | Yes | Yes                       | Yes |
| RT Beams Treatment Record Storage               | 1.2.840.10008.5.1.4.1.1.481.4 | No                    | Yes | Yes                       | No  |
| RT Plan Storage                                 | 1.2.840.10008.5.1.4.1.1.481.5 | Yes                   | Yes | Yes                       | Yes |
| Radio Therapy Brachy Treatment Record Storage   | 1.2.840.10008.5.1.4.1.1.481.6 | No                    | Yes | Yes                       | No  |
| RT Treatment Summary Record Storage             | 1.2.840.10008.5.1.4.1.1.481.7 | No                    | Yes | Yes                       | No  |
| RT Ion Plan Storage                             | 1.2.840.10008.5.1.4.1.1.481.8 | No                    | Yes | Yes                       | No  |
| RT Ion Beams Treatment Record Storage           | 1.2.840.10008.5.1.4.1.1.481.9 | No                    | Yes | Yes                       | No  |
| Hanging Protocol Storage                        | 1.2.840.10008.5.1.4.38.1      | No                    | Yes | Yes                       | No  |
| <b>Transfer (Private SOP Class)</b>             |                               |                       |     |                           |     |
| Syngo Non-Image Storage                         | 1.3.12.2.1107.5.9.1           | Yes                   |     | Yes                       |     |
| <b>Storage Commitment</b>                       |                               |                       |     |                           |     |
| Storage Commitment Push Model SOP Class         | 1.2.840.10008.1.20.1          | Yes                   |     | No                        |     |
| <b>Worklist Management</b>                      |                               |                       |     |                           |     |
| Modality Worklist Information Model - FIND      | 1.2.840.10008.5.1.4.31        | Yes                   |     | No                        |     |
| Modality Performed Procedure Step SOP Class     | 1.2.840.10008.3.1.2.3.3       | Yes                   |     | No                        |     |
| <b>Query/Retrieve</b>                           |                               |                       |     |                           |     |
| Study Root Q/R - Information Model - FIND       | 1.2.840.10008.5.1.4.1.2.2.1   | Yes                   |     | Yes                       |     |
| Study Root Q/R - Information Model - MOVE       | 1.2.840.10008.5.1.4.1.2.2.2   | Yes                   |     | Yes                       |     |
| Patient/Study Only Q/R - Information Model MOVE | 1.2.840.10008.5.1.4.1.2.3.2   | No                    |     | No                        |     |
| <b>Print Management</b>                         |                               |                       |     |                           |     |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9         | Yes                   |     | No                        |     |

| SOP Classes                                 | SOP Class UID           | User of Service (SCU) | Provider of Service (SCP) |
|---------------------------------------------|-------------------------|-----------------------|---------------------------|
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18  | Yes                   | No                        |
| Basic Film Session SOP Class                | 1.2.840.10008.5.1.1.1   | Yes                   | No                        |
| Basic Film Box SOP Class                    | 1.2.840.10008.5.1.1.2   | Yes                   | No                        |
| Basic Grayscale Image Box SOP Class         | 1.2.840.10008.5.1.1.4   | Yes                   | No                        |
| Basic Color Image Box SOP Class             | 1.2.840.10008.5.1.1.4.1 | Yes                   | No                        |
| Printer SOP Class                           | 1.2.840.10008.5.1.1.16  | Yes                   | No                        |
| Print Job SOP Class                         | 1.2.840.10008.5.1.1.14  | Yes                   | No                        |
| Presentation LUT SOP Class                  | 1.2.840.10008.5.1.1.23  | Yes                   | No                        |

Table 2: Media Services

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| <b>Compact Disk - Recordable</b>  |                          |                  |
| STD-GEN-CD                        | Yes                      | Yes              |
| AUG-GEN-CD                        | Yes                      | Yes              |
| <b>DVD</b>                        |                          |                  |
| AUG-GEN-DVD                       | Yes                      | Yes              |
| AUG- GEN-DVD-J2K                  | Yes                      | Yes              |
| STD-GEN-DVD                       | Yes                      | Yes              |
| STD-GEN-DVD-J2K                   | Yes                      | Yes              |
| <b>USB</b>                        |                          |                  |
| AUG- GEN-USB-J2K                  | Yes                      | Yes              |
| STD-GEN-USB-J2K                   | Yes                      | Yes              |

Table 3: Implementation Identifying Information

| Name                        | Value                  |
|-----------------------------|------------------------|
| Application Context Name    | 1.2.840.100008.3.1.1.1 |
| Implementation Class UID    | 1.3.12.2.1107.5.2      |
| Implementation Version Name | "SYNGO_MR_XA50A"       |

|          |                                                            |           |
|----------|------------------------------------------------------------|-----------|
| <b>3</b> | <b>INTRODUCTION</b>                                        | <b>10</b> |
| 3.1      | Revision History                                           | 10        |
| 3.2      | Audience                                                   | 10        |
| 3.3      | Remarks                                                    | 10        |
| 3.4      | Definitions, Terms and Abbreviations                       | 11        |
| 3.5      | References                                                 | 11        |
| <b>4</b> | <b>NETWORKING</b>                                          | <b>12</b> |
| 4.1      | Implementation Model                                       | 12        |
| 4.1.1    | Application Data Flow                                      | 12        |
| 4.1.2    | Functional Definitions of Application Entities             | 14        |
| 4.1.3    | Sequencing of Activities                                   | 16        |
| 4.2      | AE Specifications                                          | 19        |
| 4.2.1    | Verification AE Specification                              | 19        |
| 4.2.2    | Storage AE Specification                                   | 20        |
| 4.2.3    | Storage Commitment AE Specification                        | 25        |
| 4.2.4    | Query/Retrieve AE Specification                            | 27        |
| 4.2.5    | Modality Worklist AE Specification                         | 31        |
| 4.2.6    | Modality Performed Procedure Step AE Specification         | 38        |
| 4.2.7    | Print AE Specification                                     | 41        |
| 4.3      | Network Interfaces                                         | 50        |
| 4.3.1    | Physical Network Interface                                 | 50        |
| 4.3.2    | Additional Protocols                                       | 50        |
| 4.3.3    | IPv4 and IPv6 Support                                      | 50        |
| 4.4      | Configuration                                              | 51        |
| 4.4.1    | AE Title/Presentation Address Mapping                      | 51        |
| 4.4.2    | Parameters                                                 | 52        |
| <b>5</b> | <b>MEDIA INTERCHANGE</b>                                   | <b>53</b> |
| 5.1      | Implementation Model                                       | 53        |
| 5.1.1    | Application Data Flow Diagram                              | 53        |
| 5.1.2    | Functional Definitions of AEs                              | 53        |
| 5.1.3    | Sequencing of Real-World Activities                        | 54        |
| 5.1.4    | File Meta Information for Implementation Class and Version | 54        |
| 5.2      | AE Specifications                                          | 55        |
| 5.2.1    | Media Storage AE – Specification                           | 55        |
| 5.3      | Augmented and Private Application Profiles                 | 56        |
| 5.3.1    | Augmented Application Profiles                             | 56        |
| 5.4      | Media Configuration                                        | 56        |
| <b>6</b> | <b>SUPPORT OF EXTENDED CHARACTER SETS</b>                  | <b>57</b> |

|                                                                  |            |
|------------------------------------------------------------------|------------|
| <b>7 ATTRIBUTE CONFIDENTIALITY PROFILES</b>                      | <b>60</b>  |
| <b>7.1 De-identification</b>                                     | <b>60</b>  |
| <b>8 SECURITY</b>                                                | <b>70</b>  |
| <b>8.1 Security Profiles</b>                                     | <b>70</b>  |
| 8.1.1 Time Synchronization Profiles                              | 70         |
| 8.1.2 Basic TLS Secure Transport Connection Profile              | 70         |
| <b>8.2 Association Level Security</b>                            | <b>70</b>  |
| <b>8.3 Application Level Security</b>                            | <b>70</b>  |
| <b>9 ANNEXES</b>                                                 | <b>71</b>  |
| <b>9.1 IOD Contents</b>                                          | <b>71</b>  |
| 9.1.1 Created SOP Classes                                        | 71         |
| 9.1.2 Usage of Attributes from received IODs                     | 106        |
| 9.1.3 Attribute mapping                                          | 106        |
| 9.1.4 Coerced / Modified fields                                  | 109        |
| <b>9.2 Data Dictionary of Private Attributes</b>                 | <b>109</b> |
| <b>9.3 Coded Terminology and Templates</b>                       | <b>143</b> |
| 9.3.1 Context Groups                                             | 143        |
| 9.3.2 Template Specifications                                    | 143        |
| 9.3.3 Private Code definitions                                   | 143        |
| <b>9.4 Grayscale Image Consistency</b>                           | <b>143</b> |
| <b>9.5 Standard Extended / Specialized / Private SOP Classes</b> | <b>143</b> |
| <b>9.6 Private Transfer Syntaxes</b>                             | <b>143</b> |
| <b>9.7 DICOM Print SCU – detailed status displays</b>            | <b>144</b> |
| 9.7.1 Common Status Information                                  | 144        |
| 9.7.2 Additional DICOM Execution Status Information              | 149        |
| <b>ANNEX A: INDEX OF TABLES</b>                                  | <b>150</b> |
| <b>ANNEX B: TABLE OF FIGURES</b>                                 | <b>151</b> |

## 3 Introduction

### 3.1 Revision History

| Version                 | Date       | Change                                                                                                                                                                                                                                                                                                                                                                                          |
|-------------------------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>syngo</i> ® MR XA50A | 2022-05-19 | <i>syngo</i> MR XA50A version (predecessor versions: <i>syngo</i> MR XA30A, <i>syngo</i> MR XA31A and <i>syngo</i> MR XA40A)<br><br>New system support:<br>MAGNETOM Avanto Fit, MAGNETOM Free.Max, MAGNETOM Free.Star, MAGNETOM Lumina, MAGNETOM Skyra Fit, MAGNETOM Vida, MAGNETOM Vida Fit<br><br>IOD Tables added for: Enhanced MR Image, Enhanced MR Color Image, MR Spectroscopy, Raw Data |

### 3.2 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

### 3.3 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between this product and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [1]. DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of conformance statements is the first step towards assessing interconnectivity and interoperability between *syngo*® MR XA50A and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens representative for the most recent product information.

### 3.4 Definitions, Terms and Abbreviations

Definitions, terms, and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

|       |                                                    |
|-------|----------------------------------------------------|
| AE    | DICOM Application Entity                           |
| AET   | Application Entity Title                           |
| ASCII | American Standard Code for Information Interchange |
| DCS   | DICOM Conformance Statement                        |
| DICOM | Digital Imaging and Communications in Medicine     |
| FSC   | File Set Creator                                   |
| FSR   | File Set Reader                                    |
| FSU   | File Set Updater                                   |
| GSDF  | Grayscale Standard Display Function                |
| IOD   | DICOM Information Object Definition                |
| ISO   | International Standard Organization                |
| n. a. | not applicable                                     |
| NEMA  | National Electrical Manufacturers Association      |
| O     | Optional Key Attribute                             |
| PDU   | DICOM Protocol Data Unit                           |
| R     | Required Key Attribute                             |
| SCU   | DICOM Service Class User (DICOM client)            |
| SCP   | DICOM Service Class Provider (DICOM Server)        |
| SOP   | DICOM Service-Object Pair                          |
| SR    | Structured Report                                  |
| TFT   | Thin Film Transistor (Display)                     |
| TID   | Template ID                                        |
| U     | Unique Key Attribute                               |
| UID   | Unique Identifier                                  |
| UTF-8 | Unicode Transformation Format-8                    |
| VR    | Value Representation                               |

### 3.5 References

- [1] Digital Imaging and Communications in Medicine (DICOM), <https://www.dicomstandard.org>
- [2] Integrating the Healthcare Enterprise – IHE Radiology Technical Framework – <http://www.ihe.net>

## 4 Networking

### 4.1 Implementation Model

*syngo*® MR XA50A supports storing DICOM images to remote nodes like workstations or Archiving Systems. Using the Storage Commitment Service it can request safe keeping of previously stored instances from an Archiving system. Additionally the *syngo*® MR XA50A can query remote notes, retrieve and store selected instances from that node. Using the Modality Worklist service the *syngo*® MR XA50A can query a HIS/RIS for scheduled procedures. Performed procedure status and other procedure data can be returned to the HIS/RIS using the Modality Performed Procedure Step (MPPS) Service. Furthermore printing of color and grayscale images is supported.

#### 4.1.1 Application Data Flow

The following figures provide a functional overview of the *syngo*® MR XA50A Application Entities (AE). Relationships are shown between user-invoked activities (in the circles at the left of the AEs) and the associated real-world activities provided by DICOM service providers (in the circles at the right of the AEs)

syngo® MR XA50A

Remote Services




Figure 1: syngo® MR XA50A DICOM Data Flow diagram – Acquisition Workflow




Figure 2: syngo® MR XA50A DICOM Data Flow diagram - Printing

#### 4.1.2 Functional Definitions of Application Entities

The SCP components of the Application Entities of the syngo® MR XA50A operate as background server processes. They exist as soon as the system is powered up and wait for association requests. Upon accepting an association with a negotiated Presentation Context they start to receive and process the requests described in the following sections.

The SCU components of the Application Entity are invoked upon requests from the user interface or indirectly by trigger from internal processes.

##### 4.1.2.1 Functional Definitions of Verification AE

The syngo® MR XA50A supports the Verification service as a SCP and SCU. As a SCU, Verification can be activated from the Administrator Portal during system configuration by sending a C-ECHO-Request.

As a SCP of the Verification Service the *syngo*® MR XA50A processes and responds to incoming verification requests using the C-ECHO-Response.

#### **4.1.2.2 Functional Definition of Storage of AE**

The *syngo*® MR XA50A Storage SCU is invoked either directly by the user, by an auto-archive trigger or internally by the Query/Retrieve Application Entity that is responsible for processing retrieve requests. The job consists of data describing the composite objects selected for storage and the destination Application Entity Title. An association is negotiated with the destination Application Entity and the image data is transferred using the C-STORE-Request. The transfer status is reported to the initiator of the Storage request.

The Storage SCP component of the *syngo*® MR XA50A starts to receive the Composite Objects and import them into the database after accepting an association with a negotiated Presentation Context. The system responses to the Storage Request immediately after reception of the Data.

#### **4.1.2.3 Functional Definition of the Storage Commitment AE**

If configured, the *syngo*® MR XA50A can serve as a SCU for the DICOM Storage Commitment service. Upon successful completion of a storage job, the system uses the N-ACTION-Request to request storage commitment from a remote DICOM Storage Commitment SCP. This can either be the same as the storage destination or a different system depending on the system configuration. Storage Commitment Requests are sent after a configurable delay after storing the objects. The *syngo*® MR XA50A can receive the N-EVENT-REPORT-Request on the same or a different association.

#### **4.1.2.4 Functional Specification of Query/Retrieve AE**

The *syngo*® MR XA50A supports DICOM Query/Retrieve as a SCU: The user can initiate a query to a remote node using the C-FIND-Request. After matching the specified keys, the remote Query /Retrieve SCP uses the C-FIND-Response to return the results of its search, which will be displayed to the user. Depending on user action the *syngo*® MR XA50A Query/Retrieve DICOM SCU sends a C-MOVE-Request to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query/Retrieve SCP) to the system's Storage SCP.

The *syngo*® MR XA50A supports the following query model:

- Study Root Query Model

Furthermore, the SCU services may issue relational queries, if supported by the remote Query/Retrieve SCP node and required by the querying Application.

The *syngo*® MR XA50A DICOM Query/Retrieve SCP accepts C-FIND Request, queries the local database based on the provided matching keys and returns the matches using the C-FIND-Response. Depending on further request from the remote Query/Retrieve SCU, the *syngo*® MR XA50A responds to C-MOVE-Requests by initiating a C-STORE sub-operation to send image objects to the Storage SCP of the querying system.

#### **4.1.2.5 Functional Definition of Modality Worklist AE**

The *syngo*® MR XA50A Modality Worklist SCU issues DICOM Modality Worklist requests using C-FIND-Requests. The results in the C-FIND-Response are stored in internal database. The provided Patient and Procedure information is used for patient registration prior to starting an exam.

#### **4.1.2.6 Functional Definition of Modality Performed Procedure Step SCU AE**

The *syngo*® MR XA50A MPPS SCU uses the N-CREATE-Request to inform an Information System that a procedure step is IN PROGRESS.

The *syngo*® MR XA50A MPPS SCU uses the N-SET-Request to inform the Information System about the finalization of the Procedure Step, using either a status of COMPLETED or DISCONTINUED.

#### 4.1.2.7 Functional Definition of Print AE

The Print SCU of the *syngo*® MR XA50A is invoked by the user interface to setup film-sheet layout and whenever an image is ready to be printed on film. The Print SCU will hold and maintain all data needed to compile a complete film-sheet from the data (images, layout, configuration) received. Whenever a film-sheet is ready to print the related data is used to supply the Information to the SOP Classes of the Print Management Service Class. A queue is maintained in order to intermediately store several film-sheets in case of resource problems on printer. The SCU will only supply and require the mandatory SOP Classes of the Print Management Service Class.

#### 4.1.3 Sequencing of Activities

This section describes the sequencing of Real-World Activities performed by the *syngo*® MR XA50A Entities using a UML sequence diagram. Real-World Activities are depicted as vertical bars and arrows show the events exchanged between them

##### 4.1.3.1 System Configuration




Figure 3: Sequence Diagram for Real World Activities - System Configuration

4.1.3.2 Acquisition Workflow




Figure 4: Sequence Diagram for Real World Activities -Acquisition workflow

## 4.1.3.3 Printing Workflow




Figure 5: Sequence Diagram for Real World Activities - Printing

## 4.2 AE Specifications

This section outlines the specifications for each of the Application Entities that are part of the *syngo*® MR XA50A.

### 4.2.1 Verification AE Specification

#### 4.2.1.1 SOP Classes

The Verification AE of the *syngo*® MR XA50A provides standard conformance to the Verification SOP Class listed in Table 1: Network Services section "Verification" in the "[Conformance Statement Overview](#)".

#### 4.2.1.2 Association Policy

The *syngo*® MR XA50A Admin Portal attempts to open an association for verification request whenever the Verification function is activated.

**Table 4: Association Policies**

|                                                                         |                       |
|-------------------------------------------------------------------------|-----------------------|
| Application Context Name                                                | 1.2.840.10008.3.1.1.1 |
| PDU size                                                                | 32 kB <sup>1</sup>    |
| Maximum number of simultaneous associations as an association acceptor  | 12 <sup>1</sup>       |
| Maximum number of simultaneous associations as an association initiator | unlimited             |

##### 4.2.1.2.1 Asynchronous Nature

The *syngo*® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any size is supported.

**Table 5: Asynchronous Nature as an Association Initiator**

|                                                         |    |
|---------------------------------------------------------|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|---------------------------------------------------------|----|

##### 4.2.1.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to "**

Table 3: Implementation Identifying Information" in the "[Conformance Statement Overview](#)".

#### 4.2.1.3 Association Initiation Policy

##### 4.2.1.3.1 Activity – "Send Activity"

##### 4.2.1.3.1.1 Description and Sequencing of Activities

The *syngo*® MR XA50A serves as a SCU of the Verification Service Class. A C-ECHO-Request is initiated by the Administrator Portal whenever "Verification" is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than "Success" this will be indicated to the user and the association is closed.

<sup>1</sup> Default, the value is configurable

#### 4.2.1.3.1.2 Proposed Presentation Contexts

Table 6: Presentation Context Table "Verification" below lists the supported presentation contexts for verification requests.

Table 6: Presentation Context Table "Verification"

| Abstract Syntax |                   | Transfer Syntax           |                   | Role | Extended Negotiation |
|-----------------|-------------------|---------------------------|-------------------|------|----------------------|
| Name            | UID               | Name List                 | UID List          |      |                      |
| Verification    | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

#### 4.2.1.3.1.3 SOP specific Conformance for SOP classes

The ECHO-SCU provides standard conformance to the Verification Service Class.

#### 4.2.1.4 Association Acceptance Policy

##### 4.2.1.4.1 Activity – "Receive Verification Request"

##### 4.2.1.4.1.1 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCP of the Verification Service Class. If the Verification SCP accepts an association, it will respond to C-ECHO-Requests. If the Called AE Title does not match any pre-configured AE Title shared by SCP, the association will be rejected.

##### 4.2.1.4.1.2 Accepted Presentation Contexts

The syngo® MR XA50A DICOM application will accept Presentation Contexts as shown in the following table:

Table 7: Presentation Context Table "Verification"

| Abstract Syntax |                   | Transfer Syntax           |                   | Role | Extended Negotiation |
|-----------------|-------------------|---------------------------|-------------------|------|----------------------|
| Name            | UID               | Name List                 | UID List          |      |                      |
| Verification    | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP  | None                 |

##### 4.2.1.4.1.3 SOP specific Conformance – Verification SCP

The ECHO-SCP provides standard conformance to the Verification Service Class.

#### 4.2.2 Storage AE Specification

##### 4.2.2.1 SOP Classes

The Storage AE provides Standard Conformance to the SOP Classes listed in "Table 1: Network Services" section "SOP Classes Created by the syngo® MR XA50A" and "SOP Classes Managed by the syngo® MR XA50A" in the ["Conformance Statement Overview"](#).

#### 4.2.2.2 Association Policy

**Table 8: Association Policy**

|                                                                                |                       |
|--------------------------------------------------------------------------------|-----------------------|
| <b>Application Context Name</b>                                                | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>                                                                | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

##### 4.2.2.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 9: Asynchronous Nature as an Association Initiator**

|                                                                |    |
|----------------------------------------------------------------|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|----------------------------------------------------------------|----|

##### 4.2.2.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to "**

Table 3: Implementation Identifying Information" in the ["Conformance Statement Overview"](#).

#### 4.2.2.3 Association Initiation Policy

##### 4.2.2.3.1 Activity- "Send Storage Request"

##### 4.2.2.3.1.1 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCU of the Storage Service Class. The Storage SCU is triggered by the transfer job queue or by an external retrieve request. An association request is sent to the destination AE. Upon successful negotiation of a Presentation Context, the transfer is started. Objects will be transferred sequentially on the same open association.

##### 4.2.2.3.1.2 Proposed Presentation Contexts

For all Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with "yes" in the Image Objects Column of the table below are supported.

For all Non-Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with "yes" in the Non-Image Objects Column of the table below are supported.

<sup>1</sup> Default, the value is configurable

For a distinction between Image and Non-Image Objects please refer to the DICOM Standard PS3.3 Section A.1.4 "Overview of the Composite IOD Module Content [1]".

**Table 10: Proposed Presentation Contexts for Storage**

| UID value              | Transfer Syntax                                                                          | Image Objects | Non-Image Objects |
|------------------------|------------------------------------------------------------------------------------------|---------------|-------------------|
| 1.2.840.10008.1.2      | Implicit VR Little Endian                                                                | Yes           | Yes               |
| 1.2.840.10008.1.2.1    | Explicit VR Little Endian                                                                | Yes           | Yes               |
| 1.2.840.10008.1.2.2    | Explicit VR Big Endian                                                                   | Yes           | Yes               |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed                                               | Yes           | No                |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed                                           | Yes           | No                |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | Yes           | No                |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed                                   | Yes           | No                |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed                                             | Yes           | No                |
| 1.2.840.10008.1.2.5    | RLE Lossless compressed                                                                  | Yes           | No                |

Depending on the configuration, the Storage SCU will choose a compressed or uncompressed Transfer Syntax among those accepted by the SCP. The Transfer Syntax chosen is the preferred one among the compressed and uncompressed ones. The preference order is the order of occurrence in the configuration. For each node it is possible to select Transfer Syntaxes, that can be used, and Transfer Syntaxes, that are excluded. The configuration can even be extended, based on the combination of SOP Classes and supported Transfer Syntaxes.

An instance will be JPEG lossless compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is either MONOCHROME1, MONOCHROME2, RGB, YBR\_FULL or YBR\_FULL\_422
- Bits Allocated (0028,0100) equal to '16' or '8'
- Bits Stored (0028,0101) equal to '12' or '8'
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to '0'

An instance will be JPEG lossy compressed during transfer only if the following criteria is fulfilled:

- Is an image
- Photometric Interpretation (0028,0004) is either MONOCHROME1, MONOCHROME2 or RGB
- Bits Allocated (0028,0100) equal to '16' or '8'
- Bits Stored (0028,0101) equal to '12' or '8'
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to '0'
- Only lossy transfer syntaxes are supported (Implicit VR Little Endian is not supported) at the remote side

An instance will be JPEG 2000 lossless compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is either MONOCHROME, RGB, YBR\_FULL or YBR\_FULL\_422
- Bits Allocated (0028,0100) equal to '16' or '8'

An instance will be RLE lossless compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME, RGB, YBR\_FULL or RGB
- Bits Allocated (0028,0100) neither '16' nor '8'

An instance will be JPEG 2000 lossy compressed during transfer only if the following criteria is fulfilled:

- Is an image
- Photometric Interpretation (0028,0004) is MONOCHROME or RGB
- Bits Stored (0028,0101) equal to '12' or '8'
- Only lossy transfer syntaxes are supported (Implicit VR Little Endian is not supported) at the remote side

There is no extended negotiation as an SCU.

#### 4.2.2.3.1.3 SOP Specific Conformance for SOP Classes

The syngo® MR XA50A does not add or change private attributes by default, even in case objects are compressed or the image header is updated according to IHE [2] Patient Information Reconciliation Profile.

The behavior of syngo® MR XA50A when encountering status codes in a C-STORE response is summarized in Table 11

**Table 11: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning              | Error Code         | Behavior                                                                         |
|----------------|------------------------------|--------------------|----------------------------------------------------------------------------------|
| Error          | Any other DIMSE Error Status | Any none null Code | Send is continued till the end. Log message is created.                          |
| Success        | Image is successfully stored | 0000               | If configured, Storage Commitment is requested for successfully stored instances |

Table 12 below indicates the behavior if exceptions occur:

**Table 12: DICOM Command Communication Failure Behavior**

| Exception           | Behavior                                                   |
|---------------------|------------------------------------------------------------|
| Timeout             | Log message is created (Timeout configurable; default 30s) |
| Association Aborted | Send is failed. Log message is created.                    |

#### 4.2.2.3.1.4 Correction and Rearrangement

When a Study is moved to a different:

- Procedure received through a DICOM Modality Worklist, the Study Instance UID is overwritten with the Study Instance UID and Accession Number of the Procedure.
- Patient, the system generates a new Study Instance UID.

The system will not update references to the changed Study Instance UIDs, therefore it is possible that there will be broken links between Studies after such move operations.

In case of Patient Merge and Correction no UIDs are changed, therefore it is advised to delete any corrected or rearranged objects from the PACS before attempting to archive them again, to ensure that the PACS system can store them successfully.

When the Patient Position (0018,5100) attribute is corrected, the following attributes are recalculated by the system (no UIDs are changed):

1. Image Position (0020,0032)
2. Image Orientation (0020,0037)
3. Patient Orientation (0020,0020)
4. Data Collection Center (Patient) (0018,9313) (CT only)
5. Reconstruction Target Center (Patient) (0018,9318) (CT only)
6. Positioner Primary Angle (0018,1510) (XA only)
7. Positioner Secondary Angle (0018,1511) (XA only)

Also, the value of the Slice Location (0020,1041) attribute is emptied and a new Frame of Reference UID (0020,0052) is generated for the corrected series.

When the Patient Birth Date or the Study Date is corrected, the system recalculates the Patient Age.

A new item containing attributes that were removed or replaced by other values is added to the Original Attribute Sequence (0400,0561).

#### 4.2.2.4 Association Acceptance Policy

##### 4.2.2.4.1 Activity – “Receive Storage Request”

##### 4.2.2.4.1.1 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCP of the Storage Service Class. The storage SCP accepts incoming C-Store Request from any configured AE Title, receives supported objects transmitted on that association and stores them in the local database.

##### 4.2.2.4.1.2 Accepted Presentation Contexts

For all supported Transfer Objects (see “Table 1: Network Services” section “SOP Classes Created by the syngo® MR XA50A” and “SOP Classes Managed by the syngo® MR XA50A” in the [“Conformance Statement Overview”](#).) the Transfer Syntaxes described in Table 10.

Generally, all Presentation Contexts are accepted as long as they contain at least one suitable Transfer Syntax. All other Presentation Contexts are rejected.

There is no Extended Negotiation as an SCP.

##### 4.2.2.4.1.3 SOP-specific Conformance Statement for Storage SOP classes

The syngo® MR XA50A conforms to the Full Storage Class at Level 2.

In case of a successful C-STORE operation, the image has successfully been received in the negotiated transfer syntax.

The Storage AE of the syngo® MR XA50A returns the status “success” when the data is received, and a minimal image header validation has been performed.

The following header attributes must be available and filled:

- SOP Class UID,
- Study Instance UID,
- Series Instance UID and
- SOP Instance UID.

Table 13 below list the status codes that the syngo® MR XA50A can return:

**Table 13: Storage C-STORE Response Status**

| Service Status | Further Meaning   | Error Code | Reason                                                                                                  |
|----------------|-------------------|------------|---------------------------------------------------------------------------------------------------------|
| Success        | Success           | 0000       | Image received correctly<br>(success notification is done after receiving, before indexing and storing) |
| Failure        | Out-of-resource   | A700       | No resource left in the Short Term Storage                                                              |
| Failure        | Unable to Process | Cxxx       | Error during instance reception                                                                         |

|         |                                   |      |                                                                         |
|---------|-----------------------------------|------|-------------------------------------------------------------------------|
| Failure | Data set does not match SOP Class | A9xx | The data set is not conform to the SOP Class contained in the resource. |
|---------|-----------------------------------|------|-------------------------------------------------------------------------|

*Restriction: Depending on response configuration successful operation does not guarantee storage on disk and storage of header data in the database.*

#### 4.2.2.4.1.4 Other SOP specific behavior

If an image is received that is already stored in the database - identified by the SOP Instance UID - the new image will be ignored. The existing instance is not superseded.

### 4.2.3 Storage Commitment AE Specification

#### 4.2.3.1 SOP Classes

The Storage Commitment AE of the syngo® MR XA50A provides standard conformance to the SOP Class listed in "Table 1: Network Services" section "Storage Commitment" in the ["Conformance Statement Overview"](#).

#### 4.2.3.2 Association Policy

**Table 14: Association Policies**

|                                                                         |                       |
|-------------------------------------------------------------------------|-----------------------|
| Application Context Name                                                | 1.2.840.10008.3.1.1.1 |
| PDU size                                                                | 32 kB <sup>1</sup>    |
| Maximum number of simultaneous associations as an association acceptor  | 12 <sup>1</sup>       |
| Maximum number of simultaneous associations as an association initiator | unlimited             |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.3.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 15: Asynchronous Nature as an Association Initiator**

|                                                         |    |
|---------------------------------------------------------|----|
| Maximum number of outstanding asynchronous transactions | 10 |
|---------------------------------------------------------|----|

#### 4.2.3.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to "**

**Table 3: Implementation Identifying Information" in the ["Conformance Statement Overview"](#).**

<sup>1</sup> Default, the value is configurable

### 4.2.3.3 Association Initiation Policy

#### 4.2.3.3.1 Activity "Send Initial Storage Commitment"

##### 4.2.3.3.1.1 Description and Sequencing of Activities

The *syngo*® MR XA50A serves as a SCU of the Storage Commitment Service Class. After successful transfer of Imaging Objects to a configured Archive, the Storage Commitment SCU initiates an N-ACTION Request, if Storage Commitment is configured. This request will be sent on a different association than the storage request.

The Storage Commitment Request will be sent out with a delay, in order to ensure that the remote node properly indexes received instances. The delay time is configurable with a default delay of 10 minutes.

The system may issue one N-ACTION Request for a complete set (bundle) of instances or issue one N-ACTION-Request per instance. This behavior is configurable; the default value is "bundled".

The *syngo*® MR XA50A will accept the N-Event-Report-Request on the same association if sent immediately after the N-ACTION-Response. However, it will not wait for it. The association is closed after three seconds.

##### 4.2.3.3.1.2 Proposed Presentation Contexts

The *syngo*® MR XA50A DICOM application supports the presentation contexts listed in the following table for the Storage Commitment Service Class.

**Table 16: Proposed Presentation Contexts for Storage Commitment**

| Abstract Syntax               |                      | Transfer Syntax           |                     | Role | Ext. Neg. |
|-------------------------------|----------------------|---------------------------|---------------------|------|-----------|
| Name                          | UID                  | Name List                 | UID List            |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                               |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                               |                      | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

##### 4.2.3.3.1.3 SOP specific Conformance for SOP classes

The behavior of *syngo*® MR XA50A when encountering status codes in an N-ACTION response is summarized in Table 17:

**Table 17: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning                                | Error Code         | Behavior                                                                                                          |
|----------------|------------------------------------------------|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Error          | Any failure that occurs                        | Any none null Code | Failure reported to user; corresponding object(s) will be marked as "Archived failed"                             |
| Success        | All Instances are available on the remote node | 0000               | Success reported to user; in case failures exist, the corresponding instances will be marked as "Archived failed" |

Table 18 below indicates the behavior if exceptions occur:

Table 18: DICOM Command Communication Failure Behavior

| Exception           | Behavior                                                                                  |
|---------------------|-------------------------------------------------------------------------------------------|
| Timeout             | Failure reported to user (Timeout configurable; default 30s); the request will be retried |
| Association Aborted | Failure reported to user; the request will be retried                                     |

#### 4.2.3.4 Association Acceptance Policy

##### 4.2.3.4.1 Activity "Receive Reply to Initial Storage Commitment"

##### 4.2.3.4.1.1 Description and Sequencing of Activities

The syngo® MR XA50A supports the reverse role negotiation of the Storage Commitment Service Class as the SCU. It accepts incoming N-EVENT-REPORT-Request, if they do not arrive on the same association as the N-ACTION-Request.

##### 4.2.3.4.1.2 Accepted Presentation Contexts

The syngo® MR XA50A DICOM application supports the presentation contexts listed in the following table for the Storage Commitment Service Class.

Table 19: Presentation Context Table "Update Flag Information"

| Abstract Syntax |                      | Transfer Syntax           |                     | Role | Ext. Neg. |
|-----------------|----------------------|---------------------------|---------------------|------|-----------|
| Name            | UID                  | Name List                 | UID List            |      |           |
| Storage         | 1.2.840.10008.1.20.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None      |
| Commitment      |                      | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |
| Push Model      |                      | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |           |

##### 4.2.3.4.1.3 SOP-specific Conformance Statement for Storage Commitment SOP classes

The Storage Commitment SCU provides standard conformance to the Storage Commitment SOP Class.

#### 4.2.4 Query/Retrieve AE Specification

##### 4.2.4.1 SOP Classes

The Query/Retrieve AE provides Standard Conformance to the SOP Classes listed in "Table 1: Network Services" section "Query/Retrieve" in the "[Conformance Statement Overview](#)".

##### 4.2.4.2 Association Policy

Table 20: Association Policies

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size                 | 32 kB <sup>1</sup>    |

<sup>1</sup> Default, the value is configurable

|                                                                                |                 |
|--------------------------------------------------------------------------------|-----------------|
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup> |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited       |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.4.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 21: Asynchronous Nature as an Association Initiator**

|                                                                |    |
|----------------------------------------------------------------|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|----------------------------------------------------------------|----|

#### 4.2.4.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to “**

Table 3: Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

#### 4.2.4.3 Association Initiation Policy

##### 4.2.4.3.1 Activity “Querying a Remote Node” for Instances

##### 4.2.4.3.1.1 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCU for the following SOP Class

- Study Root Q/R Information Model – FIND SOP Class

Using the attributes specified by the user as Query Keys (in accordance with the query model) the Query SCU initiates a C-FIND Request and displays the responses to the user.

##### 4.2.4.3.1.2 Proposed Presentation Contexts

The syngo® MR XA50A will propose Presentation Contexts as shown in the following table:

**Table 22: Proposed Presentation Contexts for Query**

| Abstract Syntax                                     |                             | Transfer Syntax           |                     | Role | Ext. Neg. |
|-----------------------------------------------------|-----------------------------|---------------------------|---------------------|------|-----------|
| Name                                                | UID                         | Name List                 | UID List            |      |           |
| Study Root Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | Yes       |
|                                                     |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                                                     |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

Table 23: Extended Negotiation as an SCU

| Name                                                | UID                         | Extended Negotiation                                          |
|-----------------------------------------------------|-----------------------------|---------------------------------------------------------------|
| Study Root Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Relational Query will be negotiated as defined in DICOM PS3.4 |

#### 4.2.4.3.1.3 SOP Specific Conformance Statement to Query SOP classes

The syngo® MR XA50A checks for the following status codes in the Query SCP's C-FIND-Response:

Table 24: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning                                                                                              | Error Code         | Behavior                           |
|----------------|--------------------------------------------------------------------------------------------------------------|--------------------|------------------------------------|
| Failure        | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process                  | Any none null Code | Failure reported to user           |
| Pending        | All optional keys are supported the same manner as Required Keys.                                            | FF00               | Pending state is indicated to user |
|                | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01               | Pending state is indicated to user |
| Success        | Query has been performed successfully.                                                                       | 0000               | Success reported to user           |

Table 25 below indicates the behavior if exceptions occur:

Table 25: DICOM Command Communication Failure Behavior

| Exception           | Behavior                                                     |
|---------------------|--------------------------------------------------------------|
| Timeout             | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user                                     |

The syngo® MR XA50A supports the following query levels:

- Study
- Series

Matching Keys on Instance Level is not supported by the syngo® MR XA50A as SCU.

The following table lists the various attributes at Study and Series level, which can be used for hierarchical queries as well as return values for display. The display capabilities are highly configurable and "yes" indicates that it is possible to configure display of the data:

Table 26: Attributes supported for Study/Series Query - SCU

| Attribute Name       | Tag         | Type | User input  | UI  |
|----------------------|-------------|------|-------------|-----|
| <b>Study Level</b>   |             |      |             |     |
| Patient's Name       | (0010,0010) | O    | enter value | Yes |
| Patient ID           | (0010,0020) | O    | enter value | Yes |
| Issuer of Patient ID | (0010,0021) | O    | enter value | Yes |
| Patient's Birth Date | (0010,0030) | O    | enter value | Yes |

| Attribute Name                                            | Tag                       | Type | User input  | UI  |
|-----------------------------------------------------------|---------------------------|------|-------------|-----|
| Patient's Birth Time                                      | (0010,0032)               | O    | enter value | Yes |
| Patient's Sex                                             | (0010,0040)               | O    | enter value | Yes |
| Accession Number                                          | (0008,0050)               | O    | enter value | Yes |
| Study ID                                                  | (0020,0010)               | O    | enter value | Yes |
| Study Instance UID                                        | (0020,000D)               | U    | enter value | Yes |
| Study Date                                                | (0008,0020)               | O    | enter value | Yes |
| Study Time                                                | (0008,0030)               | O    | enter value | Yes |
| Referring Physician's Name                                | (0008,0090)               | O    | enter value | Yes |
| Study Description                                         | (0008,1030)               | O    | enter value | Yes |
| Number of Study related Instances                         | (0020,1208)               | O    | -           | Yes |
| Modalities in Study                                       | (0008,0061)               | O    | enter value | Yes |
| Number of Study Related Series                            | (0020,1206)               | O    | -           | Yes |
| <b>Series Level</b>                                       |                           |      |             |     |
| Modality                                                  | (0008,0060)               | O    | enter value | Yes |
| Series Date                                               | (0008,0021)               | O    | enter value | Yes |
| Series Time                                               | (0008,0031)               | O    | enter value | Yes |
| Number of Series related Instances                        | (0020,1209)               | O    | -           | Yes |
| Series Number                                             | (0020,0011)               | O    | enter value | Yes |
| Series Description                                        | (0008,103E)               | O    | enter value | Yes |
| Request Attributes Sequence \ Requested Procedure ID      | (0040,0275) \ (0040,1001) | O    | enter value | Yes |
| Request Attributes Sequence \ Scheduled Procedure Step ID | (0040,0275) \ (0040,0009) | O    | enter value | Yes |
| Performed Procedure Step Start Date                       | (0040,0244)               | O    | enter value | Yes |
| Performed Procedure Step Start Time                       | (0040,0245)               | O    | enter value | Yes |
| Series Instance UID                                       | (0020,000E)               | U    | -           | Yes |

- Activity "Retrieve Instances from a remote node"

#### 4.2.4.3.1.4 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCU for the following SOP Classes

- Study Root Q/R Information Model – MOVE SOP Class

The C-MOVE-Request is used to retrieve the selected imaging objects. The Retrieve AE supports the query model Patient Root, Study Root, Patient/Study Root.

#### 4.2.4.3.1.5 Proposed Presentation Contexts

The syngo® MR XA50A proposes Presentation Contexts shown in the following table:

Table 27: Proposed Presentation Contexts for Retrieve and Activity "MOVE SCU"

| Abstract Syntax                        |                             | Transfer Syntax           |                     | Role | Ext. Neg. |
|----------------------------------------|-----------------------------|---------------------------|---------------------|------|-----------|
| Name                                   | UID                         | Name List                 | UID List            |      |           |
| Study Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                                        |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                                        |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.4.3.1.6 SOP Specific Conformance Statement for Move SCU Classes

The behavior of *syngo*® MR XA50A when encountering status codes in a C-MOVE response is summarized in Table 28

Table 28: DICOM Command Response Status Handling Behavior

| Service Status | Further Meaning                                                                                                       | Error Code         | Behavior                          |
|----------------|-----------------------------------------------------------------------------------------------------------------------|--------------------|-----------------------------------|
| Error          | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process; Move destination unknown | Any none null Code | Failure reported to user          |
| Pending        | Move Operation continues                                                                                              | FF00               | Operation continues in background |
| Success        | Move has been performed successfully.                                                                                 | 0000               | Success reported to user          |

Table 29 below indicates the behavior if exceptions occur:

Table 29: DICOM Command Communication Failure Behavior

| Exception           | Behavior                                                     |
|---------------------|--------------------------------------------------------------|
| Timeout             | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user                                     |

#### 4.2.4.4 Association Acceptance Policy

The *syngo*® MR XA50A provides SCU and SCP functionality.

#### 4.2.5 Modality Worklist AE Specification

##### 4.2.5.1 SOP Classes

The Modality Worklist AE provides Standard Conformance to the SOP Class "Modality Worklist Information Model – FIND" as it is listed in "Table 1: Network Services" in the "[Conformance Statement Overview](#)".

##### 4.2.5.2 Association Policy

Table 30: Association Policies

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size                 | 32 kB <sup>1</sup>    |

<sup>1</sup> Default, the value is configurable

|                                                                                |                 |
|--------------------------------------------------------------------------------|-----------------|
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup> |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited       |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.5.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 31: Asynchronous Nature as an Association Initiator**

|                                                                |    |
|----------------------------------------------------------------|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|----------------------------------------------------------------|----|

#### 4.2.5.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to “**

Table 3: Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

#### 4.2.5.3 Association Initiation Policy

- Activity “Querying a Remote Node” for Modality Worklist

##### 4.2.5.3.1.1 Description and Sequencing of Activities

The syngo® MR XA50A serves as a SCU of the Modality Worklist service. It performs worklist queries by issuing a C-FIND request at regular intervals. In addition a worklist request can be triggered manually.

##### 4.2.5.3.1.2 Proposed Presentation Contexts

The syngo® MR XA50A will propose Presentation Contexts as shown in the following table:

**Table 32: Proposed Presentation Contexts for Worklist**

| Abstract Syntax         |                        | Transfer Syntax           |                     | Role | Ext. Neg. |
|-------------------------|------------------------|---------------------------|---------------------|------|-----------|
| Name                    | UID                    | Name List                 | UID List            |      |           |
| Modality Worklist- FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                         |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                         |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.5.3.1.3 SOP Specific Conformance for SOP Classes

##### Search Key Attributes of the Worklist C-FIND

The *syngo*® MR XA50A Modality Worklist SCU supports "broad worklist queries" with all required search keys. The following tables describe the "broad query" search keys that the SCU supports. The list is configurable in 'DICOM Modality Worklist Query'.

**Table 33: Broad Query search keys**

| Attribute Name                         | Tag          | Matching Key Type | Query Value                     |
|----------------------------------------|--------------|-------------------|---------------------------------|
| <b>Scheduled Procedure Step</b>        |              |                   |                                 |
| Scheduled Procedure Step Sequence      | (0040,0100)  | R                 |                                 |
| >Modality                              | (0008,0060)  | R                 | <configured by Modality> or "*" |
| >Scheduled Station AE Title            | (0040,0001)  | R                 | <own AET> or "*"1               |
| >Scheduled Procedure Step Start Date   | (0040,0002)  | R                 | Range from UI 2                 |
| >Scheduled Station Name                | (0040,0010)  | O                 |                                 |
| >Scheduled Procedure Step Location     | (0040,0011)  | O                 |                                 |
| >Scheduled Procedure Step Status       | (0040,0020)  | O                 |                                 |
| >Scheduled Performing Physician's Name | (0040,0006)  | O                 |                                 |
| Requested Procedure Description        | (0032,1060)  | O                 |                                 |
| Requested Procedure Priority           | (0040,1003)  | O                 |                                 |
| Patient Transport Arrangements         | (0040,1004)  | O                 |                                 |
| Requested Procedure Comments           | (0040,1400)  | O                 |                                 |
| Requested Procedure Code Sequence      | (0032,1064)  | O                 |                                 |
| >Code Value                            | (0008,0100)  | O                 |                                 |
| Requesting Physician                   | (0032,1032)  | O                 |                                 |
| Referring Physicians Name              | (0008,0090)  | O                 |                                 |
| Current Patient Location               | (0038,0300)  | O                 |                                 |
| Pregnancy Status                       | (0010, 21C0) | O                 |                                 |
| Medical Alerts                         | (0010,2000)  | O                 |                                 |
| Allergies                              | (0010,2110)  | O                 |                                 |

##### Return Key Attributes of the Modality Worklist C-FIND

The *syngo*® MR XA50A Modality Worklist SCU supports worklist queries with return key attributes of all types. The following tables describe the return keys that the SCU supports.

An "yes" in the **UI** column indicates that the attribute may be visualized when browsing the Worklist results with the Browser. The Browser display is additionally influenced by the related Browser configuration.

<sup>1</sup> This depends on user configuration (Administration Portal->Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist) if the "own AET" is provided or not.

<sup>2</sup> A time window can be configured by defining how many days to look into the past and into the future (Administration Portal-> Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist)

Table 34: Modality Worklist C-Find Return keys

| Attribute Name                         | Tag         | Return Key Type | UI  | Notes                                                                                                                          |
|----------------------------------------|-------------|-----------------|-----|--------------------------------------------------------------------------------------------------------------------------------|
| <b>SOP Common</b>                      |             |                 |     |                                                                                                                                |
| Specific Character Set                 | (0008,0005) | 1C              | No  |                                                                                                                                |
| <b>Scheduled Procedure Step</b>        |             |                 |     |                                                                                                                                |
| Scheduled Procedure Step Sequence      | (0040,0100) | 1               | No  |                                                                                                                                |
| >Modality                              | (0008,0060) | 1               | Yes |                                                                                                                                |
| >Scheduled Station AE Title            | (0040,0001) | 1               |     | "Scheduled Station AE Title" is taken as default for "Performed Station AE Title"                                              |
| >Scheduled Procedure Step Start Date   | (0040,0002) | 1               | Yes |                                                                                                                                |
| >Scheduled Procedure Step Start Time   | (0040,0003) | 1               | Yes |                                                                                                                                |
| >Scheduled Procedure Step End Date     | (0040,0004) | 3               | No  |                                                                                                                                |
| >Scheduled Procedure Step End Time     | (0040,0005) | 3               | No  |                                                                                                                                |
| >Scheduled Performing Physician's Name | (0040,0006) | 1               | Yes | "Scheduled Performing Physician's Name" is taken as default for "Performing Physician's Name"                                  |
| >Scheduled Procedure Step Description  | (0040,0007) | 1C              | Yes | "Scheduled Procedure Step Description" is taken as default for "Performed Procedure Step Description"                          |
| >Scheduled Protocol Code Sequence **   | (0040,0008) | 1C              | No  | Uses universal sequence match<br>"Scheduled Protocol Code Sequence" is taken as default for "Performed Protocol Code Sequence" |
| >>Code Value                           | (0008,0100) | 1C              | No  |                                                                                                                                |
| >>Coding Scheme Designator             | (0008,0102) | 1C              | No  |                                                                                                                                |
| >>Coding Scheme Version                | (0008,0103) | 3               | No  |                                                                                                                                |
| >>Code Meaning                         | (0008,0104) | 3               | No  |                                                                                                                                |
| >>Mapping Resource                     | (0008,0105) | 3               | No  |                                                                                                                                |
| >>Context Group Version                | (0008,0106) | 3               | No  |                                                                                                                                |
| >>Context Group Local Version          | (0008,0107) | 3               | No  |                                                                                                                                |
| >>Context Group Extension Flag         | (0008,010B) | 3               | No  |                                                                                                                                |
| >>Context Group Extension Creator UID  | (0008,010D) | 3               | No  |                                                                                                                                |
| >>Context Identifier                   | (0008,010F) | 3               | No  |                                                                                                                                |
| >Scheduled Procedure Step ID           | (0040,0009) | 1               | Yes | "Scheduled Procedure Step ID" is taken as default for "Performed Procedure Step ID"                                            |
| >Scheduled Station Name                | (0040,0010) | 2               | Yes |                                                                                                                                |
| >Scheduled Procedure Step Location     | (0040,0011) | 2               | No  | "Scheduled Procedure Step Location" is taken as default for "Performed Location"                                               |

| Attribute Name                            | Tag         | Return Key Type | UI  | Notes                                                                                 |
|-------------------------------------------|-------------|-----------------|-----|---------------------------------------------------------------------------------------|
| >Scheduled Procedure Step Status          | (0040,0020) | 3               | No  |                                                                                       |
| >Comments on the Scheduled Procedure Step | (0040,0400) | 3               | No  |                                                                                       |
| <b>Requested Procedure</b>                |             |                 |     |                                                                                       |
| Study Date                                | (0008,0020) | 3               | No  |                                                                                       |
| Study Time                                | (0008,0030) | 3               | No  |                                                                                       |
| Referenced Study Sequence                 | (0008,1110) | 2               | No  | Uses universal sequence match                                                         |
| >Referenced SOP Class UID                 | (0008,1150) | 1C              | No  |                                                                                       |
| >Referenced SOP Instance UID              | (0008,1155) | 1C              | No  |                                                                                       |
| Study Instance UID                        | (0020,000D) | 1               | No  |                                                                                       |
| Requested Procedure Description           | (0032,1060) | 1C              | Yes | Uses universal sequence match                                                         |
| Requested Procedure Code Sequence         | (0032,1064) | 1C              | No  | "Requested Procedure Code Sequence" is taken as default for "Procedure Code Sequence" |
| >Code Value                               | (0008,0100) | 1C              | No  |                                                                                       |
| >Coding Scheme Designator                 | (0008,0102) | 1C              | No  |                                                                                       |
| >Coding Scheme Version                    | (0008,0103) | 3               | No  |                                                                                       |
| >Code Meaning                             | (0008,0104) | 3               | No  |                                                                                       |
| >>Mapping Resource                        | (0008,0105) | 3               | No  |                                                                                       |
| >>Context Group Version                   | (0008,0106) | 3               | No  |                                                                                       |
| >>Context Group Local Version             | (0008,0107) | 3               | No  |                                                                                       |
| >>Context Group Extension Flag            | (0008,010B) | 3               | No  |                                                                                       |
| >>Context Group Extension Creator UID     | (0008,010D) | 3               | No  |                                                                                       |
| >>Context Identifier                      | (0008,010F) | 3               | No  |                                                                                       |
| Requested Procedure ID                    | (0040,1001) | 1               | Yes | "Requested Procedure ID" is taken as default for "Study ID"                           |
| Reason for the Requested Procedure        | (0040,1002) | 3               | No  |                                                                                       |
| Requested Procedure Priority              | (0040,1003) | 2               | Yes |                                                                                       |
| Patient Transport Arrangements            | (0040,1004) | 2               | No  |                                                                                       |
| Confidentiality Code                      | (0040,1008) | 3               | No  |                                                                                       |
| Reporting Priority                        | (0040,1009) | 3               | Yes |                                                                                       |
| Names of intended Recipients of Results   | (0040,1010) | 3               | No  |                                                                                       |
| Requested Procedure Comments              | (0040,1400) | 3               | No  |                                                                                       |
| <b>Imaging Service Request</b>            |             |                 |     |                                                                                       |
| Accession Number                          | (0008,0050) | 2               | Yes |                                                                                       |
| Referring Physician's Name                | (0008,0090) | 2               | Yes |                                                                                       |

| Attribute Name                                | Tag         | Return Key Type | UI  | Notes                                        |
|-----------------------------------------------|-------------|-----------------|-----|----------------------------------------------|
| Requesting Physician                          | (0032,1032) | 2               | Yes |                                              |
| Requesting Service                            | (0032,1033) | 3               | No  |                                              |
| Issuing Date of Imaging Service Request       | (0040,2004) | 3               | No  |                                              |
| Issuing Time of Imaging Service Request       | (0040,2005) | 3               | No  |                                              |
| Placer Order Number / Imaging Service Request | (0040,2016) | 3               | No  | Old tag (0040,2006) is retired and not used. |
| Filler Order Number / Imaging Service Request | (0040,2017) | 3               | No  | Old tag (0040,2007) is retired and not used. |
| Order entered by ...                          | (0040,2008) | 3               | No  |                                              |
| Order Enterer's location                      | (0040,2009) | 3               | No  |                                              |
| Order Callback Phone Number                   | (0040,2010) | 3               | No  |                                              |
| Imaging Service Request Comments              | (0040,2400) | 3               | No  |                                              |
| <b>Visit Identification</b>                   |             |                 |     |                                              |
| Admission ID                                  | (0038,0010) | 2               | Yes |                                              |
| Issuer of Admission ID                        | (0038,0011) | 3               | No  |                                              |
| Institution Name                              | (0008,0080) | 3               | No  |                                              |
| Institution Address                           | (0008,0081) | 3               | No  |                                              |
| <b>Visit Status</b>                           |             |                 |     |                                              |
| Current Patient Location                      | (0038,0300) | 2               | Yes |                                              |
| <b>Visit Admission</b>                        |             |                 |     |                                              |
| Admitting Diagnosis Description               | (0008,1080) | 3               | Yes |                                              |
| Admitting Date                                | (0038,0020) | 3               | No  |                                              |
| <b>Patient Identification</b>                 |             |                 |     |                                              |
| Patient's Name                                | (0010,0010) | 1               | Yes |                                              |
| Patient ID                                    | (0010,0020) | 1               | Yes |                                              |
| Issuer of Patient ID                          | (0010,0021) | 3               | No  |                                              |
| Other Patient IDs                             | (0010,1000) | 3               | Yes |                                              |
| Other Patient Birth Names                     | (0010,1001) | 3               | Yes |                                              |
| Patient's Birth Name                          | (0010,1005) | 3               | No  |                                              |
| <b>Patient Demographic</b>                    |             |                 |     |                                              |
| Patient's Birth Date                          | (0010,0030) | 2               | Yes |                                              |
| Patient's Birth Time                          | (0010,0032) | 3               | No  |                                              |
| Patient's Sex                                 | (0010,0040) | 2               | Yes |                                              |
| Patient's Insurance Plan Code Sequence        | (0010,0050) | 3               | No  | Uses universal sequence match                |
| >Code Value                                   | (0008,0100) | 1C              | No  |                                              |
| >Coding Scheme Designator                     | (0008,0102) | 1C              | No  |                                              |
| >Coding Scheme Version                        | (0008,0103) | 3               | No  |                                              |

| Attribute Name             | Tag         | Return Key Type | UI  | Notes |
|----------------------------|-------------|-----------------|-----|-------|
| >Code Meaning              | (0008,0104) | 3               | No  |       |
| Patient's Age              | (0010,1010) | 3               | No  |       |
| Patient's Size             | (0010,1020) | 3               | Yes |       |
| Patient's Weight           | (0010,1030) | 2               | Yes |       |
| Patient's Address          | (0010,1040) | 3               | Yes |       |
| Military Rank              | (0010,1080) | 3               | Yes |       |
| Branch of Service          | (0010,1081) | 3               | No  |       |
| Ethnic Group               | (0010,2160) | 3               | Yes |       |
| Patient Comments           | (0010,4000) | 3               | Yes |       |
| <b>Patient Medical</b>     |             |                 |     |       |
| Medical Alerts             | (0010,2000) | 2               | Yes |       |
| Allergies                  | (0010,2110) | 2               | Yes |       |
| Pregnancy Status           | (0010,21C0) | 2               | Yes |       |
| Smoking Status             | (0010,21A0) | 3               | Yes |       |
| Last Menstrual Date        | (0010,21D0) | 3               | Yes |       |
| Additional Patient History | (0010,21B0) | 3               | Yes |       |
| Special Needs              | (0038,0050) | 2               | Yes |       |

The behavior of the *syngo*® MR XA50A when encountering status codes in a C-FIND response is summarized in Table 35.

**Table 35: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning                                                                                              | Error Code         | Behavior                           |
|----------------|--------------------------------------------------------------------------------------------------------------|--------------------|------------------------------------|
| Error          | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process                  | Any none null Code | Failure reported to user           |
| Pending        | All optional keys are supported the same manner as Required Keys.                                            | FF00               | Pending state is indicated to user |
|                | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01               | Pending state is indicated to user |
| Success        | Query has been performed successfully.                                                                       | 0000               | Success reported to user           |

Table 36 below indicates the behavior if exceptions occur:

**Table 36: DICOM Command Communication Failure Behavior**

| Exception | Behavior                                                     |
|-----------|--------------------------------------------------------------|
| Timeout   | Failure reported to user (Timeout configurable; default 30s) |

Association Aborted Failure reported to user

#### 4.2.5.4 Association Acceptance Policy

The syngo® MR XA50A does not provide the functionality of a SCP of the Modality Worklist – Find SOP Class.

#### 4.2.6 Modality Performed Procedure Step AE Specification

##### 4.2.6.1 SOP Classes

The Modality Performed Procedure Step AE provides Standard Conformance to the SOP Class “Modality Performed Procedure Step” as it is listed in “Table 1: Network Services” section “Worklist Management” in the [“Conformance Statement Overview”](#).

##### 4.2.6.2 Association Policy

**Table 37: Association Policies**

|                                                                                |                       |
|--------------------------------------------------------------------------------|-----------------------|
| <b>Application Context Name</b>                                                | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>                                                                | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

##### 4.2.6.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 38: Asynchronous Nature as an Association Initiator**

|                                                                |    |
|----------------------------------------------------------------|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|----------------------------------------------------------------|----|

##### 4.2.6.2.2 Implementation Identifying Information

**For Implementation Identifying Information please refer to “**

Table 3: Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

<sup>1</sup> Default, the value is configurable

#### 4.2.6.3 Association Initiation Policy

##### 4.2.6.3.1 Activity "Create Modality Performed Procedure Step"

##### 4.2.6.3.1.1 Description and Sequencing of Activities

The *syngo*® MR XA50A serves as a SCU of the Modality Performed Procedure Step SOP Class. It sends N-CREATE request to inform the Information System that a Procedure Step has been started.

##### 4.2.6.3.1.2 Accepted Presentation Contexts

The *syngo*® MR XA50A proposes Presentation Contexts as shown in the following table:

**Table 39: Acceptable Presentation Contexts Activity "Create MPPS"**

| Abstract Syntax                         |                         | Transfer Syntax           |                     | Role | Ext.<br>Neg. |
|-----------------------------------------|-------------------------|---------------------------|---------------------|------|--------------|
| Name                                    | UID                     | Name List                 | UID List            |      |              |
| Modality<br>Performed<br>Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None         |
|                                         |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |              |
|                                         |                         | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |              |

#### 4.2.6.3.1.3 SOP specific Conformance for MPPS SOP class Table 40

The behavior of syngo® MR XA50A when encountering status codes in an N-CREATE-RSP response is summarized in Table 40:

**Table 40: MPPS N-CREATE Response Status Handling Behavior**

| Service Status | Further Meaning                               | Error Code         | Behavior             |
|----------------|-----------------------------------------------|--------------------|----------------------|
| Error          | MPPS creation request could not be processed. | Any none null Code | MPPS is not created. |
| Success        | MPPS creation request processed successfully. | 0000               | MPPS is created.     |

#### 4.2.6.3.2 Activity "Update Modality Performed Procedure Step"

##### 4.2.6.3.2.1 Description and Sequencing of Activities

When the procedure step has been finished, the syngo® MR XA50A sends N-SET request to inform the Information System about the finalization of the procedure step (COMPLETED or DISCONTINUED).

##### 4.2.6.3.2.2 Proposed Presentation Contexts

The syngo® MR XA50A proposes Presentation Contexts as shown in the following table:

**Table 41: Acceptable Presentation Contexts Activity "Update MPPS"**

| Abstract Syntax                   |                         | Transfer Syntax           |                     | Role | Ext. Neg. |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|-----------|
| Name                              | UID                     | Name List                 | UID List            |      |           |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                                   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                                   |                         | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.6.3.2.3 SOP specific Conformance for MPPS SOP class

The behavior of syngo® MR XA50A when encountering status codes in an N-SET-RSP response is summarized in Table 42:

**Table 42: MPPS N-SET Response Status Handling Behavior**

| Service Status | Further Meaning                                   | Error Code         | Behavior             |
|----------------|---------------------------------------------------|--------------------|----------------------|
| Error          | MPPS update request could not be processed.       | Any none null Code | MPPS is not updated. |
| Success        | MPPS update request could processed successfully. | 0000               | MPPS is updated.     |

#### 4.2.6.4 Association Acceptance Policy

The syngo® MR XA50A does not provide the functionality of a SCP of the Modality Performed Procedure Step SOP Class.

## 4.2.7 Print AE Specification

### 4.2.7.1 SOP Classes

The Print AE provides Standard Conformance to the SOP Classes listed in "Table 1: Network Services" section "Print Management" in the [Conformance Statement Overview](#).

### 4.2.7.2 Association Policy

**Table 43: Association Policies**

|                                                                                |                       |
|--------------------------------------------------------------------------------|-----------------------|
| <b>Application Context Name</b>                                                | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>                                                                | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

The syngo® MR XA50A contains a limitation of 512 kB for the maximum PDU size. By default, the PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.7.2.1 Asynchronous Nature

The syngo® MR XA50A supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 44: Asynchronous Nature as an Association Initiator**

|                                                                |    |
|----------------------------------------------------------------|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|----------------------------------------------------------------|----|

### 4.2.7.3 Association Initiation Policy

#### 4.2.7.3.1 Activity Print Film

##### 4.2.7.3.1.1 Description and Sequencing of Activities

Whenever a film-sheet is prepared by the user, it is forwarded to the Printer Job queue. As soon as the associated Printer device is available the job is activated and an association is established.

After the film sheet is internally processed, converted to a Standard/1,1 layout and the page image is sent to the printer, the status is controlled by awaiting any N-EVENT-REPORT message throughout the transfer until the last image or film-sheet is sent.

If the response from the remote application contains a status other than Success or Warning the printing is stopped and the job status is set to Aborted.

<sup>1</sup> Default, the value is configurable

#### 4.2.7.3.1.2 Proposed Presentation Context

The syngo® MR XA50A proposes Presentation Contexts as shown in the following table:

**Table 45: Presentation Contexts for the Activity “Print Film”**

| Abstract Syntax                                 |                         | Transfer Syntax           |                   | Role | Ext. Neg. |
|-------------------------------------------------|-------------------------|---------------------------|-------------------|------|-----------|
| Name                                            | UID                     | Name List                 | UID List          |      |           |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Film Session SOP Class                    | 1.2.840.10008.5.1.1.1   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Film Box SOP Class                        | 1.2.840.10008.5.1.1.2   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Grayscale Image Box SOP Class             | 1.2.840.10008.5.1.1.4   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Color Image Box SOP Class                 | 1.2.840.10008.5.1.1.4.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Printer SOP Class                               | 1.2.840.10008.5.1.1.16  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Print Job SOP Class                             | 1.2.840.10008.5.1.1.14  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

#### 4.2.7.3.1.3 SOP Specific Conformance

The syngo® MR XA50A Print SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class and the Basic Color Print Management Meta SOP Class.

The application uses a configuration platform to define the properties of the connected DICOM SCP, e.g.:

- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP

The printing is only suspended in the case of a failure return status of the SCP.

##### Basic Film Session SOP Class

The Basic Film Session information object definition describes all the user-defined parameters, which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes that are printed on one hardcopy printer.

The syngo® MR XA50A Print Management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the attributes listed in the table below:

**Table 46: Attributes for the N-CREATE-RQ of the Basic Film Session**

| Attribute Name   | Tag         | Usage SCU | Supported Values |
|------------------|-------------|-----------|------------------|
| Number of Copies | (2000,0010) | U         | 1                |
| Medium Type      | (2000,0030) | U         | BLUE FILM        |
|                  |             |           | CLEAR FILM       |
|                  |             |           | PAPER            |
|                  |             |           | MAMMO BLUE FILM  |
|                  |             |           | MAMMO CLEAR FILM |

The Affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session.

The Basic Film Session SOP class interprets the status codes (from N-CREATE-RSP messages) listed in the table below:

**Table 47: N-CREATE-RSP Status Handling Behavior for the Basic Film Session**

| Service Status | Further Meaning                   | Error Codes | Behavior                               |
|----------------|-----------------------------------|-------------|----------------------------------------|
| Warning        | Memory Allocation not supported   | B600        | Print job continues, warning is logged |
| Success        | Film session successfully created | 0000        | Print job continues                    |

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

#### Basic Film Box SOP Class

The Basic Film Box information object definition describes all user-defined parameters of one film of the film session including presentation parameters, which are common for all images on a given film sheet.

The Basic Film Box refers to one or more Image Boxes.

The *syngo*® MR XA50A Print Management SCU supports the following DIMSE Service elements for the Basic Film Box SOP Class as SCU:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP Class N-CREATE-RQ message uses the attributes listed below. The actual values for each attribute depend on DICOM printer configuration within the *syngo*® MR XA50A DICOM Print Management SCU:

**Table 48: Attributes for the N-CREATE-RQ of the Basic Film Session**

| Attribute Name                   | Tag         | Usage SCU | Supported Values      |
|----------------------------------|-------------|-----------|-----------------------|
| Image Display Format             | (2010,0010) | M         | STANDARD\1,1          |
| Referenced Film Session Sequence | (2010,0500) | M         |                       |
| > Referenced SOP Class UID       | (0008,1150) | M         | 1.2.840.10008.5.1.1.1 |
| > Referenced SOP Instance UID    | (0008,1155) | M         |                       |
| Film Orientation                 | (2010,0040) | M         | PORTRAIT, LANDSCAPE   |

| Attribute Name                                 | Tag         | Usage SCU | Supported Values                                                                      |
|------------------------------------------------|-------------|-----------|---------------------------------------------------------------------------------------|
| Film Size ID                                   | (2010,0050) | M         | 8INX10IN, 10INX12IN, 10INX14IN, 11INX14IN, 14INX14IN, 14INX17IN, 24CMX24CM, 24CMX30CM |
| Magnification Type                             | (2010,0060) | M         | BILINEAR, CUBIC, NONE, REPLICATE                                                      |
| Border Density                                 | (2010,0100) | U         | BLACK, WHITE                                                                          |
| Max Density                                    | (2010,0130) | U         | 0 < Value                                                                             |
| Min Density                                    | (2010,0120) | U         | 0 < Value < 50                                                                        |
| <b>Required if Presentation LUT is present</b> |             |           |                                                                                       |
| Reflective Ambient Light                       | (2010,0160) | U         | 0 < Value                                                                             |
| Illumination                                   | (2010,015E) | U         | 0 < Value                                                                             |
| Referenced Presentation LUT Sequence           | (2050,0500) | U         |                                                                                       |

For Page Mode printing, the Image Display format used is Standard\1,1.

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally to be further used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the syngo® MR XA50A print manager will issue an N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g., N-DELETE-RQ) on the Basic Film Box.

The Basic Film Box SOP class interprets the status codes listed in the tables below:

**Table 49: N-CREATE-RSP Status Handling Behavior for Basic Film Box**

| Service Status | Meaning                                                                                                                                                                                           | Error Codes | Behavior                                               |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|--------------------------------------------------------|
| Failure        | There is an existing Film Box that has not been printed and N-ACTION at the Film Session level is not supported. A new Film Box will not be created when a previous Film Box has not been printed | C616        | Print job is marked as failed and the reason is logged |
| Warning        | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead                                           | B605        | Print job continues and warning is logged              |
| Success        | Film Box successfully created                                                                                                                                                                     | 0000        | Print job continues                                    |

Table 50: N-ACTION-RSP Status Handling Behavior for Basic Film Box

| Service Status | Meaning                                                                                                                             | Error Codes | Behavior                                               |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------|-------------|--------------------------------------------------------|
| Failure        | Unable to create print job, print queue is full                                                                                     | C602        | Print job is marked as failed and the reason is logged |
|                | Image size is larger than images box size                                                                                           | C603        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size                                                                         | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Film box does not contain image box (empty page)                                                                                    | B603        | Print job continues and warning is logged              |
|                | Image size is larger than image box size, the image has been demagnified                                                            | B604        | Print job continues and warning is logged              |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.                                                    | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Film accepted for printing                                                                                                          | 0000        | Print job continues                                    |

#### Basic Grayscale Image Box SOP Class

The Basic Grayscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the attributes listed in the table below:

Table 51: Attributes for N-SET-RQ of Basic Grayscale Image Box

| Attribute Name                 | Tag         | Usage SCU | Supported Values |
|--------------------------------|-------------|-----------|------------------|
| Image Position                 | (2020,0010) | M         | 1                |
| Basic Grayscale Image Sequence | (2020,0110) | M         |                  |
| > Samples per Pixel            | (0028,0002) | M         | 1                |
| > Photometric Interpretation   | (0028,0004) | M         | MONOCHROME2      |
| > Rows                         | (0028,0010) | M         |                  |
| > Columns                      | (0028,0011) | M         |                  |
| > Pixel Aspect Ratio           | (0028,0034) | M         |                  |
| > Bits Allocated               | (0028,0100) | M         | 8,16             |
| > Bits Stored                  | (0028,0101) | M         | 8,12             |
| > High Bit                     | (0028,0102) | M         | 7,11             |
| > Pixel Representation         | (0028,0103) | M         | 0                |
| > Pixel Data                   | (7FE0,0010) | M         |                  |

The Basic Grayscale Image Box SOP class interprets the status codes as listed below:

**Table 52: N-SET-RSP Status Handling Behavior for the Basic Grayscale Image Box SOP Class**

| Service Status | Further Meaning                                                                                                                     | Error Codes | Behavior                                               |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------|-------------|--------------------------------------------------------|
| Failure        | Image contains more pixel than printer can print in Image Box                                                                       | C603        | Print job is marked as failed and the reason is logged |
|                | Insufficient memory in printer to store the image                                                                                   | C605        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size                                                                         | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Image size is larger than image box size, the image has been demagnified.                                                           | B604        | Print job continues and the reason is logged           |
|                | Requested MinDensity or MaxDensity outside of Printer's operating range                                                             | B605        | Print job continues and the reason is logged           |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.                                                    | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Image successfully stored in Image Box                                                                                              | 0000        | Print job continues                                    |

#### Basic Color Image Box SOP Class

The Basic Color Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Color Image Box SOP Class uses only the N-SET-RQ with the attributes listed below:

**Table 53: Attributes for N-SET-RQ of Basic Color Image Box**

| Attribute Name               | Tag         | Usage SCU | Supported Values |
|------------------------------|-------------|-----------|------------------|
| Image Position               | (2020,0010) | M         | 1                |
| Basic Color Image Sequence   | (2020,0111) | M         |                  |
| > Samples per Pixel          | (0028,0002) | M         | 3                |
| > Photometric Interpretation | (0028,0004) | M         | RGB              |
| > Planar Configuration       | (0028,0006) | M         | 0                |
| > Rows                       | (0028,0010) | M         |                  |
| > Columns                    | (0028,0011) | M         |                  |
| > Pixel Aspect Ratio         | (0028,0034) | M         |                  |
| > Bits Allocated             | (0028,0100) | M         | 8                |
| > Bits Stored                | (0028,0101) | M         | 8                |
| > High Bit                   | (0028,0102) | M         | 7                |
| > Pixel Representation       | (0028,0103) | M         | 0                |
| > Pixel Data                 | (7FE0,0010) | M         |                  |

The Color Image Box SOP class interprets the status codes listed below:

**Table 54: N-SET-RSP Status Handling Behavior for the Color Grayscale Image Box**

| Service Status | Meaning                                                                                                                             | Error Codes | Behavior                                               |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------|-------------|--------------------------------------------------------|
| Failure        | Image contains more pixel than printer can print in Image Box                                                                       | C603        | Print job is marked as failed and the reason is logged |
|                | Insufficient memory in printer to store the image                                                                                   | C605        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size                                                                         | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Image size larger than image box size                                                                                               | B604        | Print job continues and the reason is logged           |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.                                                    | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Image successfully stored in Image Box                                                                                              | 0000        | Print job continues                                    |

### Presentation LUT SOP Class

The objective of the Presentation LUT is to realize image hardcopy printing tailored for specific modalities, applications and user preferences.

The output of the Presentation LUT is Presentation Values (P-Values). P-Values are approximately related to human perceptual response. They are intended to facilitate common input for hardcopy. P-Values are intended to be independent of the specific class or characteristics of the hardcopy device.

The *syngo*® MR XA50A Print Management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Presentation LUT SOP Class uses only the N-CREATE-RQ with the attributes listed below:

**Table 55: Attributes for N-CREATE-RQ of Presentation LUT SOP Class**

| Attribute Name         | Tag         | Usage SCU | Supported Values |
|------------------------|-------------|-----------|------------------|
| Presentation LUT Shape | (2050,0020) | U         | IDENTITY         |

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and is used for later requests on the Basic Film Box (N-CREATE-RQ) and on the Presentation LUT (N-DELETE-RQ).

The Presentation LUT SOP class interprets the status codes listed below:

**Table 56: N-CREATE-RSP Status Handling Behavior for the Presentation LUT SOP Class**

| Service Status | Further Meaning                                                                                                                            | Error Codes | Behavior                                     |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------|-------------|----------------------------------------------|
| Warning        | Requested MinDensity or MaxDensity outside of HCD's operating range. HCD will use its respective minimum or maximum density value instead. | B605        | Print job continues and the reason is logged |
| Success        | Presentation LUT successfully created                                                                                                      | 0000        | Print job continues                          |

### Printer SOP Class

The Printer SOP Class provides the possibility to monitor the status of the hardcopy printer in a synchronous and in an asynchronous way.

When used synchronously the *syngo*® MR XA50A Print SCU uses the N-GET-RQ to request information about the printer status. It uses the attributes listed in the table below.

**Table 57: Attributes for N-GET-RQ of the Printer SOP Class**

| Attribute Name      | Tag         | Usage SCP | Supported Values                              |
|---------------------|-------------|-----------|-----------------------------------------------|
| Printer Status      | (2110,0010) | M         | NORMAL, FAILURE, WARNING                      |
| Printer Status Info | (2110,0020) | M         | See table in chapter 9.7 for possible values. |

The command communication failure behavior listed below applies to all SOP classes used for the "Print Film" activity:

**Table 58: DICOM Command Communication Failure Behavior**

| Exception           | Behavior                                                     |
|---------------------|--------------------------------------------------------------|
| Timeout             | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user                                     |

#### 4.2.7.4 Association Acceptance Policy

##### 4.2.7.4.1 Activity Print Film

##### 4.2.7.4.1.1 Description and Sequencing of Activities

The *syngo*® MR XA50A supports the reverse role negotiation of the Printer SOP Class. Receiving the N-EVENT-REPORT-RQ from a printer the *syngo*® MR XA50A is asynchronously informed about changes of the printer status.

##### 4.2.7.4.1.2 Accepted Presentation Context

The *syngo*® MR XA50A accepts Presentation Contexts as shown in the following table:

**Table 59: Presentation Contexts for the Activity "Print Film"**

| Abstract Syntax   |                        | Transfer Syntax           |                   | Role | Ext. Neg. |
|-------------------|------------------------|---------------------------|-------------------|------|-----------|
| Name              | UID                    | Name List                 | UID List          |      |           |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

##### 4.2.7.4.1.3 SOP Specific Conformance

The arguments of the N-EVENT-REPORT-RQ are defined in the table below:

**Table 60: Attributes for the N-EVENT-REPORT-RQ of the Printer SOP Class**

| Event-type Name | Event | Attributes          | Tag         | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal          | 1     | -                   | -           | -         |
| Warning         | 2     | Printer Status Info | (2110,0020) | U         |
| Failure         | 3     | Printer Status Info | (2110,0020) | U         |

##### 4.2.7.4.2 Activity Print Management

##### 4.2.7.4.2.1 Description and Sequencing of Activities

The *syngo*® MR XA50A supports the reverse role negotiation of the Print Job SOP Class. Receiving the N-EVENT-REPORT-RQ from a printer the *syngo*® MR XA50A is asynchronously informed about the status of a print job for monitoring its progress.

##### 4.2.7.4.2.2 Accepted Presentation Context

The *syngo*® MR XA50A accepts Presentation Contexts as shown in the following table:

**Table 61: Presentation Contexts for the Activity "Print Management"**

| Abstract Syntax     |                        | Transfer Syntax           |                   | Role | Ext. Neg. |
|---------------------|------------------------|---------------------------|-------------------|------|-----------|
| Name                | UID                    | Name List                 | UID List          |      |           |
| Print Job SOP Class | 1.2.840.10008.5.1.1.14 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

#### 4.2.7.4.2.3 SOP Specific Conformance

Attributes that can be handled by the Print AE of the product are listed in the table below.

**Table 62: Attributes for the N-EVENT-REPORT-RQ of the Print Job SOP Class**

| Event-type Name | Event | Attributes            | Tag         | Usage SCU                                              |
|-----------------|-------|-----------------------|-------------|--------------------------------------------------------|
| Pending         | 1     | Execution Status Info | (2100,0030) | U                                                      |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U                                                      |
|                 |       | Printer Name          | (2110,0030) | U                                                      |
| Printing        | 2     | Execution Status Info | (2100,0030) | U                                                      |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U                                                      |
|                 |       | Printer Name          | (2110,0030) | U                                                      |
| Done            | 3     | Execution Status Info | (2100,0030) | U                                                      |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U                                                      |
|                 |       | Printer Name          | (2110,0030) | U                                                      |
| Failure         | 4     | Execution Status Info | (2100,0030) | U                                                      |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U                                                      |
|                 |       | Printer Name          | (2110,0030) | U                                                      |

## 4.3 Network Interfaces

### 4.3.1 Physical Network Interface

The syngo® MR XA50A provides DICOM 3.0 TCP/IP network communication support as defined in Part 8 of the DICOM Standard. The network communication is independent from the physical medium over which TCP/IP executes; it inherits this from the Windows OS system upon which it executes.

### 4.3.2 Additional Protocols

None

### 4.3.3 IPv4 and IPv6 Support

IPv4 and IPv6 are supported. Regarding IPv6 please note, that the complete networking infrastructure in the hospital (firewalls, DNS-Servers, ...) must support IPv6 to get a functioning communication.

## 4.4 Configuration

### 4.4.1 AE Title/Presentation Address Mapping

AE Titles shall be unique within the hospital. A common way to achieve that is to use the hostname as part of the AE Titles. The string can be up to 16 characters and must not contain any extended characters. Only 7-bit ASCII characters (excluding Control Characters) are allowed according to the DICOM Standard.

#### 4.4.1.1 Local AE Titles

The *syngo*® MR XA50A allows to configure AE Titles, Ports and Services in any wished way. Default delivery is that all services are using the same AE title and only one port number. In case the connected systems cannot handle this default, the customer service engineer is able to configure for each service its own AE title and Port number.

**Table 63: AE Titles**

| Parameter           | Configurable | Default Value                                              |
|---------------------|--------------|------------------------------------------------------------|
| Default AE title    | Yes          | hostname in uppercase characters; limited to 16 characters |
| Default Basic Port  | Yes          | 104                                                        |
| Default Secure Port | Yes          | 2762                                                       |

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

##### 4.4.1.2.1 Remote Association Initiators

All relevant remote applications that may setup DICOM associations towards *syngo*® MR XA50A need to be configured in *syngo*® MR XA50A, before the association can be established. This behavior is configurable, but it is recommended, not to change this behavior.

The mapping of external AE Titles to TCP/IP addresses and ports is configurable and initially set at the time of installation by Installation Personnel. Changes can later on also be performed by the local system administrator. The Application Entity Titles and supported transfer syntaxes need to be known for configuration.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

##### 4.4.1.2.2 Remote Association Acceptors

For remote applications that shall be able to accept DICOM associations from *syngo*® MR XA50A, the following information needs to be available:

- Application Entity Title
- Host Name / IP address on which the remote application service runs
- Port number on which the remote application accepts association requests.

The remote system will be indicated in the UI of *syngo*® MR XA50A with a logical name, that is also entered when configuring the node in the administration UI.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

#### 4.4.1.3 Secure DICOM Communication

The system supports configuring the DICOM communication to use secure channel (TLS) between *syngo*® MR XA50A and configured remote nodes. As a security measure the certificate thumbprint or certificate trust chain of the remote nodes shall be added (pinned) to the *syngo*® MR XA50A system to authorize the incoming connection.

Detailed instructions how to set up secure DICOM communication are available in the Administrator Online Help.

Note: The default DICOM port will change to 2762.

If the certificate of remote node contains Enhanced Key Usage (Extended Key Usage) field, then:

- If the remote node acts as DICOM SCP it shall contain Server Authentication (1.3.6.1.5.5.7.3.1)
- If the remote node acts as a DICOM SCU it shall contain Client Authentication (1.3.6.1.5.5.7.3.2)

Otherwise *syngo*® MR XA50A will not accept the certificate.

#### 4.4.2 Parameters

The next table lists configuration parameters, which are true for all Application Entities.


**Table 64: Parameter List**

| Parameter                                                                                                  | Configurable | Default Value |
|------------------------------------------------------------------------------------------------------------|--------------|---------------|
| PDU size                                                                                                   | Yes          | 32768 Bytes   |
| time-out for accepting/rejecting an association request                                                    | Yes          | 30 s          |
| time-out for responding to an association open/close request                                               | Yes          | 30 s          |
| time-out for accepting a message over network                                                              | Yes          | 30 s          |
| time-out for waiting for data between TCP/IP-packets                                                       | Yes          | 5 s           |
| time-outs for waiting for a Service Request/Response message from the remote node (Storage SCP/SCU)        | Yes          | 30 s          |
| time-outs for waiting for a Service Request/Response message from the remote node (Query/Retrieve SCP/SCU) | Yes          | 30 s          |
| time-out for waiting for a C-MOVE-RSP                                                                      | No           | 1200 s        |
| number of image collection before saving to database                                                       | Yes          | 20            |
| max matches query limit                                                                                    | Yes          | 100           |
| max number of parallel receiving associations                                                              | Yes          | 12            |

## 5 Media Interchange

### 5.1 Implementation Model

#### 5.1.1 Application Data Flow Diagram



**Figure 6: Media Interchange Application Data Flow Diagram**

The *syngo*® MR XA50A provides the functionality to Import or Export DICOM Instances from and to the File System. During export, a DICOMDIR may also be generated (user selection). A complete ISO Image ready-to-burn can also be generated. All SOP Classes defined in Table 1: Network Services are supported for the Import/Export functionality.

#### 5.1.2 Functional Definitions of AEs

The *syngo*® MR XA50A application is capable of

- creating a new File-set in the File System (Export to ...)
- importing SOP Instances from the medium onto local storage
- writing the File-sets DICOMDIR information into the file system and joining it to an ISO image.

### 5.1.3 Sequencing of Real-World Activities

Whenever data is written to an external media, *syngo*® MR XA50A creates a DICOMDIR from the selected data and creates an ISO image of the selected data on the local hard disk. Depending on the selected data and options (selected media size, with or without compression) either General Purpose CD profile, DVD-J2K, BD-J2K profile is used.




Figure 7: Sequence diagram – Media creation

### 5.1.4 File Meta Information for Implementation Class and Version

This section describes the values assigned to the File Meta Information attributes (see part PS 3.10) that pertain to the Implementation Class and Version. The implementation Class UID and the Implementation Version name in the File Meta Header are the same as the values specified for networking.

Table 65: Implementation Class/Version Name - Media Interchange

|                               |                     |
|-------------------------------|---------------------|
| File Meta Information Version | 0001                |
| Implementation Class UID      | 1.3.12.2.1107.5.8.2 |
| Implementation Version Name   | SYNGO_MR_XA50A      |

## 5.2 AE Specifications

### 5.2.1 Media Storage AE – Specification

The *syngo*® MR XA50A provides conformance to the following Application Profiles as an FSC as well as an FSR. The FSU role is only supported on a non-optical storage device (e.g. USB stick).

In addition, augmented conformance is provided to store extra data objects important for the full feature support of the *syngo*®-based products. Details are listed below:

**Table 66: Media - Application Profiles and Real-World Activities**

| Application Profiles Supported | Real-World Activity                                                                      | Role          | Service Class Option |
|--------------------------------|------------------------------------------------------------------------------------------|---------------|----------------------|
| AUG-GEN-CD                     | Browse Directory Information<br>Import into Application<br>Export to local Archive Media | FSR, FSC      | Interchange          |
| AUG-GEN-DVD                    |                                                                                          |               |                      |
| AUG- GEN-DVD-J2K               |                                                                                          |               |                      |
| AUG-GEN-BD                     |                                                                                          |               |                      |
| AUG-GEN-BD-J2K                 |                                                                                          |               |                      |
| AUG- GEN-USB-J2K               | Browse Directory Information<br>Import into Application<br>Export to local Archive Media | FSR, FSC, FSU | Interchange          |
| STD-GEN-CD                     |                                                                                          |               |                      |
| STD-GEN-DVD                    |                                                                                          |               |                      |
| STD-GEN-DVD-J2K                |                                                                                          |               |                      |
| STD-GEN-BD                     |                                                                                          |               |                      |
| STD-GEN-BD-J2K                 | Browse Directory Information<br>Import into Application<br>Export to local Archive Media | FSR, FSC      | Interchange          |
| STD-GEN-USB-J2K                |                                                                                          |               |                      |

#### 5.2.1.1 Real-World Activities

##### 5.2.1.1.1 Activity “Browse Directory Information”

The *syngo*® MR XA50A acts as FSR using the interchange option when requested to read the media directory.

The *syngo*® MR XA50A will read the DICOMDIR and insert those directory entries that are valid for the application profiles supported, into a local database. The database then is used for browsing media contents.

**Note:** The “Icon Image Sequence” is also supported in DICOMDIR. But only those Icon Images with “Bits Allocated” (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the Browser.

##### 5.2.1.1.2 Real World Activity “Import into Application”

The *syngo*® MR XA50A application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the application.

The SOP Instance selected from the media directory will be copied into the running Application. Only SOP Instances, that are valid for the application profile supported and supported by *syngo*® MR XA50A can be retrieved from media.

### 5.2.1.1.3 Real-World Activity "Export to local Archive Media"

The *syngo*® MR XA50A application acts as FSU (for media with existing DICOM file-set) or FSC (media not initialized) using the interchange option when requested to copy SOP Instances from the local storage to local Archive Medium. The activity as FSU is only possible as long as the local burning SW of *syngo*® MR XA50A has not already processed the generated ISO file.

The *syngo*® MR XA50A application will receive a list of SOP Instances to be copied to the local archive medium. Depending on the profile selected (Standard: uncompressed, with DICOMDIR; Patient: compressed with DICOMDIR) the SOP Instances will be taken, and an ISO file is being generated that includes the DICOMDIR and the corresponding objects.

It is then up to *syngo*® MR XA50A local configuration (if equipped with a local media burner) to burn the ISO file to the appropriate media.

### 5.2.1.2 SOP Classes and Transfer Syntaxes

These Application Profiles are based on the Media Storage Service Class with the Interchange Option. The *syngo*® MR XA50A provides Standard Conformance to the SOP Classes listed in "Table 1: Network Services" section "SOP Classes Created by the *syngo*® MR XA50A" and "SOP Classes Managed by the *syngo*® MR XA50A" in the "[Conformance Statement Overview](#)".

Using the Application Profiles supporting compression (AUG- GEN-DVD-J2K, AUG-GEN-BD-J2K, AUG- GEN-USB-J2K, STD-GEN-DVD-J2K, STD-GEN-USB-J2K, STD-GEN-BD-J2K) the following Transfer Syntaxes are supported:

**Table 67: Transfer Syntaxes for STD-GEN-DVD-J2K, STD-GEN-BD-J2K and STD-GEN-USB-J2K**

| UID value              | Transfer Syntax                                                                          | Image Objects | Non-Image Objects |
|------------------------|------------------------------------------------------------------------------------------|---------------|-------------------|
| 1.2.840.10008.1.2.1    | Explicit VR Little Endian                                                                | Yes           | Yes               |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed                                               | Yes           | No                |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed                                           | Yes           | No                |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | Yes           | No                |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed                                   | Yes           | No                |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed                                             | Yes           | No                |
| 1.2.840.10008.1.2.5    | RLE Lossless compressed                                                                  | Yes           | No                |

Using the Application Profiles that do not support compression (AUG- GEN-DVD, AUG-GEN-BD, AUG- GEN-USB, STD-GEN-DVD, STD-GEN-BD, STD-GEN-USB) only Explicit VR Little Endian (1.2.840.10008.1.2.1) is supported.

## 5.3 Augmented and Private Application Profiles

### 5.3.1 Augmented Application Profiles

N/A

## 5.4 Media Configuration

None

## 6 Support of Extended Character Sets

The syngo® MR XA50A DICOM application supports the following character sets as defined in the four tables below:

**Table 68: Single-Byte Character Sets without Code Extension**

| Character Set Description | Defined Term | ISO registration number | Character Set        |
|---------------------------|--------------|-------------------------|----------------------|
| Default repertoire        | None         | ISO_IR 6                | ISO 646:             |
| Latin alphabet No. 1      | ISO_IR 100   | ISO_IR 100              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646:             |
| Latin alphabet No. 2      | ISO_IR 101   | ISO_IR 101              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 3      | ISO_IR 109   | ISO_IR 109              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 4      | ISO_IR 110   | ISO_IR 110              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Cyrillic                  | ISO_IR 144   | ISO_IR 144              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Arabic                    | ISO_IR 127   | ISO_IR 127              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Greek                     | ISO_IR 126   | ISO_IR 126              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Hebrew                    | ISO_IR 138   | ISO_IR 138              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 5      | ISO_IR 148   | ISO_IR 148              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Japanese                  | ISO_IR 13    | ISO_IR 13               | JIS X 0201: Katakana |
|                           |              | ISO_IR 14               | JIS X 0201: Romaji   |
| Thai                      | ISO_IR166    | ISO_IR166               | TIS 620-253 (1990)   |
|                           |              | ISO_IR 6                | ISO 646              |

Table 69: Single-Byte Characters Sets with Code Extension

| Character Set Description | Defined Term    | Standard for Code Extension | ESC sequence    | ISO registration number | Character Set     |
|---------------------------|-----------------|-----------------------------|-----------------|-------------------------|-------------------|
| Default repertoire        | ISO 2022 IR 6   | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.1       | ISO 2022 IR 100 | ISO 2022                    | ESC 02/13 04/01 | ISO-IR 100              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.2       | ISO 2022 IR 101 | ISO 2022                    | ESC 02/13 04/02 | ISO-IR 101              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.3       | ISO 2022 IR 109 | ISO 2022                    | ESC 02/13 04/03 | ISO-IR 109              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.4       | ISO 2022 IR 110 | ISO 2022                    | ESC 02/13 04/04 | ISO-IR 110              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |

Multi-Byte Character Sets without Code Extension

Table 70: Multi-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set                                         |
|---------------------------|--------------|-------------------------|-------------------------------------------------------|
| Unicode                   | ISO_IR 192   | ISO 10646               | Unicode in UTF-8                                      |
| Chinese                   | GB18030      | GB18030                 | GB 18030-2000 (China Association for Standardization) |

Table 71: Multi-Byte Character Sets with Code Extension

| Character Set Description | Defined Term    | Standard for Code Extension | ESC sequence          | ISO registration number | Character Set                       |
|---------------------------|-----------------|-----------------------------|-----------------------|-------------------------|-------------------------------------|
| Japanese                  | ISO 2022 IR 159 | ISO 2022                    | ESC 02/04 02/08 04/04 | ISO-IR 159              | JIS X 0212: Supplementary Kanji set |
|                           |                 |                             |                       | ISO-IR 87               |                                     |
| Korean                    | ISO 2022 IR 149 | ISO 2022                    | ESC 02/04 02/09 04/03 | ISO-IR 149              | KS X 1001: Hangul and Hanja         |

All Special Character Sets (SCS) listed above are supported for incoming Data. When creating new Instances, the system will use the default SCS (or SCS List) configured on the machine.

When there is a mismatch between the given character set in attribute (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM conform:

- Convert each illegal character to '?'

There are three categories of character sets which have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO\_IR 6, ISO 2022 IR 6, ISO\_IR 100, etc.  
→ encoded in ISO 2022
- ISO\_IR 192 → encoded in UTF-8
- GB18030 → encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO\_IR 192 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains ISO\_IR 192
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies is not DICOM conform. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

The syngo® MR XA50A supports Kanji characters in the byte zones after 74 (79, 7A, 7B and 7C).

## 7 Attribute confidentiality profiles

### 7.1 De-identification

The *syngo*® MR XA50A application can de-identify attributes, when exporting to Media. Three different levels of de-identification are supported:

- Full de-identification
- Reduced de-identification
- Service de-identification

The user needs to select the appropriate de-identification level during export.

Handling public attributes during de-identification:

- Attributes listed in Table 72: Application Level Confidentiality Profile attributes (standard tags) will be de-identified as specified for the different levels. Attributes not listed in the table are not PII / PHI relevant and will not be de-identified.

Handling private attributes during de-identification:

- Full de-identification: private attributes are not included in anonymized studies
- Reduced de-identification: private attributes are not included in anonymized studies except for the tags listed in Table 73
- Service de-identification: all private attributes are included in anonymized studies except for the tags listed in Table 73

In the following table for attributes marked with:

- 'Yes' - data are anonymized
- 'No' - data are kept

**Table 72: Application Level Confidentiality Profile attributes (standard tags)**

| DICOM Tag   | Attribute Name                      | Full | Reduced | Service |
|-------------|-------------------------------------|------|---------|---------|
| (0000,1000) | Affected SOP Instance UID           | Yes  | No      | No      |
| (0000,1001) | Requested SOP Instance UID          | Yes  | No      | No      |
| (0002,0003) | Media Storage SOP Instance UID      | Yes  | No      | No      |
| (0004,1511) | Referenced SOP Instance UID in File | Yes  | No      | No      |
| (0008,0014) | Instance Creator UID                | Yes  | No      | No      |
| (0008,0015) | Instance Coercion DateTime          | Yes  | No      | No      |
| (0008,0018) | SOP Instance UID                    | Yes  | No      | No      |
| (0008,0020) | Study Date                          | Yes  | No      | No      |
| (0008,0021) | Series Date                         | Yes  | No      | No      |
| (0008,0022) | Acquisition Date                    | Yes  | No      | No      |
| (0008,0023) | Content Date                        | Yes  | No      | No      |
| (0008,0024) | Overlay Date                        | Yes  | No      | No      |
| (0008,0025) | Curve Date                          | Yes  | No      | No      |
| (0008,002A) | Acquisition DateTime                | Yes  | No      | No      |
| (0008,0030) | Study Time                          | Yes  | No      | No      |
| (0008,0031) | Series Time                         | Yes  | No      | No      |

| DICOM Tag   | Attribute Name                                  | Full | Reduced | Service |
|-------------|-------------------------------------------------|------|---------|---------|
| (0008,0032) | Acquisition Time                                | Yes  | No      | No      |
| (0008,0033) | Content Time                                    | Yes  | No      | No      |
| (0008,0034) | Overlay Time                                    | Yes  | No      | No      |
| (0008,0035) | Curve Time                                      | Yes  | No      | No      |
| (0008,0050) | Accession Number                                | Yes  | Yes     | No      |
| (0008,0058) | Failed SOP Instance UID List                    | Yes  | No      | No      |
| (0008,0080) | Institution Name                                | Yes  | Yes     | No      |
| (0008,0081) | Institution Address                             | Yes  | Yes     | No      |
| (0008,0082) | Institution Code Sequence                       | Yes  | Yes     | No      |
| (0008,0090) | Referring Physician's Name                      | Yes  | Yes     | Yes     |
| (0008,0092) | Referring Physician's Address                   | Yes  | Yes     | Yes     |
| (0008,0094) | Referring Physician's Telephone Numbers         | Yes  | Yes     | Yes     |
| (0008,0096) | Referring Physician's Identification Sequence   | Yes  | Yes     | No      |
| (0008,010D) | Context Group Extension Creator UID             | Yes  | No      | No      |
| (0008,0201) | Timezone Offset From UTC                        | Yes  | No      | No      |
| (0008,1010) | Station Name                                    | Yes  | Yes     | Yes     |
| (0008,1030) | Study Description                               | Yes  | Yes     | No      |
| (0008,103E) | Series Description                              | Yes  | Yes     | No      |
| (0008,1040) | Institutional Department Name                   | Yes  | Yes     | No      |
| (0008,1048) | Physician(s) of Record                          | Yes  | Yes     | Yes     |
| (0008,1049) | Physician(s) of Record Identification Sequence  | Yes  | Yes     | No      |
| (0008,1050) | Performing Physicians' Name                     | Yes  | Yes     | Yes     |
| (0008,1052) | Performing Physicians' Identification Sequence  | Yes  | Yes     | No      |
| (0008,1060) | Name of Physician(s) Reading Study              | Yes  | Yes     | Yes     |
| (0008,1062) | Physician Reading Study Identification Sequence | Yes  | Yes     | No      |
| (0008,1070) | Operators' Name                                 | Yes  | Yes     | Yes     |
| (0008,1072) | Operators' Identification Sequence              | Yes  | Yes     | No      |
| (0008,1080) | Admitting Diagnoses Description                 | Yes  | Yes     | No      |
| (0008,1084) | Admitting Diagnoses Code Sequence               | Yes  | Yes     | No      |
| (0008,1110) | Referenced Study Sequence                       | Yes  | No      | No      |
| (0008,1111) | Referenced Performed Procedure Step Sequence    | Yes  | No      | No      |
| (0008,1120) | Referenced Patient Sequence                     | Yes  | Yes     | No      |
| (0008,1140) | Referenced Image Sequence                       | Yes  | No      | No      |
| (0008,1155) | Referenced SOP Instance UID                     | Yes  | No      | No      |
| (0008,1195) | Transaction UID                                 | Yes  | No      | No      |
| (0008,2111) | Derivation Description                          | Yes  | No      | No      |
| (0008,2112) | Source Image Sequence                           | Yes  | No      | No      |

| DICOM Tag   | Attribute Name                                    | Full | Reduced | Service |
|-------------|---------------------------------------------------|------|---------|---------|
| (0008,3010) | Irradiation Event UID                             | Yes  | No      | No      |
| (0008,4000) | Identifying Comments                              | Yes  | Yes     | No      |
| (0008,9123) | Creator Version UID                               | Yes  | No      | No      |
| (0010,0010) | Patient's Name                                    | Yes  | Yes     | Yes     |
| (0010,0020) | Patient ID                                        | Yes  | Yes     | Yes     |
| (0010,0021) | Issuer of Patient ID                              | Yes  | Yes     | No      |
| (0010,0030) | Patient's Birth Date                              | Yes  | Yes     | Yes     |
| (0010,0032) | Patient's Birth Time                              | Yes  | Yes     | No      |
| (0010,0040) | Patient's Sex                                     | Yes  | No      | No      |
| (0010,0050) | Patient's Insurance Plan Code Sequence            | Yes  | Yes     | Yes     |
| (0010,0101) | Patient's Primary Language Code Sequence          | Yes  | Yes     | Yes     |
| (0010,0102) | Patient's Primary Language Modifier Code Sequence | Yes  | Yes     | Yes     |
| (0010,1000) | Other Patient IDs                                 | Yes  | Yes     | Yes     |
| (0010,1001) | Other Patient Names                               | Yes  | Yes     | Yes     |
| (0010,1002) | Other Patient IDs Sequence                        | Yes  | Yes     | Yes     |
| (0010,1005) | Patient's Birth Name                              | Yes  | Yes     | Yes     |
| (0010,1010) | Patient's Age                                     | Yes  | No      | No      |
| (0010,1020) | Patient's Size                                    | Yes  | No      | No      |
| (0010,1030) | Patient's Weight                                  | Yes  | No      | No      |
| (0010,1040) | Patient Address                                   | Yes  | Yes     | Yes     |
| (0010,1050) | Insurance Plan Identification                     | Yes  | Yes     | No      |
| (0010,1060) | Patient's Mother's Birth Name                     | Yes  | Yes     | Yes     |
| (0010,1080) | Military Rank                                     | Yes  | Yes     | No      |
| (0010,1081) | Branch of Service                                 | Yes  | Yes     | No      |
| (0010,1090) | Medical Record Locator                            | Yes  | Yes     | No      |
| (0010,1100) | Referenced Patient Photo Sequence                 | Yes  | Yes     | No      |
| (0010,2000) | Medical Alerts                                    | Yes  | Yes     | No      |
| (0010,2110) | Allergies                                         | Yes  | Yes     | No      |
| (0010,2150) | Country of Residence                              | Yes  | Yes     | No      |
| (0010,2152) | Region of Residence                               | Yes  | Yes     | No      |
| (0010,2154) | Patient's Telephone Number                        | Yes  | Yes     | Yes     |
| (0010,2160) | Ethnic Group                                      | Yes  | No      | No      |
| (0010,2180) | Occupation                                        | Yes  | Yes     | No      |
| (0010,21A0) | Smoking Status                                    | Yes  | No      | No      |
| (0010,21B0) | Additional Patient's History                      | Yes  | Yes     | Yes     |
| (0010,21C0) | Pregnancy Status                                  | Yes  | No      | No      |
| (0010,21D0) | Last Menstrual Date                               | Yes  | No      | No      |

| DICOM Tag   | Attribute Name                            | Full | Reduced | Service |
|-------------|-------------------------------------------|------|---------|---------|
| (0010,21F0) | Patient's Religious Preference            | Yes  | Yes     | No      |
| (0010,2203) | Patient Sex Neutered                      | Yes  | No      | No      |
| (0010,2297) | Responsible Person                        | Yes  | Yes     | No      |
| (0010,2299) | Responsible Organization                  | Yes  | Yes     | No      |
| (0010,4000) | Patient Comments                          | Yes  | Yes     | Yes     |
| (0018,0010) | Contrast Bolus Agent                      | Yes  | Yes     | No      |
| (0018,1000) | Device Serial Number                      | Yes  | Yes     | No      |
| (0018,1002) | Device UID                                | Yes  | No      | No      |
| (0018,1004) | Plate ID                                  | Yes  | Yes     | No      |
| (0018,1005) | Generator ID                              | Yes  | Yes     | No      |
| (0018,1007) | Cassette ID                               | Yes  | Yes     | No      |
| (0018,1008) | Gantry ID                                 | Yes  | Yes     | No      |
| (0018,1030) | Protocol Name                             | Yes  | Yes     | No      |
| (0018,1400) | Acquisition Device Processing Description | Yes  | Yes     | No      |
| (0018,2042) | Target UID                                | Yes  | No      | No      |
| (0018,4000) | Acquisition Comments                      | Yes  | Yes     | No      |
| (0018,700A) | Detector ID                               | Yes  | Yes     | No      |
| (0018,9424) | Acquisition Protocol Description          | Yes  | Yes     | No      |
| (0018,9516) | Start Acquisition DateTime                | Yes  | No      | No      |
| (0018,9517) | End Acquisition DateTime                  | Yes  | No      | No      |
| (0018,A003) | Contribution Description                  | Yes  | Yes     | No      |
| (0020,000D) | Study Instance UID                        | Yes  | No      | No      |
| (0020,000E) | Series Instance UID                       | Yes  | No      | No      |
| (0020,0010) | Study ID                                  | Yes  | Yes     | No      |
| (0020,0052) | Frame of Reference UID                    | Yes  | No      | No      |
| (0020,0200) | Synchronization Frame of Reference UID    | Yes  | No      | No      |
| (0020,3401) | Modifying Device ID                       | Yes  | Yes     | No      |
| (0020,3404) | Modifying Device Manufacturer             | Yes  | Yes     | No      |
| (0020,3406) | Modified Image Description                | Yes  | Yes     | No      |
| (0020,4000) | Image Comments                            | Yes  | Yes     | No      |
| (0020,9158) | Frame Comments                            | Yes  | Yes     | No      |
| (0020,9161) | Concatenation UID                         | Yes  | No      | No      |
| (0020,9164) | Dimension Organization UID                | Yes  | No      | No      |
| (0028,1199) | Palette Color Lookup Table UID            | Yes  | No      | No      |
| (0028,1214) | Large Palette Color Lookup Table UID      | Yes  | No      | No      |
| (0028,4000) | Image Presentation Comments               | Yes  | Yes     | No      |
| (0032,0012) | Study ID Issuer                           | Yes  | Yes     | No      |

| DICOM Tag   | Attribute Name                                         | Full | Reduced | Service |
|-------------|--------------------------------------------------------|------|---------|---------|
| (0032,1020) | Scheduled Study Location                               | Yes  | Yes     | No      |
| (0032,1021) | Scheduled Study Location AE Title                      | Yes  | Yes     | No      |
| (0032,1030) | Reason for Study                                       | Yes  | Yes     | No      |
| (0032,1032) | Requesting Physician                                   | Yes  | Yes     | No      |
| (0032,1033) | Requesting Service                                     | Yes  | Yes     | No      |
| (0032,1060) | Requested Procedure Description                        | Yes  | Yes     | No      |
| (0032,1070) | Requested Contrast Agent                               | Yes  | Yes     | No      |
| (0032,4000) | Study Comments                                         | Yes  | Yes     | No      |
| (0038,0004) | Referenced Patient Alias Sequence                      | Yes  | Yes     | No      |
| (0038,0010) | Admission ID                                           | Yes  | Yes     | No      |
| (0038,0011) | Issuer of Admission ID                                 | Yes  | Yes     | No      |
| (0038,001E) | Scheduled Patient Institution Residence                | Yes  | Yes     | No      |
| (0038,0020) | Admitting Date                                         | Yes  | No      | No      |
| (0038,0021) | Admitting Time                                         | Yes  | No      | No      |
| (0038,0040) | Discharge Diagnosis Description                        | Yes  | Yes     | No      |
| (0038,0050) | Special Needs                                          | Yes  | Yes     | No      |
| (0038,0060) | Service Episode ID                                     | Yes  | Yes     | No      |
| (0038,0061) | Issuer of Service Episode ID                           | Yes  | Yes     | No      |
| (0038,0062) | Service Episode Description                            | Yes  | Yes     | No      |
| (0038,0300) | Current Patient Location                               | Yes  | Yes     | No      |
| (0038,0400) | Patient's Institution Residence                        | Yes  | Yes     | No      |
| (0038,0500) | Patient State                                          | Yes  | Yes     | No      |
| (0038,4000) | Visit Comments                                         | Yes  | Yes     | No      |
| (0040,0001) | Scheduled Station AE Title                             | Yes  | Yes     | No      |
| (0040,0002) | Scheduled Procedure Step Start Date                    | Yes  | No      | No      |
| (0040,0003) | Scheduled Procedure Step Start Time                    | Yes  | No      | No      |
| (0040,0004) | Scheduled Procedure Step End Date                      | Yes  | No      | No      |
| (0040,0005) | Scheduled Procedure Step End Time                      | Yes  | No      | No      |
| (0040,0006) | Scheduled Performing Physician Name                    | Yes  | Yes     | No      |
| (0040,0007) | Scheduled Procedure Step Description                   | Yes  | Yes     | No      |
| (0040,000B) | Scheduled Performing Physician Identification Sequence | Yes  | Yes     | No      |
| (0040,0010) | Scheduled Station Name                                 | Yes  | Yes     | No      |
| (0040,0011) | Scheduled Procedure Step Location                      | Yes  | Yes     | No      |
| (0040,0012) | Pre-Medication                                         | Yes  | Yes     | No      |
| (0040,0241) | Performed Station AE Title                             | Yes  | Yes     | No      |
| (0040,0242) | Performed Station Name                                 | Yes  | Yes     | No      |
| (0040,0243) | Performed Location                                     | Yes  | Yes     | No      |

| DICOM Tag   | Attribute Name                                                      | Full | Reduced | Service |
|-------------|---------------------------------------------------------------------|------|---------|---------|
| (0040,0244) | Performed Procedure Step Start Date                                 | Yes  | No      | No      |
| (0040,0245) | Performed Procedure Step Start Time                                 | Yes  | No      | No      |
| (0040,0250) | Performed Procedure Step End Date                                   | Yes  | No      | No      |
| (0040,0251) | Performed Procedure Step End Time                                   | Yes  | No      | No      |
| (0040,0253) | Performed Procedure Step ID                                         | Yes  | Yes     | No      |
| (0040,0254) | Performed Procedure Step Description                                | Yes  | Yes     | No      |
| (0040,0275) | Request Attributes Sequence                                         | Yes  | Yes     | No      |
| (0040,0280) | Comments on Performed Procedure Step                                | Yes  | Yes     | No      |
| (0040,0555) | Acquisition Context Sequence                                        | Yes  | Yes     | No      |
| (0040,1001) | Requested Procedure ID                                              | Yes  | Yes     | No      |
| (0040,1004) | Patient Transport Arrangements                                      | Yes  | Yes     | No      |
| (0040,1005) | Requested Procedure Location                                        | Yes  | Yes     | No      |
| (0040,1010) | Names of Intended Recipient of Results                              | Yes  | Yes     | No      |
| (0040,1011) | Intended Recipients of Results Identification Sequence              | Yes  | Yes     | No      |
| (0040,1101) | Person Identification Code Sequence                                 | Yes  | Yes     | No      |
| (0040,1102) | Person Address                                                      | Yes  | Yes     | No      |
| (0040,1103) | Person Telephone Numbers                                            | Yes  | Yes     | No      |
| (0040,1400) | Requested Procedure Comments                                        | Yes  | Yes     | No      |
| (0040,2001) | Reason for Imaging Service Request                                  | Yes  | Yes     | No      |
| (0040,2008) | Order Entered By                                                    | Yes  | Yes     | No      |
| (0040,2009) | Order Enterer Location                                              | Yes  | Yes     | No      |
| (0040,2010) | Order Callback Phone Number                                         | Yes  | Yes     | No      |
| (0040,2016) | Placer Order Number of Imaging Service Request                      | Yes  | Yes     | No      |
| (0040,2017) | Filler Order Number of Imaging Service Request                      | Yes  | Yes     | No      |
| (0040,2400) | Imaging Service Request Comments                                    | Yes  | Yes     | No      |
| (0040,3001) | Confidentiality Constraint on Patient Data Description              | Yes  | Yes     | No      |
| (0040,4005) | Scheduled Procedure Step Start DateTime                             | Yes  | No      | No      |
| (0040,4010) | Scheduled Procedure Step Modification DateTime                      | Yes  | No      | No      |
| (0040,4011) | Expected Completion Date Time                                       | Yes  | No      | No      |
| (0040,4023) | Referenced General Purpose Scheduled Procedure Step Transaction UID | Yes  | No      | No      |
| (0040,4025) | Scheduled Station Name Code Sequence                                | Yes  | Yes     | No      |
| (0040,4027) | Scheduled Station Geographic Location Code Sequence                 | Yes  | Yes     | No      |
| (0040,4028) | Performed Station Name Code Sequence                                | Yes  | Yes     | No      |
| (0040,4030) | Performed Station Geographic Location Code Sequence                 | Yes  | Yes     | No      |
| (0040,4034) | Scheduled Human Performers Sequence                                 | Yes  | Yes     | No      |
| (0040,4035) | Actual Human Performers Sequence                                    | Yes  | Yes     | No      |
| (0040,4036) | Human Performers Organization                                       | Yes  | Yes     | No      |

| DICOM Tag   | Attribute Name                                  | Full | Reduced | Service |
|-------------|-------------------------------------------------|------|---------|---------|
| (0040,4037) | Human Performers Name                           | Yes  | Yes     | No      |
| (0040,4050) | Performed Procedure Step Start DateTime         | Yes  | No      | No      |
| (0040,4051) | Performed Procedure Step End DateTime           | Yes  | No      | No      |
| (0040,4052) | Procedure Step Cancellation DateTime            | Yes  | No      | No      |
| (0040,A027) | Verifying Organization                          | Yes  | Yes     | No      |
| (0040,A073) | Verifying Observer Sequence                     | Yes  | Yes     | No      |
| (0040,A075) | Verifying Observer Name                         | Yes  | Yes     | No      |
| (0040,A078) | Author Observer Sequence                        | Yes  | Yes     | No      |
| (0040,A07A) | Participant Sequence                            | Yes  | Yes     | No      |
| (0040,A07C) | Custodial Organization Sequence                 | Yes  | Yes     | No      |
| (0040,A088) | Verifying Observer Identification Code Sequence | Yes  | Yes     | No      |
| (0040,A123) | Person Name                                     | Yes  | Yes     | No      |
| (0040,A124) | UID                                             | Yes  | Yes     | No      |
| (0040,A171) | Observation UID                                 | Yes  | No      | No      |
| (0040,A172) | Referenced Observation UID (Trial)              | Yes  | No      | No      |
| (0040,A192) | Observation Date (Trial)                        | Yes  | No      | No      |
| (0040,A193) | Observation Time (Trial)                        | Yes  | No      | No      |
| (0040,A307) | Current Observer (Trial)                        | Yes  | Yes     | No      |
| (0040,A352) | Verbal Source (Trial)                           | Yes  | Yes     | No      |
| (0040,A353) | Address (Trial)                                 | Yes  | Yes     | No      |
| (0040,A354) | Telephone Number (Trial)                        | Yes  | Yes     | Yes     |
| (0040,A358) | Verbal Source Identifier Code Sequence (Trial)  | Yes  | Yes     | No      |
| (0040,A402) | Observation Subject UID (Trial)                 | Yes  | No      | No      |
| (0040,A730) | Content Sequence                                | Yes  | Yes     | No      |
| (0040,DB0C) | Template Extension Organization UID             | Yes  | No      | No      |
| (0040,DB0D) | Template Extension Creator UID                  | Yes  | No      | No      |
| (0070,0001) | Graphic Annotation Sequence                     | Yes  | Yes     | No      |
| (0070,0084) | Content Creator's Name                          | Yes  | Yes     | No      |
| (0070,0086) | Content Creator's Identification Code Sequence  | Yes  | Yes     | No      |
| (0070,031A) | Fiducial UID                                    | Yes  | No      | No      |
| (0088,0140) | Storage Media Fileset UID                       | Yes  | No      | No      |
| (0088,0200) | Icon Image Sequence                             | Yes  | Yes     | No      |
| (0088,0904) | Topic Title                                     | Yes  | Yes     | No      |
| (0088,0906) | Topic Subject                                   | Yes  | Yes     | No      |
| (0088,0910) | Topic Author                                    | Yes  | Yes     | No      |
| (0088,0912) | Topic Keywords                                  | Yes  | Yes     | No      |
| (0400,0100) | Digital Signature UID                           | Yes  | Yes     | No      |

| DICOM Tag   | Attribute Name                        | Full | Reduced | Service |
|-------------|---------------------------------------|------|---------|---------|
| (0400,0402) | Referenced Digital Signature Sequence | Yes  | Yes     | No      |
| (0400,0403) | Referenced SOP Instance MAC Sequence  | Yes  | Yes     | No      |
| (0400,0404) | MAC                                   | Yes  | Yes     | No      |
| (0400,0550) | Modified Attributes Sequence          | Yes  | Yes     | No      |
| (0400,0561) | Original Attributes Sequence          | Yes  | Yes     | Yes     |
| (2030,0020) | Text String                           | Yes  | Yes     | No      |
| (3006,0024) | Referenced Frame of Reference UID     | Yes  | No      | No      |
| (3006,00C2) | Related Frame of Reference UID        | Yes  | No      | No      |
| (3008,0105) | Source Serial Number                  | Yes  | No      | No      |
| (300A,0013) | Dose Reference UID                    | Yes  | No      | No      |
| (300E,0008) | Reviewer Name                         | Yes  | Yes     | No      |
| (4000,0010) | Arbitrary                             | Yes  | Yes     | No      |
| (4000,4000) | Text Comments                         | Yes  | Yes     | No      |
| (4008,0042) | Results ID Issuer                     | Yes  | Yes     | No      |
| (4008,0102) | Interpretation Recorder               | Yes  | Yes     | No      |
| (4008,010A) | Interpretation Transcriber            | Yes  | Yes     | No      |
| (4008,010B) | Interpretation Text                   | Yes  | Yes     | No      |
| (4008,010C) | Interpretation Author                 | Yes  | Yes     | No      |
| (4008,0111) | Interpretation Approver Sequence      | Yes  | Yes     | No      |
| (4008,0114) | Physician Approving Interpretation    | Yes  | Yes     | No      |
| (4008,0115) | Interpretation Diagnosis Description  | Yes  | Yes     | No      |
| (4008,0118) | Results Distribution List Sequence    | Yes  | Yes     | No      |
| (4008,0119) | Distribution Name                     | Yes  | Yes     | No      |
| (4008,011A) | Distribution Address                  | Yes  | Yes     | No      |
| (4008,0202) | Interpretation ID Issuer              | Yes  | Yes     | No      |
| (4008,0300) | Impressions                           | Yes  | Yes     | No      |
| (4008,4000) | Results Comments                      | Yes  | Yes     | No      |
| (50xx,xxxx) | Curve Data                            | Yes  | Yes     | No      |
| (60xx,0100) | Overlay Bits Allocated                | Yes  | Yes     | No      |
| (60xx,0102) | Overlay Bit Position                  | Yes  | Yes     | No      |
| (60xx,3000) | Overlay Data                          | Yes  | Yes     | No      |
| (60xx,4000) | Overlay Comments                      | Yes  | Yes     | No      |
| (FFFA,FFFA) | Digital Signatures Sequence           | Yes  | Yes     | Yes     |
| (FFFC,FFFC) | Data Set Trailing Padding             | Yes  | Yes     | Yes     |

Table 73: Application Level Confidentiality Profile Attributes (private tags)

| DICOM Tag                         | Attribute Name                     | Full | Reduced | Service |
|-----------------------------------|------------------------------------|------|---------|---------|
| (0019, SIEMENS CT VAO COAD, 90)   | Osteo offset                       | Yes  | No      | No      |
| (0019, SIEMENS CT VAO COAD, 92)   | Osteo Regression Line Slope        | Yes  | No      | No      |
| (0019, SIEMENS CT VAO COAD, 93)   | Osteo Regression Line Intercept    | Yes  | No      | No      |
| (0019, SIEMENS CT VAO COAD, 96)   | Osteo Phantom Number               | Yes  | No      | No      |
| (0019, SIEMENS MED NM, 93)        | Phase start time                   | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A1)        | Number of Phases                   | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A5)        | Number of repeats / phase          | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A6)        | Cycles Per Repeat                  | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A7)        | Repeat Start time                  | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A8)        | Repeat Stop time                   | Yes  | No      | No      |
| (0019, SIEMENS MED NM, A9)        | Effective Repeat Time              | Yes  | No      | No      |
| (0019, SIEMENS MED NM, AA)        | Acquired Cycles Per Repeat         | Yes  | No      | No      |
| (0033, SIEMENS MED NM, 29)        | Crystal Thickness                  | Yes  | No      | No      |
| (0033, SIEMENS MED NM, 30)        | Preset Name Used for Acquisition   | Yes  | No      | No      |
| (0033, SIEMENS MED NM, 38)        | Pixel Scale factor                 | Yes  | No      | No      |
| (0035, SIEMENS MED NM, 00)        | Specialized TOMO Type              | Yes  | No      | No      |
| (0035, SIEMENS MED NM, 04)        | Repeat ID                          | Yes  | No      | No      |
| (0035, SIEMENS MED NM, 05)        | Phase ID                           | Yes  | No      | No      |
| (0041, SIEMENS MED NM, 01)        | WholeBody Tomo Position Index      | Yes  | No      | No      |
| (0041, SIEMENS MED NM, 02)        | WholeBody Tomo Number of Positions | Yes  | No      | No      |
| (0041, SIEMENS MED NM, 10)        | Effective Emission Energy          | Yes  | No      | No      |
| (0057, SIEMENS MED NM, 03)        | NM Pixel Units                     | Yes  | No      | No      |
| (0061, SIEMENS MED NM, 62)        | Recon Output Type                  | Yes  | No      | No      |
| (0061, SIEMENS MED NM, 70)        | NM Reconstruction Algorithm        | Yes  | No      | No      |
| (0061, SIEMENS MED NM, 8D)        | QSPECT Flag                        | Yes  | No      | No      |
| (0065, SIEMENS MED NM, 01)        | Original Detector Index            | Yes  | No      | No      |
| (0065, SIEMENS MED NM, 02)        | Siemens Planar Data Organization   | Yes  | No      | No      |
| (7FE3, SIEMENS MED NM, 14)        | Minimum pixel value in frame       | Yes  | No      | No      |
| (7FE3, SIEMENS MED NM, 15)        | Maximum pixel value in frame       | Yes  | No      | No      |
| (7FE3, SIEMENS MED NM, 29)        | Number of Rwaves in a frame        | Yes  | No      | No      |
| (0021, SIEMENS MR SDS 01, 19)     | MR Phoenix Protocol                | Yes  | No      | No      |
| (0029, SIEMENS CT EXAM IMAGE, 49) | Metal Artifact Reduction Type      | Yes  | No      | No      |
| (0029, SIEMENS CSA ENVELOPE, 10)  | Syngo Report Data                  | Yes  | No      | No      |
| (0029, SIEMENS CSA ENVELOPE, 11)  | Syngo Report Presentation          | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 08)    | Modality Image Header Type         | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 09)    | Modality Image Header Version      | Yes  | No      | No      |

| DICOM Tag                              | Attribute Name                           | Full | Reduced | Service |
|----------------------------------------|------------------------------------------|------|---------|---------|
| (0029, SIEMENS CSA HEADER, 10)         | Modality Image Header Info               | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 18)         | Modality Series Header Type              | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 19)         | Modality Series Header Version           | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 20)         | Modality Series Header Info              | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 40)      | Application Header Sequence              | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 41)      | Application Header Type                  | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 42)      | Application Header ID                    | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 43)      | Application Header Version               | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 44)      | Application Header Info                  | Yes  | No      | No      |
| (0029, SIEMENS CT APPL DATASET, 00)    | Dual Energy Algorithm Parameters         | Yes  | No      | No      |
| (0029, SIEMENS CT APPL ALG PARAMS, 20) | Perfusion Result Set Id                  | Yes  | No      | No      |
| (0029, SIEMENS CSA REPORT, 08)         | <i>syngo</i> Report Type                 | Yes  | No      | No      |
| (0029, SIEMENS CSA REPORT, 09)         | <i>syngo</i> Report Version              | Yes  | No      | No      |
| (0029, SIEMENS CSA REPORT, 15)         | SR Variant                               | Yes  | No      | No      |
| (0029, SIEMENS CSA REPORT, 17)         | SC SOP Instance UID                      | Yes  | No      | No      |
| (0043, GEMS_PARM_01, 1E)               | GE Delta Start Time                      | Yes  | No      | No      |
| (0049, SIEMENS CT SPP HEADER, 10)      | Raw Data Container                       | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 02)             | Scanner Console Generation               | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 03)             | Recon Parameters                         | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 05)             | Device IVK                               | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 14)             | Raw Data Description                     | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 16)             | Raw Data Series Instance UIDs            | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 17)             | Raw Data Referenced Series Instance UIDs | Yes  | No      | No      |
| (0067, SIEMENS MED MI, 18)             | Raw Data Blob Sequence                   | Yes  | No      | No      |
| (0071, SIEMENS MED PT, 22)             | Decay Correction DateTime                | Yes  | No      | No      |
| (0071, SIEMENS MED PT, 23)             | Registration Matrix                      | Yes  | No      | No      |
| (0071, SIEMENS MED PT, 24)             | Table Motion                             | Yes  | No      | No      |
| (0071, SIEMENS MED PT, 25)             | Lumped Constant                          | Yes  | No      | No      |
| (0071, SIEMENS MED PT, 26)             | Histogramming Method                     | Yes  | No      | No      |
| (0071, SIEMENS MED PT MU MAP, 01)      | SOP Class of Source                      | Yes  | No      | No      |
| (0071, SIEMENS MED PT MU MAP, 02)      | Related Mu Map Series                    | Yes  | No      | No      |

## 8 Security

### 8.1 Security Profiles

#### 8.1.1 Time Synchronization Profiles

Time Synchronization Profiles: The *syngo*® MR XA50A acts as an NTP Client in the Maintain Time Transaction.

#### 8.1.2 Basic TLS Secure Transport Connection Profile

Basic TLS Secure Transport Connection Profile supports TLS version 1.0, 1.1 and 1.2 protocols with the following features:

**Table 74: Secure Transport Connection Profiles**

| Supported TLS Feature      | Mechanism              |
|----------------------------|------------------------|
| Entity Authentication      | RSA based certificates |
| Exchange of Master Secrets | RSA                    |
| Data Integrity             | SHA                    |
| Privacy                    | Triple DES EDE, CBC    |

The default secure DICOM port is 2762 (can be reconfigured).

### 8.2 Association Level Security

It is possible to configure whether the SCP will only answer to known AETs or to any AET.

### 8.3 Application Level Security

- User must login with own password
- For configuration and Maintenance, Service Technician must login with a separate password.

## 9 Annexes

### 9.1 IOD Contents

#### 9.1.1 Created SOP Classes

##### 9.1.1.1 Enhanced MR Image IOD

This chapter describes the DICOM attributes of Enhanced MR Image Instances performed by the MR acquisition.

| IE                 | Module                                  | Reference              | Presence of Module                |
|--------------------|-----------------------------------------|------------------------|-----------------------------------|
| Patient            | Patient                                 | see chapter 9.1.1.1.1  | always                            |
| Study              | General Study                           | see chapter 9.1.1.1.2  | always                            |
|                    | Patient Study                           | see chapter 9.1.1.1.3  | always                            |
| Series             | General Series                          | see chapter 9.1.1.1.4  | always                            |
|                    | MR Series                               | see chapter 9.1.1.1.5  | always                            |
| Frame of Reference | Frame of Reference                      | see chapter 9.1.1.1.6  | always                            |
| Equipment          | General Equipment                       | see chapter 9.1.1.1.7  | always                            |
|                    | Enhanced General Equipment              | see chapter 9.1.1.1.8  | always                            |
| Image              | Image Pixel                             | see chapter 9.1.1.1.9  | always                            |
|                    | Enhanced Contrast/Bolus                 | see chapter 9.1.1.1.10 | always                            |
|                    | Multi-frame Functional Groups           | see chapter 9.1.1.1.11 | always                            |
|                    | Multi-frame Dimension                   | see chapter 9.1.1.1.12 | always                            |
|                    | Cardiac Synchronization                 | see chapter 9.1.1.1.14 | in case of cardiac triggering     |
|                    | Respiratory Synchronization             | see chapter 9.1.1.1.15 | in case of respiratory triggering |
|                    | Bulk Motion Synchronization             | see chapter 9.1.1.1.16 | in case of bulk motion triggering |
|                    | Supplemental Palette Color Lookup Table | see chapter 9.1.1.1.17 | in case of Color LUTs             |
|                    | Acquisition Context                     | see chapter 9.1.1.1.18 | always                            |
|                    | Enhanced MR Image                       | see chapter 9.1.1.1.19 | always                            |
|                    | MR Pulse Sequence                       | see chapter 9.1.1.1.20 | always                            |
|                    | SOP Common                              | see chapter 9.1.1.1.22 | always                            |
|                    | Common Instance Reference               | see chapter 9.1.1.1.23 | always                            |

##### 9.1.1.1.1 Patient Module

| Attribute Name | Tag         | Supported Values              |
|----------------|-------------|-------------------------------|
| Patient's Name | (0010,0010) | RIS defined or set by creator |

| Attribute Name           | Tag         | Supported Values              |
|--------------------------|-------------|-------------------------------|
| Patient ID               | (0010,0020) | RIS defined or set by creator |
| Issuer of Patient ID     | (0010,0021) | RIS defined                   |
| Patient's Birth Date     | (0010,0030) | RIS defined or set by creator |
| Patient's Birth Time     | (0010,0032) | RIS defined or set by creator |
| Patient's Sex            | (0010,0040) | RIS defined or set by creator |
| Other Patient IDs        | (0010,1000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Ethnic Group             | (0010,2160) | RIS defined or set by creator |
| Patient Comments         | (0010,4000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Patient Identity Removed | (0012,0062) | set by creator                |

#### 9.1.1.1.2 General Study Module

| Attribute Name             | Tag         | Supported Values              |
|----------------------------|-------------|-------------------------------|
| Study Instance UID         | (0020,000D) | RIS defined or set by creator |
| Study Date                 | (0008,0020) | set by creator                |
| Study Time                 | (0008,0030) | set by creator                |
| Referring Physician's Name | (0008,0090) | RIS defined or set by creator |
| Requesting Physician       | (0032,1032) | RIS defined or set by creator |
| Study ID                   | (0020,0010) | set by creator                |
| Accession Number           | (0008,0050) | RIS defined or set by creator |
| Study Description          | (0008,1030) | set by creator                |
| Procedure Code Sequence    | (0008,1032) | RIS defined                   |
| >Code Value                | (0008,0100) | RIS defined                   |
| >Code Scheme Designator    | (0008,0102) | RIS defined                   |
| >Code Meaning              | (0008,0104) | RIS defined                   |

#### 9.1.1.1.3 Patient Study Module

| Attribute Name   | Tag         | Supported Values              |
|------------------|-------------|-------------------------------|
| Patient's Age    | (0010,1010) | RIS defined or set by creator |
| Patient's Size   | (0010,1020) | RIS defined or set by creator |
| Patient's Weight | (0010,1030) | RIS defined or set by creator |

#### 9.1.1.1.4 General Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | MR               |

| Attribute Name                               | Tag         | Supported Values              |
|----------------------------------------------|-------------|-------------------------------|
| Series Instance UID                          | (0020,000E) | set by creator                |
| Series Number                                | (0020,0011) | set by creator                |
| Series Date                                  | (0008,0021) | set by creator                |
| Series Time                                  | (0008,0031) | set by creator                |
| Performing Physicians' Name                  | (0008,1050) | RIS defined or set by creator |
| Protocol Name                                | (0018,1030) | set by creator                |
| Series Description                           | (0008,103E) | set by creator                |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator                |
| >Referenced SOP Class UID                    | (0008,1150) | set by creator                |
| >Referenced SOP Instance UID                 | (0008,1155) | set by creator                |
| Body Part Examined                           | (0018,0015) | set by creator                |
| Patient Position                             | (0018,5100) | set by creator                |
| Request Attributes Sequence                  | (0040,0275) | RIS defined                   |
| >Requested Procedure ID                      | (0040,1001) | RIS defined                   |
| >Accession Number                            | (0008,0050) | RIS defined                   |
| >Study Instance UID                          | (0020,000D) | RIS defined                   |
| >Requested Procedure Description             | (0032,1060) | RIS defined                   |
| >Scheduled Procedure Step ID                 | (0040,0009) | RIS defined                   |
| >Scheduled Procedure Step Description        | (0040,0007) | RIS defined                   |
| > Referenced Study Sequence                  | (0008,1110) | RIS defined                   |
| >> Referenced SOP Class UID                  | (0008,1150) | RIS defined                   |
| >> Referenced SOP Instance UID               | (0008,1155) | RIS defined                   |
| > Requested Procedure Code Sequence          | (0032,1064) | RIS defined                   |
| >> Code Value                                | (0008,0100) | RIS defined                   |
| >> Code Scheme Designator                    | (0008,0102) | RIS defined                   |
| >> Code Meaning                              | (0008,0104) | RIS defined                   |
| > Scheduled Protocol Code Sequence           | (0040,0008) | RIS defined                   |
| >> Code Value                                | (0008,0100) | RIS defined                   |
| >> Code Scheme Designator                    | (0008,0102) | RIS defined                   |
| >> Code Meaning                              | (0008,0104) | RIS defined                   |
| Performed Procedure Step ID                  | (0040,0253) | set by creator                |
| Performed Procedure Step Start Date          | (0040,0244) | set by creator                |
| Performed Procedure Step Start Time          | (0040,0245) | set by creator                |
| Performed Procedure Step Description         | (0040,0254) | set by creator                |

#### 9.1.1.1.5 MR Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | MR               |

## 9.1.1.1.6 Frame of Reference Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Frame of Reference UID       | (0020,0052) | set by creator   |
| Position Reference Indicator | (0020,1040) | empty            |

## 9.1.1.1.7 General Equipment Module

| Attribute Name                | Tag         | Supported Values              |
|-------------------------------|-------------|-------------------------------|
| Manufacturer                  | (0008,0070) | Siemens Healthineers          |
| Institution Name              | (0008,0080) | RIS defined or set by creator |
| Station Name                  | (0008,1010) | set by creator                |
| Institution Address           | (0008,0081) | RIS defined set by creator    |
| Institutional Department Name | (0008,1040) | set by creator                |
| Manufacturer's Model Name     | (0008,1090) | set by creator                |
| Device Serial Number          | (0018,1000) | set by creator                |
| Software Versions             | (0018,1020) | set by creator                |

## 9.1.1.1.8 Enhanced General Equipment Module

| Attribute Name | Tag         | Supported Values     |
|----------------|-------------|----------------------|
| Manufacturer   | (0008,0070) | Siemens Healthineers |

## 9.1.1.1.9 Image Pixel Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel          | (0028,0002) | 1                |
| Photometric Interpretation | (0028,0004) | MONOCHROME2      |
| Rows                       | (0028,0010) | set by creator   |
| Columns                    | (0028,0011) | set by creator   |
| Bits Allocated             | (0028,0100) | 16               |
| Bits Stored                | (0028,0101) | 12 or 16         |
| High Bit                   | (0028,0102) | 11 or 15         |
| Pixel Representation       | (0028,0103) | 0                |
| Pixel Data                 | (7FE0,0010) | set by creator   |
| Smallest Image Pixel Value | (0028,0106) | set by creator   |
| Largest Image Pixel Value  | (0028,0107) | set by creator   |
| Pixel Data                 | (7FE0,0010) | set by creator   |

## 9.1.1.1.10 Enhanced Contrast/Bolus Module

| Attribute Name                                | Tag         | Supported Values |
|-----------------------------------------------|-------------|------------------|
| Contrast/Bolus Agent Sequence                 | (0018,0012) | set by creator   |
| >Code Value                                   | (0008,0100) | set by creator   |
| >Code Scheme Designator                       | (0008,0102) | set by creator   |
| >Code Meaning                                 | (0008,0104) | set by creator   |
| >Contrast/Bolus Agent Number                  | (0018,9337) | 1                |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Ingredient Code Sequence      | (0018,9338) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Volume                        | (0018,1041) | set by creator   |
| >Contrast/Bolus Ingredient Concentration      | (0018,1049) | set by creator   |

## 9.1.1.1.11 Multi-frame Functional Groups Module

| Attribute Name   | Tag         | Supported Values |
|------------------|-------------|------------------|
| Instance Number  | (0020,0013) | set by creator   |
| Content Date     | (0008,0023) | set by creator   |
| Content Time     | (0008,0033) | set by creator   |
| Number of Frames | (0028,0008) | set by creator   |

## 9.1.1.1.12 Multi-frame Dimension Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Dimension Organization Sequence | (0020,9221) | set by creator   |
| Dimension Index Sequence        | (0020,9222) | set by creator   |

## 9.1.1.1.13 Enhanced MR Image Functional Groups

| Functional Group<br>Macro | Attribute Name            | Tag         | Supported Values |
|---------------------------|---------------------------|-------------|------------------|
| Pixel Measures            | Pixel Measures Sequence   | (0028,9110) | set by creator   |
|                           | >Pixel Spacing            | (0018,0030) | set by creator   |
|                           | >Slice Thickness          | (0028,0050) | set by creator   |
| Frame Content             | Frame Content Sequence    | (0020,9111) | set by creator   |
|                           | >Frame Acquisition Number | (0020,9156) | set by creator   |

| Functional Group<br>Macro  | Attribute Name                                | Tag         | Supported Values |
|----------------------------|-----------------------------------------------|-------------|------------------|
|                            | >Frame Reference Date Time                    | (0018,9151) | set by creator   |
|                            | >Frame Acquisition Date Time                  | (0018,9074) | set by creator   |
|                            | >Frame Acquisition Duration                   | (0018,9220) | set by creator   |
|                            | >Dimension Index Values                       | (0020,9157) | set by creator   |
|                            | >Temporal Position Index                      | (0020,9128) | set by creator   |
|                            | >Stack ID                                     | (0020,9056) | set by creator   |
|                            | >In-Stack Position Number                     | (0020,9057) | set by creator   |
|                            | >Frame Comments                               | (0020,9158) | set by creator   |
| Plane Position             | Plane Position Sequence                       | (0020,9113) | set by creator   |
|                            | >Image Position (Patient)                     | (0020,0032) | set by creator   |
| Plane Orientation          | Plane Orientation Sequence                    | (0020,9116) | set by creator   |
|                            | >Image Orientation (Patient)                  | (0020,0037) | set by creator   |
| Referenced Image           | Referenced Image Sequence                     | (0008,1140) | set by creator   |
|                            | >Referenced SOP Class UID                     | (0008,1150) | set by creator   |
|                            | >Referenced SOP Instance UID                  | (0008,1155) | set by creator   |
|                            | >Referenced Frame Number                      | (0008,1160) | set by creator   |
|                            | >Purpose of Referenced Code Sequence          | (0040,A170) | set by creator   |
|                            | >>Code Value                                  | (0008,0100) | e.g. 121311      |
|                            | >>Coding Scheme Designator                    | (0008,0102) | e.g. DCM         |
|                            | >>Code Meaning                                | (0008,0104) | e.g. Localizer   |
| Derivation Image           | Derivation Image Sequence                     | (0008,9124) | set by creator   |
|                            | >Derivation Code Sequence                     | (0008,9215) | set by creator   |
|                            | >>Code Value                                  | (0008,0100) | set by creator   |
|                            | >>Coding Scheme Designator                    | (0008,0102) | set by creator   |
|                            | >>Code Meaning                                | (0008,0104) | set by creator   |
|                            | >Source Image Sequence                        | (0008,2112) | set by creator   |
|                            | >>Referenced SOP Class UID                    | (0008,1150) | set by creator   |
|                            | >>Referenced SOP Instance UID                 | (0008,1155) | set by creator   |
|                            | >>Referenced Frame Number                     | (0008,1160) | set by creator   |
|                            | >>Purpose of Referenced Code Sequence         | (0040,A170) | set by creator   |
|                            | >>>Code Value                                 | (0008,0100) | set by creator   |
|                            | >>>Coding Scheme Designator                   | (0008,0102) | set by creator   |
|                            | >>>Code Meaning                               | (0008,0104) | set by creator   |
| Cardiac<br>Synchronization | Cardiac Synchronization Sequence              | (0018,9118) | set by creator   |
|                            | >Nominal Percentage of Cardiac Phase          | (0020,9241) | set by creator   |
|                            | >Nominal Cardiac Trigger Delay Time           | (0020,9153) | set by creator   |
|                            | >Actual Cardiac Trigger Delay Time            | (0020,9252) | set by creator   |
|                            | >Nominal Cardiac Trigger Time Prior to R-peak | (0020,9154) | set by creator   |
|                            | >Actual Cardiac Trigger Time Prior to R-peak  | (0020,9155) | set by creator   |
|                            | >Intervals Acquired                           | (0018,1083) | set by creator   |
|                            | >Intervals Rejected                           | (0018,1084) | set by creator   |
|                            | >Heart Rate                                   | (0018,1088) | set by creator   |
|                            | >R-R Interval Time Nominal                    | (0020,9251) | set by creator   |
|                            | >Low R-R Value                                | (0018,1081) | set by creator   |
|                            | >High R-R Value                               | (0018,1082) | set by creator   |
| Frame Anatomy              | Frame Anatomy Sequence                        | (0020,9071) | set by creator   |
|                            | >Frame Laterality                             | (0020,9072) | set by creator   |
|                            | >Anatomic Region Sequence                     | (0008,2218) | set by creator   |
|                            | >>Code Value                                  | (0008,0100) | set by creator   |

| Functional Group<br>Macro        | Attribute Name                                  | Tag         | Supported Values |
|----------------------------------|-------------------------------------------------|-------------|------------------|
| Pixel Value Transformation       | >>Coding Scheme Designator                      | (0008,0102) | set by creator   |
|                                  | >>Code Meaning                                  | (0008,0104) | set by creator   |
|                                  | Pixel Value Transformation Sequence             | (0028,9145) | set by creator   |
|                                  | >Rescale Intercept                              | (0028,1052) | set by creator   |
| Frame VOI LUT                    | >Rescale Slope                                  | (0028,1053) | set by creator   |
|                                  | >Rescale Type                                   | (0028,1054) | set by creator   |
|                                  | Frame VOI LUT Sequence                          | (0028,9132) | set by creator   |
|                                  | >Window Center                                  | (0028,1050) | set by creator   |
| Real World Value Mapping         | >Window Width                                   | (0028,1051) | set by creator   |
|                                  | Real World Value Mapping Sequence               | (0040,9096) | set by creator   |
|                                  | >Real World Value Intercept                     | (0040,9224) | set by creator   |
|                                  | >Real World Value Slope                         | (0040,9225) | set by creator   |
|                                  | >Measurement Units Code Sequence                | (0040,08EA) | set by creator   |
|                                  | >>Code Value                                    | (0008,0100) | set by creator   |
|                                  | >>Coding Scheme Designator                      | (0008,0102) | set by creator   |
| Contrast/Bolus Usage             | >>Code Meaning                                  | (0008,0104) | set by creator   |
|                                  | Contrast/Bolus Usage Sequence                   | (0018,9341) | set by creator   |
|                                  | >Contrast/Bolus Agent Number                    | (0018,9337) | 1                |
|                                  | >Contrast/Bolus Agent Administered              | (0018,9342) | set by creator   |
|                                  | >Contrast/Bolus Agent Detected                  | (0018,9343) | set by creator   |
| MR Image Frame Type              | >Contrast/Bolus Agent Phase                     | (0018,9344) | set by creator   |
|                                  | MR Image Frame Type Sequence                    | (0018,9226) | set by creator   |
|                                  | >Frame Type                                     | (0008,9007) | set by creator   |
|                                  | >Pixel Presentation                             | (0008,9205) | set by creator   |
|                                  | >Volumetric Properties                          | (0008,9206) | set by creator   |
|                                  | >Volume Based Calculation Technique             | (0008,9207) | set by creator   |
|                                  | >Complex Image Component                        | (0008,9208) | set by creator   |
|                                  | >Acquisition Contrast                           | (0008,9209) | set by creator   |
|                                  | >Functional Settling Phase Frames Present       | (0018,9622) | set by creator   |
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence       | (0018,9112) | set by creator   |
|                                  | >Repetition Time                                | (0018,0080) | set by creator   |
|                                  | >Flip Angle                                     | (0018,1314) | set by creator   |
|                                  | >Echo Train Length                              | (0018,0091) | set by creator   |
|                                  | >RF Echo Train Length                           | (0018,9240) | set by creator   |
|                                  | >Gradient Echo Train Length                     | (0018,9241) | set by creator   |
|                                  | >Specific Absorption Rate Sequence <sup>3</sup> | (0018,9239) | set by creator   |
|                                  | >>Specific Absorption Rate Definition           | (0018,9179) | set by creator   |
|                                  | >>Specific Absorption Rate Value <sup>3</sup>   | (0018,9181) | set by creator   |
|                                  | >Gradient Output Type <sup>3</sup>              | (0018,9180) | set by creator   |
|                                  | >Gradient Output <sup>3</sup>                   | (0018,9182) | set by creator   |
|                                  | >Operation Mode Sequence <sup>3</sup>           | (0018,9176) | set by creator   |
|                                  | >>Operating Mode Type                           | (0018,9177) | set by creator   |
| MR FOV/Geometry                  | >>Operating Mode <sup>3</sup>                   | (0018,9178) | set by creator   |
|                                  | MR FOV/Geometry Sequence                        | (0018,9125) | set by creator   |
|                                  | >In-plane Phase Encoding Direction              | (0018,1312) | set by creator   |
|                                  | >MR Acquisition Frequency Encoding Steps        | (0018,9058) | set by creator   |

<sup>3</sup> Specific Absorption Rate Sequence, Gradient Output Type, Gradient Output, and Operating Mode Sequence are not set for GRASP reconstructed images due to technical reasons. The values can be found in the related preview images of the GRASP acquisition.

| Functional Group<br>Macro | Attribute Name                                    | Tag         | Supported Values |
|---------------------------|---------------------------------------------------|-------------|------------------|
|                           | >MR Acquisition Phase Encoding Steps in-plane     | (0018,9231) | set by creator   |
|                           | >MR Acquisition Phase Encoding Steps out-of-plane | (0018,9232) | set by creator   |
|                           | >Percent Sampling                                 | (0018,0093) | set by creator   |
|                           | >Percent Phase Field of View                      | (0018,0094) | set by creator   |
| MR Echo                   | MR Echo Sequence                                  | (0018,9114) | set by creator   |
|                           | >Effective Echo Time                              | (0018,9082) | set by creator   |
| MR Modifier               | MR Modifier Sequence                              | (0018,9115) | set by creator   |
|                           | >Inversion Recovery                               | (0018,9009) | set by creator   |
|                           | >Inversion Times                                  | (0018,9079) | set by creator   |
|                           | >Flow Compensation                                | (0018,9010) | set by creator   |
|                           | >Flow Compensation Direction                      | (0018,9183) | set by creator   |
|                           | >Spoiling                                         | (0018,9016) | set by creator   |
|                           | >T2 Preparation                                   | (0018,9021) | set by creator   |
|                           | >Spectrally Selected Excitation                   | (0018,9026) | set by creator   |
|                           | >Spatial Pre-saturation                           | (0018,9027) | set by creator   |
|                           | >Partial Fourier                                  | (0018,9081) | set by creator   |
|                           | >Partial Fourier Direction                        | (0018,9036) | set by creator   |
|                           | >Parallel Acquisition                             | (0018,9077) | set by creator   |
|                           | >Parallel Acquisition Technique                   | (0018,9078) | set by creator   |
|                           | >Parallel Reduction Factor In-plane               | (0018,9069) | set by creator   |
|                           | >Parallel Reduction Factor out-of-plane           | (0018,9155) | set by creator   |
| MR Imaging Modifier       | MR Imaging Modifier Sequence                      | (0018,9006) | set by creator   |
|                           | >Magnetization Transfer                           | (0018,9020) | set by creator   |
|                           | >Blood Signal Nulling                             | (0018,9022) | set by creator   |
|                           | >Tagging                                          | (0018,9028) | set by creator   |
|                           | >Tag Spacing First Dimension                      | (0018,9030) | set by creator   |
|                           | >Tag Spacing Second Dimension                     | (0018,9218) | set by creator   |
|                           | >Tag Angle First Axis                             | (0018,9019) | set by creator   |
|                           | >Tag Angle Second Axis                            | (0018,9219) | set by creator   |
|                           | >Tag Thickness                                    | (0018,9035) | set by creator   |
|                           | >Tagging Delay                                    | (0018,9184) | set by creator   |
|                           | >Transmitter Frequency                            | (0018,9098) | set by creator   |
|                           | >Pixel Bandwidth                                  | (0018,0095) | set by creator   |
| MR Receive Coil           | MR Receive Coil Sequence                          | (0018,9042) | set by creator   |
|                           | >Receive Coil Name                                | (0018,1250) | set by creator   |
|                           | >Receive Coil Manufacturer Name                   | (0018,9041) | set by creator   |
|                           | >Receive Coil Type                                | (0018,9043) | set by creator   |
|                           | >Quadrature Receive Coil                          | (0018,9044) | set by creator   |
|                           | >Multi-Coil Definition Sequence                   | (0018,9045) | set by creator   |
|                           | >>Multi-Coil Element Name                         | (0018,9047) | set by creator   |
|                           | >>Multi-Coil Element Used                         | (0018,9048) | set by creator   |
| MR Transmit Coil          | MR Transmit Coil Sequence                         | (0018,9049) | set by creator   |
|                           | >Transmit Coil Name                               | (0018,1251) | set by creator   |
|                           | >Transmit Coil Manufacturer Name                  | (0018,9050) | set by creator   |
|                           | >Transmit Coil Type                               | (0018,9051) | set by creator   |
| MR Diffusion              | MR Diffusion Sequence                             | (0018,9117) | set by creator   |
|                           | >Diffusion b-value                                | (0018,9087) | set by creator   |
|                           | >Diffusion Directionality                         | (0018,9075) | set by creator   |
|                           | >Diffusion Gradient Direction Sequence            | (0018,9076) | set by creator   |
|                           | >>Diffusion Gradient Orientation                  | (0018,9089) | set by creator   |

| Functional Group<br>Macro | Attribute Name                     | Tag         | Supported Values |
|---------------------------|------------------------------------|-------------|------------------|
|                           | >Diffusion b-matrix Sequence       | (0018,9601) | set by creator   |
|                           | >>Diffusion b-value XX             | (0018,9602) | set by creator   |
|                           | >>Diffusion b-value XY             | (0018,9603) | set by creator   |
|                           | >>Diffusion b-value XZ             | (0018,9604) | set by creator   |
|                           | >>Diffusion b-value YY             | (0018,9605) | set by creator   |
|                           | >>Diffusion b-value YZ             | (0018,9606) | set by creator   |
|                           | >>Diffusion b-value ZZ             | (0018,9607) | set by creator   |
| MR Averages               | MR Averages Sequence               | (0018,9119) | set by creator   |
|                           | >Number of Averages                | (0018,0083) | set by creator   |
| MR Velocity Encoding      | MR Velocity Encoding Sequence      | (0018,9197) | set by creator   |
|                           | >Velocity Encoding Direction       | (0018,9090) | set by creator   |
|                           | >Velocity Encoding Minimum Value   | (0018,9091) | set by creator   |
|                           | >Velocity Encoding Maximum Value   | (0018,9217) | set by creator   |
| MR Arterial Spin Labeling | MR Arterial Spin Labeling Sequence | (0018,9251) | set by creator   |
|                           | >ASL Technique Description         | (0018,9252) | set by creator   |
|                           | >ASL Context                       | (0018,9257) | set by creator   |
|                           | >ASL Slab Sequence                 | (0018,9260) | set by creator   |
|                           | >>ASL Slab Number                  | (0018,9253) | set by creator   |
|                           | >>ASL Slab Thickness               | (0018,9254) | set by creator   |
|                           | >>ASL Slab Orientation             | (0018,9255) | set by creator   |
|                           | >>ASL Mid Slab Position            | (0018,9256) | set by creator   |
|                           | >>ASL Pulse Train Duration         | (0018,9258) | set by creator   |
|                           | >ASL Crusher Flag                  | (0018,9259) | set by creator   |
|                           | >ASL Crusher Flow Limit            | (0018,925A) | set by creator   |
|                           | >ASL Crusher Description           | (0018,925B) | set by creator   |
|                           | >ASLBolus Cut-off Flag             | (0018,925C) | set by creator   |
|                           | >ASLBolus Cut-off Timing Sequence  | (0018,925D) | set by creator   |
|                           | >>ASL Bolus Cut-off Delay Time     | (0018,925F) | set by creator   |
|                           | >>ASL Bolus Cut-off Technique      | (0018,925E) | set by creator   |

#### 9.1.1.1.14 Cardiac Synchronization Module

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| Cardiac Synchronization Technique | (0018,9037) | set by creator   |
| Cardiac Signal Source             | (0018,9085) | set by creator   |
| Cardiac RR Interval Specified     | (0018,9070) | set by creator   |
| Low R-R Value                     | (0018,1081) | set by creator   |
| High R-R Value                    | (0018,1082) | set by creator   |
| Intervals Acquired                | (0018,1083) | set by creator   |
| Intervals Rejected                | (0018,1084) | set by creator   |

**9.1.1.1.15 Respiratory Synchronization Module**

| Attribute Name                            | Tag         | Supported Values |
|-------------------------------------------|-------------|------------------|
| Respiratory Motion Compensation Technique | (0018,9170) | set by creator   |

**9.1.1.1.16 Bulk Motion Synchronization Module**

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Bulk Motion Compensation Technique | (0018,9172) | set by creator   |
| Bulk Motion Signal Source          | (0018,9173) | set by creator   |

**9.1.1.1.17 Supplemental Palette Color Lookup Table Module**

| Attribute Name                              | Tag         | Supported Values |
|---------------------------------------------|-------------|------------------|
| Red Palette Color Lookup Table Descriptor   | (0028,1101) | set by creator   |
| Green Palette Color Lookup Table Descriptor | (0028,1102) | set by creator   |
| Blue Palette Color Lookup Table Descriptor  | (0028,1103) | set by creator   |
| Red Palette Color Lookup Table Data         | (0028,1201) | set by creator   |
| Green Palette Color Lookup Table Data       | (0028,1202) | set by creator   |
| Blue Palette Color Lookup Table Data        | (0028,1203) | set by creator   |

**9.1.1.1.18 Acquisition Context Module**

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Acquisition Context Sequence | (0040,0555) | TID 99_3100      |

**9.1.1.1.18.1 MR Acquisition Context TID 99\_3100**

The MR Acquisition Context TID 99\_3100 is used in the Acquisition Context Sequence (0040,0555) of Enhanced MR Images.

Type: Extensible

| VT   | Concept Name                          | VM | Req Type | Value Set Constraint |
|------|---------------------------------------|----|----------|----------------------|
| CODE | EV (A-52, 99SMS_CTMR, "MR Technique") | 1  | M        | DCID(A-200)          |

**9.1.1.1.18.2 MR Technique Context Group A-200**

Type: Extensible

| Coding Scheme Designator<br>(0008,0102) | Coding Scheme<br>Version (0008,0103) | Code Value (0008,0100) | Code Meaning (0008,0104) |
|-----------------------------------------|--------------------------------------|------------------------|--------------------------|
| 99SMS_CTMR                              | 1.0                                  | A-200                  | GRASP                    |
| 99SMS_CTMR                              | 1.0                                  | A-202                  | FastView                 |
| 99SMS_CTMR                              | 1.0                                  | A-204                  | Angio                    |
| 99SMS_CTMR                              | 1.0                                  | A-206                  | BOLD                     |
| 99SMS_CTMR                              | 1.0                                  | A-208                  | Diffusion                |
| 99SMS_CTMR                              | 1.0                                  | A-210                  | DTI                      |
| 99SMS_CTMR                              | 1.0                                  | A-212                  | Perfusion                |
| 99SMS_CTMR                              | 1.0                                  | A-214                  | TWIST                    |

## 9.1.1.1.19 Enhanced MR Image Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Acquisition Number                 | (0020,0012) | set by creator   |
| Acquisition Date Time              | (0008,002A) | set by creator   |
| Acquisition Duration               | (0018,9073) | set by creator   |
| Content Qualification              | (0018,9004) | set by creator   |
| Resonant Nucleus                   | (0018,9100) | set by creator   |
| k-space Filtering                  | (0018,9064) | set by creator   |
| Magnetic Field Strength            | (0018,0087) | set by creator   |
| Applicable Safety Standard Agency  | (0018,9174) | set by creator   |
| Image Comments                     | (0020,4000) | set by creator   |
| Image Type                         | (0008,0008) | set by creator   |
| Pixel Presentation                 | (0008,9205) | set by creator   |
| Volumetric Properties              | (0008,9206) | set by creator   |
| Volume Based Calculation Technique | (0008,9207) | set by creator   |
| Complex Image Component            | (0008,9208) | set by creator   |
| Acquisition Contrast               | (0008,9209) | set by creator   |
| Samples per Pixel                  | (0028,0002) | 1                |
| Photometric Interpretation         | (0028,0004) | MONOCHROME2      |
| Bits Allocated                     | (0028,0100) | 16               |
| Bits Stored                        | (0028,0101) | 12 or 16         |
| High Bit                           | (0028,0102) | 11 or 15         |
| Pixel Representation               | (0028,0103) | 0                |
| Burned In Annotation               | (0028,0301) | NO               |
| Lossy Image Compression            | (0028,2110) | 00               |
| Presentation LUT Shape             | (2050,0020) | IDENTITY         |

## 9.1.1.1.20 MR Pulse Sequence Module

| Attribute Name                         | Tag         | Supported Values |
|----------------------------------------|-------------|------------------|
| Pulse Sequence Name                    | (0018,9005) | set by creator   |
| MR Acquisition Type                    | (0018,0023) | set by creator   |
| Echo Pulse Sequence                    | (0018,9008) | set by creator   |
| Multiple Spin Echo                     | (0018,9011) | set by creator   |
| Multi-planar Excitation                | (0018,9012) | set by creator   |
| Phase Contrast                         | (0018,9014) | set by creator   |
| Velocity Encoding Acquisition Sequence | (0018,9092) | set by creator   |
| > Velocity Encoding Direction          | (0018,9090) | set by creator   |
| Time of Flight Contrast                | (0018,9015) | set by creator   |
| Arterial Spin Labeling Contrast        | (0018,9250) | set by creator   |
| Steady State Pulse Sequence            | (0018,9017) | set by creator   |
| Echo Planar Pulse Sequence             | (0018,9018) | set by creator   |
| Saturation Recovery                    | (0018,9024) | set by creator   |
| Spectral Selected Suppression          | (0018,9025) | set by creator   |
| Oversampling Phase                     | (0018,9029) | set by creator   |
| Geometry of k-Space Traversal          | (0018,9032) | set by creator   |
| Rectilinear Phase Encode Reordering    | (0018,9034) | set by creator   |
| Segmented k-Space Traversal            | (0018,9033) | set by creator   |
| Coverage of k-Space                    | (0018,9094) | set by creator   |
| Number of k-Space Trajectories         | (0018,9093) | set by creator   |

## 9.1.1.1.21 ICC Profile Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| ICC Profile    | (0028,2000) | set by creator   |
| Color Space    | (0028,2002) | SRGB             |

## 9.1.1.1.22 SOP Common Module

| Attribute Name         | Tag         | Supported Values                  |
|------------------------|-------------|-----------------------------------|
| SOP Class UID          | (0008,0016) | Enhanced MR Storage SOP Class UID |
| SOP Instance UID       | (0008,0018) | set by creator                    |
| Specific Character Set | (0008,0005) | set by creator                    |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

## 9.1.1.1.23 Common Instance Reference Module

| Attribute Name                 | Tag         | Supported Values |
|--------------------------------|-------------|------------------|
| Referenced Series Sequence     | (0008,1115) | set by creator   |
| >Series Instance UID           | (0020,000E) | set by creator   |
| > Referenced Instance Sequence | (0008,114A) | set by creator   |
| >>Referenced SOP Class UID     | (0008,1150) | set by creator   |
| >> Referenced SOP Instance UID | (0008,1155) | set by creator   |

## 9.1.1.2 Enhanced MR Color Image IOD

The Enhanced MR Color Image IOD supports the same DICOM attributes as the Enhanced MR Image IOD except for the Photometric Interpretation which is RGB. Also, Supplemental Palette Color Lookup Table module, Pixel Value Transformation Macro, Frame VOI LUT Macro, and Real World Value Mapping Macro are not part of this IOD.

| IE                 | Module                        | Reference                            | Usage                             |
|--------------------|-------------------------------|--------------------------------------|-----------------------------------|
| Patient            | Patient                       | see chapter 9.1.1.1.1                | always                            |
| Study              | General Study                 | see chapter 9.1.1.1.2                | always                            |
|                    | Patient Study                 | see chapter 9.1.1.1.3                | always                            |
| Series             | General Series                | see chapter 9.1.1.1.4                | always                            |
|                    | MR Series                     | see chapter 9.1.1.1.5                | always                            |
| Frame of Reference | Frame of Reference            | see chapter 9.1.1.1.6                | always                            |
| Equipment          | General Equipment             | see chapter 9.1.1.1.7                | always                            |
|                    | Enhanced General Equipment    | see chapter 9.1.1.1.8                | always                            |
| Image              | Image Pixel                   | see chapter 9.1.1.2.19.1.1.9         | always                            |
|                    | Enhanced Contrast/Bolus       | see chapter 9.1.1.1.10               | always                            |
|                    | Multi-frame Functional Groups | see chapter 9.1.1.1.11               | always                            |
|                    | Multi-frame Dimension         | see chapter 9.1.1.1.12               | always                            |
|                    | Cardiac Synchronization       | see chapter 9.1.1.1.14               | in case of cardiac triggering     |
|                    | Respiratory Synchronization   | see chapter 9.1.1.1.15               | in case of respiratory triggering |
|                    | Bulk Motion Synchronization   | see chapter 9.1.1.1.16               | in case of bulk motion triggering |
|                    | Acquisition Context           | see chapter 9.1.1.1.18               | always                            |
|                    | Enhanced MR Image             | see chapter 9.1.1.1.19 and 9.1.1.2.2 | always                            |
|                    | MR Pulse Sequence             | see chapter 9.1.1.1.20               | always                            |
|                    | ICC Profile                   | see chapter 9.1.1.1.21               | always                            |
|                    | SOP Common                    | see chapter 9.1.1.1.22               | always                            |
|                    | Common Instance Reference     | see chapter 9.1.1.1.23               | always                            |

## 9.1.1.2.1 Image Pixel Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel          | (0028,0002) | 3                |
| Photometric Interpretation | (0028,0004) | RGB              |
| Rows                       | (0028,0010) | set by creator   |
| Columns                    | (0028,0011) | set by creator   |
| Bits Allocated             | (0028,0100) | 8                |
| Bits Stored                | (0028,0101) | 8                |
| High Bit                   | (0028,0102) | 7                |
| Pixel Representation       | (0028,0103) | 0                |
| Planar Configuration       | (0028,0006) | 0                |
| Pixel Data                 | (7FE0,0010) | set by creator   |

## 9.1.1.2.2 Enhanced MR Image Module

| Attribute Name     | Tag         | Supported Values |
|--------------------|-------------|------------------|
| Pixel Presentation | (0008,9205) | TRUE_COLOR       |

## 9.1.1.3 MR Spectroscopy IOD

This chapter describes the DICOM attributes of MR Spectroscopy Instances performed by the MR acquisition.

| IE                 | Module                        | Reference              | Usage                             |
|--------------------|-------------------------------|------------------------|-----------------------------------|
| Patient            | Patient                       | see chapter 9.1.1.1.1  | always                            |
| Study              | General Study                 | see chapter 9.1.1.1.2  | always                            |
|                    | Patient Study                 | see chapter 9.1.1.1.3  | always                            |
| Series             | General Series                | see chapter 9.1.1.1.4  | always                            |
|                    | MR Series                     | see chapter 9.1.1.1.5  | always                            |
| Frame of Reference | Frame of Reference            | see chapter 9.1.1.1.6  | always                            |
| Equipment          | General Equipment             | see chapter 9.1.1.1.7  | always                            |
|                    | Enhanced General Equipment    | see chapter 9.1.1.1.8  | always                            |
| MR Spectroscopy    | Enhanced Contrast/Bolus       | see chapter 9.1.1.1.10 | always                            |
|                    | Multi-frame Functional Groups | see chapter 9.1.1.1.11 | always                            |
|                    | Multi-frame Dimension         | see chapter 9.1.1.1.12 | always                            |
|                    | Cardiac Synchronization       | see chapter 9.1.1.1.14 | in case of cardiac triggering     |
|                    | Respiratory Synchronization   | see chapter 9.1.1.1.15 | in case of respiratory triggering |

| IE | Module                         | Reference                       | Usage                             |
|----|--------------------------------|---------------------------------|-----------------------------------|
|    | Bulk Motion Synchronization    | see chapter 9.1.1.1.16          | in case of bulk motion triggering |
|    | Acquisition Context            | see chapter 9.1.1.1.18          | always                            |
|    | MR Spectroscopy                | see chapter 9.1.1.3.29.1.1.1.19 | always                            |
|    | MR Spectroscopy Pulse Sequence | see chapter 9.1.1.3.3           | always                            |
|    | MR Spectroscopy Data           | see chapter 9.1.1.3.4           | always                            |
|    | SOP Common                     | see chapter 9.1.1.1.22          | always                            |
|    | Common Instance Reference      | see chapter 9.1.1.1.23          | always                            |

### 9.1.1.3.1 MR Spectroscopy Functional Groups

#### Functional

| Group             | Attribute Name                       | Tag         | Supported Values |
|-------------------|--------------------------------------|-------------|------------------|
| <b>Macro</b>      |                                      |             |                  |
| Pixel Measures    | Pixel Measures Sequence              | (0028,9110) | set by creator   |
|                   | >Pixel Spacing                       | (0018,0030) | set by creator   |
|                   | >Slice Thickness                     | (0028,0050) | set by creator   |
| Frame Content     | Frame Content Sequence               | (0020,9111) | set by creator   |
|                   | >Frame Acquisition Number            | (0020,9156) | set by creator   |
|                   | >Frame Reference Date Time           | (0018,9151) | set by creator   |
|                   | >Frame Acquisition Date Time         | (0018,9074) | set by creator   |
|                   | >Frame Acquisition Duration          | (0018,9220) | set by creator   |
|                   | >Dimension Index Values              | (0020,9157) | set by creator   |
|                   | >Temporal Position Index             | (0020,9128) | set by creator   |
|                   | >Stack ID                            | (0020,9056) | set by creator   |
|                   | >In-Stack Position Number            | (0020,9057) | set by creator   |
|                   | >Frame Comments                      | (0020,9158) | set by creator   |
| Plane Position    | Plane Position Sequence              | (0020,9113) | set by creator   |
|                   | >Image Position (Patient)            | (0020,0032) | set by creator   |
| Plane Orientation | Plane Orientation Sequence           | (0020,9116) | set by creator   |
|                   | >Image Orientation (Patient)         | (0020,0037) | set by creator   |
| Referenced Image  | Referenced Image Sequence            | (0008,1140) | set by creator   |
|                   | >Referenced SOP Class UID            | (0008,1150) | set by creator   |
|                   | >Referenced SOP Instance UID         | (0008,1155) | set by creator   |
|                   | >Referenced Frame Number             | (0008,1160) | set by creator   |
|                   | >Purpose of Referenced Code Sequence | (0040,A170) | set by creator   |

| Functional Group Macro     | Attribute Name                        | Tag         | Supported Values |
|----------------------------|---------------------------------------|-------------|------------------|
|                            | >>Code Value                          | (0008,0100) | set by creator   |
|                            | >>Coding Scheme Designator            | (0008,0102) | set by creator   |
|                            | >>Code Meaning                        | (0008,0104) | set by creator   |
| Derivation Image           | Derivation Image Sequence             | (0008,9124) | set by creator   |
|                            | >Derivation Code Sequence             | (0008,9215) | set by creator   |
|                            | >>Code Value                          | (0008,0100) | set by creator   |
|                            | >>Coding Scheme Designator            | (0008,0102) | set by creator   |
|                            | >>Code Meaning                        | (0008,0104) | set by creator   |
|                            | >Source Image Sequence                | (0008,2112) | set by creator   |
|                            | >>Referenced SOP Class UID            | (0008,1150) | set by creator   |
|                            | >>Referenced SOP Instance UID         | (0008,1155) | set by creator   |
|                            | >>Referenced Frame Number             | (0008,1160) | set by creator   |
|                            | >>Purpose of Referenced Code Sequence | (0040,A170) | set by creator   |
|                            | >>>Code Value                         | (0008,0100) | set by creator   |
|                            | >>>Coding Scheme Designator           | (0008,0102) | set by creator   |
|                            | >>>Code Meaning                       | (0008,0104) | set by creator   |
| Frame Anatomy              | Frame Anatomy Sequence                | (0020,9071) | set by creator   |
|                            | >Frame Laterality                     | (0020,9072) | set by creator   |
|                            | >Anatomic Region Sequence             | (0008,2218) | set by creator   |
|                            | >>Code Value                          | (0008,0100) | set by creator   |
|                            | >>Coding Scheme Designator            | (0008,0102) | set by creator   |
|                            | >>Code Meaning                        | (0008,0104) | set by creator   |
| Contrast/Bolus Usage       | Contrast/Bolus Usage Sequence         | (0018,9341) | set by creator   |
|                            | >Contrast/Bolus Agent Number          | (0018,9337) | 1                |
|                            | >Contrast/Bolus Agent Administered    | (0018,9342) | set by creator   |
|                            | >Contrast/Bolus Agent Detected        | (0018,9343) | set by creator   |
|                            | >Contrast/Bolus Agent Phase           | (0018,9344) | set by creator   |
| MR Spectroscopy Frame Type | MR Spectroscopy Frame Type Sequence   | (0018,9227) | set by creator   |
|                            | >Frame Type                           | (0008,9007) | set by creator   |
|                            | >Volumetric Properties                | (0008,9206) | set by creator   |
|                            | >Volume Based Calculation Technique   | (0008,9207) | set by creator   |
|                            | >Complex Image Component              | (0008,9208) | set by creator   |
|                            | >Acquisition Contrast                 | (0008,9209) | set by creator   |

| Functional Group Macro           | Attribute Name                                     | Tag         | Supported Values |
|----------------------------------|----------------------------------------------------|-------------|------------------|
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence          | (0018,9112) | set by creator   |
|                                  | >Repetition Time                                   | (0018,0080) | set by creator   |
|                                  | >Flip Angle                                        | (0018,1314) | set by creator   |
|                                  | >Echo Train Length                                 | (0018,0091) | set by creator   |
|                                  | >RF Echo Train Length                              | (0018,9240) | set by creator   |
|                                  | >Gradient Echo Train Length                        | (0018,9241) | set by creator   |
|                                  | >Specific Absorption Rate Sequence                 | (0018,9239) | set by creator   |
|                                  | >>Specific Absorption Rate Definition              | (0018,9179) | set by creator   |
|                                  | >>Specific Absorption Rate Value                   | (0018,9181) | set by creator   |
|                                  | >Gradient Output Type                              | (0018,9180) | set by creator   |
|                                  | >Gradient Output                                   | (0018,9182) | set by creator   |
|                                  | >Operation Mode Sequence                           | (0018,9176) | set by creator   |
|                                  | >>Operating Mode Type                              | (0018,9177) | set by creator   |
|                                  | >>Operating Mode                                   | (0018,9178) | set by creator   |
| MR Spectroscopy FOV/Geometry     | MR Spectroscopy FOV/Geometry Sequence              | (0018,9103) | set by creator   |
|                                  | >Spectroscopy Acquisition Data Columns             | (0018,9127) | set by creator   |
|                                  | >Spectroscopy Acquisition Phase Rows               | (0018,9095) | set by creator   |
|                                  | >Spectroscopy Acquisition Phase Columns            | (0018,9234) | set by creator   |
|                                  | >Spectroscopy Acquisition Out-of-Plane Phase Steps | (0018,9159) | set by creator   |
|                                  | >Percent Sampling                                  | (0018,0093) | set by creator   |
|                                  | >Percent Phase Field of View                       | (0018,0094) | set by creator   |
| MR Echo                          | MR Echo Sequence                                   | (0018,9114) | set by creator   |
|                                  | >Effective Echo Time                               | (0018,9082) | set by creator   |
| MR Modifier                      | MR Modifier Sequence                               | (0018,9115) | set by creator   |
|                                  | >Inversion Recovery                                | (0018,9009) | set by creator   |
|                                  | >Inversion Times                                   | (0018,9079) | set by creator   |
|                                  | >Flow Compensation                                 | (0018,9010) | set by creator   |
|                                  | >Flow Compensation Direction                       | (0018,9183) | set by creator   |
|                                  | >Spoiling                                          | (0018,9016) | set by creator   |
|                                  | >T2 Preparation                                    | (0018,9021) | set by creator   |
|                                  | >Spectrally Selected Excitation                    | (0018,9026) | set by creator   |
|                                  | >Spatial Pre-saturation                            | (0018,9027) | set by creator   |
|                                  | >Partial Fourier                                   | (0018,9081) | set by creator   |
|                                  | >Partial Fourier Direction                         | (0018,9036) | set by creator   |

**Functional**

| <b>Group<br/>Macro</b> | <b>Attribute Name</b>                   | <b>Tag</b>  | <b>Supported Values</b> |
|------------------------|-----------------------------------------|-------------|-------------------------|
|                        | >Parallel Acquisition                   | (0018,9077) | set by creator          |
|                        | >Parallel Acquisition Technique         | (0018,9078) | set by creator          |
|                        | >Parallel Reduction Factor In-plane     | (0018,9069) | set by creator          |
|                        | >Parallel Reduction Factor out-of-plane | (0018,9155) | set by creator          |
| MR Receive Coil        | MR Receive Coil Sequence                | (0018,9042) | set by creator          |
|                        | >Receive Coil Name                      | (0018,1250) | set by creator          |
|                        | >Receive Coil Manufacturer Name         | (0018,9041) | set by creator          |
|                        | >Receive Coil Type                      | (0018,9043) | set by creator          |
|                        | >Quadrature Receive Coil                | (0018,9044) | set by creator          |
|                        | >Multi-Coil Definition Sequence         | (0018,9045) | set by creator          |
|                        | >>Multi-Coil Element Name               | (0018,9047) | set by creator          |
|                        | >>Multi-Coil Element Used               | (0018,9048) | set by creator          |
| MR Transmit Coil       | MR Transmit Coil Sequence               | (0018,9049) | set by creator          |
|                        | >Transmit Coil Name                     | (0018,1251) | set by creator          |
|                        | >Transmit Coil Manufacturer Name        | (0018,9050) | set by creator          |
|                        | >Transmit Coil Type                     | (0018,9051) | set by creator          |
| MR Averages            | MR Averages Sequence                    | (0018,9119) | set by creator          |
|                        | >Number of Averages                     | (0018,0083) | set by creator          |

**9.1.1.3.2 MR Spectroscopy Module**

| <b>Attribute Name</b>              | <b>Tag</b>  | <b>Supported Values</b> |
|------------------------------------|-------------|-------------------------|
| Acquisition Number                 | (0020,0012) | set by creator          |
| Acquisition Date Time              | (0018,002A) | set by creator          |
| Acquisition Duration               | (0018,9073) | set by creator          |
| Referenced Raw Data Sequence       | (0008,9121) | set by creator          |
| Content Qualification              | (0018,9004) | set by creator          |
| Resonant Nucleus                   | (0018,9100) | set by creator          |
| k-space Filtering                  | (0018,9064) | set by creator          |
| Magnetic Field Strength            | (0018,0087) | set by creator          |
| Applicable Safety Standard Agency  | (0018,9174) | set by creator          |
| Image Comments                     | (0020,4000) | set by creator          |
| Image Type                         | (0008,0008) | set by creator          |
| Volumetric Properties              | (0008,9206) | set by creator          |
| Volume Based Calculation Technique | (0008,9207) | set by creator          |

| Attribute Name                       | Tag         | Supported Values |
|--------------------------------------|-------------|------------------|
| Complex Image Component              | (0008,9208) | set by creator   |
| Acquisition Contrast                 | (0008,9209) | set by creator   |
| Transmitter Frequency                | (0018,9098) | set by creator   |
| Spectral Width                       | (0018,9052) | set by creator   |
| Chemical Shift Reference             | (0018,9053) | set by creator   |
| Volume Localisation Technique        | (0018,9054) | set by creator   |
| Volume Localization Sequence         | (0018,9126) | set by creator   |
| >Slab Thickness                      | (0018,9104) | set by creator   |
| >Slab Orientation                    | (0018,9105) | set by creator   |
| >Mid SlabPosition                    | (0018,9106) | set by creator   |
| De-coupling                          | (0018,9059) | set by creator   |
| De-coupling Nucleus                  | (0018,9060) | set by creator   |
| De-coupling Frequency                | (0018,9061) | set by creator   |
| De-coupling Method                   | (0018,9062) | set by creator   |
| De-coupling Chemical Shift Reference | (0018,9063) | set by creator   |
| Time Domain Filtering                | (0018,9065) | set by creator   |
| Number Of Zero Fills                 | (0018,9066) | set by creator   |
| Baseline Correction                  | (0018,9067) | set by creator   |
| Frequency Correction                 | (0018,9101) | set by creator   |
| First Order Phase Correction         | (0018,9198) | set by creator   |
| Water Referenced Phase Correction    | (0018,9199) | set by creator   |
| Water Reference Acquisition          | (0018,9297) | set by creator   |
| Referenced Instance Sequence         | (0008,114A) | set by creator   |
| >Purpose of Reference Code Sequence  | (0040,A170) | set by creator   |

#### 9.1.1.3.3 MR Spectroscopy Pulse Sequence Module

| Attribute Name                      | Tag         | Supported Values |
|-------------------------------------|-------------|------------------|
| Pulse Sequence Name                 | (0018,9005) | set by creator   |
| MR Spectroscopy Acquisition Type    | (0018,9200) | set by creator   |
| Echo Pulse Sequence                 | (0018,9008) | set by creator   |
| Multi Spin Echo                     | (0018,9011) | set by creator   |
| Multi-planar Excitation             | (0018,9012) | set by creator   |
| Steady State Pulse Sequence         | (0018,9017) | set by creator   |
| Echo Planar Pulse Sequence          | (0018,9018) | set by creator   |
| Spectrally Selected Suppression     | (0018,9025) | set by creator   |
| Geometry of k-Space Traversal       | (0018,9032) | set by creator   |
| Rectilinear Phase Encode Reordering | (0018,9034) | set by creator   |
| Segmented k-Space Traversal         | (0018,9033) | set by creator   |
| Coverage of k-Space                 | (0018,9094) | set by creator   |

| Attribute Name                 | Tag         | Supported Values |
|--------------------------------|-------------|------------------|
| Number of k-Space Trajectories | (0018,9093) | set by creator   |
| Echo Peak Position             | (0018,9298) | set by creator   |

#### 9.1.1.3.4 MR Spectroscopy Data Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Rows                               | (0028,0010) | set by creator   |
| Columns                            | (0028,0011) | set by creator   |
| Data Point Rows                    | (0028,9001) | set by creator   |
| Data Point Columns                 | (0028,9002) | set by creator   |
| Data Representation                | (0028,9108) | set by creator   |
| Signal Domain Columns              | (0028,9003) | set by creator   |
| Signal Domains Rows                | (0028,9235) | set by creator   |
| First Order Phase Correction Angle | (5600,0010) | set by creator   |
| Spectroscopy Data                  | (5600,0020) | set by creator   |

#### 9.1.1.4 MR Image

This chapter describes the DICOM attributes of MR Image instances created in interoperability mode. The following table lists the modules created by syngo® MR.

| IE                 | Module                    | Reference              | Usage  |
|--------------------|---------------------------|------------------------|--------|
| Patient            | Patient                   | see chapter 9.1.1.1.1  | always |
| Study              | General Study             | see chapter 9.1.1.1.2  | always |
|                    | Patient Study             | see chapter 9.1.1.1.3  | always |
| Series             | General Series            | see chapter 9.1.1.1.4  | always |
| Frame of Reference | Frame of Reference        | see chapter 9.1.1.1.6  | always |
| Equipment          | General Equipment         | see chapter 9.1.1.1.7  | always |
| Image              | General Image             | see chapter 9.1.1.4.1  | always |
|                    | General Reference         | see chapter 9.1.1.4.2  | always |
|                    | Image Plane               | see chapter 9.1.1.4.3  | always |
|                    | Image Pixel               | see chapter 9.1.1.1.9  | always |
|                    | Contrast/Bolus            | see chapter 9.1.1.4.4  | always |
|                    | MR Image                  | see chapter 9.1.1.4.5  | always |
|                    | Overlay Plane             | see chapter 9.1.1.4.6  | always |
|                    | VOI LUT                   | see chapter 9.1.1.4.7  | always |
|                    | SOP Common                | see chapter 9.1.1.1.22 | always |
|                    | Common Instance Reference | see chapter 9.1.1.1.23 | always |

## 9.1.1.4.1 General Image Module

| Attribute Name         | Tag         | Supported Values |
|------------------------|-------------|------------------|
| Instance Number        | (0020,0013) | set by creator   |
| Patient Orientation    | (0020,0020) | set by creator   |
| Content Date           | (0008,0023) | set by creator   |
| Content Time           | (0008,0033) | set by creator   |
| Image Type             | (0008,0008) | set by creator   |
| Acquisition Number     | (0020,0012) | set by creator   |
| Acquisition Date       | (0008,0022) | set by creator   |
| Acquisition Time       | (0008,0032) | set by creator   |
| Acquisition DateTime   | (0008,002A) | set by creator   |
| Image Comments         | (0020,4000) | set by creator   |
| Presentation LUT Shape | (2050,0020) | set by creator   |

## 9.1.1.4.2 General Reference Module

| Attribute Name               | Tag         | Supported Values                                  |
|------------------------------|-------------|---------------------------------------------------|
| Referenced Image Sequence    | (0008,1140) | set by creator                                    |
| >Referenced SOP Class UID    | (0008,1150) | set by creator                                    |
| >Referenced SOP Instance UID | (0008,1155) | set by creator                                    |
| >Referenced Frame Number     | (0008,1160) | If referenced SOP Instance is a multi-frame image |

## 9.1.1.4.3 Image Plane Module

| Attribute Name              | Tag         | Supported Values |
|-----------------------------|-------------|------------------|
| Pixel Spacing               | (0028,0030) | set by creator   |
| Image Orientation (Patient) | (0020,0037) | set by creator   |
| Image Position (Patient)    | (0020,0032) | set by creator   |
| Slice Thickness             | (0018,0050) | set by creator   |
| Slice Location              | (0020,1041) | set by creator   |

## 9.1.1.4.4 Contrast/Bolus Module

| Attribute Name                | Tag         | Supported Values |
|-------------------------------|-------------|------------------|
| Contrast/Bolus Agent          | (0018,0010) | set by creator   |
| Contrast/Bolus Agent Sequence | (0018,0012) | set by creator   |
| >Code Value                   | (0008,0100) | set by creator   |
| >Coding Scheme Designator     | (0008,0102) | set by creator   |

| Attribute Name                               | Tag         | Supported Values |
|----------------------------------------------|-------------|------------------|
| >Coding Scheme Version                       | (0008,0103) | set by creator   |
| >Code Meaning                                | (0008,0104) | set by creator   |
| Contrast/Bolus Route                         | (0018,1040) | set by creator   |
| Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator   |
| >Code Value                                  | (0008,0100) | set by creator   |
| >Coding Scheme Designator                    | (0008,0102) | set by creator   |
| >Coding Scheme Version                       | (0008,0103) | set by creator   |
| >Code Meaning                                | (0008,0104) | set by creator   |
| Contrast/Bolus Volume                        | (0018,1041) | set by creator   |
| Contrast/Bolus Start Time                    | (0018,1042) | set by creator   |
| Contrast/Bolus Stop Time                     | (0018,1043) | set by creator   |
| Contrast/Bolus Total Dose                    | (0018,1044) | set by creator   |
| Contrast Flow Rate                           | (0018,1046) | set by creator   |
| Contrast Flow Duration                       | (0018,1047) | set by creator   |
| Contrast/Bolus Ingredient                    | (0018,1048) | set by creator   |
| Contrast/Bolus Ingredient Concentration      | (0018,1049) | set by creator   |

#### 9.1.1.4.5 MR Image Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Image Type                 | (0008,0008) | set by creator   |
| Samples per Pixel          | (0028,0002) | set by creator   |
| Photometric Interpretation | (0028,0004) | set by creator   |
| Bits Allocated             | (0028,0100) | set by creator   |
| Bits Stored                | (0028,0101) | set by creator   |
| High Bit                   | (0028,0102) | set by creator   |
| Scanning Sequence          | (0018,0020) | set by creator   |
| Sequence Variant           | (0018,0021) | set by creator   |
| Scan Options               | (0018,0022) | set by creator   |
| MR Acquisition Type        | (0018,0023) | set by creator   |
| Repetition Time            | (0018,0080) | set by creator   |
| Echo Time                  | (0018,0081) | set by creator   |
| Echo Train Length          | (0018,0091) | set by creator   |
| Inversion Time             | (0018,0082) | set by creator   |
| Trigger Time               | (0018,1060) | set by creator   |
| Sequence Name              | (0018,0024) | set by creator   |
| Angio Flag                 | (0018,0025) | set by creator   |
| Number of Averages         | (0018,0083) | set by creator   |
| Imaging Frequency          | (0018,0084) | set by creator   |
| Imaged Nucleus             | (0018,0085) | set by creator   |

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| Echo Number(s)                    | (0018,0086) | set by creator   |
| Magnetic Field Strength           | (0018,0087) | set by creator   |
| Spacing Between Slices            | (0018,0088) | set by creator   |
| Number of Phase Encoding Steps    | (0018,0089) | set by creator   |
| Percent Sampling                  | (0018,0093) | set by creator   |
| Percent Phase Field of View       | (0018,0094) | set by creator   |
| Pixel Bandwidth                   | (0018,0095) | set by creator   |
| Nominal Interval                  | (0018,1062) | set by creator   |
| Beat Rejection Flag               | (0018,1080) | set by creator   |
| Low R-R Value                     | (0018,1081) | set by creator   |
| High R-R Value                    | (0018,1082) | set by creator   |
| Intervals Acquired                | (0018,1083) | set by creator   |
| Intervals Rejected                | (0018,1084) | set by creator   |
| PVC Rejection                     | (0018,1085) | set by creator   |
| Skip Beats                        | (0018,1086) | set by creator   |
| Heart Rate                        | (0018,1088) | set by creator   |
| Cardiac Number of Images          | (0018,1090) | set by creator   |
| Trigger Window                    | (0018,1094) | set by creator   |
| Reconstruction Diameter           | (0018,1100) | set by creator   |
| Receive Coil Name                 | (0018,1250) | set by creator   |
| Transmit Coil Name                | (0018,1251) | set by creator   |
| Acquisition Matrix                | (0018,1310) | set by creator   |
| In-plane Phase Encoding Direction | (0018,1312) | set by creator   |
| Flip Angle                        | (0018,1314) | set by creator   |
| SAR                               | (0018,1316) | set by creator   |
| Variable Flip Angle Flag          | (0018,1315) | set by creator   |
| dB/dt                             | (0018,1318) | set by creator   |
| Temporal Position Identifier      | (0020,0100) | set by creator   |
| Number of Temporal Positions      | (0020,0105) | set by creator   |
| Temporal Resolution               | (0020,0110) | set by creator   |
| B1rms                             | (0018,1320) | set by creator   |

#### 9.1.1.4.6 Overlay Plane Module

| Attribute Name         | Tag         | Supported Values |
|------------------------|-------------|------------------|
| Overlay Rows           | (60xx,0010) | set by creator   |
| Overlay Columns        | (60xx,0011) | set by creator   |
| Overlay Type           | (60xx,0040) | set by creator   |
| Overlay Origin         | (60xx,0050) | set by creator   |
| Overlay Bits Allocated | (60xx,0100) | set by creator   |

| Attribute Name         | Tag         | Supported Values |
|------------------------|-------------|------------------|
| Overlay Bit Position   | (60xx,0102) | set by creator   |
| Overlay Data           | (60xx,3000) | set by creator   |
| Overlay Description    | (60xx,0022) | set by creator   |
| Overlay Subtype        | (60xx,0045) | set by creator   |
| Overlay Label          | (60xx,1500) | set by creator   |
| ROI Area               | (60xx,1301) | set by creator   |
| ROI Mean               | (60xx,1302) | set by creator   |
| ROI Standard Deviation | (60xx,1303) | set by creator   |

#### 9.1.1.4.7 VOI LUT Module

| Attribute Name                      | Tag         | Supported Values |
|-------------------------------------|-------------|------------------|
| VOI LUT Sequence                    | (0028,3010) | set by creator   |
| >LUT Descriptor                     | (0028,3002) | set by creator   |
| >LUT Explanation                    | (0028,3003) | set by creator   |
| >LUT Data                           | (0028,3006) | set by creator   |
| Window Center                       | (0028,1050) | set by creator   |
| Window Width                        | (0028,1051) | set by creator   |
| Window Center and Width Explanation | (0028,1055) | set by creator   |
| VOI LUT Function                    | (0028,1056) | set by creator   |

#### 9.1.1.5 Raw Data

Neuro Diffusion Tensor Imaging data are encoded in the Raw Data SOP Class.

| IE                 | Module              | Reference              | Usage  |
|--------------------|---------------------|------------------------|--------|
| Patient            | Patient             | see chapter 9.1.1.1.1  | always |
| Study              | General Study       | see chapter 9.1.1.1.2  | always |
|                    | Patient Study       | see chapter 9.1.1.1.3  | always |
| Series             | General Series      | see chapter 9.1.1.1.4  | always |
| Frame of Reference | Frame of Reference  | see chapter 9.1.1.1.6  | always |
| Equipment          | General Equipment   | see chapter 9.1.1.1.7  | always |
| Image              | Acquisition Context | see chapter 9.1.1.1.18 | always |
|                    | Raw Data            | see chapter 9.1.1.5.1  | always |
|                    | SOP Common          | see chapter 9.1.1.1.22 | always |

## 9.1.1.5.1 Raw Data Module

| Attribute Name       | Tag                        | Supported Values |
|----------------------|----------------------------|------------------|
| Instance Number      | (0020,0013)                | set by creator   |
| Content Date         | (0008,0023)                | set by creator   |
| Content Time         | (0008,0033)                | set by creator   |
| Image Type           | (0008,0008)                | set by creator   |
| Acquisition DateTime | (0008,002A)                | set by creator   |
| Creator_Version UID  | (0008,9123)                | set by creator   |
| Raw Data             | (7FE1,SIEMENS MR IMA, 10)) | set by creator   |

## 9.1.1.6 Multi-frame Grayscale Word Secondary Capture Image

For some use cases Multi-frame Grayscale Word SC Image instances are created for documentation and reporting purposes. The following table lists the modules created by *syngo*® MR XA50A.

| IE                 | Module                        | Reference                       | Presence of Module |
|--------------------|-------------------------------|---------------------------------|--------------------|
| Patient            | Patient                       | see chapter 9.1.1.1.1           | always             |
| Study              | General Study                 | see chapter 9.1.1.1.2           | always             |
|                    | Patient Study                 | see chapter 9.1.1.1.3           | always             |
| Series             | General Series                | see chapter 9.1.1.1.4           | always             |
| Frame of Reference | Frame of Reference            | see chapter 9.1.1.1.6           | always             |
| Equipment          | General Equipment             | see chapter 9.1.1.1.7           | always             |
|                    | SC Equipment                  | see chapter 9.1.1.6.1           | always             |
| Image              | General Image                 | see chapter 9.1.1.4.1           | always             |
|                    | Image Pixel                   | see chapter 9.1.1.1.9           | always             |
|                    | Cine                          | n. a.                           | not used           |
|                    | Multi-frame                   | see chapter 9.1.1.6.2           | always             |
|                    | Frame Pointers                | see chapter 9.1.1.6.39.1.1.1.12 | always             |
|                    | Multi-frame Functional Groups | see chapter 9.1.1.6.4           | always             |
|                    | SC Image                      | see chapter 9.1.1.6.6           | always             |
|                    | SC Multi-frame Image          | see chapter 9.1.1.6.7           | always             |
|                    | SC Multi-frame Vector         | see chapter 9.1.1.6.8           | always             |
|                    | VOI LUT                       | see chapter 9.1.1.4.6           | always             |
|                    | SOP Common                    | see chapter 9.1.1.1.22          | always             |
|                    | Common Instance Reference     | see chapter 9.1.1.1.23          | always             |

## 9.1.1.6.1 SC Equipment Module

| Attribute Name  | Tag         | Supported Values |
|-----------------|-------------|------------------|
| Conversion Type | (0008,0064) | set by creator   |
| Modality        | (0008,0060) | MR               |

## 9.1.1.6.2 Multi-frame Module

| Attribute Name          | Tag         | Supported Values |
|-------------------------|-------------|------------------|
| Number of Frames        | (0028,0008) | set by creator   |
| Frame Increment Pointer | (0028,0009) | set by creator   |

## 9.1.1.6.3 Frame Pointers Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Representative Frame Number     | (0028,6010) | set by creator   |
| Frame Numbers of Interest (FOI) | (0028,6020) | set by creator   |
| Frame of Interest Description   | (0028,6022) | set by creator   |
| Frame of Interest Type          | (0028,6023) | set by creator   |

## 9.1.1.6.4 Multi-frame Functional Groups Module

| Attribute Name                       | Tag         | Supported Values |
|--------------------------------------|-------------|------------------|
| Shared Functional Groups Sequence    | (5200,9229) | set by creator   |
| Per-frame Functional Groups Sequence | (5200,9230) | set by creator   |
| Instance Number                      | (0020,0013) | set by creator   |
| Content Date                         | (0008,0023) | set by creator   |
| Content Time                         | (0008,0033) | set by creator   |
| Number of Frames                     | (0028,0008) | set by creator   |
| Representative Frame Number          | (0028,6010) | set by creator   |

## 9.1.1.6.5 Multi-frame Grayscale Word SC Image Functional Group Macros

| Functional Group Macro                                                                        | Attribute Name             | Tag         | Supported Values |
|-----------------------------------------------------------------------------------------------|----------------------------|-------------|------------------|
| Pixel Measures<br>(if Plane Position (Patient) or Plane Orientation (Patient) Macros Present) | Pixel Measures Sequence    | (0028,9110) | set by creator   |
|                                                                                               | >Pixel Spacing             | (0018,0030) | set by creator   |
|                                                                                               | >Slice Thickness           | (0028,0050) | set by creator   |
| Plane Position (Patient)<br>(if Pixel Measures or Plane Orientation (Patient) Macros Present) | Plane Position Sequence    | (0020,9113) | set by creator   |
|                                                                                               | >Image Position (Patient)  | (0020,0032) | set by creator   |
| Plane Orientation (Patient)                                                                   | Plane Orientation Sequence | (0020,9116) | set by creator   |

| Functional Group Macro                                         | Attribute Name               | Tag         | Supported Values |
|----------------------------------------------------------------|------------------------------|-------------|------------------|
| (if Pixel Measures or Plane Position (Patient) Macros Present) | >Image Orientation (Patient) | (0020,0037) | set by creator   |

#### 9.1.1.6.6 SC Image Module

| Attribute Name                        | Tag         | Supported Values |
|---------------------------------------|-------------|------------------|
| Date of Secondary Capture             | (0018,1012) | set by creator   |
| Time of Secondary Capture             | (0018,1014) | set by creator   |
| Nominal Scanned Pixel Spa             | (0018,2010) | set by creator   |
| Pixel Spacing                         | (0028,0030) | set by creator   |
| Pixel Spacing Calibration Type        | (0028,0A02) | set by creator   |
| Pixel Spacing Calibration Description | (0018,1012) | set by creator   |

#### 9.1.1.6.7 SC Multi-frame Image Module

| Attribute Name                        | Tag         | Supported Values |
|---------------------------------------|-------------|------------------|
| Burned In Annotation                  | (0028,0301) | set by creator   |
| Recognizable Visual Features          | (0028,0302) | set by creator   |
| Presentation LUT Shape                | (2050,0020) | set by creator   |
| Illumination                          | (2010,015E) | set by creator   |
| Reflected Ambient Light               | (2010,0160) | set by creator   |
| Rescale Intercept                     | (0028,1052) | set by creator   |
| Rescale Slope                         | (0028,1053) | set by creator   |
| Rescale Type                          | (0028,1054) | set by creator   |
| Frame Increment Pointer               | (0028,0009) | set by creator   |
| Nominal Scanned Pixel Spacing         | (0018,2010) | set by creator   |
| Pixel Spacing                         | (0028,0030) | set by creator   |
| Pixel Spacing Calibration Type        | (0028,0A02) | set by creator   |
| Pixel Spacing Calibration Description | (0028,0A04) | set by creator   |
| Digitizing Device Transport Direction | (0018,2020) | set by creator   |
| Rotation of Scanned Film              | (0018,2030) | set by creator   |

#### 9.1.1.6.8 SC Multi-frame Vector Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Frame Time Vector          | (0018,1065) | set by creator   |
| Page Number Vector         | (0018,2001) | set by creator   |
| Frame Label Vector         | (0018,2002) | set by creator   |
| Frame Primary Angle Vector | (0018,2003) | set by creator   |

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Frame Secondary Angle Vector | (0018,2004) | set by creator   |
| Slice Location Vector        | (0018,2005) | set by creator   |
| Display Window Label Vector  | (0018,2006) | set by creator   |

#### 9.1.1.7 Secondary Capture Image

This chapter describes the DICOM attributes of Secondary Capture Image instances created in interoperability mode. The following table lists the modules created by *syngo*® MR XA50A.

| IE                 | Module                    | Reference              | Usage  |
|--------------------|---------------------------|------------------------|--------|
| Patient            | Patient                   | see chapter 9.1.1.1.1  | always |
| Study              | General Study             | see chapter 9.1.1.1.2  | always |
|                    | Patient Study             | see chapter 9.1.1.1.3  | always |
| Series             | General Series            | see chapter 9.1.1.1.4  | always |
| Frame of Reference | Frame of Reference        | see chapter 9.1.1.1.6  | always |
| Equipment          | General Equipment         | see chapter 9.1.1.1.7  | always |
|                    | SC Equipment              | see chapter 9.1.1.6.1  | always |
| Image              | General Image             | see chapter 9.1.1.4.1  | always |
|                    | Image Pixel               | see chapter 9.1.1.1.9  | always |
|                    | SC Image                  | see chapter 9.1.1.6.6  | always |
|                    | Overlay Plane             | see chapter 9.1.1.4.6  | always |
|                    | Modality LUT              | see chapter 9.1.1.7.1  | always |
|                    | VOI LUT                   | see chapter 9.1.1.4.7  | always |
|                    | SOP Common                | see chapter 9.1.1.1.22 | always |
|                    | Common Instance Reference | see chapter 9.1.1.1.23 | always |

##### 9.1.1.7.1 Modality LUT Module

| Attribute Name        | Tag         | Supported Values |
|-----------------------|-------------|------------------|
| Modality LUT Sequence | (0028,3000) | set by creator   |
| >LUT Descriptor       | (0028,3002) | set by creator   |
| >LUT Explanation      | (0028,3003) | set by creator   |
| >Modality LUT Type    | (0028,3004) | set by creator   |
| >LUT Data             | (0028,3006) | set by creator   |
| Rescale Intercept     | (0028,1052) | set by creator   |
| Rescale Slope         | (0028,1053) | set by creator   |
| Rescale Type          | (0028,1054) | set by creator   |

### 9.1.1.8 Grayscale Softcopy Presentation State

For some usecases Grayscale Softcopy Presentation State instances are created referencing Enhanced MR Image or Multi-frame Grayscale Word SC Image instances. The following table lists the modules created by the *syngo*® MR XA50A.

| IE                 | Module                            | Reference              | Usage  |
|--------------------|-----------------------------------|------------------------|--------|
| Patient            | Patient                           | see chapter 9.1.1.1.1  | always |
| Study              | General Study                     | see chapter 9.1.1.1.2  | always |
|                    | Patient Study                     | see chapter 9.1.1.1.3  | always |
| Series             | General Series                    | see chapter 9.1.1.1.4  | always |
|                    | Presentation Series               | see chapter 9.1.1.8.1  | always |
| Equipment          | General Equipment                 | see chapter 9.1.1.1.7  | always |
| Presentation State | Presentation State Identification | see chapter 9.1.1.8.2  | always |
|                    | Presentation State Relationship   | see chapter 9.1.1.8.3  | always |
|                    | Displayed Area                    | see chapter 9.1.1.8.4  | always |
|                    | Graphic Annotation                | see chapter 9.1.1.8.5  | always |
|                    | Spatial Transformation            | see chapter 9.1.1.8.6  | always |
|                    | Graphic Layer                     | see chapter 9.1.1.8.7  | always |
|                    | Graphic Group                     | see chapter 9.1.1.8.8  | always |
|                    | Modality LUT                      | see chapter 9.1.1.7.1  | always |
|                    | Softcopy VOI LUT                  | see chapter 9.1.1.8.9  | always |
|                    | SOP Common                        | see chapter 9.1.1.8.10 | always |
|                    | Common Instance Reference         | see chapter 9.1.1.1.23 | always |

#### 9.1.1.8.1 Presentation Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | PR               |

#### 9.1.1.8.2 Presentation State Identification Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Presentation Creation Date | (0070,0082) | set by creator   |
| Presentation Creation Time | (0070,0083) | set by creator   |
| Instance Number            | (0020,0013) | set by creator   |
| Content Label              | (0070,0080) | set by creator   |
| Content Description        | (0070,0081) | set by creator   |
| Content Creator's Name     | (0070,0084) | set by creator   |

## 9.1.1.8.3 Presentation State Relationship Module

| Attribute Name                | Tag         | Supported Values |
|-------------------------------|-------------|------------------|
| Referenced Series Sequence    | (0008,1115) | set by creator   |
| >Series Instance UID          | (0020,000E) | set by creator   |
| >Referenced Image Sequence    | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID    | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID | (0008,1155) | set by creator   |
| >>Referenced Frame Number     | (0008,1160) | set by creator   |

## 9.1.1.8.4 Displayed Area Module

| Attribute Name                           | Tag         | Supported Values |
|------------------------------------------|-------------|------------------|
| Displayed Area Selection Sequence        | (0070,005A) | set by creator   |
| >Referenced Image Sequence               | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID               | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID            | (0008,1155) | set by creator   |
| >>Referenced Frame Number                | (0008,1160) | set by creator   |
| >Pixel Origin Interpretation             | (0048,0301) | set by creator   |
| >Displayed Area Top Left Hand Corner     | (0070,0052) | set by creator   |
| >Displayed Area Bottom Right Hand Corner | (0070,0053) | set by creator   |
| >Presentation Size Mode                  | (0070,0100) | set by creator   |
| >Presentation Pixel Spacing              | (0070,0101) | set by creator   |
| >Presentation Pixel Aspect Ratio         | (0070,0102) | set by creator   |
| >Presentation Pixel Magnification Ratio  | (0070,0103) | set by creator   |

## 9.1.1.8.5 Graphic Annotation Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Graphic Annotation Sequence     | (0070,0001) | set by creator   |
| >Referenced Image Sequence      | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID      | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID   | (0008,1155) | set by creator   |
| >>Referenced Frame Number       | (0008,1160) | set by creator   |
| >Graphic Layer                  | (0070,0002) | set by creator   |
| >Text Object Sequence           | (0070,0008) | set by creator   |
| >>Bounding Box Annotation Units | (0070,0003) | set by creator   |
| >>Anchor Point Annotation Units | (0070,0004) | set by creator   |
| >>Unformatted Text Value        | (0070,0006) | set by creator   |
| >>Text Style Sequence           | (0070,0231) | set by creator   |
| >>>Font Name                    | (0070,0227) | set by creator   |

| Attribute Name                               | Tag         | Supported Values |
|----------------------------------------------|-------------|------------------|
| >>>Font Name Type                            | (0070,0228) | set by creator   |
| >>>CSS Font Name                             | (0070,0229) | set by creator   |
| >>>Text Color CIELab Value                   | (0070,0241) | set by creator   |
| >>>Horizontal Alignment                      | (0070,0242) | set by creator   |
| >>>Vertical Alignment                        | (0070,0243) | set by creator   |
| >>>Shadow Style                              | (0070,0244) | set by creator   |
| >>>Shadow Offset X                           | (0070,0245) | set by creator   |
| >>>Shadow Offset Y                           | (0070,0246) | set by creator   |
| >>>Shadow Color CIELab Val                   | (0070,0247) | set by creator   |
| >>>Shadow Opacity                            | (0070,0258) | set by creator   |
| >>>Underlined                                | (0070,0248) | set by creator   |
| >>>Bold                                      | (0070,0249) | set by creator   |
| >>>Italic                                    | (0070,0250) | set by creator   |
| >>Bounding Box Top Left Hand Corner          | (0070,0010) | set by creator   |
| >>Bounding Box Bottom Right Hand Corner      | (0070,0011) | set by creator   |
| >>Bounding Box Text Horizontal Justification | (0070,0012) | set by creator   |
| >>Anchor Point                               | (0070,0014) | set by creator   |
| >>Anchor Point Visibility                    | (0070,0015) | set by creator   |
| >>Compound Graphic Instance ID               | (0070,0226) | set by creator   |
| >>Graphic Group ID                           | (0070,0295) | set by creator   |
| >Graphic Object Sequence                     | (0070,0009) | set by creator   |
| >>Graphic Annotation Units                   | (0070,0005) | set by creator   |
| >>Graphic Dimensions                         | (0070,0020) | set by creator   |
| >>Number of Graphic Points                   | (0070,0021) | set by creator   |
| >>Graphic Data                               | (0070,0022) | set by creator   |
| >>Graphic Type                               | (0070,0023) | set by creator   |
| >>Line Style Sequence                        | (0070,0232) | set by creator   |
| >>>Pattern On Color CIELab Value;            | (0070,0251) | set by creator   |
| >>>Pattern Off Color CIELab Value            | (0070,0252) | set by creator   |
| >>>Pattern On Opacity                        | (0070,0284) | set by creator   |
| >>>Pattern Off Opacity                       | (0070,0285) | set by creator   |
| >>>Line Thickness                            | (0070,0253) | set by creator   |
| >>>Line Dashing Style                        | (0070,0254) | set by creator   |
| >>>Line Pattern                              | (0070,0255) | set by creator   |
| >>>Shadow Style                              | (0070,0244) | set by creator   |
| >>>Shadow Offset X                           | (0070,0245) | set by creator   |
| >>>Shadow Offset Y                           | (0070,0246) | set by creator   |
| >>>Shadow Color CIELab Value                 | (0070,0247) | set by creator   |
| >>>Shadow Opacity                            | (0070,0258) | set by creator   |
| >>Graphic Filled                             | (0070,0024) | set by creator   |

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| >>Fill Style Sequence             | (0070,0233) | set by creator   |
| >>>Pattern On Color CIELab Value; | (0070,0251) | set by creator   |
| >>>Pattern Off Color CIELab Value | (0070,0252) | set by creator   |
| >>>Pattern On Opacity             | (0070,0284) | set by creator   |
| >>>Pattern Off Opacity            | (0070,0285) | set by creator   |
| >>>Fill Mode                      | (0070,0257) | set by creator   |
| >>>Fill Pattern                   | (0070,0256) | set by creator   |
| >>Compound Graphic Instance ID    | (0070,0226) | set by creator   |
| >>Graphic Group ID                | (0070,0295) | set by creator   |
| >Compound Graphic Sequence        | (0070,0209) | set by creator   |
| >>Compound Graphic Instance       | (0070,0226) | set by creator   |
| >>Compound Graphic Units          | (0070,0282) | set by creator   |
| >>Graphic Dimensions              | (0070,0020) | set by creator   |
| >>Number of Graphic Points        | (0070,0021) | set by creator   |
| >>Graphic Data                    | (0070,0022) | set by creator   |
| >>Compound Graphic Type           | (0070,0294) | set by creator   |
| >Text Style Sequence              | (0070,0231) | set by creator   |
| >>>Font Name                      | (0070,0227) | set by creator   |
| >>>Font Name Type                 | (0070,0228) | set by creator   |
| >>>CSS Font Name                  | (0070,0229) | set by creator   |
| >>>Text Color CIELab Value        | (0070,0241) | set by creator   |
| >>>Horizontal Alignment           | (0070,0242) | set by creator   |
| >>>Vertical Alignment             | (0070,0243) | set by creator   |
| >>>Shadow Style                   | (0070,0244) | set by creator   |
| >>>Shadow Offset X                | (0070,0245) | set by creator   |
| >>>Shadow Offset Y                | (0070,0246) | set by creator   |
| >>>Shadow Color CIELab Value      | (0070,0247) | set by creator   |
| >>>Shadow Opacity                 | (0070,0258) | set by creator   |
| >>>Underlined                     | (0070,0248) | set by creator   |
| >>>Bold                           | (0070,0249) | set by creator   |
| >>>Italic                         | (0070,0250) | set by creator   |
| >>Line Style Sequence             | (0070,0232) | set by creator   |
| >>>Pattern On Color CIELab Value; | (0070,0251) | set by creator   |
| >>>Pattern Off Color CIELab Value | (0070,0252) | set by creator   |
| >>>Pattern On Opacity             | (0070,0284) | set by creator   |
| >>>Pattern Off Opacity            | (0070,0285) | set by creator   |
| >>>Line Thickness                 | (0070,0253) | set by creator   |
| >>>Line Dashing Style             | (0070,0254) | set by creator   |
| >>>Line Pattern                   | (0070,0255) | set by creator   |
| >>>Shadow Style                   | (0070,0244) | set by creator   |

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| >>>Shadow Offset X                | (0070,0245) | set by creator   |
| >>>Shadow Offset Y                | (0070,0246) | set by creator   |
| >>>Shadow Color CIELab Value      | (0070,0247) | set by creator   |
| >>>Shadow Opacity                 | (0070,0258) | set by creator   |
| >>Rotation Angle                  | (0070,0230) | set by creator   |
| >>Rotation Point                  | (0070,0273) | set by creator   |
| >>Gap Length                      | (0070,0261) | set by creator   |
| >>Diameter of Visibility          | (0070,0262) | set by creator   |
| >>Major Ticks Sequence            | (0070,0287) | set by creator   |
| >>>Tick Position                  | (0070,0288) | set by creator   |
| >>>Tick Label                     | (0070,0289) | set by creator   |
| >>Tick Alignment                  | (0070,0274) | set by creator   |
| >>Tick Label Alignment            | (0070,0279) | set by creator   |
| >>Show Tick Label                 | (0070,0278) | set by creator   |
| >>Graphic Filled                  | (0070,0024) | set by creator   |
| >>Fill Style Sequence             | (0070,0233) | set by creator   |
| >>>Pattern On Color CIELab Value; | (0070,0251) | set by creator   |
| >>>Pattern Off Color CIELab Value | (0070,0252) | set by creator   |
| >>>Pattern On Opacity             | (0070,0284) | set by creator   |
| >>>Pattern Off Opacity            | (0070,0285) | set by creator   |
| >>>Fill Mode                      | (0070,0257) | set by creator   |
| >>>Fill Pattern                   | (0070,0256) | set by creator   |
| >>Graphic Group ID                | (0070,0295) | set by creator   |

#### 9.1.1.8.6 Spatial Transformation Module

| Attribute Name        | Tag         | Supported Values |
|-----------------------|-------------|------------------|
| Image Rotation        | (0070,0042) | set by creator   |
| Image Horizontal Flip | (0070,0041) | set by creator   |

#### 9.1.1.8.7 Graphic Layer Module

| Attribute Name                                     | Tag         | Supported Values |
|----------------------------------------------------|-------------|------------------|
| Graphic Layer Sequence                             | (0070,0060) | set by creator   |
| >Graphic Layer                                     | (0070,0002) | set by creator   |
| >Graphic Layer Order                               | (0070,0062) | set by creator   |
| >Graphic Layer Recommended Display Grayscale Value | (0070,0066) | set by creator   |
| >Graphic Layer Recommended Display CIELab Value    | (0070,0401) | set by creator   |
| >Graphic Layer Description                         | (0070,0068) | set by creator   |

## 9.1.1.8.8 Graphic Group Module

| Attribute Name              | Tag         | Supported Values |
|-----------------------------|-------------|------------------|
| Graphic Group Sequence      | (0070,0234) | set by creator   |
| > Graphic Group ID          | (0070,0295) | set by creator   |
| > Graphic Group Label       | (0070,0207) | set by creator   |
| > Graphic Group Description | (0070,0208) | set by creator   |

## 9.1.1.8.9 Softcopy VOI LUT Module

| Attribute Name                       | Tag         | Supported Values |
|--------------------------------------|-------------|------------------|
| Softcopy VOI LUT Sequence            | (0028,3110) | set by creator   |
| >Referenced Image Sequence           | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID           | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID        | (0008,1155) | set by creator   |
| >>Referenced Frame Number            | (0008,1160) | set by creator   |
| >VOI LUT Sequence                    | (0028,3010) | set by creator   |
| >>LUT Descriptor                     | (0028,3002) | set by creator   |
| >>LUT Explanation                    | (0028,3003) | set by creator   |
| >>LUT Data                           | (0028,3006) | set by creator   |
| >Window Center                       | (0028,1050) | set by creator   |
| >Window Width                        | (0028,1051) | set by creator   |
| >Window Center and Width Explanation | (0028,1055) | set by creator   |
| >VOI LUT Function                    | (0028,1056) | set by creator   |

## 9.1.1.8.10 Softcopy Presentation LUT Module

| Attribute Name            | Tag         | Supported Values |
|---------------------------|-------------|------------------|
| Presentation LUT Sequence | (2050,0010) | set by creator   |
| >LUT Descriptor           | (0028,3002) | set by creator   |
| >LUT Explanation          | (0028,3003) | set by creator   |
| >LUT Data                 | (0028,3006) | set by creator   |
| Presentation LUT Shape    | (2050,0020) | set by creator   |

## 9.1.1.9 Enhanced and Comprehensive SR IOD

Enhanced and Comprehensive SR SOP Classes are used by applications e. g. cardiac analysis to store evaluated results in an Evidence Document.

| IE      | Module        | Reference             | Usage  |
|---------|---------------|-----------------------|--------|
| Patient | Patient       | see chapter 9.1.1.1.1 | always |
| Study   | General Study | see chapter 9.1.1.1.2 | always |

| IE                 | Module              | Reference              | Usage  |
|--------------------|---------------------|------------------------|--------|
|                    | Patient Study       | see chapter 9.1.1.1.3  | always |
| Series             | SR Document Series  | see chapter 9.1.1.9.1  | always |
| Frame of Reference | Frame of Reference  | see chapter 9.1.1.1.6  | always |
| Equipment          | General Equipment   | see chapter 9.1.1.1.7  | always |
| Image              | SR Document General | see chapter 9.1.1.9.2  | always |
|                    | SR Document Content | see chapter 9.1.1.9.3  | always |
|                    | SOP Common          | see chapter 9.1.1.1.22 | always |

#### 9.1.1.9.1 SR Document Series Module

| Attribute Name                               | Tag         | Supported Values |
|----------------------------------------------|-------------|------------------|
| Modality                                     | (0008,0060) | SR               |
| Series Instance UID                          | (0020,000E) | set by creator   |
| Series Number                                | (0020,0011) | set by creator   |
| Series Date                                  | (0008,0021) | set by creator   |
| Series Time                                  | (0008,0031) | set by creator   |
| Protocol Name                                | (0018,1030) | set by creator   |
| Series Description                           | (0008,103E) | set by creator   |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator   |
| >Referenced SOP Class UID                    | (0008,1150) | set by creator   |
| >Referenced SOP Instance UID                 | (0008,1155) | set by creator   |

#### 9.1.1.9.2 SR Document General Module

| Attribute Name              | Tag         | Supported Values                   |
|-----------------------------|-------------|------------------------------------|
| Series Number               | (0020,0011) | set by creator                     |
| Preliminary Flag            | (0040,A496) | set by creator (PRELIMINARY/FINAL) |
| Completion Flag             | (0040,A491) | set by creator (PARTIAL/FINAL)     |
| Completion Flag Description | (0040,A492) | set by creator                     |
| Verification Flag           | (0040,A493) | set by creator (UNVERIFIED)        |
| Content Date                | (0008,0023) | set by creator                     |
| Content Time                | (0008,0033) | set by creator                     |

## 9.1.1.9.3 SR Document Content Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Value Type                 | (0040,A040) | CONTAINER        |
| Concept Name Code Sequence | (0040,A043) | set by creator   |
| >Code Value                | (0008,0100) | set by creator   |
| >Coding Scheme Designator  | (0008,0102) | set by creator   |
| >Code Meaning              | (0018,0104) | set by creator   |
| Continuity of Content      | (0008,103E) | set by creator   |
| Content Template Sequence  | (0040,A504) | set by creator   |
| >Mapping Resource          | (0008,0105) | set by creator   |
| >Mapping Resource UID      | (0008,0118) | set by creator   |
| >Template Identifier       | (0040,DB00) | set by creator   |
| Content Sequence           | (0040,A730) | set by creator   |
| >Relationship Type         | (0040,A010) | set by creator   |

## 9.1.1.9.4 Evidence Document Templates

The finding related results of applications are stored into the Content Sequence of a Structured Evidence Document. The Content Items of a Content Sequence are specified in Structured Reporting Templates.

Examples of Structured Reporting Templates are:

- BI-RADS reporting
- Cardiac reporting
- PI-RADS 2 reporting
- PhoenixZip documentation.

## 9.1.2 Usage of Attributes from received IODs

N/A

## 9.1.3 Attribute mapping

The following table describes the mapping from attributes received in DICOM Modality Worklist to other attributes:

| Attribute Name                    | Tag         | Mapped attribute / Notes                |
|-----------------------------------|-------------|-----------------------------------------|
| <b>SOP Common</b>                 |             |                                         |
| Specific Character Set            | (0008,0005) |                                         |
| <b>Scheduled Procedure Step</b>   |             |                                         |
| Scheduled Procedure Step Sequence | (0040,0100) | Request Attributes Sequence (0040,0275) |

| Attribute Name                            | Tag         | Mapped attribute / Notes                                                                                                                                                                         |
|-------------------------------------------|-------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| >Modality                                 | (0008,0060) | not mapped                                                                                                                                                                                       |
| >Scheduled Station AE Title               | (0040,0001) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step Start Date      | (0040,0002) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step Start Time      | (0040,0003) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step End Date        | (0040,0004) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step End Time        | (0040,0005) | not mapped                                                                                                                                                                                       |
| >Scheduled Performing Physician's Name    | (0040,0006) | default for "Performing Physician's Name" (0008,1050)                                                                                                                                            |
| >Scheduled Procedure Step Description     | (0040,0007) | Scheduled Procedure Step Description, default for "Performed Procedure Step Description" (0040,0254)                                                                                             |
| >Scheduled Protocol Code Sequence **      | (0040,0008) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step ID              | (0040,0009) | "Scheduled Procedure Step ID" default for "Performed Procedure Step ID" (0040,0253)                                                                                                              |
| >Scheduled Station Name                   | (0040,0010) | not mapped                                                                                                                                                                                       |
| >Scheduled Procedure Step Location        | (0040,0011) | "Scheduled Procedure Step Location" is taken as default for "Performed Location"                                                                                                                 |
| >Scheduled Procedure Step Status          | (0040,0020) | not mapped                                                                                                                                                                                       |
| >Comments on the Scheduled Procedure Step | (0040,0400) | not mapped                                                                                                                                                                                       |
| <b>Requested Procedure</b>                |             |                                                                                                                                                                                                  |
| Study Date                                | (0008,0020) | not mapped                                                                                                                                                                                       |
| Study Time                                | (0008,0030) | not mapped                                                                                                                                                                                       |
| Referenced Study Sequence                 | (0008,1110) | Referenced Study Sequence (0008,1110)                                                                                                                                                            |
| >Referenced SOP Class UID                 | (0008,1150) | >Referenced SOP Class UID (0008,1150)                                                                                                                                                            |
| >Referenced SOP Instance UID              | (0008,1155) | not mapped                                                                                                                                                                                       |
| Study Instance UID                        | (0020,000D) | Study Instance UID                                                                                                                                                                               |
| Requested Procedure Description           | (0032,1060) | Requested Procedure Description in Request Attributes Sequence (0040,0275).<br>Used to preselect examination program in scheduler (favourite list).                                              |
| Requested Procedure Code Sequence         | (0032,1064) | "Requested Procedure Code Sequence" in Request Attributes Sequence (0040,0275).<br>Default for "Procedure Code Sequence"<br>Used to preselect examination program in scheduler (favourite list). |
| >Code Value                               | (0008,0100) | mapped                                                                                                                                                                                           |
| >Coding Scheme Designator                 | (0008,0102) | mapped                                                                                                                                                                                           |
| >Coding Scheme Version                    | (0008,0103) | mapped                                                                                                                                                                                           |
| >Code Meaning                             | (0008,0104) | mapped                                                                                                                                                                                           |

| Attribute Name                                | Tag         | Mapped attribute / Notes                                                                   |
|-----------------------------------------------|-------------|--------------------------------------------------------------------------------------------|
| Requested Procedure ID                        | (0040,1001) | "Requested Procedure ID" in Request Attributes Sequence (0040,0275) default for "Study ID" |
| Reason for the Requested Procedure            | (0040,1002) | Reason for the Requested Procedure                                                         |
| Requested Procedure Priority                  | (0040,1003) | not mapped                                                                                 |
| Patient Transport Arrangements                | (0040,1004) | not mapped                                                                                 |
| Confidentiality Code                          | (0040,1008) | not mapped                                                                                 |
| Reporting Priority                            | (0040,1009) | Reporting Priority                                                                         |
| Names of intended Recipients of Results       | (0040,1010) | not mapped                                                                                 |
| Requested Procedure Comments                  | (0040,1400) | Requested Procedure Comments                                                               |
| <b>Imaging Service Request</b>                |             |                                                                                            |
| Accession Number                              | (0008,0050) | mapped                                                                                     |
| Referring Physician's Name                    | (0008,0090) | mapped                                                                                     |
| Requesting Physician                          | (0032,1032) | mapped                                                                                     |
| Requesting Service                            | (0032,1033) | not mapped                                                                                 |
| Issuing Date of Imaging Service Request       | (0040,2004) | not mapped                                                                                 |
| Issuing Time of Imaging Service Request       | (0040,2005) | not mapped                                                                                 |
| Placer Order Number / Imaging Service Request | (0040,2016) | mapped                                                                                     |
| Filler Order Number / Imaging Service Request | (0040,2017) | mapped                                                                                     |
| Order entered by ...                          | (0040,2008) | not mapped                                                                                 |
| Order Enterer's location                      | (0040,2009) | not mapped                                                                                 |
| Order Callback Phone Number                   | (0040,2010) | not mapped                                                                                 |
| Imaging Service Request Comments              | (0040,2400) | not mapped                                                                                 |
| <b>Visit Identification</b>                   |             |                                                                                            |
| Admission ID                                  | (0038,0010) | mapped                                                                                     |
| Issuer of Admission ID                        | (0038,0011) | not mapped                                                                                 |
| Institution Name                              | (0008,0080) | default for Institution Name                                                               |
| Institution Address                           | (0008,0081) | default for Institution Address                                                            |
| <b>Visit Status</b>                           |             |                                                                                            |
| Visit Status ID                               | (0038,0008) | mapped                                                                                     |
| Institution Address                           | (0008,0081) | mapped                                                                                     |
| Current Patient Location                      | (0038,0300) | mapped                                                                                     |
| Patient's Institution Residence               | (0038,0400) | mapped                                                                                     |
| <b>Visit Admission</b>                        |             |                                                                                            |
| Admitting Diagnosis Description               | (0008,1080) | mapped                                                                                     |
| Admitting Date                                | (0038,0020) | mapped                                                                                     |
| <b>Patient Identification</b>                 |             |                                                                                            |
| Patient's Name                                | (0010,0010) | mapped                                                                                     |

| Attribute Name                         | Tag         | Mapped attribute / Notes |
|----------------------------------------|-------------|--------------------------|
| Patient ID                             | (0010,0020) | mapped                   |
| Issuer of Patient ID                   | (0010,0021) | mapped                   |
| Other Patient IDs                      | (0010,1000) | mapped                   |
| Other Patient Names                    | (0010,1001) | mapped                   |
| Patient's Birth Name                   | (0010,1005) | not mapped               |
| <b>Patient Demographic</b>             |             |                          |
| Patient's Birth Date                   | (0010,0030) | mapped                   |
| Patient's Birth Time                   | (0010,0032) | mapped                   |
| Patient's Sex                          | (0010,0040) | mapped                   |
| Patient's Insurance Plan Code Sequence | (0010,0050) | not mapped               |
| Patient's Age                          | (0010,1010) | not mapped               |
| Patient's Size                         | (0010,1020) | mapped                   |
| Patient's Weight                       | (0010,1030) | mapped                   |
| Patient's Address                      | (0010,1040) | mapped                   |
| Military Rank                          | (0010,1080) | mapped                   |
| Branch of Service                      | (0010,1081) | not mapped               |
| Ethnic Group                           | (0010,2160) | mapped                   |
| Patient Comments                       | (0010,4000) | mapped                   |
| <b>Patient Medical</b>                 |             |                          |
| Medical Alerts                         | (0010,2000) | mapped                   |
| Allergies                              | (0010,2110) | mapped                   |
| Pregnancy Status                       | (0010,21C0) | mapped                   |
| Smoking Status                         | (0010,21A0) | mapped                   |
| Last Menstrual Date                    | (0010,21D0) | not mapped               |
| Additional Patient History             | (0010,21B0) | not mapped               |
| Special Needs                          | (0038,0050) | mapped                   |

#### 9.1.4 Coerced / Modified fields

N/A

## 9.2 Data Dictionary of Private Attributes

Table 75 lists private attributes created by which may be included in instances generated by the syngo® MR XA50A. These private attributes may be deprecated or replaced with standard attributes in the future.

Table 75: Private Data Element Dictionary

| DICOM Tag                                     | Name                                  | VR | VM  |
|-----------------------------------------------|---------------------------------------|----|-----|
| (0029,SIEMENS CSA ENVELOPE,10)                | syngo Report Data                     | OB | 1   |
| (0029,SIEMENS CSA ENVELOPE,11)                | syngo Report Presentation             | OB | 1   |
| (0029,SIEMENS CSA HEADER,08)                  | Modality Image Header Type            | CS | 1   |
| (0029,SIEMENS CSA HEADER,09)                  | Modality Image Header Version         | LO | 1   |
| (0029,SIEMENS CSA HEADER,10)                  | Modality Image Header Info            | OB | 1   |
| (0029,SIEMENS CSA HEADER,18)                  | Modality Series Header Type           | CS | 1   |
| (0029,SIEMENS CSA HEADER,19)                  | Modality Series Header Version        | LO | 1   |
| (0029,SIEMENS CSA HEADER,20)                  | Modality Series Header Info           | OB | 1   |
| (0029,SIEMENS CSA NON-IMAGE,08)               | Modality Data Header Type             | CS | 1   |
| (0029,SIEMENS CSA NON-IMAGE,09)               | Modality Data Header Version          | LO | 1   |
| (0029,SIEMENS CSA NON-IMAGE,10)               | Modality Data Header Info             | OB | 1   |
| (7FE1,SIEMENS CSA NON-IMAGE,10)               | Modality Data                         | OB | 1   |
| (0029,SIEMENS CSA REPORT,08)                  | syngo Report Type                     | CS | 1   |
| (0029,SIEMENS CSA REPORT,09)                  | syngo Report Version                  | LO | 1   |
| (0029,SIEMENS CSA REPORT,15)                  | SR Variant                            | US | 1   |
| (0029,SIEMENS CSA REPORT,17)                  | SC SOP Instance UID                   | UI | 1   |
| (0029,SIEMENS MEDCOM HEADER,40)               | Application Header Sequence           | SQ | 1   |
| (0029,SIEMENS MEDCOM HEADER,41)               | Application Header Type               | CS | 1   |
| (0029,SIEMENS MEDCOM HEADER,42)               | Application Header ID                 | LO | 1   |
| (0029,SIEMENS MEDCOM HEADER,43)               | Application Header Version            | LO | 1   |
| (0029,SIEMENS MEDCOM HEADER,44)               | Application Header Info               | OB | 1   |
| (0029,SIEMENS MEDCOM HEADER,70)               | Siemens Link Sequence                 | SQ | 1-n |
| (0029,SIEMENS MEDCOM HEADER,71)               | Referenced Tag                        | AT | 1   |
| (0029,SIEMENS MEDCOM HEADER,72)               | Referenced Tag Type                   | CS | 1   |
| (0029,SIEMENS MEDCOM HEADER,73)               | Referenced Value Length               | UL | 1   |
| (0029,SIEMENS MEDCOM HEADER,74)               | Referenced Object Device Type         | CS | 1   |
| (0029,SIEMENS MEDCOM HEADER,75)               | Referenced Object Device Location     | OB | 1   |
| (0029,SIEMENS MEDCOM HEADER,76)               | Referenced Object ID                  | OB | 1   |
| (0029,SIEMENS MEDCOM HEADER,77)               | Referenced Object Offset              | UL | 1   |
| (0029,SIEMENS MEDCOM HEADER2,60)              | Series Work Flow Status               | LO | 1   |
| (0029,SIEMENS SYNGO 3D FUSION MATRIX,08)      | Object Series Instance UID            | UI | 1   |
| (0029,SIEMENS SYNGO 3D FUSION MATRIX,09)      | Model Series Instance UID             | UI | 1   |
| (0029,SIEMENS SYNGO 3D FUSION MATRIX,10)      | Matrix Referenced Series Instance UID | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,00) | Presentation Name                     | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,01) | Presentation Type                     | LO | 1   |

| DICOM Tag                                     | Name                                       | VR | VM  |
|-----------------------------------------------|--------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,02) | Advanced Presentation Sequence             | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,03) | Time Point Sequence                        | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,04) | Base Image Sequence                        | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,05) | Overlay Image Sequence                     | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,06) | Registration Instance Sequence             | SQ | 0-1 |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,07) | Real World Value Mapping Instance Sequence | SQ | 0-1 |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,08) | Measurement Sequence                       | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,09) | Measurement UID                            | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,10) | Segmentation Sequence                      | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,11) | Segmentation UID                           | SH | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,12) | Navigation Sequence                        | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,13) | Navigation Name                            | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,14) | Auto Navigation Direction                  | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,15) | Auto Navigation Frame Rate                 | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,16) | Auto Navigation Mode                       | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,17) | Auto Navigation Realtime Speed             | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,18) | Auto Navigation Strategy                   | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,19) | Auto Navigation Realtime Flag              | SH | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,20) | Index Navigation Current Index             | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,21) | Index Auto Navigation Skipping Degree      | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,22) | Volume Navigation Minimum Pixel Spacing    | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,23) | Volume Navigation Scroll Unit              | CS | 1   |

| DICOM Tag                                     | Name                                       | VR | VM  |
|-----------------------------------------------|--------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,24) | Volume Navigation Step Size                | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,25) | Volume Navigation Jump Size                | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,26) | Referenced Registration Number             | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,27) | Real World Value Mapping UID               | SH | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,28) | Channel Alpha Value                        | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,30) | Measurement Application Name               | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,31) | Measurement Data Sequence                  | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,32) | Measurement Type                           | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,33) | Measurement Frame of Reference UID         | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,34) | Measurement Uid                            | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,35) | Measurement Application Number             | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,36) | Measurement Application Number Prefix Text | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,37) | Measurement Graphic Is Visible Flag        | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,38) | Referenced Syngo Uid                       | UI | 4   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,39) | Clinical Finding Uid                       | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3A) | Measurement Evaluation String Value        | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3B) | Measurement Evaluation Integer Value       | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3C) | Measurement Evaluation Decimal Value       | FL | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3D) | Measurement Line Show Center Flag          | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3E) | Measurement Angle Show ArrowTip Flag       | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,3F) | Camera Home Settings Sequence              | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,40) | Camera Zoom                                | DS | 1   |

| DICOM Tag                                     | Name                                         | VR | VM  |
|-----------------------------------------------|----------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,41) | Camera Position                              | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,42) | Camera Orientation                           | DS | 4   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,43) | Camera Far Clip Plane                        | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,44) | Camera Near Clip Plane                       | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,46) | Camera ViewPort Size                         | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,47) | Camera Aspect Ratio                          | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,48) | Camera Projection Type                       | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,49) | Camera Field of View                         | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4A) | Camera Image Plane Distance                  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4B) | Camera Image Maximum Height                  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4C) | Camera Image Minimum Height                  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4D) | ParallelShift Interval MM                    | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4E) | Measurement ArrowTip Size                    | FL | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,4F) | Measurement ArrowTip Size ScalingFactor      | FL | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,50) | Renderer Vertex Is Characteristic Flag       | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,51) | Renderer Thickness Usage Flag                | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,52) | Renderer Threshold                           | DS | 2   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,53) | Renderer Material                            | DS | 4-8 |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,54) | Renderer DirectionalLight Color              | DS | 4   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,55) | Renderer DirectionalLight Direction          | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,56) | Renderer DirectionalLight TwoSide Usage Flag | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,57) | Renderer PWL TransferFunction Sequence       | SQ | 0-n |

| DICOM Tag                                     | Name                                     | VR | VM  |
|-----------------------------------------------|------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,58) | Renderer PWL Vertex Index                | IS | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,59) | Renderer PWL Vertex Color                | DS | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5A) | Renderer Is Camera Required Flag         | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5B) | Renderer Do Depth Test Flag              | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5C) | Renderer DirectionalLight Usage Flag     | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5D) | Renderer Thickness Sequence              | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5E) | Renderer Slice Spacing Sequence          | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,5F) | Renderer Sampling Distance               | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,60) | Renderer Initial Sampling Distance       | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,61) | Segmentation Display Data Sequence       | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,62) | Segmentation Display Data UID            | ST | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,63) | Segmentation Display Parameter Sequence  | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,64) | Segmentation Display Parameter Type      | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,65) | Segmentation Display Visibility          | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,66) | Segmentation Display Color               | DS | 4   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,67) | Segmentation Display Is Selected Flag    | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,68) | Segmentation Additional Information Blob | OB | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,69) | Hash Code Value                          | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,6A) | Segmentation Version Identifier          | LO | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,6B) | Segmentation Lock Mode                   | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,70) | Segmentation Volume Size                 | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7A) | Segmentation Volume StorageDataType      | DS | 1   |

| DICOM Tag                                     | Name                                       | VR | VM  |
|-----------------------------------------------|--------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7B) | Segmentation Volume Model Matrix           | FL | 16  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7C) | Segmentation Display Is Applied Flag       | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7D) | Display Representation Instance Identifier | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,71) | Registration Referenced Frames             | UI | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,72) | Registration Referenced Registrations      | UI | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,73) | Registration Creation Algorithm Name       | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,74) | ECG Graphic Type                           | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,75) | Hidden Pixel Spacing                       | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,76) | Renderer Edge Enhancement Param            | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,77) | Renderer Gradient Modulated Opacity Param  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,78) | Renderer Volume Smoothing Param            | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7E) | Fused Presentation LUT Shape               | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,7F) | Overlay Graphic VisibleFlag                | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,80) | Camera Rotation Axis                       | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,81) | Overlay Hidden Display Attributes          | SL | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,82) | Presentation State Group Identifier        | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,83) | Temporary Smallest Image Pixel Value       | US | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,84) | Camera Rotation Center                     | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,85) | Camera Rotation Center Usage Flag          | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,86) | Camera Parallel Epiped                     | DS | 12  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,87) | Camera Max Zoom In Factor                  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,88) | Camera Min Zoom In Factor                  | DS | 1   |

| DICOM Tag                                     | Name                                     | VR | VM  |
|-----------------------------------------------|------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,89) | Temporary Largest Image Pixel Value      | US | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,8A) | Camera Rotation Axis Usage Flag          | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,8B) | Measurement Surface Normal               | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,8C) | Measurement Ellipsoid Model Matrix       | FL | 16  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,8D) | Measurement Evaluation DataRole ID       | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,8E) | Measurement Algorithm Type               | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,91) | Measurement Evaluation DataRole Sequence | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,93) | Measurement Evaluation Sequence          | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,94) | Measurement Evaluation Value             | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,95) | Measurement Evaluation ID                | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,96) | Measurement Data Points                  | FL | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,97) | Measurement Data Angles                  | FL | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,98) | Measurement Data Slice                   | FD | 9   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,99) | Measurement Data Slice Thickness         | FL | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9A) | Measurement Referenced Frames Sequence   | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9B) | Measurement Evaluation Longest Distance  | DS | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9C) | Measurement Evaluation Centroid          | DS | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9D) | Measurement Data Bounding Box            | FL | 6   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9E) | Measurement Text                         | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,9F) | Measurement Diameter                     | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A0) | Image Rotation Fractional                | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A1) | Preset Name                              | LO | 1   |

| DICOM Tag                                     | Name                                               | VR | VM  |
|-----------------------------------------------|----------------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A2) | Fusion LUT Sequence                                | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A3) | Fusion LUT Is Active Flag                          | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A4) | Scale To Fit Type                                  | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A5) | Syngo UID                                          | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A6) | Presentation Version Identifier                    | UI | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A7) | Presentation Module Sequence                       | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A8) | Presentation Module Type                           | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,A9) | Presentation State Sequence                        | SQ | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AA) | LUT Inverted Flag                                  | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AB) | Softcopy Voi Lut Viewing Index                     | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AC) | Displayed Area Bottom Right Hand Corner Fractional | FD | 2   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AD) | Displayed Area Top Left Hand Corner Fractional     | FD | 2   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AE) | Measurement Alpha                                  | FL | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,AF) | Measurement Bitmap                                 | OB | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B0) | Current Frame Number                               | US | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B1) | ImageText View Name                                | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B2) | ImageText View Content Sequence                    | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B3) | ImageText Line Names                               | LO | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B4) | ImageText Line Values                              | LO | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B5) | Measurement Evaluation Text Position Sequence      | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B6) | Measurement Link Evaluation Text Flag              | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B7) | Measurement Evaluation Text Position Vector        | DS | 3   |

| DICOM Tag                                     | Name                                | VR | VM  |
|-----------------------------------------------|-------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B8) | ImageText Alpha Blending Line Value | OB | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,B9) | ImageText Visibility                | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,BA) | ImageText Reduced                   | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,BB) | ImageText Minimal                   | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,BC) | ImageText ViewID                    | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,BD) | Unlinked Image Text View Name       | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C1) | Task Data Sequence                  | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C2) | Task Data Type                      | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C3) | Task Data Version                   | LO | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C4) | Task Data Description               | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C5) | Task Data                           | OB | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C6) | Task Data Size                      | IS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C7) | True Size Type                      | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C8) | Image Graphics Visibility           | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,C9) | Clip Plane Sequence                 | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CA) | Clip Plane Center                   | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CB) | Clip Plane Normal                   | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CC) | Clip Plane Scale                    | DS | 2   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CD) | Clip Plane Use Thickness Flag       | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CE) | Clip Plane Thickness                | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,CF) | Image Sequence                      | SQ | 1-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,DO) | Clip Plane Enable Clip              | CS | 1   |

| DICOM Tag                                     | Name                                       | VR | VM  |
|-----------------------------------------------|--------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D1) | Clip Plane Handle Ratio                    | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D2) | Clip Plane Bounding Points                 | DS | 24  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D3) | Clip Plane Motion Matrix                   | DS | 16  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D4) | Clip Plane Shift Velocity                  | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D5) | Clip Plane Enabled Flag                    | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D6) | Clip Plane Rotate Velocity                 | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,D7) | Clip Plane Show Graphics Flag              | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E0) | Crop Box Size                              | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E1) | Crop Box Enabled Flag                      | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E2) | Crop Box Absolute Origin                   | DS | 3   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E3) | Crop Box Show Graphics Flag                | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E4) | Renderer Filter Settings Sequence          | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E5) | Renderer Energy Conversion Lut             | OF | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E6) | Apply Fusion                               | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,E7) | RepresentationId                           | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,EE) | Advanced Display Representation Sequence   | SQ | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,EF) | Measurement Visibility In View             | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F0) | Measurement Label Text                     | ST | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F1) | Curved Camera Coordinates                  | DS | 0-n |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F2) | Curved Camera Point Of Interest            | DS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F3) | Curved Camera Point Of Interest Usage Flag | CS | 1   |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F4) | Curved Camera Thickness                    | DS | 1   |

| DICOM Tag                                          | Name                                                            | VR | VM |
|----------------------------------------------------|-----------------------------------------------------------------|----|----|
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F5)      | Curved Camera Extrusion Length                                  | DS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F6)      | Curved Camera Rotation Axis Mode                                | CS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F7)      | Curved Camera Roll Rotation Axis                                | DS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F8)      | Curved Camera View Port Height                                  | DS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,F9)      | Curved Camera Cut Direction                                     | DS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FA)      | Curved Camera Pan Vector                                        | DS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FB)      | Clinical Finding ID                                             | LO | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FC)      | Measurement Is Circle Flag                                      | CS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FD)      | Measurement Application TypeID                                  | LO | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FE)      | Measurement Default Evaluation Text Visibility                  | CS | 1  |
| (0029,SIEMENS SYNGO ADVANCED PRESENTATION,FF)      | Measurement Font Height                                         | SL | 1  |
| (002B,SIEMENS SYNGO ADVANCED PRESENTATION,00)      | Measurement Min Threshold Value                                 | SL | 1  |
| (002B,SIEMENS SYNGO ADVANCED PRESENTATION,01)      | Measurement Max Threshold Value                                 | SL | 1  |
| 002B,SIEMENS SYNGO ADVANCED PRESENTATION,06        | Measurement Roi Binary Mask                                     | UT | 1  |
| 002B,SIEMENS SYNGO ADVANCED PRESENTATION,08        | Segmentation Double Threshold                                   | DS | 1  |
| 002B,SIEMENS SYNGO ADVANCED PRESENTATION,09        | Serialized Custom Attributes                                    | UT | 1  |
| (0029,SIEMENS SYNGO ALPHA CAD,10)                  | Automatic Landmarking and Parsing of Human Anatomy Version      | CS | 1  |
| (0029,SIEMENS SYNGO ALPHA CAD,12)                  | Automatic Landmarking and Parsing of Human Anatomy Body Regions | SQ | 1  |
| (7FDF,SIEMENS SYNGO DATA PADDING,FC)               | Pixel Data Leading Padding                                      | OB | 1  |
| (0087,SIEMENS SYNGO ENCAPSULATED DOCUMENT DATA,20) | Study Model                                                     | OB | 1  |
| (0087,SIEMENS SYNGO ENCAPSULATED DOCUMENT DATA,30) | Report XML Schema                                               | OB | 1  |
| (0087,SIEMENS SYNGO ENCAPSULATED DOCUMENT DATA,40) | Report Identifier                                               | OB | 1  |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,01)      | Business Unit Code                                              | CS | 1  |

| DICOM Tag                                      | Name                                        | VR | VM  |
|------------------------------------------------|---------------------------------------------|----|-----|
| (0027,SIEMENS SYNGO ENHANCED IDASET API,02)    | Application Type                            | LO | 1   |
| (0027,SIEMENS SYNGO ENHANCED IDASET API,03)    | Application Attributes Sequence             | SQ | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,10) | Evidence Document Template ID               | LO | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,11) | Evidence Document Template Version          | DS | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,20) | Evidence Document Clinical Finding Data     | OB | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,21) | Evidence Document Metadata                  | OB | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,30) | Evidence Document Implementation Version    | DS | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,40) | Evidence Document Predecessor               | OB | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,50) | Evidence Document Logical ID                | LO | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,60) | Evidence Document Application Data          | OB | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,70) | Evidence Document Owner Clinical Task Name  | LO | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,71) | Evidence Document Owner Task Name           | LO | 1   |
| (0077,SIEMENS SYNGO EVIDENCE DOCUMENT DATA,80) | Evidence Document Owner Supported Templates | OB | 1   |
| (0029,SIEMENS SYNGO FRAME SET,10)              | Image Frame Sequence                        | SQ | 1-n |
| (0029,SIEMENS SYNGO FRAME SET,12)              | Type of Progression                         | CS | 1   |
| (0029,SIEMENS SYNGO FRAME SET,14)              | Representation Level                        | IS | 1   |
| (0029,SIEMENS SYNGO FRAME SET,16)              | Representation Information Sequence         | SQ | 1-n |
| (0029,SIEMENS SYNGO FRAME SET,18)              | Number of Representations                   | IS | 1   |
| (0029,SIEMENS SYNGO FRAME SET,20)              | Representation Pixel Offset                 | IS | 1   |
| (0029,SIEMENS SYNGO FUNCTION ASSIGNMENT,01)    | Data Reference                              | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,20)          | Object Insertion Date                       | DA | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,30)          | Instance Object States                      | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,31)          | Series Object States                        | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,40)          | Last Access Time                            | DT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,41)          | Delete Protected Status                     | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,42)          | Received from Archive Status                | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,43)          | Archive Status                              | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,44)          | Location                                    | AE | 1   |

| DICOM Tag                             | Name                                                  | VR | VM  |
|---------------------------------------|-------------------------------------------------------|----|-----|
| (0009,SIEMENS SYNGO INDEX SERVICE,45) | Logical Deleted Status                                | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,46) | Insert Time                                           | DT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,47) | Visible Instances on Series Level                     | IS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,48) | Series Archived                                       | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,49) | Visible Instances on Study Level                      | IS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4A) | Series Completed                                      | IS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4B) | Delete Protection Time                                | DT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4C) | Delete Protection User                                | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4D) | Study Marked Time                                     | DT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4E) | Study Marked User                                     | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,4F) | Series Delete Protected                               | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,50) | Hidden Instance                                       | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,51) | Frame Number                                          | UL | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,60) | Data Handling Status                                  | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,61) | Result Status                                         | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,62) | Clinical Result Attribute                             | US | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,63) | Instance Origin                                       | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,64) | DM file relative location for store                   | LT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,70) | Header Offset                                         | UL | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,71) | Study Archived                                        | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,72) | Study Exported                                        | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,73) | Series Exported                                       | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,74) | Instance Exported                                     | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,75) | Study Corrected                                       | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,76) | Study Marked                                          | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,77) | Review                                                | SL | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,78) | Isotope Module Level                                  | SL | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,79) | Radionuclide Total Dose on Series level               | DS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,7A) | Radionuclide Total Dose on Instance level             | DS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,7B) | Series Radiopharmaceutical Information Sequence       | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,7C) | Instance Radiopharmaceutical Information Sequence     | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,80) | Workflow Attribute Sequence                           | SQ | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,81) | Patient's Name Attribute Sequence                     | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,82) | Other Patient Names Attribute Sequence                | SQ | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,83) | Name of Physician(s) Reading Study Attribute Sequence | SQ | 1-n |

| DICOM Tag                             | Name                                            | VR | VM  |
|---------------------------------------|-------------------------------------------------|----|-----|
| (0009,SIEMENS SYNGO INDEX SERVICE,84) | Requesting Physician Attribute Sequence         | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,85) | Referring Physician's Name Attribute Sequence   | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,86) | Performing Physician's Name Attribute Sequence  | SQ | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,87) | Family Name                                     | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,88) | Given Name                                      | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,89) | Middle Name                                     | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8A) | Prefix                                          | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8B) | Suffix                                          | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8C) | Degree                                          | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8D) | Ideographic Family Name                         | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8E) | Ideographic Given Name                          | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,8F) | Ideographic Middle Name                         | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,90) | Ideographic Prefix                              | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,91) | Ideographic Suffix                              | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,92) | Ideographic Degree                              | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,93) | Phonetic Family Name                            | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,94) | Phonetic Given Name                             | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,95) | Phonetic Middle Name                            | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,96) | Phonetic Prefix                                 | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,97) | Phonetic Suffix                                 | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,98) | Phonetic Degree                                 | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,99) | Study Printed                                   | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9A) | Series Printed                                  | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9B) | Instance Printed                                | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9C) | User Specified Worklist                         | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9D) | User Specified Worklist Study Registration Time | DT | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9E) | User Specified Worklist Study Registration User | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,9F) | User Specified Worklist Sequence                | SQ | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,A0) | Sender System Device Name                       | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A1) | Instance Repetition Time                        | DS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A2) | Instance Effective Echo Time                    | FD | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A3) | Instance Inversion Times                        | FD | 1-n |
| (0009,SIEMENS SYNGO INDEX SERVICE,A4) | Instance Nominal Cardiac Trigger Delay Time     | FD | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A5) | Instance Diffusion b-value                      | FD | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A6) | Instance Creator Entity                         | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A7) | Data Source                                     | AE | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,A8) | Study Availability                              | SL | 1   |

| DICOM Tag                                 | Name                                          | VR | VM  |
|-------------------------------------------|-----------------------------------------------|----|-----|
| (0009,SIEMENS SYNGO INDEX SERVICE,AC)     | Detector View Code Sequence                   | SQ | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,AD)     | Detector Information Code Value               | SH | 1-2 |
| (0009,SIEMENS SYNGO INDEX SERVICE,AE)     | Detector Information Code Meaning             | LO | 1-2 |
| (0009,SIEMENS SYNGO INDEX SERVICE,AF)     | Detector Information Coding Scheme Designator | SH | 1-2 |
| (0009,SIEMENS SYNGO INDEX SERVICE,B0)     | Relevant Body Part Examined                   | CS | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,B1)     | Relevant Protocol Name                        | LO | 1   |
| (0009,SIEMENS SYNGO INDEX SERVICE,B2)     | Relevant Modality                             | CS | 1   |
| (0009,SIEMENS SYNGO INSTANCE MANIFEST,00) | Temporary Original Header Sequence            | SQ | 1   |
| (0009,SIEMENS SYNGO INSTANCE MANIFEST,10) | syngo Index Source AE Title                   | AE | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,02)   | Hanging Protocol Excellence Rank              | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,06)   | Data Sharing Flag                             | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,08)   | Bagging Operations Sequence                   | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,10)   | Synchronization Type                          | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,12)   | Custom Filter Type                            | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,14)   | Custom Sorter Type                            | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,20)   | Selector DT Value                             | DT | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,22)   | Selector DA Value                             | DA | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,24)   | Selector TM Value                             | TM | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,26)   | Selector UI Value                             | UI | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,30)   | Custom Property Sequence                      | SQ | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,32)   | Custom Property Type                          | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,34)   | Custom Property Name                          | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,36)   | Custom Property Value                         | LT | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,38)   | Layout Property Sequence                      | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,40)   | Synchronization Sequence                      | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,60)   | Viewport Definitions Sequence                 | SQ | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,62)   | Protocol Type                                 | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,70)   | Display Protocol Name                         | SH | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,72)   | Display Protocol Description                  | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,74)   | Display Protocol Level                        | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,76)   | Display Protocol Creator                      | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,78)   | Display Protocol Creation Datetime            | DT | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,7A)   | Referenced Data Protocol                      | UI | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,7B)   | Original Display Protocol                     | UI | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,7C)   | Display Protocol Excellence Rank              | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,7E)   | Layout Sequence                               | SQ | 1-n |

| DICOM Tag                               | Name                              | VR | VM  |
|-----------------------------------------|-----------------------------------|----|-----|
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,80) | Layout Number                     | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,82) | Layout Description                | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,84) | Segment Sequence                  | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,86) | Segment Number                    | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,88) | Segment Description               | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,8A) | Segment Type                      | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,8C) | Tile Horizontal Dimension         | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,8E) | Tile Vertical Dimension           | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,90) | Fill Order                        | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,92) | Segment Small Scroll Type         | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,94) | Segment Small Scroll Amount       | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,96) | Segment Large Scroll Type         | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,98) | Segment Large Scroll Amount       | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,9A) | Segment Overlap Priority          | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,9C) | Data Role View Sequence           | SQ | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,9E) | Data Role View Number             | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,A2) | Referenced Data Role              | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,A4) | Sharing Enabled                   | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,A8) | Referenced Data Role Views        | US | 2-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,B0) | Data Protocol Name                | SH | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,B2) | Data Protocol Description         | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,B4) | Data Protocol Level               | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,B6) | Data Protocol Creator             | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,B8) | Data Protocol Creation Datetime   | DT | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,BA) | Data Protocol Excellence Rank     | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,BC) | Data Protocol Definition Sequence | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,BE) | Data Role Sequence                | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,C0) | Data Role Number                  | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,C2) | Data Role Name                    | SH | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,C4) | IsSystem                          | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,C6) | Selector Operations Sequence      | SQ | 0-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,C8) | Selector Usage Flag               | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,CA) | Select by Attribute Presence      | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,CC) | Select by Category                | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,CE) | Select by Operator                | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,D0) | Custom Selector Type              | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,D2) | Selector Operator                 | CS | 1   |

| DICOM Tag                               | Name                              | VR | VM  |
|-----------------------------------------|-----------------------------------|----|-----|
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,D4) | Reformatting Required             | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,D6) | Registration Data Sequence        | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,D8) | Reference Data Role Number        | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,DA) | Model Data Sequence               | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,DC) | Model Data Role Number            | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,DE) | Fusion Display Sequence           | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,E0) | Transparency                      | FD | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,E2) | Time Point                        | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,E4) | First Time Point Token            | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,E6) | Last Time Point Token             | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,E8) | Intermediate Time Point Token     | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,EA) | Data Processor Sequence           | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,EC) | Data Processor Type               | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,EE) | Template Data Role Sequence       | SQ | 1-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,F0) | View Sequence                     | SQ | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,F4) | View Type                         | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,F6) | Custom Bagging Type               | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,F8) | Referenced Display Segment Number | US | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,FA) | Data Role Type                    | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,46) | Internal Flag                     | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,FC) | Unassigned Flag                   | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,FE) | Initial Display Scroll Position   | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,04) | Template Data Role ID             | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,16) | Reference Template Data Role ID   | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,18) | Model Template Data Role ID       | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,28) | Referenced Template Data Role     | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,42) | Presentation Creator Type         | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,44) | Cine Navigation Type              | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,48) | Semantic Naming Strategy          | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,50) | Parameter String                  | LO | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,52) | Sorting Order                     | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,54) | syngo Template Type               | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,56) | Sorter Type                       | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,58) | Data Display Protocol Version     | SH | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,64) | TemplateSelectorSequence          | SQ | 0-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,66) | DefaultTemplate                   | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,68) | IsPreferred                       | CS | 1   |

| DICOM Tag                               | Name                            | VR | VM  |
|-----------------------------------------|---------------------------------|----|-----|
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,6A) | TimepointInitialValueSequence   | SQ | 0-n |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,6C) | TimepointVariable               | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,5A) | TimepointValue                  | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,5B) | SharingGroupSequence            | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,5C) | TemplateSelectorOperator        | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,5D) | SharingType                     | CS | 1   |
| (0073,SIEMENS SYNGO LAYOUT PROTOCOL,FF) | VRT Preset                      | LO | 1   |
| (0029,SIEMENS SYNGO MODULES,51)         | Incomplete Frame                | CS | 1   |
| (0029,SIEMENS SYNGO MODULES,52)         | Presentation UserData           | UT | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,00) | syngo Graphic Object Sequence   | SQ | 1-n |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,01) | Fill Style Version              | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,02) | Fill Background Color           | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,03) | Fill Foreground Color           | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,04) | Siemens Fill Mode               | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,05) | Siemens Fill Pattern            | OB | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,06) | Line Style Version              | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,07) | Line Background Color           | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,08) | Line Foreground Color           | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,09) | Line Type                       | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,10) | Siemens Line Thickness          | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,11) | Line Shadow X Offset            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,12) | Line Shadow Y Offset            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,13) | Siemens Shadow Style            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,14) | Shadow Color                    | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,15) | Stipple Pattern                 | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,16) | Line Anti Aliasing              | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,17) | Line-Z-Blend Flag               | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,18) | Text Style Version              | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,19) | Text Color                      | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,1A) | Connection Line Width Selected  | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,1B) | Connection Line Stipple Pattern | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,1C) | Connection Line Width           | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,20) | Text Horizontal Align           | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,21) | Text Vertical Align             | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,22) | Text Shadow X Offset            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,23) | Text Shadow Y Offset            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,24) | Text Shadow Style               | SL | 1   |

| DICOM Tag                               | Name                                  | VR | VM  |
|-----------------------------------------|---------------------------------------|----|-----|
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,25) | Text Shadow Color                     | FL | 4   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,26) | Text Log Font                         | CS | 1-n |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,27) | Text-Z-Blend Flag                     | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,28) | syngo Graphic Bit Mask                | OB | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,29) | Visibility Flag                       | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,30) | syngo Graphic Sensitivity             | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,31) | syngo Graphic Pick Mode Type          | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,32) | syngo Graphic Layer                   | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,33) | syngo Graphic Object Version          | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,34) | syngo Graphic Coordinate System       | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,35) | syngo Graphic Custom Attributes       | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,36) | syngo Graphic Custom Attributes Key   | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,37) | syngo Graphic Custom Attributes Value | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,38) | syngo Graphic View Name               | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,39) | syngo Graphic Data                    | DS | 3   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,40) | syngo Graphic Type                    | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,41) | Number of syngo Graphic Points        | US | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,42) | Axis Main Tick Length                 | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,43) | Axis Detail Tick Length               | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,44) | Axis Main Tick Spacing                | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,45) | Axis Detail Tick Spacing              | DS | 1-n |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,46) | Axis Main Tick Count                  | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,47) | Axis Detail Tick Count                | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,48) | Axis Tick Behavior                    | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,49) | Axis Tick Alignment                   | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,50) | Axis Step                             | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,51) | Axis Step Index                       | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,52) | Axis Text Format                      | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,53) | Axis Show Center Text Flag            | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,54) | Axis Show Tick Text Flag              | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,55) | Bitmap X Orientation                  | DS | 3   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,56) | Bitmap Y Orientation                  | DS | 3   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,57) | syngo Graphic Blob                    | OB | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,58) | syngo Graphic Interpolation           | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,59) | syngo Graphic Angle                   | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,60) | syngo Graphic Size                    | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,61) | Cut Line Side                         | CS | 1   |

| DICOM Tag                               | Name                           | VR | VM |
|-----------------------------------------|--------------------------------|----|----|
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,62) | syngo Graphic Tip Length       | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,63) | Cut Line Arrow Length          | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,64) | Line Gap Length                | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,65) | syngo Graphic Circle Radius    | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,66) | Line Distance Move             | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,67) | Line Marker Length             | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,68) | syngo Graphic Center           | DS | 3  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,69) | Range Center Area Top Left     | DS | 3  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,70) | Range Center Area Bottom Right | DS | 3  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,71) | Range Tilt                     | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,72) | Range Minimum Tilt             | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,73) | Range Maximum Tilt             | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,74) | syngo Graphic Width            | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,75) | Range Minimum Width            | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,76) | Range Maximum Width            | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,77) | syngo Graphic Height           | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,78) | Range Feed                     | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,79) | Range Direction                | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,80) | Range Show Scans               | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,81) | Range Minimum Scan Distance    | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,82) | Range Orthogonal Height        | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,83) | syngo Graphic Position         | DS | 3  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,84) | syngo Graphic Closed Flag      | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,85) | Range Line Tip Mode            | SL | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,86) | syngo Graphic List Count       | SL | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,87) | Axis Flip Text Flag            | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,88) | Curve Diagram Type             | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,89) | syngo Graphic Start Angle      | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,90) | syngo Graphic End Angle        | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,91) | Live Wire Smoothness           | IS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,92) | Live Wire Spline Flag          | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,93) | Ellipse Circle Flag            | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,94) | syngo Graphic Square Flag      | CS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,95) | Curve Section Start Index      | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,96) | Curve Section End Index        | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,97) | Marker Alpha                   | DS | 1  |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,98) | Table Row Count                | IS | 1  |

| DICOM Tag                                                   | Name                                          | VR | VM  |
|-------------------------------------------------------------|-----------------------------------------------|----|-----|
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,99)                     | Table Column Count                            | IS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,9A)                     | Table Row Height                              | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,9B)                     | Table Column Width                            | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,9C)                     | Rectangle Selection Segment Offset            | IS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,9D)                     | syngo Graphic Text                            | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A0)                     | Axis Tick Spacing Coordinate System           | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A1)                     | Axis Diagram Grid Type                        | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A2)                     | Polar Plot Circle Count                       | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A3)                     | Polar Plot Lines-per-Circle                   | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A4)                     | Polar Plot Compartment Count                  | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A5)                     | Polar Plot Radius Weight                      | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A6)                     | Circle Segment Outer Radius                   | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A7)                     | Circle Segment Clockwise Flag                 | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A8)                     | Axis Diagram Auto Resize Flag                 | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,A9)                     | Axis Diagram Step Start                       | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B0)                     | Group Root                                    | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B1)                     | Group Name                                    | ST | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B2)                     | syngo Graphic Annotation Sequence             | SQ | 1-n |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B3)                     | Text Minimum Height                           | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B4)                     | Text Font Scaling Factor                      | DS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B5)                     | Text Maximum Extensions                       | SL | 2   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B6)                     | Text Segment Size                             | CS | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,B7)                     | Graphic Object Reference Label                | SL | 1   |
| (0071,SIEMENS SYNGO OBJECT GRAPHICS,C4)                     | Show Grid                                     | CS | 1   |
| (0029,SIEMENS SYNGO PRINT SERVICE,10)                       | Sheet Number                                  | IS | 1   |
| (0071,SIEMENS SYNGO REGISTRATION,20)                        | Registered Image Sequence                     | SQ | 1   |
| (0071,SIEMENS SYNGO REGISTRATION,21)                        | Registration Is Validated Flag                | CS | 1   |
| (0071,SIEMENS SYNGO REGISTRATION,22)                        | Deformable Registration Inverse Grid Sequence | SQ | 1   |
| (0029,SIEMENS SYNGO RWVM,50)                                | RWVM Version                                  | LO | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION,10)                        | Segmentation Inverted                         | CS | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION,11)                        | Segmentation Uneditable by User               | CS | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,12) | Surface Mesh Private Data                     | SQ | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,13) | Surface Mesh Private Data Task name           | CS | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,14) | Surface Mesh Private Data Version Number      | LO | 1   |

| DICOM Tag                                                   | Name                                                        | VR | VM  |
|-------------------------------------------------------------|-------------------------------------------------------------|----|-----|
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,15) | Surface Mesh ImagePlane Sequence                            | SQ | 0-N |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,16) | ImagePlane Sequence Identifier                              | LO | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,17) | ImagePlane Sequence Identification Code Sequence            | SQ | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,18) | ImagePlane Sequence Identification Code Value               | SH | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,19) | ImagePlane Sequence Identification Coding Scheme Designator | SH | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,20) | ImagePlane Sequence Identification Code Meaning             | LO | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,21) | Surface Mesh ImagePlane                                     | SQ | 1-N |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,22) | ImagePlane                                                  | FL | 9   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,23) | ImagePlane Type                                             | CS | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,24) | ImagePlane Description                                      | LO | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,25) | Surface Mesh Point List Sequence                            | SQ | 0-N |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,26) | Point List Identifier                                       | LO | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,27) | Point List Identification Code Sequence                     | SQ | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,28) | Point List Identification Code Value                        | SH | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,29) | Point List Identification Coding Scheme Designator          | SH | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,30) | Point List Identification Code Meaning                      | LO | 1   |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,31) | Point List Description                                      | LO | 0-1 |
| (0029,SIEMENS SYNGO SEGMENTATION DICOM MESH PRIVATEDATA,32) | Point List                                                  | FL | 1-N |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,10)                   | SOP Class Packing Sequence                                  | SQ | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,20)                   | Packing Version                                             | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,21)                   | Packing Originator                                          | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,30)                   | Original SOP Class UID                                      | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,31)                   | Original Study Instance UID                                 | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,32)                   | Original Series Instance UID                                | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,33)                   | Original SOP Instance UID                                   | UI | 1   |

| DICOM Tag                                  | Name                               | VR | VM  |
|--------------------------------------------|------------------------------------|----|-----|
| (0031,SIEMENS SYNGO SOP CLASS PACKING,34)  | Original Transfer Syntax UID       | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,40)  | Attributes to Set to zero-length   | AT | 1-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,41)  | Attributes to Remove               | AT | 1-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,50)  | Original Rows                      | US | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,51)  | Original Columns                   | US | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,58)  | Original Image Type                | CS | 2-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,60)  | Original Modality                  | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,70)  | Sequence of original stream chunks | SQ | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,71)  | Start tag of a stream chunk        | AT | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,72)  | End tag of a stream chunk          | AT | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,73)  | Stream chunk is a PAYLOAD          | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,80)  | Stream chunk                       | OB | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,10)  | SOP Class Packing Sequence         | SQ | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,20)  | Packing Version                    | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,21)  | Packing Originator                 | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,30)  | Original SOP Class UID             | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,31)  | Original Study Instance UID        | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,32)  | Original Series Instance UID       | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,33)  | Original SOP Instance UID          | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,34)  | Original Transfer Syntax UID       | UI | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,40)  | Attributes to Set to zero-length   | AT | 1-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,41)  | Attributes to Remove               | AT | 1-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,50)  | Original Rows                      | US | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,51)  | Original Columns                   | US | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,58)  | Original Image Type                | CS | 2-n |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,60)  | Original Modality                  | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,70)  | Sequence of original stream chunks | SQ | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,71)  | Start tag of a stream chunk        | AT | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,72)  | End tag of a stream chunk          | AT | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,73)  | Stream chunk is a PAYLOAD          | CS | 1   |
| (0031,SIEMENS SYNGO SOP CLASS PACKING,80)  | Stream chunk                       | OB | 1   |
| (0029,SIEMENS SYNGO TIME POINT SERVICE,01) | Time Point ID                      | LO | 1   |
| (0029,SIEMENS SYNGO TIME POINT SERVICE,02) | Time Point Information             | LO | 1   |
| (0029,SIEMENS SYNGO TIME POINT SERVICE,50) | Studies in Time Point Sequence     | SQ | 1   |

| DICOM Tag                                                | Name                                     | VR | VM  |
|----------------------------------------------------------|------------------------------------------|----|-----|
| (7FD1,SIEMENS SYNGO ULTRA-SOUND TOYON DATA STREAMING,01) | Padding                                  | OB | 1   |
| (7FD1,SIEMENS SYNGO ULTRA-SOUND TOYON DATA STREAMING,09) | Version ID                               | OB | 1   |
| (7FD1,SIEMENS SYNGO ULTRA-SOUND TOYON DATA STREAMING,10) | Payload                                  | LO | 1   |
| (7FD1,SIEMENS SYNGO ULTRA-SOUND TOYON DATA STREAMING,11) | After Payload                            | LO | 1   |
| (0029,SIEMENS SYNGO VOLUME,12)                           | Slices                                   | US | 1   |
| (0029,SIEMENS SYNGO VOLUME,18)                           | Volume Level                             | IS | 1   |
| (0029,SIEMENS SYNGO VOLUME,30)                           | Voxel Spacing                            | DS | 3   |
| (0029,SIEMENS SYNGO VOLUME,32)                           | Volume Position (Patient)                | DS | 3   |
| (0029,SIEMENS SYNGO VOLUME,37)                           | Volume Orientation (Patient)             | DS | 9   |
| (0029,SIEMENS SYNGO VOLUME,40)                           | Resampling Flag                          | CS | 1   |
| (0029,SIEMENS SYNGO VOLUME,42)                           | Normalization Flag                       | CS | 1   |
| (0029,SIEMENS SYNGO VOLUME,44)                           | SubVolume Sequence                       | SQ | 1-n |
| (0029,SIEMENS SYNGO VOLUME,46)                           | Histogram Number Of Bins UL              | UL | 1   |
| (0029,SIEMENS SYNGO VOLUME,47)                           | Volume Histogram Data                    | OB | 1   |
| (0029,SIEMENS SYNGO VOLUME,48)                           | Volume Histogram BinBase                 | SL | 1   |
| (0029,SIEMENS SYNGO VOLUME,50)                           | Volume Version                           | LO | 1   |
| (0029,SIEMENS SYNGO VOLUME,60)                           | Total Frame Count Of Referenced Instance | IS | 0-1 |
| (0031,SIEMENS SYNGO WORKFLOW,10)                         | WF Internal Patient UID                  | UI | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,11)                         | WF Patients Death Indicator              | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,12)                         | WF Patients Death Date                   | DA | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,13)                         | WF Patients Death Time                   | TM | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,14)                         | WF VIP Indicator                         | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,15)                         | WF Emergency Flag                        | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,16)                         | Source                                   | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,20)                         | WF Internal Visit UID                    | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,25)                         | WF Internal ISR UID                      | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,32)                         | WF Control State                         | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,34)                         | WF Local Flag                            | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,36)                         | WF Referenced Studies                    | UI | 1-n |
| (0031,SIEMENS SYNGO WORKFLOW,40)                         | WF Workflow ID                           | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,41)                         | WF Workflow Description                  | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,42)                         | WF Workflow Control State                | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,43)                         | WF Workflow Ad Hoc Flag                  | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,44)                         | WF Hybrid Flag                           | US | 1   |

| DICOM Tag                        | Name                           | VR | VM  |
|----------------------------------|--------------------------------|----|-----|
| (0031,SIEMENS SYNGO WORKFLOW,45) | WF Workflow Read State         | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,50) | WF Workitem ID                 | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,51) | WF Workitem Name               | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,52) | WF Workitem Type               | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,53) | WF Workitem Roles              | LO | 1-n |
| (0031,SIEMENS SYNGO WORKFLOW,54) | WF Workitem Description        | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,55) | WF Workitem Control State      | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,56) | WF Claiming User               | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,57) | WF Claiming Host               | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,58) | WF Taskflow ID                 | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,59) | WF Taskflow Name               | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,5A) | WF Failed Flag                 | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,5B) | WF Scheduled Time              | DT | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,5C) | WF Workitem Ad Hoc Flag        | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,5D) | WF Patient Update Pending Flag | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,5E) | WF Patient Mixup Flag          | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,60) | WF Client ID                   | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,61) | WF Template ID                 | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,70) | WF E-mail                      | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,81) | WF Institution Name            | LO | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,82) | WF Institution Address         | ST | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,83) | WF Institution Code Sequence   | SQ | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,90) | WF Issuer of Master Patient ID | UI | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A0) | Patient History Sequence       | SQ | 1-n |
| (0031,SIEMENS SYNGO WORKFLOW,A1) | Procedure Marked Deleted Flag  | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A2) | ContextfolderId                | LT | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A3) | Visible DataRoles              | UT | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A4) | Preloaded                      | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A5) | TaskflowTechnology             | SL | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A6) | WfWorkflowType                 | SH | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A7) | NumberOfRequestedProcedures    | SL | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A8) | External Master Patient ID     | LT | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,A9) | Not To Be Archived             | US | 1   |
| (0031,SIEMENS SYNGO WORKFLOW,B0) | WorklistItemChecksum           | ST | 1   |
| (0019,SIEMENS MED NM,0F)         | ICON Data Type                 | SL | 1   |
| (0019,SIEMENS MED NM,16)         | Number of views                | SS | 1   |
| (0019,SIEMENS MED NM,93)         | Phase start time               | SL | 1-n |

| DICOM Tag                | Name                                                   | VR | VM  |
|--------------------------|--------------------------------------------------------|----|-----|
| (0019,SIEMENS MED NM,A1) | Number of Phases                                       | SS | 1   |
| (0019,SIEMENS MED NM,A5) | Number of repeats / phase                              | SS | 1-n |
| (0019,SIEMENS MED NM,A6) | Cycles Per Repeat                                      | SS | 1-n |
| (0019,SIEMENS MED NM,A7) | Repeat Start time                                      | SL | 1-n |
| (0019,SIEMENS MED NM,A8) | Repeat Stop time                                       | SL | 1-n |
| (0019,SIEMENS MED NM,A9) | Effective Repeat Time                                  | SL | 1-n |
| (0019,SIEMENS MED NM,AA) | Acquired Cycles Per Repeat                             | SS | 1-n |
| (0023,SIEMENS MED NM,01) | DICOM Reader Flag                                      | US | 1   |
| (0033,SIEMENS MED NM,00) | Flood correction matrix Detector 1 ESoft To SR2_0      | FL | n   |
| (0033,SIEMENS MED NM,01) | Flood correction matrix Detector 2 ESoft To SR2_0      | FL | n   |
| (0033,SIEMENS MED NM,10) | COR Data for Detector 1                                | FL | n   |
| (0033,SIEMENS MED NM,11) | COR Data for Detector 2                                | FL | n   |
| (0033,SIEMENS MED NM,14) | MHR Y Shift data for detector 1                        | FL | 1   |
| (0033,SIEMENS MED NM,15) | MHR Y Shift data for detector 2                        | FL | 1   |
| (0033,SIEMENS MED NM,18) | NCO Data for detector 1                                | FL | n   |
| (0033,SIEMENS MED NM,19) | NCO Data for detector 2                                | FL | n   |
| (0033,SIEMENS MED NM,20) | Bed correction angle                                   | FL | 1   |
| (0033,SIEMENS MED NM,21) | Gantry correction angle                                | FL | 1   |
| (0033,SIEMENS MED NM,22) | Bed U/D correction data                                | SS | n   |
| (0033,SIEMENS MED NM,23) | Gantry L/R correction data                             | SS | n   |
| (0033,SIEMENS MED NM,24) | BackProjection Correction angle head 1                 | FL | 1   |
| (0033,SIEMENS MED NM,25) | BackProjection Correction angle head 2                 | FL | 1   |
| (0033,SIEMENS MED NM,28) | Number point sources used for NCO and MHR Calibrations | SL | 1   |
| (0033,SIEMENS MED NM,29) | Crystal thickness                                      | FL | n   |
| (0033,SIEMENS MED NM,30) | Preset Name Used for Acquisition                       | LO | 1   |
| (0033,SIEMENS MED NM,31) | Camera Configuration Angle                             | FL | 1   |
| (0033,SIEMENS MED NM,32) | Crystal Type Starburst Or Not                          | LO | 1   |
| (0033,SIEMENS MED NM,33) | Gantry step for COIN acquisitions                      | SL | 1   |
| (0033,SIEMENS MED NM,34) | Bed step for wholebody or Coin acquisitions            | FL | 1   |
| (0033,SIEMENS MED NM,35) | Transaxial acceptance width For coincidence            | FL | 1   |
| (0033,SIEMENS MED NM,36) | Coincidence weight factor table                        | FL | n   |
| (0033,SIEMENS MED NM,37) | Starburst Flags                                        | SL | 1   |
| (0033,SIEMENS MED NM,38) | Pixel Scale factor                                     | FL | 1   |
| (0035,SIEMENS MED NM,00) | Specialized TOMO Type                                  | LO | 1   |
| (0035,SIEMENS MED NM,01) | Energy window type                                     | LO | 1   |
| (0035,SIEMENS MED NM,02) | Wing position                                          | SS | 2   |

| DICOM Tag                | Name                                         | VR | VM  |
|--------------------------|----------------------------------------------|----|-----|
| (0035,SIEMENS MED NM,03) | Loid of BlankScan Image                      | LO | 1   |
| (0035,SIEMENS MED NM,04) | Repeat ID                                    | SS | 1   |
| (0035,SIEMENS MED NM,05) | Phase ID                                     | SS | 1   |
| (0035,SIEMENS MED NM,06) | Siemens Profile 2 Image Sub type             | LO | 1   |
| (0037,SIEMENS MED NM,00) | Flood correction matrix for Detector 1       | OW | n   |
| (0037,SIEMENS MED NM,80) | Flood correction matrix for Detector 2       | OW | n   |
| (0039,SIEMENS MED NM,00) | Toshiba CBF activity results                 | LT | 1   |
| (0039,SIEMENS MED NM,01) | Related CT Series Instance UID               | LT | 1   |
| (0041,SIEMENS MED NM,01) | Whole Body Tomo Position Index               | SL | 1   |
| (0041,SIEMENS MED NM,02) | Whole Body Tomo Number of Positions          | SL | 1   |
| (0041,SIEMENS MED NM,03) | Horizontal Table Position of CT Scan         | FL | 1   |
| (0041,SIEMENS MED NM,04) | Effective Energy of CT Scan                  | FL | 1   |
| (0041,SIEMENS MED NM,05) | Long Linear Drive Information for Detector 1 | FD | 1-n |
| (0041,SIEMENS MED NM,06) | Long Linear Drive Information for Detector 2 | FD | 1-n |
| (0041,SIEMENS MED NM,07) | Trunnion Information for Detector 1          | FD | 1-n |
| (0041,SIEMENS MED NM,08) | Trunnion Information for Detector 2          | FD | 1-n |
| (0041,SIEMENS MED NM,09) | Broad beam factor                            | FL | 1   |
| (0041,SIEMENS MED NM,0A) | Original Wholebody Position                  | FD | 1   |
| (0041,SIEMENS MED NM,0B) | Wholebody Scan Range                         | FD | 1   |
| (0041,SIEMENS MED NM,10) | Effective Emission Energy                    | FL | n   |
| (0041,SIEMENS MED NM,11) | Tomo Frame View Duration                     | FL | 1-n |
| (0043,SIEMENS MED NM,01) | Detector View Angle                          | FL | 1-n |
| (0043,SIEMENS MED NM,02) | FOV Transformation Matrix                    | FD | 16  |
| (0043,SIEMENS MED NM,03) | View Dependent Y Shift MHR for Detector 1    | FL | 1-n |
| (0043,SIEMENS MED NM,04) | View Dependent Y Shift MHR for Detector 2    | FL | 1-n |
| (0045,SIEMENS MED NM,01) | Planar Processing String                     | LO | 1-n |
| (0045,SIEMENS MED NM,03) | Framing Parameters                           | UT | 1   |
| (0045,SIEMENS MED NM,04) | Scan Parameters                              | UT | 1   |
| (0045,SIEMENS MED NM,05) | Total Counts                                 | SL | 1   |
| (0045,SIEMENS MED NM,21) | Collimator Type ID                           | US | 1-n |
| (0045,SIEMENS MED NM,22) | Point Source Number                          | US | 1   |
| (0045,SIEMENS MED NM,23) | Calibration Type                             | LO | 1   |
| (0045,SIEMENS MED NM,24) | Camera Configuration                         | LO | 1   |
| (0055,SIEMENS MED NM,03) | View start time                              | SL | 1   |
| (0055,SIEMENS MED NM,04) | Prompt Window Width for Coincidence          | SS | 1   |
| (0055,SIEMENS MED NM,05) | Random Window Width for Coincidence          | SS | 1   |
| (0055,SIEMENS MED NM,7E) | Collimator Thickness                         | FL | 1-n |

| DICOM Tag                | Name                                                                | VR | VM  |
|--------------------------|---------------------------------------------------------------------|----|-----|
| (0055,SIEMENS MED NM,7F) | Collimator Angular Resolution                                       | FL | 1-n |
| (0055,SIEMENS MED NM,C0) | Useful Field of View                                                | SS | 1-n |
| (0057,SIEMENS MED NM,01) | Original ImageType                                                  | LO | 1   |
| (0057,SIEMENS MED NM,02) | Dose calibration factor                                             | FL | 1   |
| (0057,SIEMENS MED NM,03) | NM Pixel Units                                                      | LO | 1   |
| (0057,SIEMENS MED NM,04) | Decay Correction                                                    | LO | 1   |
| (0057,SIEMENS MED NM,05) | Radio Nuclide Half Life                                             | FL | n   |
| (0057,SIEMENS MED NM,06) | Rescale Intercept                                                   | FL | 1   |
| (0057,SIEMENS MED NM,07) | Rescale Slope                                                       | FL | 1   |
| (0057,SIEMENS MED NM,08) | Frame Reference Time                                                | FL | n   |
| (0057,SIEMENS MED NM,09) | Number original PET data Radiopharmaceutical<br>Info Sequence items | SL | 1   |
| (0057,SIEMENS MED NM,0A) | Decay Factor                                                        | FL | n   |
| (0057,SIEMENS MED NM,0B) | Counts Source                                                       | LO | 1   |
| (0057,SIEMENS MED NM,0C) | Radionuclide Positron Fraction                                      | FL | n   |
| (0057,SIEMENS MED NM,0E) | Trigger Time of CT Slice                                            | FD | 1-n |
| (0057,SIEMENS MED NM,0F) | QSPECT Compliant Flag                                               | SS | 1   |
| (0061,SIEMENS MED NM,01) | X Principal Ray Offset Detector 1                                   | FL | 1-n |
| (0061,SIEMENS MED NM,02) | X Principal Ray Offset Detector 2                                   | FL | 1-n |
| (0061,SIEMENS MED NM,05) | Y Principal Ray Offset Detector 1                                   | FL | 1-n |
| (0061,SIEMENS MED NM,06) | Y Principal Ray Offset Detector 2                                   | FL | 1-n |
| (0061,SIEMENS MED NM,09) | X Principal Ray Angle                                               | FL | 1-n |
| (0061,SIEMENS MED NM,0A) | Y Principal Ray Angle                                               | FL | 1-n |
| (0061,SIEMENS MED NM,0B) | X Short Focal Length                                                | FL | 1-n |
| (0061,SIEMENS MED NM,0C) | Y Short Focal Length                                                | FL | 1-n |
| (0061,SIEMENS MED NM,0D) | X Long Focal Length                                                 | FL | 1-n |
| (0061,SIEMENS MED NM,0E) | Y Long Focal Length                                                 | FL | 1-n |
| (0061,SIEMENS MED NM,0F) | X Focal Scaling                                                     | FL | 1-n |
| (0061,SIEMENS MED NM,10) | Y Focal Scaling                                                     | FL | 1-n |
| (0061,SIEMENS MED NM,11) | X Motion Correction Shift Detector 1                                | FL | 1-n |
| (0061,SIEMENS MED NM,12) | X Motion Correction Shift Detector 2                                | FL | 1-n |
| (0061,SIEMENS MED NM,15) | Y Motion Correction Shift Detector 1                                | FL | 1-n |
| (0061,SIEMENS MED NM,16) | Y Motion Correction Shift Detector 2                                | FL | 1-n |
| (0061,SIEMENS MED NM,19) | X Heart Center                                                      | FL | 1   |
| (0061,SIEMENS MED NM,1A) | Y Heart Center                                                      | FL | 1   |
| (0061,SIEMENS MED NM,1B) | Z Heart Center                                                      | FL | 1   |
| (0061,SIEMENS MED NM,1C) | Image Pixel Content Type                                            | LO | 1   |

| DICOM Tag                | Name                                    | VR | VM  |
|--------------------------|-----------------------------------------|----|-----|
| (0061,SIEMENS MED NM,1D) | Auto Save Corrected Series              | SS | 1   |
| (0061,SIEMENS MED NM,1E) | Distorted Series Instance UID           | LT | 1   |
| (0061,SIEMENS MED NM,21) | Recon Range                             | SS | 1-n |
| (0061,SIEMENS MED NM,22) | Recon Orientation                       | LO | 1   |
| (0061,SIEMENS MED NM,23) | Recon Selected Angular Range            | FL | 1-n |
| (0061,SIEMENS MED NM,24) | Recon Transverse Angle                  | FL | 1   |
| (0061,SIEMENS MED NM,25) | Recon Sagittal Angle                    | FL | 1   |
| (0061,SIEMENS MED NM,26) | Recon X Mask Size                       | FL | 1   |
| (0061,SIEMENS MED NM,27) | Recon Y Mask Size                       | FL | 1   |
| (0061,SIEMENS MED NM,28) | Recon X Image Center                    | FL | 1   |
| (0061,SIEMENS MED NM,29) | Recon Y Image Center                    | FL | 1   |
| (0061,SIEMENS MED NM,2A) | Recon Z Image Center                    | FL | 1   |
| (0061,SIEMENS MED NM,2B) | Recon X Zoom                            | FL | 1   |
| (0061,SIEMENS MED NM,2C) | Recon Y Zoom                            | FL | 1   |
| (0061,SIEMENS MED NM,2D) | Recon Threshold                         | FL | 1   |
| (0061,SIEMENS MED NM,2E) | Recon Output Pixel Size                 | FL | 1   |
| (0061,SIEMENS MED NM,2F) | Scatter Estimation Method               | LO | 1-n |
| (0061,SIEMENS MED NM,30) | Scatter Estimation Method Mode          | LO | 1-n |
| (0061,SIEMENS MED NM,31) | Scatter Estimation Lower Window Weights | FL | 1-n |
| (0061,SIEMENS MED NM,32) | Scatter Estimation Upper Window Weights | FL | 1-n |
| (0061,SIEMENS MED NM,33) | Scatter Estimation Window Mode          | LO | 1-n |
| (0061,SIEMENS MED NM,34) | Scatter Estimation Filter               | LO | 1-n |
| (0061,SIEMENS MED NM,35) | Recon RawTomo Input Uid                 | LO | 1-n |
| (0061,SIEMENS MED NM,36) | Recon CT Input Uid                      | LO | 1   |
| (0061,SIEMENS MED NM,37) | Recon Z Mask Size                       | FL | 1   |
| (0061,SIEMENS MED NM,38) | Recon X Mask Center                     | FL | 1   |
| (0061,SIEMENS MED NM,39) | Recon Y Mask Center                     | FL | 1   |
| (0061,SIEMENS MED NM,3A) | Recon Z Mask Center                     | FL | 1   |
| (0061,SIEMENS MED NM,3B) | Recon First slice index                 | SL | 1   |
| (0061,SIEMENS MED NM,3C) | Non Image UID                           | LT | 1   |
| (0061,SIEMENS MED NM,3D) | Non Image Series UID                    | LT | 1   |
| (0061,SIEMENS MED NM,3E) | Non Image Associated Parent Series UID  | LT | 1-n |
| (0061,SIEMENS MED NM,3F) | Raw Bin Times                           | FL | 1-n |
| (0061,SIEMENS MED NM,51) | Raw Tomo Series UID                     | LT | 1   |
| (0061,SIEMENS MED NM,52) | LowRes CT Series UID                    | LT | 1   |
| (0061,SIEMENS MED NM,53) | HighRes CT Series UID                   | LT | 1   |
| (0061,SIEMENS MED NM,54) | Vector Map Offset                       | FL | 1-4 |

| DICOM Tag                | Name                                               | VR | VM  |
|--------------------------|----------------------------------------------------|----|-----|
| (0061,SIEMENS MED NM,55) | Collimator Hole Length                             | FL | 1-n |
| (0061,SIEMENS MED NM,56) | Collimator Entry Hole Diameter                     | FL | 1-n |
| (0061,SIEMENS MED NM,57) | Collimator Exit Hole Diameter                      | FL | 1-n |
| (0061,SIEMENS MED NM,58) | Collimator Front Padding Distance                  | FL | 1-n |
| (0061,SIEMENS MED NM,59) | Collimator Back Spacing Distance                   | FL | 1-n |
| (0061,SIEMENS MED NM,5A) | Collimator Mean Hole Area                          | FL | 1-n |
| (0061,SIEMENS MED NM,5B) | Collimator Field of view                           | FL | 1-n |
| (0061,SIEMENS MED NM,5C) | Collimator Septal Penetration                      | FL | 1-n |
| (0061,SIEMENS MED NM,5D) | Collimator Sensitivity                             | FL | 1-n |
| (0061,SIEMENS MED NM,5E) | Crystal Distance to the Interaction Plane          | FL | 1-2 |
| (0061,SIEMENS MED NM,5F) | Crystal Intrinsic Resolution                       | FL | 1-n |
| (0061,SIEMENS MED NM,60) | Detector 1 Focus Point - IQSPECT Heart Offset      | FL | 1-n |
| (0061,SIEMENS MED NM,61) | Detector 2 Focus Point - IQSPECT Heart Offset      | FL | 1-n |
| (0061,SIEMENS MED NM,62) | Recon Output Type                                  | LT | 1   |
| (0061,SIEMENS MED NM,63) | NM Recon Input CT UID                              | LT | 1   |
| (0061,SIEMENS MED NM,64) | NM Recon Input NM UID                              | LT | 1   |
| (0061,SIEMENS MED NM,65) | NM Recon Input Raw UID                             | LT | 1   |
| (0061,SIEMENS MED NM,66) | Recon Decay Corrected Date and Time                | LT | 1   |
| (0061,SIEMENS MED NM,67) | Recon Attenuation Correction Temporal Relationship | LT | 1   |
| (0061,SIEMENS MED NM,68) | Recon Attenuation Correction Source                | LT | 1   |
| (0061,SIEMENS MED NM,69) | Detector Normalization Correction                  | SS | 1   |
| (0061,SIEMENS MED NM,6A) | Recon Sensitivity Calibrated                       | SS | 1   |
| (0061,SIEMENS MED NM,6B) | Recon Image Filter                                 | LT | 1   |
| (0061,SIEMENS MED NM,6C) | Number of Subsets                                  | SL | 1   |
| (0061,SIEMENS MED NM,6D) | Number of Iterations                               | SL | 1   |
| (0061,SIEMENS MED NM,6E) | Recon Iterative Method                             | LT | 1   |
| (0061,SIEMENS MED NM,6F) | Reconstruction Angle                               | FL | 2   |
| (0061,SIEMENS MED NM,70) | NM Reconstruction Algorithm                        | LT | 1   |
| (0061,SIEMENS MED NM,71) | CT Transformation Matrix                           | FD | 16  |
| (0061,SIEMENS MED NM,72) | Adv Acquisition Zone Paramaters                    | FD | n   |
| (0061,SIEMENS MED NM,73) | Adv Acquisition Zone Paramaters DataType           | LT | 1   |
| (0061,SIEMENS MED NM,74) | Advanced Acq Frame Start Time                      | SL | 1-n |
| (0061,SIEMENS MED NM,75) | Advanced Acq Frame End Time                        | SL | 1-n |
| (0061,SIEMENS MED NM,76) | QC Motion Vector                                   | FL | 1-n |
| (0061,SIEMENS MED NM,77) | Acquisition Uncorrected Raw Frame Series UID       | LT | 1   |
| (0061,SIEMENS MED NM,78) | Acquisition Uncorrected Raw Frame Image UID        | LT | 1   |

| DICOM Tag                | Name                                          | VR | VM  |
|--------------------------|-----------------------------------------------|----|-----|
| (0061,SIEMENS MED NM,79) | NM Recon Recon Type                           | LT | 1   |
| (0061,SIEMENS MED NM,7A) | Initial Assay Dose                            | FD | 1   |
| (0061,SIEMENS MED NM,7B) | Assay Dose Date Time                          | DT | 1   |
| (0061,SIEMENS MED NM,7C) | Effective Dose - OBSOLETE                     | FD | 1   |
| (0061,SIEMENS MED NM,7D) | Residual Dose                                 | FD | 1   |
| (0061,SIEMENS MED NM,7E) | Residual Dose Date Time                       | DT | 1   |
| (0061,SIEMENS MED NM,7F) | Recon Volume Sensitivity                      | FD | 1   |
| (0061,SIEMENS MED NM,80) | Recon Parameter Block                         | LT | 1   |
| (0061,SIEMENS MED NM,81) | Legacy Corrected Series UID                   | LT | 1   |
| (0061,SIEMENS MED NM,82) | Legacy Corrected Image UID                    | LT | 1   |
| (0061,SIEMENS MED NM,83) | Collimator Septal Thickness                   | FL | 1-n |
| (0061,SIEMENS MED NM,85) | View Start Times                              | DT | 1-n |
| (0061,SIEMENS MED NM,86) | View Pause Durations                          | SL | 1-n |
| (0061,SIEMENS MED NM,87) | Reconstruction Performance Range              | SL | 1   |
| (0061,SIEMENS MED NM,88) | Injection Dose Date Time                      | DT | 1   |
| (0061,SIEMENS MED NM,89) | Effective Dose Date Time                      | DT | 1   |
| (0061,SIEMENS MED NM,8A) | Sensitivity Calibration Distance (Detector 1) | FD | 1   |
| (0061,SIEMENS MED NM,8B) | Sensitivity Calibration Distance (Detector 2) | FD | 1   |
| (0061,SIEMENS MED NM,8C) | QSPECT UTC Time zone offset                   | LO | 1   |
| (0061,SIEMENS MED NM,8D) | QSPECT Flag                                   | SS | 1   |
| (0063,SIEMENS MED NM,01) | System Sensitivity Detector 1                 | FL | 1-n |
| (0063,SIEMENS MED NM,02) | System Sensitivity Detector 2                 | FL | 1-n |
| (0063,SIEMENS MED NM,03) | Assay Dose Detector 1                         | FL | 1-n |
| (0063,SIEMENS MED NM,04) | Assay Dose Detector 2                         | FL | 1-n |
| (0063,SIEMENS MED NM,05) | Residual Dose Detector 1                      | FL | 1-n |
| (0063,SIEMENS MED NM,06) | Residual Dose Detector 2                      | FL | 1-n |
| (0063,SIEMENS MED NM,07) | Count Pass Detector 1                         | FL | 1-n |
| (0063,SIEMENS MED NM,08) | Count Pass Detector 2                         | FL | 1-n |
| (0063,SIEMENS MED NM,09) | Volume Sensitivity Factor                     | FD | 1-n |
| (0063,SIEMENS MED NM,0A) | Volume Sensitivity Factor Version             | LO | 1-n |
| (0063,SIEMENS MED NM,0B) | Volume Sensitivity Factor Volume              | FD | 1-n |
| (0063,SIEMENS MED NM,0C) | xSpect Quant Source Labels                    | LO | 1-n |
| (0063,SIEMENS MED NM,0D) | Broadquant Volume Array                       | FL | 1-n |
| (0063,SIEMENS MED NM,0E) | Broadquant Organ Names                        | LO | 1-n |
| (0063,SIEMENS MED NM,0F) | Broadquant Zone Values                        | SS | 1-n |
| (0063,SIEMENS MED NM,10) | Broadquant Zone Statistics                    | FL | 1-n |
| (0063,SIEMENS MED NM,11) | Isotope Half Life                             | FL | 1-n |

| DICOM Tag                | Name                                | VR | VM  |
|--------------------------|-------------------------------------|----|-----|
| (0063,SIEMENS MED NM,12) | System Sensitivity Date Time        | DT | 1   |
| (0063,SIEMENS MED NM,13) | Assay Sensitivity Date Time         | DT | 1   |
| (0063,SIEMENS MED NM,14) | Residual Sensitivity Date Time      | DT | 1   |
| (0063,SIEMENS MED NM,15) | Front End Count Rate for Detector 1 | FL | 1-n |
| (0063,SIEMENS MED NM,16) | Front End Count Rate for Detector 2 | FL | 1-n |
| (0063,SIEMENS MED NM,80) | Assay Dose List                     | FL | 1-n |
| (0063,SIEMENS MED NM,81) | Assay Dose Date Time List           | DT | 1-n |
| (0065,SIEMENS MED NM,01) | Original Detector Index             | SS | 1   |
| (0065,SIEMENS MED NM,02) | Siemens Planar Data Organization    | LO | 1   |
| (6001,SIEMENS MED NM,00) | Data 1 for ROI 1                    | US | 1-n |
| (6001,SIEMENS MED NM,01) | Data 1 for ROI 2                    | US | 1-n |
| (6001,SIEMENS MED NM,02) | Data 1 for ROI 3                    | US | 1-n |
| (6001,SIEMENS MED NM,03) | Data 1 for ROI 4                    | US | 1-n |
| (6001,SIEMENS MED NM,04) | Data 1 for ROI 5                    | US | 1-n |
| (6001,SIEMENS MED NM,05) | Data 1 for ROI 6                    | US | 1-n |
| (6001,SIEMENS MED NM,06) | Data 1 for ROI 7                    | US | 1-n |
| (6001,SIEMENS MED NM,07) | Data 1 for ROI 8                    | US | 1-n |
| (6001,SIEMENS MED NM,08) | Data 1 for ROI 9                    | US | 1-n |
| (6001,SIEMENS MED NM,09) | Data 1 for ROI 10                   | US | 1-n |
| (6001,SIEMENS MED NM,0A) | Data 1 for ROI 11                   | US | 1-n |
| (6001,SIEMENS MED NM,0B) | Data 1 for ROI 12                   | US | 1-n |
| (6001,SIEMENS MED NM,0C) | Data 1 for ROI 13                   | US | 1-n |
| (6001,SIEMENS MED NM,0D) | Data 1 for ROI 14                   | US | 1-n |
| (6001,SIEMENS MED NM,0E) | Data 1 for ROI 15                   | US | 1-n |
| (6001,SIEMENS MED NM,0F) | Data 1 for ROI 16                   | US | 1-n |
| (6001,SIEMENS MED NM,10) | Data 2 for ROI 1                    | FL | 1-n |
| (6001,SIEMENS MED NM,11) | Data 2 for ROI 2                    | FL | 1-n |
| (6001,SIEMENS MED NM,12) | Data 2 for ROI 3                    | FL | 1-n |
| (6001,SIEMENS MED NM,13) | Data 2 for ROI 4                    | FL | 1-n |
| (6001,SIEMENS MED NM,14) | Data 2 for ROI 5                    | FL | 1-n |
| (6001,SIEMENS MED NM,15) | Data 2 for ROI 6                    | FL | 1-n |
| (6001,SIEMENS MED NM,16) | Data 2 for ROI 7                    | FL | 1-n |
| (6001,SIEMENS MED NM,17) | Data 2 for ROI 8                    | FL | 1-n |
| (6001,SIEMENS MED NM,18) | Data 2 for ROI 9                    | FL | 1-n |
| (6001,SIEMENS MED NM,19) | Data 2 for ROI 10                   | FL | 1-n |
| (6001,SIEMENS MED NM,1A) | Data 2 for ROI 11                   | FL | 1-n |
| (6001,SIEMENS MED NM,1B) | Data 2 for ROI 12                   | FL | 1-n |

| DICOM Tag                         | Name                                     | VR | VM  |
|-----------------------------------|------------------------------------------|----|-----|
| (6001,SIEMENS MED NM,1C)          | Data 2 for ROI 13                        | FL | 1-n |
| (6001,SIEMENS MED NM,1D)          | Data 2 for ROI 14                        | FL | 1-n |
| (6001,SIEMENS MED NM,1E)          | Data 2 for ROI 15                        | FL | 1-n |
| (6001,SIEMENS MED NM,1F)          | Data 2 for ROI 16                        | FL | 1-n |
| (7FE3,SIEMENS MED NM,14)          | Minimum pixel value in frame             | OW | n   |
| (7FE3,SIEMENS MED NM,15)          | Maximum pixel value in frame             | OW | n   |
| (7FE3,SIEMENS MED NM,16)          | Total counts in frame                    | OW | n   |
| (7FE3,SIEMENS MED NM,17)          | View Duration                            | OW | n   |
| (7FE3,SIEMENS MED NM,1B)          | Pixel Overflow flag                      | OW | n   |
| (7FE3,SIEMENS MED NM,1C)          | Buffer Overflow                          | OW | n   |
| (7FE3,SIEMENS MED NM,1E)          | Frame start time                         | OW | n   |
| (7FE3,SIEMENS MED NM,29)          | Number of Rwaves in a frame              | OW | n   |
| (0067,SIEMENS MED MI,01)          | MI Scan ID                               | LT | 1   |
| (0067,SIEMENS MED MI,02)          | Scanner Console Generation               | LO | 1   |
| (0067,SIEMENS MED MI,03)          | Recon Parameters                         | OB | 1   |
| (0067,SIEMENS MED MI,04)          | Group Reconstruction ID                  | LO | 1   |
| (0067,SIEMENS MED MI,05)          | Device IVK                               | ST | 1   |
| (0067,SIEMENS MED MI,14)          | Raw Data Description                     | LO | 1   |
| (0067,SIEMENS MED MI,16)          | Raw Data Series instance UIDs            | UI | 1   |
| (0067,SIEMENS MED MI,17)          | Raw Data Referenced Series instance UIDs | UI | 1   |
| (0067,SIEMENS MED MI,18)          | Raw Data Blob Sequence                   | SQ | 1   |
| (0071,SIEMENS MED PT,21)          | Reference To Registration                | LO | 1   |
| (0071,SIEMENS MED PT,22)          | Decay Correction DateTime                | DT | 1   |
| (0071,SIEMENS MED PT,23)          | Registration Matrix                      | FD | 16  |
| (0071,SIEMENS MED PT,24)          | Table Motion                             | CS | 1   |
| (0071,SIEMENS MED PT,25)          | Lumped Constant                          | FD | 1   |
| (0071,SIEMENS MED PT,26)          | Histogramming Method                     | CS | 1   |
| (0071, SIEMENS MED PT MU MAP, 01) | SOP Class of Source                      | UI | 1   |
| (0071, SIEMENS MED PT MU MAP, 02) | Related Mu Map Series                    | UI | 1   |
| (0021, SIEMENS MR SDS 01, 0C)     | Positive PCS Directions                  | SH | 1   |
| (0021, SIEMENS MR SDS 01, 5E)     | Field Of View Text                       | LO | 1   |
| (0021, SIEMENS MR SDS 01, 5F)     | Relative Table Position Text             | SH | 1   |
| (0021, SIEMENS MR SDS 01, FE)     | Series Data Sequence                     | SQ | 1   |
| (0021, SIEMENS MR SDI 02, 4F)     | Coil String                              | LO | 1   |
| (0021, SIEMENS MR SDI 02, 56)     | PAT Mode Text                            | LO | 1   |
| (0021, SIEMENS MR SDI 02, 58)     | Acquisition Matrix Text                  | SH | 1   |
| (0021, SIEMENS MR SDI 02, 88)     | Slice Position                           | DS | 1   |

| DICOM Tag                     | Name                    | VR | VM |
|-------------------------------|-------------------------|----|----|
| (0021, SIEMENS MR SDI 02, 89) | Slice Position Text     | SH | 1  |
| (0021, SIEMENS MR SDI 02, FE) | Image Data Sequence     | SQ | 1  |
| (0021, SIEMENS MR SDR 01, 01) | Creator Identifier      | LO | 1  |
| (0021, SIEMENS MR SDR 01, 02) | Application Identifier  | LO | 1  |
| (0021, SIEMENS MR SDR 01, 03) | Cause Identifier        | LO | 1  |
| (0051, SIEMENS MR HEADER, 0A) | Meas Duration           | LO | 1  |
| (0051, SIEMENS MR HEADER, 0C) | Field Of View           | LO | 1  |
| (0051, SIEMENS MR HEADER, 0D) | Slice Position          | SH | 1  |
| (0051, SIEMENS MR HEADER, 12) | Rel Table Position      | SH | 1  |
| (0051, SIEMENS MR HEADER, 13) | Positive PCS Directions | SH | 1  |
| (7FE1, SIEMENS MR IMA, 10)    | Raw Data                | OB | 1  |

*Interpretation of the DICOM Tags from the above table:*

*(gggg, pp,ee) -> (gggg, ppee)*

*gggg - odd group number*

*pp - private creator identification code*

*ee - private element*

## 9.3 Coded Terminology and Templates

### 9.3.1 Context Groups

N/A

### 9.3.2 Template Specifications

N/A

### 9.3.3 Private Code definitions

N/A

## 9.4 Grayscale Image Consistency

The high resolution TFT display monitor option of syngo® MR XA50A comes with a DICOM Grayscale Standard Display Function (GSDF) compliant factory pre-setting. A typical working environment setup is assumed for ambient light.

## 9.5 Standard Extended / Specialized / Private SOP Classes

N/A

## 9.6 Private Transfer Syntaxes

No private Transfer Syntaxes are defined for or requested by syngo® MR XA50A DICOM application.

## 9.7 DICOM Print SCU – detailed status displays

The following tables document the behavior of the *syngo*® MR XA50A DICOM Print AE in response to messages received for the printer SOP class and the print job SOP class.

Definitions of camera symbols:

- Idle: Camera is installed and ready; idle icon is displayed.
- Interact: The user has to react in near future, but not immediately.  
Example: A camera was low in 8x10 clear sheets: LOW 8x10 CLR was sent by N-EVENT-REPORT.
- Queue Stopped: The user has to react immediately. Either the camera needs immediate interaction or a job has been aborted.  
Example: A camera is out of 8x10 clear sheets, or camera is down, or a film job is aborted.

Note: different camera symbols are displayed according to the Printer Status Info.

### 9.7.1 Common Status Information

Table 76: Common Printer Status Information

| Printer Status Info / Execution Status Info | Description                                                                                                                        | Message string visible in 'Status Bar'                       | Other action for UI / 'camera symbol'                   |
|---------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------|---------------------------------------------------------|
| NORMAL                                      | Camera is ready                                                                                                                    | Camera is ready                                              | <None>/idle                                             |
| BAD RECEIVE MGZ                             | There is a problem with the film receive magazine. Films from the printer cannot be transported into the magazine.                 | Problem with receive magazine.                               | <None>/interact                                         |
| BAD SUPPLY MGZ                              | There is a problem with the film supply magazine. Films from this magazine cannot be transported into the printer.                 | Problem with supply magazine.                                | <None>/interact                                         |
| CALIBRATING                                 | Printer is performing self-calibration, it is expected to be available for normal operation shortly.                               | Self-calibration. Please wait.                               | <None>/idle                                             |
| CALIBRATION ERR                             | An error in the printer calibration has been detected, quality of processed films may not be optimal.                              | Problem in calibration. Film quality may not be optimal.     | <None>/interact                                         |
| CHECK CHEMISTRY                             | A problem with the processor chemicals has been detected, quality of processed films may not be optimal.                           | Problem with chemistry.<br>Film quality may not be optimal.  | <None>/interact                                         |
| CHECK SORTER                                | There is an error in the film sorter                                                                                               | Error in film sorter.                                        | <None>/interact                                         |
| CHEMICALS EMPTY                             | There are no processing chemicals in the processor, films will not be printed and processed until the processor is back to normal. | Camera chemistry empty.<br>Please check.                     | <None>/interact                                         |
| CHEMICALS LOW                               | The chemical level in the processor is low, if not corrected, it will probably shut down soon.                                     | Camera chemistry low.<br>Please check.                       | <None>/interact                                         |
| COVER OPEN                                  | One or more printer or processor covers, drawers, doors are open.                                                                  | Camera cover, drawer or door open.                           | <None>/interact                                         |
| ELEC CONFIG ERR                             | Printer configured improperly for this job.                                                                                        | Camera configured improperly for this job.<br>Queue stopped. | Queue for this camera will be STOPPED/<br>Queue stopped |

| Printer Status Info / Execution Status Info | Description                                                                   | Message string visible in 'Status Bar'  | Other action for UI / 'camera symbol'                |
|---------------------------------------------|-------------------------------------------------------------------------------|-----------------------------------------|------------------------------------------------------|
| ELEC DOWN                                   | Printer is not operating due to some unspecified electrical hardware problem. | Camera electrical hardware Problem.     | <None>/interact                                      |
| ELEC SW ERROR                               | Printer not operating for some unspecified software error.                    | Camera software problem. Queue stopped. | Queue for this camera will be STOPPED/ Queue stopped |
| EMPTY 8X10                                  | The 8x10 inch film supply magazine is empty.                                  | 8x10 film supply empty.                 | <None>/interact                                      |
| EMPTY 8X10 BLUE                             | The 8x10 inch blue film supply magazine is empty.                             | 8x10 blue film supply empty.            | <None>/interact                                      |
| EMPTY 8X10 CLR                              | The 8x10 inch clear film supply magazine is empty.                            | 8x10 clear film supply empty.           | <None>/interact                                      |
| EMPTY 8X10 PAPR                             | The 8x10 inch paper supply magazine is empty.                                 | 8x10 paper supply empty.                | <None>/interact                                      |
| EMPTY 10X12                                 | The 10x12 inch film supply magazine is empty.                                 | 10x12 film supply empty.                | <None>/interact                                      |
| EMPTY 10X12 BLUE                            | The 10x12 inch blue film supply magazine is empty.                            | 10x12 blue film supply empty.           | <None>/interact                                      |
| EMPTY 10X12 CLR                             | The 10x12 inch clear film supply magazine is empty.                           | 10x12 clear film supply empty.          | <None>/interact                                      |
| EMPTY 10X12 PAPR                            | The 10x12 inch paper supply magazine is empty.                                | 10x12 paper supply empty.               | <None>/interact                                      |
| EMPTY 10X14                                 | The 10x14 inch film supply magazine is empty.                                 | 10x14 film supply empty.                | <None>/interact                                      |
| EMPTY 10X14 BLUE                            | The 10x14 inch blue film supply magazine is empty.                            | 10x14 blue film supply empty.           | <None>/interact                                      |
| EMPTY 10X14 CLR                             | The 10x14 inch clear film supply magazine is empty.                           | 10x14 clear film supply empty.          | <None>/interact                                      |
| EMPTY 10X14 PAPR                            | The 10x14 inch paper supply magazine is empty.                                | 10x14 paper supply empty.               | <None>/interact                                      |
| EMPTY 11X14                                 | The 11x14 inch film supply magazine is empty.                                 | 11x14 film supply empty.                | <None>/interact                                      |
| EMPTY 11X14 BLUE                            | The 11x14 inch blue film supply magazine is empty.                            | 11x14 blue film supply empty.           | <None>/interact                                      |
| EMPTY 11X14 CLR                             | The 11x14 inch clear film supply magazine is empty.                           | 11x14 clear film supply empty.          | <None>/interact                                      |
| EMPTY 11X14 PAPR                            | The 11x14 inch paper supply magazine is empty.                                | 11x14 paper supply empty.               | <None>/interact                                      |
| EMPTY 14X14                                 | The 14x14 inch film supply magazine is empty.                                 | 14x14 film supply empty.                | <None>/interact                                      |
| EMPTY 14X14 BLUE                            | The 14x14 inch blue film supply magazine is empty.                            | 14x14 blue film supply empty.           | <None>/interact                                      |
| EMPTY 14X14 CLR                             | The 14x14 inch clear film supply magazine is empty.                           | 14x14 clear film supply empty.          | <None>/interact                                      |

| Printer Status Info /<br>Execution Status<br>Info | Description                                                                           | Message string visible in<br>'Status Bar' | Other action for UI /<br>'camera symbol' |
|---------------------------------------------------|---------------------------------------------------------------------------------------|-------------------------------------------|------------------------------------------|
| EMPTY 14X14 PAPR                                  | The 14x14 inch paper supply magazine is empty.                                        | 14x14 paper supply empty.                 | <None>/interact                          |
| EMPTY 14X17                                       | The 14x17 inch film supply magazine is empty.                                         | 14x17 film supply empty.                  | <None>/interact                          |
| EMPTY 14X17 BLUE                                  | The 14x17 inch blue film supply magazine is empty.                                    | 14x17 blue film supply empty.             | <None>/interact                          |
| EMPTY 14X17 CLR                                   | The 14x17 inch clear film supply magazine is empty.                                   | 14x17 clear film supply empty.            | <None>/interact                          |
| EMPTY 14X17 PAPR                                  | The 14x17 inch paper supply magazine is empty.                                        | 14x17 paper supply empty.                 | <None>/interact                          |
| EMPTY 24X24                                       | The 24x24 inch film supply magazine is empty.                                         | 24x24 film supply empty.                  | <None>/interact                          |
| EMPTY 24X24 BLUE                                  | The 24x24 inch blue film supply magazine is empty.                                    | 24x24 blue film supply empty.             | <None>/interact                          |
| EMPTY 24X24 CLR                                   | The 24x24 inch clear film supply magazine is empty.                                   | 24x24 clear film supply empty.            | <None>/interact                          |
| EMPTY 24X24 PAPR                                  | The 24x24 inch paper supply magazine is empty.                                        | 24x24 paper supply empty                  | <None>/interact                          |
| EMPTY 24X30                                       | The 24x30 inch film supply magazine is empty.                                         | 24x30 film supply empty.                  | <None>/interact                          |
| EMPTY 24X30 BLUE                                  | The 24x30 inch blue film supply magazine is empty.                                    | 24x30 blue film supply empty.             | <None>/interact                          |
| EMPTY 24X30 CLR                                   | The 24x30 inch clear film supply magazine is empty.                                   | 24x30 clear film supply empty.            | <None>/interact                          |
| EMPTY 24X30 PAPR                                  | The 24x30 inch paper supply magazine is empty.                                        | 24x30 paper supply empty.                 | <None>/interact                          |
| EMPTY A4 PAPR                                     | The A4 paper supply magazine is empty.                                                | A4 paper supply empty                     | <None>/interact                          |
| EMPTY A4 TRANS                                    | The A4 transparency supply magazine is empty.                                         | A4 transparency supply empty.             | <None>/interact                          |
| EXPOSURE FAILURE                                  | The exposure device has failed due to some unspecified reason.                        | Exposure device has failed.               | <None>/interact                          |
| FILM JAM                                          | A film transport error has occurred and a film is jammed in the printer or processor. | Film jam.                                 | <None>/interact                          |
| FILM TRANSP ERR                                   | There is a malfunction with the film transport, there may or may not be a film jam.   | Film transport problem.                   | <None>/interact                          |
| FINISHER EMPTY                                    | The finisher is empty.                                                                | Finisher is empty.                        | <None>/interact                          |
| FINISHER ERROR                                    | The finisher is not operating due to some unspecified reason                          | Finisher problem.                         | <None>/interact                          |
| FINISHER LOW                                      | The finisher is low on supplies.                                                      | Finisher low.                             | <None>/interact                          |
| LOW 8X10                                          | The 8x10 inch film supply magazine is low.                                            | 8x10 film supply low.                     | <None>/interact                          |
| LOW 8X10 BLUE                                     | The 8x10 inch blue film supply magazine is low.                                       | 8x10 blue film supply low.                | <None>/interact                          |

| Printer Status Info /<br>Execution Status<br>Info | Description                                       | Message string visible in<br>'Status Bar' | Other action for UI /<br>'camera symbol' |
|---------------------------------------------------|---------------------------------------------------|-------------------------------------------|------------------------------------------|
| LOW 8X10 CLR                                      | The 8x10 inch clear film supply magazine is low.  | 8x10 clear film supply low.               | <None>/interact                          |
| LOW 8X10 PAPR                                     | The 8x10 inch paper supply magazine is low.       | 8x10 paper supply low.                    | <None>/interact                          |
| LOW 10X12                                         | The 10x12 inch film supply magazine is low.       | 10x12 film supply low.                    | <None>/interact                          |
| LOW 10X12 BLUE                                    | The 10x12 inch blue film supply magazine is low.  | 10x12 blue film supply low.               | <None>/interact                          |
| LOW 10X12 CLR                                     | The 10x12 inch clear film supply magazine is low. | 10x12 clear film supply low.              | <None>/interact                          |
| LOW 10X12 PAPR                                    | The 10x12 inch paper supply magazine is low.      | 10x12 paper supply low.                   | <None>/interact                          |
| LOW 10X14                                         | The 10x14 inch film supply magazine is low.       | 10x14 film supply low.                    | <None>/interact                          |
| LOW 10X14 BLUE                                    | The 10x14 inch blue film supply magazine is low.  | 10x14 blue film supply low.               | <None>/interact                          |
| LOW 10X14 CLR                                     | The 10x14 inch clear film supply magazine is low. | 10x14 clear film supply low.              | <None>/interact                          |
| LOW 10X14 PAPR                                    | The 10x14 inch paper supply magazine is low.      | 10x14 paper supply low.                   | <None>/interact                          |
| LOW 11X14                                         | The 11x14 inch film supply magazine is low.       | 11x14 film supply low.                    | <None>/interact                          |
| LOW 11X14 BLUE                                    | The 11x14 inch blue film supply magazine is low.  | 11x14 blue film supply low.               | <None>/interact                          |
| LOW 11X14 CLR                                     | The 11x14 inch clear film supply magazine is low. | 11x14 clear film supply low.              | <None>/interact                          |
| LOW 11X14 PAPR                                    | The 11x14 inch paper supply magazine is low.      | 11x14 paper supply low.                   | <None>/interact                          |
| LOW 14X14                                         | The 14x14 inch film supply magazine is low.       | 14x14 film supply low.                    | <None>/interact                          |
| LOW 14X14 BLUE                                    | The 14x14 inch blue film supply magazine is low.  | 14x14 blue film supply low.               | <None>/interact                          |
| LOW 14X14 CLR                                     | The 14x14 inch clear film supply magazine is low. | 14x14 clear film supply low.              | <None>/interact                          |
| LOW 14X14 PAPR                                    | The 14x14 inch paper supply magazine is low.      | 14x14 paper supply low.                   | <None>/interact                          |
| LOW 14X17                                         | The 14x17 inch film supply magazine is low.       | 14x17 film supply low.                    | <None>/interact                          |
| LOW 14X17 BLUE                                    | The 14x17 inch blue film supply magazine is low.  | 14x17 blue film supply low.               | <None>/interact                          |
| LOW 14X17 CLR                                     | The 14x17 inch clear film supply magazine is low. | 14x17 clear film supply low.              | <None>/interact                          |
| LOW 14X17 PAPR                                    | The 14x17 inch paper supply magazine is low.      | 14x17 paper supply low.                   | <None>/interact                          |

| Printer Status Info /<br>Execution Status<br>Info | Description                                                                                                                                         | Message string visible in<br>'Status Bar' | Other action for UI /<br>'camera symbol' |
|---------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|------------------------------------------|
| LOW 24X24                                         | The 24x24 inch film supply magazine is low.                                                                                                         | 24x24 film supply low.                    | <None>/interact                          |
| LOW 24X24 BLUE                                    | The 24x24 inch blue film supply magazine is low.                                                                                                    | 24x24 blue film supply low.               | <None>/interact                          |
| LOW 24X24 CLR                                     | The 24x24 inch clear film supply magazine is low.                                                                                                   | 24x24 clear film supply low.              | <None>/interact                          |
| LOW 24X24 PAPR                                    | The 24x24 inch paper supply magazine is low.                                                                                                        | 24x24 paper supply low.                   | <None>/interact                          |
| LOW 24X30                                         | The 24x30 inch film supply magazine is low.                                                                                                         | 24x30 film supply low.                    | <None>/interact                          |
| LOW 24X30 BLUE                                    | The 24x30 inch blue film supply magazine is low.                                                                                                    | 24x30 blue film supply low.               | <None>/interact                          |
| LOW 24X30 CLR                                     | The 24x30 inch clear film supply magazine is low.                                                                                                   | 24x30 clear film supply low.              | <None>/interact                          |
| LOW 24X30 PAPR                                    | The 24x30 inch paper supply magazine is low.                                                                                                        | 24x30 paper supply low.                   | <None>/interact                          |
| LOW A4 PAPR                                       | The A4 paper supply magazine is low.                                                                                                                | A4 paper supply low.                      | <None>/interact                          |
| LOW A4 TRANS                                      | The A4 transparency supply magazine is low.                                                                                                         | A4 transparency supply low.               | <None>/interact                          |
| NO RECEIVE MGZ                                    | The film receive magazine is not available.                                                                                                         | Film receiver not available.              | <None>/interact                          |
| NO RIBBON                                         | The ribbon cartridge needs to be replaced.                                                                                                          | Replace ribbon cartridge.                 | <None>/interact                          |
| NO SUPPLY MGZ                                     | The film supply magazine is not available.                                                                                                          | Film supply not available.                | <None>/interact                          |
| CHECK PRINTER                                     | The printer is not ready at this time, operator intervention is required to make the printer available.                                             | Check camera.                             | <None>/interact                          |
| CHECK PROC                                        | The processor is not ready at this time, operator intervention is required to make the printer available.                                           | Check processor.                          | <None>/interact                          |
| PRINTER DOWN                                      | The printer is not operating due to some unspecified reason.                                                                                        | Camera down.                              | <None>/interact                          |
| PRINTER INIT                                      | The printer is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state.   | Camera initializing.                      | <None>/Idle                              |
| PRINTER OFFLINE                                   | The printer has been disabled by an operator or service person.                                                                                     | Camera off-line.                          | <None>/interact                          |
| PROC DOWN                                         | The processor is not operating due to some unspecified reason.                                                                                      | Processor down.                           | <None>/interact                          |
| PROC INIT                                         | The processor is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state. | Processor initializing.                   | <None>/Idle                              |
| PROC OVERFLOW FL                                  | Processor chemicals are approaching the overflow full mark.                                                                                         | Processor chemicals near overflow.        | <None>/interact                          |

| Printer Status Info / Execution Status Info | Description                                                                                                                            | Message string visible in 'Status Bar'                                   | Other action for UI / 'camera symbol'                |
|---------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|------------------------------------------------------|
| PROC OVERFLOW HI                            | Processor chemicals have reached the overflow full mark.                                                                               | Processor chemicals overflow.                                            | <None>/interact                                      |
| QUEUED                                      | Print job in Queue                                                                                                                     | --                                                                       | <None>/Idle                                          |
| RECEIVER FULL                               | The film receive magazine is full.                                                                                                     | Receiver full.                                                           | <None>/interact                                      |
| REQ MED NOT INST                            | The requested film, paper, or other media supply magazine is installed in the printer but may be available with operator intervention. | Install media supply.                                                    | <None>/interact                                      |
| REQ MED NOT AVAI                            | The requested film, paper, or other media requested is not available on this printer.                                                  | Media supply not available on this camera. Queue stopped. Change camera. | Queue for this camera will be STOPPED/ Queue stopped |
| RIBBON ERROR                                | There is an unspecified problem with the print ribbon.                                                                                 | Error with print ribbon.                                                 | <None>/interact                                      |
| SUPPLY EMPTY                                | The printer is out of film.                                                                                                            | Camera out of film.                                                      | <None>/interact                                      |
| SUPPLY LOW                                  | The film supply is low.                                                                                                                | Film supply low.                                                         | <None>/interact                                      |
| UNKNOWN                                     | There is an unspecified problem.                                                                                                       | Unspecified problem with camera.                                         | <None>/interact                                      |

### 9.7.2 Additional DICOM Execution Status Information

Printer Status Info and Execution Status Info are defined terms and can therefore be extended or reduced by camera manufacturers. Therefore *syngo*® MR XA50A shall be flexible.

If any other printer status info or execution status info is received (as described in Table Table 77, *syngo*® MR XA50A will react as shown in the following table:

**Table 77: Additional Printer Status Information**

| Printer Status / Execution | Printer / Execution Status Info | Description               | Message string visible in the Job status bar | Other action for syngo / camera symbol              |
|----------------------------|---------------------------------|---------------------------|----------------------------------------------|-----------------------------------------------------|
| WARNING                    | <any other>                     | <not defined status info> | Camera info: <status info>                   | <None>/Interact                                     |
| FAILURE                    | <any other>                     | <not defined status info> | Camera info: <status info><br>Queue stopped. | Queue for this camera will be STOPPED/Queue stopped |

## Annex A: Index of Tables

|                                                                                                |    |
|------------------------------------------------------------------------------------------------|----|
| Table 1: Network Services .....                                                                | 2  |
| Table 2: Media Services.....                                                                   | 7  |
| Table 3: Implementation Identifying Information .....                                          | 7  |
| Table 4: Association Policies.....                                                             | 19 |
| Table 5: Asynchronous Nature as an Association Initiator .....                                 | 19 |
| Table 6: Presentation Context Table "Verification".....                                        | 20 |
| Table 7: Presentation Context Table "Verification".....                                        | 20 |
| Table 8: Association Policy .....                                                              | 21 |
| Table 9: Asynchronous Nature as an Association Initiator .....                                 | 21 |
| Table 10: Proposed Presentation Contexts for Storage .....                                     | 22 |
| Table 11: DICOM Command Response Status Handling Behavior .....                                | 23 |
| Table 12: DICOM Command Communication Failure Behavior .....                                   | 23 |
| Table 13: Storage C-STORE Response Status .....                                                | 24 |
| Table 14: Association Policies.....                                                            | 25 |
| Table 15: Asynchronous Nature as an Association Initiator .....                                | 25 |
| Table 16: Proposed Presentation Contexts for Storage Commitment.....                           | 26 |
| Table 17: DICOM Command Response Status Handling Behavior .....                                | 26 |
| Table 18: DICOM Command Communication Failure Behavior .....                                   | 27 |
| Table 19: Presentation Context Table "Update Flag Information".....                            | 27 |
| Table 20: Association Policies.....                                                            | 27 |
| Table 21: Asynchronous Nature as an Association Initiator .....                                | 28 |
| Table 22: Proposed Presentation Contexts for Query .....                                       | 28 |
| Table 23: Extended Negotiation as an SCU .....                                                 | 29 |
| Table 24: DICOM Command Response Status Handling Behavior .....                                | 29 |
| Table 25: DICOM Command Communication Failure Behavior .....                                   | 29 |
| Table 26: Attributes supported for Study/Series Query - SCU.....                               | 29 |
| Table 27: Proposed Presentation Contexts for Retrieve and Activity "MOVE SCU" .....            | 31 |
| Table 28: DICOM Command Response Status Handling Behavior .....                                | 31 |
| Table 29: DICOM Command Communication Failure Behavior .....                                   | 31 |
| Table 30: Association Policies.....                                                            | 31 |
| Table 31: Asynchronous Nature as an Association Initiator .....                                | 32 |
| Table 32: Proposed Presentation Contexts for Worklist.....                                     | 32 |
| Table 33: Broad Query search keys .....                                                        | 33 |
| Table 34: Modality Worklist C-Find Return keys .....                                           | 34 |
| Table 35: DICOM Command Response Status Handling Behavior .....                                | 37 |
| Table 36: DICOM Command Communication Failure Behavior .....                                   | 37 |
| Table 37: Association Policies.....                                                            | 38 |
| Table 38: Asynchronous Nature as an Association Initiator .....                                | 38 |
| Table 39: Acceptable Presentation Contexts Activity "Create MPPS".....                         | 39 |
| Table 40: MPPS N-CREATE Response Status Handling Behavior .....                                | 40 |
| Table 41: Acceptable Presentation Contexts Activity "Update MPPS".....                         | 40 |
| Table 42: MPPS N-SET Response Status Handling Behavior .....                                   | 40 |
| Table 43: Association Policies.....                                                            | 41 |
| Table 44: Asynchronous Nature as an Association Initiator .....                                | 41 |
| Table 45: Presentation Contexts for the Activity "Print Film".....                             | 42 |
| Table 46: Attributes for the N-CREATE-RQ of the Basic Film Session .....                       | 43 |
| Table 47: N-CREATE-RSP Status Handling Behavior for the Basic Film Session .....               | 43 |
| Table 48: Attributes for the N-CREATE-RQ of the Basic Film Session .....                       | 43 |
| Table 49: N-CREATE-RSP Status Handling Behavior for Basic Film Box .....                       | 44 |
| Table 50: N-ACTION-RSP Status Handling Behavior for Basic Film Box .....                       | 45 |
| Table 51: Attributes for N-SET-RQ of Basic Grayscale Image Box .....                           | 45 |
| Table 52: N-SET-RSP Status Handling Behavior for the Basic Grayscale Image Box SOP Class ..... | 46 |
| Table 53: Attributes for N-SET-RQ of Basic Color Image Box.....                                | 47 |
| Table 54: N-SET-RSP Status Handling Behavior for the Color Grayscale Image Box .....           | 47 |
| Table 55: Attributes for N-CREATE-RQ of Presentation LUT SOP Class .....                       | 48 |
| Table 56: N-CREATE-RSP Status Handling Behavior for the Presentation LUT SOP Class .....       | 48 |
| Table 57: Attributes for N-GET-RQ of the Printer SOP Class.....                                | 48 |
| Table 58: DICOM Command Communication Failure Behavior .....                                   | 48 |
| Table 59: Presentation Contexts for the Activity "Print Film".....                             | 49 |

|                                                                                           |     |
|-------------------------------------------------------------------------------------------|-----|
| Table 60: Attributes for the N-EVENT-REPORT-RQ of the Printer SOP Class .....             | 49  |
| Table 61: Presentation Contexts for the Activity "Print Management" .....                 | 49  |
| Table 62: Attributes for the N-EVENT-REPORT-RQ of the Print Job SOP Class .....           | 50  |
| Table 63: AE Titles .....                                                                 | 51  |
| Table 64: Parameter List .....                                                            | 52  |
| Table 65: Implementation Class/Version Name - Media Interchange .....                     | 54  |
| Table 66: Media - Application Profiles and Real-World Activities .....                    | 55  |
| Table 67: Transfer Syntaxes for STD-GEN-DVD-J2K, STD-GEN-BD-J2K and STD-GEN-USB-J2K ..... | 56  |
| Table 68: Single-Byte Character Sets without Code Extension .....                         | 57  |
| Table 69: Single-Byte Characters Sets with Code Extension .....                           | 58  |
| Table 70: Multi-Byte Character Sets without Code Extension .....                          | 58  |
| Table 71: Multi-Byte Character Sets with Code Extension .....                             | 58  |
| Table 72: Application Level Confidentiality Profile attributes (standard tags) .....      | 60  |
| Table 73: Application Level Confidentiality Profile Attributes (private tags) .....       | 68  |
| Table 74: Secure Transport Connection Profiles .....                                      | 70  |
| Table 75: Private Data Element Dictionary .....                                           | 110 |
| Table 76: Common Printer Status Information .....                                         | 144 |
| Table 77: Additional Printer Status Information .....                                     | 149 |

## Annex B: Table of Figures

|                                                                                        |    |
|----------------------------------------------------------------------------------------|----|
| Figure 1: <i>syngo</i> ® MR XA50A DICOM Data Flow diagram – Acquisition Workflow ..... | 13 |
| Figure 2: <i>syngo</i> ® MR XA50A DICOM Data Flow diagram - Printing .....             | 14 |
| Figure 3: Sequence Diagram for Real World Activities - System Configuration .....      | 16 |
| Figure 4: Sequence Diagram for Real World Activities -Acquisition workflow .....       | 17 |
| Figure 5: Sequence Diagram for Real World Activities - Printing .....                  | 18 |
| Figure 6: Media Interchange Application Data Flow Diagram .....                        | 53 |
| Figure 7: Sequence diagram – Media creation .....                                      | 54 |

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some / All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory- new components, we guarantee the quality of these recycled components.

---

Siemens Healthineers Headquarters  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

---

HRT HOOD05162003239645| © Siemens Healthcare GmbH, 05.2022