



Valid from
October 2020

DICOM Conformance Statement

teamplay Receiver

1 DICOM Conformance Statement Overview

The **teamplay Receiver** is designed to be integrated into an environment of medical, DICOM-based devices. The **teamplay Receiver** supports Storage of images utilizing the DICOM “Storage Service Class” and the retrieval of images from DICOM Archives utilizing the DICOM “Query/Retrieve Service Class”.

The **teamplay Receiver** provides an AET for the analytics offering functionality.

The default AET is TEAMPLAY.

It is possible to change the configuration of the AET.

The **teamplay** analytics offerings AE **teamplay** provides the following functionality:

- Receive DICOM objects from service class users
- Query remote nodes for objects based on Study Root Model and retrieve them.
- Receive data and optionally deidentifies it.

The **teamplay Receiver** conforms to the DICOM Standard and supports the network services as described in Table 1: Network Services and the media services as described in Table 2 - Media Services.

Table 1: Network Services

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|--|-------------------------------|-----------------------|------|---------------------------|---------|
| Verification | | | | | |
| Verification | 1.2.840.10008.1.1 | Yes | | Yes | |
| SOP Classes created by <i>teamplay Receiver</i> | | | | | |
| | | Create | Send | Store | Display |
| None | | | | | |
| SOP Classes managed by <i>teamplay Receiver</i> | | | | | |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | No | No | Yes | No |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | No | No | Yes | No |
| Digital X-Ray Image Storage – For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | No | No | Yes | No |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2 | No | No | Yes | No |
| Digital Mammography X-Ray Image Storage – For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | No | No | Yes | No |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | No | No | Yes | No |
| Enhanced CT Image Storage | 1.2.840.10008.5.1.4.1.1.2.1 | No | No | Yes | No |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | No | No | Yes | No |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | No | No | Yes | No |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | No | No | Yes | No |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | No | No | Yes | No |
| Enhanced MR Color Image Storage | 1.2.840.10008.5.1.4.1.1.4.3 | No | No | Yes | No |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | No | No | Yes | No |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | No | No | Yes | No |

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|--|--------------------------------|-----------------------|----|---------------------------|----|
| Multi-frame Single Bit Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.1 | No | No | Yes | No |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2 | No | No | Yes | No |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3 | No | No | Yes | No |
| Multi-frame True Color Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | No | No | Yes | No |
| 12-lead ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.1 | No | No | Yes | No |
| General ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.2 | No | No | Yes | No |
| Ambulatory ECG Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.1.3 | No | No | Yes | No |
| Hemodynamic Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.2.1 | No | No | Yes | No |
| Cardiac Electrophysiology Waveform Storage | 1.2.840.10008.5.1.4.1.1.9.3.1 | No | No | Yes | No |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 | No | No | Yes | No |
| Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.2 | No | No | Yes | No |
| Pseudo-Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.3 | No | No | Yes | No |
| Blending Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.4 | No | No | Yes | No |
| X-Ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | No | No | Yes | No |
| Enhanced XA Image Storage | 1.2.840.10008.5.1.4.1.1.12.1.1 | No | No | Yes | No |
| X-Ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | No | No | Yes | No |
| Enhanced XRF Image Storage | 1.2.840.10008.5.1.4.1.1.12.2.1 | No | No | Yes | No |
| X-Ray 3D Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.1 | No | No | Yes | No |
| X-Ray 3D Craniofacial Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.2 | No | No | Yes | No |
| Breast Tomosynthesis Image Storage | 1.2.840.10008.5.1.4.1.1.13.1.3 | No | No | Yes | No |
| Raw Data Storage | 1.2.840.10008.5.1.4.1.1.66 | No | No | Yes | No |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | No | No | Yes | No |
| Spatial Fiducials Storage | 1.2.840.10008.5.1.4.1.1.66.2 | No | No | Yes | No |
| Deformable Spatial Registration SOP Class | 1.2.840.10008.5.1.4.1.1.66.3 | No | No | Yes | No |
| Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.4 | No | No | Yes | No |
| Surface Segmentation Storage | 1.2.840.10008.5.1.4.1.1.66.5 | No | No | Yes | No |
| Real World Value Mapping Storage | 1.2.840.10008.5.1.4.1.1.67 | No | No | Yes | No |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | No | No | Yes | No |

| SOP Classes | SOP Class UID | User of Service (SCU) | | Provider of Service (SCP) | |
|---|-------------------------------|-----------------------|----|---------------------------|----|
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 | No | No | Yes | No |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | No | No | Yes | No |
| Procedure Log Storage Storage | 1.2.840.10008.5.1.4.1.1.88.40 | No | No | Yes | No |
| Mammography CAD SR Storage | 1.2.840.10008.5.1.4.1.1.88.50 | No | No | Yes | No |
| Key Object Selection Document Storage | 1.2.840.10008.5.1.4.1.1.88.59 | No | No | Yes | No |
| X-Ray Radiation Dose SR Storage | 1.2.840.10008.5.1.4.1.1.88.67 | No | No | Yes | No |
| Encapsulated PDF Storage SOP Class | 1.2.840.10008.5.1.4.1.1.104.1 | No | No | Yes | No |
| Positron Emission Tomography Image Storage | 1.2.840.10008.5.1.4.1.1.128 | No | No | Yes | No |
| RT Image Storage | 1.2.840.10008.5.1.4.1.1.481.1 | No | No | Yes | No |
| RT Dose Storage | 1.2.840.10008.5.1.4.1.1.481.2 | No | No | Yes | No |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | No | No | Yes | No |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | No | No | Yes | No |
| RT Plan Storage | 1.2.840.10008.5.1.4.1.1.481.5 | No | No | Yes | No |
| RT Treatment Summary Record Storage | 1.2.840.10008.5.1.4.1.1.481.7 | No | No | Yes | No |
| Digital Intra-Oral X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.3 | No | No | Yes | No |
| Digital Intra-Oral X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.3.1 | No | No | Yes | No |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | No | No | Yes | No |
| Enhanced PET Image Storage | 1.2.840.10008.5.1.4.1.1.130 | No | No | Yes | No |
| Enhanced US Volume Storage | 1.2.840.10008.5.1.4.1.1.6.2 | No | No | Yes | No |
| Transfer (Private SOP Class) | | | | | |
| Syngo Non-Image Storage | 1.3.12.2.1107.5.9.1 | Yes | | Yes | |
| Storage Commitment | | | | | |
| Storage Commitment Push Model SOP Class | 1.2.840.10008.1.20.1 | No | | No | |
| Worklist Management | | | | | |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | No | | No | |
| Modality Performed Procedure Step SOP Class | 1.2.840.10008.3.1.2.3.3 | No | | No | |
| Query/Retrieve | | | | | |
| Patient Root Q/R Information Model - FIND | 1.2.840.10008.3.1.2.3.3 | No | | No | |
| Patient Root Q/R - Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.1.2 | No | | No | |
| Study Root Q/R - Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | | No | |
| Study Root Q/R - Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | | No | |
| Patient/Study Only Q/R - Information Model FIND | 1.2.840.10008.5.1.4.1.2.2.1 | No | | No | |

| SOP Classes | SOP Class UID | User of Service (SCU) | Provider of Service (SCP) |
|---|-----------------------------|-----------------------|---------------------------|
| Patient/Study Only Q/R - Information Model MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | No | No |
| Print Management | | | |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | No | No |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 | No | No |
| Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | No | No |
| Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | No | No |
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | No | No |
| Basic Color Image Box SOP Class | 1.2.840.10008.5.1.1.4.1 | No | No |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | No | No |
| Print Job SOP Class | 1.2.840.10008.5.1.1.14 | No | No |
| Presentation LUT SOP Class | 1.2.840.10008.5.1.1.23 | No | No |

Table 2 - Media Services

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| Compact Disk - Recordable | | |
| STD-GEN-CD | No | No |
| AUG-GEN-CD | No | No |
| | | |
| DVD | | |
| AUG-GEN-DVD | No | No |
| AUG- GEN-DVD-J2K | No | No |
| STD-GEN-DVD | No | No |
| STD-GEN-DVD-J2K | No | No |
| | | |
| USB | | |
| AUG- GEN-USB-J2K | No | No |
| STD-GEN-USB-J2K | No | No |
| | | |

Table 3 - Implementation Identifying Information

| Name | Value |
|--------------------------|------------------------|
| Application Context Name | 1.2.840.100008.3.1.1.1 |

2 Table of Contents

| | | |
|------------|--|-----------|
| 1 | <i>DICOM Conformance Statement Overview</i> | 2 |
| 2 | <i>Table of Contents.....</i> | 6 |
| 3 | <i>Introduction.....</i> | 8 |
| 3.1 | Audience..... | 8 |
| 3.2 | Remarks..... | 8 |
| 3.3 | Definitions, Terms and Abbreviations | 9 |
| 3.4 | References | 9 |
| 4 | <i>Networking</i> | 10 |
| 4.1 | Implementation Model | 10 |
| 4.1.1 | Application Data Flow | 11 |
| 4.1.2 | Functional Definitions of Application Entities | 12 |
| 4.1.3 | Sequencing of Activities | 13 |
| 4.2 | Application Entity Specification | 16 |
| 4.2.1 | teamplay AE Specification | 16 |
| 4.3 | Network Interfaces | 21 |
| 4.3.1 | Physical Network Interface | 21 |
| 4.3.2 | Additional Protocols | 21 |
| 4.3.3 | IPv4 and IPv6 Support | 21 |
| 4.4 | Configuration..... | 21 |
| 4.4.1 | AE Title/Presentation Address Mapping | 21 |
| 4.4.2 | Parameters | 22 |
| 5 | <i>Media Interchange</i> | 23 |
| 5.1 | Implementation Model | 23 |
| 5.1.1 | Application Data Flow Diagram..... | 23 |
| 5.1.2 | Functional definitions of AEs..... | 23 |
| 5.1.3 | Sequencing of Real-World Activities | 23 |
| 5.1.4 | File Meta Information for Implementation Class and Version | 23 |
| 5.2 | AE SPECIFICATIONS | 23 |
| 5.2.1 | Media Storage AE – Specification | 23 |
| 5.3 | AUGMENTED AND PRIVATE APPLICATION PROFILES..... | 23 |
| 5.3.1 | Augmented Application Profiles | 23 |
| 5.4 | MEDIA CONFIGURATION..... | 24 |
| 6 | <i>Support of Extended Character Sets.....</i> | 25 |
| 6.1 | Character sets for teamwork Receiver..... | 25 |
| 7 | <i>Attribute confidentiality profiles.....</i> | 27 |
| 7.1 | De-identification | 27 |

| | | |
|----------|---|-----------|
| 8 | <i>Security</i> | 27 |
| 8.1 | Security Profiles | 27 |
| 8.2 | Association Level Security | 27 |
| 8.3 | Application Level Security | 27 |
| 9 | <i>Annexes</i> | 27 |
| 9.1 | IOD Contents | 27 |
| 9.2 | Data Dictionary of Private Attributes | 28 |
| 9.3 | Coded Terminology and Templates | 28 |
| 9.3.1 | Context Groups | 28 |
| 9.3.2 | Template Specifications | 28 |
| 9.3.3 | Private Code definitions | 28 |
| 9.4 | Grayscale Image Consistency | 28 |
| 9.5 | Standard Extended / Specialized / Private SOP Classes | 28 |
| 9.6 | Private Transfer Syntaxes..... | 28 |

3 Introduction

3.1 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

3.2 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between **teamplay Receiver** and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard (Digital Imaging and Communication in Medicine (DICOM)). DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of conformance statements is the first step towards assessing interconnectivity and interoperability between **teamplay Receiver** and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens Healthineers representative for the most recent product information.

3.3 Definitions, Terms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

| | |
|-------|--|
| AE | DICOM Application Entity |
| AET | Application Entity Title |
| ASCII | American Standard Code for Information Interchange |
| DCS | DICOM Conformance Statement |
| DICOM | Digital Imaging and Communications in Medicine |
| FSC | File Set Creator |
| FSR | File Set Reader |
| FSU | File Set Updater |
| GSDF | Grayscale Standard Display Function |
| IOD | DICOM Information Object Definition |
| ISO | International Standard Organization |
| n. a. | not applicable |
| NEMA | National Electrical Manufacturers Association |
| O | Optional Key Attribute |
| PDU | DICOM Protocol Data Unit |
| R | Required Key Attribute |
| SCU | DICOM Service Class User (DICOM client) |
| SCP | DICOM Service Class Provider (DICOM Server) |
| SOP | DICOM Service-Object Pair |
| SR | Structured Report |
| TFT | Thin Film Transistor (Display) |
| TID | Template ID |
| U | Unique Key Attribute |
| UID | Unique Identifier |
| UTF-8 | Unicode Transformation Format-8 |
| VR | Value Representation |

3.4 References

[1] NEMA PS3 / ISO 12052, Digital Imaging and Communications in Medicine (DICOM) Standard, National Electrical Manufacturers Association, Rosslyn, VA, USA (available free at <https://www.dicomstandard.org/>)

4 Networking

4.1 Implementation Model

- **Verification**

Verification requests will be processed and responded by the **teampplay** AE. The AE can also initiate an association and request verification to a remote AE.

- **Storage**

The **teampplay** AE Storage SCP starts to receive the Information Objects and store them on file system after accepting an association with a negotiated Presentation Context.

- **Query**

The **teampplay** AE supports the Query/Retrieve services in a SCU role. The C-FIND request to the remote SCP is invoked by teampplay.

- **Retrieve**

The **teampplay** AE initiates a C-MOVE request to the remote Retrieve SCP. The remote Retrieve SCP in turn starts C-STORE sub operations to the teampplay Storage SCP. The **teampplay** AE supports Study Root - Query/Retrieve Information model as SCU.

4.1.1 Application Data Flow

The following figures provide a functional overview of the **teampplay** Application Entity (AE). Relationships are shown between user-invoked activities (in the circles at the left of the AEs) and the associated real-world activities provided by DICOM service providers (in the circles at the right of the AEs)




Figure 4.1-1: Application Data Flow Diagram teampplay AE

- The teampplay DICOM service opens an association when a "verification" of a remote application is requested during a configuration session. This can be done when entering new data to configure a remote application or to verify existing configuration data.
- The Storage SCP of the **teampplay** AE can receive incoming DICOM images and add them to the local database. It can respond to external Storage and Verification Requests as a Service Class Provider (SCP) for C-STORE and C-ECHO requests.
- The Query part of the Query SCU **teampplay** AE uses C-FIND to search a remote SCP for Patient Study and Series information.
- The **teampplay** AE initiates a C-MOVE request to the remote Retrieve SCP. The remote Retrieve SCP in turn starts C-STORE sub operations to the teampplay Storage SCP.

4.1.2 Functional Definitions of Application Entities

The **teampplay Receiver** operate as background server processes. Upon accepting an association with a negotiated Presentation Context it starts to receive and process requests.

4.1.2.1 Functional Definition of teampplay AE

The **teampplay** AE can initiate verification requests to remote AE title as well as respond to verification requests from remote AEs.

The **teampplay** AE Storage SCP receives the Information Objects and stores them on file system after accepting an association with a negotiated Presentation Context. Once stored on the file system a file uploader will upload the same to cloud storage and the local copy on file system will be deleted.

The **teampplay** AE supports the Study Root Query/Retrieve Information Model - FIND as an SCU.

The C-FIND request to a remote SCP is invoked directly by **teampplay** AE for analytics offerings.

The remote SCP returns a list of responses with defined data (which can be displayed to the user).

The **teampplay** AE initiates a C-MOVE request to the remote Retrieve SCP. The remote Retrieve SCP in turn starts C-STORE sub operations to the **teampplay** AE Storage SCP. The **teampplay** AE supports the Study Root Query/Retrieve Information Model - MOVE as an SCU.

For supported SOP Classes or Data Objects refer to the "[Conformance Statement Overview](#)".

4.1.3 Sequencing of Activities

4.1.3.1 Verification

The communication between **teamplay** AE and an external DICOM node in case of Verify is depicted in Figure 4.1-2 in more detail.




Figure 4.1-2: Sequence diagram – Verify

4.1.3.2 Storage

The communication between **teampay** AE and an external DICOM node in case of triggering the transfer of objects to the teampay Receiver from the external node is depicted in Figure 4.1-3 in more detail.




Figure 4.1-3: Sequence diagram – Storage

4.1.3.3 Query / Retrieve

The communication between the **teampplay** AE and an external DICOM node in case of querying for objects from a remote DICOM node and retrieving to the teampplay Receiver is depicted in Figure 4.1-4 in more detail.




Figure 4.1-4: Sequence diagram – Query/Retrieve

4.2 Application Entity Specification

This section outlines the specifications for the Application Entity that is part of the *teampay* Receiver.

4.2.1 *teampay* AE Specification

4.2.1.1 SOP Classes

The **teampay** AE provides standard conformance to the SOP Class listed in “Table 1: Network Services” section “Verification” in the [“Conformance Statement Overview”](#).

4.2.1.2 Association Policy

The association policies for the **teampay** AE listed in table 4-1.

Table 4-1: Association Policies

| | |
|--|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
| PDU size | 64 kB |
| Maximum number of simultaneous associations as an association acceptor | Unlimited |
| Maximum number of simultaneous associations as an association initiator | Unlimited |

4.2.1.2.1 Asynchronous Nature

The **teampay Receiver** DICOM software does not support asynchronous communication (multiple outstanding transactions over a single association).

4.2.1.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

4.2.1.3 Association Initiation Policy

4.2.1.3.1 Activity – “Send Verification” Request

4.2.1.3.1.1 Description and Sequencing of Activities

The Verification SCU C-ECHO request is initiated by Service and Configuration SW whenever "Verification" is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than "Success", will be indicated in the service environment and the association is closed.

4.2.1.3.1.2 Proposed Presentation Contexts

The **teampplay Receiver** DICOM application will propose Presentation Contexts as shown in the following table:

Table 4-2 - Presentation Context Table "Verification"

| Presentation Context Table – “Verification” | | | | | |
|---|-------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.3.1.3 SOP specific Conformance for SOP classes – Verification SCU

The Application conforms to the definition of the Verification SCU in accordance to the DICOM Standard.

4.2.1.3.2 Activities “Query a Remote Node for Instances” / “Retrieve Instances from a remote node”

4.2.1.3.2.1 Description and Sequencing of Activities

The **teampplay** AE serves as a SCU for the following SOP Classes

- Study Root Q/R Information Model –FIND SOP Class
The associated Real-World activity is a C-Find request initiated by the user **teampplay** AE. The user specifies some attributes and will send a C-Find request (according to the query model) and will then return the results to the initiating application.
- Study Root Q/R Information Model –MOVE SOP Class
The Retrieve request(C-MOVE) is initiated for the C-FIND results.

4.2.1.3.2.2 Proposed Presentation Contexts

The **teamplay** AE will propose Presentation Contexts as shown in the following table:

Table 4-3: Proposed Presentation Contexts for Query/Retrieve

| Presentation Context Table | | | | | |
|---|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax | | Transfer Syntax | | Role | Ext. Neg. |
| Name | UID | Name List | UID List | | |
| Study Root Query/ Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Study Root Query/ Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.2.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | No |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.3.2.3 SOP Specific Conformance Statement to Query / Retrieve SOP classes

The **teamplay** AE supports the following query levels:

- Study

Matching Keys on Series and Instance Levels are not supported by **teamplay** Receiver as SCU.

The **teamplay** AE checks for the following status codes in the Query SCP's C-FIND-Response:

Table 4-4: C-FIND DICOM status codes

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|-----------------------|--------------------------------|
| Error | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process | Any; none; null; Code | Log message is created |
| Pending | All optional keys are supported the same manner as Required Keys. | FF00 | Log message is created |
| | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01 | Log message is created |
| Success | Query has been performed successfully. | 0000 | Query request is sent success- |

The **teamplay** AE checks for the following status codes in the Retrieve SCP's C-MOVE-Response:

Table 4-5: C-MOVE DICOM status codes

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|------------|---|
| Refused | Out of Resources - Unable to calculate number of matches | A701 | Retrieve will be retried again after some time. |
| | Out of Resources - Unable to perform sub operations | A702 | Retrieve will be retried again after some time. |
| | Move destination unknown | A801 | Retrieve will be retried again after some time. |
| Failed | Identifier does not match SOP Class | A900 | N.A. |
| | Unable to process | Cxxx | Retrieve will be retried again after some time. |
| Cancel | Sub-operations terminated due to Cancel Indication | FE00 | Retrieve will be retried again after some time. |
| Warning | Sub-operations Complete - One or more Failures of Warnings | B000 | Data is processed for teamplay . |
| Success | Sub-operations Complete - No Failures or Warning | 0000 | Data is processed for teamplay . |
| Pending | Sub-operations are continuing | FF00 | Wait for completion. |

Table 4-6: Query/Retrieve DICOM Command Communication Failure Behaviour

| Exception | Behavior |
|---------------------|--|
| Timeout | Log message is created. (default timeout value is 30 seconds) |
| Association Aborted | Log message is created. |

4.2.1.4 Association Acceptance Policy

4.2.1.4.1 Activity – “Receive Storage Request”

4.2.1.4.1.1 Description and Sequencing of Activities

The **teamplay** AE receiving process will accept an association, receive any objects transmitted on that association and store the objects on disk.

4.2.1.4.1.2 Accepted Presentation Contexts

For all supported Transfer objects (see SOP Classes in the [“Conformance Statement Overview”](#)).the Transfer Syntaxes described in Table 4-7 are supported.

Generally all Presentation Contexts are accepted as long as they contain at least one suitable Transfer syntax specified in table 4.3.

Table 4-7: Available Transfer Syntax

| Presentation Context | Transfer Syntax UID |
|--|------------------------|
| Explicit Value Representation Little Endian | 1.2.840.10008.1.2.1 |
| Implicit Value Representation Little Endian | 1.2.840.10008.1.2 |
| Explicit Value Representation Big Endian | 1.2.840.10008.1.2.2 |
| JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) | 1.2.840.10008.1.2.4.70 |
| JPEG 2000 Image Compression (Lossless Only) | 1.2.840.10008.1.2.4.90 |
| RLE Lossless | 1.2.840.10008.1.2.5 |
| JPEG Extended (Process 2 & 4) | 1.2.840.10008.1.2.4.51 |
| JPEG Baseline (Process 1) | 1.2.840.10008.1.2.4.50 |
| JPEG 2000 Image Compression | 1.2.840.10008.1.2.4.91 |

4.2.1.4.1.3 SOP specific Conformance for SOP classes

In case of a successful C-STORE operation, the image has successfully been written to disk. The Storage **teamplay** AE return the status “success” when the data is stored to disk and a minimal image header validation has been performed. After a successful storage to disk, the upload functionality will upload the same data into the **teamplay** cloud blob storage using non-DICOM communication. The copy on the disk will be deleted upon successful upload.

The following header attributes must be available and filled:

- Patient Name (0010,0010)
- Study Instance UID (0020,000D)
- Series Instance UID (0020,000E) and
- SOP Instance UID (0008, 0018).

The C-STORE response codes that are available for the Storage are described in the table 4-8.

Table 4-8: Storage C-STORE Response Status

| Service Status | Further Meaning | Error Code | Reason |
|----------------|----------------------------------|------------|--|
| Success | Success | 0000 | Image received correctly (success notification is done after receiving, before indexing and storing) |
| Failure | Out-of-resource | A700 | No resource left in the storage |
| Failure | Unable to Process | Cxxx | Error during instance reception |
| Failure | Dataset does not match SOP Class | A9xx | The Dataset is not conform to the SOP Class contained in the resource. |

4.3 Network Interfaces

4.3.1 Physical Network Interface

The **teamplay Receiver** provides DICOM 3.0 TCP/IP network communication support as defined in Part 8 of [1]. The network communication is independent from the physical medium over which TCP/IP executes; it inherits this from the Windows OS system upon which it executes.

4.3.2 Additional Protocols

N/A

4.3.3 IPv4 and IPv6 Support

Only IPv4 support is provided in this version.

4.4 Configuration

4.4.1 AE Title/Presentation Address Mapping

AE Titles shall be unique within the hospital. A common way to achieve that is to use the hostname as part of the AE Titles. The string can be up to 16 characters and must not contain any extended characters. Only 7-bit ASCII characters (excluding Control Characters) are allowed according to DICOM [1].

4.4.1.1 Local AE Titles

The **teamplay Receiver** allows configuring the AETitle and the related Port.

Table 4-9: Default Local Analytics teamwork AE Title

| Parameter | Configurable | Default Value |
|------------------|--------------|--|
| Default AE title | Yes | hostname in uppercase characters; limited to 16 characters |
| Default Port | Yes | 104 |

4.4.1.2 Remote AE Title / Presentation Address Mapping

4.4.1.2.1 Remote SCU's

All relevant remote applications that may initiate DICOM associations with teamplay Receiver need to be configured in the **teamplay Receiver**, before the association can be established.

The mapping of external AE Titles to TCP/IP addresses and ports is configurable and initially set at the time of installation. Changes can also be performed later on. The Application Entity Title and supported transfer syntaxes need to be known for configuration.

4.4.1.2.2 Remote SCP's

Remote applications can accept DICOM associations from the **teamplay Receiver**, the following information needs to be available:

- Application Entity Title
- Host Name / IP address on which the remote application service runs
- Port number on which the remote application accepts association requests.

4.4.2 Parameters

There are no additional parameters to be configured.

5 Media Interchange

5.1 Implementation Model

The *teamplay Receiver* does not support Media Storage.

5.1.1 Application Data Flow Diagram

N/A

5.1.2 Functional definitions of AEs

N/A

5.1.3 Sequencing of Real-World Activities

N/A

5.1.4 File Meta Information for Implementation Class and Version

N/A

5.2 AE SPECIFICATIONS

5.2.1 Media Storage AE – Specification

N/A

5.2.1.1 Real-World Activities

N/A

5.2.1.2 SOP Classes and Transfer Syntaxes

N/A

5.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

N/A

5.3.1 Augmented Application Profiles

N/A

5.4 MEDIA CONFIGURATION

N/A

6 Support of Extended Character Sets

6.1 Character sets for teamplay Receiver

The **teamplay Receiver** DICOM application supports the following character sets as defined in the four tables below.

Table 6-1: Single-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|----------------------|
| Default repertoire | None | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 1 | ISO_IR 100 | ISO_IR 100 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 2 | ISO_IR 101 | ISO_IR 101 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 3 | ISO_IR 109 | ISO_IR 109 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 4 | ISO_IR 110 | ISO_IR 110 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Cyrillic | ISO_IR 144 | ISO_IR 144 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Arabic | ISO_IR 127 | ISO_IR 127 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Greek | ISO_IR 126 | ISO_IR 126 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Hebrew | ISO_IR 138 | ISO_IR 138 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Latin alphabet No. 5 | ISO_IR 148 | ISO_IR 148 | Supplementary set |
| | | ISO_IR 6 | ISO 646 |
| Japanese | ISO_IR 13 | ISO_IR 13 | JIS X 0201: Katakana |
| | | ISO_IR 14 | JIS X 0201: Romaji |
| Thai | ISO_IR166 | ISO_IR166 | TIS 620-253 (1990) |
| | | ISO_IR 6 | ISO 646 |

Table 6-2: Single-Byte Characters Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------|-------------------------|-------------------|
| Default repertoire | ISO 2022 IR 6 | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.1 | ISO 2022 IR 100 | ISO 2022 | ESC 02/13 04/01 | ISO-IR 100 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.2 | ISO 2022 IR 101 | ISO 2022 | ESC 02/13 04/02 | ISO-IR 101 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.3 | ISO 2022 IR 109 | ISO 2022 | ESC 02/13 04/03 | ISO-IR 109 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |
| Latin alphabet No.4 | ISO 2022 IR 110 | ISO 2022 | ESC 02/13 04/04 | ISO-IR 110 | Supplementary set |
| | | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6 | ISO 646 |

Multi-Byte Character Sets without Code Extension

Table 6-3: Multi-Byte Character Sets without Code Extension

| Character Set Description | Defined Term | ISO registration number | Character Set |
|---------------------------|--------------|-------------------------|---|
| Unicode | ISO_IR 192 | ISO 10646 | Unicode in UTF-8 |
| Chinese | GB18030 | GB18030 | GB 18030-2000 (China Association for Standardization) |

Table 6-4: Multi-Byte Character Sets with Code Extension

| Character Set Description | Defined Term | Standard for Code Extension | ESC sequence | ISO registration number | Character Set |
|---------------------------|-----------------|-----------------------------|-----------------------|-------------------------|-------------------------------------|
| Japanese | ISO 2022 IR 159 | ISO 2022 | ESC 02/04 02/08 04/04 | ISO-IR 159 | JIS X 0212: Supplementary Kanji set |
| Korean | ISO 2022 IR 149 | ISO 2022 | ESC 02/04 02/09 04/03 | ISO-IR 149 | KS X 1001: Hangul and Hanja |

All SCS (Special Character Sets) listed above are supported for incoming Data.

When there is a mismatch between the given character set in attribute (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM conform:

- Convert each illegal character to a '?’.

There are now three categories of character sets which have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO_IR 6, ISO 2022 IR 6, ISO_IR 100, etc. → encoded in ISO 2022
- ISO_IR 192 → encoded in UTF-8
- GB18030 → encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO_IR 192 \leftrightarrow (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030 \leftrightarrow (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022 \leftrightarrow (0008,0005) contains ISO_IR 192
- An attribute value is encoded in ISO 2022 \leftrightarrow (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies does not conform to DICOM. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

The *teamplay Receiver* supports Kanji characters in the byte zone after 74 (79, 7A, 7B and 7C).

7 Attribute confidentiality profiles

7.1 De-identification

The *teamplay Receiver* can de-identify attributes using proprietary mechanisms.

8 Security

8.1 Security Profiles

The *teamplay Receiver* does not support any specific security measures.

8.2 Association Level Security

It is possible to configure whether the SCP will only answer to known AETs or to any AET.

8.3 Application Level Security

N/A

9 Annexes

9.1 IOD Contents

N/A

9.2 Data Dictionary of Private Attributes

N/A

9.3 Coded Terminology and Templates

N/A

9.3.1 Context Groups

N/A

9.3.2 Template Specifications

N/A

9.3.3 Private Code definitions

| Attribute Name | Attribute Tag | VR | VM | Private Creator Code |
|--------------------|---------------|----|----|---|
| Private Creator ID | (0031, 00xx) | LO | 1 | Constant value "SIEMENS SY CLOUD". Reserve elements in block (0031, xx00-xxFF) |
| Snapshot | (0031, xx01) | LT | 1 | The value present in this field represent the deidentification rules applied. |
| Toolkit version | (0031, xx02) | LO | 1 | The value present in the field represent the version of the toolkit that performed the deidentification |

9.4 Grayscale Image Consistency

N/A

9.5 Standard Extended / Specialized / Private SOP Classes

N/A

9.6 Private Transfer Syntaxes

None

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/ All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 80
[siemens.com/healthineers](https://www.siemens.com/healthineers)

© Siemens Healthcare GmbH, 04.2020