

# One view or two views for wide-angle tomosynthesis with synthetic mammography in the assessment setting?



**Authors** Clauser P, Baltzer PAT, Kapetas P, Woitek R, Weber M, Leone F, Bernathova M, Helbich TH.

## Objectives

To evaluate the diagnostic performance in the assessment setting of three protocols: one-view wide-angle digital breast tomosynthesis (WA-DBT) with synthetic mammography (SM), two-view WA-DBT/SM, and two-view digital mammography (DM).

## Methods

Included in this retrospective study were patients who underwent bilateral two-view DM and WA-DBT. SM were reconstructed from the WA-DBT data. The standard of reference was histology and/or 2 years follow-up. Included were 205 women with 179 lesions (89 malignant, 90 benign). Four blinded readers randomly evaluated images to assess density, lesion type, and level of suspicion according to BI-RADS. Three protocols were evaluated: two-view DM, one-view (mediolateral oblique) WA-DBT/SM, and two-view WA-DBT/SM. Detection rate, sensitivity, specificity, and accuracy were calculated and compared using multivariate analysis. Reading time was assessed.

## Results

The detection rate was higher with two-view WA-DBT/SM ( $p = 0.063$ ). Sensitivity was higher for two-view WA-DBT/SM compared to two-view DM ( $p = 0.001$ ) and one-view WA-DBT/SM ( $p = 0.058$ ). No significant differences in specificity were found. Accuracy was higher with both one-view WA-DBT/SM and two-view WA-DBT/SM compared to DM ( $p = 0.003$  and  $> 0.001$ , respectively). Accuracy did not differ between one- and two-view WA-DBT/SM. Two-view WA-DBT/SM performed better for masses and asymmetries. Reading times were significantly longer when WA-DBT was evaluated. could be up to 0.63% and 0.56%, if the interval cancers selected for consensus were detected at screening. In the former scenario, screen-reading volume would be reduced by 50%, while the latter would reduce the volume by 90%.



## Conclusion

One-view and two-view WA-DBT/SM can achieve a higher diagnostic performance compared to two-view DM. The detection rate and sensitivity were highest with two-view WA-DBT/SM. Two-view WA-DBT/SM appears to be the most appropriate tool for the assessment of breast lesions.

> [Click here](#) to read the article

Clauser P, Baltzer PAT, Kapetas P, Woitek R, Weber M, Leone F, Bernathova M, Helbich TH. One view or two views for wide-angle tomosynthesis with synthetic mammography in the assessment setting? *Eur Radiol.* 2022;32:661-670. doi: 10.1007/s00330-021-08079-2